

JC06 Rec'd PCT/PTO 19 MAY 2005

Product and Method

The present invention relates to oligonucleotide probes, for use in assessing gene transcript levels in a cell, which may be used in analytical techniques, particularly diagnostic techniques. Conveniently the probes are provided in kit form. Different sets of probes may be used in techniques to prepare gene expression patterns and identify, diagnose or monitor different states, such as diseases, conditions or stages thereof. Also provided are methods of identifying suitable probes and their use in methods of the invention.

The identification of quick and easy methods of sample analysis for, for example, diagnostic applications, remains the goal of many researchers. End users seek methods which are cost effective, produce statistically significant results and which may be implemented routinely without the need for highly skilled individuals.

The analysis of gene expression within cells has been used to provide information on the state of those cells and importantly the state of the individual from which the cells are derived. The relative expression of various genes in a cell has been identified as reflecting a particular state within a body. For example, cancer cells are known to exhibit altered expression of various proteins and the transcripts or the expressed proteins may therefore be used as markers of that disease state.

Thus biopsy tissue may be analysed for the presence of these markers and cells originating from the site of the disease may be identified in other tissues or fluids of the body by the presence of the markers.

Furthermore, products of the altered expression may be

- 2 -

released into the blood stream and these products may be analysed. In addition cells which have contacted disease cells may be affected by their direct contact with those cells resulting in altered gene expression and their expression or products of expression may be similarly analysed.

However, there are some limitations with these methods. For example, the use of specific tumour markers for identifying cancer suffers from a variety of defects, such as lack of specificity or sensitivity, association of the marker with disease states besides the specific type of cancer, and difficulty of detection in asymptomatic individuals.

In addition to the analysis of one or two marker transcripts or proteins, more recently, gene expression patterns have been analysed. Most of the work involving large-scale gene expression analysis with implications in disease diagnosis has involved clinical samples originating from diseased tissues or cells. For example, several recent publications, which demonstrate that gene expression data can be used to distinguish between similar cancer types, have used clinical samples from diseased tissues or cells (Alon et al. 1999, PNAS, 96, p6745-6750; Golub et al. 1999, Science, 286, p531-537; Alizadeh et al, 2000, Nature, 403, p503-511; Bittner et al., 2000, Nature, 406, p536-540).

However, these methods have relied on analysis of a sample containing diseased cells or products of those cells or cells which have been contacted by disease cells. Analysis of such samples relies on knowledge of the presence of a disease and its location, which may be difficult in asymptomatic patients. Furthermore, samples can not always be taken from the disease site, e.g. in diseases of the brain.

In a finding of great significance, the present inventors identified the previously untapped potential of all cells within a body to provide information

- 3 -

relating to the state of the organism from which the cells were derived. WO98/49342 describes the analysis of the gene expression of cells distant from the site of disease, e.g. peripheral blood collected distant from a cancer site.

This finding is based on the premise that the different parts of an organism's body exist in dynamic interaction with each other. When a disease affects one part of the body, other parts of the body are also affected. The interaction results from a wide spectrum of biochemical signals that are released from the diseased area, affecting other areas in the body. Although, the nature of the biochemical and physiological changes induced by the released signals can vary in the different body parts, the changes can be measured at the level of gene expression and used for diagnostic purposes.

The physiological state of a cell in an organism is determined by the pattern with which genes are expressed in it. The pattern depends upon the internal and external biological stimuli to which said cell is exposed, and any change either in the extent or in the nature of these stimuli can lead to a change in the pattern with which the different genes are expressed in the cell. There is a growing understanding that by analysing the systemic changes in gene expression patterns in cells in biological samples, it is possible to provide information on the type and nature of the biological stimuli that are acting on them. Thus, for example, by monitoring the expression of a large number of genes in cells in a test sample, it is possible to determine whether their genes are expressed with a pattern characteristic for a particular disease, condition or stage thereof. Measuring changes in gene activities in cells, e.g. from tissue or body fluids is therefore emerging as a powerful tool for disease diagnosis.

- 4 -

Such methods have various advantages. Often, obtaining clinical samples from certain areas in the body that is diseased can be difficult and may involve undesirable invasions in the body, for example biopsy is often used to obtain samples for cancer. In some cases, such as in Alzheimer's disease the diseased brain specimen can only be obtained post-mortem. Furthermore, the tissue specimens which are obtained are often heterogeneous and may contain a mixture of both diseased and non-diseased cells, making the analysis of generated gene expression data both complex and difficult.

It has been suggested that a pool of tumour tissues that appear to be pathogenetically homogeneous with respect to morphological appearances of the tumour may well be highly heterogeneous at the molecular level (Alizadeh, 2000, supra), and in fact might contain tumours representing essentially different diseases (Alizadeh, 2000, supra; Golub, 1999, supra). For the purpose of identifying a disease, condition, or a stage thereof, any method that does not require clinical samples to originate directly from diseased tissues or cells is highly desirable since clinical samples representing a homogeneous mixture of cell types can be obtained from an easily accessible region in the body.

We have now identified a set of probes of surprising utility for identifying one or more diseases. Thus, we now describe probes and sets of probes derived from cells which are not disease cells and which have not contacted disease cells, which correspond to genes which exhibit altered expression in normal versus disease individuals, for use in methods of identifying, diagnosing or monitoring certain conditions, particularly diseases or stages thereof.

Thus the invention provides a set of oligonucleotide probes which correspond to genes in a cell whose expression is affected in a pattern characteristic of a particular disease, condition or

- 5 -

stage thereof, wherein said genes are systemically affected by said disease, condition or stage thereof. Preferably said genes are metabolic or house-keeping genes and preferably are constitutively moderately or highly expressed. Preferably the genes are moderately or highly expressed in the cells of the sample but not in cells from disease cells or in cells having contacted such disease cells.

Such probes, particularly when isolated from cells distant to the site of disease, do not rely on the development of disease to clinically recognizable levels and allow detection of a disease or condition or stage thereof very early after the onset of said disease or condition, even years before other subjective or objective symptoms appear.

As used herein "systemically" affected genes refers to genes whose expression is affected in the body without direct contact with a disease cell or disease site and the cells under investigation are not disease cells.

"Contact" as referred to herein refers to cells coming into close proximity with one another such that the direct effect of one cell on the other may be observed, e.g. an immune response, wherein these responses are not mediated by secondary molecules released from the first cell over a large distance to affect the second cell. Preferably contact refers to physical contact, or contact that is as close as is sterically possible, conveniently, cells which contact one another are found in the same unit volume, for example within 1cm^3 .

A "disease cell" is a cell manifesting phenotypic changes and is present at the disease site at some time during its life-span, e.g. a tumour cell at the tumour site or which has disseminated from the tumour, or a brain cell in the case of brain disorders such as Alzheimer's disease.

- 6 -

"Metabolic" or "house-keeping" genes refer to those genes responsible for expressing products involved in cell division and maintenance, e.g. non-immune function related genes.

5 "Moderately or highly" expressed genes refers to those present in resting cells in a copy number of more than 30-100 copies/cell (assuming an average 3×10^5 mRNA molecules in a cell).

10 Specific probes having the above described properties are provided herein.

Thus in one aspect, the present invention provides a set of oligonucleotide probes, wherein said set comprises at least 10 oligonucleotides selected from:

15 an oligonucleotide as described in Table 1 or derived from a sequence described in Table 1, or an oligonucleotide with a complementary sequence, or a functionally equivalent oligonucleotide.

20 "Table 1" as referred to herein refers to Table 1a and/or Table 1b. Table 1b contains reference to additional clones and sequences as disclosed herein. Similarly Tables 2 and 4 comprise 2 parts, a and b.

25 The invention also provides one or more oligonucleotide probes, wherein each oligonucleotide probe is selected from the oligonucleotides listed in Table 1, or derived from a sequence described in Table 1, or a complementary sequence thereof. The use of such probes in products and methods of the invention, form further aspects of the invention.

30 As referred to herein an "oligonucleotide" is a nucleic acid molecule having at least 6 monomers in the polymeric structure, ie. nucleotides or modified forms thereof. The nucleic acid molecule may be DNA, RNA or PNA (peptide nucleic acid) or hybrids thereof or modified versions thereof, e.g. chemically modified
35 forms, e.g. LNA (Locked Nucleic acid), by methylation or made up of modified or non-natural bases during synthesis, providing they retain their ability to bind

- 7 -

to complementary sequences. Such oligonucleotides are used in accordance with the invention to probe target sequences and are thus referred to herein also as oligonucleotide probes or simply as probes.

5 An "oligonucleotide derived from a sequence described in Table 1" (or any other table) refers to a part of a sequence disclosed in that Table (e.g. Table 1-4), which satisfies the requirements of the oligonucleotide probes as described herein, e.g. in
10 length and function. Preferably said parts have the size described hereinafter.

Preferably the oligonucleotide probes forming said set are at least 15 bases in length to allow binding of target molecules. Especially preferably said
15 oligonucleotide probes are from 20 to 200 bases in length, e.g. from 30 to 150 bases, preferably 50-100 bases in length.

As referred to herein the term "complementary sequences" refers to sequences with consecutive
20 complementary bases (ie. T:A, G:C) and which complementary sequences are therefore able to bind to one another through their complementarity.

Reference to "10 oligonucleotides" refers to 10 different oligonucleotides. Whilst a Table 1
25 oligonucleotide, a Table 1 derived oligonucleotide and their functional equivalent are considered different oligonucleotides, complementary oligonucleotides are not considered different. Preferably however, the at least 10 oligonucleotides are 10 different Table 1
30 oligonucleotides (or Table 1 derived oligonucleotides or their functional equivalents). Thus said 10 different oligonucleotides are preferably able to bind to 10 different transcripts.

Preferably said oligonucleotides are as described
35 in Table 1 or are derived from a sequence described in Table 1. Especially preferably said oligonucleotides are as described in Table 2 or Table 4 or are derived

- 8 -

from a sequence described in either of those tables. Especially preferably the oligonucleotide (or the oligonucleotide derived therefrom) has a high occurrence as defined in Table 3, especially preferably >40%, e.g. >80 or >90, e.g. 100%.

A "set" as described refers to a collection of unique oligonucleotide probes (ie. having a distinct sequence) and preferably consists of less than 1000 oligonucleotide probes, especially less than 500 probes, e.g. preferably from 10 to 500, e.g. 10 to 100, 200 or 300, especially preferably 20 to 100, e.g. 30 to 100 probes. In some cases less than 10 probes may be used, e.g. from 2 to 9 probes, e.g. 5 to 9 probes.

It will be appreciated that increasing the number of probes will prevent the possibility of poor analysis, e.g. misdiagnosis by comparison to other diseases which could similarly alter the expression of the particular genes in question. Other oligonucleotide probes not described herein may also be present, particularly if they aid the ultimate use of the set of oligonucleotide probes. However, preferably said set consists only of said Table 1 oligonucleotides, Table 1 derived oligonucleotides, complementary sequences or functionally equivalent oligonucleotides, or a sub-set thereof (e.g. of the size as described above), preferably a sub-set for which sequences are provided herein (see Table 1 and its footnote). Especially preferably said set consists only of said Table 1 oligonucleotides, Table 1 derived oligonucleotides, or complementary sequences thereof, or a sub-set thereof.

Multiple copies of each unique oligonucleotide probe, e.g. 10 or more copies, may be present in each set, but constitute only a single probe.

A set of oligonucleotide probes, which may preferably be immobilized on a solid support or have means for such immobilization, comprises the at least 10 oligonucleotide probes selected from those described

- 9 -

hereinbefore. Especially preferably said probes are selected from those having high occurrence as described in Table 3 and as mentioned above. As mentioned above, these 10 probes must be unique and have different sequences. Having said this however, two separate probes may be used which recognize the same gene but reflect different splicing events. However oligonucleotide probes which are complementary to, and bind to distinct genes are preferred.

As described herein a "functionally equivalent" oligonucleotide to those described in Table 1 or derived therefrom refers to an oligonucleotide which is capable of identifying the same gene as an oligonucleotide of Table 1 or derived therefrom, ie. it can bind to the same mRNA molecule (or DNA) transcribed from a gene (target nucleic acid molecule) as the Table 1 oligonucleotide or the Table 1 derived oligonucleotide (or its complementary sequence). Preferably said functionally equivalent oligonucleotide is capable of recognizing, ie. binding to the same splicing product as a Table 1 oligonucleotide or a Table 1 derived oligonucleotide. Preferably said mRNA molecule is the full length mRNA molecule which corresponds to the Table 1 oligonucleotide or the Table 1 derived oligonucleotide.

As referred to herein "capable of binding" or "binding" refers to the ability to hybridize under conditions described hereinafter.

Alternatively expressed, functionally equivalent oligonucleotides (or complementary sequences) have sequence identity or will hybridize, as described hereinafter, to a region of the target molecule to which molecule a Table 1 oligonucleotide or a Table 1 derived oligonucleotide or a complementary oligonucleotide binds. Preferably, functionally equivalent oligonucleotides (or their complementary sequences) hybridize to one of the mRNA sequences which corresponds

- 10 -

to a Table 1 oligonucleotide or a Table 1 derived oligonucleotide under the conditions described hereinafter or has sequence identity to a part of one of the mRNA sequences which corresponds to a Table 1
5 oligonucleotide or a Table 1 derived oligonucleotide. A "part" in this context refers to a stretch of at least 5, e.g. at least 10 or 20 bases, such as from 5 to 100, e.g. 10 to 50 or 15 to 30 bases.

In a particularly preferred aspect, the
10 functionally equivalent oligonucleotide binds to all or a part of the region of a target nucleic acid molecule (mRNA or cDNA) to which the Table 1 oligonucleotide or Table 1 derived oligonucleotide binds. A "target" nucleic acid molecule is the gene transcript or related
15 product e.g. mRNA, or cDNA, or amplified product thereof. Said "region" of said target molecule to which said Table 1 oligonucleotide or Table 1 derived oligonucleotide binds is the stretch over which complementarity exists. At its largest this region is
20 the whole length of the Table 1 oligonucleotide or Table 1 derived oligonucleotide, but may be shorter if the entire Table 1 sequence or Table 1 derived oligonucleotide is not complementary to a region of the target sequence.

25 Preferably said part of said region of said target molecule is a stretch of at least 5, e.g. at least 10 or 20 bases, such as from 5 to 100, e.g. 10 to 50 or 15 to 30 bases. This may for example be achieved by said functionally equivalent oligonucleotide having several
30 identical bases to the bases of the Table 1 oligonucleotide or the Table 1 derived oligonucleotide. These bases may be identical over consecutive stretches, e.g. in a part of the functionally equivalent oligonucleotide, or may be present non-consecutively,
35 but provide sufficient complementarity to allow binding to the target sequence.

Thus in a preferred feature, said functionally

- 11 -

equivalent oligonucleotide hybridizes under conditions of high stringency to a Table 1 oligonucleotide or a Table 1 derived oligonucleotide or the complementary sequence thereof. Alternatively expressed, said functionally equivalent oligonucleotide exhibits high sequence identity to all or part of a Table 1 oligonucleotide. Preferably said functionally equivalent oligonucleotide has at least 70% sequence identity, preferably at least 80%, e.g. at least 90, 95, 98 or 99%, to all of a Table 1 oligonucleotide or a part thereof. As used in this context, a "part" refers to a stretch of at least 5, e.g. at least 10 or 20 bases, such as from 5 to 100, e.g. 10 to 50 or 15 to 30 bases, in said Table 1 oligonucleotide. Especially preferably when sequence identity to only a part of said Table 1 oligonucleotide is present, the sequence identity is high, e.g. at least 80% as described above.

Functionally equivalent oligonucleotides which satisfy the above stated functional requirements include those which are derived from the Table 1 oligonucleotides and also those which have been modified by single or multiple nucleotide base (or equivalent) substitution, addition and/or deletion, but which nonetheless retain functional activity, e.g. bind to the same target molecule as the Table 1 oligonucleotide or the Table 1 derived oligonucleotide from which they are further derived or modified. Preferably said modification is of from 1 to 50, e.g. from 10 to 30, preferably from 1 to 5 bases. Especially preferably only minor modifications are present, e.g. variations in less than 10 bases, e.g. less than 5 base changes.

Within the meaning of "addition" equivalents are included oligonucleotides containing additional sequences which are complementary to the consecutive stretch of bases on the target molecule to which the Table 1 oligonucleotide or the Table 1 derived oligonucleotide binds. Alternatively the addition may

- 12 -

comprise a different, unrelated sequence, which may for example confer a further property, e.g. to provide a means for immobilization such as a linker to bind the oligonucleotide probe to a solid support.

5 Particularly preferred are naturally occurring equivalents such as biological variants, e.g. allelic, geographical or allotypic variants, e.g. oligonucleotides which correspond to a genetic variant, for example as present in a different species.

10 Functional equivalents include oligonucleotides with modified bases, e.g. using non-naturally occurring bases. Such derivatives may be prepared during synthesis or by post production modification.

15 "Hybridizing" sequences which bind under conditions of low stringency are those which bind under non-stringent conditions (for example, 6x SSC/50% formamide at room temperature) and remain bound when washed under conditions of low stringency (2 X SSC, room temperature, more preferably 2 X SSC, 42°C). Hybridizing under high
20 stringency refers to the above conditions in which washing is performed at 2 X SSC, 65°C (where SSC = 0.15M NaCl, 0.015M sodium citrate, pH 7.2).

25 "Sequence identity" as referred to herein refers to the value obtained when assessed using ClustalW (Thompson et al., 1994, Nucl. Acids Res., 22, p4673-4680) with the following parameters:

Pairwise alignment parameters - Method: accurate, Matrix: IUB, Gap open penalty: 15.00, Gap extension penalty: 6.66;

30 Multiple alignment parameters - Matrix: IUB, Gap open penalty: 15.00, % identity for delay: 30, Negative matrix: no, Gap extension penalty: 6.66, DNA transitions weighting: 0.5.

35 Sequence identity at a particular base is intended to include identical bases which have simply been derivatized.

The invention also extends to polypeptides encoded

- 13 -

by the mRNA sequence to which a Table 1 oligonucleotide or a Table 1 derived oligonucleotide binds. The invention further extends to antibodies which bind to any of said polypeptides.

5 As described above, conveniently said set of oligonucleotide probes may be immobilized on one or more solid supports. Single or preferably multiple copies of each unique probe are attached to said solid supports, e.g. 10 or more, e.g. at least 100 copies of each unique
10 probe are present.

One or more unique oligonucleotide probes may be associated with separate solid supports which together form a set of probes immobilized on multiple solid support, e.g. one or more unique probes may be
15 immobilized on multiple beads, membranes, filters, biochips etc. which together form a set of probes, which together form modules of the kit described hereinafter. The solid support of the different modules are conveniently physically associated although the signals
20 associated with each probe (generated as described hereinafter) must be separately determinable.

Alternatively, the probes may be immobilized on discrete portions of the same solid support, e.g. each unique oligonucleotide probe, e.g. in multiple copies,
25 may be immobilized to a distinct and discrete portion or region of a single filter or membrane, e.g. to generate an array.

A combination of such techniques may also be used, e.g. several solid supports may be used which each
30 immobilize several unique probes.

The expression "solid support" shall mean any solid material able to bind oligonucleotides by hydrophobic, ionic or covalent bridges.

"Immobilization" as used herein refers to
35 reversible or irreversible association of the probes to said solid support by virtue of such binding. If reversible, the probes remain associated with the solid

- 14 -

support for a time sufficient for methods of the invention to be carried out.

5 Numerous solid supports suitable as immobilizing moieties according to the invention, are well known in the art and widely described in the literature and generally speaking, the solid support may be any of the well-known supports or matrices which are currently widely used or proposed for immobilization, separation etc. in chemical or biochemical procedures. Such 10 materials include, but are not limited to, any synthetic organic polymer such as polystyrene, polyvinylchloride, polyethylene; or nitrocellulose and cellulose acetate; or tosyl activated surfaces; or glass or nylon or any surface carrying a group suited for covalent coupling of 15 nucleic acids. The immobilizing moieties may take the form of particles, sheets, gels, filters, membranes, microfibre strips, tubes or plates, fibres or capillaries, made for example of a polymeric material e.g. agarose, cellulose, alginate, teflon, latex or 20 polystyrene or magnetic beads. Solid supports allowing the presentation of an array, preferably in a single dimension are preferred, e.g. sheets, filters, membranes, plates or biochips.

Attachment of the nucleic acid molecules to the 25 solid support may be performed directly or indirectly. For example if a filter is used, attachment may be performed by UV-induced crosslinking. Alternatively, attachment may be performed indirectly by the use of an attachment moiety carried on the oligonucleotide probes and/or solid support. Thus for example, a pair of 30 affinity binding partners may be used, such as avidin, streptavidin or biotin, DNA or DNA binding protein (e.g. either the lac I repressor protein or the lac operator sequence to which it binds), antibodies (which may be 35 mono- or polyclonal), antibody fragments or the epitopes or haptens of antibodies. In these cases, one partner of the binding pair is attached to (or is inherently

- 15 -

part of) the solid support and the other partner is attached to (or is inherently part of) the nucleic acid molecules.

5 As used herein an "affinity binding pair" refers to two components which recognize and bind to one another specifically (ie. in preference to binding to other molecules). Such binding pairs when bound together form a complex.

10 Attachment of appropriate functional groups to the solid support may be performed by methods well known in the art, which include for example, attachment through hydroxyl, carboxyl, aldehyde or amino groups which may be provided by treating the solid support to provide suitable surface coatings. Solid supports presenting
15 appropriate moieties for attachment of the binding partner may be produced by routine methods known in the art.

Attachment of appropriate functional groups to the oligonucleotide probes of the invention may be performed
20 by ligation or introduced during synthesis or amplification, for example using primers carrying an appropriate moiety, such as biotin or a particular sequence for capture.

Conveniently, the set of probes described
25 hereinbefore is provided in kit form.

Thus viewed from a further aspect the present invention provides a kit comprising a set of oligonucleotide probes as described hereinbefore immobilized on one or more solid supports.

30 Preferably, said probes are immobilized on a single solid support and each unique probe is attached to a different region of said solid support. However, when attached to multiple solid supports, said multiple solid supports form the modules which make up the kit.
35 Especially preferably said solid support is a sheet, filter, membrane, plate or biochip.

Optionally the kit may also contain information

- 16 -

relating to the signals generated by normal or diseased samples (as discussed in more detail hereinafter in relation to the use of the kits), standardizing materials, e.g. mRNA or cDNA from normal and/or diseased samples for comparative purposes, labels for incorporation into cDNA, adapters for introducing nucleic acid sequences for amplification purposes, primers for amplification and/or appropriate enzymes, buffers and solutions. Optionally said kit may also contain a package insert describing how the method of the invention should be performed, optionally providing standard graphs, data or software for interpretation of results obtained when performing the invention.

The use of such kits to prepare a standard diagnostic gene transcript pattern as described hereinafter forms a further aspect of the invention.

The set of probes as described herein have various uses. Principally however they are used to assess the gene expression state of a test cell to provide information relating to the organism from which said cell is derived. Thus the probes are useful in diagnosing, identifying or monitoring a disease or condition or stage thereof in an organism.

Thus in a further aspect the invention provides the use of a set of oligonucleotide probes or a kit as described hereinbefore to determine the gene expression pattern of a cell which pattern reflects the level of gene expression of genes to which said oligonucleotide probes bind, comprising at least the steps of:

- a) isolating mRNA from said cell, which may optionally be reverse transcribed to cDNA;
- b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotide probes or a kit as defined herein; and
- c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce said pattern.

The mRNA and cDNA as referred to in this method,

- 17 -

and the methods hereinafter, encompass derivatives or copies of said molecules, e.g. copies of such molecules such as those produced by amplification or the preparation of complementary strands, but which retain the identity of the mRNA sequence, ie. would hybridize to the direct transcript (or its complementary sequence) by virtue of precise complementarity, or sequence identity, over at least a region of said molecule. It will be appreciated that complementarity will not exist over the entire region where techniques have been used which may truncate the transcript or introduce new sequences, e.g. by primer amplification. For convenience, said mRNA or cDNA is preferably amplified prior to step b). As with the oligonucleotides described herein said molecules may be modified, e.g. by using non-natural bases during synthesis providing complementarity remains. Such molecules may also carry additional moieties such as signalling or immobilizing means.

The various steps involved in the method of preparing such a pattern are described in more detail hereinafter.

As used herein "gene expression" refers to transcription of a particular gene to produce a specific mRNA product (ie. a particular splicing product). The level of gene expression may be determined by assessing the level of transcribed mRNA molecules or cDNA molecules reverse transcribed from the mRNA molecules or products derived from those molecules, e.g. by amplification.

The "pattern" created by this technique refers to information which, for example, may be represented in tabular or graphical form and conveys information about the signal associated with two or more oligonucleotides. Preferably said pattern is expressed as an array of numbers relating to the expression level associated with each probe.

- 18 -

Preferably, said pattern is established using the following linear model:

$$y = Xb + f$$

Equation 1

wherein, X is the matrix of gene expression data and y is the response variable, b is the regression coefficient vector and f the estimated residual vector. Although many different methods can be used to establish the relationship provided in equation 1, especially preferably the partial Least Squares Regression (PLSR) method is used for establishing the relationship in equation 1.

The probes are thus used to generate a pattern which reflects the gene expression of a cell at the time of its isolation. The pattern of expression is characteristic of the circumstances under which that cells finds itself and depends on the influences to which the cell has been exposed. Thus, a characteristic gene transcript pattern standard or fingerprint (standard probe pattern) for cells from an individual with a particular disease or condition may be prepared and used for comparison to transcript patterns of test cells. This has clear applications in diagnosing, monitoring or identifying whether an organism is suffering from a particular disease, condition or stage thereof.

The standard pattern is prepared by determining the extent of binding of total mRNA (or cDNA or related product), from cells from a sample of one or more organisms with the disease or condition or stage thereof, to the probes. This reflects the level of transcripts which are present which correspond to each unique probe. The amount of nucleic acid material which binds to the different probes is assessed and this information together forms the gene transcript pattern standard of that disease or condition or stage thereof. Each such standard pattern is characteristic of the disease, condition or stage thereof.

- 19 -

In a further aspect therefore, the present invention provides a method of preparing a standard gene transcript pattern characteristic of a disease or condition or stage thereof in an organism comprising at least the steps of:

- 5 a) isolating mRNA from the cells of a sample of one or more organisms having the disease or condition or stage thereof, which may optionally be reverse transcribed to cDNA;
- 10 b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotides or a kit as described hereinbefore specific for said disease or condition or stage thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and
- 15 c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce a characteristic pattern reflecting the level of gene expression of genes to which said oligonucleotides bind, in the sample with the disease, condition or stage thereof.

For convenience, said oligonucleotides are preferably immobilized on one or more solid supports.

The standard pattern for a great number of diseases or conditions and different stages thereof using particular probes may be accumulated in databases and be made available to laboratories on request.

"Disease" samples and organisms as referred to herein refer to organisms (or samples from the same) with an underlying pathological disturbance relative to a normal organism (or sample), in a symptomatic or asymptomatic organism, which may result, for example, from infection or an acquired or congenital genetic imperfection. Such organisms are known to have, or which exhibit, the disease or condition or stage thereof under study.

A "condition" refers to a state of the mind or body of an organism which has not occurred through disease,

- 20 -

e.g. the presence of an agent in the body such as a toxin, drug or pollutant, or pregnancy.

"Stages" thereof refer to different stages of the disease or condition which may or may not exhibit particular physiological or metabolic changes, but do exhibit changes at the genetic level which may be detected as altered gene expression. It will be appreciated that during the course of a disease or condition the expression of different transcripts may vary. Thus at different stages, altered expression may not be exhibited for particular transcripts compared to "normal" samples. However, combining information from several transcripts which exhibit altered expression at one or more stages through the course of the disease or condition can be used to provide a characteristic pattern which is indicative of a particular stage of the disease or condition. Thus for example different stages in cancer, e.g. pre-stage I, stage I, stage II, III or IV can be identified.

"Normal" as used herein refers to organisms or samples which are used for comparative purposes. Preferably, these are "normal" in the sense that they do not exhibit any indication of, or are not believed to have, any disease or condition that would affect gene expression, particularly in respect of the disease for which they are to be used as the normal standard. However, it will be appreciated that different stages of a disease or condition may be compared and in such cases, the "normal" sample may correspond to the earlier stage of the disease or condition.

As used herein a "sample" refers to any material obtained from the organism, e.g. human or non-human animal under investigation which contains cells and includes, tissues, body fluid or body waste or in the case of prokaryotic organisms, the organism itself. "Body fluids" include blood, saliva, spinal fluid, semen, lymph. "Body waste" includes urine, expectorated

- 21 -

matter (pulmonary patients), faeces etc. "Tissue samples" include tissue obtained by biopsy, by surgical interventions or by other means e.g. placenta.

5 Preferably however, the samples which are examined are from areas of the body not apparently affected by the disease or condition. The cells in such samples are not disease cells, e.g. cancer cells, have not been in contact with such disease cells and do not originate from the site of the disease or condition. The "site of
10 disease" is considered to be that area of the body which manifests the disease in a way which may be objectively determined, e.g. a tumour or area of inflammation. Thus for example peripheral blood may be used for the diagnosis of non-haematopoietic cancers, and the blood
15 does not require the presence of malignant or disseminated cells from the cancer in the blood. Similarly in diseases of the brain, in which no diseased cells are found in the blood due to the blood:brain barrier, peripheral blood may still be used in the
20 methods of the invention.

It will however be appreciated that the method of preparing the standard transcription pattern and other methods of the invention are also applicable for use on living parts of eukaryotic organisms such as cell lines
25 and organ cultures and explants.

As used herein, reference to "corresponding" sample etc. refers to cells preferably from the same tissue, body fluid or body waste, but also includes cells from tissue, body fluid or body waste which are sufficiently
30 similar for the purposes of preparing the standard or test pattern. When used in reference to genes "corresponding" to the probes, this refers to genes which are related by sequence (which may be complementary) to the probes although the probes may
35 reflect different splicing products of expression.

"Assessing" as used herein refers to both quantitative and qualitative assessment which may be

- 22 -

determined in absolute or relative terms.

The invention may be put into practice as follows. To prepare a standard transcript pattern for a particular disease, condition or stage thereof, sample mRNA is extracted from the cells of tissues, body fluid or body waste according to known techniques (see for example Sambrook et. al. (1989), Molecular Cloning : A laboratory manual, 2nd Ed., Cold Spring Harbor Laboratory Press, Cold Spring Harbor, N.Y.) from a diseased individual or organism.

Owing to the difficulties in working with RNA, the RNA is preferably reverse transcribed at this stage to form first strand cDNA. Cloning of the cDNA or selection from, or using, a cDNA library is not however necessary in this or other methods of the invention. Preferably, the complementary strands of the first strand cDNAs are synthesized, ie. second strand cDNAs, but this will depend on which relative strands are present in the oligonucleotide probes. The RNA may however alternatively be used directly without reverse transcription and may be labelled if so required.

Preferably the cDNA strands are amplified by known amplification techniques such as the polymerase chain reaction (PCR) by the use of appropriate primers. Alternatively, the cDNA strands may be cloned with a vector, used to transform a bacteria such as E. coli which may then be grown to multiply the nucleic acid molecules. When the sequence of the cDNAs are not known, primers may be directed to regions of the nucleic acid molecules which have been introduced. Thus for example, adapters may be ligated to the cDNA molecules and primers directed to these portions for amplification of the cDNA molecules. Alternatively, in the case of eukaryotic samples, advantage may be taken of the polyA tail and cap of the RNA to prepare appropriate primers.

To produce the standard diagnostic gene transcript pattern or fingerprint for a particular disease or

- 23 -

condition or stage thereof, the above described oligonucleotide probes are used to probe mRNA or cDNA of the diseased sample to produce a signal for hybridization to each particular oligonucleotide probe species, ie. each unique probe. A standard control gene transcript pattern may also be prepared if desired using mRNA or cDNA from a normal sample. Thus, mRNA or cDNA is brought into contact with the oligonucleotide probe under appropriate conditions to allow hybridization.

When multiple samples are probed, this may be performed consecutively using the same probes, e.g. on one or more solid supports, ie. on probe kit modules, or by simultaneously hybridizing to corresponding probes, e.g. the modules of a corresponding probe kit.

To identify when hybridization occurs and obtain an indication of the number of transcripts/cDNA molecules which become bound to the oligonucleotide probes, it is necessary to identify a signal produced when the transcripts (or related molecules) hybridize (e.g. by detection of double stranded nucleic acid molecules or detection of the number of molecules which become bound, after removing unbound molecules, e.g. by washing).

In order to achieve a signal, either or both components which hybridize (ie. the probe and the transcript) carry or form a signalling means or a part thereof. This "signalling means" is any moiety capable of direct or indirect detection by the generation or presence of a signal. The signal may be any detectable physical characteristic such as conferred by radiation emission, scattering or absorption properties, magnetic properties, or other physical properties such as charge, size or binding properties of existing molecules (e.g. labels) or molecules which may be generated (e.g. gas emission etc.). Techniques are preferred which allow signal amplification, e.g. which produce multiple signal events from a single active binding site, e.g. by the catalytic action of enzymes to produce multiple

- 24 -

detectable products.

Conveniently the signalling means may be a label which itself provides a detectable signal. Conveniently this may be achieved by the use of a radioactive or
5 other label which may be incorporated during cDNA production, the preparation of complementary cDNA strands, during amplification of the target mRNA/cDNA or added directly to target nucleic acid molecules.

Appropriate labels are those which directly or
10 indirectly allow detection or measurement of the presence of the transcripts/cDNA. Such labels include for example radiolabels, chemical labels, for example chromophores or fluorophores (e.g. dyes such as fluorescein and rhodamine), or reagents of high electron
15 density such as ferritin, haemocyanin or colloidal gold. Alternatively, the label may be an enzyme, for example peroxidase or alkaline phosphatase, wherein the presence of the enzyme is visualized by its interaction with a suitable entity, for example a substrate. The label may
20 also form part of a signalling pair wherein the other member of the pair is found on, or in close proximity to, the oligonucleotide probe to which the transcript/cDNA binds, for example, a fluorescent compound and a quench fluorescent substrate may be used.
25 A label may also be provided on a different entity, such as an antibody, which recognizes a peptide moiety attached to the transcripts/cDNA, for example attached to a base used during synthesis or amplification.

A signal may be achieved by the introduction of a
30 label before, during or after the hybridization step. Alternatively, the presence of hybridizing transcripts may be identified by other physical properties, such as their absorbance, and in which case the signalling means is the complex itself.

35 The amount of signal associated with each oligonucleotide probe is then assessed. The assessment may be quantitative or qualitative and may be based on

- 25 -

binding of a single transcript species (or related cDNA or other products) to each probe, or binding of multiple transcript species to multiple copies of each unique probe. It will be appreciated that quantitative results will provide further information for the transcript fingerprint of the disease which is compiled. This data may be expressed as absolute values (in the case of macroarrays) or may be determined relative to a particular standard or reference e.g. a normal control sample.

Furthermore it will be appreciated that the standard diagnostic gene pattern transcript may be prepared using one or more disease samples (and normal samples if used) to perform the hybridization step to obtain patterns not biased towards a particular individual's variations in gene expression.

The use of the probes to prepare standard patterns and the standard diagnostic gene transcript patterns thus produced for the purpose of identification or diagnosis or monitoring of a particular disease or condition or stage thereof in a particular organism forms a further aspect of the invention.

Once a standard diagnostic fingerprint or pattern has been determined for a particular disease or condition using the selected oligonucleotide probes, this information can be used to identify the presence, absence or extent or stage of that disease or condition in a different test organism or individual.

To examine the gene expression pattern of a test sample, a test sample of tissue, body fluid or body waste containing cells, corresponding to the sample used for the preparation of the standard pattern, is obtained from a patient or the organism to be studied. A test gene transcript pattern is then prepared as described hereinbefore as for the standard pattern.

In a further aspect therefore, the present invention provides a method of preparing a test gene

- 26 -

transcript pattern comprising at least the steps of:

a) isolating mRNA from the cells of a sample of said test organism, which may optionally be reverse transcribed to cDNA;

5 b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotides or a kit as described hereinbefore specific for a disease or condition or stage thereof in an organism and sample thereof corresponding to the organism and sample thereof under
10 investigation; and

c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce said pattern reflecting the level of gene expression of genes to which said oligonucleotides bind, in said test sample.

15 This test pattern may then be compared to one or more standard patterns to assess whether the sample contains cells having the disease, condition or stage thereof.

20 Thus viewed from a further aspect the present invention provides a method of diagnosing or identifying or monitoring a disease or condition or stage thereof in an organism, comprising the steps of:

a) isolating mRNA from the cells of a sample of said organism, which may optionally be reverse
25 transcribed to cDNA;

b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotides or a kit as described hereinbefore specific for said disease or
30 condition or stage thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation;

c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce a characteristic pattern reflecting the level
35 of gene expression of genes to which said oligonucleotides bind, in said sample; and

d) comparing said pattern to a standard

- 27 -

5 diagnostic pattern prepared according to the method of the invention using a sample from an organism corresponding to the organism and sample under investigation to determine the presence of said disease or condition or a stage thereof in the organism under investigation.

The method up to and including step c) is the preparation of a test pattern as described above.

10 As referred to herein, "diagnosis" refers to determination of the presence or existence of a disease or condition or stage thereof in an organism. "Monitoring" refers to establishing the extent of a disease or condition, particularly when an individual is
15 known to be suffering from a disease or condition, for example to monitor the effects of treatment or the development of a disease or condition, e.g. to determine the suitability of a treatment or provide a prognosis.

The presence of the disease or condition or stage thereof may be determined by determining the degree of
20 correlation between the standard and test samples' patterns. This necessarily takes into account the range of values which are obtained for normal and diseased samples. Although this can be established by obtaining
25 standard deviations for several representative samples binding to the probes to develop the standard, it will be appreciated that single samples may be sufficient to generate the standard pattern to identify a disease if the test sample exhibits close enough correlation to
30 that standard. Conveniently, the presence, absence, or extent of a disease or condition or stage thereof in a test sample can be predicted by inserting the data relating to the expression level of informative probes in test sample into the standard diagnostic probe
35 pattern established according to equation 1.

Data generated using the above mentioned methods may be analysed using various techniques from the most

- 28 -

basic visual representation (e.g. relating to intensity) to more complex data manipulation to identify underlying patterns which reflect the interrelationship of the level of expression of each gene to which the various probes bind, which may be quantified and expressed mathematically. Conveniently, the raw data thus generated may be manipulated by the data processing and statistical methods described hereinafter, particularly normalizing and standardizing the data and fitting the data to a classification model to determine whether said test data reflects the pattern of a particular disease, condition or stage thereof.

The methods described herein may be used to identify, monitor or diagnose a disease, condition or ailment or its stage or progression, for which the oligonucleotide probes are informative. "Informative" probes as described herein, are those which reflect genes which have altered expression in the diseases or conditions in question, or particular stages thereof. Probes of the invention may not be sufficiently informative for diagnostic purposes when used alone, but are informative when used as one of several probes to provide a characteristic pattern, e.g. in a set as described hereinbefore.

Preferably said probes correspond to genes which are systemically affected by said disease, condition or stage thereof. Especially preferably said genes, from which transcripts are derived which bind to probes of the invention, are metabolic or house-keeping genes and preferably are moderately or highly expressed. The advantage of using probes directed to moderately or highly expressed genes is that smaller clinical samples are required for generating the necessary gene expression data set, e.g. less than 1ml blood samples.

Furthermore, it has been found that such genes which are already being actively transcribed tend to be more prone to being influenced, in a positive or

- 29 -

negative way, by new stimuli. In addition, since transcripts are already being produced at levels which are generally detectable, small changes in those levels are readily detectable as for example, a certain
5 detectable threshold does not need to be reached.

In preferred methods of the invention, the set of probes of the invention are informative for a variety of different diseases, conditions or stages thereof. A sub-set of the probes disclosed herein may be used for
10 diagnosis, identification or monitoring a particular disease, condition or stage thereof.

Thus the probes may be used to diagnose or identify or monitor any condition, ailment, disease or reaction that leads to the relative increase or decrease in the activity of informative genes of any or all eukaryotic or prokaryotic organisms regardless of whether these
15 changes have been caused by the influence of bacteria, virus, prions, parasites, fungi, radiation, natural or artificial toxins, drugs or allergens, including mental conditions due to stress, neurosis, psychosis or
20 deteriorations due to the ageing of the organism, and conditions or diseases of unknown cause, providing a sub-set of the probes as described herein are informative for said disease or condition or stage
25 thereof.

Such diseases include those which result in metabolic or physiological changes, such as fever-associated diseases such as influenza or malaria. Other diseases which may be detected include for example
30 yellow fever, sexually transmitted diseases such as gonorrhea, fibromyalgia, candida-related complex, cancer (for example of the stomach, lung, breast, prostate gland, bowel, skin, colon, ovary etc), Alzheimer's disease, disease caused by retroviruses such as HIV,
35 senile dementia, multiple sclerosis and Creutzfeldt-Jakob disease to mention a few.

The invention may also be used to identify patients

- 30 -

with psychiatric or psychosomatic diseases such as schizophrenia and eating disorders. Of particular importance is the use of this method to detect diseases, conditions, or stages thereof, which are not readily detectable by known diagnostic methods, such as HIV which is generally not detectable using known techniques 1 to 4 months following infection. Conditions which may be identified include for example drug abuse, such as the use of narcotics, alcohol, steroids or performance enhancing drugs.

Preferably said disease to be identified or monitored is a cancer or a degenerative brain disorder (such as Alzheimer's or Parkinson's disease).

In particular, a set of oligonucleotide probes, wherein said set comprises at least 10 oligonucleotides selected from:

an oligonucleotide as described in Table 4 or an oligonucleotide derived therefrom or an oligonucleotide with a complementary sequence, or a functionally equivalent oligonucleotide, may be used for diagnosis or identification or monitoring the progression of Alzheimer's disease. Similarly Table 2 probes and Table 2 derived probes and their functional equivalents may be used to diagnose, identify or monitor the progression of breast cancer. Especially preferably the probes used for breast cancer analysis are selected based on their occurrence as set forth in Table 3 and as described hereinbefore.

The diagnostic method may be used alone as an alternative to other diagnostic techniques or in addition to such techniques. For example, methods of the invention may be used as an alternative or additive diagnostic measure to diagnosis using imaging techniques such as Magnetic Resonance Imaging (MRI), ultrasound imaging, nuclear imaging or X-ray imaging, for example in the identification and/or diagnosis of tumours.

The methods of the invention may be performed on

- 31 -

cells from prokaryotic or eukaryotic organisms which may be any eukaryotic organisms such as human beings, other mammals and animals, birds, insects, fish and plants, and any prokaryotic organism such as a bacteria.

5 Preferred non-human animals on which the methods of the invention may be conducted include, but are not limited to mammals, particularly primates, domestic animals, livestock and laboratory animals. Thus preferred animals for diagnosis include mice, rats,
10 guinea pigs, cats, dogs, pigs, cows, goats, sheep, horses. Particularly preferably the disease state or condition of humans is diagnosed, identified or monitored.

As described above, the sample under study may be
15 any convenient sample which may be obtained from an organism. Preferably however, as mentioned above, the sample is obtained from a site distant to the site of disease and the cells in such samples are not disease cells, have not been in contact with such cells and do
20 not originate from the site of the disease or condition. In such cases, although preferably absent, the sample may contain cells which do not fulfil these criteria. However, since the probes of the invention are concerned with transcripts whose expression is altered in cells
25 which do satisfy these criteria, the probes are specifically directed to detecting changes in transcript levels in those cells even if in the presence of other, background cells.

It has been found that the cells from such samples
30 show significant and informative variations in the gene expression of a large number of genes. Thus, the same probe (or several probes) may be found to be informative in determinations regarding two or more diseases, conditions or stages thereof by virtue of the particular
35 level of transcripts binding to that probe or the interrelationship of the extent of binding to that probe relative to other probes. As a consequence, it is

- 32 -

possible to use a relatively small number of probes for screening for multiple disorders or diseases. This has consequences with regard to the selection of probes, discussed in relation to random identification of probes hereinafter, but also for the use of a single set of probes for more than one diagnosis. Table 9 which represents preferred probes of the invention discloses probes which are informative for both Alzheimer's and breast cancer.

Thus, the present invention also provides sets of probes for diagnosing, identifying or monitoring two or more diseases, conditions or stages thereof, wherein at least one of said probes is suitable for said diagnosing, identifying or monitoring at least two of said diseases, conditions or stages thereof, and kits and methods of using the same. Preferably at least 5 probes, e.g. from 5 to 15 probes, are used in at least two diagnoses.

Thus, in a further preferred aspect, the present invention provides a method of diagnosis or identification or monitoring as described hereinbefore for the diagnosis, identification or monitoring of two or more diseases, conditions or stages thereof in an organism, wherein said test pattern produced in step c) of the diagnostic method is compared in step d) to at least two standard diagnostic patterns prepared as described previously, wherein each standard diagnostic pattern is a pattern generated for a different disease or condition or stage thereof.

Whilst in a preferred aspect the methods of assessment concern the development of a gene transcript pattern from a test sample and comparison of the same to a standard pattern, the elevation or depression of expression of certain markers may also be examined by examining the products of expression and the level of those products. Thus a standard pattern in relation to the expressed product may be generated.

- 33 -

In such methods the levels of expression of a set of polypeptides encoded by the gene to which an oligonucleotide of Table 1 or a Table 1 derived oligonucleotide, binds, are analysed.

5 Various diagnostic methods may be used to assess the amount of polypeptides (or fragments thereof) which are present. The presence or concentration of polypeptides may be examined, for example by the use of a binding partner to said polypeptide (e.g. an
10 antibody), which may be immobilized, to separate said polypeptide from the sample and the amount of polypeptide may then be determined.

"Fragments" of the polypeptides refers to a domain or region of said polypeptide, e.g. an antigenic
15 fragment, which is recognizable as being derived from said polypeptide to allow binding of a specific binding partner. Preferably such a fragment comprises a significant portion of said polypeptide and corresponds to a product of normal post-synthesis processing.

20 Thus in a further aspect the present invention provides a method of preparing a standard gene transcript pattern characteristic of a disease or condition or stage thereof in an organism comprising at least the steps of:

25 a) releasing target polypeptides from a sample of one or more organisms having the disease or condition or stage thereof;

b) contacting said target polypeptides with one or more binding partners, wherein each binding partner is
30 specific to a marker polypeptide (or a fragment thereof) encoded by the gene to which an oligonucleotide of Table 1 (or derived from a sequence described in Table 1) binds, to allow binding of said binding partners to said target polypeptides, wherein said marker polypeptides
35 are specific for said disease or condition thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and

- 34 -

c) assessing the target polypeptide binding to said binding partners to produce a characteristic pattern reflecting the level of gene expression of genes which express said marker polypeptides, in the sample with the disease, condition or stage thereof.

As used herein "target polypeptides" refer to those polypeptides present in a sample which are to be detected and "marker polypeptides" are polypeptides which are encoded by the genes to which Table 1 oligonucleotides or Table 1 derived oligonucleotides bind. The target and marker polypeptides are identical or at least have areas of high similarity, e.g. epitopic regions to allow recognition and binding of the binding partner.

"Release" of the target polypeptides refers to appropriate treatment of a sample to provide the polypeptides in a form accessible for binding of the binding partners, e.g. by lysis of cells where these are present. The samples used in this case need not necessarily comprise cells as the target polypeptides may be released from cells into the surrounding tissue or fluid, and this tissue or fluid may be analysed, e.g. urine or blood. Preferably however the preferred samples as described herein are used. "Binding partners" comprise the separate entities which together make an affinity binding pair as described above, wherein one partner of the binding pair is the target or marker polypeptide and the other partner binds specifically to that polypeptide, e.g. an antibody.

Various arrangements may be envisaged for detecting the amount of binding pairs which form. In its simplest form, a sandwich type assay e.g. an immunoassay such as an ELISA, may be used in which an antibody specific to the polypeptide and carrying a label (as described elsewhere herein) may be bound to the binding pair (e.g. the first antibody:polypeptide pair) and the amount of label detected.

- 35 -

Other methods as described herein may be similarly modified for analysis of the protein product of expression rather than the gene transcript and related nucleic acid molecules.

5 Thus a further aspect of the invention provides a method of preparing a test gene transcript pattern comprising at least the steps of:

- a) releasing target polypeptides from a sample of said test organism;
- 10 b) contacting said target polypeptides with one or more binding partners, wherein each binding partner is specific to a marker polypeptide (or a fragment thereof) encoded by the gene to which an oligonucleotide of Table 1 (or derived from a sequence described in Table 1)
- 15 binds, to allow binding of said binding partners to said target polypeptides, wherein said marker polypeptides are specific for said disease or condition thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and
- 20 c) assessing the target polypeptide binding to said binding partners to produce a characteristic pattern reflecting the level of gene expression of genes which express said marker polypeptides, in said test sample.

25 A yet further aspect of the invention provides a method of diagnosing or identifying or monitoring a disease or condition or stage thereof in an organism comprising the steps of:

- a) releasing target polypeptides from a sample of said organism;
- 30 b) contacting said target polypeptides with one or more binding partners, wherein each binding partner is specific to a marker polypeptide (or a fragment thereof) encoded by the gene to which an oligonucleotide of Table 1 (or derived from a sequence described in Table 1)
- 35 binds, to allow binding of said binding partners to said target polypeptides, wherein said marker polypeptides are specific for said disease or condition thereof in an

- 36 -

organism and sample thereof corresponding to the organism and sample thereof under investigation; and

c) assessing the target polypeptide binding to said binding partners to produce a characteristic pattern reflecting the level of gene expression of genes which express said marker polypeptides in said sample; and

d) comparing said pattern to a standard diagnostic pattern prepared as described hereinbefore using a sample from an organism corresponding to the organism and sample under investigation to determine the degree of correlation indicative of the presence of said disease or condition or a stage thereof in the organism under investigation.

The methods of generating standard and test patterns and diagnostic techniques rely on the use of informative oligonucleotide probes to generate the gene expression data. In some cases it will be necessary to select these informative probes for a particular method, e.g. to diagnose a particular disease, from a selection of available probes, e.g. the probes described hereinbefore (the Table 1 oligonucleotides, the Table 1 derived oligonucleotides, their complementary sequences and functionally equivalent oligonucleotides). The following methodology describes a convenient method for identifying such informative probes, or more particularly how to select a suitable sub-set of probes from the probes described herein.

Probes for the analysis of a particular disease or condition or stage thereof, may be identified in a number of ways known in the prior art, including by differential expression or by library subtraction (see for example WO98/49342). As described hereinafter, in view of the high information content of most transcripts, as a starting point one may also simply analyse a random sub-set of mRNA or cDNA species and pick the most informative probes from that sub-set. The following method describes the use of immobilized

- 37 -

oligonucleotide probes (e.g. the probes of the invention) to which mRNA (or related molecules) from different samples is bound to identify which probes are the most informative to identify a particular type of sample, e.g. a disease sample.

The immobilized probes can be derived from various unrelated or related organisms; the only requirement is that the immobilized probes should bind specifically to their homologous counterparts in test organisms. Probes can also be derived from commercially available or public databases and immobilized on solid supports or, as mentioned above, they can be randomly picked and isolated from a cDNA library and immobilized on a solid support.

The length of the probes immobilised on the solid support should be long enough to allow for specific binding to the target sequences. The immobilised probes can be in the form of DNA, RNA or their modified products or PNAs (peptide nucleic acids). Preferably, the probes immobilised should bind specifically to their homologous counterparts representing highly and moderately expressed genes in test organisms. Conveniently the probes which are used are the probes described herein.

The gene expression pattern of cells in biological samples can be generated using prior art techniques such as microarray or macroarray as described below or using methods described herein. Several technologies have now been developed for monitoring the expression level of a large number of genes simultaneously in biological samples, such as, high-density oligoarrays (Lockhart et al., 1996, Nat. Biotech., 14, p1675-1680), cDNA microarrays (Schena et al, 1995, Science, 270, p467-470) and cDNA macroarrays (Maier E et al., 1994, Nucl. Acids Res., 22, p3423-3424; Bernard et al., 1996, Nucl. Acids Res., 24, p1435-1442).

In high-density oligoarrays and cDNA microarrays,

- 38 -

hundreds and thousands of probe oligonucleotides or cDNAs, are spotted onto glass slides or nylon membranes, or synthesized on biochips. The mRNA isolated from the test and reference samples are labelled by reverse transcription with a red or green fluorescent dye, mixed, and hybridised to the microarray. After washing, the bound fluorescent dyes are detected by a laser, producing two images, one for each dye. The resulting ratio of the red and green spots on the two images provides the information about the changes in expression levels of genes in the test and reference samples. Alternatively, single channel or multiple channel microarray studies can also be performed.

In cDNA macroarray, different cDNAs are spotted on a solid support such as nylon membranes in excess in relation to the amount of test mRNA that can hybridise to each spot. mRNA isolated from test samples is radio-labelled by reverse transcription and hybridised to the immobilised probe cDNA. After washing, the signals associated with labels hybridising specifically to immobilised probe cDNA are detected and quantified. The data obtained in macroarray contains information about the relative levels of transcripts present in the test samples. Whilst macroarrays are only suitable to monitor the expression of a limited number of genes, microarrays can be used to monitor the expression of several thousand genes simultaneously and is, therefore, a preferred choice for large-scale gene expression studies.

A macroarray technique for generating the gene expression data set has been used to illustrate the probe identification method described herein. For this purpose, mRNA is isolated from samples of interest and used to prepare labelled target molecules, e.g. mRNA or cDNA as described above. The labelled target molecules are then hybridised to probes immobilised on the solid support. Various solid supports can be used for the

- 39 -

purpose, as described previously. Following hybridization, unbound target molecules are removed and signals from target molecules hybridizing to immobilised probes quantified. If radio labelling is performed, 5 PhosphoImager can be used to generate an image file that can be used to generate a raw data set. Depending on the nature of label chosen for labelling the target molecules, other instruments can also be used, for example, when fluorescence is used for labelling, a 10 FluoroImager can be used to generate an image file from the hybridised target molecules.

The raw data corresponding to mean intensity, median intensity, or volume of the signals in each spot can be acquired from the image file using commercially 15 available software for image analysis. However, the acquired data needs to be corrected for background signals and normalized prior to analysis, since, several factors can affect the quality and quantity of the hybridising signals. For example, variations in the 20 quality and quantity of mRNA isolated from sample to sample, subtle variations in the efficiency of labelling target molecules during each reaction, and variations in the amount of unspecific binding between different macroarrays can all contribute to noise in the acquired 25 data set that must be corrected for prior to analysis.

Background correction can be performed in several ways. The lowest pixel intensity within a spot can be used for background subtraction or the mean or median of the line of pixels around the spots' outline can be used 30 for the purpose. One can also define an area representing the background intensity based on the signals generated from negative controls and use the average intensity of this area for background subtraction.

35 The background corrected data can then be transformed for stabilizing the variance in the data structure and normalized for the differences in probe

- 40 -

intensity. Several transformation techniques have been described in the literature and a brief overview can be found in Cui, Kerr and Churchill
[http://www.jax.org/research/churchill/research/](http://www.jax.org/research/churchill/research/expression/Cui-Transform.pdf)
5 [expression/Cui-Transform.pdf](http://www.jax.org/research/churchill/research/expression/Cui-Transform.pdf)). Normalization can be performed by dividing the intensity of each spot with the collective intensity, average intensity or median intensity of all the spots in a macroarray or a group of spots in a macroarray in order to obtain the relative
10 intensity of signals hybridising to immobilised probes in a macroarray. Several methods have been described for normalizing gene expression data (Richmond and Somerville, 2000, *Current Opin. Plant Biol.*, 3, p108-116; Finkelstein et al., 2001, In "Methods of Microarray
15 Data Analysis. Papers from CAMDA, Eds. Lin & Johnson, Kluwer Academic, p57-68; Yang et al., 2001, In "Optical Technologies and Informatics", Eds. Bittner, Chen, Dorsel & Dougherty, *Proceedings of SPIE*, 4266, p141-152; Dudoit et al, 2000, *J. Am. Stat. Ass.*, 97, p77-87; Alter
20 et al 2000, *supra*; Newton et al., 2001, *J. Comp. Biol.*, 8, p37-52). Generally, a scaling factor or function is first calculated to correct the intensity effect and then used for normalising the intensities. The use of external controls has also been suggested for improved
25 normalization.

One other major challenge encountered in large-scale gene expression analysis is that of standardization of data collected from experiments performed at different times. We have observed that
30 gene expression data for samples acquired in the same experiment can be efficiently compared following background correction and normalization. However, the data from samples acquired in experiments performed at different times requires further standardization prior
35 to analysis. This is because subtle differences in experimental parameters between different experiments, for example, differences in the quality and quantity of

- 41 -

mRNA extracted at different times, differences in time used for target molecule labelling, hybridization time or exposure time, can affect the measured values. Also, factors such as the nature of the sequence of
5 transcripts under investigation (their GC content) and their amount in relation to the each other determines how they are affected by subtle variations in the experimental processes. They determine, for example, how efficiently first strand cDNAs, corresponding to a
10 particular transcript, are transcribed and labelled during first strand synthesis, or how efficiently the corresponding labelled target molecules bind to their complementary sequences during hybridization. Batch to batch difference in the printing process is also a major
15 factor for variation in the generated expression data.

Failure to properly address and rectify for these influences leads to situations where the differences between the experimental series may overshadow the main information of interest contained in the gene expression
20 data set, i.e. the differences within the combined data from the different experimental series. Figure 1 provides one such example showing a classification based on Principal Component Analysis (PCA) of combined data from two experimental series where the main goal is to
25 distinguish between Alzheimer/non-Alzheimer patients.

PCA (also known as singular value decomposition) is a technique for studying interdependencies and underlying relationships of a set of variables. The data are modelled in terms of a few significant factors
30 or principal components (PC's), plus residuals. The PC's contain the main phenomena and define the systematic variability present in the data, while the residuals represent the variability interpreted as noise. Details on PCA can be found in Jolliffe (1986,
35 Principal Component Analysis, Springer-Verlag, NY), and Jackson (1991, A User's Guide to Principal Components, Wiley, NY). The results of Figure 1 show that two

- 42 -

clusters are formed representing the data from two experimental series rather than the Alzheimer/non-Alzheimer differentiation. There were eight samples in common between the two series of experiments, which ideally should have fallen on top of, or in near proximity to, each other if appropriately standardized.

We have now found that gene expression data between different experiments can be efficiently standardized by including a subset of samples from one experimental series in the next experimental series and using a direct standardization method (DS), originally described by Wang and Kowalski (Anal. Chem., 1991, 63, p2750 and J. Chemometrics, 1991, 5, p129-145). Although the method of DS is well known in the field of analytical chemistry, it remains undescribed and unused in the field of gene expression data analysis.

In DS, the secondary data representing for example experimental series 2 (secondary measurements, R_2) are corrected to match the data measured on the primary measurements representing data from series 1 (R_1), while the calibration model remains unchanged. In DS, response matrices for both experimental series are related to each other by a transformation matrix F , i.e.

$$R_1 = R_2 F \quad (1)$$

Where F is a square matrix dimensioned gene by gene. From (1), the transformation matrix is calculated as:

$$F = R_2^+ R_1 \quad (2)$$

The transformation matrix F in equation (2) is calculated using a relatively small subset of samples which are measured on both the master primary and the secondary series of data.

Finally, the response of the unknown sample measured on the secondary series $r_{2,un}^T$, is standardized

- 43 -

to the response vector $\hat{r}_{1,un}^T$ expected from the primary series

5

$$\hat{r}_{1,un}^T = rT_{21,un}\hat{F} \quad (3)$$

From the preceding equation it can be seen that the column i of the transformation matrix contains the multiplication factors for a set of genes measured in the secondary series to obtain the intensity at spot i of the corrected series.

The number of samples that are repeated in the experimental series, R_1 and R_2 , should be equal to their ranks, which in this case is equal to the number of principal components retained for explaining the variation in the R_1 and R_2 . For example, if three principal components are retained for explaining the variation in the data set, a minimum of three samples should be repeated between R_1 and R_2 . The samples that should be repeated between different series should ideally be those that exhibit high leverages in the gene expression pattern. At times, two samples may suffice, while at other times, more than two samples should be ideally be included for good representativity. In some cases, the samples selected can be the same in all the experimental series to be compared (reference samples), while in other cases, representative samples can be selected sequentially by analyzing the expression pattern after each experiment. The selected samples with high leverages are then included in the next experimental series. The results of using Direct Standardization are shown in Figure 1.

Another approach for normalizing and standardizing the gene expression data set is to hybridize each DNA array with target molecules prepared from a test sample and an equal amount of labelled target molecules prepared from representative reference samples. In order to measure the intensity of labelled target

- 44 -

molecules hybridizing to the immobilized probes it is necessary that the labelled molecules are prepared from test and reference samples using different labels, for example, different fluorescent dyes can be used for
5 preparing the labelled material. The labelled molecules prepared from reference samples can be added to the hybridization solution together with the labelled material prepared from test samples. A data file from each array representing the expression pattern of
10 different genes in the test sample and reference samples can then be obtained, normalized and standardized by the direct standardization method as described above. An instant advantage of including the differentially labelled target molecules from reference samples during
15 hybridization is that it enables an efficient comparison of new test samples to the data sets already stored in a database.

Monitoring the expression of a large number of genes in several samples leads to the generation of a
20 large amount of data that is too complex to be easily interpreted. Several unsupervised and supervised multivariate data analysis techniques have already been shown to be useful in extracting meaningful biological information from these large data sets. Cluster
25 analysis is by far the most commonly used technique for gene expression analysis, and has been performed to identify genes that are regulated in a similar manner, and or identifying new/unknown tumour classes using gene expression profiles (Eisen et al., 1998, PNAS, 95,
30 p14863-14868, Alizadeh et al. 2000, supra, Perou et al. 2000, Nature, 406, p747-752; Ross et al, 2000, Nature Genetics, 24(3), p227-235; Herwig et al., 1999, Genome Res., 9, p1093-1105; Tamayo et al, 1999, Science, PNAS, 96, p2907-2912).

35 In the clustering method, genes are grouped into functional categories (clusters) based on their expression profile, satisfying two criteria: *homogeneity*

- 45 -

- the genes in the same cluster are highly similar in expression to each other; and *separation* - genes in different clusters have low similarity in expression to each other.

5 Examples of various clustering techniques that have been used for gene expression analysis include hierarchical clustering (Eisen et al., 1998, *supra*; Alizadeh et al. 2000, *supra*; Perou et al. 2000, *supra*; Ross et al, 2000, *supra*), K-means clustering (Herwig et
10 al., 1999, *supra*; Tavazoie et al, 1999, *Nature Genetics*, 22(3), p. 281-285), gene shaving (Hastie et al., 2000, *Genome Biology*, 1(2), research 0003.1-0003.21), block clustering (Tibshirani et al., 1999, Tech report Univ Stanford.) Plaid model (Lazzeroni, 2002, *Stat. Sinica*,
15 12, p61-86), and self-organizing maps (Tamayo et al. 1999, *supra*). Also, related methods of multivariate statistical analysis, such as those using the singular value decomposition (Alter et al., 2000, *PNAS*, 97(18), p10101-10106; Ross et al. 2000, *supra*) or
20 multidimensional scaling can be effective at reducing the dimensions of the objects under study.

 However, methods such as cluster analysis and singular value decomposition are purely exploratory and only provide a broad overview of the internal structure
25 present in the data. They are unsupervised approaches in which the available information concerning the nature of the class under investigation is not used in the analysis. Often, the nature of the biological
perturbation to which a particular sample has been
30 subjected is known. For example, it is sometimes known whether the sample whose gene expression pattern is being analysed derives from a diseased or healthy individual. In such instances, discriminant analysis can be used for classifying samples into various groups
35 based on their gene expression data.

 In such an analysis one builds the classifier by training the data that is capable of discriminating

- 46 -

between member and non-members of a given class. The trained classifier can then be used to predict the class of unknown samples. Examples of discrimination methods that have been described in the literature include

5 Support Vector Machines (Brown et al, 2000, PNAS, 97, p262-267), Nearest Neighbour (Dudoit et al., 2000, supra), Classification trees (Dudoit et al., 2000, supra), Voted classification (Dudoit et al., 2000, supra), Weighted Gene voting (Golub et al. 1999, supra),

10 and Bayesian classification (Keller et al. 2000, Tec report Univ of Washington). Also a technique in which PLS (Partial Least Square) regression analysis is first used to reduce the dimensions in the gene expression data set followed by classification using logistic

15 discriminant analysis and quadratic discriminant analysis (LD and QDA) has recently been described (Nguyen & Rocke, 2002, Bioinformatics, 18, p39-50 and 1216-1226).

A challenge that gene expression data poses to

20 classical discriminatory methods is that the number of genes whose expression are being analysed is very large compared to the number of samples being analysed. However in most cases only a small fraction of these genes are informative in discriminant analysis problems.

25 Moreover, there is a danger that the noise from irrelevant genes can mask or distort the information from the informative genes. Several methods have been suggested in literature to identify and select genes that are informative in microarray studies, for example,

30 t-statistics (Dudoit et al, 2002, J. Am. Stat. Ass., 97, p77-87), analysis of variance (Kerr et al., 2000, PNAS, 98, p8961-8965), Neighbourhood analysis (Golub et al, 1999, supra), Ratio of between groups to within groups sum of squares (Dudoit et al., 2002, supra), Non

35 parametric scoring (Park et al., 2002, Pacific Symposium on Biocomputing, p52-63) and Likelihood selection (Keller et al., 2000, supra).

- 47 -

In the methods described herein the gene expression data that has been normalized and standardized is analysed by using Partial Least Squares Regression (PLSR). Although PLSR is primarily a method used for regression analysis of continuous data (see Appendix A), it can also be utilized as a method for model building and discriminant analysis using a dummy response matrix based on a binary coding. The class assignment is based on a simple dichotomous distinction such as breast cancer (class 1) / healthy (class 2), or a multiple distinction based on multiple disease diagnosis such as breast cancer (class 1) / Alzheimer (class 2) / healthy (class 3). The list of diseases for classification can be increased depending upon the samples available corresponding to other diseases or conditions or stages thereof.

PLSR applied as a classification method is referred to as PLS-DA (DA standing for Discriminant analysis). PLS-DA is an extension of the PLSR algorithm in which the Y-matrix is a dummy matrix containing n rows (corresponding to the number of samples) and K columns (corresponding to the number of classes). The Y-matrix is constructed by inserting 1 in the k th column and -1 in all the other columns if the corresponding i th object of X belongs to class k . By regressing Y onto X , classification of a new sample is achieved by selecting the group corresponding to the largest component of the fitted, $\hat{y}(x) = (\hat{y}_1(x), \hat{y}_2(x), \dots, \hat{y}_k(x))$. Thus, in a -1/1 response matrix, a prediction value below 0 means that the sample belongs to the class designated as -1, while a prediction value above 0 implies that the sample belongs to the class designated as 1.

An advantage of PLSR-DA is that the results obtained can be easily represented in the form of two different plots, the score and loading plots. Score plots represent a projection of the samples onto the principal components and shows the distribution of the

- 48 -

samples in the classification model and their relationship to one another. Loading plots display correlations between the variables present in the data set.

5 It is usually recommended to use PLS-DA as a starting point for the classification problem due to its ability to handle collinear data, and the property of PLSR as a dimension reduction technique. Once this purpose has been satisfied, it is possible to use other
10 methods such as Linear discriminant analysis, LDA, that has been shown to be effective in extracting further information, Indahl et al. (1999, Chem. and Intell. Lab. Syst., 49, p19-31). This approach is based on first decomposing the data using PLS-DA, and then using the
15 scores vectors (instead of the original variables) as input to LDA. Further details on LDA can be found in Duda and Hart (Classification and Scene Analysis, 1973, Wiley, USA).

 The next step following model building is of model
20 validation. This step is considered to be amongst the most important aspects of multivariate analysis, and tests the "goodness" of the calibration model which has been built. In this work, a cross validation approach has been used for validation. In this approach, one or
25 a few samples are kept out in each segment while the model is built using a full cross-validation on the basis of the remaining data. The samples left out are then used for prediction/classification. Repeating the simple cross-validation process several times holding
30 different samples out for each cross-validation leads to a so-called double cross-validation procedure. This approach has been shown to work well with a limited amount of data, as is the case in some of the Examples described here. Also, since the cross validation step
35 is repeated several times the dangers of model bias and overfitting are reduced.

 Once a calibration model has been built and

- 49 -

validated, genes exhibiting an expression pattern that is most relevant for describing the desired information in the model can be selected by techniques described in the prior art for variable selection, as mentioned
 5 elsewhere. Variable selection will help in reducing the final model complexity, provide a parsimonious model, and thus lead to a reliable model that can be used for prediction. Moreover, use of fewer genes for the purpose of providing diagnosis will reduce the cost of
 10 the diagnostic product. In this way informative probes which would bind to the genes of relevance may be identified.

We have found that after a calibration model has been built, statistical techniques like Jackknife
 15 (Effron, 1982, The Jackknife, the Bootstrap and other resampling plans. Society for Industrial and Applied mathematics, Philadelphia, USA), based on resampling methodology, can be efficiently used to select or confirm significant variables (informative probes).

20 The approximate uncertainty variance of the PLS regression coefficients B can be estimated by:

$$S^2B = \sum_{m=1}^M ((B-B_m)g)^2$$

25

where

S^2B = estimated uncertainty variance of B ;

B = the regression coefficient at the cross validated
 30 rank A using all the N objects;

B_m = the regression coefficient at the rank A using all objects except the object(s) left out in cross validation segment m ; and

g = scaling coefficient (here: $g=1$).

35 In our approach, Jackknife has been implemented together with cross-validation. For each variable the difference between the B -coefficients B_i in a

- 50 -

cross-validated sub-model and B_{tot} for the total model is first calculated. The sum of the squares of the differences is then calculated in all sub-models to obtain an expression of the variance of the B_i estimate for a variable. The significance of the estimate of B_i is calculated using the t-test. Thus, the resulting regression coefficients can be presented with uncertainty limits that correspond to 2 Standard Deviations, and from that significant variables are detected.

No further details as to the implementation or use of this step are provided here since this has been implemented in commercially available software, The Unscrambler, CAMO ASA, Norway. Also, details on variable selection using Jackknife can be found in Westad & Martens (2000, J. Near Inf. Spectr., 8, p117-124).

The following approach can be used to select informative probes from a gene expression data set:

- a) keep out one unique sample (including its repetitions if present in the data set) per cross validation segment;
- b) build a calibration model (cross validated segment) on the remaining samples using PLSR-DA;
- c) select the significant genes for the model in step b) using the Jackknife criterion;
- d) repeat the above 3 steps until all the unique samples in the data set are kept out once (as described in step a). For example, if 75 unique samples are present in the data set, 75 different calibration models are built resulting in a collection of 75 different sets of significant probes;
- e) select the most significant variables using the frequency of occurrence criterion in the generated sets of significant probes in step d). For example, a set of probes appearing in all sets (100%) are more informative than probes appearing in only 50% of the

- 51 -

generated sets in step d).

Once the informative probes for a disease have been selected, a final model is made and validated. The two most commonly used ways of validating the model are

5 cross-validation (CV) and test set validation. In cross-validation, the data is divided into k subsets. The model is then trained k times, each time leaving out one of the subsets from training, but using only the omitted subset to compute error criterion, RMSEP (Root

10 Mean Square Error of Prediction). If k equals the sample size, this is called "leave-one-out" cross-validation. The idea of leaving one or a few samples out per validation segment is valid only in cases where the covariance between the various experiments is zero.

15 Thus, one sample at-a-time approach can not be justified in situations containing replicates since keeping only one of the replicates out will introduce a systematic bias in our analysis. The correct approach in this case will be to leave out all replicates of the same samples

20 at a time since that would satisfy assumptions of zero covariance between the CV-segments.

The second approach for model validation is to use a separate test-set for validating the calibration model. This requires running a separate set of

25 experiments to be used as a test set. This is the preferred approach given that real test data are available.

The final model is then used to identify a disease, condition or stage thereof in test samples. For this

30 purpose, expression data of selected informative genes is generated from test samples and then the final model is used to determine whether a sample belongs to a diseased or non-diseased class or has a condition or stage thereof.

35 Thus viewed from a yet further aspect the present invention provides a method of identifying probes useful for diagnosing or identifying or monitoring a disease or

- 52 -

condition or stage thereof in an organism, comprising the steps of:

- 5 a) immobilizing a set of oligonucleotide probes, preferably as described hereinbefore, on a solid support;
- b) isolating mRNA from a sample of a normal organism (normal sample), which may optionally be reverse transcribed to cDNA;
- 10 c) isolating mRNA from a sample from an organism, corresponding to the sample and organism of step (b), which is known to have said disease or condition or a stage thereof (diseased sample), which may optionally be reverse transcribed to cDNA;
- 15 d) hybridizing the mRNA or cDNA of steps (b) and (c) to said set of immobilized oligonucleotide probes of step (a); and
- e) assessing the amount of mRNA or cDNA hybridizing to each of said oligonucleotide probes to determine the level of gene expression of genes to which said oligonucleotide probes bind in said normal and diseased samples to generate a gene expression data set for each sample;
- 20 f) normalizing and standardizing said data set of step (e);
- g) constructing a calibration model for classification, preferably using the statistical techniques Partial Least Squares Discriminant Analysis (PLS-DA) and Linear Discriminant Analysis (LDA);
- 30 h) performing JackKnife analysis and identifying those oligonucleotide probes which are required for classification of said disease and normal samples into their respective groups.
- 35

Preferably a model for classification purposes is

- 53 -

generated by using the data relating to the probes identified according to the above described method. Preferably the sample is as described previously. Preferably the oligonucleotides which are immobilized in
5 step (a) are randomly selected as described below or are the probes as described hereinbefore. Such oligonucleotides may be of considerable length, e.g. if using cDNA (which is encompassed within the scope of the term "oligonucleotide"). The identification of such
10 cDNA molecules as useful probes allows the development of shorter oligonucleotides which reflect the specificity of the cDNA molecules but are easier to manufacture and manipulate.

The above described model may then be used to
15 generate and analyse data of test samples and thus may be used for the diagnostic methods of the invention. In such methods the data generated from the test sample provides the gene expression data set and this is normalized and standardized as described above. This is
20 then fitted to the calibration model described above to provide classification.

The method described herein can also be used to simultaneously select informative probes for several related and unrelated diseases or conditions. Depending
25 upon which diseases or conditions have been included in the calibration or training set, informative probes can be selected for the said diseases or conditions. The informative probes selected for one disease or condition may or may not be similar to the informative probes
30 selected for another disease or condition of interest. It is the pattern with which the selected genes are expressed in relation to each other during a disease, condition, or stage thereof, that determines whether or not they are informative for the disease, condition or
35 stage thereof.

In other words, informative genes are selected based on how their expression correlates with the

- 54 -

expression of other selected informative genes under the influence of responses generated by the disease, condition or stage thereof under investigation. In examples 1 and 2 provided hereinafter, 139 informative probes were selected for breast cancer diagnosis and 182 probes were selected for Alzheimer's disease diagnosis by training the gene expression data set of genes representing 1435 or 758 randomly picked cDNA clones for breast cancer/non breast cancer samples, or Alzheimer/non-Alzheimer samples, respectively. Among the probes selected for breast cancer and Alzheimer, about 10 probes were informative both for breast cancer and Alzheimer disease diagnosis.

For the purpose of isolating informative probes or identifying several related and unrelated diseases, conditions and stages thereof simultaneously, the gene expression data set must contain the information on how genes are expressed when the subject has a particular disease, condition or stage thereof under investigation. The data set is generated from a set of healthy or diseased samples, where a particular sample may contain the information of only one disease, condition or stages thereof or may also contain information about multiple diseases, conditions or stages thereof. For example, if the isolation of informative probes for Alzheimer disease, breast cancer and diabetes is sought, whole blood samples can be obtained from an Alzheimer patient who has breast cancer and diabetes. Hence, the method also teaches an efficient experimental design to reduce the number of samples required for isolating informative probes by selecting samples representing more than one disease, condition or stage thereof.

As mentioned previously, in view of the high information content of most transcripts, the identification and selection of informative probes for use in diagnosing, monitoring or identifying a particular disease, condition or stage thereof may be

- 55 -

dramatically simplified. Thus the pool of genes from which a selection may be made to identify informative probes may be radically reduced.

Unlike, in prior art technologies where informative probes are selected from a population of thousands of genes that are being expressed in a cell, like in microarray, in the method described herein, the informative probes are selected from a limited number of randomly obtained genes. For example, from a population of 1435 cDNA clones, randomly picked from a human whole blood cDNA library, we were able to select 139 informative probes for breast cancer diagnosis (see Example 1 and Table 2).

Thus in a preferred aspect of the above mentioned method of identifying probes useful for diagnosing or identifying or monitoring a disease or condition or stage thereof in an organism, said set of oligonucleotides which are immobilized in step (a) are randomly selected from a larger set of oligonucleotides, e.g. from a cDNA library or other oligonucleotide pool, which may be, but is preferably not selected from the set provided herein. Preferably said larger set comprises oligonucleotides which correspond to moderately or highly expressed genes. Thus preferably in methods of the invention, the set of oligonucleotides according to the invention are replaced with a set of oligonucleotides which are randomly selected, e.g. from commercially available oligonucleotide or cDNA libraries.

As referred to herein "random" refers to selection which is not biased based on the extent of information carried by the transcripts in relation to the disease, condition or organism under study, ie. without bias towards their likely utility as informative probes. Whilst a random selection may be made from a pool of transcripts (or related products) which have been biased, e.g. to highly or moderately expressed

- 56 -

transcripts, preferably random selection is made from a pool of transcripts not biased or selected by a sequence-based criterion. The larger set may therefore contain oligonucleotides corresponding to highly and moderately expressed genes, or alternatively, may be enriched for those corresponding to the highly and moderately expressed genes.

Random selection from highly and moderately expressed genes can be achieved in a wide variety of ways. A strategy used in this work, but not limiting in itself involves randomly picking a significant number of cDNA clones from a cDNA library constructed from a biological specimen under investigation. Since, in a cDNA library, the cDNA clones corresponding to transcripts present in high or moderate amount are more frequently present than transcripts corresponding to cDNA present in low amount, the former will tend to be picked up more frequently than the latter. A pool of cDNA enriched for those corresponding to highly and moderately expressed genes can be isolated by this approach.

To identify genes that are expressed in high or moderate amount among the isolated population for use in methods of the invention, the information about the relative level of their transcripts in samples of interest can be generated using several prior art techniques. Both non-sequence based methods, such as differential display or RNA fingerprinting, and sequence-based methods such as microarrays or macroarrays can be used for the purpose. Alternatively, specific primer sequences for highly and moderately expressed genes can be designed and methods such as quantitative RT-PCR can be used to determine the levels of highly and moderately expressed genes. Hence, a skilled practitioner may use a variety of techniques which are known in the art for determining the relative level of mRNA in a biological sample.

- 57 -

Especially preferably the sample for the isolation of mRNA in the above described method is as described previously and is preferably not from the site of disease and the cells in said sample are not disease cells and have not contacted disease cells.

The following examples are given by way of illustration only in which the Figures referred to are as follows:

Figure 1 shows the effect of Direct Standardization (DS) on the Alzheimer data measured in two different series of experiments in which AD denotes Alzheimer's samples and A,B are non-Alzheimer's samples. The samples in both series have been labelled systematically as (xx_7/xx_8), whereas the corrected samples from series 8 (in b,c,d) have been labelled as (xx_c), thus, for example, AD2-7 denotes Alzheimer disease sample number 2 in experiment series 7. The circled spots represent the samples chosen as the transfer samples. The connecting lines in figures b,c,d show the proximity of the replicated samples after applying DS. The dashed lines in figures a,c,d represent the decision boundary separating the classes. These lines have not been drawn on the basis of any statistical criteria, but serve the purpose of visually separating the classes. All the four figures show scores plot (PC1-PC2) from PCA analysis based on (a) non-standardized data, (b) scores plot after direct standardization using 3 transfer samples, (c) scores plot after direct standardization using 4 transfer sample, (d) scores plot after direct standardization using 8 transfer samples;

Figure 2 shows the projection of normal (including benign) and breast cancer samples onto a classification model generated by PLSR-DA using the data of 44 informative genes, in which PC is the principal components and N and C are normal and breast cancer samples, respectively;

Figure 3 shows the projection of individuals with

- 58 -

and without Alzheimer's disease onto a classification model generated by PLSR-DA using 182 informative genes;

Figures 4, 6 and 8 show projection plots as Figure 2 in which the classification model is generated using 719, 111 and 345 cDNAs, respectively, wherein PC is the principal components, N denotes normal and B denotes breast cancer samples;

Figures 5, 7 and 9 show prediction plots based on 3 principal components using the data of 719, 111 and 345 cDNAs, respectively;

Figure 10 shows a projection plot as Figure 3 in which the classification model is generated using 520 cDNAs; and

Figure 11 is the prediction plot corresponding to Figure 10.

Example 1: Diagnosis of Breast Cancer

Methods

Whole blood was obtained from the arms of breast cancer patients and patients with benign tumours (Ullevål and Haukland hospitals in Norway). All of the patients with breast cancer had a malignant tumour of the breast (disease samples). Healthy blood was collected from the above two hospitals, or collected at a Health station at Ås, Norway or at DiaGenic AS, Norway, from the arms of female donors with no reported signs of breast cancer. The blood from healthy individuals or with benign tumours comprise the normal samples. The blood was either collected in tubes containing EDTA and stored immediately at -80°C or was collected in PAXgene tubes and stored for 12-24 hours at room temperature before finally storing them at -80°C before use. Further details of the breast cancer and benign tumour patients from which blood was taken is provided in Table 5.

- 59 -

mRNA was isolated from the blood of the 29 breast cancer patients and 46 normal donors and used to prepare labelled probes by reverse transcribing in the presence of $\alpha^{33}\text{P}$ -dATP. The first strand cDNA of the normal and diseased samples was bound, separately to 1435 cDNA clones immobilized on a solid support (nylon membrane). These cDNA clones were randomly picked, without any prior knowledge of their gene sequences, from a cDNA library constructed using whole blood of 550 healthy individuals (Clontech, Palo Alto, USA). These methods were conducted as follows.

For amplification of inserts, bacterial clones were grown in microtiter plates containing 150 μl LB with 50 $\mu\text{g/ml}$ carbenicillin, and incubated overnight with agitation at 37°C. To lyse the cells, 5 μl of each culture were diluted with 50 μl H₂O and incubated for 12 min. at 95°C. Of this mixture, 2 μl were subjected to a PCR reaction using 20 pmoles of M13 forward and reverse primer in presence of 1.5 mM MgCl₂. PCR reactions were performed with the following cycling protocol: 4 min. at 95°C, followed by 25 cycles of 1 min. at 94°C, 1 min. at 60°C and 3 min. at 72°C either in a RoboCycler® Temperature Cycler (Stratagene, La Jolla, USA) or DNA Engine Dyad Peltier Thermal Cycler (MJ Research Inc., Waltham, USA). The amplified products were denatured by incubating with NaOH (0.2 M, final concentration) for 30 min. and spotted onto Hybond-N+ membranes (Amersham Pharmacia Biotech, Little Chalfont, UK), using MicroGrid II workstation according to the manufacturer's instructions (BioRobotics Ltd, Cambridge England). The immobilized cDNAs were fixed using a UV cross-linker (Hoefer Scientific Instruments, San Francisco, USA).

In addition to the 1435 cDNAs, the printed arrays also contained controls for assessing background level, consistency and sensitivity of the assay. These were

- 60 -

spotted at multiple positions and included controls such as PCR mix (without any insert); positive and negative controls of SpotReport™ 10 array validation system (Stratagene, La Jolla, USA) and cDNAs corresponding to constitutively expressed genes such as b-actin, g-actin, GAPDH, HOD and cyclophilin. Also, oligonucleotides corresponding to SIX1, b-tubulin, TRP-2, MDM2, Myosin Light C, CD44, Maspin, Laminin, and SRP 19 were included to detect disseminated cancer cells.

The total RNA from blood collected in EDTA tubes was purified using Trizol LS Reagent protocol (Invitrogen/Life Technologies). From blood contained in PAXgene tubes, the total RNA was purified according to the supplier's instructions (PreAnalytiX, Hombrechtikon, Switzerland). Contaminating DNA was removed from the isolated RNA by DNAase I treatment using DNA-free kit (Ambion, Inc. Austin, USA). RNA quality was determined visually by inspecting the integrity of 28S and 18S ribosomal bands following agarose gel electrophoresis. The concentration and purity of extracted RNA was determined by measuring the absorbance at 260 nm and 280 nm. mRNA was isolated from the total RNA using Dynabeads as per the supplier's instructions (DynaL AS, Oslo, Norway).

Labelling and hybridization experiments were performed in batches. The number of samples assayed in each batch varied from six to nine. In the case of samples that were assayed more than once (replicates), aliquots derived from the same mRNA pool were used for probe synthesis. For probe synthesis, aliquots of mRNA corresponding to 4-5 µg of total RNA were mixed together with oligodT_{25NV} (0.5 µg/ml) and mRNA spikes of SpotReport™ 10 array validation system (10 pg; Spike 2, 1 pg), heated to 70°C to remove secondary structures, and then chilled on ice. Probes were prepared in 35µl

- 61 -

reaction mixes by reverse transcription in the presence of 50 μ Ci [α^{33} P] dATP, 3.5 μ M dATP, 0.6 mM each of dCTP, dTTP, dGTP, 200 units of SuperScript reverse transcriptase (Invitrogen, LifeTechnologies) and 0.1 M DTT, labelling for 1.5 hr at 42°C. Following synthesis, the enzyme was deactivated for 10 min. at 70°C and mRNA removed by incubating the reaction mix for 20 min. at 37°C in 4 units of Ribo H (Promega, Madison USA). Unincorporated nucleotides were removed using ProbeQuant G 50 Columns (Amersham Biosciences, Piscataway, USA).

Prior to hybridization, the membranes were equilibrated in 4 x SSC for 2 hr at room temperature and prehybridized overnight at 65°C in 10 ml prehybridisation solution (4 x SSC, 0.1 M NaH₂PO₄, 1 mM EDTA, 8% dextran sulphate, 10 x denhardt's solution, 1% SDS). Freshly prepared probes were added to 5 ml of the same prehybridisation solution, and hybridization continued overnight at 65°C. The membranes were washed at 65°C at increasing stringency (2 x 30 min. each in 2 x SSC, 0.1% SDS; 1 x SSC, 0.1% SDS; 0.1 x SSC, 0.1% SDS) to remove unspecific signals.

The amount of labelled first strand cDNA binding to each spot was assessed and quantified using a PhosphorImager to generate a gene expression data set. The data was generated using Phoretix software version 3 (Non Linear Dynamics, England). Background subtraction was performed on the generated data by subtracting the median of the line of pixels around each spot outline from the total intensity obtained from the respective spots.

The background-subtracted data was then normalized and transformed by selecting out 50 lowest and 50 maximum signals from each membrane. This step was to exclude genes that were expressed with a high degree of

- 62 -

variance. Since the genes varied from membrane to membrane, the expression data from 497 genes were removed from the data set. The values for the remaining 938 genes were then normalised by using different approaches such as external controls, dividing each spot by the median intensity of the observed signal in the respective membrane, range normalizing the data from each membrane, and then log transforming the data obtained.

10

The processed data obtained above was then used to isolate the informative probes by:

- a) keeping one unique sample (including all repetitions of the selected sample) out per cross validation segment;
- b) building a calibration model (cross validated) on the remaining samples using PLSR-DA;
- c) selecting the set of significant genes for the model in step b using the Jackknife criterion;
- d) repeating steps a), b) and c) until all the unique samples were kept out once (hence, in all 75 different calibration models were built (after repeating step b) 75 times), resulting in 75 different sets of significant probes (after repeating step c) 75 times));
- e) selecting significant variables using the frequency of occurrence criterion amongst the 75 different sets of significant probes.

25

The selected informative probes based on occurrence criterion were used to construct a classification model. The result of the classification model based on probes appearing in at least 90% of the generated sets after the step of isolating informative probes as described above is shown in Figure 2 in which it is seen that the expression pattern of these genes was able to classify most women with breast cancer and women with no breast cancer into distinct groups. In this figure PC1 and PC2

30

35

- 63 -

indicate the two principal components statistically derived from the data which best define the systemic variability present in the data. This allows each sample, and the data from each of the informative probes to which the sample's labelled first strand cDNA was bound, to be represented on the classification model as a single point which is a projection of the sample onto the principal components - the score plot.

The ability of the generated model, based on isolated informative probes, to predict future samples was determined by the double cross-validation approach. The performance of the diagnostic test for breast cancer based on the occurrence criterion is presented in Table 6.

Correct prediction of most breast cancer cells was achieved. These included all three samples obtained from women with ductal carcinoma in situ (DCIS), 11/15 samples obtained from women with stage I breast cancer, all five samples obtained from women with stage II breast cancer, and one of two samples obtained from women with stage III breast cancer. Interestingly, two correctly predicted stage I samples were obtained from women having a tumour size of <5 mm in diameter.

The model also correctly predicted the class of most non-cancer samples (41/46), including those that were obtained from women with non-cancerous breast abnormalities.

Confirmation that the gene transcripts are not from cells which are disseminated disease cells has been confirmed by several lines of evidences. Firstly, the informative genes were expressed constitutively at high or moderate levels in blood cells of women irrespective of whether they had cancer or not. Secondly, in the

- 64 -

assay described in this Example, in order to identify transcripts, at least 720 disseminated cells in blood samples would be required. Since, the average number of disseminated cells present in blood during different stages of breast cancer is much lower (organ confined breast cancer, 0.8 cells per ml; invasive breast cancer spread to lymph nodes only, 2.4 cells per ml; and metastatic breast cancer, 6 cells per ml; SD>100%) (29), we believe that the signals being detected originated from peripheral blood cells and could not have originated from disseminated cells. Thirdly, we were not able to detect any signal from the eight cancer markers known to have elevated expression in malignant cancer cells, including cancer cells that are disseminated in the blood.

Example 2: Diagnosis of Alzheimer's disease

Similar experiments were conducted with samples from Alzheimer's patients. In this method 7 patients diagnosed with Alzheimer's Disease at the Memory Clinic at Ullevål University Hospital were used in the trial. The patients were confirmed as having Alzheimer's disease based on the following criteria:

- * A standardized interview with a care-giver using IQCODE, an ADL scale and a scale measuring behaviour of the patient (Green scale).
- * Neuropsychological evaluation using MMSE, Clock drawing test, Trailmaking test A and B (TMT A and B), Kendrick object learning test (visual memory test), part of the Wechsler battery and Benton test.
- * A psychiatric evaluation using scales for detection of depression, MADRS for interviewing the patient and Cornell scale for interviewing the care-giver.
- * A physical examination.
- * Laboratory tests of blood samples to rule out other

- 65 -

diseases.

* CT scan of the brain.

* SPECT of the brain.

5 The mean age of the patients was 72.3 with an age range of 69-76. The mean MMSE score was 22.0 (the maximum score attainable being 30).

Six age-matched individuals without diagnosed
10 Alzheimer's disease were used as a control. All had been tested with MMSE and had a minimum score of 28 (mean: 28.4). The mean age of the normal control group was 73.0 and the age range 66-81. A sample from a 16-
15 year old individual, with a consequent minimal chance of having Alzheimer's disease, was also included as an additional control.

Using the methods described above (except that
hybridization to 758 rather than 1435 cDNA clones was
20 performed), informative probes were selected based on occurrence criterion and used to construct a classification model. The results of the classification model based on probes appearing at least once in the generated sets after the method to isolate informative
25 probes as described above is shown in Figure 3 in which it will be seen that the expression pattern of these genes was able to classify individuals with or without Alzheimer's disease into distinct groups. In this Figure PC1 and PC2 indicate the 2 principal components
30 statistically derived from the data which define the systematic variability present in the data. This allows each sample, and the data from each of the informative probes to which the samples' cDNA was bound, to be represented on the classification model as a single
35 point which is a projection of the sample onto the principal components - the score plot.

- 66 -

The ability of the generated model, based on isolated informative probes, to predict future samples was determined by the double cross-validation. The performance of the diagnostic test for Alzheimer's disease is presented in Table 7.

5

Appendix A

Partial Least Squares regression (PLSR)

5 Let a multivariate regression model be defined as:

$$Y = XB + F$$

where

10 X a $N \times P$ matrix with N predictor variables (genes);
 Y ($N \times J$) being the J predicted variables. In our case Y
 represents a matrix containing dummy variables;
 B is a matrix of regression coefficients; and
 F is a $N \times J$ matrix of residuals.

15

The structure of the PLSR model can be written as:

$$X = TP^T + E_A, \text{ and}$$

$$Y = TQ^T + F_A, \text{ where}$$

20

where

T ($N \times A$) is a matrix of score vectors which are linear
 combinations of the x-variables;

25 P ($P \times A$) is a matrix with the x-loading vectors p_a as
 columns;

Q ($J \times A$) is a matrix with the y-loading vectors q_a as
 columns;

E_a ($N \times P$) is the matrix for X after A factors; and

F_a ($N \times J$) is the matrix for Y after A factors.

30

The criterion in PLSR is to maximize the explained
 covariance of $[X, Y]$. This is achieved by the loading
 weights vector w_{a+1} , which is the first eigenvector of
 $E_a^T F_a F_a^T E_a$ (E_a and F_a are the deflated X and Y after a
 35 factors or PLS components).

The regression coefficients are given by:

- 68 -

$$B = W(P^TW)^{-1}Q^T$$

A PLSR model with full rank, i.e. maximum number of components, is equivalent to the MLR solutions. Further
5 details on PLSR can be found in Martens & Naes, 1989, Multivariate Calibration, John Wiley & Sons, Inc., USA and Kowalski & Seasholtz, 1991, supra.

- 69 -

Example 3: Validation of Example 1, diagnosis of breast cancer

The results in Example 1 were validated by using the
5 informative probes identified in Example 1 on new breast
cancer and control samples.

Methods

The methods, essentially as described in Example 1, were
10 used. Blood was taken from patients as described in
Table 8. However, blood was collected in PAXgene tubes
and the first strand labelled cDNAs were hybridized to
719 cDNAs spotted on nylon membranes along with other
controls as described in Example 1. After background
15 subtraction using control spots, the data of each
membrane was normalized using the inter quantile range.
The data was analysed as described in Example 1 and the
model validated by cross validation.

20 The 719 cDNAs which were spotted are a subset of the
cDNAs spotted in Example 1 and include 111 cDNAs
described in Table 2 and which were found to be
informative in Example 1.

Results

25 The results are shown in Figures 4 to 9. Figures 4, 6
and 8 are projection plots similar to Figure 2 and show
the projection of normal and breast cancer patients'
samples onto a classification model generated using all
30 719 cDNA. Figure 6 is similar but uses a classification
model generated with the 111 probes common to Example 1.
Figure 8 uses the 345 sequences of the 719 for which
sequence information is provided herein. In each case
classification of normal and breast cancer groups was
35 possible. Figures 5, 7 and 9 show prediction plots
which reflect the ability of the generated models to
correctly diagnose breast cancer. In the 3 prediction

- 70 -

plots shown, the disease samples appear on the x axis at +1 and the non-disease samples appear at -1. The y axis represents the predicted class membership. During prediction, if the prediction is correct, disease samples should fall above zero and non-disease samples should fall below zero. In each case almost all samples are correctly predicted.

Example 4: Validation of Example 2, diagnosis of
Alzheimers

The results in Example 2 were validated by using the informative probes identified in Example 2 on new Alzheimer's patient samples.

Methods

The methods, essentially as described in Example 2, were used. Twelve female patients diagnosed with Alzheimer's disease at the Memory Clinic at Ullevål University Hospital who were confirmed as having Alzheimer's disease based on the criteria of Example 2 were used in the trial. The mean age of the patients was 72.3 with an age range of 66-83. The mean MMSE score was 22.0 (the maximum score attainable being 30).

Sixteen age-matched female individuals without diagnosed Alzheimer's disease were used as the normal control group. All had been tested with MMSE and had a minimum score of 29. The mean age of the normal control group was 74.0 and the age range 66-86.

After transfer of the blood to PAXgene tubes, total mRNA was isolated from the blood of the Alzheimer's disease and from the control group donors according to the manufacturers's instructions (PreAnalytiX, Hombrechtikon, Switzerland). The isolated mRNA was labelled during reverse transcription in the presence of

- 71 -

$\alpha^{33}\text{P}$ -dATP, yielding a labelled first strand cDNA.
Hybridization was performed as described previously onto
730 cDNA clones picked from a cDNA library from whole
blood of 550 healthy individuals without knowledge of
5 the gene sequence of the random cDNA clones.

Results

The results are shown in Figures 10 and 11. Figure 10
is a projection plot generated using 520 probes which
10 have been sequenced. Figure 11 is a prediction plot and
shows correct prediction of almost all samples.

Table 1a

List of probes informative for disease diagnosis

	Clone ID	Sequence ID	No. of nucleotides
1	I-01	-	-
2	I-02	-	-
3	I-13	-	-
4	I-21	-	-
5	I-24	308	373
6	I-28	310	564
7	I-30	1180	622
8	I-34	313	554
9	I-37	-	-
10	I-42	-	-
11	I-52	-	-
12	I-54	1181	155
13	I-58	326	554
14	I-71	-	-
15	I-72	-	-
16	I-86	-	-
17	I-95	-	-
18	II-03	361	622
19	II-05	363	628
20	II-06	364	528
21	II-10	368	329
22	II-24	381	534
23	II-25	382	444
24	II-26	383	566
25	II-33	390	523
26	II-34	391	566
27	II-41	397	534
28	II-42	398	512
29	II-47	-	-
30	II-57	411	505
31	II-61	415	596
32	II-69	423	387
33	II-70	424	420
34	II-75	429	535
35	II-83	-	-
36	II-84	438	577
37	II-87	441	552
38	II-88	442	606
39	II-90	-	-
40	II-94	448	329
41	III-02	453	747
42	III-05	-	-
43	III-06	458	682

44	III-08	460	536
45	III-10	-	-
46	III-13	464	615
47	III-15	-	-
48	III-17	-	-
49	III-20	1183	479
50	III-23	473	694
51	III-26	476	476
52	III-35	485	551
53	III-39	487	224
54	III-40	488	349
55	III-43	490	382
56	III-44	491	382
57	III-53	500	390
58	III-56	503	109
59	III-57	504	374
60	III-60	-	-
61	III-60	-	-
62	III-61	507	521
63	III-63	509	575
64	III-68	-	-
65	III-74	518	502
66	III-80	523	585
67	III-82	-	-
68	III-85	526	516
69	III-89	530	660
70	III-92	-	-
71	III-96	-	-
72	IV-14	684	545
73	IV-15	1185	628
74	IV-23	-	-
76	IV-26	1186	494
75	IV-26	-	-
77	IV-29	-	-
78	IV-31	687	268
79	IV-32	688	569
80	IV-34	-	-
81	IV-35	-	-
82	IV-41	-	-
83	IV-45	-	-
84	IV-53	61	362
85	IV-62	-	-
86	IV-69	192	286
87	IV-80	701	579
88	IV-82	-	-
89	IV-93	-	-
90	IX-10	736	641
91	IX-12	-	-
92	IX-38	757	583
93	IX-39	758	424
94	IX-42	-	-
95	IX-48	764	626
96	IX-77	785	556
97	V-01	-	-
98	V-02	-	-
99	V-03	706	496

100	V-04	707	397
101	V-06	-	-
102	V-07	708	293
103	V-11	1188	599
104	V-12	711	498
105	V-15	-	-
106	V-17	-	-
107	V-21	-	-
108	V-25	-	-
109	V-32	-	-
110	V-35	-	-
111	V-39	-	-
112	V-42	-	-
113	V-43	-	-
114	V-47	-	-
115	V-49	-	-
116	V-52	-	-
117	V-54	-	-
118	V-55	77	412
119	V-58	-	-
120	V-59	-	-
121	V-65	-	-
122	V-68	-	-
123	V-71	-	-
124	V-75	-	-
125	V-79	-	-
126	V-80	726	260
127	V-90	-	-
128	V-91	-	-
129	V-92	-	-
130	V-94	-	-
131	VI-02	-	-
132	VI-04	865	122
133	VI-07	93	405
134	VI-09	-	-
135	VI-10	-	-
136	VI-12	869	667
137	VI-14	871	642
138	VI-17	-	-
139	VI-20	876	115
140	VI-21	-	-
141	VI-23	878	634
142	VI-34	-	-
143	VI-41	-	-
144	VI-42	-	-
145	VI-43	-	-
146	VI-44	-	-
147	VI-48	891	626
148	VI-49	-	-
149	VI-50	893	585
150	VI-52	-	-
151	VI-53	895	560
152	VI-55	897	509
153	VI-65	-	-
154	VI-70	108	550
155	VI-71	-	-

156	VI-72	-	-
157	VI-74	905	655
158	VI-76	907	582
159	VI-78	-	-
160	VI-79	-	-
161	VI-84	-	-
162	VI-87	911	595
163	VI-88	912	651
164	VI-90	-	-
165	VI-93	-	-
166	VI-95	915	230
167	VI-96	-	-
168	VII-02	-	-
169	VII-03	1196	412
170	VII-06	-	-
171	VII-10	-	-
172	VII-11	-	-
173	VII-15	1199	439
174	VII-19	562	580
175	VII-21	564	671
176	VII-25	-	-
177	VII-32	571	457
178	VII-36	575	209
179	VII-39	576	541
180	VII-42	579	502
181	VII-43	580	316
182	VII-46	583	631
183	VII-47	1200	526
184	VII-48	1201	613
185	VII-59	593	565
186	VII-60	-	-
187	VII-63	595	98
188	VII-66	598	362
189	VII-67	-	-
190	VII-72	600	595
191	VII-73	601	522
192	VII-75	-	-
193	VII-76	603	624
194	VII-77	1203	692
195	VII-80	605	338
196	VII-81	606	556
197	VII-83	-	-
198	VII-86	-	-
199	VII-88	-	-
200	VII-90	612	576
201	VII-91	613	341
202	VII-93	615	379
203	VIII-01	-	-
204	VIII-02	-	-
205	VIII-03	-	-
206	VIII-06	-	-
207	VIII-09	618	598
208	VIII-10	-	-
209	VIII-15	-	-
210	VIII-20	628	419
211	VIII-22	-	-

212	VIII-26	-	-
213	VIII-28	634	511
214	VIII-29	635	592
215	VIII-30	636	572
216	VIII-31	637	482
217	VIII-32	638	545
218	VIII-33	639	624
219	VIII-39	-	-
220	VIII-41	645	649
221	VIII-42	646	600
222	VIII-44	-	-
223	VIII-46	649	425
224	VIII-48	651	251
225	VIII-58	-	-
226	VIII-64	663	627
227	VIII-65	-	-
228	VIII-66	665	345
229	VIII-67	666	252
230	VIII-74	-	-
231	VIII-76	675	591
232	VIII-78	-	-
233	VIII-82	-	-
234	VIII-83	-	-
235	VIII-85	-	-
236	VIII-87	-	-
237	VIII-91	-	-
238	VIII-92	-	-
239	VIII-93	-	-
240	VIII-95	-	-
241	X-04	-	-
242	X-07	808	641
243	X-15	814	132
244	X-29	821	370
245	X-34	-	-
246	X-35	-	-
247	X-54	837	603
248	X-56	839	71
249	X-68	1207	642
250	X-72	849	622
251	X-94	860	501
252	XI-07	-	-
253	XI-13	1209	620
254	XI-50	-	-
255	XI-58	-	-
256	XI-81	1212	374
257	XII-07	1213	567
258	XII-17	-	-
259	XII-26	-	-
260	XII-27	-	-
261	XII-31	-	-
262	XII-32	-	-
263	XII-35	1214	620
264	XII-36	-	-
265	XII-52	-	-
266	XII-59	1216	484
267	XIII-19	1219	559

268	XIII-29	-	-
269	XIII-52	939	513
270	XIII-62	-	-
271	XIII-84	-	-
272	XIII-92	1221	741
273	XV-18	-	-
274	XV-22	1099	561
275	XV-24	-	-
276	XV-25	1224	485
277	XV-28	-	-
278	XV-34	-	-
279	XV-42	-	-
280	XV-68	-	-
281	XV-74	-	-
282	XV-93	-	-
283	XV-94	-	-
284	XV-96	-	-
285	XVI-36	1056	435
286	XVI-53	1230	741
287	XVI-59	-	-
288	XVI-66	1074	689
289	XVI-76	1083	198
290	XVI-77	1084	198
291	XVII-07	-	-
292	XVII-08	-	-
293	XVII-17	-	-
294	XVII-28	-	-
295	XVII-29	-	-
296	XVII-31	1139	503
297	XVII-36	-	-
298	XVII-39	-	-
299	XVII-40	1231	203
300	XVII-48	1148	587
301	XVII-55	-	-
302	XVII-58	-	-
303	XVII-67	-	-
304	XVII-72	-	-
305	XVII-76	1160	650
306	XVII-82	-	-
307	XVII-87	1165	502
308	XVII-95	1172	648

Table 1 b**List of sequences of probes informative for disease diagnosis**

Please see the note at the bottom

Clone ID	Sequence ID
I-09	298
I-10	299
I-13	1331
I-14	1178
I-15	300
I-16	301
I-17	302
I-19	304
I-20	305
I-22	306
I-23	307
I-24	308
I-25	309
I-28	310
I-30	1180
I-31	311
I-32	312
I-34	313
I-37	1440
I-38	314
I-39	315
I-40	316
I-42	1332
I-44	317
I-45	318
I-46	319
I-47	320
I-48	321
I-49	322
I-53	323
I-54	1181
I-56	324
I-57	325
I-58	326
I-60	327
I-64	328
I-67	330
I-69	331
I-71	332
I-72	333
I-73	334
I-77	335
I-79	336
I-80	337

I-81	338
I-82	339
I-86	1336
I-88	1182
I-95	1337
II-02	360
II-03	361
II-05	363
II-06	364
II-07	365
II-08	366
II-09	367
II-10	368
II-11	369
II-12	370
II-13	371
II-14	372
II-15	373
II-16	374
II-17	375
II-18	376
II-20	377
II-21	378
II-22	379
II-23	380
II-24	381
II-25	382
II-26	383
II-27	384
II-28	385
II-29	386
II-30	387
II-31	388
II-32	389
II-33	390
II-34	391
II-35	392
II-37	393
II-38	394
II-39	395
II-40	396
II-41	397
II-42	398
II-43	399
II-44	400
II-46	401
II-47	402
II-48	403
II-49	404
II-50	405
II-52	406

II-53	407
II-54	408
II-55	409
II-56	410
II-57	411
II-58	412
II-59	413
II-60	414
II-61	415
II-62	416
II-63	417
II-64	418
II-65	419
II-66	420
II-67	421
II-68	422
II-69	423
II-70	424
II-71	425
II-72	426
II-73	427
II-74	428
II-75	429
II-76	430
II-77	431
II-78	432
II-79	433
II-80	434
II-81	435
II-82	436
II-83	437
II-84	438
II-85	439
II-86	440
II-87	441
II-88	442
II-89	443
II-90	444
II-91	445
II-92	446
II-93	447
II-94	448
II-95	449
II-96	450
III-01	452
III-02	453
III-03	454
III-04	455
III-05	457
III-06	458
III-07	459

III-08	460
III-09	461
III-11	462
III-12	463
III-13	464
III-14	465
III-15	466
III-16	467
III-17	468
III-18	469
III-19	470
III-20	1183
III-21	471
III-22	472
III-23	473
III-24	474
III-25	475
III-26	476
III-27	477
III-28	478
III-29	479
III-31	481
III-32	482
III-33	483
III-34	484
III-35	485
III-37	486
III-39	487
III-40	488
III-42	489
III-43	490
III-44	491
III-45	492
III-46	493
III-47	494
III-48	495
III-49	496
III-50	497
III-51	498
III-52	499
III-53	500
III-54	501
III-55	502
III-56	503
III-57	504
III-58	505
III-59	506
III-61	507
III-62	508
III-63	509
III-64	510

III-65	511
III-66	512
III-67	513
III-69	514
III-70	515
III-71	516
III-73	517
III-74	518
III-75	519
III-77	520
III-78	521
III-79	522
III-80	523
III-81	524
III-82	1348
III-83	525
III-85	526
III-86	527
III-87	528
III-88	529
III-89	530
III-91	531
III-92	1351
III-93	532
III-94	533
III-95	534
III-96	535
IV-02	681
IV-04	682
IV-13	683
IV-14	684
IV-15	1185
IV-17	685
IV-23	1353
IV-26	1186
IV-28	686
IV-31	687
IV-32	688
IV-35	1355
IV-37	g6
IV-38	689
IV-40	690
IV-42	691
IV-43	1239
IV-44	692
IV-47	693
IV-53	61
IV-55	694
IV-56	695
IV-61	696
IV-64	697

IV-65	698
IV-69	192
IV-72	699
IV-73	700
IV-80	701
IV-82	196
IV-85	702
IV-93	703
IV-95	704
IV-96	705
IX-10	736
IX-12	738
IX-13	739
IX-24	747
IX-38	757
IX-39	758
IX-48	764
IX-50	766
IX-56	768
IX-62	773
IX-65	776
IX-72	782
IX-77	785
IX-91	796
IX-96	801
V-01	1361
V-03	706
V-04	707
V-07	708
V-08	709
V-09	710
V-11	1188
V1-16	873
V1-19	875
V-12	711
V-17	1364
V-18	712
V-20	713
V-24	714
V-25	1365
V-28	1189
V-35	1366
V-37	716
V-38	1190
V-39	1109
V-40	717
V-41	718
V-47	1368
V-48	719
V-49	1369
V-55	77

V-57	720
V-58	1370
V-61	721
V-64	722
V-65	723
V-68	1448
V-71	1495
V-74	724
V-75	1372
V-80	726
V-81	727
V-87	728
V-90	1374
VI-02	340
VI-03	341
VI-04	342
VI-06	343
VI-07	344
VI-08	345
VI-09	346
VI-11	347
VI-12	869
VI-13	870
VI-14	871
VI-16	873
VI-18	348
VI-19	349
VI-20	350
VI-21	351
VI-22	352
VI-23	878
VI-24	879
VI-25	353
VI-26	354
VI-27	355
VI-31	356
VI-32	885
VI-33	357
VI-35	358
VI-39	887
VI-43	1382
VI-44	1193
VI-45	889
VI-48	359
VI-49	892
VI-50	893
VI-53	895
VI-55	897
VI-58	899
VI-66	903
VI-67	904

VI-70	108
VI-71	1387
VI-74	905
VI-75	906
VI-76	907
VI-77	110
VI-79	1389
VI-80	908
VI-85	910
VI-87	911
VI-88	912
VI-90	1390
VI-93	1391
VI-95	915
VI-96	1392
VII-02	547
VII-03	548
VII-04	549
VII-05	550
VII-06	551
VII-07	552
VII-08	553
VII-09	554
VII-10	555
VII-11	556
VII-12	557
VII-14	558
VII-15	559
VII-17	560
VII-18	561
VII-19	562
VII-20	563
VII-21	564
VII-22	565
VII-23	566
VII-24	567
VII-25	1397
VII-26	250
VII-27	568
VII-28	569
VII-29	570
VII-32	571
VII-33	572
VII-34	573
VII-35	574
VII-36	575
VII-39	576
VII-40	577
VII-41	578
VII-42	579
VII-43	580

VII-44	581
VII-45	582
VII-46	583
VII-47	1200
VII-48	584
VII-49	585
VII-50	586
VII-52	587
VII-53	588
VII-54	589
VII-55	590
VII-57	591
VII-58	592
VII-59	593
VII-62	594
VII-63	595
VII-64	596
VII-65	597
VII-66	598
VII-67	1399
VII-71	599
VII-72	600
VII-73	601
VII-74	602
VII-76	603
VII-77	604
VII-80	605
VII-81	606
VII-82	607
VII-83	608
VII-84	609
VII-86	1453
VII-87	610
VII-89	611
VII-90	612
VII-91	613
VII-92	614
VII-93	615
VII-94	616
VII-96	617
VIII-09	618
VIII-10	619
VIII-11	620
VIII-12	621
VIII-13	622
VIII-15	623
VIII-16	624
VIII-17	625
VIII-18	626
VIII-19	627
VIII-20	628

VIII-21	629
VIII-22	1455
VIII-23	630
VIII-24	631
VIII-25	632
VIII-26	1456
VIII-27	633
VIII-28	634
VIII-29	635
VIII-30	636
VIII-31	637
VIII-32	638
VIII-33	639
VIII-34	640
VIII-36	641
VIII-37	642
VIII-38	643
VIII-40	644
VIII-41	645
VIII-42	646
VIII-43	647
VIII-45	648
VIII-46	649
VIII-47	650
VIII-48	651
VIII-50	652
VIII-51	653
VIII-53	654
VIII-54	655
VIII-55	656
VIII-56	657
VIII-57	658
VIII-58	659
VIII-59	660
VIII-60	661
VIII-61	662
VIII-64	663
VIII-65	664
VIII-66	665
VIII-67	666
VIII-68	667
VIII-69	668
VIII-70	669
VIII-71	670
VIII-72	671
VIII-73	672
VIII-74	673
VIII-75	674
VIII-76	675
VIII-77	676
VIII-78	677

VIII-79	678
VIII-80	679
X-07	808
X-15	814
X-20	817
X-29	821
X-34	825
X-46	833
X-54	837
X-56	839
X-68	1207
X-72	849
X-73	1208
X-94	860
XI-13	1209
XI-37	1460
XI-43	1210
XI-67	1211
XI-81	1212
XII-07	1213
XII-35	1214
XII-36	1215
XII-59	1216
XII-65	1028
XII-92	1217
XIII-03	917
XIII-04	1218
XIII-19	1219
XIII-24	926
XIII-51	938
XIII-52	939
XIII-67	947
XIII-69	949
XIII-88	1220
XIII-92	1221
XV-22	1099
XV-24	1101
XV-25	1224
XV-42	1108
XV-62	1226
XV-64	1118
XV-84	1125
XVI-19	1228
XVI-36	1056
XVI-53	1230
XVI-60	1071
XVI-66	1074
XVI-74	1081
XVI-76	1083
XVI-77	1084
XVII-31	1139

XVII-40	1231
XVII-48	1148
XVII-76	1160
XVII-87	1165
XVII-95	1172

Note

Sequences not available for sequence IDs in Table 1, and corresponding sequence IDs in Table 2 and 4.

298,301,305,307,312,317,318,319,320,332,333,334,336,340,341,342,343,344,345,346,347,348,349,350,351,352,353,354,355,356,357,358,359,367,372,375,376,377,379,385,392,393,404,437,439,440,443,444,445,449,455,457,465,466,467,468,470,486,498,501,511,514,516,517,520,522,528,531,535,547,548,549,550,551,552,553,554,555,556,557,558,559,573,584,604,608,616,620,623,640,659,662,664,667,668,673,677,678,679,681,695,702,712,716,825,886,894,902,909,916,1101,1108,1109,1177,1187,1193,1204,1220,1239,1255,1256,1342,1347,1354,1357,1362,1363,1364,1373,1375,1379,1403,1404,1405,1406,1413

Table 2a

List of informative probes for diagnosis of breast cancer

Clone ID	Sequence ID
I-24	308
I-28	310
I-30	1180
I-52	-
I-54	1181
II-41	397
II-70	424
II-87	441
III-06	458
III-20	1183
III-40	488
III-57	504
III-60	-
III-61	507
III-89	530
IV-14	684
IV-15	1185
IV-26	1186
IV-32	688
IV-41	-
IV-53	61
IV-62	-
IV-69	192
IV-80	701
IV-82	196
IX-10	736
IX-12	-
IX-38	757
IX-39	758
IX-42	-
IX-48	764
IX-77	785
V-11	1188
V-32	-
V-39	-
V-55	77
V-80	726
V-94	-
VI-07	93
VI-34	-
VI-41	-
VI-48	891
VI-49	-
VI-52	-
VI-55	897
VI-65	-
VI-70	108

Clone ID	Sequence ID
VI-72	-
VI-78	-
VI-84	-
VII-03	1196
VII-15	1199
VII-32	571
VII-39	576
VII-47	1200
VII-48	1201
VII-60	-
VII-73	601
VII-77	1203
VII-90	612
VIII-20	628
VIII-29	635
VIII-30	636
VIII-31	637
VIII-39	-
VIII-44	-
VIII-46	649
VIII-48	651
VIII-66	665
VIII-74	-
VIII-76	675
X-04	-
X-07	808
X-15	814
X-29	821
X-34	-
X-35	-
X-54	837
X-56	839
X-68	1207
X-72	849
X-94	860
XI-07	-
XI-13	1209
XI-50	-
XI-58	-
XI-81	1212
XII-07	1213
XII-17	-
XII-26	-
XII-27	-
XII-31	-
XII-32	-
XII-35	1214

Clone ID	Sequence ID
XII-36	-
XII-52	-
XII-59	1216
XIII-19	1219
XIII-29	-
XIII-52	939
XIII-62	-
XIII-84	-
XIII-92	1221
XV-18	-
XV-22	1099
XV-24	-
XV-25	1224
XV-28	-
XV-34	-
XV-42	-
XV-68	-
XV-74	-
XV-93	-
XV-94	-
XV-96	-
XVI-36	1056
XVI-53	1230
XVI-59	-
XVI-66	1074
XVI-76	1083
XVI-77	1084
XVII-07	-
XVII-08	-
XVII-17	-
XVII-28	-
XVII-29	-
XVII-31	1139
XVII-36	-
XVII-39	-
XVII-40	1231
XVII-48	1148
XVII-55	-
XVII-58	-
XVII-67	-
XVII-72	-
XVII-76	1160
XVII-82	-
XVII-87	1165
XVII-95	1172

Table 2b**List of sequences of probes informative for breast cancer**

Please see the note at the bottom of Table 1. Some sequences are missing.

Clone ID	Sequence ID
I-13	1331
I-14	1178
I-24	308
I-25	309
I-28	310
I-30	1180
I-37	1440
I-42	1332
I-48	321
I-54	1181
I-60	327
I-72	1335
I-81	338
I-82	339
I-86	1336
I-88	1182
I-95	1337
II-02	360
II-03	361
II-06	364
II-07	365
II-10	368
II-21	378
II-23	380
II-24	381
II-25	382
II-27	384
II-33	390
II-34	391
II-41	397
II-42	398
II-46	401
II-47	1338
II-48	403
II-52	406
II-57	411
II-58	412
II-59	413
II-60	414
II-61	415
II-62	416
II-64	418

II-67	421
II-69	423
II-70	424
II-74	428
II-80	434
II-82	436
II-84	438
II-87	441
II-88	442
II-96	450
III-01	452
III-02	453
III-06	458
III-08	460
III-12	463
III-13	464
III-17	1344
III-18	469
III-20	1183
III-21	471
III-23	473
III-24	474
III-25	475
III-26	476
III-27	477
III-28	478
III-29	479
III-32	482
III-33	483
III-35	485
III-39	487
III-40	488
III-42	489
III-45	492
III-46	493
III-47	494
III-48	495
III-56	503
III-57	504
III-58	505
III-59	506
III-61	507
III-62	508
III-63	509
III-64	510
III-66	512
III-67	513
III-70	515
III-74	518
III-75	519
III-78	521

III-80	523
III-81	524
III-82	1348
III-85	526
III-86	527
III-88	529
III-89	530
III-92	1351
III-93	532
III-95	534
III-96	1352
IV-04	682
IV-13	683
IV-14	684
IV-15	1185
IV-17	685
IV-23	1353
IV-26	1186
IV-31	687
IV-32	688
IV-35	1355
IV-37	96
IV-38	689
IV-42	691
IV-43	1239
IV-47	693
IV-53	61
IV-61	696
IV-64	697
IV-69	192
IV-72	699
IV-80	701
IV-82	196
IV-85	702
IV-93	1360
IV-96	705
IX-10	736
IX-12	738
IX-13	739
IX-24	747
IX-38	757
IX-39	758
IX-48	764
IX-50	766
IX-56	768
IX-62	773
IX-65	776
IX-72	782
IX-77	785
IX-91	796
IX-96	801

V-01	1361
V-03	706
V-04	707
V-07	708
V-08	709
V-11	1188
V-12	711
V-17	1364
V-24	714
V-25	1365
V-28	1189
V-35	1366
V-38	1190
V-39	1109
V-41	718
V-47	1368
V-49	1369
V-55	77
V-57	720
V-58	1370
V-61	721
V-64	722
V-65	1371
V-68	1448
V-71	1495
V-74	724
V-75	1372
V-80	726
V-90	1374
VI-03	864
VI-04	865
VI-07	93
VI-08	867
VI-09	1378
VI-12	869
VI-13	870
VI-14	871
VI-16	873
VI-19	875
VI-20	876
VI-21	1380
VI-23	878
VI-24	879
VI-25	1192
VI-26	881
VI-32	885
VI-39	887
VI-43	1382
VI-44	1193
VI-45	889
VI-48	891

VI-49	892
VI-50	893
VI-53	895
VI-55	897
VI-58	899
VI-66	903
VI-67	904
VI-70	108
VI-71	1387
VI-74	905
VI-75	906
VI-76	907
VI-77	110
VI-79	1389
VI-80	908
VI-85	910
VI-87	911
VI-88	912
VI-90	1390
VI-93	1391
VI-95	915
VI-96	1392
VII-02	1195
VII-03	1196
VII-06	1394
VII-08	1197
VII-09	1198
VII-10	1395
VII-11	1396
VII-15	1199
VII-17	560
VII-19	562
VII-21	564
VII-22	565
VII-23	566
VII-24	567
VII-25	1397
VII-26	250
VII-27	568
VII-29	570
VII-32	571
VII-33	572
VII-36	575
VII-39	576
VII-41	578
VII-42	579
VII-43	580
VII-46	583
VII-47	1200
VII-48	1201
VII-49	585

VII-54	589
VII-57	591
VII-58	592
VII-59	593
VII-62	594
VII-63	1202
VII-64	596
VII-66	598
VII-67	1399
VII-72	600
VII-73	601
VII-77	1203
VII-80	605
VII-82	607
VII-86	1453
VII-87	610
VII-90	612
VII-91	613
VII-92	614
VII-93	615
VII-96	617
VIII-09	618
VIII-10	619
VIII-13	622
VIII-16	624
VIII-20	628
VIII-21	629
VIII-22	1455
VIII-23	630
VIII-24	631
VIII-25	632
VIII-26	1456
VIII-27	633
VIII-28	634
VIII-29	635
VIII-30	636
VIII-31	637
VIII-32	638
VIII-33	639
VIII-34	1204
VIII-38	643
VIII-40	644
VIII-41	645
VIII-46	649
VIII-48	651
VIII-55	656
VIII-57	658
VIII-59	660
VIII-60	661
VIII-61	1205
VIII-64	663

VIII-66	665
VIII-73	672
VIII-74	673
VIII-76	675
VIII-80	679
X-07	808
X-15	814
X-20	817
X-29	821
X-34	825
X-46	833
X-54	837
X-56	839
X-68	1207
X-72	849
X-73	1208
X-94	860
XI-13	1209
XI-37	1460
XI-43	1210
XI-67	1211
XI-81	1212
XII-07	1213
XII-35	1214
XII-36	1215
XII-59	1216
XII-65	1028
XII-92	1217
XIII-03	917
XIII-04	1218
XIII-19	1219
XIII-24	926
XIII-51	938
XIII-52	939
XIII-67	947
XIII-69	949
XIII-88	1220
XIII-92	1221
XV-22	1099
XV-24	1101
XV-25	1224
XV-42	1108
XV-62	1226
XV-64	1118
XV-84	1125
XVI-19	1228
XVI-36	1056
XVI-53	1230
XVI-60	1071
XVI-66	1074
XVI-74	1081

XVI-76	1083
XVI-77	1084
XVII-31	1139
XVII-40	1231
XVII-48	1148
XVII-76	1160
XVII-87	1165
XVII-95	1172

Table 3

List of informative probes (Clone ID) selected for breast cancer diagnosis based on their occurrence criterion during variable selection.

Occurrence*	Clone ID
100%	XI-8,XVI-66,VIII-66,XVI-59,VII-03,XIII-19,XII-35,X-35,XI-50,XII-26,IV-53,XIII-29,XIII-62,I-30,III-06,XV-22,XV-94,VII-15,VII-39,IX-39,XVII-39,III-40,VII-32
90%	I-52,VI-65,VI-34,IV-62,XV-34,XVII-58,V-11,VI-78,XII-36,XIII-92,VIII-29,XVI-53,XVI-77,XI-13,XIII-84,IV-14,XII-31,V-80,VII-48,XVII-29,XVII-72
80%	III-60,VIII-74,IX-12,X-04,XIII-52,VIII-30,IX-38
70%	VI-49,X-29,VIII-48
60%	IV-82,IX-10,VI-52,X-68,VII-77
50%	IV-15
40%	XV-28,II-70,V-55
30%	XVII-17,XVII-67
20%	XI-58,XVI-36,VIII-39,VIII-44,III-61,IV-69,XV-68,X-72
10%	IX-42,IX-77,X-94,XV-96,XVII-55
5%	XII-59,XVI-76,I-54,XV-18,V-94,X-54,VI-07,VII-47,XVII-31,XVII-87,XVII-48
In at least one model	II-41,VI-41,III-57,III-89,VII-73,XV-25,IV-26,X-34,IV-41,VII-90,XV-42,XVII-82,XII-27,VIII-20,I-28,VII-60,VIII-76,III-20,VI-84,XI-07,XVII-28,XII-17,XVII-36,XII-52,XVII-76,VIII-46,VI-70,XV-74,XV-93,VIII-31,II-87,V-39,VI-55,X-07,X-15,XII-07,XVII-07,XVII-08,XVII-95,I-24,IV-32,V-32,VI-48,VI-72,IV-80,IX-48,X-56,XV-24,XII-32,XVII-40

*100% = Genes appearing in all the 75 cross validated models; 90% = Additional genes appearing in at least 68 out of 75 cross validated models; 5% = Additional genes appearing in at least 4 out of 75 cross validated models and so on.

Table 4a

List of informative probes for diagnosis of Alzheimer disease

Clone ID	Sequence ID
I-01	-
I-02	-
I-13	-
I-21	-
I-34	313
I-37	-
I-42	-
I-58	326
I-71	-
I-72	-
I-86	-
I-95	-
II-03	361
II-05	363
II-06	364
II-10	368
II-24	381
II-25	382
II-26	383
II-33	390
II-34	391
II-42	398
II-47	-
II-57	411
II-61	415
II-69	423
II-75	429
II-83	-
II-84	438
II-88	442
II-90	-
II-94	448
III-02	453
III-05	-
III-06	458
III-08	460
III-10	-
III-13	464
III-15	-
III-17	-
III-23	473
III-26	476
III-35	485
III-39	487
III-43	490
III-44	491
III-53	500
III-56	503

Clone ID	Sequence ID
III-60	-
III-63	509
III-68	-
III-74	518
III-80	523
III-82	-
III-85	526
III-92	-
III-96	-
IV-23	-
IV-26	-
IV-29	-
IV-31	687
IV-34	-
IV-35	-
IV-45	-
IV-80	701
IV-82	-
IV-93	-
V-01	-
V-02	-
V-03	706
V-04	707
V-06	-
V-07	708
V-12	711
V-15	-
V-17	-
V-21	-
V-25	-
V-35	-
V-42	-
V-43	-
V-47	-
V-49	-
V-52	-
V-54	-
V-58	-
V-59	-
V-65	-
V-68	-
V-71	-
V-75	-
V-79	-
V-80	726
V-90	-
V-91	-
V-92	-

Clone ID	Sequence ID
VI-02	-
VI-04	865
VI-09	-
VI-10	-
VI-12	869
VI-14	871
VI-17	-
VI-20	876
VI-21	-
VI-23	878
VI-41	-
VI-42	-
VI-43	-
VI-44	-
VI-48	891
VI-49	-
VI-50	893
VI-53	895
VI-71	-
VI-74	905
VI-76	907
VI-78	-
VI-79	-
VI-87	911
VI-88	912
VI-90	-
VI-93	-
VI-95	915
VI-96	-
VII-02	-
VII-03	-
VII-06	-
VII-10	-
VII-11	-
VII-19	562
VII-21	564
VII-25	-
VII-36	575
VII-42	579
VII-43	580
VII-46	583
VII-59	593
VII-63	595
VII-66	598
VII-67	-
VII-72	600
VII-73	601
VII-75	-
VI-02	-
VI-04	866
VI-09	-
VI-10	-
VI-12	873
VI-14	875
VI-17	-

Clone ID	Sequence ID
VII-91	613
VII-93	615
VIII-01	-
VIII-02	-
VIII-03	-
VIII-06	-
VIII-09	618
VIII-10	-
VIII-15	-
VIII-22	-
VIII-26	-
VIII-28	634
VIII-30	636
VIII-32	638
VIII-33	639
VIII-41	645
VIII-42	646
VIII-48	651
VIII-58	-
VIII-64	663
VIII-65	-
VIII-67	666
VIII-78	-
VIII-82	-
VIII-83	-
VIII-85	-
VIII-87	-
VIII-91	-
VIII-92	-
VIII-93	-
VIII-95	-

Table 4b**List of sequences of probes informative for Alzheimer disease**

Please see note to Table 1

Clone ID	Sequence ID
I-09	298
I-10	299
I-15	300
I-16	301
I-17	302
I-19	304
I-20	305
I-22	306
I-23	307
I-24	308
I-25	309
I-28	310
I-31	311
I-32	312
I-34	313
I-38	314
I-39	315
I-40	316
I-44	317
I-45	318
I-46	319
I-47	320
I-48	321
I-49	322
I-53	323
I-56	324
I-57	325
I-58	326
I-60	327
I-64	328
I-67	330
I-69	331
I-71	332
I-72	333
I-73	334
I-77	335
I-79	336
I-80	337
I-81	338
I-82	339
VI-02	340

VI-03	341
VI-04	342
VI-06	343
VI-07	344
VI-08	345
VI-09	346
VI-11	347
VI-18	348
VI-19	349
VI-20	350
VI-21	351
VI-22	352
VI-25	353
VI-26	354
VI-27	355
VI-31	356
VI-33	357
VI-35	358
VI-48	359
II-02	360
II-03	361
II-05	363
II-06	364
II-07	365
II-08	366
II-09	367
II-10	368
II-11	369
II-12	370
II-13	371
II-14	372
II-15	373
II-16	374
II-17	375
II-18	376
II-20	377
II-21	378
II-22	379
II-23	380
II-24	381
II-25	382
II-26	383
II-27	384
II-28	385
II-29	386
II-30	387
II-31	388
II-32	389
II-33	390
II-34	391
II-35	392

II-37	393
II-38	394
II-39	395
II-40	396
II-41	397
II-42	398
II-43	399
II-44	400
II-46	401
II-47	402
II-48	403
II-49	404
II-50	405
II-52	406
II-53	407
II-54	408
II-55	409
II-56	410
II-57	411
II-58	412
II-59	413
II-60	414
II-61	415
II-62	416
II-63	417
II-64	418
II-65	419
II-66	420
II-67	421
II-68	422
II-69	423
II-70	424
II-71	425
II-72	426
II-73	427
II-74	428
II-75	429
II-76	430
II-77	431
II-78	432
II-79	433
II-80	434
II-81	435
II-82	436
II-83	437
II-84	438
II-85	439
II-86	440
II-87	441
II-88	442
II-89	443

II-90	444
II-91	445
II-92	446
II-93	447
II-94	448
II-95	449
II-96	450
III-01	452
III-02	453
III-03	454
III-04	455
III-05	457
III-06	458
III-07	459
III-08	460
III-09	461
III-11	462
III-12	463
III-13	464
III-14	465
III-15	466
III-16	467
III-17	468
III-18	469
III-19	470
III-21	471
III-22	472
III-23	473
III-24	474
III-25	475
III-26	476
III-27	477
III-28	478
III-29	479
III-31	481
III-32	482
III-33	483
III-34	484
III-35	485
III-37	486
III-39	487
III-40	488
III-42	489
III-43	490
III-44	491
III-45	492
III-46	493
III-47	494
III-48	495
III-49	496
III-50	497

III-51	498
III-52	499
III-53	500
III-54	501
III-55	502
III-56	503
III-57	504
III-58	505
III-59	506
III-61	507
III-62	508
III-63	509
III-64	510
III-65	511
III-66	512
III-67	513
III-69	514
III-70	515
III-71	516
III-73	517
III-74	518
III-75	519
III-77	520
III-78	521
III-79	522
III-80	523
III-81	524
III-83	525
III-85	526
III-86	527
III-87	528
III-88	529
III-89	530
III-91	531
III-93	532
III-94	533
III-95	534
III-96	535
VII-02	547
VII-03	548
VII-04	549
VII-05	550
VII-06	551
VII-07	552
VII-08	553
VII-09	554
VII-10	555
VII-11	556
VII-12	557
VII-14	558
VII-15	559

VII-17	560
VII-18	561
VII-19	562
VII-20	563
VII-21	564
VII-22	565
VII-23	566
VII-24	567
VII-27	568
VII-28	569
VII-29	570
VII-32	571
VII-33	572
VII-34	573
VII-35	574
VII-36	575
VII-39	576
VII-40	577
VII-41	578
VII-42	579
VII-43	580
VII-44	581
VII-45	582
VII-46	583
VII-48	584
VII-49	585
VII-50	586
VII-52	587
VII-53	588
VII-54	589
VII-55	590
VII-57	591
VII-58	592
VII-59	593
VII-62	594
VII-63	595
VII-64	596
VII-65	597
VII-66	598
VII-71	599
VII-72	600
VII-73	601
VII-74	602
VII-76	603
VII-77	604
VII-80	605
VII-81	606
VII-82	607
VII-83	608
VII-84	609
VII-87	610

VII-89	611
VII-90	612
VII-91	613
VII-92	614
VII-93	615
VII-94	616
VII-96	617
VIII-09	618
VIII-10	619
VIII-11	620
VIII-12	621
VIII-13	622
VIII-15	623
VIII-16	624
VIII-17	625
VIII-18	626
VIII-19	627
VIII-20	628
VIII-21	629
VIII-23	630
VIII-24	631
VIII-25	632
VIII-28	634
VIII-29	635
VIII-30	636
VIII-31	637
VIII-32	638
VIII-33	639
VIII-34	640
VIII-36	641
VIII-37	642
VIII-38	643
VIII-40	644
VIII-41	645
VIII-42	646
VIII-43	647
VIII-45	648
VIII-46	649
VIII-47	650
VIII-48	651
VIII-50	652
VIII-51	653
VIII-53	654
VIII-54	655
VIII-55	656
VIII-56	657
VIII-57	658
VIII-58	659
VIII-59	660
VIII-60	661
VIII-61	662

VIII-64	663
VIII-65	664
VIII-66	665
VIII-67	666
VIII-68	667
VIII-69	668
VIII-70	669
VIII-71	670
VIII-72	671
VIII-73	672
VIII-74	673
VIII-75	674
VIII-76	675
VIII-77	676
VIII-78	677
VIII-79	678
VIII-80	679
IV-02	681
IV-04	682
IV-13	683
IV-14	684
IV-17	685
IV-28	686
IV-31	687
IV-32	688
IV-38	689
IV-40	690
IV-42	691
IV-44	692
IV-47	693
IV-55	694
IV-56	695
IV-61	696
IV-64	697
IV-65	698
IV-72	699
IV-73	700
IV-80	701
IV-85	702
IV-93	703
IV-95	704
IV-96	705
V-03	706
V-04	707
V-07	708
V-08	709
V-09	710
V-12	711
V-18	712
V-20	713
V-24	714

V-37	716
V-40	717
V-41	718
V-48	719
V-57	720
V-61	721
V-64	722
V-65	723
V-74	724
V-80	726
V-81	727
V-87	728
VI-13	870
VI-14	871
VI-16	873
VI-23	878
VI-24	879
VI-28	883
VI-32	885
VI-38	886
VI-39	887
VI-45	889
VI-46	890
VI-49	892
VI-50	893
VI-52	894
VI-53	895
VI-54	896
VI-55	897
VI-57	898
VI-58	899
VI-63	900
VI-65	902
VI-66	903
VI-67	904
VI-74	905
VI-75	906
VI-76	907
VI-80	908
VI-81	909
VI-85	910
VI-87	911
VI-88	912
VI-91	913
VI-94	914
VI-95	915
VI-96	916
I-13	1177
I-14	1178
I-30	1180

I-54	1181
I-88	1182
III-20	1183
IV-15	1185
IV-26	1186
IV-62	1187
V-11	1188
V-28	1189
V-38	1190
V-45	1191
VI-44	1193
VII-47	1200
I-42	1332
I-52	1333
I-86	1336
I-95	1337
III-10	1342
III-60	1347
III-82	1348
III-92	1351
IV-23	1353
IV-34	1354
IV-35	1355
IV-41	1356
IV-45	1357
IV-82	1359
V-01	1361
V-02	1362
V-06	1363
V-17	1364
V-25	1365
V-35	1366
V-42	1367
V-47	1368
V-49	1369
V-58	1370
V-75	1372
V-79	1373
V-90	1374
V-91	1375
V-94	1376
VI-10	1379
VI-41	1381
VI-43	1382
VI-71	1387
VI-72	1388
VI-79	1389
VI-90	1390
VI-93	1391
VII-25	1397
VII-60	1398

VII-67	1399
VIII-22	1403
VIII-26	1404
VIII-39	1405
VIII-44	1406
I-37	1440
V-32	1445
V-52	1447
V-68	1448
V-92	1449
VI-42	1450
VI-78	1452
VII-86	1453
VII-88	1454
IV-29	1490
V-15	1491
V-39	1492
V-54	1493
V-59	1494
V-71	1495

Table 5

Samples

Diagnosis	No. of women
Normal /Benign	42*
DCIS	3
Invasive cancer	26

* From one woman, whole blood was collected at weeks 1,2,3,4,5 following menstruation. Hence, the number of unique normal/benign samples tested in the experiment is 75

Information about women with breast cancer

Sample	AGE	Stage	Cancer type	Size hist. (mm)	Nodes
1	51	II	IDC	20	1/7
2	84	II	IDC	22	2/2
3	50	I	DCIS+ 1 IDC	>50 DCIS; 5 x 14	0/7
4	47	I	IDC	15	0
5	69	III	ILC g.2 + tubular adenocarcinoma	50 + 3	1 av 12 + 1 av 7
6	50	II	IDC	24	0
7	65	I	IDC	15	0
8	63	II	IDC	23	0
9	55	I	IDC + DCIS	4	0 av 1
10	52	0	DCIS + small colloid carcinoma foci	50 + 3	0
11	60	II	IDC	24	0
12	54	I	IDC	11	0
13		0	DCIS	20	0
14	49	0	DCIS	9	0
15	48	I	IDC	4	0
16	56	I	IDC	4	0
17	68	I	IDC	14	0
18	68	I	IDC	7	0
19	63	I	IDC	10	0
20	45	I	IDC	19	1
21	57	III	IDC	60	8/20
22	55	II	IDC/DCIS	35 + 55	0
23	71	I	IDC/extensive DCIS	8	0
24	56	I	IDC	9	?

25	66	II	IDC	26	0
26	66	I	IDC	15	?
27	61	I	IDC	9	?
28	?	?	?	?	?
29	65	I	IDC	11	0

Other diseases /conditions present in the women tested

Other diseases /conditions present in the women tested

Disease/condition
Diabetes
Asthma
Ulcerous colitis
Hemochromatose
Crohn's disease
Fibromyalgia
Psoriasis
Atopic eczema
Rheumatism
Allergies

Prior history of cancer in the women tested

Cancer type	No. of women
Breast	3
Colon	2
Stomach	1
Skin	1

Table 6

Number of samples tested by double cross validation and success of the diagnostic test for breast cancer based on selected informative genes

Number of samples tested by double cross validation

Number of unique samples tested	75
Number of unique non cancer samples tested	48
Number of cancer samples tested	29

Success of the diagnostic test for breast cancer based on selected informative genes

Occurrence in percentage*	Number of informative probes	Specificity	Sensitivity	Accuracy	False Positive rate	False negative rate	Total error rate
100.00	23	84.78	75.86	81.33	15.22	24.14	18.67
90.00	44	91.30	79.31	86.67	8.70	20.69	13.33
80.00	51	86.96	79.31	84.00	13.04	20.69	16.00
70.00	54	88.13	75.86	84.00	10.87	24.14	16.00
60.00	58	89.13	75.86	84.00	10.87	24.14	16.00
50.00	59	89.13	75.86	84.00	10.87	24.14	16.00
40.00	63	88.13	75.86	84.00	10.87	24.14	16.00
30.00	66	86.96	75.86	82.67	13.04	24.14	17.33
20.00	74	89.13	75.86	84.00	10.87	24.14	16.00
10.00	79	89.13	75.86	84.00	10.87	24.14	16.00
5.00	90	86.96	79.31	84.00	13.04	20.69	16.00
1.33	139	84.78	72.41	80.00	15.22	27.59	20.00

*100% = Genes appearing in all the 75 cross validated models; 90% = Genes appearing in at least 68 out of 75 cross validated models;

5% = Genes appearing in at least 4 out of 75 cross validated models; and so on.

Table 7

Double cross-validation and details of the success of the diagnostic test for Alzheimer disease based on the expression 182 informative genes

Validation Result

Total number of samples tested	14
Number of Alzheimer's disease samples tested	7
Number of Alzheimer's disease samples incorrectly predicted	1
Number of non-Alzheimer's disease samples tested	7
Number of non-Alzheimer's disease samples incorrectly predicted	0

Success of diagnostic test

Performance	Description	%
Accuracy	Percentage of the total number of predictions that were correct	92.9
Sensitivity	Percentage of positive cases that were correctly identified	85.7
Specificity	Percentage of negatives cases that were correctly predicted	100
False positive rate	Percentage of negatives cases that were incorrectly classified as positive	0.0
False negative rate	Percentage of positives cases that were incorrectly classified as negative	14.3
Total error rate	Percentage of the total cases incorrectly predicted	7.1

Table 3

Some relevant features of the blood donors. B, Female donors with breast cancer; N, Female donors with suspected mammogram but no breast cancer; IDC, invasive ductal carcinoma; DCIS, ductal carcinoma in situ; na, not available nd, not determined; ++, no degradation of mRNA and no ribosomal contamination in the sample, +, no degradation of mRNA but ribosomal contamination in the sample.

		AGE	Cancer type /breast abnormality	Size Hist. (mm)	mRNA Quality
1	B1	na	IDC	5	++
2	B2	49	DCIS	8	nd
3	B3	54	IDC	18	++
4	B4	59	IDC	12	+
5	B5	61	DCIS+micro invasive cancer	15+1.5	++
6	B6	55	IDC	12+17	nd
7	B6		IDC	12+17	nd
8	N1	45	Fibroadenoma	-	nd
9	N2	52	na	-	+
10	N3	55	Cyst	-	++
11	N4	54	na	-	++
12	N5	51	Benign ductal epithelium	-	nd
13	N6	57	Benign	-	nd
14	N7	50	na	-	++
15	N8	52	na	-	+

Table 9

List of sequence of probes informative for both alzheimer and breast cancer disease

Clone ID	Sequence ID
I-24	308
I-25	309
I-28	310
I-48	321
I-60	327
I-72	333
I-81	338
I-82	339
II-02	360
II-03	361
II-06	364
II-07	365
II-10	368
II-21	378
II-23	380
II-24	381
II-25	382
II-27	384
II-33	390
II-34	391
II-41	397
II-42	398
II-46	401
II-47	402
II-48	403
II-52	406
II-57	411
II-58	412
II-59	413
II-60	414
II-61	415
II-62	416
II-64	418
II-67	421
II-69	423
II-70	424
II-74	428
II-80	434
II-82	436
II-84	438

II-87	441
II-88	442
II-96	450
III-01	452
III-02	453
III-06	458
III-08	460
III-12	463
III-13	464
III-17	468
III-18	469
III-21	471
III-23	473
III-24	474
III-25	475
III-26	476
III-27	477
III-28	478
III-29	479
III-32	482
III-33	483
III-35	485
III-39	487
III-40	488
III-42	489
III-45	492
III-46	493
III-47	494
III-48	495
III-56	503
III-57	504
III-58	505
III-59	506
III-61	507
III-62	508
III-63	509
III-64	510
III-66	512
III-67	513
III-70	515
III-74	518
III-75	519
III-78	521
III-80	523
III-81	524
III-85	526
III-86	527
III-88	529
III-89	530
III-93	532
III-95	534

III-96	535
IV-04	682
IV-13	683
IV-14	684
IV-17	685
IV-31	687
IV-32	688
IV-38	689
IV-42	691
IV-47	693
IV-61	696
IV-64	697
IV-72	699
IV-80	701
IV-85	702
IV-93	703
IV-96	705
V-03	706
V-04	707
V-07	708
V-08	709
V-12	711
V-24	714
V-41	718
V-57	720
V-61	721
V-64	722
V-65	723
V-74	724
V-80	726
VI-03	341
VI-04	342
VI-07	344
VI-08	345
VI-09	346
VI-12	869
VI-14	871
VI-19	349
VI-20	350
VI-21	351
VI-23	878
VI-25	353
VI-26	354
VI-48	359
VI-50	893
VI-53	895
VI-74	905
VI-76	907
VI-87	911
VI-88	912
VI-95	915

VII-02	547
VII-03	548
VII-06	551
VII-08	553
VII-09	554
VII-10	555
VII-11	556
VII-15	559
VII-17	560
VII-19	562
VII-21	564
VII-22	565
VII-23	566
VII-24	567
VII-27	568
VII-29	570
VII-32	571
VII-33	572
VII-36	575
VII-39	576
VII-41	578
VII-42	579
VII-43	580
VII-46	583
VII-48	584
VII-49	585
VII-54	589
VII-57	591
VII-58	592
VII-59	593
VII-62	594
VII-63	595
VII-64	596
VII-66	598
VII-72	600
VII-73	601
VII-77	604
VII-80	605
VII-82	607
VII-87	610
VII-90	612
VII-91	613
VII-92	614
VII-93	615
VII-96	617
VIII-09	618
VIII-10	619
VIII-13	622
VIII-16	624
VIII-20	628
VIII-21	629

VIII-23	630
VIII-24	631
VIII-25	632
VIII-28	634
VIII-29	635
VIII-30	636
VIII-31	637
VIII-32	638
VIII-33	639
VIII-34	640
VIII-38	643
VIII-40	644
VIII-41	645
VIII-46	649
VIII-48	651
VIII-55	656
VIII-57	658
VIII-59	660
VIII-60	661
VIII-61	662
VIII-64	663
VIII-66	665
VIII-73	672
VIII-74	673
VIII-76	675
VIII-80	679

- 124 -

Nucleotide sequences

Sequence ID - 93

nt: 405

5 GGATCCTGTGGCCACAGAGCTGCCCCAGCAGACGCTCCGCCCCACCCGGTGATGG
AGCCCCGGGGGACAATCGTGCCTGGGGAGGAGCAGGGTACAGCCCATTCCCCCAG
CCCTGGCTGACCTGGCCTAGCAGTTTGGCCCTGCTGGCCTTAGCAGGGAGACAGGG
GAGCAAAGAACGCCAAGCCGGAGGCCCGAGGCCAGCCGGCCTCTCGAGAGCCAGAG
CAGCAGTTGAATGTAATGCTGGGGACAGGCATGCTGCCGCCAGTAGGGCGGGGACC
CGGACAGCCAGGTGACTACCAGTCCTGGGGACACACTCACCATAAACACATCCCCA
10 GGCAGGACAGATCGGGGAAGGGGTGTGTACCAGGCTATGATTTCTCTTGCAATTA
ATGTATTATTATT

Sequence ID - 108

nt: 550

15 GGCTTTGACAGAGTGCAAGACGATGACTTGCAAAATGTCGCATCTGGAACGCAACA
TAGANACCATCATCAACACCTTCCACCAATACTCTGTGAAGCTGGGGCACCAGAC
ACCCTGAACCAGGGGGAATTCAAAGAGCTGGTGCGAAAAGATCTGCAAAATTTTCT
CAAGAAGGAGAATAAGAATGAAAAGGTCATAGAACACATCATGGAGGACCTGGACA
CAAATGCAGACAAGCAGCTGAGCTTCGAGGAGTTCATCATGCTGATGGCGAGGCTA
ACCTGGGCCTCCACGAGAAGATGCACGAGGGTGACGAGGGCCCTGGCCACCACCA
20 TAAGCCAGGCCTCGGGGAGGGCACCCCTAAGACCACAGTGGCCAAGATCACAGTG
GCCACGGCCACGGCCACAGTCATGGTGGCCACGGCCACAGCCACTAATCAGGAGGC
CAGGCCACCCTGCCTNTACCCAACCAGGGCCCCGGGGCCTGTTATGTCAAAGTGT
TTGGCTGTGGGGCTAGGGGCTGGGGCCAAATAAAGTCTCTTTCTCC

25 Sequence ID 110

ACGAAGACAGACATCTGTGGAATGATTCACATCCTCTCAAGTTAGGAGGATGGAGG
CCTGCTTCATTAAGAAGCTGGGGGTAGGGTGGGGGTGGGGAGAACACTTAACAACA
TGGGGACCAGTCAGGGGAATCCCCTTATTTCTGTTTTGCATATGAGGAACCCTAGA
GCAGCCAGGTGAGGCTCTCTAGTTTAATAAAAAATCATGGAAAGACTCTTAATGCAG
30 ACTCTTCTTAAGTGTTAATAGGGATTTTTTCAGCTTATTTTGGTTGCAGTTTCCAA
TTTTTA AAAATGTTGAGGTAATCTTTCCACCTTCCCAAACCTAATTCTTGTAGAT
GCATTAGTGTTGAACCAATGCTTTCTCATGTCTCAATTCTTTGTATATGCATTCTT
TTCAGATGTATTAAACAAACAAAACCCCTTC

35 Sequence ID - 192

nt: 286

CCGGTAATAGAATAGAAAAGGGAGAGTGTCTTCATGCAATGTGGCATCCTGGATTG
GGTCTCGNNACAAAAACAGGACATTAGTGGGAAAATTGGAAATCTGAAAAAAGTCT

- 125 -

GAATTTTAGTTAATATAACCAATTTTCAGTCTCTTGGTTTTGACAGATGTACCATGGT
GATGTAAGATGTTGACCTTGGGGTAGGCTGGGTGAAGGGTATACAGGAACTCTTTG
TACTATCTCTGCAACTTCTCTGTAAATCTAGTATCATTCCAAAATAAAAGTTTATT
TAATTT

5

Sequence ID 250

GTGGAAGTGACATCGTCTTTAAACCCTGCGTGGCAATCCCTGACGCACCGCCGTGA
TGCCCAGGGAAGACAGGGCGACCTGGAAGTCCAACACTTCTTAAAGATCATCCAA
CTATTGGATGATTATCCGAAATGTTTCATTGTGGGAGCAGACAATGTGGGCTCCAA
10 GCAGATGCAGCAGATCCGCATGTCCCTTCGCGGGAAGGCTGTGGTGCTGATGGGCA
AGAACACCATGATGCGCAAGGCCATCCGAGGGCACCTGGAAAACAACCCAGCTCTG
GAGAACTGCTGCCTCATATCCGGGGGAATGTGGGCTTTGTGTTACCAAGGAGGA
CCTCACTGAGATCAGGGACATGTTGCTGGCCAATAAGGTGCCAGCTGCTGCCCCGTG
CTGGTGCCATTGCCCCATGTGAAGTCACTGTGCCAGCCAGAACACTGGTCTCGGG
15 CCCGAGAAGACCTCCTTTTTCCAGGCTTTAGGTATCACCATAAAATCTCCAGGGG
CACCATTGAAATCCTGAGTGATGTGCACTGATCAAGACTGG

Sequence ID 299

CAGCGCAGGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCGCAGATAAG
20 TTTTTTTCTCTTTGAAAGATAGAGATTGNTACAACTACTTAAAAAATATAGTCAAT
AGGTTACTAAGATATTGCTTAGCGTTAAGTTTTTAACGTAATTTTAATAGCTTAAG
ATTTTAAGAGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAAA
AGGTTTCTAAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAATG
TATTTAAAAGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACCAATAGA
25 AGGGCAATGCTTTTAGATTAAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTTT
AAAAGTTGTAGGTGATTAAAATAATTTGAAGGCGATCTTTTAAAAGAGATTAAAC
CGAAGGTGATTAAAAGACCTTGAAATCCATGACGCANGGAGAATTGCGCATTTAAA
GCCTAGTTACGCATTTACTAAACGCAGACGAAAATGGGAAGATTAATTGGGAGTG
TAGGATGAAACAATTTTGGAGAAGATAGAAG

30

Sequence ID 300

CTCAAAGGAGAAAAAAACCTTGTAAAAAAAGCAAAAATGACAACAGAAAAACAAT
CTTATTCCGAGCATTCCAGTAACCTTTTTTGTGTATGTACTTAGCTGTACTATAAGT
AGTTGGTTTGTATGAGATGGTTAAAAAGGCCAAAGATAAAAGGTTTCTTTTTTTTTT
35 CCTTTTTTGTCTATGAAGTTGCTGTTTATTTTTTTTTTGGCCTGTTTGATGTATGTGT
GAAACAATGTTGTCCAACAATAAACAGGAATTTTATTTTGCTGAGTTGTTCTAAAA
AAAAAAAAAAAAAAAAAA

- 126 -

Sequence ID 302

AGTAGAGACGGGGTTTCACTGTGTTAGCCAGGATGGTCTCGATCTCCTGACCTCGT
GATCCGGCCACCTCGGCCTCCCGAAAGTGCTGGGATTACAGGCGTGAGCCACGGCG
CCCAGCCCCAGCCTGTCACTTAAACTGATAAACGACAGATTAACAGTAGAAAAATT
5 TTATTTTGCATACATAATGAGGCTTCACAAAAGAGAAGTGAAAACCCAAGTAGGAG
TTTAGGGCTGGGGGCTTATATACCATTTAACAAGGGGTGATAAATTGTAAGAGAAT
AG

Sequence ID 304

10 TCCTTGGTTTTCGATTTGTGGCAACAATCCAGTCTTTTTGTTTTTTTCAGGGATACC
ATATGTAACAGGTGCCATTGTTACTGTAACCTTTTCACACATGCCTTCAGTTTGATG
TCAAAGTCATCATTTAGTGTAACAGCAAGTTATCTGTTAGGCTGCACATCATGAA
CTTTACTTTTAGAAAGTCTTATCTTTTATGCCACAGAAATAGCATTGCTATTAG
TCATGGATGGCAAAGAAATTAATTTTGAGTTGTTTGGATAAAAATGTTTCAGTTGA
15 CTGTAGTGTGTATTGAGAGACACTGCCAGTAAACAAACTCTCTTGGTAGGTGGAAA
TCCCCTAGAAGTTACAGAAAATTGGGAGGAGGTGAACTTAATTAAATAACTTGAAT
TGTTTAGACATATTCAGAGCTTCTTATGACCTTGAAGAAATCACCCAACCTCAAAA
GACCTCGGTTTCTTCATTTGTAAAATTAGGGAGTTTGACTAGATGTGTAAATCTAG
TTGTTAGTTAACTTCTAAGATGTAAAAACCTCTTGTTTAAACAAAACCTACAAGA
20 TCAAGTTGCTTATCTGAAATCTTTATGAATCAACACTAGTCACTAAGTCTAGCTCG
ACC

Sequence ID 306

CTTTTCCTCCCGCTGTCCCCCACGGAGGGGACTGCTCTCCCCCGCTGCATCCTTTC
25 TGTGAGGTACCTTACCCACCTCAGCACCTGAGAGGGTGAAATAGAATTCTAACCTC
GACATTCGGGAAGTGTTTTTGTGAGAAGTCTCGGTCCGTAAGGGAAGTCTTCCAAGTC
CGTGCAGCACTAACGTATTGGCACCTGCCTCCTCTTCGGCCACCCCCCAGATGAGG
CAGCTGTGACTGTGTCAAGGGAAGCCACGACTCTGACCATAGTCTTCTCTCAGCTT
CCACTGCCGTCTCCACAGGAAACCCAGAAGTTCTGTGAACAAGTCCATGCTGCCAT
30 CAAGGCATTTATTGCAGTGTACTATTTGCTTCCAAAGGATCAGGCCCTGAGAACAA
TGACCTTATTTCTTACAACAGTGTCTGGGTTGCGTGCCAGCAGATGCCTCAGATAC
CAAGAGATAACAAAGCTGCAGCTCTTTTGATGCTGACCAAGAATGTGGATTTTGTG
AAGGATGCNCATGAANAAATGGACNAGCTGTG

35 Sequence ID - 308

nt: 373

AAGTGGGTCTTGCCATCCCTGAACTGNAATCATCCCTAACATATTCATACCTGTTT
TCATTTTAAAAGTTGGGTCAGTTTTTTTATTAGTACATGTATTTCTATCCTACTGA

- 127 -

TTTATTTGCTATATCATCTAATTTAGTTTGAATATTCCATAATTTACTTAATTAGT
CCTGTATGGAGACCTAGCTCTTCTCAGTGTCTACTATTATAAACAATGCTACAGTG
AATATTGGTGNATAAATCCATACNCACCACGTACATATCTTAAGTTCTGGAAGAGA
TATTGCTAAACCAGAAGATAACCTGCATTTAAAATTTGACTGCTAGGGNCAGGGNC
5 ACATTTAATTAAATTAGAACAANGAATGCATAATGNC

Sequence ID 309

CCGGAATCGCGGCCGCGTCGACGAAAATATGTGCCCTGGCCAACTCCACAGGACTA
GTTCTAGGCAATCTGAAGGAAACCAGAAAATGTGAATTTCTCTTCCCTCAAAAAGC
10 TATACTGAAGTAGTATTTAATATTCAAGTACTTGTAATTTGCAGAACAGTACTTT
TTAATTTGACCCATGAATTTCTATTTAAATTTGTCACTTAATATTTAGCCAAGAAGC
AAACCATCTAAAAAGATTTCTGGTTTATTTCTCCAACCTCCTAATAAATAGGGTCAC
ATATTTTTTTAACTTTTTTTCTAATTTGAAAAGTAATACAGGCATATGGTATTTTTAAA
AATGAAACAACACAAAGGGATATGTTTTGAAAAGTGGTCTTGCCATCCCTGAACTG
15 TAATCATCCCTAACATATTCATACCTGTTTTCATTTTTAAAAGTTGGGTCAGTTTTT
TTATTAGTACATGTATTTCTATCCTACTGATTTATTTGCTATATCATCTAATTTAG
TTTGAATATTCCATAATTTACTTAATTAGTCCTGTATGGAGACCTAGCTCTTCTCA
GTGTCTACTATTATAAACAATGCTACAGTGAATATTGGTGNATAAATCCTACACAC
CACGTAACATATCTTAAGTTCCTGGAAGAGATATTGCTAAACCAGAAGATAACCTG
20 CATTTAAAATTTGACTGCTAGGGTCAGGGTCACATTTAAATTAAATTAGAACAAGG
AATGCATAATGTCTTCGATAGCAATCTATTCAAGGTGCACCGTGGTCACAAAGGAA
AGCAAAACTGTC

Sequence ID - 310 nt:564

CCTGGNCAGAGGCCTCTATCCTGTANTGATAATTGCCATCAAAATTGTCAAAAANG
ATTTAATTTCTATGGGNAATAGTCCTTTTCTTAGCTTCTGCCNNTCACTTGCTTAT
TTTTTGTGTGGGAATGGGGTTGGATAAACCAATGAACTTTATTATAAACAAATCCC
ACCTATATCTANCAAATTTATATTTTCGGTGAAATACAGATATTTGCCTTTCTGGGA
GTANTATAGAAGCTGTCAATATGTATCTACTGTACAGTACTAAATAGTATTCATTT
30 ATGAAATGAGTAGTGTTTGGGTGGCTGGGGTTAAGGAAAAATGAGACTTGGAATTG
TAGCTTTTATCCAAGTTTTTGAGTATAAATAGGGTTTTGTTTTGTTTTTTTAACTT
AAAACTGAAATGCCATATAGAAAAACAGCATTGTTTTTACAGTTTGTAGTAAGTA
ACTTTTTAAAGATTTTATCAAAAAGAATTTTGTCTATNGTGAGTAAAAGAAGTTCT
AATAATGGCCTAATCACTGCATTTTTTAAAAAACAAAGTTCAACACAAATGACATTT
35 GTTT

Sequence ID 311

- 128 -

5 CCTCTCCTCCATCTAAAGGCAACATTCCTTACCCATTAGTCTCAGAAATTGTCTTA
AGCAACAGCCCCAAATGCTGGCTGCCCCCGGCCAAGCATTGGGGCCGCCATCCTGC
CTGGCACTGGCTGATGGGCACCTCTGTTGGTTCCATCAGCCAGAGCTCTGCCAAAG
GCCCCGCAGTCCCTCTCCCAGGAGGACCCTAGAGGCAATTAAATGATGTCCTGTTT
CATTGG

Sequence ID - 313

nt: 554

10 CCCGGAATCGCGGCCCGCTCGACAACAAACCTGCATGTTCTGCACATGTATCCAG
GAACTTAAAAAAAAAAAAAGATAGTTTGTGTGTCTTAATTGAATAATAGTAGATTT
ATAGATTAAAGATCTATGGGTTTTTAATATGGATTANAAATCTGTGGGTTTTTGAT
ATGGATTANAAATCTGTGGGTTTTTAATATGGATTGGAAATCTGTGGGTTTTTAAT
ATGGATTAAAAAACATCTGTGGGTTTTTAATATGGATTAAACATCTGTGGGTTTTT
AATATGGATTAAACATCTGGGTTTTTAATATGGATTAAACATCTGTGGGTTTTTAA
TATGGGTTAAAAATCAAAGAAAATGAACTATTTGCTCCAGTGCAGGAAAATACAG
15 GCAATACTGGATAACAATTAGATGGTCAGGAGCGATAACCCGGTTGCCATTGTTTGA
AGAAGAGAATAAGGNGCTAGCATTCCTATCCGTAGATAATTTGACAGCTAGGAAAT
AGGGGGAGTCTTCTATGTAGTTAGTGAAGGCTAAATGAACTATTATATGC

Sequence ID 314

20 CTTTTCTCCCGCTGTCCCCACGGAGGGGACTGCTCTCCCCCGCTGCATCCTTTC
TGTGAGGTACCTTACCCACCTCAGCACCTGAGAGGGTGAAATAGAATTCTAACCTC
GACATTTCGGGAAGTGTTTTTGAGAAGTCTCGGTTCGGTAAGGGAAGTCTTCCAAGTC
CGTGCAGCACTAACGTATTGGCACCTGCCTCCTCTTCGGCCACCCCCCAGATGAGG
CAGCTGTGACTGTGTCAAGGGAAGCCACGACTCTGACCATAGTCTTCTCTCAGCTT
25 CCACTGCCGTCTCCACAGGAAACCCAGAAGTTCTGTGAACAAGTCCATGCTGCCAT
CAAGGCATTTATTGCAGTGTAATTTGCTTCCAAAGGATCAGGCCCTGAGAACAA
TGACCTTATTTCTACAACAGTGTCTGGGTTGCGTGCCAGCAGATGCCTCAGATAC
CAAGAGATAACAAAGCTGCAGCTCTTTTGATGCTGACCAAGAATGTGGATTTTGTG
AAGGATGCACATGAAGAAATGGAGCAGGCTGTGGAAGAATGTGACCCTTACTCTGG
30 CCTCTTGAATGATACTGAGGAGAACAACCTCTGACAACCACAATCATGAGG

Sequence ID 315

35 TGGTACAGATACAAACTGGACTCTCAGGACAAAACGACACCAGCCAAACCAGCAGC
CCCTCAGCATCCAGCAGCATGAGCGGAGGCATTTTCTTTTCTTCGTGGCCAATGC
CATAATCCACCTCTTCTGCTTCAGTTGAGGTGACACGTCTCAGCCTTAGCCCTGTG
CCCCCTGAAACAGCTGCCACCATCACTCGCAAGAGAATCCCCTCCATCTTTGGGAG
GGGTTGATGCCAGACATCACCAGGTTGTAGAAGTTGACAGGCAGTGCCATGGGGGC

- 129 -

AACAGCCAAAATAGGGGGGTAATGATGTACGGGCCAAGCACTGCCCAGCTGGGGGT
CAATAAAGTTACCCTTGTA CTTG

Sequence ID 316

5 CGCCACTTATCCAGTGAACCACTATCACGAAAAAACTCTACCTCTCTATACTAAT
CTCCCTACAAATCTCCTTAATTATAACATTACAGCCACAGAACTAATCATATTAA
AA

Sequence ID 321

10 CAGAACAGTACTTTTTAATTTGACCCATGAATTCTATTTAAATTTGTCACTTAATA
TTTAGCCAAGAAGCAAACCATCTAAAAAGATTTCTGGTTTATTTCTCCAACCTCCTA
ATAAATAGGGTCACATATTTTTTAACTTTTTTCTAATTTGAAAAGTAATACAGGCA
TATGGTATTTTTAAAAATGAAACAACACAAAGGGATATGTTTTGAAAAGTGTTCTT
GCCATCCCTGAACTGTAATCATCCCTAACATATTCATACCTGTTTTCATTTTTAAAA
15 GTTGGGTCAGTTTTTTTTATTAGTACATGTATTTCTATCCTACTGATTTATTTGCTA
TATCATCTAATTTAGTTTGAATATTCCATAATTTACTTAATTAGTCCTGTATGGAG
ACCTAGCTCTTCTCAGTGTCTACTATTATAAACAATGCTACAGTGAATATTGGTGT
ATAAATCCATACACACCACGTAACATATCTTAAGTTCCTGGAAGAGATATTGCTAA
ACCAGAAGATAACCTGCATTTTAAAATTTTGA CTGCTAGGGTCAGGGTCACATTTAA
20 ATTAAATTAGAACAAGGAATGCATAATGTCTTCGATAGCAATCTATTCCAGGTGCA
CCGTGGTCACAAAGGAAAGCAAACTGTCAATAACTTTCTTCTCA

Sequence ID 322

TAGCATTTGGCCTTTTAAACATTTGTTTATTTTTTTTTCTGAGAATGGCTAACACA
25 CTTTATTGAGGTTGAAATTAATAAAGAAAATAAAAGAAATGTATCTTCATTCATT
CTGTATGTTAGTGTTTTAATTACCCTTAGAATATATGGATAAAAAATACTATTCTT
TGTCTTGGAGAAGGTAAGAGTCTAGTTAGATGAATAAGGGTTATCTATGTAGAACA
ACTAGAGAATGAGAAGAGAGCTTATGAGATTGAGTACTACGTTATGCAGTAGAGTA
GCACGTCATCTGCTACTGAGTATGGTGTGATAACATTGTGTAACAGGAAAGTATGA
30 TCAATATCTACTTAAATTAAGGACAATATTAGCACTACATTGCTTTATTTTAAAG
TAAAAATTAGAGAACTAAACACAAGCATTGTAAGTACAATAAAAGCTGATCTTTCT
AGTTAAGCAGAATAATACATGTTCAAGCATCTGCTAAATCATTAATATAAGAATA
TAGGGGTTTTCTATAATCTTATTTTCTTTGGAAGAGTACCTCATTTTCAAGANGAG
AAGTTTCTAATTGCCACTTCTTTAAAAATAAAACAGGGTTTTAATGTTCCCAGCAC
35 AAAAATTAATATCTCTTCAAAAAGTCTCTTGTGATTAAAGTTTGAATCCCTTGTCAT
ACTGCTTCTAATATTGACACTGACCTCCTTAGGTATTTTTT CAGGGGTTATAATCTT
TTCTTAAGGTATCTTTTTTCAAGAATTGGATACCTTGGGCTT

- 130 -

Sequence ID 323

CGCGTCGACTTTTAAAGTCATCTCTATAGGAAGGTGCTGGGCAGGGATCCCAGAGA
AAGAAAGGGTCCAAGACTCCATTAAGTCCCTGGATGAAGGGCACTGCTACAGCAG
CTAGTACCAGAGACTCTCCTATCTCACGGTTGAGGCAGACCCAGGATAGAATAGAG
5 AATAAAAGGAATGCTTATAGGAAACAATTTTGTATGGAATGCTAGATGGCCAAGCC
TCAGCCTTTGGTCCAGTGCAACCCCTTGCCTCGCTTGTCAACAGTGAAAAATTAGTT
TGGTTAGAAGAACCATCTGGAAACACACCAGCTTCTGCTACCTTCATGCTCATTGT
TAAAAAAGATTAAACAGTGTGAACATTCTGATCTGTTAATTCCAGGGACTGTTTT
CTTTCCAATGGACTGTTTGTGGTAGAATAACCCCAAAGCTCAAAGCTAAAATG
10 CATCATCAGTCCTAGTCGGCAGTTCCTTAAGAATGGACTGGCGGCGTGGTTGAGCT
GATATGGAAAAGCTGCACCTTCCTGCAGAAGATCAACTGACCTGCTATCCCCACCC
AAATTCAACCTGAGGTATATTTAGTGAAGCAGGTAGCTGTGCTTCTCAAAGCAGA
GAAGCAGTTTTTAAGAACCAAAAAGGTAGAGGAAATCTA

15 Sequence ID 324

GTTTGTTACAGGCAGAATTGGATAGATACAGCCCTACAAATGTATATGCCCTCCCC
TGAAAAAATTGGATGAAAATCTGCACAGCAAAGTGAAACACACAGATAATAGGAA
CAAAATGTAGTTCCCATGTGCCAAACAAAATAAATGAAATCTCTGCATGTTTGCAG
CATATCTGCCTTTTGGGAATGTAATCAAGGNATAATCTTTGGCTAGTGTTATGTGC
20 CTGTATTTTTTTTAAAATGGTACACCAGAAAAGGACTGGCAGTCTACTTCTACCATA
GTTAAACTTCACCCTCTTTAATTTCAACAACATATTCTTTGGAAGCAGGAAGAAATG
CTCATAAAGAGGATCAGACCTTCTTTCCCGTGAAACCAGTATTTGGCGCCATATAT
AAGCCTGGTTAAATTGGTCATCTAAAGCTGTCAAATAAGACATTCTGTGAAAGGTA
AACATCGAACTGGTTATAAGTAAAACCATCAAGCCAACAACAGGGTCTTGAGATA
25 ACCTTTGAAGCTTATTGTCTGGCCTGCACCAGAAGATGTCTGCATTACTCATTGCT
AAAAATGTGTACACAGAACTGCACTAGGATTAATTGGTTCAAGAAGAAATTTAAAC
TTACGTTTGGGTTTCCATACAGCACTCTATTGAATACATGCATCTGAATTTAAGTT
GCAA

30 Sequence ID 325

GACCAGTAATGGCTTTTAAAGAGTCCATTTTGTCAATTGTCTCCCTAGTTAATTACAG
GTGGGGGATCTTTTGCCTCTATTCTCTTCATATTGAAATGAATCATACTCATGTTT
TGTGGAACCTCCTTAAAGTTGTAGCTGTCAATGATCAGATTTTTTTTATATTTCTCA
GCTTAACTCTGCTACTTGATTTACAGTGACCCATAACCTACTCATCCTTGGTTTAT
35 AGTGACACATAATCTTATCTCTTTATAGAACCCTTAAATTTTATCATTATTTTCGCT
TAGAATACAGCATTTCTTTGCTTCTGTTGCTGGTTTGACTTAAGAAATAAGGCAGT
AACTCTGATCAATCAATTATCCATAAGGAAGGGCTTTTCATGGGTTCTATTAATTT

- 131 -

GTTAGTACCCCTAAGTATATCTGAAAAATATGTCTATTGAGAGAAGATTTTGGCATT
CCAGATGGTATAGTCTATATATATATTTAAAGTTTTGAATTTGCTTATATATACTCAG
CTTTCTTTTTCTAGCATTTTTGCATTTACCTGTTAATTGAAGTATACCCCCACAT
ATAAAAGTTCCTCTTAAAGACACTGGACTCTTTCTGGGGGGCTAAAATA

5

Sequence ID - 326

nt: 554

CCCGGAATCGCGGCCCCGCGTCGACAACAAACCTGCATGTTCTGCACATGTATCCAG
GAACTTAAAAAAAAAAAAAGATAGTTTTGTGTGTCTTAATTGAATAATAGTAGATTT
ATAGATTAAAGATCTATGGGTTTTTAATATGGATTANAAATCTGTGGGTTTTTGGAT
10 ATGGATTANAAATCTGTGGGTTTTTAATATGGATTGGAAATCTGTGGGTTTTTAAT
ATGGATTAAAAACATCTGTGGGTTTTTAATATGGATTAAACATCTGTGGGTTTTT
AATATGGATTAAACATCTGGGTTTTTAATATGGATTAAACATCTGTGGGTTTTTAA
TATGGGTTAAAAATCAAAGAAAATGAACTATTTGCTCCAGTGCAGGAAAATACAG
GCAATACTGGATAACAATTAGATGGTCAGGAGCGATAACCCGGTTGCCATTGTTTTGA
15 AGAAGAGAATAAGGNGCTAGCATTCCCTATCCGTAGATAATTTGACAGCTAGGAAAT
AGGGGGAGTCTTCTATGTAGTTAGTGAAGGCTAAATGAACTATTATATGC

Sequence ID 327

CGGCTACCGACAGAAGGACTATTTTCATCGCCACCCAGGGGCCACTGGCACACACGG
20 TTGAGGACTTCTGGAGGATGATCTGGGAGGGGAAGTCCCACACTATCGTGATGCTG
ACGGAGGTGCAGGAGAGAGAGCAGGATAAATGCTACCAGTATTGGCCAACCGAGGG
CTCAGTTACTCATGGAGAAATAACGATTGAGATAAAGAATGATACCCTTTTCAGAAG
CCATCAGTATACGAGACTTTCTGGTCACTCTCAATCAGCCCCAGGCCCGCCAGGAG
GAGCAGGTCCGAGTAGTGCGCCAGTTTCACTTCCACGGCTGGCCTGAGATCGGGAT
25 TCCCGCCGAGGGCAAAGGCATGATTGACCTCATCGCAGCCGTGCAGAAGCANCAGC
AGCAGACAGGCAACCACCCCATCACCGTGCAGTGCAGTGCCGGAGCTGGGCGAACA
GGTACATTCATAGCCCTCAGCAACATTTTGGAGCGAGTAAAAGCCGAGGGACTTTT
ANATGTATTTCAAGCTGTGAAGAGTTTACGACTTCAGAGACCACATATGGTGCAAC
CCTGGAACAGTATGAAATGTGCTACAAAGTGGTACAAGATTTATTGATATATTTCT
30 GATTATGCTAATTTCAATGAAGATCCTGCCTTAAATATTTTTTAATTTAATGGCAN
AT

Sequence ID 328

CAAGACTCCATCTCAAAAAAAAAAAAAAATCTACAGTGCTGAGTATATAAAATTAT
35 TAACACATTTTACAACAATATGTGTTTGTGGAGTTAAATATTTTTTGTCTTTAAAA
CAGGTAATTTTAGTGCATACTTAATTTGATGATTAAATATGGTAGAATTAAGCATT
TTAAATGTTAATGTTTGTACATTGTTCAAGAAATAAGTAGAAATATATTCCTTTG

- 132 -

TTTTTTATTTAAATTTTGTTCCTCTGTAACTAAAAGAACACGAAGTAATTGGTC
ACAATTACTGGTGTTTAACTGCCAAATATGGGTAAATAAGGGAAAATTTTGTTTAA
TATTTAGTCCTTCTGAGATGGCTTGAATATTTGAATTTTGTGTACGTCTATACTG
GGTAGTCACAAGTCTTATAAACACTTTAGAGGAAAGATGGATTTTCACTGTATTT
5 TTAACATCATTTATTTTAAATCTGGTGCTGAAAAATAAGAAAAAAATTAACCTGC
ATTCTGCTGTTCTTCTTTANAAGCATTCCTGCGTAAATACTGCTGTAATACTGTCA
TGCAAAGTGTATCCTTTCTTGTTCGTATCCTTTTTGGGGCAGTGGTTTTT

Sequence ID 330

10 GCGGGAATCGCGGCCCGCGTCGACCTCAAAGGAGAAAAAAACCTTGTAaaaaaag
CAAAAATGACAACAGAAAAACAATCTTATTCCGAGCATTCCAGTAACTTTTTGTG
TATGTACTTAGCTGTACTATAAGTAGTTGGTTTGTATGAGATGGTTAAAAAGGCCA
AAGATAAAAGGTTTCTTTTTTTTTCTTTTTTGTCTATGAAGTTGCTGTTTATTTT
TTTTGGCCTGTTTGATGTATGTGTGAAACAATGTTGTCCAACAATAAACAGGAATT
15 TTATTTTGCTGAGTTGTTCTAAAAAaaaaaaaaaaaaaaaaaaaaaaaaaaaaa
AAAAAAAAAAAAAATTTTAAAATTTTAAAATAAAACCCTTGGTTAT

Sequence ID 331

GCCGCGTCGACCTGCATGAGCCACAGTTTCTTGACTGGAGGCCATCAACCCTCTTG
20 GTTGAGGCCTTGTTCTGAGCCCTGACATGTGCTTGGGCACTGGTGGGCCTGGGCTT
CTGAGGTGGCCTCCTGCCCTGATCAGGGACCCTCCCCGCTTTCCTGGGCCTCTCAG
TTGAACAAAGCAGCAAAACAAAGGCAGTTTATATGAAAGATTANAAGCCTGGAAT
AATCAGGCTTTTTAAATGATGTAATTCCCACTGTAATAGCATAGGGATTTTGGAAG
CAGCTGCTGGTGGCTTGGGACATCANTGGGGCCAAGGGTTCTCTGTCCCTGGTTCA
25 ACTGTGATTTGGCTTTCCCGTGTCTTTCCTGGTGATGCCTTGTTTGGGGTTCTGTG
GGTTTGGGTGGGAAGAGGGCCATCTGCCTGAATGTAACTGCTAGCTCTCCGAAGC
CCTGCGGGCCTGGCTTGTGTGAGCGTGTGGACAGTGGTGGCCGCGCTGTGCCTGCT
CGTGTTGCCTACATGTCCCTGGCTTGTTGAGGCGCTGCTTCAACCTGCACCCCTCC
TTGTCTCATAGATGCTCCTTTTGACCTTTTCAAATTAATATGGATGGGAAAGCTC
30 CTATGCCTTTTGGCTTCTGTTAGAAGGCGGGATGCCCAAGGGTCTGCCTGGGTGT
GGATTGGATGCTTGGGGTGTGGGGGTGGAACTGTCTTGTGGCCCACTTGGGCCC
C

Sequence ID 335

35 CCCGCGTCGACTTTTAAAGTCATCTCTATAGGAAGGTGCTGGGCAGGGATCCCAGA
GAAAGAAAGGGTCCAAGACTCCATTAACCTGCCCTGGATGAAGGGCACTGCTACAGC
AGCTAGTACCAGAGACTCTCCTATCTCACGGTTGAGGCAGACCCAGGATAGAATAG

- 133 -

AGAATAAAAGGAATGCTTATAGGAAACAATTTTGTATGGAATGCTAGATGGCCAAG
CCTCAGCCTTTGGTCCAGTGCAACCCTTGCCTCGCTTGTCAACAGTGAAAAATTAG
TTTGGTTAGAAGAACCATCTGGAAACACACCAGCTTCTGCTACCTTCATGCTCATT
GTTAAAAAAGATTAACCAGTGTGAACATTCTGATCTGTTAATTCAGGGACTGTT
5 TTCTTTCCAATGGACTGTTTGTGGTAGAATAACCCCCAAAAGCTCAAAGCTAAAA
TGCATCATCAGTCCTAGTCGGCAGTTCCTTAAGAATGGACTGGCGGCGTGGGTGAG
CTGATTTGGAAAAC TGCCCTTCTGCAAAAACACTGGCCTGCTTTCCA

Sequence ID 337

10 CAAGACTCCATCTCAAAAAAAAAAAAAATCTACAGTGCTGAGTATATAAAATTAT
TAACACATTTTCAACAACATATGTGTTTGTGGAGTTAAATATTTTTTGTCTTTAAAA
CAGGTAATTTTAGTGCACTTAATTTGATGATTAAATATGGTAGAATTAAGCATT
TTAAATGTTAATGTTTGTACATTGTTCAAGAAATAAGTAGAAATATATTCCTTTG
TTTTTTATTTAAATTTTGTTCCTCTGTAAACTAAAAGAACACGAAGTAATTGGTC
15 ACAATTACTGGTGTTTAACTGCCAAATATGGGTAAATAAGGGAAAATTTTGTTTAA
TATTTAGTCCTTCTGAGATGGCTTGAATATTTGAATTTTGTGTACGTCTATACTG
GGTAGTCACAAGTCTTATAAACACTTTAGAGGAAAGATGGATTTTCAGTCTGTATTT
TTAAACATCATTTATTTTAAATCTGGTGCTGAAAAATAAGAAAAAATTAACTGC
ATTCTGCTGTTCTTCTTTAGAACATTCTGCGTAAATACTGCTGTAATACTGTCA
20 TGCAAAGTGTATCCTTTCTTGTGCTATCCTTTTTGGGGCAGTGTT

Sequence ID 338

CTGGACTGCATGACCAGATCTGATGGGTGAGACTCAGGTGGCATGGAAGAGCCGAA
AGAGGATACCATATGTGGGTGCCGGGGGGGATAGGTGAGAAGTACTAGAAGGCGGA
25 ATGGAAGGACACTTCTGCTCAGCTCTGTGACACGGGCAGGGACCTGCAGGGCTCA
GGTCCTTTAACACAGCAGCTTCATTCTAACACCAGCAGCGTTGGAACACACGTACA
AGTATGCAGACTAAGCTCTTGCTTGGCTGATACGGCTTTTTGGGTTTTTAGAGAAC
ATGCATATATGTTCTCATTCATGGTACATGAACTCAGAAGCCTTACTGCCTATTTT
TGTTAATACTTCTGGGCAAACATTACCACCTTACAACCTCACACCAGTTAGAAATCAT
30 TTGTAAAATGTTATTTAATAAAGCCAAAGAACTAAATCATATTTATTTCCAAGGN
TTTCTAAGATCTCTGAACTAATGAGGTTTTTTAAATCCCCATTAAGTACTCATCA
CTGCTAGTAAAAGCAGTTGTCTTTACCTTTAATTCCAGTGAGTCCCCTTAAATTTA
TTTTTTATTATCTTTGGCTACATTGCCTTAGACAAAATGTGGTCACCCTAATTTAA
NGGATAAAATTCACATCCTCACAGATTTCTTATTAAGAGGGTCTAANCCTTGAATA
35 ATCANCAGTGGAATGGAAGTCTTCTTTACTGGNTTTNATCCTTTCCCTTTTTTAT
CCCATG

- 134 -

Sequence ID 339

TTTTTTTTTAAATAAAGCTGTCGGCACTCAAGGGTAATTTTCATATCAGTGTGNTCT
ACAAGCTGGGGGAAAATGAGTTCTAATTGTTCANAGCTACCAAATCCTTCACCTTTA
GCATAAAGGTTTAAAGATATCACAAAGATGCCAAGTGATTAATAATGTTTTAAACC
5 ACCCCTTTTTCTGTCTGAAAAACAACATAAACAATATTACAACAGTATAGTTACA
GAAGGGTTCTATTTTCATATGTTTTATGCACACTGTGCCTCAAAGGTAATTTAA
ATATATATACTTTTGAGGGGGTGGCTAATGCAGAAACACCCAAGACCTAAGGAAGA
TACAACCCCATTTCTAGGTGTGAGGTCTAAATGCTTCACACACCCACTTGTGACCT
TTTTTCATGAAGAATCATAACACTGTGCAGTGAGAAACAGTGGCAAAGCAATACTG
10 AAAGCATTTTAAATTATTTACTAGGTAAAAGGGTGAAGTACTTTAAATACAT
CAAATTTTCATCAT

Sequence ID 360

GCAAGTGAGAGCCGGACGGGCACTGGGCGACTCTGTGCCTCGCTGAGGAAAAATAA
15 CTAACATGGGCAAAGGAGATCCTAAGAAGCCGAGAGGCAAATGTCATCATATGC
ATTTTTTGTGCAAACCTTGTCTGGGAGGAGCATAAGAAGAAGCACCCAGATGCTTCAG
TCAACTTCTCAGAGTTTTCTAAGAAGTGCTCAGAGAGGTGGAAGACCATGTCTGCT
AAAGAGAAAGGAAAATTTGAAGATATGGCAAAGCGGACAAGGCCCGTTATGAAAG
AGAAATGAAAACCTATATCCCTCCCAAAGGGGAGACAAAAAGAAGTTCAAGGATC
20 CCAATGCACCCAAGAGGCCTCCTTCGGCCTTCTTCCTCTTCTGCTCTGAGTATCGC
CCAAAAATCAAAGGAGAACATCCTGGCCTGTCCATTGGTGATGTTGCGAAGAACT
GGGAGAGATGTGGAATAACACTGCTGCAGATGACAAGCAGCCTTATGAAAAGAAGG
CTGCGAAGCTGAAGGAAAAATACGAAAAGGTA

25 Sequence ID - 361

nt: 622

CTGTNATNGAATCTGCTTGTNACTNAAATGCTAAACTCAATTCTGTAATTCAATAG
GTGCACCTNTCTGAGAAACATANNAGACAATGAGGAAAAGGATTCANCATTCCGTG
GAATTTGTACCATGATCAGTGTGAATCCCANTGGCGTAATCCAAGTAAGATGTTCA
CAAAGATTTGTTTTTAATGTCTAATTAATAAAATTTTAAAGGAAGAAACATTCTAA
30 TACTTTAATTATAAAAAGTTAACTATTTTCAAAGGTATCAAATACAGTTAAACCT
TTAAATGTATATTTCTTAATATCTTGAAATTGTAATGCCTTTTTTTTTTCTAAA
TTTTTTTTTGTTCATGAAATGAGATAGTAACAGCAGATTGGGACAACAAGGTTATATT
CTTGTCTTGAATCAGGCCATGGCTTCTTTCATCCAAATTTTACAGACCTCATTTATTT
ACTTTGTCCCTGCCTCCCATCCCTGGATATCANGTTTGTGGATATCTACAGTTAAT
35 AGAGTGACCAAATAGTAGGAATACTGTCTCTCTATTCTGAATAAAATACTTTGAAT
CAGATTTAGAAATAATGAATAAAATACAAATCACCATTGAAATTGCTCTAATTTTG
AGAGCT

- 135 -

Sequence ID - 363

nt: 628

ATCACNTGAGGCAAGAGTTTGAGCCAGCCTAGCTAACATGGTGAAACCCCATCTCT
ACAAAAATATAAAAATTAGCCTGGGTGGTGATGGGCACCTGTAACCCAGCTACTC
GGGAGGCTGAGGTAGGAGAATCACTTGAACCCGGGAGATGGAGGTTGCAGTGAGCC
5 AAGATCGTGCCACTGCACTCCAGCCTGTGTGACAGAACAAGACTCTGTCTCAAAAA
AAAATAATAATAATAATAATAAAAAAGGAATAACATAGCTAGGAATAAATTTAA
TCAAAGAGGTGAAAGACTTATACACTTAAAACTACAAAAAAAATCACTGAAGGA
ATTATAGACCCAAATAAAAATAAAAAGACATTCTGTGTTTTAGGGAAAGAAG
ACTTAATATTGTTAAGATGTCAATACTACCCAAAGTGATCTACAGATTCAACATAA
10 TCCCTATCAAAATTCCAACAGCCTACTTTGTAGAAATGGAAAAGCCAATTTTCAAA
TTCAGATGGAATTGCGAGGGGTCTGAATAACAAAAACAATCTTGGGGAAAAAAA
CAAAAAACAAAGTCAAAGAACTCACACTTCTCTATTTATAAATTTACTACAAAGTT
ATAGTAATCAAA

15 Sequence ID - 364

nt: 528

TGAACATCCAGCCATGTCATTTCTTCCATTCCCTGCCCTGGAGTAAAGTAGATTTAC
TGAGCTGATGACTTGTGTGCATTTGTACATTGCAACCTTAGCTTACCTCTTGAAGC
ATGTAGAGCATTCATCACCCACCATTCATTCACTGCCTACTCCCACCACAGCTGTT
TCGTGGTCTGTCTGCTCCCTGTGCCACCCCCACCCCATCAGGTGGGCCTTTTGCAA
20 GTGATGAAGTCACCTGTGGGGGAAGAGCTTTCCTTTCTCTCTCAACTCAGAAGG
CCTCTTCTCTTGCTCAAGAGGGTGCTGCTGCTTTCTGCCTCCTTCCCCGGCCGGC
CTCCATCCCAGTTCACCTTTTCAGAAATGGCCCCTCAGTCAACTCTTCCCTTTTCT
CCTGGCTTTTTATTCTCTCCAGTCTCTTAAGAGTATCCTTAGCTTTAAAAACAATA
ACACAGAGGATGGGTGCAGTGGCTCATGCCTGTAATCCCAGCACTTTGGAGCCTGG
25 GGCGGGCGGATCACTTGAGGNCA

Sequence ID 365

GTCCCGGAATCGCGGCCGCGTCGACCTTTTCTATGCCTGCTATATAAACAGTACCT
TGCAAGATGTCCTGTCTGATATCCACAAAGGGGTATTGTCAACCCCAAGTTCAGAC
30 AGCTTTGTATTCTTCTGTCCCTGGATACATGAATTACTGCCATCTTTACACAGCGC
CCTAAAATACCAACGCGAAGTTACCTGCTCAGCTTGAAGCTGCGCTGTACCCTGGA
ACCAGCACTTCTGCTGAATGACTCAGGATGAAGCCTCGACTTCTCCTTCCCATCCC
ATGCCAGACCCAGTGGCTCCTTTCCCAATCTGATCCAGTGACTTTAAGTCCAGC
TGTTGCAACCTGGGCATGAGGAGGAGTGCAAGATGGCTTTGTCCTACCTGGAAAGA
35 GGCTTTCTGGA

Sequence ID 366

- 136 -

CACCATTTACACACAGTGGGTCTTGAATAGCATCGTTTTATTCAATGTCATTTTG
TTATAACATTGAGAAAAAATTGATTCCCGGCTGGGGCCACTGTCTGTGCACCGT

Sequence ID - 368

nt: 329

5 GAAAGATCTAAAATCGACACCCTAACATCACAATTAAAAGAACTAGAGAAGCAAGA
GCAAATTCAAAAGCTAGCAGAAGGCAAGAAATAACTAAGATCAGAGCAGAGCTGAA
AGAGATAGAGACACAAAAAACCATTCAAAAAAACAATGAATCCAGGAGTTTTTT
TTTTAAAAGATCAACAGAATTGACAGACTGCTAGCAAGACTAATAAAGAAGAGAG
AAGCATCAAATAGACTCAATAAAAAATGATAAAGGGGATATCACCACCAATCCAC
10 AGAAATACAACTACCATCAGAGAACACTATAAACACCTCTATGCAAAT

Sequence ID 369

GAAAGATCTAAAATCGACACCCTAACATCACAATTAAAAGAACTAGAGAAGCAAGA
GCAAATTCAAAAGCTAGCAGAAGGCAAGAAATAACTAAGATCAGAGCAGAGCTGAA
15 AGAGATAGAGACACAAAAAACCATTCAAAAAAACAATGAATCCAGGAGTTTTTT
TTTTAAAAGATCAACAGAATTGACAGACTGCTAGCAAGACTAATAAAGAAGAGAG
AAGCATCAAATAGACTCAATAAAAAATGATAAAGGGGATATCACCACCAATCCAC
AGAAATACAACTACCATCAGAGAACACTATAAACACCTCTATGCAAATAAAGTAG
AAAAT

20

Sequence ID 370

GAAAGATCTAAAATCGACACCCTAACATCACAATTAAAAGAACTAGAGAAGCAAGA
GCAAATTCAAAAGCTAGCAGAAGGCAAGAAATAACTAAGATCAGAGCAGAGCTGAA
AGAGATAGAGACACAAAAAACCATTCAAAAAAACAATGAATCCAGGAGTTTTTT
25 TTTTAAAAGATCAACA

Sequence ID 371

GCCCGGAATCGCGGCCGCTCGACGTAAGCTCGGCTGAATCCACGGTTCAAGAACA
GGAAAGAAGGCCAAGGCATAGGGAGTGGGGCAGTTGGGTGAATATTAGTACCTTTC
30 CCTCAGNTNCATTAATTACCCCTGCCTACTCTGCACAAAAGGATNTAACAACAGTT
TCCTTTTTAATGGCCAGGTACAGCTGCTTATATGGANGGGCATTTNTNAATGATAT
CCTTNATCACTGTCTTAATCATCACATNCTTAAACAATCACTTTATTGTGTTAAG
GAAGATAAAAATGGCTGGGTTCATTTCCGTTCTGGAAGAAATCGANTNAAAAGGT
AACCATTTAATAATGCANAGGGCANTTTCACTGCAGACCCTAATACTGGAAATTTT
35 TAAAAACAAATGAAAACTTCTACTTTTTCTTCTAAGCTTACTTAACCACCCAAAT
TTTCCAGCCACATATCTTCCTAGTCTACAACCTGCCTTTAAGAGATGCTCA
AAAAAATGTAAATTCTCAAATACATTCTTATTACAATTACTGCTAACCT

- 137 -

Sequence ID 373

CCAGTGTGCTGGGATTACAGGCATGAGCCCTGCACCCAGCCTCTTAAACTGATCAT
ATGATATTGGTTCTCAACCAAGGGTGACTTTGCCCCCAGAGGATACTTGGCAATGT
CTGGAGATACTCAGTTGTTCATGACTTGGACAGGTGCTACTGTACCCAGTGGGTAG
5 AGGTCAGGGATGGTGCTAAACATAGGACAGCTGTCAAGAGAAAAGAATGTACCCAG
CCCCAAATGTCAGTAGGGCTGAGGTTGAGAAACCCAGCTGTAGCTGACGTGTGAAG
GACAGACTGGCCTGGAAGTGTGTTTTCTGCCCTTTCCACCCCTGCATATTAGTTA
AGGCCAAAGGAAAAAAGGAATGCAGGAAATGCCCGTTAAAAATCTTCAAAACAATA
TAAATGATCAATTCCACTAAAACCCCTTTACACATTTAAGTATAAAGGTATTGGTA
10 GGAAAATTTGTTATTCACTGCTTTTCTCAGTGTTCATGAAATAATTATTTCTGCTGT
CAGTTT

Sequence ID 374

AAAAAAAAAATCACTGAAGGAATTATAGACCCAAATAAAAATAAATAAAAAGACAT
15 TCTGTGTTTTAGGGAAAGAAGACTTAATATTGTTAAGATGTCAATACTACCCAAAG
TGATCTACAGATTCAACATAATCCCTATCAAAATCCAACAGCCTACTTTGTAGAA
ATGGAAGGCAATTTTCAAATTCAGATGGAATTGCGAGGGGTTCTGAATAACAAA
AACAATCTTGGGGAAAAAAAACAAAAACAAAGTCAAAGAACTCACACTTCTCTAT
TTATAATTTACTACAAAGTTATAGTAATCAAAGTCGACGCGGCCGCGATTCCGGG
20

Sequence ID 378

CGACTGCGGCTCTTCCTCGGGCAGCGGAAGCGGCGGCGGTCGGAGAAGTGGCCT
AAAACCTTCGGCGTTGGGTGAAAGAAAATGGCCCGAACCAAGCAGACTGCTCGTAAG
TCCACCGGTGGGAAAGCCCCCGCAAACAGCTGGCCACGAAAGCCGCCAGGAAAAG
25 CGCTCCCTCTACCGGCGGGGTGAAGAAGCCTCATCGCTACAGGCCCGGGACCGTGG
CGCTTCGAGAGATTGTCGTTATCAGAAGTCGACCGAGCTGCTCATCCGGAAGCTG
CCCTTCAGAGGTTGGTGAGGGANATCGCCAGG

30 Sequence ID 380

GCAATTTAATTTTTAATAACAAAGATACTGTATTTTAACATGGTGAAATATACTTG
GCTAAGTCCAGATTAAAAAAAAGTATCTAGCCCAACAGTACAATTATACAGCT
TTGTACAGAACATTCCATAGATCAACAGAAAATACATTTGAGCGCAAAAATAAAAA
ATATTTAAGGAGAATCTCTAAGCAGCATTTTATTTCTGCAAAAGACATATCTTGTC
35 TGATTAAATATCTACAAGTGCTTTTCCTTTCAAAAATACATATATTCTTAATAGAC
TAAGTCATTAACAATGACCTGGTAATTCTTTCACCTTCAATTTGAATGATTTATAAG
CTAAATCTTCAACCACAAAAGGTTTTTATTTGTATTAAAGATGTTACCACTTTTGA

- 138 -

CAAAAAGCTTAAAATATTTTATATTTCAAAGGAAAATTAGCAACATAACTTTACAA
TATATTTCTATGATATTTTGATTGTGAGGGCTACTCTATTTAAAACTGATGATCTCT
GTTGTGTTGCTCAGATGCAGGAAAGCAGCAAAA

```

5      Sequence ID - 381                               nt:      534
      GACTTANATCTAAATGGACCACATTCTCTACTTAAAAAATGCTATTAACCATGTG
      ATCTTCTCAGTCATGAGGTAATCTGGTGACTACCCTTCCTCAAAGCCAGTTGGGAT
      ATTCTTTGAATAGAGTAAACAGTGTTTCTAGGCTGGGAGACACCAGACATAGTTG
      AGGACAGAGGTGCTAGAAAATAGGAAGTTTAAAAGCATGTGCGGTGATGCTCAGAG
10     GAGGTAAACCCCAACCTCATGCTCATAGCTTCCAATCATTTTCTCTAGTTCTTAAC
      TCTTAAATGTGAGAAATGCTTGAAGATTCTAGTCATCTGAAGAAAGTCTCTTTATT
      AAAGATTTTTCATAAAAGAGACCAAAGCAGACAAACAGAAAAAGACATCTTGGGGAA
      AAAAAAAGGATAATGGGAAGAGAAGGAAAGTTTTTAAAAATTATCAATATCCTCAG
      GGGGACAAAATATTATATCCTATAAAGACAGATTTTTTATTTTTTAAAAAAATAGAA
15     AGCAAAACAAGCTCCTAAAAATAAAGTTTG

```

```

Sequence ID - 382                               nt:      444
GTTAAGGAAGTCAGCACTTACATTAAGAAAATTGGCTACAACCCCGACACAGTAGC
ATTTGTGCCAATTTCTGGTTGGAATGGTGACAACATGCTGGAGCCAAGTGCTAACA
20  TGCCTTGGTTCAAGGGATGGAAAGTCACCCGTAAGGATGGCAATGCCAGTGGAAAC
ACGCTGCTTGAGGCTCTGGACTGCATCCTACCACCAACTCGTCCAACTGACAAGCC
CTTGCGCCTGCCCTCTCCAGGATGTCTACAAAATTGGTGGTATTGGTACTGTTCTTG
TTGGCCGAGTGGAGACTGGTGTCTCAAACCCGGTATGGTGGTCACCTTTGCTCCA
GTCAACGTTACAACGGAAGTAAAATCTGTGCAAATGCACCATGAAGCTTTGAGTGA
25  AGCTTTTCTTGGGGACAATGTGGGGCTTCAATGTCAAGAATGTGTCTGTCAAG

```

Sequence ID - 383 nt: 566

CTTTGAAGAACTTTGCCAAATACTTTCTTACCAATCTCATGAGGAGAGGGGAACATG
CTGAGAAACTGATGAAGCTGCAGAACCAACGAGGTGGCCGAATCTTCCTTCAGGAT
30 ATCAAGAAACCAGACTGTGATGACTGGGAGAGCGGGCTGAATGCAATGGAGTGTGC
ATTACATTTGGAAAAAATGTGAATCAGTCACTACTGGAAGTGCACAACTGGCCA
CTGACAAAAATGACCCCCATTTGTGTGACTTCATTGAGACACATTACCTGAATGAG
CAGGTGAAAGCCATCAAAGAATTGGGTGACCACGTGACCAACTTGCGCAAGATGGG
AGCGCCCGAATCTGGCTTGGCGGAATATCTCTTTGACAAGCACACCTGGGAGACA
35 GTGATAATGAAAGCTAAGCCTCGGGCTAATTTCCCCATAGCCGTGGGGTGACTTCC
CTGGTCACCAAGGCAGTGCATGCATGTTGGGGTTCCTTTTACCTTTTCTATAAGTT
GTACCAAAACATCCACTTAAGTTCTTTGATTTGTCCATTTCCTTCAAATAAAGAAAT

- 139 -

TTGGTA

Sequence ID 384

TTTTGGGGTTTATATATAAGCCTGGTTCTTGCTGAAACTGCTTATGTTGATAACCA
5 GTTAGTGAGTTCCTCTCTATTGACTTGCTGGGAAGTTTATAGAGACATTTTTTATG
CATTTCAGAGATTTTCAGTACAAATCTTGAAAAAGGGACATTTAGGCCGGGCGCGGTG
GCTCACATCTGTAACCCTAGCACTCTGGGAGGCTGAGGTGGGTGGATCATGAAGTC
AAGAGATAGAGACCATCCTGGCAAAAATTAGCTGGGCGTGGTGGGGTGCGCCCGTA
GTCCCAGCTACTCGGGAGGCTGAGGCAGGAGAATTGCTTGAGCCCGGGAGGCGGAG
10 GTTTCATTGAGCCGAGATAGTGCCACTGCACTCCAGCCTGGACAACAGAGCGAGAC
TGTGTCTT

Sequence ID 386

CTAAGGGTTTAAAGATGGAAAGAGGCATTGATGAACAGCTGGGGAAGGAGTAGTTT
15 GAGGTAGATGTGCAGATGGAATGAAGAGAAGGTCTCAAGAAGAGGGTGGAGCCAAA
GAGGGCTGCAGATTTAGAAGGCTAAAGTCTTTAGATGGCTTTGGATAGCCTGTTGT
ATCTTGACCATGCAGGTTACAGTGGAGCATGGAGTGGGGACAGAAGTGGAGGAAG
GAACCAGGGAACATGGAGTGAGAAGCTAAAGGAAAGTGATGCAGTAGATACATGGC
TCTAAAGTACTCAGGACTTTCAGAGGCTTAAACATAGGGTGACCAACTATCCCCT
20 ATGCCTGATACTAAGGGCATTCCCTGGATGTGGACCTTTCATTCCCCAAATTAGGA
AAGTCTTGGGCATACCAAGACAAGTTGGCCACCCTACTCAAAGTATGTAAGCTAA
CATATCTGTTCTCTAAGAGGTTAAAGCTGGATGGGGATACCAGATGTATGTACGTG
ATGCAGTTAAACAGCAATACAAGGGGGCAAGTCTACCTGATCGGCCAATTCAATGG
GA

25

Sequence ID 387

GAAGCCAAACCAAAGGAGCTTCTACTTCATGATGCCATTTATGTAAAGTTCAGGCA
GAGAAAATCAGTGGTTTAAAGAAGTTAGAATAATGATTATCTTTGGAGGGATTGCAA
CTGGAAGAAGTCATGATTGGGATTTCTGGGTCCTAATAGTGCTCTGTGTCTTGATC
30 TGAGTGCCGACTACATGAGTGGTTAGGTTTGCAAAATTCATTGAGTTATGCACTTA
ATGGTGTTGTCTTATTAGAGCTGATGGAGGAGAGAGGGCTTCAATTTGCACAACCTG
AGTAATCAGCTAGGCCCAGTCACTAGGTGAACAACTTACTGCTCCAATCAGCCTTA
GAGCAGGAATCAAACCTCATGTCTCAGAAAAGTTATTAATTCAGCTTGTCTTGGGAC
TTCCTTCAGAGTCACTCTTGAATAGCTGAAATAGTAAATGTTAAATCTGTGGATGC
35 AAGTGTGTAAATTATTTTAGTCATCAGCTCTAATAAGATGGCCTTTGGGGAAATGA
GTATAAGGTCACGAAAATGAAATGGCAAGAAGGAGGTCTACTATTTCTTCTGTAAT
ACTGATTTTTACCCCATCAGGGTCAGTCCCAGAGGTGTAAATGTGAAGCTTG-T

- 140 -

CTTTTCTTTAATAA

Sequence ID 388

CTTTGGACACTAGGAAAAACCTTGTAGAGAGAGTAAAAATTTAACACCCATAGT
5 AGGCCCTAAAAGCAGCCACCAATTAAGAAAGCGTTCAAGCTCAACACCCACTACCTA
AAAAATCCCAAACATATAACTGAACTCCTCACACCCAATTGGACCAATCTATCACC
CTATAGAAGAACTAATGTTAGTATAAGTAACATGAAAACATTCTCCTCCGCATAAG
CCTGCGTCAGATTAAAACACTGAACTGACAATTAACAGCCCAATATCTACAATCAA
CCAACAAGTCATTATTACCCTCACTGTCAACCCAACACAGGCATGCTCATAAGGAA
10 AGGTTAAAAAAGTAAAAGGAACTCGGCAAATCTTACCCCGCCTGTTTACCAAAAA
CATCACCTCTAGCATCACCAGTATTAGAGGCACCGCCTGCCAGTGACACATGTTT
AACGGCCGCGGTACCCTAACCGTGC AAAGGTAGCATAATCACTTGTTCTTAATTA
GGGACCTGTATGAATGGCTCCACGAGGGTTCAGCTGTCTCTTACTTTTAACCAGTG
AAATTGACCTGCCCCGTGAAGAGGCGGGCATAACACAGCAAGACGAGAAGACCCTAT
15 GGAGCTTTAATTTATTAATGCAAACAGTCCTAACAAACCCAGGTCCTAAACTCCA
AACCTGCATTAAA

Sequence ID 389

CGACCCGGAATTTCGCGGCCGCGTCGACTGAGTTCTTGACAAGAGTGTTTTTCCCTT
20 CCCGTCACAGAGTGGGCCCCAACGACCTACGGCACTTTGACCCCGAGTTTACCGAAG
AGCCTGTCCCCAACTCCATTGGCAAGTCCCCTGACAGCGTCCTCGTCACAGCCAGC
GTCAAGGAAGCTGCCGAGGCTTTCCCTAGGCTTTTCCTATGCGCCTCCCACGGACTC
TTTCCTCTGAACCCTGTTAGGGCTTGTTTTAAAGGATTTTATGTGTGTTTTCCGAA
TGTTTTAGTTAGCCTTTTGGTGGAGCCGCCAGCTGACAGGACATCTTACAAGAGAA
25 TTTGCACATCTCTGGAAGCTTAGCAATCTTATTGCACACTGTTGCTGGAAGCTTT
TTGAAGAGCACATTCTCCTCAGTGAGCTCATGAGGTTTTCATTTTTATTCTTCCTT
CCAACGTGGTGCTATCTCTGAAACGAGCGTTAGAGTGCCGCCTTAGACGGAGGCAG
GAGTTTCGTTAGAAAGCGGACGCTGTTCT

30 Sequence ID - 390

nt: 523

GAATCCCTAGAAAAAGAGAATTCCCAACTTGATGAGGAAAACCTTAGAACTGCGAAG
GAATGTAGAATCTTTGAAGTGTGCAAGCATGAAAATGGCTCAGCTACAGCTAGAAA
ACAAAGAACTGGAAAGTGAAAAAGAGCAACTTAAGAAGGGTTTGGAGCTCCTGAAA
GCATCTTTCAAGAAAACAGAACGCTTAGAAGTTAGCTACCAGGGTTTAGATATAGA
35 AAATCAAAGACTGCAAAAACTTTAGAGAACAGCAATAAAAAAATCCAGCAATTAG
AGAGTGAAC TACAAGACTTAGAGATGGAAAATCAAACATTGCAGAAAAACCTAGAA
GAACTAAAAATATCTAGCAAAAGACTAGAACAGCTGGAAAAAGAAAATAAATCATT

- 141 -

AGAGCAAGAGACTTCTCAACTGGAAAAGGATAAGAAACAATTGGAGAAGGAAAATA
AGAGACTCCGACANCAAGCAGAAATTAAAGATCCACATTTGAAGAAAATAATGTGA
AGATTGGAAATTTGGAAAA

[illegible]

Sequence ID 394

20 GACCCGGAATCGCGGCCGCGTCGACCATTTTAGCCAAGGTGCCTCTATAGGGGTCA
AGACATCATGTGCCAGACCTAAGGTCAGGAATGTCATATTTTTCTGTAAAATCA
TTTTATTTCTGTGTATCTTACCTTTAAATCATTGTGGTTTACTCTGAGATTCTGTA
GTCCTAATATTGTATCATTGTGCTGTCTGCAAAACAACTTGAATCTATTTTGTTTG
CATCTTTTGTTACATGTAACGCAGCTGTACTTTATGTTCTTTGCAACTGTTTCCAT
25 TATGAGAACGCTGTGCTATTTACAAGGTTACATTTTTCTTGGCCAGGCGAGGTGGT
CATGCCTGTGATCCCAGCACTTTGGGAGGCCAAGGTGGGCGGATCACTTGAGGTAA
AGAGTTGAGACCAGCCTGGCTAGCATGGCGAAGCCCAGTCTCTACTAAAAATACAA
AAATTGGCCGGGTGAAATTAGCCGGGCGTGGTGGTGTGTGCTTGTAATCCCAGCTA
CTCGGGAGGCTGAGGCAGGAGAATCGCTTGAATCCGGGAGGCAGAGGTTGCAGTGA
30 GCCAAGATCANGCCACTGCACTCCACCTCGGGGTCAAGAGCGAAACTCTGTCTCAA

Sequence ID 395

35 CCGTTTTAGTCAGGATGGTCTCGATCTCCTGACCTCGTGATCCGCCTGCCTCGGCC
TCCCAAAGTGCTGGGATTACAGGCGTGAGCCACCGCGCCGGCGTAAATCAGGTTT
TTTAAATGTTTGCCAAACCTTATCACTGACTTTTATAACAAAATTATTTACTATAA
TCATTAGGGAATATTTAAGTTCTGCTAATACTTAAATTGCAGAGTGCTAAAACCA
GCAGTGAGTTTAGAATCAAGCTAAGCTTTATTGTTGCTACTATTTGAGGCATATTA
GTTGACTGGTGTTTCATATGCAAGGCAGTCTACTGGGTGCAACAAGGGTTAGAAGGA

- 142 -

TATTTTAAAAAACTGACCCTATTCTCAGGATGAAAATAATACACTAGTAATAGTC
TGCTCTGTTGGTTAACTCCTCGTAAGGAGGTCAATTAAAAATGCTGTAGTGTTGCAA
GGGAAGGAGAGGAAGAATCATATTCCTTCACTAGCAGGATCAAGAAAGCTTTTATA
GAAATATACAAAATCTTCACCTTCTTGAAGGATTGGTAAAAATTTAATAGCCAACATT
5 GGGCACTTATTCACTTCTCTGAGTAAATATTTATTGCAT

Sequence ID 396

CTTAAATCTAAATGGACCACATTCTCTACTTAAAAAAATGCTATTAACCATGTGAT
CTTCTCAGTCATGAGGTAATCTGGTGACTACCCTTCCTCAAAGCCAGTTGGGATAT
10 TCTTTGAATAGAGTAAAACAGTGTTTCTAGGCTGGGAGACACCAGACATAGTTGAG
GACAGAGGTGCTAGAAAATAGGAAGTTTAAAAGCATGTGCGGTGATGCTCAGAGGA
GGTAAACCCCAACCCTCATGCTCATAGCTTCCAATCATTTTCTCTAGTTCTTAACTC
TTAAATGTGAGAAATGCTTGAAGATTACTAGTCATCTGAAGAAAGTCTCTTTATTA
AAGATTTTCATAAAAGAGACCAAAGCAGACAAACAGAAAAAGACATCTTGGGGAAA
15 AAAACAAGGATAATGGGAAGAGAAGGAAAGTTTAAAAATTATCAATATCCTCAGG
GGGACAAAATATTATATCCTATAAAGACAGATTTTATTTTTTAAAAAAATAGAAA
GCAAAACAAGCTCCTAAAAA

Sequence ID - 397

nt: 534

20 GACCCGGAATCGCGGCCGCGTCGACGGAAGCTCCTGCCCCCTCCTAAAGCTGAAGCC
AAAGCGAAGGCTTTAAAGGCCAAGAAGGCAGTGTTGAAAGGTGTCCACAGCCACAA
AAAGAAGGAGATCCGCACGTCACCCACCTTCCGGCGGCCGAAGACACTGCGACTCC
GGAGACAGCCCAAATATCCTCGGAAGAGCGCTCCAGGAGAAACAAGCTTGACCAC
TATGCTATCATCAAGTTTCCGCTGACCACTGAGTCTGCCATGAAGAAGATAGAAGA
25 CAACAACACACTTGTGTTTCAATTGTGGATGTTAAAGCCAACAAGCACCAGATTAAAC
AGGCTGTGAAGAAGCTGTATGACATTGATGTGGCCAAGGTCAACACCCTGATTTCGG
CCTGATGGAGAGAAGAAGGCATATGTTTCGACTGGCTCCTGATTACGATGCTTTGGA
TGTTGCCAACAAAATTGGGATCATTTAAACTGAGTCCAGCTGCCTAATTCTGAATA
TATATATATATATATATCTTTTCACCATAA

30

Sequence ID - 398

nt: 512

GGGGAGCCCCCTCTTCCCTCAGTTGTTTCTACTCAGACTGTTGCACTCTAAACCTA
GGGAGGTTGAAGAATGAGACCCTTAGGTTTAAACACGAATCCTGACACCACCATCT
ATAGGGTCCCAACTTGGTTATTGTAGGCAACCTTCCCTCTCTCCTTGGTGAAGAAC
35 ATCCCAAGCCAGAAAGAAGTTAACTACAGTGTTTTCTTTGCACCGATCCCCACCC
CAATTCAATCCCGGAAGGGACTTACTTAGGAAACCTTCTTTACTAGATATCCTGG
CCCCCTGGGCTTGTGAACACCTCCTAGCCACATCACTACAGTACAGTGAGTGACCC

- 143 -

CAGCCTCCTGCCTACCCCAAGATGCCCCCTCCCCACCCTGACCGTGCTAACTGTGTG
TACATATATATTCTACATATATGTATATTAAACTGCACTGCCATGTCTGCCCTTT
TTTGTGGTGTCTAGCATTAACCTATTGTCTAGGCCAAAGCGGGGGTGGGAGGGGAA
TGCCACAG

5

Sequence ID 399

TTTTGGCATTACTTAATCCAATTATAAAAACTGAATTTTTTAAAAACAGCACTTGT
TTTTTCTTCCAAGATTAAATTTGAATTTTTTTATGGACATTAGAAAACATTGCAGTT
TAGTCATAATCAAAAATAAATCTTGAGGCTGGTAGAGCAGCTTTGTGTGCTGTTTAT
10 ATTTTTATTGCTTACTGGATTTTCAGTGTTACCTAGTGCCATCAGTTTGGTATTTTG
CCACCTTGACATTTCAGTGATGTTTGATTTTTCTTTTTCTTTTTTCATATTACT
TTTAAATCCTGAATAGTTTGTGGCAGCTGGAGATCACCTAGTCCACCACTGTCCAA
CATGGCAATGGTAAGTAATATTGAGTAAAGAATAGAAAATTAGTAAAATGCATGGC
TTCAGAAATTATAGCAATTTGCAAAATAGGTTAATGGATGAAAATTAGAATGACCAG
15 TTTAACTTTCCCCCAGCAGATTCTTCTGTAAACAATGCCCCCTTCAAAATAAAGG
AAGAACAAGTGGGTGTTATACCTATGTTATTTGGCTATGTTAGCACAATATGATGG
ACTAATTTGAGAAAAGCATTCTTCTTTACTATTACTTCTTTTCTTTATAGGG
CTAAGTCTGCCTTCTGGGTCTTTGAA

20

Sequence ID 400

GAAGAAGCGCGAAGAGCCGTTAGTCATGCCGGTGTGGTGGCGGCGGCGGAGACTGC
GGGCCCCGTAGCTGGGCTCTGCGAGGTGCAAGAAAGCCTTTGAGGTGAAGGTGTATG
AAAGTCATCATAACAGATGTTTTCCAAAACTTGTAGAAGGTTGTGAAAAAACTAC
TAGGATCACGCGGCATGTATTGAGCATATAGGTTGCTGTAGATGAATGTTCTTAGC
25 TGTCATGTTTAAAAATACTTCTGCTTCGTTACCTCAAGTGTGGCATGCAGCATTTT
GGAAGGAAAATTGAAGACGTGTTCAAGAAAACATGAACAGAAGCAAATGATGAAAA
TGAGCATTTTACTTGATGTTGATAACATCACAATAAATTATGGAGAAAAATACATA
TTTGGCTAACTTTTAATTGCTGAACAATAAAGTGTTTTCTTTTAAATCNAAAAA

30

Sequence ID 401

GAAGCCAAACCAAGGGAGCTTCTACTTCATGATGCCATTTATGTAAAGTTCAGGC
AGAGAAAATCAGTGTTTAAAGAAGTTAGAATAATGATTATCTTTGGAGGGATTGCA
ACTGGAAGAAGTCATGATTGGGATTTCTGGGTCTTAATAGTGCTCTGTGTCTTGAT
CTGAGTGCCGACTACATGAGTGGTTAGGTTTGCAAAATTCATTGAGTTATGCACTT
35 AATGGTGTGTCTTATTAGAGCTGATGGAGGAGAGAGGGCTTCAATTTGCACAACT
GAGTAATCAGCTAGGCCAGTCACTAGGTGAACAACCTTACTGCTACCAATCAGCCT
TAGAGCAGGAATCAAACCTCATGTCTCAGAAAAGTTATTAATTCAGCTTGTCTTGGG

- 144 -

5 ACTTCCTTCAGAGTCACTCTTGAATAGCTGAAATAGTAAATGTTAAATCTGTGGAT
GCAAGTGTGTAAATTATTTTAGTCATCAGCTCTAATAAGATGGCCTTTGGGGAAAT
GAGTATAAGGTCACGAAAATGAAATGGCAAGAAGGAGGTCTACTATTTCTTCTGTA
ATACTGATTTTTTACCCCATCAGGGTCAGTCCCCAAAGGTTGTAAATGTGAAGCTTG
GTCTTTTTCTTTA

Sequence ID 402

10 GACCCTATTCTCAGGATGAAAATAATACACTAGTAATAGTCTGCTCTGTTGGTTAA
CTCCTCGTAAGGAGGTACAATTAAAATGCTGTAGTGTTGCAAGGGAAGGAGAGGAA
GAATCATATTCCTTCACTAGCAGGATCAAGAAAGCTTTTATAGAAATATACAAAAT
CTTCACTTCTTGAAGGATTGGTAAAATTTAATAGCCAACATTGGGCACTTATTCAT
TCTCTGAGTAAATATTTATTGCATGCTTATCTTGTATCAACATTGNGATGAAAGCN
CAAGAATGAAAGAGGAGGGAGAATGTTTANAGAATAAGGCTGAAACACAGATTTTG
TAGGGAGCGTAGGGGAGACTGANAAAACAG

15

Sequence ID 403

20 AAGACACCTGATAGATTGTCTTGTATTATTTTTCTTTGCCTTCTTACAATCTCAG
TGATTAGAATTGGGCTGAAAACAATACATCAAATTCTCAGCAAAATCCTTATGGGT
TGCTGGATACCGAGGGTTTTTAAGATCTTTAGACTTCACTATATAGAACAAATGTT
GAATGGGAATTTTCTTTATTTCTATANCGTTTNG

Sequence ID 405

25 CCCGGAATCGCGGCCGCGTCGACGATGAGCATTTTTTTCATGTGTCTTTTGGCTGCA
TAAATGTCTTCTTTTGAGAAGTGTGCGTTCATATCCTTTGCCCACTTTTGTATGGG
GTTGTTTTTTTTCTTGTAATTTGTTTGAGTTCATTGTAGATTCTGGATATTAGCCC
TTTGTGATGAGTAGGTTGCGAAAATTTTCTCCCATTTTGTAGGTTGCCTGTTCA
CTCTGATGGTAGTTTCATTTGCTGTGCAGAAGCTCTTTAGTTTAATTAGATCCCAT
TTGTCAATTTTGGCTTTTGTGTCATTTGCTTTTGGTGTTTGTAGACTTGAAGTCCTT
GCCCATGCCTATGTCCTGAATGGTAATGCCTAGGTTTTCTTCTAGGGTTTTGATGG
30 TTTTAGGTCTAACGTTTCAGTCTTTAATCCATCTTTTAAAAGTCTCTTCACAGTAC
ATGAGTAGTAGTGACACCAATAATGTCAGAGCAGGGAAGTCCCAGGTTCTGCCCAT
CCACAAAAACAACAAATAAGCTGGCAAAAACCTTAAGAATCAACTTTTGCAGATCT
CTGAAATCTAGTCAAACTTAAACAGAGGAAAGATTAATAAAGACNGGCTGCCTGA
GATAACACTAACACACAC

35

Sequence ID 406

CATCAAATAAATAAATAAATAAATTTTAAAAGTCACAGCATTGAATTTTTAAATGT

- 145 -

TTGGGATGATAAAGCACCTGCTTATCATGAAGCTANAGAAATTCAATGACACGTTT
GCCAGGGTCTTTGCTAGTGATGTTGGAACAAGTCTGTAATGCTGATGAAACATCAC
TGTTCTGGGCATTATTGCCCCAGAAAGACACTGACTGCAGCTGATGAAACAGCCCTT
CCAAGAATTAAGGATGCCAAAGACCAAATAAAGTGTGCTGAGATATACTTACGCAGC
5 AGGCATGCATAAGTGTAAACTTGCTGTTATAAGCAAAAGCTTGCGTTCTCACTGTT
TTCAAGGAGTGAATTTTACATACCAATCCATTATTATGCTAATAAAAAGGCATGGATC
ACCAGGGACATCTTTTCAGATTGGTTTCACAAACATTTTGTACCAGCAGCTTGTGC
T TACTGCAGGGAAGCTGACTGGATGATGACTGCAAGATTTTGTATATCTTAACAA
CTGTTGTGCTCATCCTCCAGCTGAAATTTCTCATCAAAAATAATGTTTATGGCTCAC
10 ACCTGTAATCTCAACACTTTGGGAGGATTGCCTGACCCAGGAGTTCAAGCCCACCC
TGGGCAACACAGCAAGACCCAACCTNTC

Sequence ID 407

TTTTAAAAATCATAAAACGTTTCTTACAAAAGAGCATTACATTTNTGCACACTGCTC
15 TGAACAGATGCCAGGGACATGTGGACTATTGTTACTTTTCCCTCCCTGTCCCACCCC
CCAAATGTTACAGTGACCACAAAGCAAGGTGTTTACAATAATTACATGGGGGGAAT
TTTTTAAACCACCAACAATAACGAAAAATAAAATCCACTCACTCTGCTGCTGTTTC
AAAATTTCAATGTTAGTTTTTGCACGCCCTTCCCCCCCCCAACCCTGTTTGTAAAG
AACTAAAACATTACATCTGGTGAACAGCAAAGATTTCACTACACCTCAAATGCAGA
20 ACACCTATGAAGCAGAGGAATGTTGGCTTTTTTAAACAGAAGCAGATAAAAAAAAAA
GATGCAGGACTCCTTCAGTTCTTCACTAGTCTTAGAAAACTTTCCAGAATACTGC
TTCACACTATAAAAAAGAAAAATATCTTGCATTAGAATCCTTCAACATCTGCATA
CTGCTTCACACTGTTTCGTTTCTAGGAGCACTTTGTACAGGACACTTCTGCTTATA
TTTCTTTAATCAGAACTTAGTTGGATGGGCCGGGCATGGTGGCTCACGCCTGTAAT
25 CCCAGCACTTTGGGAGGCCGAGG-GGGTGGATCACC

Sequence ID 408

CCATCTCCAAATTTAGTATTTCATTCTGTTTAGCATATTATCAGTTGCCATCTATTT
GTTTTAACTGATTACTTGAATCTGATTAAACATCACAGAAATGGGCTTTGATAAGA
30 ACAATATTGAATAAGAAATTTTAAATAACAAAACAGCTTATAGAAAAATTTCAGCAT
AACTTTTCCATCACCTTCACCACCTTGCCTTTTATTATCCTGTCCTGTATCACTG
CTTTCTGTTAGCAGTGTTGTGTGAGTTAGGATTTGGGCAGGAAAGCAAAAGCAACC
ACCCGTCATTTTCCCAGAATGAAGGGTTTGACGTAGGATGTAGACTTTGTATAGTA
GTTGGGAGAGCTGTGGGAGTGAAGGTCAGGGATGTCACCTACAGAAGTCAGGGAAT
35 CTGCCACCAGAGATCCTGCATCAGAAACAGCCAACAGCGTGCTTCTGAAGAACTAG
TGGGGAAGTGGCTATAATTCTTAGGAATCCCAGCAAGTCCGCACCACTGTCTCAGT
CTACAGCAGTGGAGAAAGGGGTTTCCAGGAGCTCTCTGGAAAGTTCTTGCCACAC

- 146 -

TTTGCAACAATCTTCAGAGGATAATGGGCTTCTCTTCCAGCTTCCACACCCAACAA
GAGTGCCTTTCATCGGCCAACTCTAACCTGGAACCCCTATGGCAGAGGGGATTTAGG
AGACAGTTTGTNATGTCTGTGGAATGCAAATGAANANGTANCAATGCTTANTTGAC
AGCGGNCATACACAAATNTNGAAA

5

Sequence ID 409
GATCCGTNGACT

Sequence ID 410

10 CTCTTCCCAGCCCCTGAGCCCAGCCCCTTCCCAAGTGGTGCCAGACAAAAAACTAC
ATGGCCCTTTCGTGTCTTGGGGGTGGAAGGGAGGGATGAATTGGGGTGATAGAAC
CCTGGTGAATTCAGAGTAATCTTCTTTAGAAAACCTGGTGTTTTCTAAAGAAACAG
GATAGGAGTTTAGAGAAGGCACCAAAGCTTTCACCTTGGTTTGGCACCAGTTTCTA
ACCATCTGTTTTTTCTACCCTAGCTATCTTTTATTGGTAAAATATAAATGTATAAT
15 TATGTTTGTAGAGCTTTACCAAGGAGTTTCCCTCCTTTTTTGTGTTGATTAGCA
AATTTTTGATTCTCCATTTTCCAAAAGTAAGAGACTCCAGCATGGCCTTCTGTTTG
CCCCGCAGTAAAGTAACTTCCATATAAAATGGTATTTGAAAGTGAGAGTTCATGAC
AACAGACCGTTTTCCATTTTCATCTGTATTTTATCTCCGTGACTCCACTTGTGGGTT
T

20

Sequence ID - 411

nt: 505

TGGAGCTGAAAAATTCCTATTACCTAGGGGCATCACAACGCATTGCATTTGCCCCG
TGTTTGGGATGATGCTGGTGTAAACCTACTATGCTGCCAGTCATGTAAAGTATAG
CACACACAATTAGTAGGTAATGCTTGCAAATAATAATGAAAGACTCTGCTACTGGT
25 TTATGTATTTACTATGCTATACTTTTTGTCACTTACTTTAGAGTGTACTCCTACTTT
TTTTTTTTTTTTTTTTTGAGATGGAGTTTCACTCTTGTCCTGTAGGCTGGAGCGAAN
TGGCGCGATCTCGGCTTACTGCAACCTCCACCTCCTGGGTTCAAGCGATTCTCCTG
CCTCANCTTCCCAGAGTAGCTGAGATTACAGGCATGCACCGCCACGCACGGGTAAT
TTTGTATTTTTGGTAGAGACAGGGTTTCACCATGTTGGCCAGGCTGGTCACCAACT
30 CCTGACCTCAGGTGACCCGCCTCCTCACCTCCAGAGTGTTGGGATTACAGGNGTGA
G

Sequence ID 412

ATAAAAATTAGCTGGGGGTGATGGGCCCTGTACCCAGCTACTCGGGAGGTGAGGT
35 AGGAGAATCACTTGAACCCGGGAGATGGAGGTTGCAGTGAGCCAAGATCGTGCCAC
TGCACTCCAGCCTGTGTGACAGAACAAGACTCTGTCTCAAAAAAAAAATAATAATA
TAATAATAATAAAAAGGAATAACATAGCTAGGAATAAATTTAATCAAAGAGGTGAA

- 147 -

AGACTTATACACTTAAACTACAAAAAAAATCACTGAAGGAATTATAGACCCAA
ATAAAAATAAATAAAAAGACATTCTGTGTTTTAGGGAAAGAAGACTTAATATTGTT
AAGATGTCAATACTACCCAAAGTGATCTACAGATTCAACATAATCCCTATCAAAAT
TCCAACAGCCTACTTTGTAGAAATGGAAAAGCCAATTTTCAAATTCAGATGGAATT
5 GCGAGGGGTTCTGAATAACAAAACACAATCTTGGGGAAAAAAAACAAAAACAAAG
TCAAAGAACTCACACTTCTCTATTTATAATTTACTACAAAGTTATAGNATCAAAGT
CGACGCGCCGCGATCCGGGC

Sequence ID 413

10 CACAGTACTCCATTTTGGGGTCCAACTGTAATGCTCAAAATAAATAAATGCTTACA
CGAAAATTATTTATTGAGAATATTCATATAAAAATTACCTAAAGCAAAGTAAAAAA
AGTAAAATCAAGGTGGTATATTTGAAGTGAATGGTGATTGGAAATTTTGTAGCTGTA
ACAAAAAGAAAGAAACAACCTTTTTTTAAAGCCTCATTCTCTTTTCTTTCAAATG
TACCTTATTTCCACACACTCTTGGGCTGACCTTTATTTTATCAATAAGCTCAATAT
15 TACTTTGTTTTAAATAAGATGCTTCAGCAAAAGTCATTCTCTCTTTAACCATATAA
TTTAAAAACTCCTCTTCACGATTGATAGCAAAATCAGAAACGTTAGGGCACCAGTG
AGTTGAAAAAACTGGTCTTAAGTTGGAAAACTATTATTAATAATATTATCCTATC
CATCCATATCTATTGAAATTGTCAGGTCCATAATTTCATTTTAATTAATTATAGGA
AAGAAGAAAAGATAATACCCATTTGTTCTAT

20

Sequence ID 414

CTCAGACTCTTTCTGCCCTAATGGCCATTACTATCCAGTCTGTATTGCTACAAGGG
ACCCACTGGTACCCCTTTTAGATTCTATCAAAAGGAACAGGGTTTTCTTAGAGGCA
GGCAGCCTGGTGGTATGGCACAGCAGAAGCTTACTGCTAATGAAATGGGAACCTCC
25 CCCTCCCTTGTGGTTTCAGCACAGAACCTGAATGCCAGGAAAAATTCCTGGGCCAA
GAAGCTAAAGCTAAAGAAACCTTCCTTTTTTTCAACGTTTTTTTTCTTTCAAACCTG
TAGGGTCACTTTTGATTGAGGCAAAGGGTCTACTGTAAGTGGAAAAGACTCACT
CCCCTAACATAAGTTTTCACTGTGGTGGGATGGTGCCGCCCCGATATGCTTGATATG
CTTTTCCTTCCACATGTTAAGCTAGGAAACCTAACAGGATGTCAGCAGGGCAGTTA
30 ACTCTGGACTIONAGCCCTCAAGGGCATGTGGCANAACTCATGGCATNCAAGACC
A

Sequence ID - 415

nt: 596

GTATAATTGATTCTTTTGAACCTAAAGTATAAGACTTCACGATTAGAAAAAATTA
35 TCCAAAGACTAATGTAATTAAGTGAGGAAAAGGTGCTGGAGGAACTGGATAACCAC
ATGGAAATGTATGAACCATGACCTCTATGTCACATACTATATATAAACTTAATTT
GAGGTGTATCACAGAGCTAACTGTGGGGGCTAAAACGTTGAAGCCTTTGGATGGCC

- 148 -

GCACAAGAGATGTCTGCATTCATAACCTTGGGGAGGGTATGAACATTTCTTGGTAA
CATGCAAAAAGCACTAACTGTAAAAGAGAACAGTTGGTCAGTTGAATTTTCATGAAA
CATTTGTAAACTTCTGCTAAACAACCTGACACCATTGAAGAATGTGGAAAAAGGCTGGG
CACAGTGGCTCATGCCCTATAATCCCAGCATTTTGGGAGGCCGGGGCGGGAGAATCA
5 CTTGAGGCCAGGAGTTTGAAACCAGCCTGGGCAACATGGCAAGACCCCGACTCTAC
AAAAATATTTTAAAAATTAGTTGGGTGTGGTGATGCACTCCTGTAGTCCTAGCTG
CCAGGANGCTAAGGNGGAAGGATCACTTAACCCCTGG

Sequence ID 416

10 CTGGTGGCGGCGGTCTGTGCGGACGCAAACATGCAGATCTTTGTGAAGACCCTCACT
GGCAAAACCATCACCCCTTGAGGTCGAGCCCAGTGACACCATTGAGAATGTCAAAGC
CAAATTCAGACAAGGAGGGTATCCACCTGACCAGCAGCGTCTGATATTTGCCG
GCAAACAGCTGGAGGATGGCCGCACTCTCTCAGACTACAACATCCAGAAAGAGTCC
ACCCTGCACCTGGTGTGCGCCTGCGAGGTGGCATTTATTGAGCCTTCTCTCCGCCA
15 GCTTGCCCAGAAATACAACCTGCGACAAGATGATCTGCCGCAAGTGCTATGCTCGCC
TTCACCCTCGTGCTGTCAACTGCCGCAAGAAGAAGTGTGGTCACACCAACAACCTG
CGTCCCAAGAAGAAGGTCAAATAAGGTTGTTCTTTCCTTGAAGGGCAGCCTCCTGC
CCAGGCCCCCGTGGCCCTGGAGCCTCAATAAAGTGTCCCTTTCATTGACTGGAGCAG

20 Sequence ID 417

GCAGGGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCCGCAGATAAGTTT
TTTTCTCTTTGAAAGATAGAGATTAATACAACCTACTTAAAAAATATAGTCAATAGG
TTACTAAGATATTGCTTAGCGTTAAGTTTTTAACGTAATTTTAATAGCTTAAGATT
TTAAGAGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAAAAGG
25 TTTCTAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAATGTAT
TTAAAAGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACCAATAGAAGG
GCAATGCTTTTAGATTAAAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTTTAAA
AGTTGTAGGTGATTAAAAATAATTTGAAGGCGATCTTTTAAAAAGAGATTAAACCGA
AGGTGATTAAAGACCTTGAAATCCATGACGCAGGGAGAATTGCGTCATTTAAAGC
30 CTAGTTAACGCATTTACTAAACGCAGACCAAAATGGAAAGATTAATTGGGAGTGGT
AGGA

Sequence ID 418

35 CCCGGAATCGCGGCCGCGTCGACGGGAGGTGATAGCATTGCTTTCGTGTAAATTAT
GTAATGCAAAATTTTTTAAATCTTCGCCTTAATACTTTTTTATTTTGTTTTATTTT
GAATGATGAGCCTTCGTGCCCCCCTTCCCCCTTTTTTGTCCCCCAACTTGAGATG
TATGAAGGCTTTTGGTCTCCCTGGGAGTGGGTGGAGGCAGCCAGGGCTTACCTGTA

- 149 -

CACTGACTTGAGACCAGTTGAATAAAAGTGCACACCTTATAAAAA

Sequence ID 419

5 CCCGGAATCGCGGCCGCGTCGACGGGAGGTGATAGCATTGCTTTCGTGTAAATTAT
GTAATGCAAAATTTTTTTAATCTTCGCCTTAATACTTTTTTATTTTGTTTTATTTT
GAATGATGAGCCTTCGTGCCCCCCTTCCCCCTTTTTTGTCCCCCAACTTGAGATG
TATGAAGGCTTTTGGTCTCCCTGGGAGTGGGTGGAGGCAGCCAGGGCTTACCTGTA
CACTGACTTGAGACCAGTTGAATAAAAGTGCACACCTTATAAAA

10 Sequence ID 420

CTTCATTTGAAATGGTTGAATCTGCTGTGTAATAAAGTGGTTCAACCATGATTAGG
AACTGAAATTTAGTAGAAGAGGGGAAAAGGAGTTAATGTAACAAATTATTTTAGCTA
CAAACCCCGGTAATAGAGCACTTGGGGGATGGGATGGGGTGGGTGGTGAGACAAT
CAGAATGGTAAATTGATTAAATGCTCCTAACCTGTAATTTTGTGCATAGAGCACC
15 CTATGCTGTGGAAATAACTGTTCTTAGATTTTCATTGTAAGTGGACTGTTTCAGGTTG
CCCAGAGGGGAAAGAACATTCTAATTCTAATAAAATAAACTTTTTATTTTGTTTA

Sequence ID 421

TGTCATTGAATCTGCTTGTTACTTAAATGCTAAACTCAATTCTGTAATTCAATAGG
20 TGCACCTCTCTGAGAAACATAAGAGACAATGAGGAAAAGGATTCAGCATTCCGTGG
AATTTGTACCATGATCAGTGTGAATCCCAGTGGCGTAATCCAAGTAAGATGTTTAC
AAAGATTTGTTTTTTAATGTCTAATTAATAAAATTTTAAAGGAAGAAACATTCTAAT
ACTTTAATTATAAAAAGTTAACTATTTTCAAAGGTATCAAAATACAGTTAAACCTT
TAAAATGTATATTTCTTAATATCTTGAAATTGTAATGCCTTTTTTTTTTCTTAAAT
25 TTTTTTTGTGCATGAAATGAGATAGTAACAGCAGATTGGGACAACAAGGTTATATTC
TTGTCTTGAATCAGGCCATGGCTTCTTTCATCCAAATTTTCAAGCTCATTATTTTA
CTTTGTCCCTGCCTCCCATCCCTGGATATCAGTTTGTGGATATCTACAGTTAATAG
AGTGACCAAATAGTAGGAATACTGTCTCTCTATTCTGAATAAAATCTTTGAATCAG
ATTTAGAAATAATGAATAAAATACAAATCAGCCATTGAAATTGCTCTAATTTTGAG
30 AGCTTATGATTTATTCATCTTTGGTTTCCAAGTTCAAGTTATATGTAGACATTTTA
ATT

Sequence ID 422

35 GCTTCCTAGGTGAGGTCACGAGGAAACCTGCTGGCCAAGTGACCTGGCAGGGTGTG
GCCAGTGTGGCCAGGGCCGCCGAGCCTGCTTTCCTTCCCTGCAGCAGGAACCCTTC
TGGGGCTGTGATCCTGCGATGGTGCCTGGGTGGGAGTGGGGGTGGGGGGCGGGATG
GTCTCCCTACCTGCCAGCTTCTTGGTTTGAGGTGAGGACAGCCCCGAAGCTCANA

- 150 -

CTTGGCTCCTGTCCATGTACTTGGGGCCATGAGCTCTGCAGGGACCTTGGAAAGAN
AGAGACGGGTGGTGTANGGCANGGGAAGGCATTGTCTTCAAACAGGAAAAAGCTGA
NAATGGAAACAGGCGAAACTTACCAAGTGTAACATCACCTGGAAGTGAAGGAGGGT
GGGAAGGTTTTAATTATTTTAAAAATAGAGATGGGGTCTCACTATGTTGCCCAGGC
5 TGGTCTCAAACACTACTGGGCTCAAGTGAACCTCCTTCT

Sequence ID - 423

nt: 387

TGTTTCTCNAGGGCGAGAGGCTGTCTTANAGCACCATTCTCTGGCCCTNGTCCCAT
GAGAAGGAACCGCACTCAGGAGCCACACTCTCCCACTNCCCTTGCCCANAAAGACTC
10 ACAGAGGGCACGGAGCTGGCTGTGGTGAGAGGAGGTCCANCAAATTCCTGTCTGCA
NAAGGGTTCTGAACACCACCGCCTGGCAGCGTGCTGGAGGAGGGATTCTCTTTTC
CTCACAGCAATTCTGACCAGAAACCTGTCAAATCAGGAATGGCTAAAAATAAGACCA
GGGTATGAATGACCATCAGCCACAGTAAACCAAGGCACAGCTCTCCTGAGCCAC
CCAAGCTGCTGTGGCCCAGACTGGTGACATCACCTCAGGGCAAAAAAAAAA

15

Sequence ID - 424

nt: 420

CGCAGAATGGCTCCCGCAAAGAAGGGTGGCGAGAAGAAAAAGGGCCGTTCTGCCAT
CAACGAAGTGGTAAACCGAGAATACACCATCAACATTCACAAGCGCATCCATGGAG
TGGGCTTCAAGAAGCGTGACCTCGGGCACTCAAAGAGATTCGGAATTTGCCATG
20 AAGGAGATGGGAACTCCAGATGTGCGCATTGACACCAGGCTCAACAAAGCTGTCTG
GGCCAAAGGAATAAGGAATGTGCCATACCGAATCCGTGTGCGGCTGTCCAGAAAC
GTAATGAGGATGAAGATTCACCAAATAAGCTATATACTTTGGTTACCTATGTACCT
GTTACCACTTTCAAAAATCTACAGACAGTCAATGTGGATGAGAACTAATCGCTGAT
CGTCAGATCAAATAAAGTTATAAAATTG

25

Sequence ID 425

GGAAACTGATGCCAGTCAGAAACTCAGATCAAATGAAGGGGTGAAGAGAACCAGAA
TTGATCTCTCTGTAGGAGAATATAAATGACTTTTTTAAAGTACATATTTTCTGTGA
AAGACAGTTTTTTGTTTAATGCAAAAATGTTAACAATGTTTATATCATGTAGAAGT
30 AAAAGATCGTGAAACAGCACAGAGAACAGTAGTAAGACAGATTGAATTGCACTGTT
GTAAGATGATGAACTTACAATATTAAGTGAAGGTAGACTGTGATAGATTAAGGATA
TATATTGTAATCCCTAGAGCAATTGTCAAAGTGGTACAGGTAAAAAGCCAATAGAG
GTGATAAAATGGAATACTAAAAAATATCAGATGAATAATAAGAAGACAGGAAATG
AGGAACAGTGGAACAGAATGAATAAAAAACAAGACCATTAACTTAATCATTAATAA
35 TTACTTTAAATGGGTAAACATTATGGTTATAAGGCAGAGATTTTCAGACTAGATA
AAAGAGCAAGCTCCACTATATACTGTCTACAAGAGATATACTTTAAAGTGTATATT
ATATTTAAATATAAAGATTTGGAATAAATAAACCTAAGAATAAGCTTACTAGGGAA

- 151 -

GTGAAAGATCTGTACAACAAGAATTACAAAACACTGCTGAACGAAATCATAGGTGA
CCA

Sequence ID 426

5 GTCCCGGAATCGCGGCCGCGTCGACGTTTCCTCAAAATTTATCTTCCTGTTAATGT
CAGGCATGTATCTCCTTAGCTTGCCACAAATAACTATATATACCACAGACCTTCCT
TTGTAGGGCTAACAGTGTTGCATTGTAAGTGGAGGCCTCATAGATACCTGGCCTTT
TCCTACCTTATTCCAAAGATGGTTGCATCTTATAAATAATGTCATTCTTCAGCAAA
TGGTATGGAAATGAGATTGTAATGTCATTATTTCTCTTTAAATAATCAGGACAAC
10 TCATGATACAAAGAGCTCTTCTCTATAAAAGGTGGGACTTTTTTTTTTTAGTAATAG
CAAAAATAAAATTGTACCTCCTTAATCTTCTACAGAAAGATGGATTTCAATTTCAA
CATTAAGAGGTAGTTTTAAGAAGCAGTAGAAGTCAGCCTGGGCAGCATGGTGAAAC
CCCGTCTCTACAAAAAAGTTAGCTGGGCTTAGTAGTTGCAATCCCAGCTACTCTGG
AGGCTGAGGTTGGAGATCATCTGANCTGGGGAGGTCNAGGCTGCAATGATACANT
15 GAGCCCTGATTGTGCCACTCCACCTGGTTGCAGA

Sequence ID 427

TTCCAATCTTCGTGTTCACTTTAAGAACACTCGTGAAACTGCTCAGGCCATCAAGG
GTATGCATATACGAAAAGCCACGAAGTATCTGAAAGATGTCACTTTACAGAAACAG
20 TGTGTACCATTCCGACGTTACAATGGTGGAGTTGGCAGGTGTGCGCAGGCCAAGCA
ATGGGGCTGGACACAAGGTCGGTGGCCCCAAAAGAGTGCTGAATTTTGTCTGCACA
TGCTTAAAAACGCAGAGAGTAATGCTGAACTTAAGGGTTTAGATGTAGATTCTCTG
GTCATTGAGCATATCCAAGTGAACAAAGCACCTAAGATGCGCCGCCGACCTACAG
AGCTCATGGTCGGATTAACCCATACATGAGCTCTCCCTGCCACATTGAGATGATCC
25 TTACGGAAAAGGAACAGATTGTTCTTAAACCAGAAGAGGAGTTGCCCAGAAGAAA
AAGATATCCCAGAAGAACTGAAGAAACAAAACTTATGGCACGGGAGTAAATTCA
GCATTAAAATAAATGTAATTAAAAGG

30 Sequence ID 428

TGCAGGATCCGTCGACTCTAGATAACATGGCTAGAAAAGAGAATGAAAAAGTTGGA
ATTTTTTAATTGCCATGGTATGGGGGGTAATCAGGTTTTCTCTTATACTGCCAACAA
AGAAATTAGAACAGATGACCTTTGCTTGGATGTTTCCAACTTAATGGCCCAGTTA
CAATGCTCAAATGCCACCACCTAAAAGGCAACCAACTCTGGGAGTATGACCCAGTG
35 AAATTAACCCTGCAGCATGTGAACAGTAATCAGTGCCTGGATAAAGCCACAGAAGA
GGATAGCCAGGTGCCCAGCATTAGAGACTGCAATGGAAGTCGGTCCCAGCAGTGGC
TTCTTCGAAACGTCACCCTGCCAGAAATATTCTGAGACCAAATTT

- 152 -

Sequence ID - 429

nt: 535

CACAGTACTCCATTTTGGGGTCCAACTGTAATGCTCAAAATAATAAATGCTTACA
CGAAAATTATTTATTGAGAATATTCATATAAAAATTACCTAAAGCAAAGTAAAAAA
AGTAAAATCAAGGTGGTATATTTGAAAGTGAATGGTGATTGGAAATTTTGTAGCTGTA
5 ACAAAAAGAAAGAAAACAACCTTTTTTTTAAAGCCTCATTCTCTTTTCTTTCAAAATG
TACCTTATTCCACACACTCTTGGGCTGACCTTTATTTTATCAATAAGCTCAATAT
TACTTTGTTTTAAAATAAGATGCTTCAGCAAAAGTCATTCTCTCTTTAACCATATAA
TTTAAAAACTCCTCTTCACGATTGATAGCAAAATCAGAAACGTTAGGGCACCAGTG
AGTTGAAAAAACTGGTCTTAAGTTGGAAAACTATTATTAATAATATTATCCTATC
10 CATCCATATCTATTGAAATTGTCAGGTCCATAATTTTCATTTTAATTAATTATAGGA
AAGAAGAAAAGATAATACCCATTTGTTCTAT

Sequence ID 430

CAGGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCGCAGATAAGTTTTT
15 TTCTCTTTGAAAGATAGAGATTAATACAACCTACTTAAAAATATAGTCAATAGGTT
ACTAAGATATTGCTTAGCGTTAAGTTTTTAACGTAATTTTAATAGCTTAAGATTTT
AAGAGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAAAAGGTT
TCTAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAATGTATTT
AAAAGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACCAATAGAAGGGC
20 AATGCTTTTATAGATTAAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTTTAAAAG
TTGTAGGTGATTAAAATAATTTGAAGGCGATCTTTTAAAAGAGATTAAACCGAAG
GTGATTAAAAGACCTTGAAATCCATGACGCAGGGAGAATTGCGTCATTTAAAGCCT
AGTTAACGCATTTACTAAACGCAGACGAAAATGGAAAGATTAATTGGGAGTGGTAG
GATGAAACAATTTGGAGAAGATAGAAGTTTGAAGTGGAAAACCTGGAAGACAGAAGT
25 ACC

Sequence ID 431

CGCTGGGTGCCTGCAGCGCCTCCCTTGTCTCATATGGTGTGTCCAGCACTCTATTG
TTGTAAACTGTTGNTTTGNCTGACCTAAATTNTCTTTACTAAACANATTTAATAGT
30 TNAAAAAANANCA

Sequence ID 432

TTTTAAAGTCATCTCTATAGGAAGGTGCTGGGCAGGGATCCCAGAGAAAGAAAGGG
TCCAAGACTCCATTAACCTGCCCTGGATGAAGGGCACTGCTACAGCAGCTAGTACCA
35 GAGACTCTCCTATCTCACGGTTGAGGCAGACCCAGGATAGAATAGAGAATAAAAGG
AATGCTTATAGGAAACAATTTTGTATGGAATGCTAGATGGCCAAGCCTCAGCCTTT
GGTCCAGTGCAACCCTTGCCCTCGCTTGTCAACAGTGAAAAATTAGTTTGGTTAGAA

- 153 -

GAACCATCTGGAAACACACCAGCTTCTGCTACCTTCATGCTCATTGTTAAAAAAG
ATTAACCAGTGTGAACATTCTGATCTGTTAATTCCAGGGACTGTTTTCTTTCCAAT
GGACTGTTTGTGTTAGTAATAACCCCCAAAAGCTCAAAGCTAAAAATGCATCATCAG
TCCTAGTCGGCAGTTCCTTAAGAATGGACTGGCGGCGTGGTTGAGCTGATATGGAA
5 AAGCTGCACCTTCCTGCAGAAGATCAACTGACCTGCTATCCCACCCCCAAATTTCAA
CCTGAGGTATATTTCAATGAAGGCAGGTAGCTGTGCTTCTCAGAGCA

Sequence ID 433

TCCCGGAATCGCGGCCGCTCGACCCGCCGCCGAGGATTGAGCAGCCTCCCCCTTG
10 AGCCCCCTCGCTTCCCGACGTTCCGTTCCCCCTGCCCGCCTTCTCCCGCCACCGC
CGCCGCCGCCTTCCGCAGGCCGTTTCCACCGAGGAAAAGGAATCGTATCGTATGTC
CGCTATCCAGAACCTCCACTCTTTCGACCCCTTTGCTGATGCAAGTAAGGGTGATG
ACCTGCTTCCTGCTGGCACTGAGGATTATATCCATATAAGAATTCAACAGAGAAAC
GGCAGGAAGACCCTTACTACTGTCCAAGGGATCGCTGATGATTACGATAAAAAGAA
15 ACTAGTGAAGGCGTTTAAGAAAAAGTTTGCCTGCAATGGTACTGTAATTGAGCATC
CGGAATATGGAGAAGTAATTCAGCTACAGGGTGACCAACGCAAGAACATATGCCAG
TTCCTCGTAGAGATTGGACTGGCTAAGGACGATCAGCTGAAGGTTTCATGGGTTTTA
AGTGCTTGTGGCTCACTGAAGCTTAAGTGAGGATTTCTTGCAATGAGTAGAATTT
CCCTTCCTCCCTTGTCACAGGTTTAAAAACCTCACAGCTTGTATAATGTAACCATT
20 TGGGGTCCGCTTTTAACTTGGACTIONCTTCATGCAATAAACTGAAAAG
ACCATGCTGCTANTC

Sequence ID 434

TTCGGACGCAAGAAGACAGCGACAGCTGTGGCGCACTGCAAACGCGGCAATGGTCT
25 CATCAAGGTGAACGGGCGGCCCTGGAGATGATTGAGCCGCGCACGCTACAGTACA
AGCTGCTGGAGCCAGTTCCTGCTTCTCGGCAAGGAGCGATTTGCTGGTGTAGACATC
CGTGTCGCTGTAAAGGGTGGTGGTACGCTGGCCCANATTTATGCTATCCGTCAGTC
CATCTCCAAAGCCCTGGTGGCCTATTACCANAAATATGTGGATGAGGCTTCCAAGA
AGGAGATCAAAGACATCCTCATCCAGTATGACCGGACCCTGCTGGTAGCTGACCCT
30 CGTCGCTGCGAGTCCAAAAAGTTTGGAGGCCCTGGTGCCCGCGCTCGCTACCAGAA
ATCCTACCGATAAGCCCATCGTGACTCAAACTCACTTGTATAATAAACAGTTTTT
GAGGGATTTTAAAA

Sequence ID 435

CTGCAATGTGCAATAGTTGCACCACTGCACTCCAGCCTGGGTGACAGAGTGAGAAC
35 CTATCTCTTAAAAAAAAAAAAAAAAAAGGAAGAAGAGACATGAGAGGGCCCAAGT
CACTTGCTCACTCACTTTCCGTGTACATGTACCAAGAAAAGGCCATGTGGGAAAGA

- 154 -

GCAAGAAGGCAGCCGCTTCAAGACAGGAAGAGAGCCCTCACCAGAACTGAGCCA
GAACCTTGAATTCCAGCCTCCANAACGTGTGAGAAAAGAATTTTCTGTTGTTTCAG
TCCCCCACACTATGGCATTGTTACGGCAGCCTGAGCTAATACTCCTACTTTGTC
CTGCATTTACTTGGTCTTCCAGTTAGTTTTTTAGACTTTGGGAATCAGAGCAGTCA
5 GTTGTGAGATTTTAGCTTACAGTTGTCCTACCTGTGCAACTGAAATTTCTTCCATT
TTAAACCAGAGCAGAGTTTTAGAGTCAAAGAAACCAGATCTTTTAGTGCAGAAGC
TTTCCACTGTATTANAAGTGAGGAAGTTGGT

Sequence ID 436

10 AAAAAAAGCTCCAGAGAAGTTTATAGAAAGAGATGACATGTAAACCCTGCTGAAAAA
TAGTTTCATTTGTTAGAATATAATTGTCTTCCACTAAAAAAGAAAAAAGCA
TTTAAGGCTCTAAGATCTCTTGAAGTACCCTTTTCTGAATCCCAGAGTTTTTAT
GTGCATTATTTTTATGCGTTTGTAGTTTGATATGTTGTATTTATAAGTAGTTTTAG
CTTTCCATTATGAATTCTTCTTTGACCCATGAGTTATTTAGGTAAGTGTAAAAA
15 TTTACAATAGTTTATATATGCAAATATTATGTTGTTAGAGTTGGTTTTTCATGTCAT
TTTTACATATACAGGGGCAGTTTCCCCAACTAAATTGTATATTCTTAAAGCAGCA
CTCTTAAATTTTATTTCTGTGTCAATTTCTTGNCTGTGTTTCCTGGCATGGAATAC
ATGGCATAAAATTTGTTATGTAATTAAATGAAATATTATTATACTTTCTATTTTTT
AGAAAAAA

20

Sequence ID - 438

nt: 577

GTCGACAGGGATGACATAACTATTAGTGGCAGGTTAGTTGTTGGTCACTTTCAACT
CTGGGTTCAAGCGATTCTCCTACCTCAGCCTCCCGAGTAGCTGGGATTACAGGCAT
25 GCACCGCCACACCTAATTTTCTATTCTTAGTAGAGACGGGGTTTCTCCCTGTTGGT
CAGGCTGGTCTCGAACTCCCGACCTCAGGTGATCTGCCTGCCTCAGTCTCCCAAAG
TCCTGGAACCACAGACATGAGCCACCACGCCTGGCCCCCTTTTAAAATATTTCTGCT
CATTGATGATGCACCCAGTCACCCAAGTGCTCTGATGGAGATGTATAAGGAGATGA
ATGCTGTTTTTCATGGCTGCTAATAACAATTATTCTGCAACCCCCAAATCAAGAA
30 GTAATTTTGACTTTCAAGTCTTATTATTTAAGAAATATATTTTGCAAGACTATAGC
TGCCATAGACCGTGATTCTCTGATGGATCAGACAACTAAAATGAAAACCTCCTG
CAACGTATTCATCATTCTAGATCCCTGAGGAATCGCCACACTGACTTNCACAATGG
GTGAACTGGGTTACAGT

35 Sequence ID - 441

nt: 552

AAACAAAATTATTCTCTGAGAGGGAAAGGACATTTGAGGGAAACATCAAATTTCCC
CATAAATAAATGAATGGAGTTTGCAGGAAGGTGAGGGTGAGCAGAGATGTGTGTGG

ACATCTCTGACCATCCATCGCTGTATTCAAATGGATTGTTTTATTCCATTCTGGTC
TCAGGCATGACCACGTCCAGTGAAGACATTTGAGGCAGCACATCTCAGGACCCAGG
CAATAGACTGGCCCCAACTCAGGCTGGACTAAGGTGTGATTAATTCTTTGTTTTTT
GTGTGGAACAGCTCACCTTGTGAGACAGCCTCAGGGCATCTCTGAGACACAGGGGC
5 AGAAAATGACATTCATCTTTTGAGTCCTCATCCATGGAGTGCTGTGTTTGGGGGGC
TGCATCTGCTGAAGCGAGAACCCCATTCCTGCCACCCACCAGGATGCCCATTCCTCC
AGGACTTCTCCAACCTTACTATTAGACTAAACCAGAACAAAGCAACAAACTGTATTTA
TGCAAGCAAAATTGATGAGAAAATTATATTCAAATAAAGCAAAAATTA

Sequence ID 446

CGGACTCCTGTGCTAATTGTCAGCTTACATATCATTGTATAGAGACTGTTTATTCT
25 GTACCAAACTGATTTCAAAGTACTACATNGAAAATAAACCGGTGACTGTTTTTCT
TCATAAAGTTCTGCGTTTGGCATCTTCACTCTTCCAAAATGTATCTGTACATCAN
AAATGTCACTATTCCAAGTGTCCTTTTAGTGTGGCTTTAGTATGGCTTCCTTTTAA
TATTGNACATACATTGNATCTTTGTTTTATGGNAATAAGTAATAAAAATGTAGACT
TCATATTTTGTACAAAATGTCCTATGTACAGAATAAAAAAGTTCATAGAAACAGCC
30 NANAA

Sequence ID 447

AGGCCGAGGCAGGCAGATCNCNTGAGGTCAAGAGTTTGAGACCAGCNTAGCTAACATGGTGA

TGAAACCCCATCTCTACAAAAATATA-AAAATTAGCCTGG-GTGGTGATGGGCACCTGT

AACCCAGCTACTCGGGAGGCTGAGGTAGGAGAATCACTTGAACCCGGGAGATGGAGGT

TGCAGTGAGCCAAGATCGTGCCACTGCACTCCAGCCTGTGTGACAGACAAGACTCTGT

CTCAAAAAAAAAATAATAATAATAATAATAAAAAAGGAATAAC

- 156 -

ATAGCTAGGAATAAATTTAATCAAAGAGGTGAAAGACTTATACACTTAAAACTACA
AAAAAAAAAATCACTGAAGGAATTATAGACCCAAATAAAAATAAAAGACATT
CTGTGTTTTAGGGAAAGAAGACTTAATATTGTTAAGATGTCAATACTACCCAAAGT
GATCTACAGATTCAACATAATCCCTATCAAAATTCCAACAGCCTACTTTGTAGAAA
5 TGGAAAAGCCAATTTTCAAATTCAGATGGAATTGCGAGGGGTTNTGAATAACAAAA
CACNATCTTGGGGAAAAAAAACAAAAACAAAGTCAAAGAACTCACACTTCTNTAT
TTATAAATTTACTACAAAGTTATAGTAATCNAA

Sequence ID - 448 nt: 329
10 TACGCACACGAGAACATGCCTCTCGCAAAGGATCTCCTTCATCCCTCTCCAGAAGA
GGAGAAGAGGAAACACAAGAAGAAACGCCTGGTGCAGAGCCCCAATTCCTACTTCA
TGGATGTGAAATGCCCAGGATGCTATAAAATCACCACGGTCTTTAGCCATGCACAA
ACGGTAGTTTTTGTGTGTTGGCTGCTCCACTGTCCTCTGCCAGCCTACAGGAGGAAA
AGCAAGGCTTACAGAAGGATGTTCCCTTCAGGAGGAAGCAGCACTAAAAGCACTCTG
15 AGTCAAGATGAGTGGGAAACCATCTCAATAAACACATTTTGGGTAAAA

Sequence ID 450
GAGCAGTGGCATGATCACACCTTACTGCGGCCTCCAACCCCTGAGCTTAAGTGATT
CTCCCGCATTATCCTCCTGAGTAGCTGAGACTACAGGTGCATGCCACCATACTACTA
20 CTAAATTTGGGTGCGGTGGTGGTGGTGAATTTTTTAATATTTTTGTAGAGACAGGGT
CTCACTGTGATGCCCAGGCTGGTCTTGAACCTCCTGGGCTCAAGCAGTCACCCACCT
CAGCCTCCCAAAGCACTGGGATTACAGGTGTGAGCCACCACACTGGCCAGCTTTGT
TTTGTTTTGATGACTAAGCTGCTCTTGCTAAAAGGGCTTCTCTCTGAACCTCCCTA
CCTTTCTTCTGTTTCCCTGGGCTAGGGCTCCATGTTGGCAGTCCCTACTCCCAATTA
25 ACCTGGGGCTGTCTGGTTAACCTTTATAAGATCTGCAGTCATTGGGAGACCCGGGG
ACCAGGAATATTGTTGTTGAGGGAGCTACCCTGGAAAGTGGATGGGTGGCCAAAGG

Sequence ID 452
TTTGGCTTTGCCTCTAGGCATTAGATGTTATCTTTGGAGGCATCCTTCTATGAGCA
30 TTCATTTTTGGACCAAGCCTGGATTTACAATTCTATTACTGGCCCAGACTTCATTT
CTATCCAATTTCAATCCACTGTGCTATAGTTTACAACATATAATTTGACTTATAAA
TAATTCCTGACTATGGGTTTAAAGACTGAAAATGGATCAATAGAACTTTGAAAAT
GTTAACATCTTGATTGCTTTTCTCAGTGTAGAAATGGACAATGTTTAGCTTAAAAA
CTGCATGTTTTTAATGAGATACGGGGTTGAAAGACTTATTCCTGGAATTTATTGTT
35 CTGGAGAAAGCCTGTTGCTATCTGCCATACCTTGGTTTACTTTGTGCAAAATGAGC
TTCTTTTTTAAGTAATGAGCTCTTTCATGTTTCAGCTTAAATTGCTGTCTTAGACAC
TTCATCAGGGTTCCTGCTCTGCCTCATTCCCCCTTTTGCTCACTTGCAGCCTTTG

- 157 -

ACATAATCCTGGGAGGCAATTGGCATCATACATATTTTGCTTTGTAATCTCCTGCT
TTGATTCTGACTGGGACCCAGC

Sequence ID - 453

nt: 747

5 GGATCTAAGACCAGCCTGGCAGCCACCAGATGGTGATTCTAGTCCTGGCTCAGTCA
GTAATAGGTCACTGACCCAGAGAAATCAATTCAGCCTCCCCAGGTCCTTGGATTT
CTTTCTGTGAAAATGAAAGCATAGGTAGGAATTTCCCATGGAACAGCTAGCAGAGG
AGAAATATTAAAAGTCAGGAGACTCATGCTATAGTTTTTCATACTTCATTACAACAA
TGTTGTTTAGGACAAGTGAGTTAACCTGTTAGCTTCCTCTATATAAAATGGAAAGT
10 CATTAAAAACCTACATAGCAGGGTTCTTGTGAAGATCAAGTGATAATGTAGGAAGC
ATGTACAAATGTCACATTCTGCCGTCACGTAATGGTCCTCACAGCTTGAGGTAGCA
TTTAGCATGTGTGTCATGATTTAGTACAAGGGTTGGCAAACCTGTTGCTCTTGGATTAA
GTCTGGCTCATTTGCCTGTTTTTCAAAGAAAAAAATTGTATATGTGTGTATATATGT
TATATATAGGTACACACACATATGTGCTATATATAGCATATATACACACATAATAT
15 ATAAACATGTACATATATAGCATTATATATATACCGTGTATAATATCTCCAGTCCT
CATGACCAGCCATGCTTGTTCAATTACATTTGCATACTCTATGATTGCTTTCATGC
AACAATGGCAGAGTTGAGTGATTGTTTTGCACAGANACTGTATGGCCCACTAAACC
TAAATATTAATCTCTGCC

20 Sequence ID 454

CTCCTGCCGGGCTCGTGGCGGCTTCTGTCCGCTCCGCGGAGGGAAGCGCCTTCCCC
ACAGGACATCAATGCAAGCTTGAATAAGAAAAACAAATTCTTCCTCCTAAGCCATG
GCATATCAGTTATACAGAAATACTACTTTGGGAAACAGTCTTCAGGAGAGCCTAGA
TGAGCTCATACAGTCTCAACAGATCACCCCCCAACTTGCCCTTCAAGTTCTACTTC
25 AGTTTGATAAGGCTATAAATGCAGCACTGGCTCAGAGGGTCAGGAACAGAGTCAAT
TTCAGGGGCTCTCTAAATACGTACAGATTCTGCGATAATGTGTGGACTTTTGTACT
GAATGATGTTGAATTCAGAGAGGTGACAGAACTTATTAAAGTGGATAAAGTGAAAA
TTGTAGCCTGTGATGGTAAAAATACTGGCTCCAATACTACAGAATGAATAGAAAA
ATATGACTTTTTTACACCATCTTCTGTTATTCAATTGCTTTTGAAGAGAAGCATAGA
30 AGAGACTTTTTTATTTATT

Sequence ID - 458

nt: 682

TGCCACTGAAGATCCTGGTGTGCGCATGGGCCGCGCCCCGCGCTTGTTACCGGT
ATTGTAAGAACAAGCCGTACCCAAAGTCTCGCTTCTGCCGAGGTGTCCCTGATGCC
35 AAGATTTCGCATTTTTTGACCTGGGGCGGAAAAAGGCAAAAGTGGATGAGTTTCCGCT
TTGTGGCCACATGGTGTGAGATGAATATGAGCAGCTGTCCTCTGAAGCCCTGGAGG
CTGCCCGAATTTGTGCCAATAAGTACATGGTAAAAAGTTGTGGCAAAGATGGCTTC

- 158 -

CATATCCGGGTGCGGCTCCACCCCTTCCACGTCATCCGCATCAACAAGATGTTGTC
CTGTGCTGGGGCTGACAGGCTCCAAACAGGCATGCGAGGTGCCTTTGGAAAGCCCC
AGGGCACTGTGGCCAGGGTTCACATTGGCCAAGTTATCATGTCCATCCGCACCAAG
CTGCAGAACAAGGAGCATGTGATTGAGGCCCTGCGCAGGGCCAAGTTCAAGTTTCT
5 GGCCGCAGAAGATCCACATCTCAAAGAAGTGGGGCTTCACCAAGTTCAATGCTGAT
GAATTTGAAGACATGGTGGCTGAAAAGCGGCTCATCCCANATGGCTGTGGGGTCAA
GTACATCCCCAATCGTGGCCCTCTGGACAAGTGGCGGCCCTGCACTCATGAAGGCT
TTCAATGTGC

10 Sequence ID 459

TCCCGGAATCGCGGCCGCGTCGACCTTGTCCTTGAGCGTCAACCTTCTTTCCCTGA
AGTGGCTGGGGTTCCTGTTTCCTTCTTTGATTGACAACTTGTGTTAACCTTCGCAC
ATCTCTGGGCCAATTTTTGCTTGTAAGTCTTTCGGGAGACCCCTGGAATTTAAATC
ATTAGCACCGCGCCCTTCCCCGAAGAGTCTTCGAAGGGTTGCCGCTTTTCGGTGGC
15 GCAGTTCTCGCGAGAAGGTGACTTTCTTCTCGGTATTTCTTGGTTTCCAGAATCC
TTAGCGCGAGGCGGAAAAAATATTTCTCCAGCTTGTGTTGATGCCGCGATTTTGA
CTGAGACTTCTTCCCACGATTTCTGTTTTTGCTTCTCCAAGGAAAATGGCAGCTCC
CGAGCAGCCGCTTGCGATATCAAGGGGATGCACGAGCTCCTCCTCGCTTTCCCCGC
CTCGGGGCGACCGAACCCTTCTGGTCAGGCACCTGCCGGCTGAGCTTACTGCTGAG
20 GAGAAAGAGGACTTGCTGAAGTACTTCGGGGCTCAGTCTGTGCGGGTCTGTCAGA
TAAGGGGCGACTGAAACATACAGCTTTTGCCACATTCCCTAATGAAAAGCAGCTN
TAAAGGCATTGACAACTNCATCAACTGAAACTTTTAGTCATACTTTAATCG

Sequence ID - 460

nt: 536

25 CAGAGATCAAAATAGGCCTTACACAGTGCGACGCGAATTTAAAAGATTACCCATT
CAGGTGTATGGATTTTGCAGTATTAAAGATGCTGCCTGGAATAGGTCATTATCTTC
TCCAAGTACTCTGTTAAGTCAATGAGTCACATAGAGTATAAGGTTTATTATCTGCT
TTTCTTTTATTAAATAAATCTTTATTGAATTTCTACTACATTAATAAAACCAAACCA
AAACAAAACAAACAAAAAAACACTTCCCTGAGCCATAAAGGAGAAGGTAGTTTTG
30 ACTGGAACCTTGAAGGATGGGTAACTTTTCCAGCAGATAAAGATTGAGAGAAGACCT
TCCAGGTAGAGAAAGCAGTGTGGGCACAGGCAAAGATGGAAGAACACACGTGGCTG
TGGGAAACACAGCTAGAAGCCAGTGCGGATAGAGAGTAGGCTATGATGTGCAAAGG
TTANACACTGGGAGAGACAGGTCCATGAGAGTAGCTTGGACTAACACAGGGAGGGT
TTGGAATCCCAACTGGGGAACCTANAAATCAA

35

Sequence ID 461

TAGGAGGCTTATTCACTGATTTCCCCTATTCTCAGGCTACACCCTAGACCAAACCT

- 159 -

ACGCCAAAATCCATTTCACTATCATATTCATCGGCGTAAATCTAACTTTCTTCCCA
CAACACTTTCTCGGCCTATCCGGAATGCCCCGACGTTACTCGGACTACCCCGATGC
ATACACCACATGAAACATCCTATCATCTGTAGGCTCATTCAATTTCTCTAACAGCAG
TAATATTAATAATTTTCATGATTTGAGAAGCCTTCGCTTCGAAGCGAAAAGTCCTA
5 ATAGTAGAAGAACCCTCCATAAACCTGGAGTGACTATATGGATGCCCCCACCCTA
CCACACATTCGAAGAACCCGTATACATAAAAT

Sequence ID 462

TCTTTATCAAGTTGAGAAAGTTCCTCCCCTCTATTCTTAGTTTGCTAAGAGTCCTT
10 CTATCCTATTTCTTAATGGTTTAGTAGATGACTCTGTGGTACTTTGAAGGTTGTTT
GCAGAATTTCCATGCCATAGGCAATTTACCTTTCCTTGACATTTGAAGGATTGATG
TTGGTGCCAAGTATAGAATCTTCACAGAGTCCTCCTGTAGCTTCTAAAGGTTTAGC
TTGAAAATGTTAATTGCTTAACGCTAGTAAGTGAGTGAAAAAGCTGGGGATAAATT
TTGTATCTTGCTTATATTTTCAAGTTCCACCTCTGTCCNGACNAAACCCCATATAT
15 AA

Sequence ID 463

TAGTTTACATATCCCAACCTTTAAAAATATTCCTCTTATTAGCTTTATATTTCACTT
TATAGAAGTTGAGTTTTAATTAAAAATTCTTGGCATCCTGAAGTATGTCACATAGCA
20 TGTGCTCCTTATAAATATGTTGATATCTCAGAAGACAGCATCCCGGTTTTTCAATTTT
ATAAAGTACCATACTTAAGAATGCTGTAATACTTATCTTTTATAACATGTTTCCTT
CGCTTTGCTTGNCTTTTATGNCATCAGTTTTAACTGTTTACTTCATTTAACAGNTT
ACATCATNCAACAGTTTACTTCATTAAACAGTAGGTGGAAAAATAGATGCCAGTCT
ATGAAAATCTTCCCATCTATATCAAATACTTTCAAGGATATACTTT

25

Sequence ID - 464

nt: 615

CGACTTTCAACCATCAAGTGAGGAATACCTTCACATAACTGAGCCTCCCTCTTTAT
CTCCTGACACAAAATTAGAACCCTCAGAAGATGATGGTAAACCTGAGTTATTAGAA
GAAATGGAAGCTTCTCCACAGAACTTATTGCTGTGGAAGGAACTGAGATTCTCCA
30 AGATTTCCAAAACAAAACCTATGGTCAAGTTTCTGGAGAAGCAATCAAGATGTTTC
CCACCATTAAAACACCTGAGGCTGGAAGTGTATTACAAGTCCCGATGAAATTGAA
TTAGAAGGTGCTACACAGTGGCCACACTCTACTTCTGCTTCTGCCACCTATGGGGT
CGAGGCAGGTGTGGTGCCTTGGCTAAGTCCACAGACTTCTGAGAGGCCACGCTTT
CTTCTTCTCCAGAAATAAACCTGAACTCAAGCAGCTTTAATCAGAGGGCAGGAT
35 TCCACGATAGCAGCATCAGAACAGCAAGTGGCAGCGAGAATTCTTGATTCCAATGA
TCAGGCAACAGTAAACCCTGTGGAATTTAATACTGAGGGTGCAACACCCCATTTTC
CCTTCTGGAGACTTCTAATGAAACANATTTCTGATTGGCATTAAATGAANAGTCA

- 160 -

Sequence ID 469

GATTTTTTAAAAATACATATAGCAAAAATATTACAGGGTCAGGGGAGACAATTAGAA
TGATATAATTCAAAGTGGATTAAAAAAAACCTGTCACCCAGAATACAATACCCAG
CAAAGTTGTCCTTCATAAATGAAAGAAAAATNAAATCTTTNCCNAACNA

5

Sequence ID 471

TCCCGGGAATCTGCAGGATCCGTCGACT

Sequence ID 472

10 GACAGTGCCCAGGGCTCTGATATGTCTNTCACANCTTGNAAGTGTGAGACAGCTG
CCTTGTGTGGGACTGAAAGGCAAGATTTGTTCTGCCCCTTCCCTTTGTGACTTGAA
GAACCCTGACTTTGTTTCTGCAAAGGCACCTGCATGTGTCTGTGTTCTTGTAGGCA
TAATGTGAGGAGGTGGGGANACCACCCACCCCATGTCCACCATGACCTCTTNC
CACNCTNACCTGTGCTCCCTCCCCAATCATNTTT

15

Sequence ID - 473

nt: 694

TGGGCTTTGGGCTGGCTGCAGTCTGTCTGAGGGCGGCCGAAGTGGCTGGCTCATTT
AAGATGAGGCTTCTGCTGCTTCTCCTAGNGGCGGCGTCTGCGATGGTCCGGAGCGA
GGCCTCGGCCAATCTGGGCGGCGTGCCAGCAAGAGATTAAAGATGCAGTACGCCA
20 CGGGGCCGCTGCTCAAGTTCCAGATTTGTGTTTCTGAGGTTATAGGCGGGTGTTT
GAGGAGTACATGCGGGTTATTAGCCAGCGGTACCCAGACATCCGCATTGAAGGAGA
GAATTACCTCCCTCAACCAATATATAGACACATAGCATCTTTCCTGTCACTCTTCA
AACTAGTATTAATAGGCTTAATAATTGTTGGCAAGGATCCTTTTGCTTTCTTTGGC
ATGCAAGCTCCTAGCATCTGGCAGTGGGGCCAAGAAAATAAGGTTTATGCATGTAT
25 GATGGTTTTTCTTCTTGAGCAACATGATTGAGAACCAGTGTATGTCAACAGGTGCAT
TTGAGATAACTTTAAATGATGTACCTGTGTGGTCTAAGCTGGAATCTGGTCACCTT
CCATCCATGCAACAACTTGTTCAAATTCTTGACAATGAAATGAACTCAATGTGCA
TATGGGATTCAATCCCCACCATCGATCATAGCACCCCTATCAGCACTGNAAACTC
TTTTGCATTAAGGGATCATTGC

30

Sequence ID 474

GGCAGCGCGGGGAGCCCGTCCGGCGCCGGCGGGCGGGCCGGTTTTCGAAGTTGATGCA
ATCGGTTTTAAACATGGCTGAACGCGTGTGTACACGGGACTGACGCAACCCACGTGT
AACTGTCAGCCGGGCCCTGAGTAATCGCTTAAAGATGTTCTACGGGCTTGTTGCT
35 GTTGATGTTTTGTTTTGTTTTGTTTTTGGTCTTTTTTTGTATTATAAAAAATAAT
CTATTTCTATGAGAAAAGAGGCGTCTGTATATTTTGGGAATCTTTTCCGTTTCAAG
CATTAAGAACACTTTTAATAAACTTTTTTTTTGATAATGGTTAAAAA

- 161 -

A

Sequence ID 475

5 CATAATAAAAAACAATCAACAAACAGGGAATGGAAAGAACTTCCTCAGCATGGTG
AAGGCCACATATGAAAATCCCACAGCTAACATCATACTCAATGATGAAAGACTGAA
AGCTTTTCTCCTGAGATCAGGAACAAGACAAAGATGTCACCTTTTGTCACTTCTAT
TCAACTCATTATTGGAAGTTTTTGCCAGAGCAATTAGGTAAG

10 Sequence ID - 476

nt: 476

CAGAATCTTTTCATAGGCTGAATGTTGCTCCACAATGTGTCTTTGACTATCTCTG
GCTAATTATTATTTTAATCTCTTCTCAGCTTTTCCAAGAACATAACGTTAACCAAA
GATCTTAGGCCATTCACAACTCTTTTGTAATAAATTAATGTGGATGTGAAACGAGGC
AACAAATCCTGAAGTAGAAAGTTATTCCTGGCCAGGCACGGTGGCTCACGCCTGTA
15 ATCCTGGCACTTTGGGAGGCCGAGGTGGGTGGATCATGAGGACAGGAGATCGAGAC
CATCCTGGCCAACATGATGAAACCCCATCTCTACTAAAATACAAAAAATTAGCTGG
GCATGGTGACGCGTGCCTGTAGTCCCAGTTACTCGGGAGGCTGAGGCAGGGGAATT
GCTTGAACCTCGGAGGTGGGAGGTTGCAGTGTGCCGAGATCACGCTACTGCACTCC
AGCCTGGCAACAGAGCAAGACTCCATCT

20

Sequence ID 477

AAACAGAAAGTTTCTTCTAAAGGCATGATTCAGTTAAGTCATTCTTAAGTGTTAAA
AAATTGTGAAAAATGTGCCTGTAATCCCAACACTTTGGGAGGCCGAGGCAGGCAGA
TCACGAGGTCAGGAGATCAAGACCATCCTGGCTAACAAGGTGAAACCCCGTCTCTA
25 CGAAAAATACCAAAAACATTAGCCGGGCGTGGTTGTGGGCGCCTGTAGTCCCAGCT
ACTTGAGAGGCTGAGGCAGGAGAATG

Sequence ID 478

TTCTTGGGATATTGATGACTACTGTCTGAGAGGTGCTGTGGGGAGATTTTCAGGAT
30 TGTGTGGTCTTTGAGGGGGGTGTTTTTTTAAGACAACATTGACCACTGTCCACTGT
CCACATGATCATTGTAAAATTGCAATGCCGCATGCTAGTTGGTTACATAAGACATA
ATTCCAGTGATTGAAGGTGGTTACACTGTATGGTGGTGTGTTCAAGATGGCACTGG
CATCTTTGAGCAGAGCCTGGCTATGCAGCATCATTGAGTTTTTTTAAACACCCTAN
AGGTCTGGTTGTTGTTGCTGTTGTCTTTTCTGTGAAAGTCACAAANANAAGTTACA
35 GTCCAGGTGAACCTGGAGTTTATAGGTTGGTTTTTGTCTGNTATATATATATATA
TATATATTTTTTTTTTTTTTTTAAACATTTACCTGTAGTGCTGTAGCTGTTGATACTA
TCACCTGCATGCTATTTCTAGTGAGTGCTAAATACAGTATGGTCCAATGACAATAA

- 162 -

CAGCCCATGGTACTGCCAG

Sequence ID 479

5 CATCAGTCTGTTATCCATGCTGACTTTCGGAAGACTTGCAGCTACTGCATTGATAT
CTTTCCTGCCAATAAGCAAAGTGTTGAACACTTCACAAAATATTTTACTGAGGCAG
GCTTGAAAGAGCTTTCAGAATATGTTTCGGAATCAGCAAACCATCGGAGCTCGTAAG
GAGCTCCAGAAAGAACTTCAAGAACAGATGTCCCGTGGTGATCCATTTAAGGATAT
AATTTTATATGTCAAGGAGGAGATGAAAAAAAAACAACATCCCAGAGCCAGTTGTCA
TCGGAATAGTCTGGTCAAGTGTAATGAGCACTGTGGAATGGAACAAAAAGAGGAG
10 CTTGTAGCAGAGCAAGCCATCAAGCACTTGAAGCAATACAGCCCTCTACTTGCTGC
CTTTACTACTCAAGGTCAGTCTGAGCTGACTCTGTTACTGAAGATTAGGGAGTATT
GCTATGACAACATTCATTTTCATGAAAGCCTTCCANAAAA

Sequence ID 481

15 CACACTTTCATGATAAAAAACAGAACCTAGGAATGAAAAGAAATTATAGCAACATAA
TAAAGACCATATATGAGAAGCCCACAGCTAACATACTGTATGGTGAAAACTGAAA
GCTCTTCTCTAAGATCAGGAACAAGGCAAGGATGCCCATTCCTTGCCACTTCTATC
GAACGTAGTACTGGAAGCCCTAGCCAGAACAAC TAGGCAATAGAAAGAAATTAAAG
GCATCCATNTCAGAAAGGAAGANCAAAATGCTGTCTGTTTAANATGACA

20

Sequence ID 482

TTTCTATANAAAAAAATTTTTTAAAATAATTGTAAAGTTAGATTTAAAATTGTAAA
ATATAAAATCACAAAGGAATGTACCCAATAAAATGTAAATGCNCCATAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAA

25

Sequence ID 483

CGNTAACGTGCAATCCGCCGCACGCCAGCAAACCTGGACAAACTCCGGGATCTCATC
GAAGCGATTGAGCACCAGTACCAGAGTAATACCGGACTGATGTAACGAGGCGAGTC
GCTCATCCAGCTTGCTGACGTGAGGCAACATCCAGGCCATCGAACGGNTCATCAAG
30 AATCAACAAGTCAGGCTCCGACATCAGCGCCTGACACAGCAGGGTTTTTTCGCGTCT
CGCCAGTGGAAGGTATTTAAAGCGTCNGTCGAGGAGGGCGGTAATACCGAACTGC
TGCGCCAGTTGCATGCAACGCGGTGCATCCTTTACTTTCATCCTGAATGATCTCAGC
CGTAGTGCGTCCGGTGCCATCTTCGCCAGGGCCGAGCATATCGGTGTTATTCCGCT
GCCATTGCTCGCTGACGAGTTTTTGCAATTGCTCGAAGGAGAGACGAGTGATGTGG
35 GAAAACTGGCTTTGCCGTTACCTTTCAAAGCGGGAAGTTCCCCCGCCAGCGCGC
GGGCCAGGGCCCGAT

- 163 -

Sequence ID 484

TTTTTTTTTTTTTATTCTATTAAAAAATGTTNNTGAAAAAAGATACTTAAATTTTAA
AGATAACTNAATTCCTAANGATTTAAAAATAATCCAAGCAGAGATGAAAGANCAAAT
GCAAATGCNTAAAAAGACCCCANAGCATTGTTAGCAAAAAGCAAATATAGTTAGCC
5 AAGCATATATATNTCATAAAAGCAATAANAAGGCNTAAAGCAAGTTTGGGGAGAGC
TTATTTAAAACTTGTA AAAAATCATTTGAATTTTAAAGTTTTCAAAC

Sequence ID - 485

nt : 551

TTTGGAACACAAAGTTCCCTTTTGTAGAAGAATAGGTATTGAGCCCTTGAGCGTGGG
10 TAGAAAGATAGAGACAGAGTGATTTGCAAAATAATGGAGGATCATATTTATATATG
AATTTTCACTTATTTGAACTTTTCAGATATCANCTTNAAAANC'TTTGGTTTAAGTAA
AGTNTNTTAATGAGACTCCTTGGATGAAAGTAACCAAAACCAGTAAAAATAAGGTA
ATAAGGATGTAATAGTTTCTTATGGACACTCAACAGCTAGAATGCAGTTAGTCTCA
GAAAAGAATTAGAACAAATAACTGGAAGGCCATCAGGAGTCCAAAACCATCACTCT
15 TTTATATTTTATATTTTATTTTCTCTCTTCANATGAGCATTCTCTTCTATGTCC
ATATGGTANAAGGCGGCAGCTCCATAGATTATGGCTTCAGATGTTACAGTTCCGCT
NAATGCAGGGACAGACTTGCTATCTTTCAGTCCCCTTACATATCCTGGGGAGAGAG
CAAATGATTGACTGGCTTGAGTCAGGTGCCCGTTCCCTTTCCAATCT

20 Sequence ID - 487 nt:224

GTTTGNTTGTGACCATCTGTACTTGTAATTTCTTTACNTTCATTGGTATGAAAAAT
ATGTTCTTAGAAGCANGAAAAAGAATTCAGNTTTGCTTTGTATACTAAATTAAATG
CTGTAATTTTGATAAAATGAAAAATCTGCTTTATTTGCAACAATTGGTTTCTTCCT
TGACGTCAGCCTCACTCTTGGACTTTGGTATTCAGCCNGNCACCCCTGGGAATTCC
25

Sequence ID - 488

nt : 349

GTGCCTCCCTGTGTGAGTAGCCTAAGGTGCATTGAAAAAGACTGGGATGTGTTTTA
TTTTTTTGTATTAGATAGCATTAACCTTACTGTTGAAGTATTTTGGTGGAGTATT
AGTGACAAGCCATTGAGTCTTAAGCCTTACGGCTTCCTATAAAATCACTAATTTTCG
30 TGTGTGTTTGTGTGTAGGTTACGTTATATATAGGATTCGTGTTCCCGTGGTGGCC
GAAAACGCCCAGTTCCTAAGGGTGCAACTTACGGCAAGCCTGTCCATCATGGTGTT
AACCAGCTAAAGTTTGCTCGAAGCCTTCAGTCCGTTGCAGAGGANCGAGCTGGACN
CCCTGGGGGGGCTC

35 Sequence ID 489

TTAACAGCTGCATAGAGTTTTTAAAGTACATTATATTTTGTGAGACAAGTAAAAATA
TCTGTTTTTTCACGCAAAAAAAGCCATGAAATACGTAATTTTTTAAAGACAAAAAAT

- 164 -

CATCTTTTGAGTTTGCTCTTTGGTTTTCTTCATTCCTTTGAGGATTGGGAAAAC
AGAAAGATTCTTTGATTTGGGTAATGAAGAGGTAATTTGGGACAGTGTGGTGGTAC
CAGGAAGAAAGAGGATTGGAAAGGCCAGTACTGTTTTAGTTGCTCGGCACTGTTGG
TTTTGTTTTAATGTGGTTGCCCTGTCCACTACATGGTTCATCAGTAGTGTAATCC
5 ATTTTCAATGTAAAGCTCTTTTAGTTTTTGTTCATAGACATAAATTAATATTTTGAG
AGGCATCCCTCACCTGTTTCATTTCTTCTGTGTTGAAATGAAGTACTTAAAATTACC
GTTATACATGAACTTTGTGGACTGTAAGATTTGTTATATATGTTCAAATGCCTTTT
AGCTGGCTTTTTAATTAATATGCCTGTTTTGAGTGCTTAATACAATGTAATGNNGA
TTGTAAATCATACCTATTTTAAATCATTCTTCCTGTATATTTGNACTCAGAGAGC
10 CTTATTTTATTCTTCCAGC

Sequence ID - 491

nt: 382

TTTTCTTAGAACTTTATTTTTTCTGGCCAGGCGCAGTGGCTCACACCTGTAATCCC
AGCACTTTGGGAGGCCAAGGCAGGTCGATCACCTGAGGTGAGGAGCTCAAGACCAG
15 CCTGGCCAACATGGTGAAACCCTGTCTCTACTAAAAATACAAAAATTAGCTGGGCG
TGGTGGCGCATGCCTGTAATCCCANCTACTCAGGAGGCTGAGGCAGGAGAATTGTT
TGAACCCGGGAGGCGGAGGTTGCANTGAGCCGAGATTGCGCCACTGCCTCCAGCC
TGGGCAACAGAGCGAACTCCATCTCAAAAAAAAAAAAAAAAAACAACCTTTATTT
TTTCTGATTTTAAAGTAATAACTAGTTTGTAGAAACATTAAAGT
20

Sequence ID 492

ACCCTAAACATAACTTAAAATTTGTTNGGAATTTGAAAGTACAGAATTTTCCTGTA
ATTGAGACTNTTTAAACTTTTGTGGTTGGAGAAGGTATTCTATTTTTTGAAAATAT
CTGTAAGTTTTATCTAAATAGTAACTCTAAGTATTCTTCCCCTTTACTTACAGCC
25 ACCCTGGGAATCTGAGACTAGAGAAAATAAGTTTGTCTCTTGTTCTAAGGAGGGT
CTGGTTTAGAAATCTGATTTAGACATAGAAAAATTGCAAGAAGCTTGAGGTGATTG
GAAGATACGATTTTGTATCAAAGNATGTTTCTGTTTTATAGATTTTATTCATCTA
CAACTCCTTATTAATATATTTAAGAAGTCATTAACCCACCATTGATTACTTGATAT
AAAAGGAGAANC GGTTGGTAAAAGGTGAAATANAATTTTTTAATTTTTTTTTTTTAA
30 GTTTAGGATTTTTTTTTTAAATTCTAAGAGTTTCTGTCATTTGGGGACAATCAGAA

Sequence ID 493

TGGGAATCATAATTNGTTAACTGAAGCTNATAAGATGAGAGCATTCANAGAGAAAA
GAACGGAAAGATTGAATATCAGTTTCCCTTCTTTAAAAAAATTGTGGATATGTGAT
35 CTAGCTTCTTGAGCATCACAGTGACTGATTGGCTCGTGGTAATTGATCGCTATGCT
GACAATCTTATCTCCACCTATGTCATTCAATTTTCTAAGAGGCAAAATCCTTAATC
AGGAGGAGAGTTTAGCTCTAGCTAAATTTCCCTTGTCAGCATGCTCCTGCTCCCC

- 165 -

CAACTTGTGGAAACAGCTAAAGGATTGGACTAGGAGCANAAGTTTGGAATGGTTAA
AATGTAGCAACATGTGTTTCCTGAAACAAAATTCCACTATAATAAAAAAAGCATTT
GAATGCTCCCTTGTAATTCTGTTGGAGCTTGTTGCCTTTTTTATGACACAACCATA
ATCAGTGATAGACAGTAGCATAAAGAAGCAAGAGCAAAGCAATTAAGTAATAATAG
5 CACTACAAAAATGTGTGCTGTACTTACCAAACACGACATTTATGAATTATTANATA
GGAATAAGGGGATGGT

Sequence ID 494

10 GACCCAGCCATCTAAATAAGTTRTACATGTTGCGTATTTTTTTTGTAGGGACTTAT
CTTCCGAAGAGGAAAGGTTTATGAAACCTAAAGTAACAATGATAGCTTGGAATCAA
AATGATAGCATTGTTGGCACAGCTGTGAATGATCATGTCCTCAAAGTGTGGAATTC
TTACACTGGACAACCTGCTTCATAACTTAATGGGACATGCTGATGAAGTATTTGTTC
TGGAGACACATCCCTTTGATTCCAGAATTATGTTATCTGCAGGACATGATGGCAGC
15 ATATTTATATGGGATATTACAAAAGGTACCAAGATGAAACATTATTTTAATATGGT
AAGTGAAGTGAGATGTACCTTGATACATGCTTGATAATTTGTTTAGAGTATTTGGG
TTATGCGGCTTACCCAGAAATTGATCTGCTTGTTTTGGCAGTTTGTTTTTACAAAT
CAACATATTCAAAGCCTGCTAAATATTAGACAGCTACATGTATATACGTACATACA
TGAA

20

Sequence ID 495

TTTC

Sequence ID 496

25 CTCGCTGGCGGGAGGCCACGGGCTTTCCACAGCGCGGGGGAACGGGAGGCTGCAGG
ATGGTCAAGCTGACGGCGGAGCTGATCGAGCAGGCGGCGCAGTACACCAACGCGGT
GCGCGACCGGGAGCTGGACCTCCGGGGGTGATCTGGACCTCTGGCATCTCTCAA
TCGCTGACTTACCTAAGTATCCTAAGAAATCCGGTAACCAATAAGAAGCATTACAG
ATTGTATGTGATTTATAAAGTTCCGCAAGTCATAGTACTGGATTTCCAGAAAGTGA
30 AACTAAAATTTTAATCCAGGTGCTGGTTTGCCAACTGACAAAAAGAAAGGTGGGCC
ATCTCCAGGGGATGTAAAAGCAATCAAGAATGCCATAGCAAATGCTTNAACTCTGG
CTGAAGTGGANAGGCTGAANGGTTGCTGCAGTCTGGTC

Sequence ID 497

35 GAAGACCTCACATCTGAGAGCTCATCTGCGTTGGCATTCTGGAGAACGCCCTTTTG
TTTGTAACCTGGATGTACTGTGGTAAAAGATTTACTCGAAGTGATGAATTACAGAGG
CACAGAAGAACACATACAGGTGAGAAGAAATTTGTTTGTCCAGAATGTTCAAACG

- 166 -

CTTTATGANAAGTGACCACCTTGCCAAACATATTAAAACACACCAGAATAAAAAAG
GTATTCACCTCTANCAGTACAGTGCTGGCATCTGTGGAAGCTGCGCGAGATGATACT
TTGATTACTGCAGGAGGAACAACGCTTATCCTTGCAAATATTCAACAAGGTTCTGT
TTCAGGGATAGGAACTGTTAATACTTCCGCCACCAGCAATCAAGATATCCTTACCA
5 AACTGAAATACCTTTACAGCTTGTACAGTTTCTGGAAATGAGACAATGGGAGTA
AATATTACACAAATACTTATTTCATTGNGGTTATTTTTATACAGTAGTGAGAAGAAT
ATTGTTCTTAAGTTCTTAGATATCTTTTTTTGGATGTGCAAAAATTTTTGGATTGA
CAGTAACTTGGGTATACATGACACTGAAATGCCTTACTTTGGATGA

10 Sequence ID 499

TGCCTGCGGGCCAGGACCTCGCCCAGCCCATGTTTCATCCAGTCAGCCAACCAGCCC
TCCGANGGGCAGGCCCCCAGGTGACCGGCGACTGAGGGCCTGAGCTGGCAAGGCC
AAGGACACCCAACAACAATTTTTGCCATACAGCCCCAGGCAATGGGCACAGCCTTCC
TCCCCANAGGACCCGGCCGACCTCAGCGCCTCCTGCAGGCTAGGACACTGGTGCAC
15 TACACCCCATGCCTGGGGGCCGAGATTCTCCAGCAGAAAGATGCAATATTTTTTGT
TTCTTTTTTTTCCATTTTTTTTCTCTAAGGAATCAATATTTCAATATGTTGAGTGTG
TGTCCAATGCTATGAAATTAAATATTAATAACATATTTATGGCATTTTCTTGAA
GAGTGTGGTTGAAGAAATATTTCTCCTTTTGTTTTTCTTTTTTTTTGNTTGNTAC
TGCCACTTCTTTTTTAGGAGCAAATCTCCCCAGGGGTGTACGGNATTTCTTGACTCT
20 GGGAACAGCTGCTACCCCCAAGACTTGCCACGTTGTTCTGCCCTCAAATGGAATTA
AGTG

Sequence ID - 500

nt: 390

GGAATATGGTCAGGATCTTCTCCATACTGTCTTCAAGAATGGCAAGGTGACAAAAA
25 GCTATTCATTTGATGAAATAAGAAAAAATGCACAGCTGAATATTGAACTGGAAGCA
GCACATCATTAGGCTTTATGACTGGGTGTGTGTTGTGTGTATGTAATACATAATGT
TTATTGTACANATGTGTGGGGTTTGTGTTTTATGATACATTACAGCCAAATTATTT
GTTGGTTNATGGACATACTGCCCTTTCATTTTTTTCTTTTCCAGTGTTTAGGTGAT
CTCAAATTAAGAAATGCATTTAACCATGTAAANATGANTGCTAAAGTCAGCTTTT
30 TAGGGCCCTTTGCCAATAGGTANTCATTCAATCTGGTATTGATCTTTTCACAAA

Sequence ID 502

ACCCGCCATCTTCCAGTAATTCGCCAAAATGACGAACACAAAGGGAAAGAGGAGAG
GCACCCGATATATGTTCTCTAGGCCTTTTANAAAACATGGAGTTGTTCCCTTTGGCC
35 ACATATATGCGAATCTATAAGAAAGGTGATATTGTAGACATCAAGGGAATGGGTAC
TGTTCAAAAAGGAATGCCCCACAAGTGTTACCATGGCAAACTGGAAGAGTCTACA
ATGTTACCCAGCATGCTGTTGGCATTGTTGTAAACAAACAAGTTAAGGGCAAGATT

- 167 -

CTTGCCAAGAGAATTAATGTGCGTATTGAGCACATTAAGCACTCTAAGAGCCGAGA
TAGCTTCCTGAAACGTGTGAAGGAAAATGATCAGAAAAAGAAAGAAGCCAAAGAGA
AAGGTACCTGGGTTCAACTAAAGCGCCAGCCTGCTCCACCCAGAGAAGCACACTTT
GTGAGAACCAATGGGAAGGAGCCTGAGCTGCTGGAACCTATTCCCTATGAATTCAT
5 GGCATAATAGGTGTTAAAAAATAAAGGACCTCTGGG

Sequence ID - 503

nt: 109

ACATTTTCCGGNCCTTTTGCCATACACAGTTACAGAGATCAGTCAAATCCATACCA
10 CCACTGAGATCTCATTTATTGCCACAGATGCACAAAATAAATAACCCAAAATC

Sequence ID - 504

nt: 374

CCAGCAACGACCCATACCTCAGACCCGACGGCCCGGAGCGGAGCGCGCCCTGCCCT
GGCGCAGCCAGAGCCGCCGGGTGCCCGCTGCAGTTTCTTGGGACATAGGAGCGCAA
AGAAGCTACAGCCTGGACTTACCACCACTAACTGCGAGAGAAGCTAAACGTGTTT
15 ATTTTCCCTTAAATTATTTTGTAAATGGTAGCTTTTCTACATCTTACTCCTGTTG
ATGCAGCTAAGGTACATTTGTAAAAAGAAAAAACCAGACTTTTCANACAAACCC
TTTGTATTGTANATAAGAGGAAAAGACTGAGCATGCTCACTTTTTTATATTAATTT
TTACAGTATTTGTAAGAATAAAGCANCATTTGAAATCG

20 Sequence ID 505

GTACAGGAGGTAAATTGGATACCCCATCTAAGGGGATCTGTGAGACCAGGTAGTTA
TTTGGAATGAAAGAGTAAGATATTAAACCAGCCAGCATGTCAACAGGTGGGTGATA
GTCTTGTTCTCACAGACAACAGATGGCCATCATCTTAAACAACATTTATGTTAAC
25 CAGCAGATAAGGGACTCCTGCATTGTGAGTGGACTTTGAGCCTGAGTTTTTCTACT
TGCATAGGTGAAAGTGGACTGCAATGCTAGTATAAATGCCGTATGATGACTAGTAC
CCCTTAGGGAGCTCCAGTTTGCCTTCCTGGGGAACCACAGACCCCAAGTGTAATTT
CCTGAGGACAGCCCGACTTCT

30 Sequence ID 506

GTTACTGTGAGCCTGTCAGTAGTGGGTACCAATCTTTTGTGACATATTGTCATGCT
GAGGTGNGACACCTGCTGCACTCATCTGATGTAAAACCATCCCANAGCTGGCGAGA
GGATGGAGCTGGGTGGAACTGCTTTGCACTATCGTTTGCTTGGTGTTTGTTTTAA
ACGCACAACCTTGCTTGTACAGTAACTGTCTTCTGTACTATTTAACTGTAAAATGG
35 AATTTTGAAGTATTGTTACAATAATAACTCTGAGATGTGTGAAAAAAAAAAAA
AAAAAAAAAAAAA

- 168 -

Sequence ID - 507

nt: 521

CTGCGGTGGAGCCGCCACCAAAATGCAGATTTTCGTGAAAACCCCTTACGGGGAAGA
CCATCACCCCTCGAGGTTGAACCCCTCGGATACGATAGAAAATGTAAAGGCCAAGATC
CAGGATAAGGAAGGAATTCCTCCTGATCAGCAGAGACTGATCTTTGCTGGCAAGCA
5 GCTGGAAGATGGACGTACTTTGTCTGACTACAATATTCAAAGGAGTCTACTCTTC
ATCTTGTGTTGAGACTTCGTGGTGGTGCTAAGAAAAGGAAGAAGAAGTCTTACACC
ACTCCCAAGAAGAATAAGCACAGAGAAAGAAGGTTAAGCTGGCTGTCCTGAAATA
TTATAAGGTGGATGAGAATGGCAAATTAGTCGCCTTCGTGAGAGTGCCCTTCTG
ATGAATGTGGTGGTGGGGTGTATGGCAAGTCACCTTGACAGACATTATTGTGGC
10 AAATGTTGTCTGACTTACTGTTTCAACAAACCAGAAGACAAGTAACTGTATGAGTT
AATAAAAGACATGAACT

Sequence ID 508

AAGCTCATGATTTTAAATGTATTTTTCTAATAAACTATACTCCCATTTAAAAATCA
15 CCAATACCTTAATGTTTCAATTATATAAGCTAATTAAAAATAAAGGCTGGGCGTGG
TGGCTCACTTTGGAAGACCGAGGCAGGCAGATCACCTGAGGTCAGGAGTTCGAGAC
CAGCCTGCCCAACATGGAGAAACCCCATCTCTACTAAAAATACAAAATTAGCCAGG
CATGGTGGCACATGCCCGTAATCCAGCTACTGGGGAAGCTGAGGCAGGAGAATCA
CTTGAACCTGGGAGGCAGGGGCTGCAGTGAGCCGAGATCATGCCATTGCACTCCAG
20 TCTGGGCAACAATAGTGGAATCCATCTCAAAAATAATAAAAAAATAAAATAAAA
ATAAAATTCAAACCTAAATAGATGCTCTACTTCAGGAGTGGGCAAATTAATCACC
TGCATCCTTTTTTTGGGCTTTC

Sequence ID - 509

nt: 575

TTTTTTTCTAAATGGNGATTACTAATATATGTGGAGACTATTAATCTCTTTCTGT
25 TGCCATTAGTTCATTTTTTCCCCAAAAGCCAATACATGTTTATTACAAAAATGAATT
ATAAAATATAAGTTAAAAGAAAAACATAAAACCCTACAATCTTACCCACCCAGACA
ACTACTATTAATACCTTAGTATTAACATATACACATCATGTATATGTATAAATTTA
TCTTAAACAAAAATAAAATTATTCTTTACATATTGTTTTAAACCTATTTATCTGG
30 CCAGGTGCCGTGGCTCACGCTTGTAATCCAGCACTTTGGGAGGCTGAGGCACGTG
GATCACCTGAGGTCAGGAATTCGAGACCAGCCCAGCCAACATGGTGAAACCCTGTC
TCTAATGGTTTAAATACCAAAAAATTAGCTGGGCATGGTGGCACATGCCTGTAATA
TCAGCTAACATGGGAGGCTGAGGCAGGAGAATCACTTGAACCANGGAGGGGGAGGT
TGCAGTGAGCCGAAATCACACCACTTCACTGCAGCCTGGGCAACAAAGCAAGACTG
35 TCTCAAAAAGAAAAA

Sequence ID 510

- 169 -

CACTGTCATTCCCAGGAGGCTTTGGAGTCAGAACTGGATTCAAATTCTGACTNTAT
GTTGTGTGACTTTGGGCCAATAGCTTCTTTNTGTGCCTCAGTTTCTTTAGCTGTAAA
TANACGGGTAGGTCACCCCTTACCCCATAGGTTATGGGGAAAGTTACAGAAAATGG
TCAGCTGGGCNCAGTGGCTCAAGCCTGTGGTCCCAGCNCCTTGGGAGGCCAAGGTG
5 AGCAGATTGCTTGAGCCCAGGAGTTTGACACCAGTNTGGCAACGTGACGAAACCCT
ATCNCTGTGAAAAATACAAAAAATTAGCCAGGCATGGTGGTGTGTGTCTGTGGTTC
CAGCTGCTTGAGAGTTTGAAGTGGGAGGATCACCTGAGCCCAGAAGGTCGAGGCTG
CAGTGAGCTGTGATCGCGTCACTGCACTCCAGCCTGGC-GACAGAGTGAGA-CCCC
T-TTTGAAAAAAAAAAAAAAAAAAT

10

Sequence ID 512

GTGAGCGGTGGTGGTTTATTCTTCCGTGGAGTTAAGGGCTCCGTGGACATCTCAGG
TCTTCAGGGTCTTCCATCTGGAACATATAAAGTTTCAGAAAACATGTCTCGAAGAT
ATGACTCCAGGACCACTATATTTTCTCCAGAAGGTCGCTTATACCAAGTTGAATAT
15 GCCATGGAAGCTATTGGACATGCAGGCACCTGTTTGGGAATTTTAGCAAATGATGG
TGTTTTTGCTTGCAGCAGAGAGACNCAACATCCACAAGCTTCTTGATGAAGTCTTTT
TTTCTGAAAAAATTTATAAACTCAATGAGGACATGGCTTGCAGTGTGGCAGGCATA
ACTTCTGATGCTAATGTTCTGACTAATGAACTAAGGCTCATTGCTCAAAGGTATTT
ATTACAGTATCAGGAGCCAATACCTTGTGAGCAGTTGGTTACAGCGCTGTGTGATA
20 TCAAACAAGCTTATACACAATTTGGAGGAAAACGTCCCTTTGGTGTTTCATTGCTG
TACATTGGCTGGGATAAGCACTATGGCTTTCAGCTCTATCAGAGTGACCCTAGTGG
AAATTCGGGGGATGGGAAGGCCACATGCATTGGAAATAATANCGCTGCAGCTGTGT
CAATGTTGAAACAAG

25

Sequence ID 513

TTTTTTTTTTTATAAACTCCAATCATTTCCAGAGCTACTTAGCTCAGCATCTTTTTT
TTCCACGCTCTTAAGTTGTGTTTATACATTTTTTGATACAGTTAGATTGTTTTTGT
ACATTCTTCAATTCTATCCTGGGATCCCCCAACCACCTAAGTGGATTTTTTGATAAT
TTGCATGCTTTAAGGATAACTCTTCATTCTGNAAAGGGCTATGGGTTTTGGCAAAT
30 GCAGAGTCATGTATCCAAGATTACAATATCGCACAGAAGAGTTTCATCACTATATA
AAACTCACCAGTCTTCCTCCTATTCAACCATCTCCATGCCTTCTTCCCAGCCCTAA
CTCCTTAAAACCACTCATATCTTTACTATTGCTATAGTATTGCCTCTTCCACCATG
TCATATAAATGGAAACATACAGTATTAGTCTTCTCAAAGTAGTTTCTTTTACCTAA
CAACATGCATTTAAGATTCATAGTGTCTTTAATGACTTGATAGATTATTTCTTTG
35 TAGCTGAATAATATTGCATCTTATAGATGTAACCGTTTGTATATCCATATTTCTC
ACAGCCTATGACTTGNCCTTTTGATTCTCTGAACAGGCCATTACAAAGCAGAAGTT
TTAATTTTTATAAAGCTAATGNATCAACTT

- 170 -

Sequence ID 515

CCTGGATGACAGCATATCTGTTTATAGCTCAGTTTACTGAATACTTTAAGCCCACT
GTTGAAACCTGCT

5 Sequence ID - 518

nt: 502

GATGCATGTCCAGCATAGGCAGGATTGCTCGGTGGTGAGAAGGTTAGGTCCGGCTC
AGACTGAATAAGAAGAGATAAAATTTGCCTTAAACTTACCTGGCAGTGGCTTTGC
TGCACGGTCTGAAACCACCTGTTCCCAACCCTCTTGACCGAAATTTCTTGTGACAC
AGAGAAGGGCAAAGGTCTGAGCCCAGAGTTGACGGAGGGAGTATTTAGGGTTTAC
10 TTCAGGGGCTCCCAAAGCGACAAGATCGTTAGGGAGAGAGGGCCAGGGTGGGGACT
GGGAATTTAAGGAGAGCTGGGAACGGATCCCTTAGGTTTCAAGGAACTTCTGTGCAA
GCTGCGAGGATGGCTTGGGCCGAAGGGTTGCTCTGCCC GCCGCTAGCTGTGAGC
TGAGCAAAGCCCTGGGCTCACAGCACCCCAAAGCCTGTGGCTTCAGTCCTGCGTC
TGCACCACACATTCAAAGGATCGTTTTGTTTTGTTTTTAAAGAAAGGTGANAT

15

Sequence ID 519

CTGCGATNGAGTTTGTGAGAGGAAGGANTAAAGTNCTCATCTCNGACGGTGAGAAAG
ATCATNACTAAGGAAACGCAGGGTTGGAAGCAGTGCTGANTGTCCAGTTGAGTTTC
ATGANCAAACATTTGCTGTGGGACCAGTTTTTCATGGNGGTTTGTCATTTTGTCCAG
20 CTGCCTGGAGCTGCTTGGTTGAAGGCACAGAATAATCAGGATTAATTGTTNAACTT
GTATGAATTTCTTTATTTTAAATAGGAATAATATCTGCCTTGGGAGCAAGTTGTA
AGAGTTAACTGAAAGCTTNAGGAAAACTTCCCTTGCTATTTAAGTAGGGCTTTA
CAAGTTACAATTCTATCACAGTTTTAAGATTATAAAC

25 Sequence ID 521

GCGGCGCANCTGCGGATCCANAAGGNCATAAACGANCNGAACCTGCCCAANNCGTG
TGATATCACCTTCTNAGATCCAGACNACCTCCTCAACTTCAAGCTGGTCATCTGTC
CTGATGAGGGCTTCNACAAGAGTGGGAAGTTTGTCTCAAAAAA

30 Sequence ID - 523

nt: 585

GATTTACTGTGGGAATTTGCTCATGCAATTATGGAAACCTAGAAGTCCCATAATAT
GCCATCTTCAAGCTGGAATCCCAGGAAAGCAGGTGGTGTAATTCTGAGATTGAAGT
CTTGAGAACCGGGGGAGTCAATGGTGTAACCTCCAATCTAGGGCTTAAGGCCCAAG
GACCAGGGCTGCTGGTGTGAGATGCAAATCCTGGAGTTCAAAGGATTGAGAACCA
35 GGAGCTCTGGTGTCTGAGGGCAGTAGAAGATGGATGTTCCAGCTCAAGAAGGGAAA
GTAAGAATCCGTCCTTCCCTCCACTTTTTTGTTCCTATTTCAGATGAGCCCTCAATGGA
CTGAACGATGCTCACCCACACTGTGAGGGCTGGTCTTCTTTATTCAATCCACTGAG

- 171 -

TTAAGTGCTGATCTCTTCTGGAAACACCTTCACAGACACACCCAGAAATAATGTTT
TACCAGCCATGGGCCTGTTACTTAGCCAGTCAAGTTGACACAGAAAATTAGCTAT
CACAACATCTGTGTGTGTATATACATATGTATTTGCATGTGTGTGTATATATGGNG
TATATATATTCATGTGTGTGTATAT

5

Sequence ID 524

CTTTTGCCAGTAGGCCCTGAGTAGGTTCTCTATCTTTTGGCATGACCCAGAA
GTCTTTGATAACTTCCTTGCTTTCTGATGTGACAAGACATCCAGGGCCAGATTGTC
CATATCCTGCCCCGGATGCACGATGCACGTGTTTCTCCAAGAATCCCTGTGTCCTTT
10 GCTGATGATGCCATGATTTTAAAGTTCTCTAATATAGTTTTATCTCTTTGTTTCAGA
TAATGCTTTTGTGTTCTCACATGTCCTGCTCTCTCTCTCTCTCTCATTTTGGTGT
GATCAGTCTTTCCATAAGATTGTTTATTTCACTAGTCCTTCATTCTTCTTTTCT
AAATTTACTCTTCTTGACTAGTATCCTGTCACTTCTGAGGACTCATATTTTGGAA
CTTGAAAATTATTCTTATTTATTTAAGTATATGTTNCTGAAACTCTCATTAGACAC
15 ATTTTG

Sequence ID 525

GTTAAAAAAGTAAAAGGAACTCGGCAAATCTTACCCCGCCTGTTTACCAAAAACA
TCACCTGGTAGCATCACCAGTATTAGAGGCACCGCCTGCCAGTGACACATGTTTA
20 ACGGCCGCGGTACCCTAACCGTGCAAAGGTAGCATAATCACTTGTTCCCTTAAATAG
GGACCTGTATGAATGGCTCCACNAGGGTTCANCTGTCTCTTACTTTTAAACCAGTGA
AATTGACCTGCCCGTGAAGAGGCGGGCATAACACAGCTGAAAAAAAAAAAAAAAAA
AAA
AAAAAATTTT

25

Sequence ID - 526

nt: 516

CTTTTCATGGTCTCTTGTTTATTAATCATCTAAAATCCAAGCNCAGAGAATTCAAT
TTTAGATGGTCTCCAGAGCAGAATTTGATGTATAATCTTAATTACAAATCATAGAT
AATTAATATTGNTTACAAAATCANAATACGATTAGAGGTAGGGATCCTGCACACAC
30 CCTATTTTCTCCTCCCCAGTGTTCTGACCGAGAGACTAATTAATAATTCAAGGAACCT
ACAGTGAATGANAACCCATGGTTTTTGCTTAATTATCAGAACAGCTAGATCTGAGAA
CAGCTGTCTCCACATGGATAGACACTTATTCACCCATTTGCAGGTAGAATAGCT
GGCAATAATAAGTCCTTCCCATTTGGATATGTTGAAAGGTGCCTGCCATGGCATAGT
TGCCACAAGAGAGGAAGAAATGGACACAAATGTAGGCTGTTTTCAGGGCANAGGGA
35 AGGTGGGAGGAAACCAANTTGCTGGTTTTTCACACACCCTCTGGGGAACACCCATGC
ACCTATGANATG

- 172 -

Sequence ID 527

GACAAAAGCTGAGAGAATTTTTTTCTTGAATATTTGCACTAAAAGATAGGTTAAAA
TTCTTCAGGCTGAAGAGAGCATACCAGGTGGAGATTTGGATCTACAAAAGGAAGG
AAGATTTGGAAATGGATTTGGCACCATTGACTCAATTTCCAGAACAAGAAAGCAGG
5 GACAGTTTTGGGAAGCTCAAGACACACTGCCCATGAGCAGCAATTTGGACCTCCTG
CTGCATCCACTGTGCATCAAACACACACTGTACAGACAAAGACTCCCAGGAAAAGA
AGTATAAACATGGACTAACACAGAGATGGGCAAACCTACAGCCTGTGACCCAGCCAC
CTGTTTATGTAGAATCCAAAGTAAGAATCTTTAACTTACACATAAACTT

10 Sequence 529; 660nt

GACAGCAGAGCACACAAGCTTNTAGGACAAGAGCCAGGAAGAAACCACCGGAAGGA
ACCATCTCACTGTGTGTAAACATGACTTCCAAGCTGGCCGTGGCTCTCTTGGCAGC
CTTCCTGATTTCTGCAGCTCTGTGTGAAGGTGCAGTTTTGCCAAGGAGTGCTAAAG
AACTTAGATGTCAGTGCATAAAGACATACTCCAAACCTTTCCACCCCAAATTTATC
15 AAAGAACTGAGAGTGATTGAGAGTGGACCACACTGCGCCAACACAGAAATTATTGT
AAAGCTTTCTGATGGAAGANAGCTCTGTCTGGACCCCAAGGAAAACCTGGGTGCANA
GGGTTGTGGANAAGTTTTTTGAAGAGGGCTGAGAATTCATAAAAAAATTCATTCTCT
GTGGTATCCAAGAATCAGTGAAGATGCCAGTGAAACTTCAAGCAAATCTACTTCAA
CACTTCATGTATTGTGTGGGTCTGTTGTAGGGTTGCCAGATGCAATACAAGATTCC
20 TGGTTAAATTTGAATTTTCAGTAAACAATGAATAGTTTTTTCATTGTACCATGAAATA
TCCAGAACATACTTATATGTAAAGTATTATTTATTTGAATCTACAAAAACAACAA
ATAATTTTTTAGATATAAGGATTTTCCTGGATATTGCACGGGAGA

Sequence ID 529

25 GACAGCAGAGCACACAAGCTTNTAGGACAAGAGCCAGGAAGAAACCACCGGAAGGA
ACCATCTCACTGTGTGTAAACATGACTTCCAAGCTGGCCGTGGCTCTCTTGGCAGC
CTTCCTGATTTCTGCAGCTCTGTGTGAAGGTGCAGTTTTGCCAAGGAGTGCTAAAG
AACTTAGATGTCAGTGCATAAAGACATACTCCAAACCTTTCCACCCCAAATTTATC
AAAGAACTGAGAGTGATTGAGAGTGGACCACACTGCGCCAACACAGAAATTATTGT
30 AAAGCTTTCTGATGGAAGANAGCTCTGTCTGGACCCCAAGGAAAACCTGGGTGCANA
GGGTTGTGGANAAGTTTTTTGAAGAGGGCTGAGAATTCATAAAAAAATTCATTCTCT
GTGGTATCCAAGAATCAGTGAAGATGCCAGTGAAACTTCAAGCAAATCTACTTCAA
CACTTCATGTATTGTGTGGGTCTGTTGTAGGGTTGCCAGATGCAATACAAGATTCC
TGGTTAAATTTGAATTTTCAGTAAACAATGAATAGTTTTTTCATTGT
35

Sequence ID - 530

nt: 660

GACAGCAGAGCACACAAGCTTNTAGGACAAGAGCCAGGAAGAAACCACCGGAAGGA

- 173 -

5 ACCATCTCACTGTGTGTAAACATGACTTCCAAGCTGGCCGTGGCTCTCTTGGCAGC
CTTCCTGATTTCTGCAGCTCTGTGTGAAGGTGCAGTTTTTGCCAAGGAGTGCTAAAG
AACTTAGATGTCAGTGCATAAAGACATACTCCAAACCTTTCCACCCCAAATTTATC
AAAGAACTGAGAGTGATTGAGAGTGGACCACACTGCGCCAACACAGAAATTATTGT
10 AAAGCTTTCTGATGGAAGANAGCTCTGTCTGGACCCCAAGGAAAACCTGGGTGCANA
GGGTGTGTGGANAAGTTTTTGAAGAGGGCTGAGAATTCATAAAAAAATTCATTCTCT
GTGGTATCCAAGAATCAGTGAAGATGCCAGTGAAACTTCAAGCAAATCTACTTCAA
CACTTCATGTATTGTGTGGGTCTGTTGTAGGGTTGCCAGATGCAATACAAGATTCC
TGGTTAAATTTGAATTTTCACTAAACAATGAATAGTTTTTTCATTGTACCATGAAATA
15 TCCAGAACATACTTATATGTAAAGTATTATTTATTTGAATCTACAAAAACAACAA
ATAATTTTTTAGATATAAGGATTTTCCTGGATATTGCACGGGAGA

Sequence ID 532

15 GAATTGTGATAGTTCAGCTTGAATGTCTCTTAGAGGGTGGGCTTTTGTGTGATGAGG
GAGGGGAAACTTTTTTTTTTTCTATAGACTTTTTTCANATAACATCTTCTGAGTCA
TAACCAGCCTGGCAGTATGATGGCCTANATGCAGAGAAAACAGCTCCTTGGTGAAT
TGATAAGTAAAGGCAGAAAAGATTATATGTCATACCTCCATTGGGGAATAAGCATA
ACCCTGAGATTCTTACTACTGATGAGAACATTATCTGCATATGCCAAAAAATTTTA
AGCAAATGAAAGCTACCAATTTAAAGTTACGGAATCTACCATTTTAAAGTTAATTG
20 CTTGTCAAGCTATAACCACAAAAATAATGAATTGATGAGAAATACAATGAAGAGGC
AATGTCCATCTCAAAATACTGCTTTTACAAAAGCAGAATAAAAGCGAAAAGAAATG
AAAATGTTACACTACATTAATCCTGGAATAAAAGAAGCCGAAATAAATGAGAGATG
AGTTGGGATCAAGTGGGATTGANGANGCTGTGCTGTGT

25 Sequence ID 533

CTTGAACCTCGGAGGCAGAGGTTGCAGTGAGCCGAGATCACGCCACTGCACTCCAG
CCTCGGGGACAGAGCAAGACTCCATCTCAAAACACACACACACACACACACACA
CACACACACACACAAAACAGATATACACTGAACACAGCACAAGTGGGACATAAGAG
ATTTAAAGGGTTAGAGATGTAAATGGATCTAGGAATGGAAACCATAAGGNGGGA
30 TTTATCAACTGGATTCTGCANAATGCTGTTAAGGCCAGATGTTAGCAGGTGTTACA
TAAAAAAGGGATACCATGAGCAAAAGTATTTGAACATGGGCAATGGTTGAAACAAG
TTTAAACAGATTATNTTTATTACCAAATCTCTCAAACCTTTAATATGCTATAAACA
TTGTGAAACAATAAAAAAACTTTCCAAAA

35 Sequence ID 534

GGGAAGGGAGCTATGAGTGTGTGTGTTGTGTATGGACTCACTCCCAGGTTACCTG
GCCACAGGTGCACCCCTTCCCACACCCTTTACATTCCCCAGAGCCAAGGGAGTTTAA

- 174 -

5 GTTTGCAGTTACAGGCCAGTTCTCCAGCTCTCCATCTTANAGAGACAGGTCACCTT
GCAGGCCTGCTTGCAGGAAATGAATCCAGCAGCCAACTCGAATCCCCCTAGGGCTC
AGGCACTGAGGGCCTGGGGACAGTGGAGCATATGGGTGGGAGACAGATGGAGGGTA
CCCTATTTACAACTGAGTCAGCCAAGCCACTGATGGGAATATACAGATTTAGGTGC
TAAACCGTTTATTTTCCACGGATGAGTCACAATCTGAAGAATCAAACCTTCCATCCT
GAAAATCTATATGTTTCAAACCACTTGCCATCCTGTTAGATTGCCAGTTCCTGGG
ACCAGGCCTCANACTGTGAAAGTA

Sequence ID 560

10 GGCGGAGGTTGCAGTGAGCTGAGATGGCGCCATTGCTCTCCCAGCCTGGGTGACAA
GAGCAAACTCCGTCTCAAAAAAAAAAAAAAAAAAAGCAATTTACTTTAAAAAC
ATACAAACACAGAGACAAGTATTTTGGAGAAACAAATACCTTTTTCATTTTTTATA
CCAATGTAACAATAATCCATTAAACACACCTTTACTAACTGTTTTCTAGGAGTCTG
ATATGATGAGGAAATAGGTAAACCTTTAATAGCCAGTACTAAATTAGAGTGGCACA
15 ACTTTCCTGAGGAAAAAGATGGGTATTTTACTTTTCTGTTTTAGAAAAAGTGGCTT
GACAACAGTATGCTTATGTCTTAGAGTTTGAAATTCAAGTTCTTGAACATTATTAA
TGGCTACAATCATTACATCCACATTGGGCTGTATTCTTGATGAATCCAAAGTGAT
TTTCACCTCAACTCTGAATTTCAATCTCCTCTTTTGAATATAATACAACCATCTCA
CTAGAGGAAGCATTTTCAGTCTTTTCTGATTGGAGATTCAATTATTGTTTTAGATAAT
20 GTTTTCATTTGCTTATGGGTATATAAAAAATTTTATCTTAAAAATATTTCTCTCA
TTTAGCTAGCAACATTGTTTTTC

Sequence ID 561

25 CTCAGGGTGATCTCTGAACCCAACTTGCCCCAAAGAAGGTTGCTCTGTCCTCTCC
ACATCCCCATCTCCTCCCTAGGGCCTTGTTGGGGAGAGGCTCCTCCATCTTTCCCA
AGTCACACCATCGTTTCCTACGTGGTCTGGACAAGAGCAAGAGCACACCTTGTCCC
CACCTTCTCCAGAGCAGCCAGAACCCACCTCAGGTGCCTTCCCCATCCGGTGCAGT
TAAGGCACTTCTGCCAGCACCATGGTATGAGCACTAGACTTGGAGTTAAGATTTGA
GAGCCCCCTCTGTCACTGTGGAAGCTTGAGCATGTTGCTTGATCTCTCTGAACCTT
30 GTGTTTCTCATCTGTGAAAGGTGATAATGTGGGGCTGCTGTGAGATTTAAAGGACA
TAATGCACCTACGGTCCAAGCACTGCCTGGAATACAGCANAAGCTCAACAGATACT
GGACAACCCATCCCCTTAGTAGAGGCACTAACCATGTGACCCAAGGCAAAAGTGCT
TAAAAAAA

35 Sequence ID - 562

nt: 580

ATTGCATGCAAGTTTGCTGAGCTGAAGGAAAAGATTGATCGCCGTTCTGGTAAAAA
GCTGGAAGATGGCCCTAAATTCCTGAAGTCTGGTGATGCTGCCATTGTTGATATGG

- 175 -

5 TTCCTGGCAAGCCCATGTGTGTTGAGAGCTTCTCAGACTATCCACCTTTGGGTCGC
TTTGCTGTTTCGTGATATGAGACAGACAGTTGCGGTGGGTGTCATCAAAGCACTGGA
CAAGAAGGCTGCTGGAGCTGGCAAGGTCACCAAGTCTGCCCAGAAAGCTCAGAAGG
CTAAATGAATATTATCCCTAATACCTGCCACCCCACTCTTAATCAGTGGTGGAAGA
ACGGTCTCAGAACTGTTTGTTCATTGGCCATTTAAGTTTAGTAGTAAAAGACTG
GTTAATGATAACAATGCATCGTAAAACCTTCAGAAGGAAAGGAGAATGTTTTGTGG
ACCACTTTGGTTTTCTTTTTTGGCGTGTGGCAGTTTTAAGTTATTAGTTTTTAAAT
CAGTACTTTTTAATGGAAACAACCTTGACCAAAAATTTGTCACAGAATTTTGAGACC
CATTA AAAAAGTTAAATGAG

10

Sequence ID 563

GCAACCTGCACAACCCCGCCCTGTTTCGAGGGCCGGAGCCCTGCCGTGTGGGAGCTG
GCCGAGGAGTATCTGGACATCGTGCGGGAGCACCCCTGCCCCCTGTCCTACGTCCG
GGCCACCTCTTCAAGCTGTGGCACCACACGCTGCAGGTGCACCAGGAGCTGCGAG
15 AGGAGCTGGCCAAGGTGAANACCCTGGAGGGCATCGCTGCTGTGAGCCAGGAGCTG
AAGCTGCGGTGTCAGGAGGAGATATCCAGGCAGGAGGGAGCGAAGCCCACCGGCGA
CTTGCCCTTCCACTGGATCTGCCAGCCCTACATCCGGCCGGGGCCAGGGAGGGGA
GCAAGGAGAAGGCAGGTGCGCGCAGCAAGCGGGCCCTGGAGGAAGAGGAGGGTGGC
ACGGAGGTCTGTCCAAGAACAAGCAAAAGAAGCAGCTGAGGAACCCCCACAAGAC
20 CTTTCGACCCCTCTCTGAACCAAAATATGCAAAGTGTGACCAGTGTGGAAACCCAAA
GGGCAACAGATGTGTGTTTCAGCCTGTGCCGCGGNTTG

Sequence ID - 564

nt: 671

25 GGAATAGAATTTTAAATAGTAATAACTGCTTGTTTTTTTTGTGCAAGTACTTTTAT
ACATAAGATAAAACAAAACCTTACCACCAAACATACCAAAATGCACCTCTTTTCATA
AGTGAGTTACTAAGATTTCTATACCTGGAATATCATGTATGTTTTCATTTACTGGAT
GTTTACATTTTATAGGAAGGAAAATAGTTTTGTATTATTAAACAACCTGAATACTTATA
AACTGTTGTTTCTGGAAGTTATTTATTCATAAAAAATTTGTTCTTTTGTGCATGAA
TTTATAATTCCTAAATGAAGACCAGAAAGTACAAATTGCTGGGAGGAAGAATAGGC
30 TTTATTAATCAACTGATGTCTTGATTTTTCTAAATGGGAAGATTGCTTTATTTTTTA
ACACTAATTATGGGAGCAGATTCTTAGCAAACCTTCTTTGGAAAAGTTAATGTTATG
ATGTGCATTAGGCTGCCCCATCGTGTATATAAATGAAGCAGATTTGATTTTTGTAT
TCTTACGTTTCTCTGCTTTGTAGTTGTGGCTGTACTTAAAGAAATACAGAATTTCA
TATATTTAAAAATGTTTAAATGTGACCCACAGACATTGTAAATGGATTNAAAAC
35 AACATGAAAAATATTCAACCTAAAAGAATTCTTAACTTCACAAGTGTTTTACTTC

Sequence ID 565

- 176 -

CTTGGTTCCGCGTTCCCTGCACAAAATGCCCGGCGAAGCCACAGAAACCGTCCCTG
CTACAGAGCAGGAGTTGCCGCAGCCCCAGGCTGAGACAGGGTCTGGAACAGAATCT
GACAGTGATGAATCAGTACCAGAGCTTGAAGAACAGGATTCACCCAGGCAACCAC
ACAACAAGCCCAGCTGGCGGCAGCAGCTGAAATCGATGAAGAACCAGTCAGTAAAG
5 CAAAACAGAGTCGGAGTGAAAAGAAGGCACGGAAGGCTATGTCCAAACTGGGTCTT
CGGCAGGTTACAGGAGTTACTAGAGTCACTATCCGGAAATCTAAGAATATCCTCTT
TGTCATCACAAAACCAGATGTCTACAAGAGCCCTGCTTCAGATACTTACATAGTTT
TTGGGGAAGCCAAGATCGAAGATTTATCCCAGCAAGCACAACTAGCAGCTGCTGAG
AAATTCAAAGTTCAAGGTGAAGCTGTCTCAAACATTCAAGAAAACACACAGACTCC
10 AACTGTACAAGAGGAGAGTGAAGAGGAAGAGGTTCGATGAAACAGGTGTAGAAGTTA
AGGACATAGAATTTGGTCATTGTTCACAAAGCAAATGTGTCTGAGAGCA

Sequence ID 566

GTCACCAAGAGCTTGTGTGTCAGGTTTTCACTTGCTATTTCGCAGAGATTTTTTTTAA
15 AGGCACTATTTGTAGTGTTAAAAGGGTGAATTTATCANAAGGCATAATAATCATAA
ATGTGTATATGCCTAATAATAGAACTTTAAAAGGCATGAAGCAACACTCAAAGGA
TTAAAGGGAGATCATCTCACCCCCTTCTTACCAATTGATAGAATGATCTGATGAAA
ACAGTAAAATAACAACAGATCTGAACACTGTCAACCATCTTGACAAATACTTATGC
CTAGTGTTCCATTATTGGAACACTAAACATGTGGAATGATTTATATCCTACTGCTC
20 AAGGTCATCACCAAGGTCTAATTGTAAAATTTCAAAAAATTGCAACCTCAGGCATA
AATGGGTTAATCGACATTTATAGCACACACATGCAACATGTACCAGAGATTCTTTC
TTTTCTATGAACATGGTACTTCCACCAAGATAGACCACATTGTGAACTATAAAACA
AATCTAAAACATTTGAAATGAAGGAAATTATATAAAATATGTTCTCTTGATCTCA
ATGAAATTAAATTAATACTATAT

25

Sequence ID 567

CTCATGGCGGCCAATGTAGGCCCAAACTTCCTCAAGTCAAACCTCTCCAGGCCAC
CTTCTGCTTCCCGGTGGCATCAACAGGCCCAGCTTTGACTTGAGAACAGCCTCTGC
AGGCCCTGCTCTTGCCCTCCAGGGGCTTTTTCCAGGCCCAGCTCTTGCCCTCATGGC
30 AGCTGCCCCAGGCCAAATTTCTGCCTGCCTGCCAGCAGCCTCAACAGGCACAGCTC
CTCCCTCACAGTGGCCCATTTAGGCCCAACTCATGACTGTGAGGCCATTTCCAGGC
CTAGTGCCCTGCCTCGTGGCTGACTCTTGAAGCCCAAACTTCCTCAAATCAGCCTT
TTGCCCAACTTCTGTCTACTGTCGGACTCTACAGGTCAGCCTCTGCCTCACAGTGG
ACCTTCCAGACCCAGATGGTGTCTNCTGTGGCATCCTCAGGCGAAGCTCCTGCCTT
35 TCGGCAGCCTCTCCAGGCCCAGCTCCTCCTGCTCCAGCCTTCTCTCCAGGCTCTGA
ACTTTCTCAGGTCTCCCTCTGTTGTCCAAGGCTGGAGTGTAGTAG

- 177 -

Sequence ID 568

TATATATGTAATGCCCTTAACCTAGTGTGTTGGCATGATCGTTGCTGAAAGGGAAGC
TTGTGGGTACAGTGTCCCCTCAGAAGCCAAAGCCCAGGGAAGGTCGCCTGCCCAGG
TCAGGCTCCCAGCGAGTTTGTCTGGGGAGGGG¹CATTTCATACCTCCAGGTCAGGAC
5 AGAGGCTCGGGCTGAGGGAACCTACACAGGTCCTGGAAGCAGATCCTTCCTGCCT
AAGCCAGCAGGACAGCTCAACAGGAAGCATCTTCAGCCACGGGAGGAGAGGCAGC
ACCTTTTTTTGGAACCATAACAGAGCTAAGAATGGTGGTACAAGTAATAGATTCTGTA
CTGGCAACCCCACTTGGTGGAGCAAGTTCTAGGAAAAGGGGGCTGTCCTTGAGTCA
GCCATGGGGTCAGCCACACAGTCACCGCAGCTGCTCTTTGGCACCGGGCGCTGGAA
10 AGACCTAGGATGACACAGCCTGGAAAGAGCTTGGGAAAAGCTCATCTTCACAGAA
CTACCTGCTATAACCAGCCAGGGCAGGTGCTTATTTCCACAACAGCCCTCTGTTGTA
GGCGGCAGTGCCATCCTGAANGTGCCGTGGTACCTTCTGAANACCCAGCTGAGGGC
CTGTAATGGCACTTGCATGCCACATGGNACACCCTTTCCCGGTAA

15 Sequence ID 570

ACCGCGGCCGCGTNAANAAAAAAAAAAAAAAAAAGAATTCCACTTGATCAACTTAATTCC
TTNTCTTTATCTTCCCTCCCTCACTTCCCTTTTCTCCCACCCTCTTTTCCAAGCTG
TTTCGCTTTGCAATATATTACTGGTAATGAGTTGCAGGATAATGCAGTCATAACTT
GTTTTCTCCTAAGTATTTGAGTTCAAACTCCTGTATCTAAAGAAATACGGTTGGG
20 GTCATTAATAAAGAAAATCTTTCTATCTTAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AA
AAA

Sequence ID - 571

nt: 457

25 TTAGAGAGGTGAGGATCTGGTATTTCTGGACTAAATTCCCCTTGGGGAAGACGAA
GGGATGCTGCAGTTCCAAAAGAGAAGGACTCTTCCAGAGTCATCTACCTGAGTCCC
AAAGCTCCCTGTCCTGAAAGCCACAGACAATATGGTCCCAAATGACTGACTGCACC
TTCTGTGCCTCAGCCGTTCTTGACATCAAGAATCTTCTGTTCCACATCCACACAGC
CAATACAATTAGTCAAACCACTGTTATTAACAGATGTAGCAACATGAGAAACGCTT
30 ATGTTACAGGTTACATGAGAGCAATCATGTAAGTCTATATGACTTCAGAAATGTTA
AAATAGACTAACCTCTAACAACAAATTAAAAGTGATTGTTTCAAGGTGATGCAATT
ATTGATGACCTATTTTATTTTCTATAATGATCATATATTACCTTTGTAATAAAC
ATTTTTCCC

35 Sequence ID 572

CGTCTATTTGNGTTTCTTCTCACAATTGGTAAGTTCTCTGTATTGATTGATGGCTA

- 178 -

AGTTTGATTAGTGTTTTCTCTAGTTGGTAATTATATTCTAGTATTTTATCATCTT
ATTGTTTACTCAACTNAAAGTGNACAGAAGAGTTGCCAGGTTTCTCTTTGATATG
AGATCTCTNNTTGATTGGAATGCAAATCANAAGTGTCATGTTTTGAATAAAGGGA
CCAGATGACTTATAGGTATTCTTTCTCTAAATATAACTAAGGTAAGATTTTTGTTT
5 TGAGGTACTTAATCTATATAAGTGGTAAAGAATTTACTTGAATTTCTCCAAATTCT
CATGTCTAAAGTCTGATTGATTAAATTCATTCTTGGTATTTCAATTTGAAAAGAAT
GTAGCTTTAGCAAACCTCTTTGTATAAATGCAGTGGGATTAAGGTCATTTAAAAAA
TTGTTATATCATTGTATTTTTTAAATTTACCAGTTTTATTTTTCTTTTTACCCTTT
AGCCCGGCCTCAGAAAGTGTGTTTGTGTCCATTTCTCCAGCGCACCTCTGCATA
10 TCTCTACCCACTTGTCTATAATTCAGCATCCAGCAGAGGAAAACAAAGTGTTGCGTA
CAGTTCCTCTACTAGCAGCATGCCTCCCCCAGGACAAGTGTA

Sequence ID 574

TTATTGCTGACATAAAAATGGTGCACATCGGCCAGGGCCCAGGATGAATCAGCCAA
15 TCTGCACCATTTATACATGGAAGTGGAGAACATTGTGCCAATAATCATTTAATATA
TGCCAAATCTTACACGTCTACTCTAAACTGCTCTAATGAAGTTTCAGTGACCTTGA
GGGCTAAAGATTGTTCTTCTGGGTAAAGAGCTCTTGGGCTGGTTTTTANAGCAGAG
TTCTTGTGTGGGTAGACTGTGACTAGGTTTACAGCCTTTGTGGAACATTCCGTAT
AACGGCATTGTGGAAGCAATAACTAGTTCCTATGAAAGAACCAGAGCTGGGAAGAT
20 GGCTGGGAAGCCAGGCCAAAGTGGGGGCAACAGCTTGCTTCTCTTTCTCTCTCAC
CCTCAGTTTGTATGGGAAAATGGAGATGTCCTCTCCACTTTATCCCACGATATCTA
AATG

Sequence ID - 575

nt: 209

25 CAGGATATCGAGACCATCCCAGACAGCATGGTGAAACTCCGTCTCTACTGGAATAC
AAAAAGTTAGCCGTGTGTGGTGGCACGCGCTCTAATCCCAGCTATTCGGGAGGCT
TAGGCAGGAGAATTACTTGAACCCGGGAGGCGAAGGTTGCAGTGAGCTGAGATCGC
ACCATTGCACTCCACCCTGG-CGACAGAGCAAGACTCCGTCT

30 Sequence ID - 576

nt: 541

CAGCCAACCCAGAAGGAGCCAGTCTACAACTATGCCTGATCCTCCTCATGGCAGGC
CACGAAGCATTGCTGCCATGTGTTGAATTATAAAACCCACATTGCTTTTTGAACCC
TGTTGCGGGTAAAAATAACCAAATTATCAGTCCTTGGAACCCAGGCAATCAAGTG
AGTACAAGGTAAAGATAAGTATGGTTTAGAGGAGAAATTATGTTCTGAACTGGTG
35 TCCTTTGATGGCAGCGTCAGCCTTGCTAAGTCAGAGTAGAGGGAGCAGTGACCTTA
ATAAGCTTTGGTGAGCATCATGTGCACGCGTGGGTGGGAGTCCCTTTCACTGATGC
TTTTAAAAGTGCTTTTGCAGACCCTGGAAGGGATCCTCCACACATATGAGGTGTGG

- 179 -

GACAGGTAGGCCAGAGAGGATTAGCCCTGCTTTTCGAGACTAGAAATCTACAGTCCT
GAAGGAGCAGTAATTAATTGGTACACCTGTCAGGGCCAGCCCCAGGTCTCCTGGC
TTTTTCCAGGTTTTCTGTCTCACATGATTTTGCTTTT

5 Sequence ID 577

CTTTAATTTTTCAAGTGTTTAAAAACAATTTTATACTTAAGCCAGCCTTGAAGAT
AAGCACAAAATTTACCAGTTTACATTTAAAAACAAACAAAAACGACAACAACCTC
AAGCACCCGCTCTGTGCATAGCACTATTCTAGGTGCAATAAAAGGGAATCTTAACC
TTAGAAATATGAGTTCACTTTCTGGAATTGTATTATCTCCTTTTCCAGAGAGTAAA
10 AATAAATAAAATCACCATTGTTTACTACAGATCTGCCCCAAACCACATCTGGTTCA
CAGAAAGGCTAATTTCTGCCAAATTAAAGATGTAATGAACTCAGTTCCTGCTTTCC
CAAAAACACGAAAGCAGAATTCCTTTTCTACTGAAAAAATAAACAGTTTTCCATGC
AAGGGCAGTTTGCTTCTAATAAGTATTTTTTAAAAAATTTTTTTTCTCTAGCTT
TTCTTTAAATTTTCTTCTCTAATATTGCCTTTTCTTGTACAAGGCAGACCAGGTA
15 TCTTTTATGCTGTTTTTCTTTTACTAAGAAAGTATTGCATCTTGAAGACAAACC
ATTTCCAGAGTAGTGATAAAAAATAACACTAAAAAACTTTTAAAGGTGAGTCACT
TCATCACCTTGATGAAGTAAAAAA

Sequence ID 578

20 GGAAAAATATTTCCACTTAGATATTTTACATGGTTTTGTTTAAAATTACCATTAC
TTGTTTTTTTAAAAACACATGACCACATATGTATATGTATATCTACCTAAACATTGT
ATCATGGTTTTCAGTATGTTATTCTATGTATTACTGGGAGATGCTACCAAGAAACCAA
CCCAAAGAAAATTCTGGAAAATACATTTCTATTTATAGAATAAATGTTTTCATTTAT
ATAAAAGCAAAGAAGACTTAGAGTTCTAATAAATGGGATGTCTAATAAATTATGAAG
25 TTAGTGATTTGAATATATTTATATTTTTATAACTTCCTTGCCAAAGTCCTGATTTAG
TACATTAGAGAACCTGTGTTTTCTCTCTCTCTACCATTCTCTCTCTCCATACA
GTCATTTGGGCTTTTTTACTCAAAGAGAATCAAGAAATAATAAGGTATAACAAGCTT
GGCAAAGTGTTGGCTTTTTTAAAAAAAATTTTTTTTAAATCTCTAGCAGTTTGGTAAT
TTAGCAGCATCATTTATTTGGGATTCTTTTATCTGATTTCAACAGTGAAAAACATC
30 CCTATGATAAAGCCTAATGACCCATTTCCAAAAGATGGAATTGCCCTTCCTAGAAA
ATATGACGGAGAAAAGT

Sequence ID - 579

nt: 502

CGAATAGCCAAGTGGTCTGACAAGATCGAGAGTAATGAGGCCCATACTTTAGTACA
35 GTCTTGAATGGCCAGATGGTGCTGGGCATACCCCAACCAGAGATATGTAAGTCTTT
ATGTTGTCAAATTTCCCAGAAACATGAATTTCCCACTAAGATTCATTAAGGAAAA
CTAGAATGAAAACAAAACGTTCTTGTATAATATTCATTANAAAGAAATGAAGAA

- 180 -

GGCCGGGCATGGTGGCTCACGCCTGTAATCCAGCACTTTGAGAGGCCAAGGTAGG
CAGATCATGAGGTCAGGAGTTTGAGACCAGCCTGGCCAACATAGTGAAATCCCGTC
TCTACCAAAAATACAAAAAATTAGCCGGGCATGGTGGCACACACCTGTCATCCCA
GCTACTCAGGAGGCTGAGGCAGGAGAATTGCTTGAACCTGGGAGGTGGAGGTTGCA
5 GTGAGCTGAGATTGCACCACTGTACTACAGCCTAGGTGACAGTGCAAGACTCTG

Sequence ID - 580

nt: 316

CCTATGCCAAACTAAAGAAAGCTTGCCTGGCCTACAGGCCTAAAGGTTCAAATGNG
GATTAAAAAACACAGTAGTCACATAAAATGTCTGCTGGCTGGCTGGAATTCCATC
10 ACCTACAATTTAÇCTGCTTTCAAAACTGTGTTCAACATTGAGAAAACAGAAAACC
ACTTATCTTGAGCTTAATATGGGCTTCTTTTTCTTAACTGTAGAACACTTACTGA
AATATCAAATCAATGGTTAGGATATGTATCCTAGGCAGGCCTAAACCATTAACT
TGGTTTAAGCAACTTTGTATAATTNACCTCCTAAAT

15 Sequence ID 581

CTTCATGAGTGCCCGGTTGCCCAAGTCAAAAACCTGGGAGTGATATAAACTCCCCA
CACATCCAGTCAGTCACTCATCAACTCTATTGATTCTG-CTGCTAAATATATCTCA
ATTGTATTAACCTAAACATATGCATAATACATCTTCTTCTTCACTGCATTTTTGTG
GGCTGCACTTACCTTTCAGGTAACAACAACACTGGCCCCCTCTTGCCCTTCTAGTCA
20 GAAGTGCCAAAATGATGAGAGCTAGCCATGACAAACCCACAGCCAACATTACACTG
AATGTGCAAAACTGGAAGGGCATCCAAACAGAGGAGG

Sequence ID 582

TAGAATTCTCGCCTGCCTTGGCTTCTCCCTCTAGTTGTTCCCTTCTCTGTCTTCTGT
25 GGGCTTCTTATTGTCTGCTCACTCCTTCTTCACTGTCTCTCATGGGCTTCCCTTCC
CTTCTCAGCTGATGCCATCACCTGGGGAATCACAGTTACTCAGCAGCACTGGGGCC
TCTCTATCTCTATGCTGGTCATGCCTATGTGTGAGCTGCAGACCCAGTGGAATTTT
CATTTGTGCATCCCATGCCCAGCCCACCCTCCACCAGCCTCGAATGCAGCTGTTCA
GCCCTACCCCAGTCCTCAGAAAAGTTCTCTCCCTGGATCCTCTTTTTCTTTCATG
30 AGTGCCCGGTTGCCCAAGTCAAAAACCTGGGAGTGATATAAACTCCCCACACATCC
AGTCAGTCACTCATCAACTCTATTGATTCTGTCTGCTAAATATATCTCAATTGTAT
TAACTTAAACATATGCATAATACATCTTCTTCTTCACTGCATTTTTGTGGGCTGCA
CTTACCTTTCAGGTAACAACAACACTGGCCCCCTCTTGCCCTTCTAGTCAGAAGTGC
CAAAATGATGAGAGCTAGCCATGACAAACCCACAGCCAACATTACACTGAATGTGC
35 AAAACTGGAAGGGCATCCAAACAGAGGA

Sequence ID - 583

nt: 631

- 181 -

CTGAGGTGGGAGGATTCCACTCTCACCCATTTCTTCTTTCATTTTCAGTTTCTCCA
GTTAGTAACTGAAGATGTTCTTTGAGTAATTAAGTGAGTGAGAAAATTTTAAAGTG
AGAAATCTATAAAAAGAACCATGTAAACATAAATATTTTCAGTCCTTACAAGTTGGT
ATTGACTTTTCTCATTGGTAATCTGACTGATTTAATACTGCTCATTCCAATATCTG
5 GTGATGTAATTCTGGTTATGAATCCTTGTATTAATAACACCTCCTGGGAGGTTTTT
TTTCCCCAACATTACATTCAGAATATTAGAGCTGAAAATACCTTTTTTAAAGGTTAT
CAGGAGGAGGGAGCTTATGTTTAAATGTGGTGGATAAACTTAACTGCTGGTTAATA
CAATTGTTATTCAGGTGAAATTCCCTAAACTTTTCACGTGCAAAGTTTGTATGTA
TACAGACATTTGGGGAAAAGTTTTATCATCCCTAAAACCGGTTACTGTCCAGAAAA
10 TGATAAGAATCCCTGGGTTCCAAATCCTTCATAAGGTATTTATTCAATTTATTTATT
CAACACATTTACTCAATGCCTCCGCTCTGCTGCAACTACACTGACATTCTGCTTCT
AATCTAACCGAAAAT

Sequence ID 585

15 TTTCAAATTGTACAATAACACAAACAACCTTTGTTAAGGCCATGTTTTATTTGCTGA
TTAATGGACAAAAGGCAATGTAATTTATTTTCAAGTATTTTCTTGAAAGTCTGTGC
TCATAAAAATCATGAAAAGTTGGAAAGACTGTTAAATCACTGAAACTTCAAATATA
TCTTACACAATCTTGTTTGTACAAAATAACAAGTTAAATATAAACATAAAGCAATC
ATGGTAATTTTATGCAAATCTGTTTATGTGATCATCAGTTATATATAAAAGTTTC
20 TCAGTTCTGTTATTTGTGAAAAGATCAATACCAGATTGAATGACTACCTATTGGCA
AAGGGCCCTAAAAGCTTACTTTAGCACTCATCTTTTACATGGTTAAATGCATTTTC
CTAATTTGAGATCACCTAAACACTGGAAAAGAAAAAAATGAAAGGGCAGTATGTC
CATAAACCAACAAATAATTTGGCTGTAATGTATCATAAAACACAAACCCACACAT
CTGTACAATAAACATTATGTATTACATACACACACACACCCAGTCATAAAGCC
25 TAATGATGTGCTGCTTCCAGTTCAATATTAGCTGTGCATTTTTTCTTATTTTCATC
AAATGAATAGCTTTTTGTACC

Sequence ID 586

GTAAACTGTTCTCTCCGAGGGAAAAAATGGAAGTTATCCTCACAGTTCACTGCCGT
30 GGTATTTCTTCTGTCCCATGCTTTGCATGACTGCCATGGTACAGCCTTGTTTCAA
CTGTTCACTGTGATCTGTGGGTCTTTGAGTTTCAGTGAGTTTGCTGAAATGTCGAA
GAAGTAGTTCCAACTTCAATGTTCAATGAAATTTTTGTTCAAGTTTGAAATGGAG
AGAGCAGCTTTAAAAGGTACTAAGCCTTTTACAAATTGGTGAGTACTGGCACATGA
GAT

35

Sequence ID 587

TTTTTTTTTTTTCCTTAAAAGGTAACCCCTAAACACAGCTAAAACCTATGCCATCAGC - - - - -

- 182 -

TGACTCCAAGGNACACACAGTCCTGTATCTGGAAC TACTGAGTGGCAGGCATCTTT
CTCTGCCTCTGACAGTGGAGTCCCCATCACTGCAGAGCATAGCCAAAGGAGTCAAA
GGTCTCAGCGGGTCACTGCCTTATCAACCCTCACCAGTCCCTTATGTTTTTTAATA
TTTTATAATCTTGACATGACACCAAGATGCTTTAATAAAAAAGCACCTCTAACTCG
5 GTCTTGATTCACTTACCTTGAGCCTGGGACTTCTCTAGGCTCCTGAGGCAAAAAC
AGGTAGAGGGGAGATGGTGGAACATAAAACACAATTTTGCTTGGCACCCACCTTGG
CGTCTGTCCCCATGACCAGGTCTTCAATTGATGATTTTGTCAATTGATGGAGGAG
CGATATCGTTTCTCAATGATATTATGGGTTGTCCGCCTTCTCCTTCTTTGGGGGG
CTCAAGCTGCTTGACTCCCCCAGGTACCTGCTTAATGGGGCACTTTCTCTTGCCCC
10 ATCATTACAGGCATTGTGGTCAGAATGGTCCCCTGCTGCCCACCAGGGTCTA

Sequence ID 588

CTAGTCTTTTCATAGTCTGCATAGAGTCTGGCCATTACCATCAGTTTTTAAGATGT
CCATATTGTGGCCGGCGCGGTGGCTCACGCCTGGTAGTCCCAGCACTTTGGGAGG
15 CTGAGGCAGGTGGATCATGAGGTGAGGAGATCGAGACCATCCTGGCTAACACGGTG
AAACCCGTCTCTACTAAAAAAATATTAATAAAATTTGGCCAGGCCTGGTGGTGGGCG
CCTGTGGTCCCGGCTGCTTGGGAGGCTGAGGCAGGANAATGGTGTGAACCCGGAAG
TCGGAGGTTGCAGTGAGCCAAGATTGCACCTGGGCAACACAGCGAGACTCCGTCTC
AAAAAAAAAAAAAA

20

Sequence ID 589

CAATTATTTATTACCTTTCCATTTGTTTCGCCTGATGATGTGACAATGCATGGTCTT
TGTGCATGCTGCTAGACACTTTTCTTTCCCAGCCGAAAAGTCTATTATGTAATTTT
TACATTATAATTTTAATGTGGATGATCAGGATTAAATCAAGATATATATCTGGAA
25 CCTCTTATAAATGGAGCACTTAGAAATTTGTTGTTCTGCACTTAACCTAGAGAGAG
AAAAAATGCTTTTCTTTGTGAAAAATCTGAATTCCTGTCCTGACCTTCTGTGATGT
GGAAACCCTAGGCTCTGAGACACACTCTCTGGTGTCTGAGACAGAACCAAGCAAT
AACGTTGTGATGCCCCACAGGCCTGGAGCCAGCTAGCGACCTTGTGCCGCCAGCTG
TCCATGGCCCGTGCAGAGCAGAGGACAGTGAGTGTCTGCACTGAGAACCCTTAAACC
30 ACAGTTGAACATACCCACACCTGTTTGTCTTAAGCTATAGTGTAACCAAGTTT
GGGCTCTGAAAATTTAACTGAAAAAGATTTCTTCTTGT

Sequence ID 590

GTGGCAGCAGGCGCAGCCCAGCCTCGAAATGCAGAACGACGCCGGCGAGTTCGTGG
35 ACCTGTACGTGCCGCGGAAATGCTCCGCTAGCAATCGCATCATCGGTGCCAAGGAC
CACGCATCCATCCAGATGAACGTGGCCGAGGTTGACAAGGTCACAGGCAGGTTTAA
TGGCCAGTTTAAACTTATGCTATCTGCGGGGCCATTCTAGGATGGGTGAGTCAG

- 183 -

ATGATTCCATTCTCCGATTGGCCAAGGCCGATGGCATCGTCTCAAAGAACTTTTGA
CTGGAGAGAATCACAGATGTGGAATATTTGTCATAAATAAATAATGAAAACCTAAA

Sequence ID 591

5 CAGCAGCAGAAATGTTTGCAAGATAGGCCAAAATGAGTACAAAAGGTCTGTCTTCC
ATCAGACCCAGTGATGCTGCGACTCACACGCTTCAATTCAAGACCTGACCGCTAGT
AGGGAGGTTTATTTCANATCGCTGGCAGCCTCGGCTGAGCAGATGCACAGAGGGGAT
CACTGTGCAGTGGGACCACCCTCACTGGCCTTCTGCAGCAGGGTTCTGGGATGTTT
TCAGTGGTCAAAATACTCTGTTTAGAGCAAGGGCTCAGAAAACAGAAATACTGTCA
10 TGGAGGTGCTGAACACAGGGAAGGTCTGGTACATATTGGAAATTATGAGCAGAACA
AATACTCAACTAAATGCACAAAGTATAAAGTGTAGCCATGT

Sequence ID 592

TACTCAATGAAAAACCATGATAATTCTTTGTATATAAAATAAACATTTGAAAAAAA
15 AAAAAAA

Sequence ID - 593

nt: 565

CAGGATCAAGGTGAAAAGGAGAACCCCATGCGGGAACTTCGCATCCGCAAACCTCTG
TCTCAACATCTGTGTTGGGGAGAGTGGAGACAGACTGACGCGAGCAGCCAAGGTGT
20 TGGAGCAGCTCACAGGGCAGACCCCTGTGTTTTCCAAAGCTAGATACACTGTCAGA
TCCTTTGGCATCCGGAGAAATGAAAAGATTGCTGTCCACTGCACAGTTCGAGGGGC
CAAGGCAGAAGAAATCTTGGAGAAGGGTCTAAAGGTGCGGGAGTATGAGTTAAGAA
AAAACAACCTTCTCAGATACTGGAACTTTGGTTTTTGGGATCCAGGAACACATCGAT
CTGGGTATCAAATATGACCCAAGCATTGGTATCTACGGCCTGGACTTCTATGTGGT
25 GCTGGGTAGGCCAGGTTTCAGCATCGCAGACAAGAAGCGCAGGACAGGCTGCATTG
GGGCCAAACACAGAATCAGCAAAGAGGAGGCCATGCGCTGGTTCCAGCAGAAGTAT
GATGGGATCATCCTTCCTGGCAAATAAATTCCCGTTTCTATCCAAAAGAGCAATAA
AAAGT

30 Sequence ID 594

CAGAAGAGTAAGCAAATCTCAAAGCAGCGAAAGGGAAGAACTAAAAAAGGTAGAG
CAGAAATAAGAGAAAATAGAGAAGAGAACAAATTGAGAAAAATAATTGAAACCAAAA
GGTGGTTCTTTGAAAAGCCTAACAAAATGGACACATCTTTAGTTAGAGTGACCAAG
AAAAAAGGGCAGTGACTCAGATTACTTCATTCAAGAGTGAAAGAGGGGCACATCACT
35 ACCAATTTACAGAAATAAAAAGGATTATGAGGAAATACTACAGATAATTGATGACA
TTAACTTAGAAGAATATATTTCAAGAAAGACACAACTACTGAAACCGACTCAAGA
AGAAACAGAAAATCTGAACAGACCTATAAAAAATAGAGATTTAATTGATATTCAGA

- 184 -

AAGTTTCCCAAAAAGAAAAGCACTGGCCAAGATGACTTCACTGGTGAATTCTATCA
AGTGTCAAAGATGAATTACTGACATTCATTACACTCCTTTAAGAAATAGAAGAGG
GGACATCACTTTTCAAAGCATCGACATTCTAATCATTAGTCCCTTGGTTTCCTGCT
CCCAAAGCCAGGTGATGTATCACAAAAAACCCCTACAGACCCACTGGGCACAATG
5 GCTTTATGCCTAT

Sequence ID - 595

nt: 98

CTTTGCTCGAATNGTCAGATAAGGATTCTGTGAANGGAGATGAGATTTCCATCCAT
GCTGACTTTGANAATACATGTTCCCGAATTGGGGNCCCCAAA
10

Sequence ID 596

CTCAAGTGTTCCCTCAGCTTAGGCTTTGTTTAAATGATCCACCCAGGGGCGATGG
TAGGGAACAACAGGGTCACTAACTATTTGGCTGGCTACAACTCTGGGAAATGGTA
AGACAGGGAAAGGCCATGTTGTTTCATTCCCTTGTGCAGATCTAGGGAGAACCAGCAG
15 AGAGAACAGTTAGCATTTCTTGTTCATGAATTATCCTATTAAGAACACTGGATGT

Sequence ID 597

CGGNCGCGGTGACGCTACTCCTACCTATCTCCCCTTTTATACTAATAATCTTATA
AAAAAAAAAAAAAAAAAAAAAAAAAAAA
20

Sequence ID - 598

nt: 362

GGCATGTGCCTGTAGTCCTAGTTGCTGAGGTAAGAGGATTGCTTGAGCCCAAGAGT
TCAAGGCTGCAACAAGCTTTGATTGCGCCACTGCACTCCANCCTTGGCGACAGACT
AAAACGCTGTCTCAAAAAAAAAAACAAAAACGACNAAAAAAAAACAAAACAGAAAAA
25 ATTAACCTTAGGCAATGACAGTCCCTGGCAAATGCTGGGAGGGAGGCAACANTGGTC
AAGGAAGGTAACCCTGAANCAGGACTTGTAAGCAAATAANATTGGGAGGCCAAGG
TGGGTGGATCACNAGGTCAGGAGTTCGAGACCAACCTGGCCAACATAGTGAAACCC
CGTCTTTCTAAAAATACAAAAAATT

30 Sequence ID 599

GACAAAAGAACCATTTGGATACATAGGTATGGTCTGAGCTATGATATCAATTGGCT
TCCTAGGGTTTATCGTGTGAGCACACCATATATTTACAGTAGGAATAGACGTAGAC
ACACGAGCATATTTACCTCCGCTACCATAATCATCGCTATCCCCACCGGCGTCAA
AGTATTTAGCTGACTCGCCACACTCCACGGAAGCAATATGAAATGATCTGCTGCAG
35 TGCTCTGAGCCCTAGGATTCATCTTTCTTTTCACCGTAGGTGGCCTGACTGGCATT
GTATTAGCAAACCTCATCACTAGACATCGTACTACACGACACGTACTACGTTGTAGC
TCACTTCCACTATGTCCTATCAATAGGAGCTGTATTTGCCATCATAGGAGGCTTCA

- 185 -

TTCAGTGAATTTCCCCTATTCTCAGGCTACACCCTAGACCAAACCTACGCCAAAATC
CATTTCACTATCATATTCATCGGCGTAAATCTAACTTTCTTCCCACAACACTTTCT
CGGCCTGTCCGGAATGCCCCGACGTTACTCGGACTACCCCGATGCATACACCACAT
GAAACATCCTATCATCTGGAG

5

Sequence ID - 600

nt: 595

TTCAAATTCCTTGNTAANAGTCTTTGTTCTGAATTTTACTTTGTCTGTTATTCCTAT
AGCCTTTCCAATTTTCTTTTCGCTTGGAATTTTACGTGATAAGTTTTTCCCCCATTT
TACTTTTANCAACTCTATATTTTTTAGTTGAGGTTGGGTTTCTTGTAACAGCATA
10 TAATTTGGGTTTTTTAATCCAATCTGAAAATTAATGTCTTAATTTTGTGTTTATA
CCATTTACACATAATGTACTCATATATAAGGTTTAACTGAAACCTACTATCTTGCT
AGTTGTGCTCTACTTGAATTTTTTTTTTAGTATTCTGTTTTAATTGACCAACATTTG
ACTGTATCTCTTTGTGTAATCTTTTACAGGTTGCTGTAGGCATGACAATATATAC
ACTTAACTTTTCTCAGTACACTGAGAGTTGAAATTGTAGTACTTCGAGGAAAACAT
15 AGAAAACCTTGCAATGATATCGGTTACATTTTACCACCTCCATATGTTGCAATTATT
AAATGTATTAGATCTGCCTACCTCGAAAACCCATCAGTCTTTTAACTTTGCTCTCA
ATGGTGATTATATTTTTTAAAAAACTTGAGGCAA

Sequence ID - 601

nt: 522

20 TCGACCGGGTTTGGAGCAGTGCCTTGTTTGCTGTGCAGCGGATACTCTACAGGTAC
ATTTCCTTTTTGGAAACCAAAGGGAGGGATTGACAATATTGATGGTAGATCTTTT
TTCTTTAGCAAGAATTAAGGATTTTGGTGGGTGGGGGAGGCTTCTGTGGGGACCA
AGACAATGTACTGTCACTCAGGATTTAAGTCGAACTACCTCATCCCTTGCCCCAGA
GAACAGTTGATCGTGTTTTAAACCAAAGGTGCGGAATGGAGAGAGGGAGGCGGTG
25 CATTCAGCTTCCGATAGAGCTTTTTATTTTTGGATATCAGGAACCAATTTTGAAG
ATTTCTTAAGAAAGTCATTTACATCAGGGACATGAAGAGCAAAGTAGGTATTTTTG
GTCAGTACTTGAATTTGATAGGCTTTATGCAAACAACCTCTCCCTCTGCTGGAGTCT
GGCAAGTTTGCTTTTCACTGGACGCTAATTCAAGTGCCATACAAAACCTAAAATAAN
AGTTTTACTTATAACACA

30

Sequence ID 602

CAGAAATCGCAATTGAAGACCAGATTTGTCAAGGTTTGAAACTGACATTTGATACT
ACCTTCTACCAAACACAGGAAAGAAAAGTGGTAAAATCAAGTCTTCTTACAAGAG
GGAGTGTATAAACCTTGGTTGTGATGTTGACTTTGATTTTGCTGGACCTGCAATCC
35 ATGGTTCAGCTGTCTTTGGTTATGAGGGCTGGCTTGCTGGCTACCAGATGACCTTT
GACAGTGCCAAATCAAAGCTGACAAGGAATAACTTTGCAGTGGGCTACAGGACTGG
GGACTTCCAGCTACACACTAATGTCAATGATGGGACAGAATTTGGAGGATCAATTT

- 186 -

ATCAGAAAGTTTGTGAAGATCTTGACACTTCAGTAAACCTTGCTTGGACATCAGGT
ACCAACTGCACTCGTTTTTGGCATTCGAGCTAAATATCAGTTGGATCCCACTGCTTC
CATTTCTGCAAAAGTCAACAACTCTAGCTTAATTGGAGTAGGCTATACTCAGACTC
TGAGGCCTGGTGTGAAGCTTACACTCTCTGCTCTGGTAGATGGGAAGAGCATTAAT
5 GCTGGAGGCCACAAGGTTGGGCTCG

Sequence ID - 603

nt: 624

GACACACGAGCATATTTACCTCCGCTACCATAATCATCGCTATCCCCACCGGCGT
CAAAGTATTTAGCTGACTCGCCACACTCCACGGAAGCAATATGAAATGATCTGCTG
10 CAGTGCTCTGAGCCCTAGGATTCATCTTTCTTTTCACCGTAGGTGGCCTGACTGGC
ATTGTATTAGCAAACCTCATCACTAGACATCGTACTACACGACACGTACTACGTTGT
AGCCCACTTCCACTATGTCCTATCAATAGGAGCTGTATTTGCCATCATAGGAGGCT
TCATTTCACTGATTTCCCCCTATTCTCAGGCTACACCCTAGACCAAACCTACGCCAAA
ATCCATTTCACTATCATATTCATCGGCGTAAATCTAACTTTCTTCCCACAACACTT
15 TCTCGGCCTATCCGGAATGCCCCGACGTTACTCGGACTACCCCGATGCATACACCA
CATGAAACATCCTATCATCTGTAGGCTCATTCAATTTCTCTAACAGCAGTAATATTA
ATAATTTTCATGATTTGAGAAGCCTTCGCTTCGAAGCGAAAAGTCCTAATAGTAGA
AGAACCCTCCATAAACCTGGAGTGACTATATGGATGCCCCCACCCTACCACACAT
TCGAAGAA

20

Sequence ID - 605

nt: 338

ACCTGAGGCCTCGGTGGGGCCAGTGCGACGCTGGCTTAAGGAGCTGGAGGGGTTCC
TAATACACATTTAATTCAGTTTCTCTTCCCTAAGAGGCTGCCGGAGTTGGGGCCTC
CTCCAGCAGAGACCCTCGGACCCCTGCAGGGCCTGGACTTGGGGTGAACAGGGCTT
25 CAGTCAGCGCAAGTATTCCATTTGCATTTGGTAATTTTTTCATGCCACCTATTTATG
AATATATAAATCTTTATACCAAATCTATTTTTTAAACATGGAAAAGTTGCCTTTA
TGGAAACTTGGCAGAGCCAGAGTGTACACATTCCTAAACCATTAAACAGATTTCTA
TA

30

Sequence ID - 606

nt: 556

GGATAATGATACCTCTGACCTTTCTTCCTTTTGGGAAGTACTTGAGTGTGCAGCTG
CATGAGGCCTCAGCAGGAGAGAGATTTTAGGTCCAAGAAGCTATACCAGTAGGACA
AGGCAGGAAAATACTACACTTTTCAGGATCAAGCCCCTCTGACTCTCATTTGGAAAC
TGGATGTTTGCTAAGCACCTGCTTCTTAAGGATGCCGAGGGATTAAATGATACTCC
35 CAGAAACCTGGAGAGATTAATGGGGCCTATGGAGAAGTGCTCTGAACTCAGTGTTG
GGACTTGAATAAAATTAACCATTGTCATGTTTTTCAGAACAACTAAGCTGTTTTATA
TTTCATGTGCATGAAAGCCCTAGAACTAAGTTGTGTTATTTCCAGAAATGAAATAG

- 187 -

ATCCCACAGTTAGATGATGTGGCCATTAGGAAGTACCAAATTTATAAAAATCACTG
GAGGTCTGTCTGAGCAGTACCTAATAAAATATAGTATACTGAAAGTGAACAGATCT
TTGTCTCTTTCTTTGGCTGCTTGATACTTTATCTGTGTCTGCCGGACAGTGC

5 Sequence ID 607

CAATAAAAGCAGGTTAACCTCAATGATAGCAGTTAAAATGTTCTATCTTATGTATT
TCTTTTAAGTATTACCATTATGGTGCTACTGAGCGTTTTCTTTTGGTAAAAAGAAA
AATGCCATGGGCTGCAGTCTTCTTCCATCACTTTTCCCTACCAGGTCCATTAATAT
GCTTATAACACTAGTGCCAGTTATTTTATTTGATAATGCTTATGGTATTTGTATAT
10 TTGTTTGCATTCCAAATTTGTTTAATAATGAGTGTGTAAACTGCATACGTAAATA
AATGTAAATACTAATGTACTGCTGC

Sequence ID 609

TTTTATTACCCAAGTTTTAACCTCTGTCTGGTGATTTGTTGTTGTTGTTGTTGTNG
15 TTGTTGTTGAAGTTCAGGCTGCATGTGGGATAGGTTTGCTCAGGCATACTTCTTAG
GAAGTAGTCACTTGCATGACTGTTTTTGGGATAACTCTTTGAGTATTTGGAGAGGT
CTATTGTAACCTTCTGAAAGGCATTGTTTTTACGTATGAATGTTCTAAAATTCATTC
TAAATGGTCATGAAAAGAAAAGGATTCACATTTTAGAATGGCAATAGTCCCTGAGG
ACTATTATGTCTTTTAGATTTTCTGTGGGTTTCTAGGAATGTTAGTGTAACCTTANA
20 TTCCACCTACCTGATTTCTGGATGTGCCTATTGGAACCTGCTGAGATCTTTTTTT
TTCCTTAACATGTTGTCCCTTGACCCGTACTTCGAAACTAAACATATTATTTTAT
TTGCTTACACTTCAGGAGGCAATTGGCAGACACCAGGCCAACAGTCT

Sequence ID 610

25 GCTCTGACCCCAAGTTGGAAATGTATCTGTACTTTGTCCGGCTTCCACTCAAGGACC
ATTTATGACATTGCTTGGTGTGCTGACAGGGGCTCTGGCCACAGCTTGTGGGGA
TGACGCGATCCGCGTGTTCAGGAGGATCCCAACTCGGATCCACAGCAGCCCACCT
TCTCCCTGACAGCCCCTTGCATCAGGCCCATTTCCAGGATGTCAACTGTGTGGCC
TGGAACCCCAAGGAGCCAGGGCTACTGGCCTCCTGCAGTGATGATGGGGAGGTGGC
30 CTTCTGGAAGTATCAGCGGCCTGAAGGCCTCTGAGCTACCTCGACTTTGGACAGAG
TAATGACTCCCCAGAAAACGTCATATAAGACTTTACCAGCCCCTGAGAGGACCAGG
AGGAGCATCCTTGACCTTCATTTAACTTGGCTCACTTCTCTTCANACTTGGGTAGA
AGTGCAGAGCCACAAAATTGCTTTCCTTCCCCGCCTTTGACATGAGGCCTTCAGTA
AAG

35

Sequence ID 611

TGCAGGATCCGTCGACT

- 188 -

Sequence ID - 612

nt: 576

5 GAGAAATATAAGATTATGTATAGATCAAATCTACCTCTATTTGGTGTCTGAAAGA
GATGAGGAGAATGGGACAAACTTGGAAAGCTTATTTCAAGATAACATTCTGAGAA
CTTCCCCAATCTTGCTAGAGAGGGCCAACATTAAAATTCAGTAAATGCTGAAAATC
10 CAGTAAGATATTTCTTAAGAAAATTATTTCCCAAGATATATACTCATCAAATTATCT
AAGGTCAAATGAAGGAAAAAATTTTATAGGCAGCTAGAGAGAAATGTCAGGTCACC
TACAAAGAGAATGGCATAAGACAAAAAGTAGAACTCCCAGCAGAAACTCTAAAAGC
CAGAAGAGATTAGGGGCCAATATTTAACATTCTGAAAGAAATTCCAACAAGGAATT
TCATATCCAGCCAAACTAAGCTTCATAATTGAAGGAGAAATAAGATATTTTCCAGA
15 CAAGCAAATGCTGATGAAATCCATCACCACCAGACCTGCCTTATAAGAGCTCCTGA
GGGAAGCACTAAATATTGAAAGGGAAGAACTTTATGAACCATTTCAAAAACACATT
TAAGTNCACAAAGCAG

Sequence ID - 613

nt: 341

15 CCTTATTTTACAGGTGAAAAACCACGAATCAGATAGATTTTTTATTTGCCCAAGTCA
CATAATATTAAGAACAGGCCAAGTGTGGTGGCTCATGTCTGTAATCTGAGCACTTT
GGGAGGCTAAGGCGGGTGGATTTCCTGAGCCTAGGAGTTTGAGATCAGCCTGGGCA
ACATGGCGAAACCTCATCTCTACAAAACATACAAAATTAGTCAGTGTGGTGGTGA
GAGCCTGTAGTCCTGGCTACTCGTGAGGCTGAGGTGGGAGCATCACCTGAGCCTGG
20 GAAGTCGAGGCTGCAGTGGCAACAGAATGGGTAACCTGGACATCAGAGTGAGACCC
TGTCT

Sequence ID 614

25 CTCACACCTGTAATTCCATTACTTTGGAAGGCTGAGAGAGGAGGATCAGTGGAGCC
CAGGAGTTTGAGACCAGCCTGGGCAATATAGGGAGACCCTGTCTCTACAAAATGA
AATAGCCAGGCGAGGTGGCATGTGCCTGTGGTCCCAGCTACTTGGGAGACTGAGGT
GGAAGGCTGCCTTGAGCCCAGGAGTTCCAGGCTGCAGTGAGCCATCATTATGCCAC
TGCACTCCAACCTGGGAGACAGAGTGAGAGAGACCCTGTCTCAAACAAACAAACCC
AAAATAGGCCAGGCACAGTGAATCATGCCTGTAATCCCAGCACTTTGGGAGGCTGA
30 AATAGGCGGATCATTGAGGTGAGGAGTTCAAATTCAGACCAGCCCGGCCAAT
GGCAAAACCACATCTCTACTACAAATAAAAAATTAGTTGGGTGTGGNGGAGCATTC
CTGTAATCACAGCTATTTCAGGAGGCTGAGGCATGANAACCGCTTCA

Sequence ID - 615

nt: 379

35 TAAATTTAAAACATTTTAAATTAGCTGGCATGATGGCATGCACCTGTAGTCCTACCT
ACTTGGGAGGCCAAGGCAGGAAGATTGCTTGAGCCCAGGAGTTTGAGCTTACTGTG
AGCTGTGATCACACCACTGCACTCCAGCCTGGGTGACAAAGGAAGACCGTATTTCT

- 189 -

AAAAAATAAAAAATACAAATACAACTACAACTAGCACTAGACCAACAGTGACTAT
GTACCATGAACTGAGGAATATTATTAATTCACCATTTCATCTGAGGTTAACAAT
ATGTCAATGACTTAAATAACATCATATCTCTGAGAGTAATTTCTCCTATATTTCCA
TGACAAATGTTAGATAATTTTCCATTTTTTCCATTCAACAAAA

5

Sequence ID 617

TTTTCAGGCATGTCAGAGAAGGGAGGACTCACTAGAATTAGCAAACAAAACCACCC
TGACATCCTCCTTCAGGAACACGGGGAGCAGAGGCCAAAGCACTAAGGGGAGGGCG
CATACCCGAGACGATTGTATGAAGAAAATATGGAGGAACGTGTACATGTTCCGTAC
10 TAAGTCATTTTTCA/GGGGATTGAAAGACTATTGCTGGATTTTCATGATGCTGACTGGC
GTTAGCTGATTAAACCCATGTAAATAGGCACTTAAATAGAAGCAGGAAAGGGAGACA
AAGACTGGCTTCTGGACTTCCTCCCTGATCCCCACTCTTACTCATCACCTGCAGTG
GCCAGAATTAGGGACTCAGAATCAAACCAGTGTAAGGCAGTGCTGGCTGCCATTGC
CTGGTCACATTGAAATTGGTGGCTTCATT

15

Sequence ID - 618

nt: 598

GATTAACCTTTCATTTTAAGCTCTTCTCTACTAATTCTGTTCGTATGTTTATTCATT
TTGCGTTGATCATATTTTGTACACCAGGCACCTCTTCTCAGTTTTATATGTGTGTTA
ATTTACTCCTTTCAAGAGCCCTATGATACATGAATTTATCTCCATTTTATAGATGA
20 GGAAATTAAGACCTAGAGTTACTGAACTTGCCCAAGGTTATACAGCTGATGGGTAG
GGCCAGAACTTTGCCTCAGAGAATCTGAATTTCCAAAAAATAACCTAAAAGAGAAA
TTTAAGTACTAATTAGTAAGCAAAGAAATGCACATTTAAGGAAGACAGTGACATT
TAAGGAAGACAGTAACCTTTTATCTATTAGAGAAAAACACACATTCTGTCTTTAAC
ACACACATAAATCTTATATTGGCAGGGATTTTCTTTATTCAGCAATTATTTATTGG
25 TTGTCTGCTTTGTGGTACACATAAATGCTGGGGATAAACACTTAATAAAATATACT
TCCTTCTCTTGAATATCTTGCACTTTAAGTGGAAGGTAAGTCAACAGAGTAGAGG
TGATATATCCAAGTGATAGACTGTTTCATTGCCAGTAG

Sequence ID 619

30 GTTGCCTGAGAGTGACCTTTGCATCTGCCTGTCCAGCCAGCATGGAACCAAAGCGG
ATCAGAGAGGGGCTACCTTGTGAAGAAGGGGAGCGTGTTCAATACGTGGAAACCCAT
GTGGGTTGTATTGTTAGAAGATGGAATTGAATTCTATAAGAAGAAAAGTGACAACA
GCCCCAAAGGAATGATCCCGCTGAAAGGGAGCACTCTGACTAGCCCTTGTCAAGAC
TTTGGCAAAGGATGTTTGTGTTTAAGATCACTATGACCAAACAGCAGGACCACTT
35 CTTCCAGGCAGCCTTCCTGGAGGAGAGAGATGCCTGGGTTCCGGGATATCAATAAGG
CCATTAAATGCATTGAAGGAGGCCAGAAATTTGCCAGGAAATCTACCAGGAGGTCC
ATTCGACTGCCAGAAACCATTGACTTAGGTGCCTTATATTTGTCCATGAAAGACAC

- 190 -

TGAAAAAGGAATAAAAGAACTGAAT

Sequence ID 621

5 TGGTACTGAACCTACGAGTACACCGACTACGGCGGACTAATCTTCAACTCCTACAT
ACTTCCCCCATTTATTCCTAGAACCAGGCGACCTGCGACTCCTTGACGTTGACAATC
GAGTAGTACTCCCGATTGAAGCCCCCATTCGTATAATAATTACATCACAAGACGTC
TTGCACTCATGAGCTGTCCCCACATTAGGCTTAAAAACAGATGCAATTCCCGGACG
TCTAAACCAAACCACTTTTACCGCTACACGACCGGGGTATACTACGGTCAATGCT
CTGAAATCTGTGGAGCAAACCACAGTTTCATGCCCATCGTCCTAGAATTAATTCCC
10 CTAAAAATCTTTGAAATAGGGCCCGTATTTACCCTATAGCACCCCCTCTACCCCCT

Sequence ID 622

TTTTTCTTGTTTTTGTGTGTCTACCTTGGCATATACTAAAGGAAGGTGTGTATTCA
TTTATTACATGATATCTCTGGGTATAATTATTTACATATATGAATTTGAAAGAAA
15 GATTGAGAGGGATATGTGTGACCTTTGTTTCATTATGATCATTACATGACTAAAG
ATAAAGATCATATGTCTGATTTTCAGTTTAATGGCAAGTTACTTAAAATAAATGAA
ATATGTTTTTIATTGTTTTTCGTGGGTTTGATGCTTTGTGTTTTATTTCAAGTAACTT
GAGAATGCATTGTGTTTGGTACTGTTTTTTATGAATATCATTAAAAATTTATTTAA
GGAGAGAGTAATTTTGCAATAATATTTTTGATTTATTTGAAAATAAAATTCAAGAT
20 AAATGAAATAATTGAAATTTTCTAAAGAAGGAATTGAATATATTTTTTACATTTGAA
TGAACTAAGGATTAAGTGAACCATTTATATATAGTACTTTCAGAACTGAATGTCTT
AAATGATAAAGCTCTAATTGGTTAAAGTGACTTTCTTTCAAGTCAAAGAACCCAGA
AACTGAATAGATGATCTAACTACTGCCACTGAGGTTTTGGATTAGTGAGTATAAAT
TT
25

Sequence ID 624

TGCAGGATCCGTGCGACT

Sequence ID 625

30 GACAATCAGAGCAGATCTTGGGCTTCTGTGGCTCATCTCAGCCCTTTATAACTGGC
CTGAGAAGAGGGTTTATCTACTTGTGCAAGTGGCCAGAAATCTCACTCGTACATG
AGGCTTTGGAACATCCTTGCAAAGGTACGCTGAAAGCAAATTGCTGTTTTCTGCT
GGTCTGCACGTTTTCTAACTTTTATCATAGTTTGTATTTTCATTATTTAAGAAAAA
ATAAAAAATCCAAAGACCATAAGATGGCATTAGATTTTTTACCATTAAATTATTA
35 TGCCTATTTGGTGCTCATAAAGATTAATCATGTACGCATGTTTCCAATCTTTCTT
TTGCAGTATATTATTTTCTAAAAATTGTTACATGCAAATTTAAACCAAGATTTATC
AGTA

- 191 -

Sequence ID 626

TTGGAAGAAATAAACCAAGGCAGAAAAATTTTAAATGGCCAAAATAAATTGTATTG
CTAACTTAGATGGCCACAGATGGGGGCAGGGGTGGAGAGAGGAGAAATTGAAAACN
CCACAAAGACCCCGCAATGGCTAGA~~A~~CTTGAAATCTCTGGATATTGCAACAATAGC
5 AGCCTCCTTAAGTCAGCAAAAAGATAAAGATTGATCCAATGTTCTATATTACAGAA
CAGAGCAGATTGTCAATATAGCAAATAAAGTTACCGTTGAGTGGACTGCGCTGTNT
AAGCTGCTTGGTTGGCCTTAAGTGCCGACAATTAAGAGATGAAGGCAATGAGAACT
GAAACAAACATTTAAGTTCAAGACCCAGTTTACTGACACTGGGACTATTACTATAT
CTCTTTGGGCCTCAGTTTACTTATCTGTAAACATTAAGAGGTTGGATTACATGATGT
10 CTCACGATTCTTTT~~T~~TTTTTATTTAGAGATGGGGTTTTGCTCTGTTGCCCAGGCTGG
AGTGCAGTGGCATGATCATAGCTCACAGCAG

Sequence ID 627

CCAGCCTGTCACTGGCCTGGCCAAGGAGGAGAGACAGGCCAGGGATTCTGGTCCTA
15 ACTCTACTGGCCACACTGTGTGGCCTGAGACCCCCCTTCCCTCCCAAGCCCCCTGC
CTCCGCATCTGCGTGGTGAAGGCCATTGGCCCTCATCGGTGGATCTGCGTTTCCTC
GGGCCTACACTGTCTAGGATTGTGCGGGGCTGGTGAGAGAACAAGATCTCTTCCGT
GTTCAAGGCAGACTTCCTGCCCCCTGCACCCTGCTCTCTCCAGGCCTTGAGGTCA
GTGTGAGCCCCAAGGGCAAGAACACTTCTGGAAGGGAGAGTGGATTTGGCTGGGCC
20 ATCTGGATGGAAGGTAAAAAAGAAAATCCCTTGAAAGGAGATTGAGGGAAGTTT

Sequence ID - 628

nt: 419

AAGAGAAAGGACTCAGTGTGTGATCCGGTTTCTTTTGTCTCGCCCCGTGTTTTTGT
AGAATCTCTTCATGCTTGACATACCTACCAGTATTATTC~~CC~~GACGACACATATACA
25 TATGAGAATATACCTTATTTATTTTTGTGTAGGTGTCTGCCTTCACAAATGTCATT
GTCTACTCCTAGAAGAACCAAATACCTCAATTTTGT~~TT~~TTTGAGTACTGTACTATC
CTGTAAATATATCTTAAGCAGGTTTGT~~TT~~TTTTCAGCACTGATGGAAAATACCAGTGTT
GGGT~~TTTTTTTT~~TTAGTTGCCAACAGTTGTATGTTTGCTGATTATTTATGACCTGAAA
TAATATATTTCTTCTTCTAAGAAGACATTTTGT~~T~~TACATAAGGATGACTTTTTTATA
30 CAATGGGAATAAATTATGGCATT~~TTTT~~T

Sequence ID 629

CTGAGAGTCACTGTGTTTTAGCCAAATCTAAGGGAGAAAATGAATATTGATAGCA
GCATGCTGTAGCCAGCTCCTTAAAGGAAGGATGGTGCCTGGTACAGAGTTAGAGTT
35 AGTGCTTCAGTAAATAATGAATGTGTGCTAGGTAGGTTCTGCTGGGTAGGCTGCAT
GCATTGACCAATTTATTCCTCCTTGT~~TT~~TTCAAACAGGATTTAAGGGCACTTATATA
TATATATTTTTTAGT~~TTTTTT~~TAATGTAAATGAGAGAATAAAGATATATATATATGT

- 192 -

CTATATATGTATATATGTATATATATGTCTATATGTCTATATGTATATATGTCTAT
ATGTATATATGTGTGTGTGTATATATATATATATATATAAGTTTTCTGTTGCTA
GCATAACAACTACCAGAACTTAGCAACTGAAACAACATGAATTTATCTTACGGT
TCTATAGTTT CAGAAGTCTAACGTGTCACTGGGATGAAATCCAGGTTTCAACAGGAC
5 TGGGTTCCCTTCTAGCTCATT CAGCTACCTGGCTCATT CAGGTTGTNGGCAGAATA
TACTTCCATGAACTGTAGGGCTGAGACCCCGTTCTTCTTCTGGCTATCATCTGAAA
ACTTTC

Sequence ID 630

10 AGGCGCAGCCCAGCCTCGAAATGCAGAACGACGCCGGCGAGTTCGTGGACCTGTAC
GTGCCGCGGAAATGCTCCGCTAGCAATCGCATCATCGGTGCCAAGGACCACGCATC
CATCCAGATGAACGTGGCCGAGGTTGACAAGGTCACAGGCAGGTTTAATGGCCAGT
TTAAAACTTATGCTATCTGCGGGGCCATTTCGTAGGATGGGTGAGTCAGATGATTCC
ATTCTCCGATTGGCCAAGGCCGATGGCATCGTCTCAAAGAACTTTTGACTGGAGAG
15 AATCACAGATGTGGAATATTTGTCATAAATAAATAATGAAAACCTAAAAAAAAAAAA
AAAAAAAAAAAAAA

Sequence ID 631

TNCACTCACACACTCCCAAACCTTAACAAACACATACATGTGCAGCCAACCCAATG
20 GGCCAGCCTCTTTTATGCTCCTCACATGTTTCCTTTAACTGGAATACCCATGACAG
CTCCCTACATAGTTACTTGTAAACTCCTCCTCTCTGTATAAGTTTTCTGAATTTT
TTTGATAAAATTAAGTTGTGCCACCCCTTTATGCTCTCTTANAACCTTGTCTGT
CTCATGGCTGTTCTGCAACGAATCTCATTGTGTTCTCCTACTCAATTACATTCCTG
CGTCTCCCACTAGATGGCAGACTCTTTGAGAGTAGGAGATTCCCTTGTTATCTCTG
25 GATCCCTGGCACTTGCGAGAAAGCCTGTTACGTAATAATTGCTCAACAATTAGTTTT
TAAATAAATGAATTATTTTAAACGCCAAAATTACAATGATTGTGCATTAGTGA
AAGATGACCATCTAAAACATAAAGCCATGCTTCATGACATTGGC

Sequence ID 632

30 GACCATT CAGGGAAATTTTATAAAAAATGCAGATACTGTCTTGAGCAGATCGAAAT
GCCGATGAGGTGGATGCAATTTCTTTTGTGCAAGCAGTGCACGGTGCCCCCCCCCT
CGGGTGTCCGTGCTGTGCCTTAGCTTCCCCAGGTGCCGGGACTCACACCTGCTAGG
GGCTGGGCAAGGCCCGGCTCTGCTTTCTCTGAAGGGCTTGTCCAAGTTCATTGCC
CTGTTACAGGTGGTCAAGACGTCCGGCCGCTTGACCCAGGCTACCCTTAGCCAAT
35 ATCTCTGCCCCCTGGGTGGTTGGTGGCTGGGCCTCAGGGTGGGCAACGTTAGGGGT
TTGGCGAAAGCCCGCCCCATGGGATTGAGGGACGGGGCTGCACTCCAACCGTCTGC
ACCTGCTCTTCCCCCACCCTGTGGGACCTCATCTTCACGTGCCATGTGTGCTGAA

- 193 -

GGCCCAGGGCCCAGCAGGGGGCAGTGGCACCTGTTGACGGAAAAGCCGAGGTGCTT
ACCAATGGACCTTCTGGCCCCGCCCTCCCCTGTACTTGTCGGGCATTGAGGGCCCCG
ACCTGTGCCTACCCGCA

5 Sequence ID 633

CAACTGTGTTCACTAGCAACCTCAAACAGACACCATGGTGCACCTGACTCCTGAGG
AGAAGTCTGCCGTTACTGCCCTGTGGGGCAAGGTGAACGTGGATGAAGTTGGTGGT
GAGGCCCTGGGCAGGCTGCTGGTGGTCTACCCTTGGACCCAGAGGTTCTTTGAGTC
CTTTGGGGATCTGTCCACTCCTGATGCTGTTATGGGCAACCCTAAGGTGAAGGCTC
10 ATGGCAAGAAAGTGCTCGGTGCCTTTAGTGATGGCCTGGCTCACCTGGACAACCTC
AAGGGCACCTTTGCCACACTGAGTGAGCTGCACTGTGACAAGCTGCACGTGGATCC
TGAGAACTTCAGGCTCCTGGGCAACGTGCTGGTCTGTGTGCTGGCCCATCACTTTG
GCAAAGAATTCACCCCAACAGTGACAGGCTGCCTATCANAAAGTGGTGGCTGGTGTG
GGCTAATGCCTGGCCCCACAAGTATCACTAAGCTCGCTTTCTTGCTGTCCAATTTG
15 TATTAAAGGTTCTTTGTTCCCTAAGTCCAACCTACTAAACTGGGGGATATTATGAA
GGGCCTTG

Sequence ID - 634

nt: 511

TTTTTTAATTTACCAAATTTGTTGACGTCCCTTGATTTGCTGATAGGGACAATA
20 ATTAAATATTTTCCACTTGTTTTTATAAAACTGTAATGGTGATTTGTTTAACAGA
TGTTGACTTAGCACCTTCTCTCTTTTTTTTTTTTTTTTTTTGAGTTGGAGTCTTGC
TCTGTCACCCAGCTGGAGTGACAGTGGCACGATTTCCGGCTCACTGCAACCTCCGCCT
CCCAGGTTCCGGGCGCTTCTCCTGCCTCAGCCTCCCANATAGTTGGGATTACAGGTG
CATGCCGCCACNCCTAGCTAATGTTTTTTGTATCTTGGGTANANATGGNGTTTCACC
25 TTGTTGCCCATGCCGCTCTTGAACCTCTTGGCCTCCCAAAGTGTTAGGATTACAGG
CGTGAGCCACTGTGCCTGGCCCCAATTTANACCTTACTGGGTGCTGAGGCTGTGA
GCCATAGTAGAATGCATGTGATCCAGGGCCTTGCTGAATTCATGGGCTAATAGGGA
GCCTGAC

30 Sequence ID - 635

nt: 592

TGAGCGTTGGGCTGTAGGTGCTGTGCTGTGTGATCCCCAGAGCCATGCCCGAGA
TAGTGGATACCTGTTTCGTTGGCCTCTCCGGCTTCCGTCTGCCGGACCAAGCACCTG
CACCTGCGCTGCAGCGTCGACTTTACTCGCCGGACGCTGACCGGGACTGCTGCTCT
CACGGTCCAGTCTCAGGAGGACAATCTGCGCAGCCTGGTTTTGGATACAAAGGACC
35 TTACAATAGAAAAAGTAGTGATCAATGGACAAGAAGTCAAATATGCTCTTGGAGAA
AGACAAAGTTACAAGGGATCGCCAATGGAAATCTCTCTTCCTATCGCTTTGAGCAA
AAATCAAGAAATTGTTATAGAAATTTCTTTTGAGACCTCTCCAAAATCTTCTGCTC

- 194 -

TCCAGTGGCTCACTCCTGAACAGACTTCTGGGAAGGAACACCCATATCTCTTTAGT
CAGTGCCAGGCCATCCACTGCAGAGCAATCCTTCCTTGTCAGGACACTCCTTCTGN
GAAATTAACCTATACTGCAGAGGTGTCTGTCCCTAAAGAACTGGTGGCACTTATGA
GTGCTATTTCGTGATGGAGAAACACCTGACCCA

5

Sequence ID - 636

nt: 572

CTTANAAGAGTTGCTCATTCACACCCACGCCCTTGCCCAAGGCTGGCCCACTCAGA
GCGAAACTTAACTTTTGTCTGGATGGGAAGAGAAGTAAGTCTACCCCGAGGTTGCC
ATGTTGAAGAGTGAGAGGTCCAAGTGATTCTGTGCATTGAAACCAAGACACCCAC
10 CCAGAACACTTCTTCCCTCCCTCAGCCCAAACCAAAGGCTGGGGTTCTCATCTCCA
AGTGGCTGTTCTCCAACCTTTCCCAAGCCGCTTGCAATCCCCAGACTGGACTACTGT
GGCGGTTAGGTTAGATTTGAAGACGGGGCCAGGCTGGGTATGAACGGGTGCAGCC
CTCTTCTCCTCTTCCCCCCCACATCTCTCATGAGAGAGGTAGTGGCATTTCCTTCT
CAGGGAGCTTCAATGGGAAAGGTCTCGAAAGCTTCAGGAGGAGCAGAATACCAACG
15 CAGGGGGATGGCTGTAACGATCTCACCGTCTCCTAACCTCAGTCCCTTTTTTGAGA
GTGAATGGTGGAGGGTGGGAAAGGGACCCAAATTTGTAGATCTCTTTGTCTGGGGG
AGGGGAANGATG

Sequence ID - 637

nt: 482

TTAAAACAGGCGCAGGGGTAAAAATGAGAATGAATCTGAAAAAAGAGAGTTGGTGT
20 TTAAAGAGGATGGACAAGAGTATGCTCAGGTAATCAAAATGTTGGGAAATGGACGA
TTGGAAGCATTGTGTTTTGATGGTGTAAAGAGGTTATGCCATATCAGAGGGAAATT
GAGAAAAAAGGTTTGGATAAATACATCAGACATTATATTGGTTGGTCTACGGGACT
ATCAGGATAACAAAGCTGATGTAATTTTAAAGTACAATGCAGATGAAGCTAGAAGC
25 CTGAAGGCATATGGCGAGCTTCCAGAACATGCTAAAATCAATGAAACAGACACATT
TGGTCTTGAGATGATGATGAAATCCAGTTTGACGATATTGGAGATGATGATGAAG
ACATTGATGATATCTAAATTGAACCAAGTGTTTTTACATGACAAGTTCTCTGAGGA
TGGTTCACAGTTGGGATTTTGGCCATCATCAAC

30 Sequence ID - 638

nt: 545

TTTGAAGGCAAAGAGGGATTAATCTGTGCTGGCATCATGTAAGGAGACTTGATAGA
TAAGAAAAAGCTTTACCTAAGTTTTGAAGAATAGGTTTTTTCATAATGGAAAATTTA
AGGGAAAAATCTCCAAAAAAGTGCTACTCAAGTTTTATCCATTTGTATTTCCAACA
CAGCCTAGGACAGTACCTGCACATAGTAGGTGATTAATAAAAATTTAGAAAGCATT
35 AATACTAAAGAGGAAAAATAGCAATGGCAAGAAAACACATGTAGGGAACACATGTA
GCCAAAAAATAATATATAATCAGAGAAATAATAGGACTTCTGGAAAAAAAAGATGA
GATCAGATTGGTTAGGATCTTTACTAACATGACAAGAGCATGAATTTTTTTTCTGT

- 195 -

AGATAATAAGTATGAAAGAATTTTAGCTTAAAAATTAGCATAATTTGGATCCACAT
ATGCAAATCAATGAATGTAATTCATAATATAAACAGAACTAAACACAAAAACCACG
TGATTATCTCAATAGACACAGAAAAGGCCTTCAAAAAAATT

5 Sequence ID - 639

nt: 624

GACACACGAGCATATTTACCTCCGCTACCATAATCATCGCTATCCCCACCGGCGT
CAAAGTATTTAGCTGACTCGCCACACTCCACGGAAGCAATATGAAATGATCTGCTG
CAGTGCTCTGAGCCCTAGGATTCATCTTTCTTTTACCGTAGGTGGCCTGACTGGC
ATTGTATTAGCAAATCATCACTAGACATCGTACTACACGACACGTACTACGTTGT
10 AGCCCACTTCCACTATGTCCTATCAATAGGAGCTGTATTTGCCATCATAGGAGGCT
TCATTCACTGATTTCCCCTATTCTCAGGCTACACCCTAGACCAAACCTACGCCAAA
ATCCATTTCACTATCATATTCATCGGCGTAAATCTAACTTTCTTCCCACAACACTT
TCTCGGCCTATCCGGAATGCCCCGACGTTACTCGGACTACCCCGATGCATACACCA
CATGAAACATCCTATCATCTGTAGGCTCATTCATTTCTCTAACAGCAGTAATATTA
15 ATAATTTTCATGATTTGAGAAGCCTTCGCTTCGAAGCGAAAAGTCCTAATAGTAGA
AGAACCCTCCATAAACCTGGAGTGACTATATGGATGCCCCCACCCTACCACACAT
TCGAAGAA

Sequence ID 641

20 CAAGATGACAAAGAAAAGAAGGAACAATGGTCTGTGCCAAAAGGGCCGCGGCCACG
TGCAGCCTATTCGCTGCACTAACTGTGCCCCGATGCGTGCCCAAGGACAAGGCCATT
AAGAAATTCGTCATTCGAAACATAGTGAGGCGCAGCAGTCAGGGACATTTCTGA
AGCGAGCGTCTTCGATGCCTATGTGCTTCCCAAGCTGTATGTGAAGCTACATTACT
GTGTGAGTTGTGCAATTCACAGCAAAGTAGTCAGGAATCGATCTCGTGAAGCCCGC
25 AAGGACCGAACACCCCCACCCCGATTTAGACCTGCGGGTGCTGCCCCACGTCCCCC
ACCAAAGCCCATGTAAGGAGCTGAGTTCTTAAAGACTGAAGACAGGCTATTCTCTG
GAGAAAAATAAAATGGAAATTGTACTTAA

Sequence ID 642

30 TGCTTGGCCCTCTACCTCCTGCCCTCTTCTGTTCATCTCCCAACCACTGCACTCT
TGATTTTATATACCACACAGAAGGTAAGAAAATTCTAGGAACCCTAAGGATCAATCC
TCTCCATTTTCACTCAAATGCCTGGGGCCCAGCTCTGCAATGACTGACTCCAGGGC
CTCTTTCCTCACTGCCAGCATAGAAGTCAGGGGAGCCAGCTGGGCCCTGCGGTCAG
GAAGGTCTCATTTTTTGGAGCATTCCCTGAGCCCAGATCATAGGAGCAGCTGTCCC
35 TGGTGGGACACAGGAGTCATGACTCCTACCCTCCACCCTCCACACCCACCAGGCAT
TTAGCAGTCTGTCTATGCAAGACAGATGAATTCTCAGCCAGGATACCTCAAGGCA
GGCAAAGGTGAGTGGAGGGAAAATTCACAAACATTTCAGGGTGTGTGGTGTGCTGGCAT

- 196 -

CACCATGGCCAAATCCAAGAGGTCTTCCTGGAAGAGGGCCCAAACCTGGAACCAAAA
GAATGCTGTCAGCAGTTGGAATAGAGCTGTGAATT

Sequence ID 643

5

CTTTCCAAGAGGAATCCTCGGCAGATAAACTGGACTGTCCTCTACAGAAGGAAGCA
CAAAAAGGGACAGTCGGAAGAAATTCAAAAGAAAAGAACCCGCCGAGCAGTCAAAT
TCCAGAGGGCCATTACTGGTGCATCTCTTGCTGATATAATGGCCAAGAGGAATCAG
AAACCTGAAGTTAGAAAGGCTCAACGAGAACAAGCTATCAGGGCTGCTAAGGAAGC
10 AAAAAAGGCTAAGCAAGCATCTAAAAAGACTGCAATGGCTGCTGCTAAGGCACCTA
CAAAGGCAGCACCTAAGCAAAAAGATTGTGAAGCCTGTGAAAGTTTCAGCTCCCCGA
GTTGGTGGAAAACGCTAAACTGGCAGATTAGATTTTAAATAAAGATTGGATTATA
ACTCT

15

Sequence ID 644

CTTTGATAGAGAAGAAAATTCTCCTAGGATACAAGAGCCTCAACATTTTAAAGATT
TTCTGCATCTCAAAAGCGTAGGCTCCTTGCTGGGCAAGGTGAGCCTCTGTGAGTCC
TCATAGGACCGAGCAAATCTGATTCACCCAGAAAATCCAATATCGAAGCTGAGCT
TTGGCCTGAGCGGGTTCCATTTCTCCCCAGATCCTATTTAGGAAGTGTCTCCTGA
20 CAACCTCCAAAAGGTGCTAACATGCAACGTTCTGAAGGGTTATTGCTCAAAAACAA
GATTTTCCTTGTGGTCAAGACTCTGCGAGCCTCGAACACGATGAATCCGCTCGAAT
GGGCTTGGGCTTTGCCCGGGTGGCGCACGCTCACACGCTGGAAGCACAGCTTTGAC
GATCTCCACACACGCACAGGCACACAGCCACAGATGATGCCGGCTCATTTCTCAGG
GGGTGTCTAAGTTCTGCTTTAAATATTTACCCCTAATTGTACAAACAATAGGGGC
25 ATGAGCCTGGTACTCGATAAATGGGGACTTNCTTAAAA

Sequence ID - 645

nt: 649

30

CTACAGCCTGGGCAGCGCGCTGCGCCCCAGCACCAGCCGAGCCTCTACGCCTCGT
CCCCGGGCGGCGTGTATGCCACGCGCTCCTCTGCCGTGCGCCTGCGGAGCAGCGTG
30 CCCGGGGTGC GGCTCCTGCAGGACTCGGTGGACTTCTCGCTGGCCGACGCCATCAA
CACCGAGTTCAAGAACACCCGCACCAACGAGAAGGTGGAGCTGCAGGAGCTGAATG
ACCGCTTCGCCAACTACATCGACAAGGTGCGCTTCCTGGAGCAGCAGAATAAGATC
CTGCTGGCCGAGCTCGAGCAGCTCAAGGGCCAAGGCAAGTCGCGCCTGGGGGACCT
CTACGAGGAGGAGATGCGGGAGCTGCGCCGGCAGGTGGACCAGCTAACCAACGACA
35 AAGCCCGCGTCGAGGTGGAGCGCGACAACCTGGCCGAGGACATCATGCGCCTCCGG
GAGAAATTGCAGGAGGAGATGCTTCAGAGAGAGGAAGCCGAAAACACCCTGCAATC
TTTCAGACAGGAAATCCAGGAGCTGCAGGCTCAGATTCAGGAACAGCATGTCCAAA

- 197 -

TCGATGTGGATGTTTCCAAGCCTGACCTCACGGCTGCCTTGCGTGACGTACGTANC
AATATGAAAGTGTGGCTGCCAAAAACCTTGACG

Sequence ID - 646

nt: 600

5 GAGATGTCTCGCTCCGTGGCCTTAGCTGTGCTCGCGCTACTCTCTCTTTCTGGCCT
GGAGGCTATCCAGCGTACTCCAAAGATTCAGGTTTACTCACGTCATCCAGCAGAGA
ATGGAAAGTCAAATTTCTGAATTGCTATGTGTCTGGGTTTCATCCATCCGACATT
GAAGTTGACTTACTGAAGAATGGAGAGAGAATTGAAAAAGTGGAGCATTACAGACTT
GTCTTTCAGCAAGGACTGGTCTTCTATCTCTTGTACTACACTGAATTCACCCCCA
10 CTGAAAAAGATGAGTATGCCTGCCGTGTGAACCATGTGACTTTGTACAGCCCAAG
ATAGTTAAGTGGGATCGAGACATGTAAGCAGCATCATGGAGGTTTGAAGATGCCGC
ATTTGGATTGGATGAATTCCAAATCTGCTTGCTTGTCTTTTAATATTGATATGCT
TATACACTTACACTTTATGCACAAAATGTAGGGTTATAATAATGTTAACATGGACA
TGATCTTCTTTATAATTCTACTTTGAGTGCTGTCTCCATGTTTGATGTATCTGAGC
15 AGGGTGCTCCACAGGTAGCTCTAGGAGGGCTGGCAACTTA

Sequence ID 647

CGAATGTGCAGGTTTGTACATAGGTATATATATGCCATGATGGAAATATTTATTT
TTTTAAGCGTAATTTTGCCAAATAATAAAACAGAAGGAAATTGAGATTAGAGGGA
20 GGTGTTTAAAGAGAGGTTATAGAGTAGAAGATTTGATGCTGGAGAGGTTAAGGTGC
AATAAGAATTTAGGGAGAAATGTTGTTTCATTATTGGAGGGTAAATGATGTGGTGCC
TGAGGTCTGTACGTTACCTCTTAACAATTTCTGTCCTTCAGATGGAACTCTTTAA
CTTCTCGTAAAAGTCATATACCTATATAATAAAGCTACTGATTTCCAAAAA

25 Sequence ID 648

AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Sequence ID - 649

nt: 425

30 CAAAAAACGAAGAAAAGTGACGACAGTCTGAGGGACTTATGGGAGATCATCAAGT
GAACCACTATATGTGTAATGTAAGTCTTGGAATGAGAAGAGAGAAGGAGAAGGAGG
AGAGAGCTTATTTGTAGAAATAATGGCTGAAAACATCCCAAACTTTCCTTTTGTG
AGGAAAGAAATAGGCATACAAGTTCAAGAACTCAAGGAACTCCAGAGAGGACAAT
TCTAAAGACACCCCTCTAACATACATTATAATCAAATTGTCAAAGTAAATACA
AAGAGAATCTTTTAAATTGACAAGAGAAAAGCAGCTGGTCACGTTCAAGGGAGTTC
35 TATAAGAATTTTCAAGCAGATTTCTCAGCAGAAACCTTGACAGGCCAACAGGCAGTGGG
ATGATACATTCAAAGTGCAAAAAAAAAAAAAA

- 198 -

Sequence ID 650

CGAGAGTTTACCAGTNGCCTAATAATGCAATAAAAAATGCTTTGAGATAGCTAACN
GCCCATAAAAACAACTCAAATTGCTTATAAAGTTTCTTCCCATGTTCCCATTTGAT
5 'GAAAAGTCTTACATCACATATAACTGGGAAGCAGGGGTCCCTCCTCAATTTTCAGA
CATTTTGAAAGGATGACAGTTCTGTTTGTAGATGAGTAAACCTCTATATTCATAA
GTTCTAAAATCCTTCATTATGAGGGATTCAAAGTATTTATAAAAAACACTGCCCTCT
AAAAATTTCTCAGATCTGAAGTATGGNCTTGGNCCTGAATATACAGTGTTATCCT
ATGTTTAAAAGGGTGATCCAGACATGAGACGCAACTAGTTGGTGATAAGAAGGCC
CCACTTGGCTATTTTCATATCTACCTACAATTGACCAAAAAAATTTTTTAGGCCAG
10 CAATTATTATTTAGCTTCGCTCTTTCTAGTGCAAGAACTGCAGGCTGGATCAGTA
GTTCAACAGCTAAACAGTCATAAAATAGTCATTGGCATGTTAAATTTCTTTCAATG
CTTCAAAGATAAATTCCAATTCTATTTACTTATTTCATTGNGACNGNATTACTAAAC
AGGTAAGGATGGGAATA

15 Sequence ID - 651

nt: 251

CTTTGGGAGGCCGAGGCGGGCGGATCACTTGAGGTCAGGGGTTCGAGACCAGTCTG
GCCAACATGGTGAAACCCCAACTCTACTAAAAATACAAAAGTTAGCCAAGTGTGGT
GGCAAGTGCCTGTAATCCCAGCTACTCGGGAGGCTGAGACAGGAGAATCACTTTGA
ACCTGGGAGGCGGAGGTTGCAGTGAGCCAAGATCGTGCCACTGCACTTCAGCCTGG
20 GCAACAGAGCAAGATTCCGTCCATCTC

Sequence ID 652

CTTCTTCAGCCTTGCAGACACCTAAACATCATGTAATTACCTAAGGAATTCCCAA
GTGCCTCTTCCAGGTTATACGTGTAAATAGCTGTTTTTATGCAAGATTAGTTAGAT
25 ACTGCTCTTTACAGGATGAGTGGTGTGTTCTTTGGCTGGGGGGGNCCTAAATGTGT
TTCTAATGTGTGTGTCAAATAATTACCTGTAAACAGACTGCCAATCTGGCTGAAG
CCAATGCTTCTGAAGAAGATAAAATTAAAGCAATGATGTCGCAATCTGGCCATGAA
TACGACCCAATCAATTACATGAAGAAACCTCTAGGTCCACCACCTCCATCTTACAC
GTGTTTCCGTTGTGGTAAACCTGGACATTATATTAAGAATTGCCCAACAAATGGGG
30 ATAAAACTTTGAATCTGGTCCTAGGATTAAAAAGAGCACTGGAATTCCCAGAAGT
TTCATGATGGAAGTGAAAGATCCTAATATGAAAGGTGCAATGCTTACCAACACTGG
AAAATATGCAATCCAACCTATAGATGCAGAAGCATATGCAATTGGGAAGAAAGAGAA
ACCTCCTTNTTACCAGAGAGCCATCTTNTTTCT

35 Sequence ID 653

GTTGTGACTCGTTGGCATGTGATCTGAAGTTCCTGCCCTGCAGCTGACGAGCCAGT
GTTTCAATAATTAAAAACAACTCAACTCACTGTCCTCCTGCCTTGAATTTGATCAT

- 199 -

TGCGCTTTGCATGTATGTATCACAATACCACATGTACCCCATAAATATGTACAAAG
ATTATGTGTCAATAAAAAACAAAATTAAAATCCCAATTTTTTA

Sequence ID 654

5 GTTGCTAGTAGCGGCAGGAAGATGTCAGGCTCACTTTCCTCTGATTCCCGAAATGG
GGGGAACCTCTAACCATAAAGGAATGGTAGAACAGTCCATTCCCTCGGATCAGAGAA
AAATGCAGACATGGTGTACCTGGATTTTTTTCTGCCCATGAATGTTGCCAGTCAG
TACCTGTCTCCTTGTCTCTATTTTTTGGTTATGAATGTTGGGGTTACCACCTGC
ATTTAGGGGAAAATTGTGTTCTG

10

Sequence ID 655

GTCCCCGGGAATCGCGGCCGCTCGACGGTTTATTTTCAGTGCTTGAAGATACATT
CACAAATACTTGGTTTGGGAAGACACCGTTTAATTTTAAGTTAACTTGCATGTTGT
AAATGCGTTTTATGTTTAAATAAAGAGGAAAATTTTTTGAAAAAAAAAAAAAAAAAA
15 AA
AAAAAAAAAAAAAAAAAATTTTT

15

Sequence ID 656

TAGAGGCCTGAATAGGTAGACAATGGCAGCAGCGTTTTTAATCACAGTCCTATTCA
20 TGCCCTAATTCGGGAGTGATGATTAAAGGACATTAGAGGGAGCACTTTGACATCTG
ATCCTTTGAACTGACGTCTGTGCAGGCTGCACTCCATAGAGCTCACTTGGCCAAAC
TGATTTCCCTTAAATAAAGTGCTGTGATTTCCAATGTAGGAAATATTACATTAGAGC
CTATTGAAATGATTAGGAATTGAGGAGCTTTTCTTTAGGTGGGAATGTGGTGTATG
CTGTATACTCACAAAAGTGAGATCATTAATATTGCATGTACTACTTTGAATATCAG
25 GGACCACAGAGAAATAGCATGAGAAACGCCTTCCTGCAGTCATGCACTTAAATGA
ATATGAACAAAATGTGGAACCTGCTGTCATAGCTCTCCG

25

Sequence ID 657

GGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCGCAGATAAGTTTTTTT
30 CTCTTTGAAAGATAGAGATTAATACAACCTTAAAAAATATAGTCAATAGGTTACT
AAGATATTGCTTAGCGTTAAGTTTTTAACGTAATTTTAATAGCTTAAGATTTTAAG
AGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAAAAGGTTTCT
AAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAATGTATTTAAA
AGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACCAATAGAAGGGCAAT
35 GCTTTTAGATTAAATGAAGGTGACTTAAACAGCTTAAAGTTTA

35

Sequence ID 658

- 200 -

GACCTTTGAGAAAATTAATTTAAATCCTAGAACTTTGGGTGAACCGAAGAAATTTA
TAATATTTGTTTGTAGTTAATAACAGATAAAAAGGAAAGATTCAAGCCTATTGGATGA
GAATTTGTACATTATTTTAGAGCTAATAATAATGGTTTTTCAGTTTAGTGAGGATTT
AAAAAATGTTTTTGAATCAAACCTTTTTTCTTTATAATCCTTTTTTAACCTAAGTCAAG
5 GAAATAAGGTATTATGAAATCCACACACTGTTACCTCCTTAAAGTATGAGGATACT
TCCCCTGTTTGGTCCACTAGTGGCTGATTATTTTGTGTGGATTATTTGTAATT
TTCTTTTTTAATTCTTCCTTAAAGAGCATGGCATTGAGTACACAGACCTATATTTG
AATCCTGTCATTTACTAGCGTTTTGACCTTGAACAATTATGCTCAGAGTCTCAGTT
TTTTCTTGTAAGTGATGATGATACTACTTAACTCACAGGGTTGTAGTGAAGATCA
10 AATGAGATCATGTCTGTANAACACCCTGCCCGGCACTCAATAAGTATTAATAGGAA
CCCATATACCTC

Sequence ID 660

TGTTTTTATTTTTTAAAAGGTATAAACACCAAAAAAAAAAATTAACATTGTATGAAG
15 ATGGAAAATAAGAAGATGCACCTTCTGTAACTTTGTCTAAGGATTTAAATTACTAA
CTTATGAACTCCAATTTGAATTGAACTTAACTATCGGCTTCTTACTGGTAAAATT
ATATGGTTTTATTTTAAATGCGTACATATTGACCAATGGCCTCTGAAAAGCACATT
TTAGATACTGAAATTGAAGGAAAGAAAATGCATCTTCAAACATTTTTTTGGAATCTC
ACCACATATACTTTGTTANATTTGTGTATTGTAGGGTGTTTGTGTTTGTATTTTTGT
20 ATTGATATATGAACTTTTTTTTAAATGTGACAGTTAAACACATCTTTAAAAGCATAGT
CACAGACAAAAGCATACAGTATAAAAATTTCTTGAAAACCTCTACAATATTATAT
TTGGAGGCAGCTTCAGACTGTTTTATTGG

Sequence ID 661

CTCTGGCACACATTAGTTCCTCTTATATTACATTGATATAAGCAAGTCATATGGAT
25 TTATCTGAGTGTAAGGAGAGCTGGAAAAATAGTTTCTAGCAGGTCAGCCACCTCC
CAGTGAGGGCTGCATACCATAGAAGGGGAGAATGAATTTTGGGAAAACAGGTAATT
ATCTCTGTCACAGAAGGGGATGAAAAGTATGGTAGTTACNCAAGTTANACATCTGT
ATGGAAAATACCACTTGGTTCTACAAATGNNG

30

Sequence ID - 663

nt: 627

GCCTCCCGGGTTCAGGGATTTCTCCTGCCTCAGCCTCCTGAGTGGCTGCATTGCAG
GCACCTGCCACCACGCCTTGCAAATTTTTGTGTTTTTAGTGAGATGGGGTTTTGC
CATGTTGGCCAGGCTGGTCTCGGACTCCTGACCTCAGGTGATCCGCCCGCCTCAGC
35 CTCCCAGAGGGCTGGGATTACAGGCGTGAGCCACTGTGCCTGGCCCCAAGTTTTGC
ATCTTTTAATGCCCTCTGAACAAATACATAGAGAAAACCTCTCAGAACAATTA AAC
CTGCAGAGCAACAGTGTCTCCATGTCTTAGGTTTCAAGTTTGCCTCTAAAATTCT

- 201 -

AATCCATATTTTTCTACTTCTCAGATAATTTATGTGTGTGTACTCTTCCTAGACGT
ACAAGAGACTTTTAAATGCTAAATATTTGTGAGTGCTTAACAAAACTCAATTTCA
CATTACTCATATTGTTTTTGTTTTAAATTGAATGTGAATTAAATTTTTATTAGTTAT
TTGATTTGGAATGTTATGTATGCCATTAACACTATTAGGGGAATCTCTAGCATTTC
5 TGTATTTTTTAAAGAATTTGATTCTTTTGTANATTCTGCCTGTGTGGCATTTTAAAC
ATGTGTGACAT

Sequence ID - 665

nt: 345

ACCGGCGACATGGCCAAACGTACCAAGAAAGTCGGGATCGTCGGTAAATACGGGAC
10 CCGCTATGGGGCCTCCCTCCGGAAAATGGTGAAGAAAATTGAAATCAGCCAGCACG
CCAAGTACACTTGCTCTTTCTGTGGCAAAACCAAGATGAAGAGACGAGCTGTGGGG
ATCTGGCACTGTGGTTCCTGCATGAAGACAGTGGCTGGCGGTGCCTGGACGTACAA
TACCACTTCCGCTGTACAGGTAAAGTCCGCCATCAGAAGACTGAAGGAGTTGAAAG
ACCAGTAGACGCTCCTCTACTCTTTGAGACATCACTGGCCTATAATAAATGGGTTA
15 ATTTATGTA

Sequence ID - 666

nt: 252

ATAATTCAGAACTTCTTCATATGCTCGAGTCTCCAGAGTCACTCCGTTCTAAGGTT
GATGAAGCTGTAGCTGTACTACAAGCCCACCAAGCTAAAGAGGCTGCCAGAAAGC
20 AGTTAACAGTGCCACCGGTGTTCCAAGTGTTTAAAATTGATCAGGGACCATGAAAA
GAAACTTGTGCTTCACCGAAGAAAAATATCTAAACATCGAAAACTTAAATATTAT
GGAAAAAAACATTGCAAAATATAAAAT

Sequence ID 669

TTACTTTTAAACCAGNGAAATTGACCTGCCCCGTGAANAGGCGGGCNTGACACAGCAA
25 GACGAGAAGACCCTATGGAGCTTTAATTTATTAATGCAAACGGTACCTAACAAACC
CACAGGTCCTAAACTACCAAACCTGCATTAAAAATTTTCGGTTGGGGCGACCTCGGA
GCAGAACCCAACCTCCGAGCAGTACATGCTAAGACTTCACCAGTCAAAGCGAACTA
CTATACTCAATTGATCCAATAACTTGACCAACGGAACAAGTTACCCTAGGGATAAC
30 AGCGCAATCCTATT

Sequence ID 670

GGCTGATTCCTGAGCTATAAAAGCATAATTGCTTTATATTTTGGATCATTTTTTAC
TGGGGGCGGACTTGGGGGGGGTGCATACAAAGATAACATATATATCCAACCTTCT
35 GAAATGAAATGTTTTTAGATTACTTTTTCAACTGTAAATAATGTACATTTAATGTC
ACAAGAAAAAATGTCTTCTGCAAATTTTCTAGTATAACAGAAATTTTTGTAGATG
AAAAAATCATTATGTTTAGAGGTCTAATGCTATGTTTTTCATATTACAGAGTGAAT

- 202 -

TTGTATTTAAACAAAAATTTAAATTTTGGAATCCTCTAAACATTTTTGTATCTTTA
ATTGGTTTATTATTAAATAAATCATATAAAAATT

Sequence ID 671

5 CAGGAAGTCACCTGGGATTGGCTGCCTCACCCACTCACAGTGCCATCCCTGCCCCA
GGCCTCCCAGTGGCAATTCCAAACCTGGGTCCCTCCCTGAGCTCTCTGCCTTCTGC
TCTGTCTTTAATGCTACCAATGGGTATTGGGGATCGAGGGGTGATGTGTGGGTTAC
CTGAAAGAACTACACCCTACCTCCACCACCTTACCCTCACCTGGAGAGCAGTTAT
TTCAGAACCATTCTACCTGGCATTTTTATCTTATTTAGCTGACAGACCACCTCCACA
10 GTACATCCACCCTAACTCTATAAATGTTGATGGTAATACAGCATTATCTATCACCA
ATAACCCTTCAGCACTA

Sequence ID 672

CAGGAAGTCACCTGGGATTGGCTGCCTCACCCACTCACAGTGCCATCCCTGCCCCA
15 GGCCTCCCAGTGGCAATTCCAAACCTGGGTCCCTCCCTGAGCTCTCTGCCTTCTGC
TCTGTCTTTAATGCTACCAATGGGTATTGGGGATCGAGGGGTGATGTGTGGGTTAC
CTGAAAGAACTACACCCTACCTCCACCACCTTACCCTCACCTGGAGAGCAGTTAT
TTCANAACCATTCTACCTGGCATTTTTATCTTATTTAGCTGACAGACCACCTCCACA
GTACATCCACCCTAACTCTATAAATGTTGATGGTAATACAGCATTATCTATCACCA
20 ATAACCCTTCAGCACTAGATCCCTATCAGTCCAATGGAAATGTTGGATTANAACCA
GGCATTGTTTCAATANACTCTCGCTCTGTGAACACACATGG

Sequence ID 673

GGGTTTTCTTTTCGGAAGCGCGCCTTGTGTTGGTACCCGGAATTCGCGGCCGCGTC
25 GACTGCTAAACAGAATACTGCTATTTTGAGAGAGTCAAGACTCTTTCTTAAGGGCC
AAGAAAGCCACNTGNNCCCTNGGNCTAATCTGGCTGAGTAGTCAGTTATAAAAGCC
NTAATNGCTTNNNTNTTTGGNNTCNTTTTNNCNGGGGNCGGNCTTGGGGGGGGTTG
CNTCCAAAGATANCATNTNTTCCAACTTTNTNAANNNAANNGTTTTAAATCCCT
TTTCCNCCNGAAAAANANNGCCCTTTAAGNGCCNCAAAAAAANNGTNTTCTGCAN
30 NTTTTCTANTATNACAAANNTTTTNGTAGAANAAAAATTTTTTTTTAGNGGCTACC
CTTTNTTTNTTANNCANNGGAGTTTNTTTTTACAAAAAANATTTGGGNCCCCT
CCACAACCTTGGGTCTNTAATNGGGGGGTTTTTAAATAAANCNTNTNTAAATCCCC
CNNNNNNNNNCNNNNNNNNNNCCNNNNNNNNNNNNNNNCCNNNNAAAAAATTTTTNC
TCCCCCNCCCTTTTTCTTCCTGCCGGCCCCAATTTAAGCCCGGGCGCTTGGGGCAA
35 ATCCCCCTTTAGNGGGGGGGTTTTANAAAAACNGGGGCGGGGNTTTAAACCNCGG
GGNNNGGGGAA

- 203 -

Sequence ID 674

ACCTCTAGCATCACCAGTATTAGAGGCACCGCCTGCCAGTGACACATG-TTTAAC
GGCCGCGGTACCCTAACCGTGCAAAGGTAGCATAATCACTTGTTCCTTAATTAGGG
ACCTGTNTGAATGGCTCCACNAGGGTTCACCTGTCTCTTACTTTTAACCAGTGAAA
5 TTGACCTGCC

Sequence ID - 675

nt: 591

GTATAGAAAATAATGTCCCCAGNGCATAGAAAAAATGAGTCTCTGGGCCAGTGAAT
ACAAAACATCATGTGCGAGAATCATTGGAAGATATACAGAGTTCGTATTTTTCAGCTTT
10 GTTTATCCTTCCTGTTAAGAGCCTCTGAGTTTTTTAGTTTTTAAAAGGATGAAAAGCT
TATGCAACATGCTCAGCAGGAGCTTCATCAACGATATATGTCAGATCTAAAGGTAT
ATTTTCATTCTGTAATTATGTTACATAAAAGCAATGTAAATCAGAATAAATATGTT
AGACCAGAATAAAATTAATTATATTCTGGTCTTCAAAGGACACACAGAACAGATAT
CAGCAGAATCACTTAATACTTCATAGAACAAAAATCACTCAAAACCTGTTTTATAAC
15 CAAAGAATTTCATGAAAAAGAAAGCCTTTGCCATTTGTCTTAGAAAGTTATTTTTTA
AAAAAAAATCATACTTACTATTAGTATCTATGGAAGTATATGTAACAATTTTTTATG
TAAAGGTCATCTTTCTGTGATAGTGAAAAAATATGTCTTTACTAAGTTGAAATGAA
TACTTTCTGNCTTTGCTAATGGATAGTTATT

20 Sequence ID 676

CTCAATTCTACTAAAAAGCCCCCAAGAAAAGCGAATGAGAAAACAGAGTCATCCT
CTGCACAGCAAGTAGCAGTGTCACGCCTTAGCGCTTCCAGCTCCAGCTCAGATTCC
AGCTCCTCCTCTTCCTCGTCGTCTTCAGACACCAGTGATTTCAGACTCAGGCTA
AGGGGTCAGGCCAGATGGGGCAGGAAGGCTNCGCAGGACCGGACCCCTAGACCACC
25 CTGCCCCACCTGCCCCCTTCCCCCTTTGCTGTGACACTTCTTCATCTCACCCCCCCC
TGCCCCCCTCTAGGAGAGCTGGCTCTGCAGTGGGGGAGGGATGCAGGGA

Sequence ID 679

GNANCNTTTCCTNTCGNAAANCGCGCCTTGTGTTGGTACCCGGGAATTCGCGGCCG
30 CGTCGACAAAAAANNTNTAGACTCGANCAAGCTT
ATGCANGCNTGCGGCCGCAATTCGAGCTCGGCCGACTTGGCCAATTCGCCCTATAG
NGAGTCGTATTACAATTCAGTGGCCGTCGTTTTACACGTCGNGACTGGGAAAACC
CTGGCGTTACCCAACTTAATCGCCTTGCAGCACATCCCCCTTTCGCCAGCTGGCGT
AATANCGAANAGGCCCGCACCGATCGCCCTTCCCAACAGTTGCGCAGCCTGAATGG
35 CGAANGGAAATTGTAAGCGTTAATATTTTGTAAATTCGCGTTAAATTTTTGTTA
AATCAGCTCATTTTTTAACCAATAGGCCGAAATCGGCAAAATCCCTTATAAATCAA
AAGAATAGACCGAGATAGGGTTGAGNGTTGTTCCAGTTTGGAAACAANAGTCCACTN

- 204 -

TTAAAGAACGNGGACTCCAACGTCAAAGGGCGAAAAACCGTCTATCAGGGCGATGG
CCCACTACGTGAACCATCNCCTAATCAAGTTTTTTGGGGTTCGAGGNGCCGTAAAG
CACTAAATCGGAACCTAAAGGGAGCCCCGATTTAAAGCTTGACGGGGAAAGCCC
GGCGAACGTGGCGAAA

5

Sequence ID 682

CACCTGCAGTCCAAGTACATCGGCACGGGCCACGCCGACACCACCAAGTGGGAGTG
GCTGGTGAACCAACACCGCGACTCGTACTGCTCCTACATGGGGCACTTCGACCTTC
TCAACTACTTCGCCATTGCGGAGAATGAGAGCAAAGCGCGAGTCCGCTTCAACTTG
10 ATGGAAAAGATGCTTCAGCCTTGTGGACCGCCAGCCGACAAGCCCCGAGGAAAAGCTG
AAACTTTGCTTAACNACCGAATGGNGGGGANCTTTTCCAACGNTTTT

Sequence ID 683

TTGGTTTCATACTGNTGGGGNTTGAATGNTCCCTNCAACACTNATGTTGANACTTA
15 ATCCCTAATGNGGCAATACTGAAAGGTGGGGCCTTTGAGATGTGATTGGATCGTAA
GGCTGTGCCTTCATTCATGGGTTAATGGATTAATGGGTTATCACAGGAATGGGACT
GGTGGCTTTATAAGAAGAGGAAAAGAGAACTGAGCTTGCATGCCC

Sequence ID - 684

nt: 545

20 GTGGAAGNGACATCGTCTTTAAACCCTGCGTGGCAATCCCTGACGCACCGCCGTGA
TGCCCANGGAAGACAGGGCGACCTGGAAGTCCAATACTTCCTTAAGATCATCCAA
CTATTGGATGATTATCCGAAATGTTTCATTGTGGGAGCAGACAATGTGGGCTCCAA
GCAGATGCAGCAGATCCGCATGTCCCTTCNCGGGAAGGCTGTGGTGCTGATGGGCA
AGAACACCATGATGCGCAAGGCCATCCGAGGGCACCTGGAAAACAACCCAGCTCTG
25 GAGAAACTGCTGCCTCATATCCGGGGGAATGTGGGCTTTGTGTTACCAAGGAGGA
CCTCACTGANATCAGGGACATGTTGCTGGCCAATAAGGTGCCAGCTGCTGCCCCGTG
CTGGTGCCATTGCCCCATGTGAAGTCACTGTGCCAGCCCAGAACACTGGTCTCGGG
CCCGATAAGACCTCCTTTTTCCAGGCTTTAGGTATCACCCTAAATCTCCAGGGG
CACCATTGAAATCCTGAGTGATGTGCACTGATCAAGACTGG

30

Sequence ID 685

GGAAAGGGCCATTTTATTGCCTAAAACCACCTGGNTTTTNAGGTAACAGTTCCAAC
ATGTCCTTTTTTTGAATAGCTGTTCTAATTATTATATATTCAGCTGATTAATAGGAG
TACTTGATAGGTGGACTGTGTCAGGTAGCCTCAGGCAATCCTACTTCAACAAGCTG
35 TCAGGGAGCCATGCCATGCTTCTTTATGACATAGGTGAATTTGATAGGCTCACTAG
CAGAACATGGGATCACAAGGTGGAACNNTTCCNTTT

- 205 -

Sequence ID 686

GACCCCTTCCTTACACCTTATACAAAAAACTGAAACTGGACCCCTTCCTTACACC
TTATACAAAAATTA ACTCAATTTTATTATGTTGTATTAAATTAAGTTGGGTTTAAT
TAAGATGGATTAAAGACTTAATTATAAGACCTAAAACCATAAAAACCCTAGAAGAA
5 AACCTAGGCCATACCATT CAGGACACGGGTATGGGCAAAGACTTCATAACTAAAC
ACCAAAGCAATGGCAACGAAGTCCAAATAGACAAATTGGACCTGATTAACTAAA
GAGCTTCAGCACAGCAGAAGAGACTATCGTCAGAGTGAACAGGCAACCCACAGAAT
GGAAGAAAATTCTTGCAATCTATCCATCTGACAAGGGGCTAATATCCAAAATCTAC
AAAGAACTTAAACAAATTTACAAGGAAAAACACAAACAACCCCATCAAAAAGTGGG
10 CTAAGGATGTGAACAGACACTTCTCAAAGAAAACATTTATGCAGCCAACAAACAT
GAAAAAAAGTTCATCATCACTGCTCATTAGAGACATGCAAATCAAACCACAATGA
GATCCCATCCACACCAGTTAGAAATGGCAATCATTA AAAATGT

Sequence ID - 687

nt: 268

15 TTTATGTGTTTTTGGCTTGGGGGGCGCTGGGCCTAGCCAGAGTAGTGCTTGCTCCC
CCTGCCTTGTCCCACCAGGGAGGCAGCAGACTCAGGCCCTCCATGGTCCCTCTTTGT
CATTTTGTGACATGCATTCCTCCTTTTGT CATCTTGTTGGGGGGAGGGGATTAAC
CAAAGGCCACCCTGACTTTGTTTTTGTGGACACACAATAAAAGCCCCGTTTATTTG
TAAA

20

Sequence ID - 688

nt: 569

CTTTAGCCAGCCTGATCAGAAAAAAACAAAAGAAGAGGAAAGACGTAGATTACCAA
CATCAAGAATGTGAGTTATGATATCACTACAGACTCTCCAGGTATTAAAAGCATAA
TTAGAGAATGATATGAGCAGCTATATGCAAATAAGTTCAACATTGGACAAATGGAC
25 AAATTTCTTGAAAGATAAATTATGAAATTTCA TTCTGAAAGAACTACATGACCTTA
ATTGTCTTACATCTATTAAATAAGTGGAATTTGTAGTTTAGAACTTTCCACAAA
GAAACTCTAGGCCCAGATGGCATCAAATAATATTCAGATGAATGAAATGGAGAA
AGGATAGCCTTTTCAACAAATGGTGGTGGAACAATTGGATTTCATATGCAAAAAA
ATAGAGATGGACGCAGAGGTGTGTGCTTAGGAGGCTGAGGTGAGAGGATTGTTTGA
30 GGCCAGCCTGGGCAACATAGCAAGACCCCATTTCAAAAACAAAATAAAGAACTTG
TAGCCTTACCTTGTGCCATATTATGAAAATGTATCATAGGCTTAAATGTGAAACGT
AAAACAAA

Sequence ID 689

35 CGCAGGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCGCAGATAAGTTT
TTTTCTCTTTGAAAGATAGAGATTAATACAAC TACTTAAAAAATATAGTCAATAGG
TTACTAAGATATTGCTTAGCGTTAAGTTTTTAACGTAATTTTAATAGCTTAAGATT

- 206 -

TTAAGAGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAAAAAGG
TTTCTAAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAATGTAT
TTAAAAGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACTAATAGAAGG
GCAATGCTTTTAGATTAAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTTTAAA
5 AGTTGTAGGTGATTAAAATAATTTGAAGGCGATCTTTTAAAAGAGATTAAACCGA
AGGTGATTAAAAGACCTTGAAATCCATGACGCGAGGGAGAATTGC

Sequence ID 690

CGAAAAGCAAATATAACTTGCCACTAACCAAGATCACCTCTGCAAAAAGAAATGAA
10 AACAACTTTTGGCAGGATTCTGTTTCATCTGACAGAATTCAGAAGCAGGAAAAAAA
GCCTTTTAAAAATACCGAGAACATTAAAAATTCGCATTTGAAGAAATCAGCATTTC
TAACTGAAGTGAGCCAAAAGGAAAATTATGCTGGGGCAAAGTTTAGTGATCCACCT
TCTCCTAGTGTTCTTCAAAGCCTCCTAGTCACTGGATGGGAAGCACTGTTGAAAA
TTCCAACCAAAACAGGGAGCTGATGGCAGTACACTTAAAACGCTCCTCAAAGTTC
15 AAACCTTAGATTTTCAGATTT

Sequence ID 691

CCGGTCTCTACACAATATATAGAAATCTGGGCATGGTGGTGCCTGGCTGTAGTCTC
AGCTACCTAGTTGGGTGAGGTGGGAGAGTCGCTTGAGTCCTGGAGGTTGAGGCTGT
20 AGTGAGCCAGGGCTGCACCACTGCATTCCAGCCTGGGTAACAGAGTGAGACCCTGT
CTCAAAAAGAAAAAAAAAAATTGCTAATTTTAACAAATCACAAAACCTGACTCAGGC
AAGTTGTCTGACTCAAAGCCCTTGAAAAACCATCAAAGACAGTAGAATGTTAACT
GGTCATTTACGTAAAATAGTGTTTCATTAAATTTTGGTTCATTTAGGATAATCATT
TTAAATGAGACTGTATTTGAGACTGTATACACATACATATACATGTTTACACACAT
25 ATACGTACAATATATGTACATTCTATCTAAAAGATCATACTGTGTGTACATATAT
GTTTTTAAAAGTCAAACCTGACATATTAATGGAAACAGTGCTTACATCTCTGGTAGT
GATTTTCTATTAGCAGCAGCCCTACATATGCTGCGTCTCTGAACAGCATGTCAGTG
CCATGACTGTCTAAACATGCAAATATGACTGACAGACTCTTGAGACAGCTTTCACC
TTG

30

Sequence ID 692

AATTCGNGGCCGCGTCNNCCTANGAGGCACCAGGAAATCCCGCGGGGTGGCCCATG
CAGACCAGGCGCACGTGGCTCATGGGGCANAATTGCCAAGGACAGCTCACGACAGT
GCCACCTTCTCACCATTCCAGCCAAGGAGAGATGTGACGTTGGAACCTGCTCTGGCA
35 CTTCTGTCAAGCCTCCCCCGCCCCAATTGCCTTGAGATCTCTGCTCTTTGTCAGAG
ATTTGCAAAGACTCACGTTTTTGTGTTTTCTCATCATTCATTGTGATACTAAGA
AACTAAGAAGCTTAATGAAAAGAAATAAAATGCCTATGTTGTTGTTCT

- 207 -

Sequence ID 693

CTAGAACCCATGACTCCTAGGTCTTATACTGCAACCACAGTATCAGCAAATAATCT
TTCATAAGGGGATTATTCTCTGATTAACAGGAAATACAGGAATTTAATTTGTGAAC
ACGCTAGGTAGAAGCAGAAACCCAAATCCAAATCCAAATTTAAACATTTAAAATTC
5 ATTCTATAACTAAGATCTAACAGTCATTTTCTTCCCAGTAAGAAATAACCAAAGCA
TGCTAAAAATCACTGGACTAAATTGGTGTCAAACTGCCACATTGCCAGGCATGGG
GGGGTCATACTTGTAATCCCAGCACTTTGGGAGGCCGAGGTGGGAAAATTGCTTGA
GGCCAGGAGTTCGAAACCAGCCTGGGCAACACAGTGAGACCCCATCTCCACAAAAA
AAAAAAATTAAAAAACAAAACAAAACATTAGCTGGGCATGGTGGTACACGCCTGTA
10 GTCCCAGCTACTCAGGAGCCTGAAGTGAGAGGATCACTGAAGCCCAGGAGGTAGAG
CTATGACTGTAGTGAGCTATGACTGTGCCACTACACTCCACCTGGGTGACAGGGGA
CTC

Sequence ID 694

15 CGACTTCCATTTGTATTAATGGAATACTAAGTCCCTCTGTGATTTCTGAACCAAGC
TATTCCTAGGCCTGAGTTTTATTTTGTGACACAGAAATAAATTANAAGGCCAAGC
GTGGTGGCATGTGCCTGTAGTCCTAGTTGCTGAGGTAAGAGGATTGCTTGAGCCCA
GGAGTTCAAGGCTGCAGCAAGCTTTGATTGCGCCACTGCACTCCAGCCTTGGCGAC
AGACTAAGACGCTGTCTCAAAAAAAAACAAAAA

20

Sequence ID 696

GGTTATCAATGAGATTAAGAGACAACTAGAGTAAAAACAAAAGAAAAGAAAAGAAA
NGAAAACAACAGAAGCTCTATTAAGTACCTCTAACCAATACAACAGGTAACTGA
TGTTCTCCATTCTGTATATAAAAATCCCAGTGGACACCCACAACACAGGCTTCAGG
25 CTTGTAGGACACTTTCTAGTTCATCTGAGCACTTTTGTTCCTCAGCAGTTGAGCTGT
ATACTTAGCAACATTTGGTGCTTCCAAACCCATTTGTGCCTGTAGCACTTACTATT
GAAATACATAATTTAATTAAATATTATATAAAGGAATGGAATACGAGTTGGACAAG
AAAAAGAGTTAAATCTGAAGGTTAGGTAAAAAGAGCAACTTCTTTTCTCTGTTTTG
CAGGTTGGCAAAATCATTTAAAAACAATTGGAAGTATTATATGTTCTGCATTAAAGT
30 TGTCAATTTTACTTAAAAACTAGGCATCAAAGATGATGCATAATAAATTTAGTGTAT
GCAAGAATGACTGCTTGGGACCTCAATATATGAATTCTTAATCCAAGGAAAGTCCT
TGGCCTTACATTTAAAAGTCGGCAAATAAGTGTACGTTCAATT

Sequence ID 697

35 GAACATTTAAAAATAATGCAAATAAGGCTGGGCGTGGGGGCTCACACCTGTAATCC
CAGCACTTTGGGAGGCCGAGGCAGGCAGATCACGAGGTCAGGAGATTGAGACCATC
CTGGCTAACACAGTGAAACCCTGTCTCTACTTAAAAAATAAAAAAATTAGCCAGGC

- 208 -

GTGGTGGTGGGCGCCTGTAGTCCCAGCTACTCAGGAGGCTGAGGCAGGAGAATGGT
GTGAACCCGGGAGGCGGAGCTTGCANTGAGCTGAGATCGTGCCACTGCACTCCAGC
CTGAGCGACAGAGCGAGACTCTGTCT

5 Sequence ID 698

TCATTAGAATCCAAGCTTTGAAAATTTCTGATTAATGCTCATGTATTTCTTTATCT
TTGTTTTTTCCTTGTGAAGAAAGACTTTCACCACTGTCTGAGTGATGATGCTGTTGA
TAAGGATGATGTCGATGACTACTATATTGCATCTCTCAGGAACAGCTGATGGGAAG
GGAGGGGCTGCTGAGTTCCTTGTCTAGCTAGCAGCACGCTCCTCANAGAGGGGG
10 CCGAGTTACAGACAGCAGCCGCAATTCTCATGCAAAATTAGTTTTAACTGCTAGTG
TGGGCATCGGTA¹CCTTTTGCCTGGGTGATACCGAAGAATTGTTGAGGATTTAGTAT
GCTCCGTAGAGACAGTTCAGCCAGTCATTTCTGCATTGGAGAGACTTCTCATACTT
TCTTTGAAGACTCATAGAAAGCTGGAT

15 Sequence ID 699

ATTAAGGTTTGTNCCCAACAAGAATAGATGTAATTAGAAAAAANTGNCTTCCTTAC
CTATTGCCTCTGATNTTTACTTGCTTAAATTTTTTTTATTGNAAATCCAGAAAAAG
NGGATTTAGAGAACAACACTAACTCCACCTAATCTATGACAGANATGTACAANAN
AGTACCTGTGAAAAATGTGAAAGNATNTGAAAAATGTAACCTTTGGCAGCCTGAGC
20 ATAGTCAACCAGAAAACTATCTGAATTAAATAAATTGGTCCATAGGTACTATTTT
ATTTGGTCCATAAGGATTATTTTTTTCAACTTTTTTTTCAAGTGTATTATTATGTCA
TTTCCCACGTAGGTTACTGATACCTGAAGACTTTTTTNCACCTTTAACCTTNCTCGT
TGAGGAGCTTTGTANTCTAATAAAAGAGAAATATAAGTAAATGTTAGATATATGGG
NGGATAATGGTAACTATGTGCTTAAAGAGGTATAAAAGAAGGGTAGGGAGCAGATA
25 AGACAAAGGAAGGGCTATATTATAANGAAGAATATTCCAAGTAGGGAAGAGAAAAA
GATATGTTATCCATATAATATTTTATGTGCAGTAGAGAACATGTTCTATAGAANAG
ACAGAAGATG

Sequence ID 700

30 CTTGAGCCCAGGAATTCCAGCCTGGGCAATATAGTAAGACTCCGTCTCTACAAAAG
ATACAAAATTAGCCAGATGTGGTGGTGCGTGCTGTAGTTCCAGATACTGGAAAG
ACTGAGGCAGGAGGATTGCTTGAGCATGGGAAGTTGAGGCTGCAATGAGCTGTGAT
TACGCCACTACACTCCAGCCTGGGCAACAGAGTAAGATCTTGTCTCAAAAAAAAAA
TTGAATTCAGCTAAAAATAATAAAATTTTAAATAATTTTAAAAAGCCCTCAACAG
35 CTTTGTTTTTCTCTCCTTGCCAGCTTCTCTGCAGCCTATAGCCTGCAGGCTGGCTG
CTGCGAGCCAGGACAAGCGGTGGGAAATGCAATCACAGCGTGAAATCTCTGTGTTC
AGAGACACGCAGGAAGCAGGTGAACCATGAAGGGCCAACACATGCCCCCAGTTAGC

- 209 -

AGGGTGTAGAGACCGGGGCAGGGCTTTCTTCTTCCTTCTGGGTTATAAATATCCAT
GTCCTGCCATTTGAAGCTGCAAGTGGCACACATGGATGCTGGACAGGCGCTCGCAC
TTTCTGGGCAGGGCANGGGGCTCAAAGGCAGGACAGCTGGGCAAAAGCACCTTGCG
TGGGCCC

5

Sequence ID - 701

nt: 579

10

CTTTGGAGCTTCTGTCTGTGCTGTGGACCTCAATGCAGATGGCTTCTCAGATCTGC
TCGTGGGAGCACCCATGCAGAGCACCATCAGAGAGGAAGGAAGAGTGTGTGTAC
ATCAACTCTGGCTCGGGAGCAGTAATGAATGCAATGGAAACAAACCTCGTTGGAAG
TGACAAATATGCTGCAAGATTTGGGGAATCTATAGTTAATCTTGGCGACATTGACA
ATGATGGCTTTGAAGGTAATTAATAATTATCAAATTGGTGCTTGATTTCTGCTTTTA
AAATGGTTTATGGAAGAAAATATGATTAAAGTTTTGTATTGTTTTCTTCCTTATAG
AAGATGGAGCCAGAATGGCATGCTAAGTTTTTCTTTCTTTAGTGTTATATATGA
CTTCTCCTCAATTGTCACCCATTGATCTTTACCACTGTTAATAATGGATGATATTC
AAAATACCTTATTTCAAGTATTCTAAGGCACCATTGATTAGAACTGCATTATTAT
TTATGTGTCCCTAAAAGCTACCTATTAAGCTGTTACACCCACCATTTTTCTGTAA
GAAAATCCTGATTTTCAGAA

15

Sequence ID 702

20

GTNNTCCTCTCGGAACGCGCCTTNTGTAGCCAGGTGCTACCAGACCNAATACACGG
TTGTTCCAGCTTGCGCATTCACCGATGGCGTAGATATCCGGATCGGAAGTCTGGCA
GGAATCATTAATGACAATACCCCCACGCGGAGCAACGTCCAGACCACACTGGGTTG
CCAGCTTATCGCGCGGACGGATACCGGTAGAGAAGACGATAAAGTCGACTTCCAGT
TCGCTGCCGTGCGCAAAACGCATGGTTTTACGCGCTTCAACACCTTCCTGCACAAT
CTCAAGGGTGTTTTTGCTGGTGTGAACGCGCACGCCATACTTTGATTTTGCGAC
GCAGCTGCTCGCCGCCCATCTGATCAAGCTGTTCTGCCATCAGCATAGGGGCAAAT
TCGATAAC

25

30

GTGGGTTTCAATACCTAAGTTTTTCAGCGCGCCTGCGGCTTCCAGACCTAACAGGC
CGCAATTTCGAGCTCGGCCGACTTGGCCAATTCGCCCTATAGTGAGTCGTATTACAA
TTCACTGGCCGTGTTTTACAACGTGCTGACTGGGAAAACCCTGGCGTTACCCAAC
TTAATCGCCTTGCAGCACATCCCCCTTCGCCAGCTGGCGTAATAGCGAAAGAGGC
CCGCACCCGATCGCCCTTTCCAACAGTTGCGCACCTGAATGGCGAATGGAAATTGT
AAGCGTTAATATTTTGTTAAATTTCGCGT

35

Sequence ID 703

CTGCGCAGACCAGACTTCGCTCGTACTCGTGCGCCTCGCTTCGCTTTTCCTCCGCA
ACCATGTCTGACAAACCCGATATGGCTGAGATCGAGAAATTCGATAAGTCGAACT

- 210 -

GAAGAAGACAGAGACGCAAGAGAAAAATCCACTGCCTTCCAAAGAAACGATTGAAC
AGGAGAAGCAAGCAGGCGAATCGTAATGAGGCGTGCGCCGCCAATATGCACTGTAC
ATTCCACAAGCATTGCCTTCTTATTTTACTTCTTTTAGCTGTTTAACTTTGTAAGA
TGCAAAGAGGTTGGATCAAGTTTAAATGACTGTGCTGCCCCTTTCACATCAAAGAA
5 C TACTGACAACGAAGGCCGCGCCTGCCTTTCCCATCTGTCTATCTATCTGGCTGGC
AGGGAAGGAAAGAACTTGCATGTTGGTGAAGGAAGAAGTGGGGTGAAGAAGTGGG
GGTGGGACGACAGTGAAATCTAA

Sequence ID 704

10 CTTGTATTCAAGAACTACTGTAATGCATTAGTGGTCTGGCTTCATTTTGTATGATG
CCAGATCCTTAATTTACCCAGCACAAATCATTTTCAGTAGTTTCCTATGGCTCCTGCA
AAAATGCAAACAGAAACCACCACAGGAACAGCCCCTTGCTGCCTCCTGTTGCTGAG
GTAGTAGTCGCTAAAGAAAATTGAAGGCTCCTTACAATCTATATTTGAAAAGTAGA
ACTTCTGTAGAAACACACAGATCCCGATCTTAGAAGTTGTACAGGACAATCTGGTA
15 AAAGTGACATAATTGTGATTTATTAACATGAATTTAAATGCCCAACAGTGCTTCA
GTGTGACAGTATATTTAAAATAAAAAAGAAATTAAAGGTCATATACTGTACTACTT
TCACAAAGATCCACAGTTTTTGCAAAGACTTGTTCATATGTACAATGCTATATATCA
AATGAGAAAAGCTGTAAGCAATTATATACGCAAAGAAATGGCAGTA

20 Sequence ID 705

TTCCAGTCCTTTTCATTTAGTATAAAAGAAATACTGAACAAGCCAGTGGGATGGAAT
TGAAAGAACTAATCATGAGGACTCTGTCTTGACACAGGTCCTCAAAGCTAGCAGAG
ATACGCAGACATTGTGGCATCTGGGTAGAAGAATACTGTATTGTGTGTGCAGTGCA
CAGTGTGTGGTGTGTGCACACTCATTCTTCTGCTCTTGGGCACAGGCAGTGGGTG
25 TAGAGGTAACCAGTAGCTTTGAGAAGCTACATGTAGCTCACCAGTGGTTTTCTCTA
AGGAATCACAAAGGTAACTACCCAACCACATGCCACGTAATATTTTCAGCCATTCA
GAGGAAACTGTTTTCTCTTTATTTGCTTATATGTTAATATGGTTTTTAAATTGGTA
ACTTTTATATAGTATGGTAACAGTATGTTAATACACACATACATATGCACACATGC
TTTGGGTCCTTCCATAATACTTTTATATTTGTAAATCAATGTTTTTGGAGCAATCC
30 CAAGTTTAAGGGAAATATTTTTGTAA

Sequence ID - 706

nt: 496

CAACCCTCTCTCCTCAGCGCTTCTTCTTTCTTGGTTTGATCCTGACTGCTGTCATG
GCGTGCCCTCTGGAGAAGGCCCTGGATGTGATGGTGTCCACCTTCCACAAGTACTC
35 GGGCAAAGAGGGTGACAAGTTCAAGCTCAACAAGTCAGAACTAAAGGAGCTGCTGA
CCCGGGAGCTGCCCAGCTTCTTGGGGAAAAGGACAGATGAAGCTGCTTTCAGAAG
CTGATGAGCAACTTGGACAGCAACAGGGACAACGAGGTGGACTTCCAAGAGTACTG

- 211 -

TGTCTTCCTGTCCTGCATCGCCATGATGTGTAACGAATTCTTTGAAGGCTTCCCAG
ATAAGCAGCCCAGGAAGAAATGAAACTCCTCTGATGTGGTTGGGGGGTCTGCCAG
CTGGGGCCCTCCCTGTCGCCAGTGGGCACTTTTTTTTTTTCCACCCTGGCTCCTTCA
ACACGTGCTTGATGCTGAGCAAAGTTCAATAAAGATTTTGGGAAGTTT

5

Sequence ID - 707

nt: 397

10

CGGATGTGGTGGCAGGCGCCTCTAGTCCCAGCTACTCGGCAGGCTGAGGTAGGAGA
ATGGCTTGAACCCAGGAGGTGGAGCTGACAGTGAGCCGAGATCGCGCCACTGCACT
CCAGCCTGGGCGGCAGAGCGAGACTCCATCTCAAAAAAAAAAAAAAAAAAATAGA
CTTTGAGACCAGCCTGACCAACATAGTGAAACCCGTCACTACTAAAATACAAAAA
TTACCCGGGCGTGGTGACGGGCGCCTGTAATCCCAGCTACTTGGGAGGCTGAGACA
GGAGAATCACTTGAACCAGGGAGGCGGAGGTTGTAGTGAAC TGAAATCGTGCCCCCT
GCACTCCAGCCTGGGTAACAAGAGCGAAACTCCGTCTCAAAAATAAATAAATAAAT
AAAAT

15

Sequence ID - 708

nt: 293

20

CCAGCTTTTTATGGTGTTTAATCTAATACACTTAAGCTGCAGTCCCAAATTAGGG
GTCCTTCAGTCTTGGAGACTATAAGGGAGCCTCTGCACCCAGGGAAAATGTTACCC
TTTACAGGGGGGAAGGGTAAACCAGTAGGGAATACAGTACAATCCCAACCCTACTG
GGAGGGGGCGGGAGGGAGGTGTTGCCGTCACTGTATTAAAGTCGATGTTGGGAAACGT
TTTAACATCTGGAGCCTTTGTGGGTGGAAATATGTCTCCAGTTACAAC TCCGCAGT
GGATGTGAAGAAG

Sequence ID 709

25

GGAAGCTACAATGATTTTGGGAATTACAACAATCAGTCTTCAAATTTTGGACCCAT
GAAGGGAGGAAATTTTGGAGGCAGAAGCTCTGGCCCCCTATGGCGGTGGAGGCCAAT
ACTTTGCAAACACGAAACCAAGGTGGCTATGGCGGTTCCAGCAGCAGCAGTAGC
TATGGCAGTGGCAGAAGATTTTAATTAGGAAACAAAGCTTANCAGGAGAGGAGAGC
CAGAGAAGTGACAGGGAAGCTACAGGTTACAACAGATTTGTGAACTCAGCCAAGCA
CAGTGGTGGCAGGGCCTAGCTGCTACAAAGAAGACATGTTTTAGACAAATACTCAT
GTGTATGGGCAAAAAACTCGAGGACTGTATTTGTGACTAATTGTATAACAGGTTAT
TTTAGTTTCTGTTCTGTGGAAAGTGTAAGCATTCCAACAAAGGGGTTTTAATGTA
NATT

30

35

Sequence ID 710

TGGATTCCCGTCGTAACCTTAAAGGGAACTTTTCACAATGTCCGGAGCCCTTGATGT
CCTGCAAATGAAGGAGGAGGATGTCCTTAAGTTCTTGCAGCAGGAACCCACTTAG

- 212 -

GTGGCACCAATCTTGACTTCCAGATGGAACAGTACATCTATAAAAGGAAAAGTGAT
GGCATCTATATCATAAATCTCAAGAGGACCTGGGAGAAGCTTCTGCTGGCAGCTCG
TGCAATTGTTGCCATTGAAAACCTGCTGATGTCAGTGTTATATCCTCCAGGAATA
CTGGCCAGAGGGCTGTGCTGAAGTTTGCTGCTGCCACTGGAGCCACTCCAATTGCT
5 GGCCGCTTCACTCCTGGAACCTTCACTAACCAGATCCAGGCAGCCTTCCGGGAGCC
ACGGCTTCTTGTGGTTACTGACCCAGGGCTGACCACCAGCCTCTCACGGAGGCAT
CTTATGTTAACCTACCTACCATTGCCCTGTGT

Sequence ID - 711

nt: 498

10 GTGGTACATATACACAAAGGAAAACCTATGTAGCCATTAAAAGAAAAGGAACTCCTA
TCATTTGTAACAACATAAATAAATCTGGAGGAGATTAGGCTAAGGTGAAATAAGCC
AGGCACAAAAAGACAACTACCATATGATCTTACTTATACGTGTGTGGAATCTAAAA
AGGTGGAATTTACAGAAGCAGAGAGTAGAATGGTGATTACCAGAGGCTGGGGAGTG
AGGGCAGGAGGTTGGAGAAATGTTGGTCAAAGGATACAAAGTTTCAGTTATACAGG
15 ATGAATAAGTTCAAGAGATCTATTGTACAACGTGGTGGCTATAGTTGATAACAATG
TATTGTGTTCTTGAAAAATGCTGAGAGAGTAGATTTTAAGTGTTCTCACCACAAAA
CATAAGTATGTGAGGTAATGCATGTGTTAATTANCTTAATTTAGACATTTCATAAT
GTATTATACATATTTCAAAACCACGTTGTACATGAGAAAGATACACAATT

20 Sequence ID 713

GCCCAGTCGACCCATGTTCTCCTTTCTACACCAGCATTAGACGCTGTCTTCACAGA
TTTGGAATCCTGGCTGCCATTTTTCAGCTGCCATCCATGACGTTGATCATCCTG
GAGTCTCCAATCAGTTTCTCATCAACACAAATTCAGAACTTGCTTTGATGTATAAT
GATGAATCTGTGTTGGAAAATCATCACCTTGCTGTGGGTTTCAAACCTGCTGCAAGA
25 AGAACACTGTGACATCTTCATGAATCTCACCAAGAAGCAGCGTCAGACACTCAGGA
AGATGGTTATTGACATGGTGTAGCAACTGATATGTCTAAACATATGAGCCTGCTG
GCAGACCTGAAGACAATGGTAGAAACGAAGAAAGTTACAAGTTCAGGCGTTCTTCT
CCTAGACAACTATACCCGATCGCATTGAGGTCCTTCGCAACATGGTCACTGTGCAG
ACCTGAGCAACCCACCAAGTCCTTG

30

Sequence ID 714

CTGTAACAGAGATTCCTTTTTTCAATAATCTTAATTCAAAGCATTATTAGACTTG
AAAGGGTTTGATAATCTCCAGTCCTTAGTAAAGATTGAGAGAGGCTGGAGCAGTT
TTCAGTTTTTAAATGAGTCTGCAGTTAATATCAAATGTGAGTTTGGGACTGCCTGGC
35 AACATTTATATTTCTTATTCAGAACCCTTGATGAGACTATTTTTTAAACATACTAGT
CTGCTGATAGAAAGCACTATACATCCTATTGTTTCTTTCTTTCCAAAATCAGCCTT
CTGTCTGTAACAAAATGTACTTTATAGAGATGGAGGAAAAGGTCTAATACTACAT

- 213 -

AGCCTTAAGTGTTCCTGTCATTGTTCAAGTGTATTTCTGTAAACAGAAACATATTT
GGAATGTTTTTCTTTTCCCCTTATAAATTGTAATTCCTGAAATACTGCTGCTTTAA
AAAGTCCCACCTGTCAGATTATATTATCTAACAATTGAATATTGNAAATATACTTGG
CTTACCTCTCAATAAAAGGGTCTTTTCTATT

5

Sequence ID 717

TCCACCCACCTTGACCTCCCAAAGTGCTGGGATTATAGGCGTGAGCCACCTCGCCC
AGCCCGATACTAGGACTTATGCAGAAAAACCTTGACATGGAGGAAAGTAAGATCT
AAATAAATACTGTATTCATAGATTAAAAGACTCAGCATAATAAATATAACCATTCT
10 CCCCAGATTGATGTACAGATTTAACACAATTCCTATCAAGATCCCAGCAAGATTTT
TGTAAGATATGTAAAAGATTATTCAAAAATGTAAAAGGAAGGACAAAGGACTAGAAT
AGATAAAACAAAATGGAGAAAGATTTAATAGGAATCACTGTAAGTGAATTTAAGAC
ATACAGAACAATAATAGAACTGCTTGTATTAGTCCATTTTCACGCTGCTGATAAA
GACATACCTGAGATTGGCAATTACAAAGGAAAGANGTTTATTGGCTTACAGTTCCC
15 ATGGCTGGGGAGGCCT

Sequence ID 718

CTCCTCTGGGTTGAAACCCGGGCGCCGCCAAGATGCCGGCTTACCACTCTTCTCTC
ATGGATCCTGATACCAAACCTCATCGGAAACATGGCACTGTTGCCTATCAGAAGTCA
20 ATTCAAAGGACCTGCCCCCAGAGAGACAAAAGATACAGATATTGTGGATGAAGCCA
TCTATTACTTCAAGGCCAATGTCTTCTTCAAAAACCTATGAAATTAAGAATGAAGCT
GATAGGACCTTGATATATATAACTCTCTACATTTCTGAATGTCTGAAGAACTGCA
AAAGTGCAATTCCAAAAGCCAAGGTGAGAAAGAAATGTATACGCTGGGAATCACTA
ATTTTCCCATTCTCTGGAGAGCCTGGTTTTCCACTTAACGCAATTTATGCCAAACCT
25 GCAAACAAACAGGAAGATGAAGTGATGAGAGCCTATTTACAACAGCTAAGGCAAGA
GACTGGACTGAGACTTTGTGAGAAAAGTTTTCGACCCTCAGAATGATAAACCAGC
AAGTGGNGGGCTTGCTTTGTGAAGAGACAGTTCATGAACAANAGTCTTTCAGGACC
TGGACAGTGAAGGGAGCCCCGGGCAGCCA

30 Sequence ID 719

CGNGGCCGCGTNAACTTTTGATCGTCAGCTGGGGCTGGCAGGCACCTAAATGGGAA
GGGTGATAGCAGTGTGTTGGGGGGAGTTTAGGGAACGGTCCTCTACCGATAGAGGC
AGCANCTCATTGGAATTTCTCCTGAAGTTGTCTTGCCCCCTGAATCCTGCAGGAA
GGCTGGCAAATGGCCATTTCCCTTCCACTTGAATAGAGACCCATAACTCAAGTATC
35 TGCCCTTAAGACACCACAGGACTGTTCTTCGCGGGCCCTGCCCTGGATTTGGGAG
AGGCAGTCCANCTCACCCAACCTAGGCTCTGCANGGGGACCANGAGGGATGGGTTGT
GTCCACAGGACCAGCCAGACTGATGAGGGATGCGGCAAGCATATTCTCACCACCTT

- 214 -

CTTTCACGTTTACAACANACCAGCNTTCCCTGTGTGGCAGGGTTACATTGGTCAC
CGAGGACCTANAATCATGGAGTGCTCTGGGGATCCGGGCTTGGA

Sequence ID 720

5 TCAGTGTGGAATTTTGTGTCAGACACTTTCTCTGCATCAATTGGTATGACCATGTGAT
TTTTTTTCTGTAGCCTGTAAATATGGTTAATTTTCAAATATTGAGCTGATTAATTT
TCAAATATTGAGCTCTCCTTGTCATCTCTGGAATAAGTACCACTTGGTCGTGGTATA
TATTTCTTTTAATATATTGCTGAATTCTGTTTGATCATGTTTTCTTAAAGACTTTC
GTGTCTGTTTTTCATGATAGATACTGGTCTATAGTTTTGTTGTAATATCTTGGTTTG
10 ATTTTGATATCAAGATAATGCTACCTTAATAGAATGAATTGGAGCCAAGTATGGTG
GCAAATGCCTATAGTCCTAGCTACTCAGGAGGCTGAGGTGGTGGGGACTGCTTGAC
CCANGAGTTCAAATCTAGCTTGGGCAATGTAGCAAGAC

Sequence ID 721

15 TAGAAGGAATGACTATTCATGTCCAAAGTGAATGGTTTTGTGCAGTGAACAACACA
TGGCGAGGTACTAACTGAGAACTTTTTTCATGCTTTATGCCTACCTCTTGTAGTTG
TTGCAGAGCAAATATAAATTGTAATAAGATAGCTAGGCCTTGCAGAAACAAACAGA
AAAACCTAAAAAAAATGATATAAGAGCTGGAGTCTAGTATTTATATGAATCTGTG
AGAGATAATTTTTTTTGGTCTCACTGCAATGAACCAAAAGCGGCTGAGTTTGGTTTT
20 TAATTGTAGCCATGTATTGAAGGCATCTTTTTGACCAACTCTTGTTGGTTCTGTCT
TGAACCATTGTTAATCACTGTGCTGTAATTAGTATAGCTAAATCTTTCCTTCCTT
GCTCCTCCCCCAGCCACCCCGTCTTCCCTTAACATTTTTTTCAGGGGGGGTTGGGA
GTGGTTTTCATTTTAATGTGAGTGGATGTTTTGATAGTTGTAAGGAAAAAATGCATT
TCAGACACATTTACACATGAGCTATTTTCTTACACAGTATGTCTTATTGGTAATA
25 AGAATGTAATTCAT

Sequence ID 722

CNTTCCNTAAGAATACAAAAAATTAGCTGGGCGTGGTGGCAGGCGCCTGTAATCCC
ATCTACTCAGGAAGCTGAGGCTGGAGAATCGCTTGAACCCGGGAGGCGGAGGTTGC
30 AGTGAGCAGAGATCACGCCACTGCAGTCCAGCCTGGGCAACAGTGCGAGACTCTGT
CTCAAAAAAAAATAAATAAATTACCTGGGTGTGGCAGCGCGTGCCTGTAATCCCA
GCTACCCAGGAGGCTGAGGCAAGAGAATGCTTGAACCCAGGAGGCAGAGGTTGCA
TGGAGCTGAGATGGCGCCACTGCACTCCAGTCTGGTGACAGAGTGAG

35 Sequence ID 724

CTCTCTACTAAAAATACAAAAATTAGCTGGGCACGGNGGTGCATGCCTGTAAACCC
AGCTACCAGGTACTCGGGAGGCTGAGGCAGGAGAATCGCTTGAACCAGGGAGTCCG

- 215 -

AGGTTGCGGCGAGCTGAGATCATGCCACTGCACTGCGGCCTGGAGACAAGAGCAAG
ACTCCGTCTCAAAAAAAAAAAAAAAAAAAAAAAAAAAGACNTCACCTAATTGCAG
NGNGNGGACCTTATTTGGCTNTTAATTCAAACCTATTAAAAATGTGAACN

5 Sequence ID - 726 nt: 260

CGGGGTCTGTACCGGGCTGGCCTGTGCCTATCACCTCTTATGCACACCTCCCACCC
CCTGTATTTCCACCCCTGGACTGGTGGCCCTGCCTTGGGGAAGGTCTCCCATGT
GCCTGCACCAGGAGACAGACAGAGAAGGCAGCAGGCGGCCTTTGTTGCTCAGCAAG
GGGCTCTGCCCTCCCTCCTTCCTTCTTGCTTCTCATAGCCCCGGTGTGCGGTGCAT
10 ACACCCCCACCTCCTGCAATAAAATAGTAGCATCGG

Sequence ID 727

CTGAGTNTAGAAATGATGCCATTAATACTGATTGCAAAAACATTACAACCTCAGTAC
TGCAGCTTTCATTCAAATAGGTTATATGTATAAACTGAGTTCAACAATATTGTATT
15 TGAGATGGTAAAGTTAAAGAAATGCAATAATGTAAATAATACTTAAGAAAATAAGA
TCTCAGGAACTGTATATACTCTGTACTTTTATGCAACTTTATCAGATCATTTTCAG
TATATGCATCAAGGATATAGTGTATATGACATGAACCTTGAGTGCAAAAACCTGTAC
TATGTACCTTTTGTATTATTTTGCTGTCAACATCTAAATAAAGGTTTTTTTGTGTTGT
TTTTTGTGTTTTTTAATTGTTTTGTTTTTAAAGATTGTTTTAATTAATTAATAAATA
20 ATTGTTTTAATTAACAATTGTTTTAATTGTTTTTAAAGTCGCCAGGCTGAGGCAGGT
GAATCACAAGCTTAGGAGTTGGAGGCTAGCCTGCCAACATGGTGAAACCCCGTCTC
TACTAAAAATACAAAAAATTAACCTGGGTGTGGG

Sequence ID 728

CCCATCTGCACCAGTACACAGGCAGGCATTATCATTCTTCACCTACTTTTTTAAATA
GTGGCAACTTGGGATTCATTCTGGTGATTCTGAACCTTGCCTCATAGCTTAAAGTA
TAAAAAAGATTCAAGAGCAGTGAGGTTTGTTCCTTCCAGTGAATGGTGGACTGAGT
GGTGCAGAGGTGGAGGGCTAACAAGAGGAAAGAACTACATTCTTCAGAATACAGTGA
TGAAAATTCATTTTGAACTCAAATATTTTCATTTTGGATATTCTCCTGTTTTTAT
30 TAAACCAGTGATTACACCTGGCCATCCCTCTAAATGTTCTAGGAAGGCATGTCTAT
TGTGATTTTGATGAAGACAGAATTATTTTTCTCTGTAGAAACACAGATACCACTTT
ATCAGGGGAAGTTAGTCAAATGAAATGGAAATTGGTAAATGGACAAAAGCTAGCTA
GTAAAAAGGACGACCCAGCAACATGCTTTAACCCCATGTATGTTTTGTGGAAAGAG
CATAGTTTAAACATCTTGAGAAATTGGGACATAAAAGTTTTTCATNGGTAGACAGTT
35 CATGGCAGTATATGAATTGACATAATGGAAATAATCTGATTTTATTTTTACAACCTA
ACATCCTTTCCCC

- 216 -

Sequence ID - 736

nt: 641

5 GGAATTCCAAGTGCTTGGGGATAATGATACCTCTGACCTTTCTTCCTTTTGGGAAG
TACTTGAGTGTGCAGCTGCATGAGGCCTCAGCAGGAGAGAGATTTTAGGTCCAAGA
AGCTATACCAGTAGGACAAGGCAGGAAAATACTACACTTTCAGGATCAAGCCCCCTC
TGACTCTCATTTGGAAACTGGATGTTTGCTAAGCACCTGCTTCTTAAGGATGCCGA
GGGATTTAATGATACTCCCAGAAACCTGGAGAGATTAATGGGGCCTATGGAGAAGT
GCTCTGAACTCAGTGTTGGGACTTGAATAAAATTAACCATTGTCATGTTTTCAGAA
CAACTAAGCTGTTTTATATTTTCATGTGCATGAAAGCCCTAGAACTAAGTTGTGTTA
TTTCCAGAAATGAAATAGATCCCACAGTTAGATGATGTGGCCATTAGGAAGTACCA
10 AATTTATAAAAAATCACTGGAGGTCTGTCTGAGCAGTACCTAATAAAATATAGTATA
CTGAAAGTGAACAGATACTTTGTCTCTTTCTTTGGCTGCTTGATCTTTATCTGTGT
CTGCCGTACAGTGCACCCTTAAAGTATTCTACACCAGTGCTTCTCAAACCTGGAAAT
GTGCATGTAAGTCACCCANGGGTCT

15 Sequence ID 739

TGCATGCCCATAGTCCCAGCTATTTGGGAGGCTGAGGCAGGAAAATCGCTTGAACC
CGGGAGCCAGAGGTTGCAGTGAGCCGAGATCGCACTCCAGCTTGGCGACAGAACAA
GACTCTGTCTCAAAAAAAAAAAAAAAAAAGAAATCTTGGGATCCTGAACCCCTTACTC
GAAGGGCTAAGGTAGCATCTCAGCATGTCTTATTCGAGACTTCGTANAACCAGACC
20 TGCTGTTTGTAGATGTTAATTAATCAAACCTTTCTCTACTCATTCTGGACCAGTTA
AGGTTTTCTCCTTCTCCGTATGAGTTTTGATTTTCGTCTCCTTGGTTGGAGATCA
CACTTTGGTCTGCTGCTAAGTTGGATGCCTCCCACTGTCTTTCCCTAAGTCTAGGG
CTTCANACCCCAGTGTGGGGAGAGGGACTTTCGTTTCCTGCCCCCTCACCACATCAG
ACACAGGCAGGCAAGAATAAGATGGCCAAAAGGCCGATGAACTTCTTGACCTAGCC
25 TGGGACATTACCTGTTACTAGGTGGACTTCACTGCCTGTGAATGGAAGCTGAAGGG
CTGTTTTTTTGGTTTGTATTTGGACAGGCCAGGCTTANAGAGGGAGAGAACTGGGC
TACTCTTCAGCAGTGATCTTTAAAATGCC

Sequence ID 747

30

CAGAGTGCAAGACGATGACTTGCAAAATGTGCGCAGCTGGAACGCAACATAGAGACC
ATCATCAACACCTTCCACCAATACTCTGTGAAGCTGGGGCACCAGACACCCTGAA
CCAGGGGGAATTCAAAGAGCTGGTGCGAAAAGATCTGCAAAATTTTCTCAAGAAGG
AGAATAAGAATGAAAAGGTCATAGAACACATCATGGAGGACCTGGACACAAATGCA
35 GACAAGCAGCTGAGCTTCGAGGAGTTCATCATGCTGATGGCGAGGCTAACCTGGGC
CTCCCACGAGAAGATGCACGAGGGTGACGAGGGCCCTGGCCACCACCATAAGCCAG
GCCTCGGGGAGGGCACCCCTAAGACCACAGTGGCCAAGATCACAGTGGCCACGGC

- 217 -

CACGGCCACAGTCATGGTGGCCACGGCCACAGCCACTAATCAGGAGGCCAGGCCAC
CCTGCCTCTACCCAACCAGGGGCCCCGGGGCCTGTTATGTCAAACGTCTTGGCTGT
GGGGCTAGGGGCTGGGGCCAAATAAAGTCTCTTTCCTC

5 Sequence ID - 757

nt: 583

GAACCCTGCGGAGGGACTTCAATCACATCAATGTAGAACTCAGCCTTCTTGGAAG
AAAAAAAAAGAGGCTCCGGGTTGACAAATGGTGGGGTAACAGAAAGGAACTGGCTAC
CGTTCGGACTATTTGTAGTCATGTACAGAACATGATCAAGGGTGTACACTGGGCT
TCCGTTACAAGATGAGGTCTGTGTATGCTCACTTCCCCATCAACGTTGTTATCCAG
10 GAGAATGGGTCTCTTGTGAAATCCGAAATTTCTTGGGTGAAAAATACATCCGCAG
GGTTCGGATGAGACCAGGTGTTGCTTGTTCAGTATCTCAAGCCCAGAAAGATGAAT
TAATCCTTGAAGGAAATGACATTGAGCTTGTTCAAATTCAGCGGCTTTGATTCAG
CAAGCCACAACAGTTAAAAACAAGGATATCAGGAAATTTTTGGATGGTATCTATGT
CTCTGAAAAGGAACTGTTTACGAGGCTGATGAATAAGATCTAAGAGTTACCTGGC
15 TACAGAAAGAAGATGCCAGATGACACTTAAGACCTACTTGTGATATTTAAATGATG
CAATAAAAGACCTATTGATTTGG

Sequence ID - 758

nt: 424

CTTGGCTCCTGTGGAGGCCTGCTGGGAACGGGACTTCTAAAAGGAACTATGTCTGG
20 AAGGCTGTGGTCCAAGGCCATTTTTGCTGGCTATAAGCGGGGTCTCCGGAACCAAA
GGGAGCACACAGCTCTTCTTAAATTTGAAGGTGTTTACGCCCAGATGAAACAGAA
TTCTATTTGGGCAAGAGATGCGCTTATGTATATAAAGCAAAGAACAACACAGTCAC
TCCTGGCGGCAAACCAACAAAACCAGAGTCATCTGGGGAAAAGTAACTCGGGCCC
ATGGAAACAGTGGCATGGTTCGTGCCAAATTCGAAGCAATCTTCTGCTAAGGCC
25 ATTGGACACAGAATCCGAGTGATGCTGTACCCCTCAAGGATTTAAACTAACGAAAA
ATCAATAAATAAATGTGGATTTGTGCTCTTGT

Sequence ID - 764

nt: 626

GATTTTTTTTTTTTTTTTTGAGATGGAGTCTTTCTCTGTCGCCCAGGCTGGAGTGCA
30 GTGGTGAAATCTCGACTCACTGCAACCTCCGTCTCCTGGGTTCAAGCAATTCTCCT
GCCTCAGCCTCCTGAGTAGCTGGGATTACAGGCACCAGCCACCACGCCCCGGCTAAT
TTTTGTATTTTTAGTAGAGACAGGTTTTTACCATGTTGGCTAGGCTGATTTTGAAC
TCATGACCCCAAGTGATCTGCCCGCCTCGGCCTCCCAAGTGCTGGAATTACAGGT
GTGAGCTACCACTCCCAGCCAATGATTACATTTTATAAGGTAAAATAACTTGTGCCA
35 ATCTGTACAAGTGAATTCAGATTTAAAATTTTAATTGTAAAAAGATATCCAGGTGA
TATTTCTCCCTGAATAATTTAGTTTCCTTTTCTATTTCTTGATATAAAAGTACTCA
GCATTGAAGTAATTGCTATCTTCACATTTCTTCCTATTTGAGCTGTCTAAATAAGT

- 218 -

AGTCCTACATATTTTCCCCCAACACAAAAAACCCAGAAAAGAATTATTTTATACT
GGATTTTTTTGGTTGTAGCAGGAACCTAAAGGNGCCAATTGTAACATGCATGTTCT
TTTTGGCAAA

5 Sequence ID 766

GTCCATCCTGCAGGCCACAAGCTCTGGATGAGGAACTTGAGGCAAGTCACCAGCCC
CTGATCATTTTCGCCTAAAAGAGCAAGGACTAGAGTTCCTGACCTCCAGGCCAGTCC
CTGATCCCTGACCTAATGTTATCGCGGAATGATGATATATGTATCTACGGGGGCCT
GGGGCTGGGCGGGCTCCTGCTTCTGGCAGTGGTCCTTCTGTCCGCCTGCCTGTGTT
10 GGCTGCATCGAAGAGTAAAGAGGCTGGAGAGGAGCTGGGCCCAGGGCTCCTCAGAG
CAGGAACTCCACTATGCATCTCTGCAGAGGCTGCCAGTGCCAGCAGTGAGGGACC
TGACCTCAGGGGCAGAGACAAGAGAGGCACCAAGGAGGATCCAAGAGCTGACTATG
CCTGCATTGCTGAGAACAAACCCACCTGAGCACCCAGACACCTTCCTCAACCCAG
GCGGGTGGACAGGGTCCCCCTGTGGTCCAGCCAGTAAAAACCATGGTCCCCCACT
15 TCTGTGTCTCAGTCCTCTCAGTCATCTCGAGCCTCCGTTCAAATGATCATCATCA
AACTTATGTGGCTTTTTGACCTTTGAATAGGGAATTTTTTAAATTTTTTAAAAA
TT

Sequence ID 768

20 CCAGCGCAGGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCGCAGATAA
GTTTTTTTCTCTTTGAAAGATAGAGATTAATACAACCTACTTAAAAAATATAGTCAA
TAGGTTACTAAGATATTGCTTAGCGTTAAGTTTTTAACGTAATTTTAATAGCTTAA
GATTTTAAGAGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAA
AAGGTTTCTAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAAT
25 GTATTTAAAGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACCAATAG
AAGGGCAATGCTTTTAGATTAAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTT
TAAAAGTTGTAGGTGATTAAAATAATTTGAAGGCGATCTTTTAAAAGAGATTAAA
CCGAAGGTGATTAAAAGACCTTGAAATCCATGACGCAGGGAGAATTGCGTCATTTA
AAGCCTAGTTAACGCATTTACTAAACGCAGACGAAAATGGAAAGATTAATTGGGAG
30 TGGTAGGATGAAACAATTTGGAGAAGATAGAAGTTT

Sequence ID 773

GAGGAAAGGGGAGTTAATATTTAGTGGACAGAATTTTCACTTTTACAGATGAAAAGA
GTTCTGGAGATAGACGGTGTTGATAGTTGCACAGCAGTGTGAATGTGCTCATTGTT
35 ACCGAACTTAAAAATGTTTAAACATAGTATTATGTGATTTTTTATTTTGCCACTTAAA
AAAAAGAATGAAGTACTGATACATGCTACAACATGGGTGAGCTTTAAATACATTC
TGCTCAGTGAAATAAGCCAGATGCAAAAGATCACATATTATATAATCCACTTATAC

- 219 -

GAGATACCTAGAATAGGCAAATTCATAGAGACAGAAAGTAGAATAGTGGTTCCCAG
GGGCTGGGGACAAGGGGGCAGTGAGAGATTGAGAGTTATTATTAATGCGTACAGAG
TTTCAGTTTGGGCTGATAAAAAAGTTCTGAAGATGGATGGTGATGATGGTTGTACA
TCAATGTGAGTGTAATTACCGCCACTGAACTGCCCTTAAAAACGTTTAAAAGAGTA
5 AATTTTATGTTGNGTATATTTTACCATAAT

Sequence ID 776

TTTTTTTTTTCATAAGAGGCAAGTACAAGAAAAAGCTTAATTACTTTAACTTCTAAG
TAGTTTGGGAATCTAAATAAATAGGAGTTACCAAATATATGCGCTTCTGTGAATAGT
10 TTTCCCCCACATGTTTATTTATATTTTGCATCTCATCAAACCTAACAGATTCTAA
AGTCTCTGGTGATAATGACAATATCTGCTACGGAGAGACTAGCCTGGGGGAAGAGG
ATCTCCCTGAACAAGGATAGCGGAGTTGCTGCAGCTTTCAAATGAAGCTGGACATT
TAGCTGCGGGGGTAGCACCTTTGATCAAGGCAGCCCAAAGATGAGTTTCAGGGAT
GGGACTGACAGAAGAGAAAAGTTCTTCCCAGCCCTTTCTACTTTTTCTCTTTGTTT
15 CTCAGGCTTCTGGCCGTCTTCAGTTTTTCAAGTTTCACTCTCAACCCTAAACAGT
ACTTCTGTGAAGTACCCTTTGGCCCCCTCGTTTTTCAGCTCCTAAACTCACCTGGAAA
TAGATGTCAATCTAATTTTGGGTCTGACTAGTGCAGTAGGCATTTTTTGGTGA

Sequence ID 782

CTCACACAGAACAAAAATGAATGAGTGTGGCTGTGTGCCACTATCACTGTGTCTAC
20 AAAACAGCCAGTGGGCCTGATTTGGCCCTTGGCTGCAGTGCGCCCGTCTCTGTTT
TTGAGGAATAAAATCGCATCATTTTCATATGGCTAATGCAATTTTTTTTCCCATCTGG
AAGCAACATCTGATTGGACTCATCTTGTATGGTGCTTGTTACAGTCTCTGTAAATG
GGAGAGGGTCCGAGAATAGCTCTTCCGTGTTTTTCATCAGGACTGTTTTTAGGGATGG
25 CAAAGAAGTCAGTGTGTCCAGCCTGTGTCCCTCCTCACCACGTGGCTGATTCTTGAA
TCTGCATGTGCANCACTGCCGTTGTCTGGGGCATGATCTGTGTGA

Sequence ID - 785

nt: 556

CTTTTCTCTGGGTATAGATTTACCCTAGCACCTATCTCATTATATTGAATTTTCCA
30 GCATATTTAAATAAACTATTAATTAGTCACACTATTTCTTAAAAGTCACACTATCA
ACTAATCGTGACCGCAATTATCTAGGGGTGATAATCTGCTGAGTCTACTCTTTAAA
TACACTGGGACCCAGCATATTGAGTTATATTGGCACAGAACTTCACTCTGGGTAT
AGATTTACCCTAGTACCTTGCCGGCAGGATCCTATTATTATCATGGTTGTACAAGCAA
GGTTCAGGGAAGAGGCTGGCACAGAGAAGGTACCTGGTAAGTGTGTTTGAGGCTG
35 AATTCAGCTCAACTCAGCTCCAGTAGAGATGGTGTCCCTTCTCTACCGTGTTGAG
ATAGTGTGCAGTCCCTTCCTAAGGGCTGTTACCCACCGCAATAGGACTTGTGAGCT
TCAACTTTTAAATTTCTCTGCTCCCGCTGGGACCCACCCGCTTCAAAAATCATCAT

- 220 -

GGNGGNNTT TAGCACCAATTTAGTAAACACAAACTGTCTGAAATATTTTGGAT

Sequence ID 796

GAACATTCAAGATAGTGAGAGGAAGAAAAAGATATGGCTGTACGGGACCGAGGTCT
5 CTTCTATTATCGCCTCCTCTTAGTTGGCATTGATGAAGTTAAGCGGATTCTGTGTA
GCCCTAAATCTGACCCTACTCTTGGA CTTTTGGAGGATCCGGCAGAAAGACCTGTG
AATAGCTGGGCCTCAGACTTCAACACACTGGTGCCAGTGTATGGCAAAGCCCACTG
GGCAACTATCTCTAAATGCCAGGGGGCAGAGCGTTGTGACCCAGAGCTTCCTAAAA
CTTCATCCTTTGCCGCATCAGGACCCTTGATTCCCTGAAGAGAACAAGGAGAGGGTA
10 CAAGAACTCCCTGATTCTGGAGCCCTCATGCTAGTCCCCAATCGCCAGCTTACTGC
TGATTATTTTGAGAAA CTTGGCTTAGCCTTAAAGTTGCTCATCAGCAAGTGTTC
CTTGGCGGGGAGAATCCATCCTGACACCCTCCAGATGGCTCTTCAAGTAGTGAAC
ATCCAGACCATCGCAATGAGTAGGGCTGGGTCTCGGCCATGGAAAGCATACCTCAG
TGCTCANGATGATACTGGCTGTCTGTTCTTAACAGAACTGCTATTGGAGCCTGGAA
15 ACTCAGAATGCAGATCTTTTGTGAACAAAATGAAGCAAGAACCGGAGACNCTGAAT
AGTTTTATTTCTGTATTAAAACTGNGATTGGAACAATTGAAGA

Sequence ID 801

CCACTCCACCTTACTACCAGACAACCTTAGCCAAACCATTTACCCAAATAAAGTAT
20 AGGCGATAGAAATTGAAACCTGGCGCAATAGATATAGTACCGCAAGGGAAAGATGA
AAAATTATAACCAAGCATAATATAGCAAGGACTAACCCCTATACCTTCTGCATAAT
GAATTAAGTAGAAATGAGGATTCTGACCTTGACTTTGATATCAGCAAATTGGAACA
GCAGAGCAAGGTGCAAAACACAGGACATGGAAAACCAAGAGAAAAGTCCATAATAG
ACGAGAAATTCTTCCAACCTCTCTGAAATGGAGGCTTATTTAGAAAACAGAGAAAAA
25 GAAGAGGAACGAAAAGATGATAATGATGATGAGTCAGGTAAAAGTTCCAGAAATGT
GAACAACAAAGATTTTTTTTGATCCAGTTGAAAGTGATGAAGACATAGCAAGTGATC
ATGATGATGAGCTGGGTTCAAACAAGATGATGAAATTGCTGAAGAAGAAGCAGAAG
AAGGAAGCATTCTGAAATATGAATGAAAAAATTACATCTTTAGAAAAGAGTTA
TTAGAAAAAAGCCTTGGCAGCCGTCNGGGGGAAGTGACGCACAGAAGAGACCAGAG
30 AATAGCTTCCTGGANGAGACCCTGCACTTTACCCATGCTGCTGGATGG

Sequence ID - 808

nt: 641

CCGGGTTTTAGTATTTAACCAAGAGCCTTTTAAATATTGAAAACCCATAGTTCAGA
AAATGTTAGTATTGCTGCCCTTCTTCACATAAATTTTTTTTTTAAATTATACTATTA
35 TTTTGCTTAATTTTATATTGGGTAAACAACCTTCAAGAAGGTTAACTAGGAAAG
AAGACCTTTTTGTTTTATTTTTACTATTTATATATAGAAGACAAATCAGCATTG
TGATAGTTTTACATGACCAGTTATCAAACGGTCATAGTATGAAGTGTGCAGTTGTT

- 221 -

CATTATTAGTAAATTATGTTTGATTTTTTAAACTATTTAGTACTAATAGTTGAGATG
AAAACCTGAAGAAAAATGCCAATGTGACGTTTGTGTATAGCTAGCCTTAAAAAACTT
CCCATGTTTTTAGGTGACTTTTTTCCCCCTCTTAGTACTCTGGAGAAACAATGAAG
ATGGGCCATCTCAATTCAGATGTAAACAAAAAGTAATTTTTTATTTCAACATTTAA
5 TGTAACCTGCTATTATTGNGGATTCTTGNCTTGNGTATTTTCTTTCCCTTATTCAAG
TAATATAGAATAACTTTCCTTAAATGATTTGATCCAAGATACGTCATTTCTGTAT
TGGCAAATGCCNCTATTAAAGTGT

Sequence ID - 814

nt: 132

10 GTTAAAGTGATACATTTTTTATACCAAATGTGTTTATTTTTTTGTGCAAGTAATCCT
TAAAATTGCAATTGTATTAGGTGTTAAAATAAAGTTTTTAAAAAATTAAAAAAA
AAAAAAAAAAAAAAAAAAAAA

Sequence ID 817

15 GACAACCTTAGCCAAACCATTTTACCCAAATAAAGTATAGGCGATAGAAATTGAAAC
CTGGCGCAATAGATATAGTACCGTAAGGGAAAGATGAAAAATTATAACCAAGCATA
ATATAGCAAGGACTAACCCTTATACCTTCTGCATAATGAATTAACTAGAAATAACT
TTGCAAGGAGAGCCAAAGCTAAGACCCCCGAAACCAGACGAGCTACCTAAGAACAG
CTAAAAGAGCACACCCGTCTATGTAGCAAAATAGTGGGAAGATTTATAGGTAGAGG
20 CGACAAACCTACCGAGCCTGGTGATAGCTGGTTGTCCAAGATAGAATCTTAGTTCA
ACTTTAAATTTGCCCACAGAACCCTCTAAATCCCCTTGTAATTTAACTGTTAGTC
CAAAGAGGAACAGCTCTTTGGACACTAGGAAAAAACCTTGTAGAGAGAGTAAAAAA
TTTAACACCCATAGTAGGCCTAAAAGCAGCCACCAATTAAGAAAGCGTTCAAGCTC
AACACCCACTACCTAAAAAAATCCCAAACATATAACTGAACTCCTCACACCCAATT
25 GGACCAATCTATCACCTATAGAAGACTAATGTTAGTATAAGTAACATGAAAACAT
TCTTCTNCGCATAAGCCTGCGTCAGATTAAAACACTGAACTGACAATTAA

Sequence ID - 821

nt: 370

30 AAAGAGCTCCCAAATGCTATATCTATTTCAGGGGCTCTCAAGAACAATGGAATATCA
TCCTGATTTANAAAATTTGGATGAAGATGGATATACTCAATTACACTTCGACTCTC
AAAGCAATACCAGGATAGCTGTTGTTTCANAGAAAGGATCGTGTGCTGCATCTCCT
CCTTGGCGCCTCATTGCTGTAATTTGGGAATCCTATGCTTGGTAATACTGGTGAT
AGCTGTGGTCCTGGGTACCATGGCTGGTTTCAAAGCTGTGGAATTCAAAGGATAAA
TTAATGAAGAAAACAAGCGGAGCTGAAGAAGAAAGTACAATATGGTGCTGTCTTCC
35 TAATGAAATAAATTCACTAAATGGACATTAAAAA

Sequence ID 825

- 222 -

AGACTCGAGCAAGCTTATGCATGCATGCGGCCGCAATTCGAGCTCGGCCACTTGGC
CAATTCGCCCTATAGTGAGTCGTATTACAATTCAGTGGCCGTCGTTTTACAACGTC
GTGACTGGGAAAACCTTGGCGTTACCCAACCTTAATCGCCTTGACAGCACATCCCCCT
TTCGCCAGCTGGCGTAATAGCGAAGAGGCCCGCACCGATCGCCCTTCCCAACAGTT
5 GCGCAGCCTGAATGGCGAATGGAAATTGTAAGCGTTAATATTTTGTAAATTCGC
GTTAAATTTTTGTAAATCAGCTCATTTTTTAACCAATAGGCCGAAATCGGCAAAA
TCCCTTATAAATCAAAGAATAGACCGAGATAGGGTTGAGTGTTGTTCCAGTTTGG
AACAGAGTCCACTATTAAAGAACGTGGACTCCAACGTCAAAGGGCGAAAAACCGT
CTATCAGGGCGATGGCCCACTACGTGAACCATCACCTAATCAAGTTTTTTGGGGT
10 CGAGGTGCCGTAAGCACTAAATCGGAACCCTAAAGGGAGCCCCGATTTAAAGCT
TGACGGGGAAAGCCGGCGAACGTGGCGAGAAAGGAAGGGAAAAAGCCAAANGGAG
CCGGCGCTAGGGCCTGGCAAGTGACGGGCACGCTGCGCGTAACCACCCACACCCC
GCCGNGCTTAATGCCCCNTTCAGGGCGCGTNCCTGATGCCGNATTTTNTCTTACNCA
TNTGTGCNGGNTT

15

Sequence ID 833

TAAAATAATGGCAAAAAACAAACAAAAACAAGTTCTCTAAACAGAAAGGAAATTA
CTAAAGAAGGAATCTTGAAATAACAGGAAAGAGGAAATACCACAGTAGGCAACATT
ATGGGTAAATAAAACAGACTTTCCTTCTTTAGTTTCTTAAATATGTTTGATGATT
20 AATGCAAAAATTACAATATTTTCTTATGTAGCACTAAAGGTATGTAGAGAAAATAT
TTAAGATAATTGTACTGTAAGCGGGAGATGACAGTGACATAAAGGCAACGTTTTTA
TACTTCACTCAAACTTTATGTATTAATGTAATCCATAAAGCAACCAAAAAAGCTAT
ACTAAGTACATTCAAAAACACAATAGATAAACCAACAAATTTCTAAAGGATGTAC
AAGTAACCCACTGGAAGCTGCAAAAAATGTAAACAGAACTAAAAACAGAGAATAA
25 ATGAAAAATTAAAAACGAAATGGCAGACTTAGGCCCTAATATACAAATTATCACAT
TAAATATAAATGGTCTAAATACACCAACTGTAAGACAGAGATTAGCAAAGTCGATT
TAAAAACATGACTCAACTACGTGCTGTCTACAAGAACTCACTTCAAATATACCAA
GATAGGAAGGTTGAAAGTAAAACGATGGAAAAAGATGTATCATGTGAACATTAATC
AAAGGAAAGCAGGGGTGGCTATATTAAACATCAGGTAAAATAAACTTT

30

Sequence ID - 837

nt: 603

TGAGGNTGGTCATGATGCANAAGCTACTCAAATGCAGTCGGCTTGTCCTGGCTCTT
GCCCTCATCCTGGTTCTGGAATCCTCAGTTCAAGGTATCCTACGCGGAGAGCCAG
GTACCAATGGGTGCGCTGCAATCCAGACAGTAATTCTGCAAACCTGCCTTGAAGAAA
35 AAGGACCAATGTTCGAACTACTTCCAGGTGAATCCAACAAGATCCCCCGTCTGAGG
ACTGACCTTTTTCCAAAGACGAGAATCCAGGACTTGAATCGTATCTTCCCACTTTC
TGAGGACTACTCTGGATCAGGCTTCGGCTCCGGCTCCGGCTCTGGATCAGGATCTG

- 223 -

GGAGTGGCTTCCTAACGGAAATGGAACAGGATTACCAACTAGTAGACGAAAGTGAT
GCTTTCCATGACAACCTTAGGTCTCTTGACAGGAATCTGCCCTCAGACAGCCAGGA
CTTGGGTCAACATGGATTAGAAGAGGATTTTATGTTATAAAAGAGGATTTTCCCAC
CTTGACACCAGGCAATGTAGTTAGCATATTTTATGTACCATGGNTATATGATTAAT
5 CTTGGGACAAAGAATTTTATAGAAATTTTAAACATCTGAAAA

Sequence ID - 839

nt: 71

ATTTATCTAATATTTGGTTTAATAAAATGTGAATAATGAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAA

10

Squence 849

nt: 622

15

TTTTTTTTTATTTTTTTGAGAATGGAGTCTTGCTCTGCCGTCCAGGCTAGAGTTCAG
TGGTGCGATCTCAGCTCACTGCCACCTCACCTCCTAGGTTCCAGAGATTCTTGTC
TTCAGCCTCCTCAGTAGTTGAGAATACAGGAACACGCCACCACGCCTAGCTAATTT
TTGTATTTTTAGTAGAGATGGGGTTTCACCATGTTGGCCAGGCTGGTCTCAAATCTC
CTGGCCTAAGTGACCCACCTGCCTCAGCCTCCCAAAGTGCTGGGATTATAGGCGTG
AGTCATTGTCCCCAGCCGGATGTTTTTCATCTTGATTTGCCTTAGTTTCTAAATCTC
ATCCTCTCCATTTTCTCCTGTTAGTAGTCACAGAGAACCAAATTCTGTCAAGTTAT
GAAACTAAAGTCTCTCTTCCACAAGTCTTCCTGTGTTCTGCCTCAAGTGAACCTGA
20 AAGAACATCAGTTTGTGGGAAGGTTGAAGACCGAATGATCTGCTGGGAAATCACTG
AGGCATTGCCATTCTCTTGAGGAATTTCATTTTCATCGAAGTTTCGGTTTATATCC
CTTTCTTGGTGAGTACTATTGCTGTTATGTAAATTAAATGAGTCGTCATCCTTCTT
NTGAGC

25

Sequence ID - 860

nt: 501

30

GTGAAATCACTTTCATGGATTATTAATGGATTTAAGAGGGCATCAATCAGCTCAAC
TCAAGATTTTCATAATCATTTTTTAGTATTTAGATTGTGCCTCAAAGTTGTAGTACCT
CACAATACCTCCACTGGTTTCCTGTTGTAAAAACCTTCAGTGAGTTTGACCATTGT
GCTCTTGGCTCTTGGGCTGGAGTACCGTGGTGAGGGAGTAAACACTAGAAGTCTTT
AGTACAAAACCTGCTCTAGGGACACCTGGTGATTCCTACACAAGTGATGTTTATATT
TCTCATAAAGAGTCTTCCCTATCCCAAGGTCTTCATGATGCCAGTAGCCATATATG
ATAAATTATGTTTCAGTGATAACTTAGTTATCAGAAATCAGCTCAGTGGTCTTCCCC
GCCATGATTCACATTTGATGAGTTTTTAAAAATCAAAGTGATTTTGAAAATCTCTA
ATGGCTCAGAAAATAAAAACATCCAGTTTGTGGATGACTATATTTAGATTTCT

35

Sequence ID 864

TTGTGTTTTTAGGACTCCTTATCTAAATTAAGGCAGAGAAGTTACAGTATTTATAT

- 224 -

CTGCATTAAATCTCAATTCCAGAAAAACCTTTTGAAAAATTATTTAATCCTCTGGA
AACTATTGATATGATACAGGAGAAATTTTCAGAAGTTTATTGAATAATTTAATATC
ATTTAATAGGACACTCTGGCTTGTATATAAGCAGATACGTTACTCAGACTTCTTGG
CTGTACTCTAAAATAATATATGTACTAGTCTCCTAAATATTACTAGCTCACCTTTC
5 AAAATGCATACTAATATTTCAATGTCTTTCTTCAATTTGAAAAGCTCTTGAATATC
TACTTGTGATAGCCCTAAGAGCTGAGATAATTATTTCCAGGAGGTTGAATCCCTGA
TTCTTAACCTGTTTCAGCAATGCATAAGCAAGAGAGAATATGACATAAGAGGACCATT
TCTACATTAGCCATTTTTTTTTTCAACAAGATACCTATGTGAATACAGGGCACCTGGGA
GGGTAAGTGGAGGACTATTTCTAACTATATTTATAAGCACATACTGATATTGGTGA
10 ATCAAAACCTACAGCAGTGCTTCTCAGATGGGAAGGGAGACAATGTGTAAGGAGAT
CAGGAATTCATTAG

Sequence ID - 865

nt: 122

CCANAATCCACTCTCCAGTCTCCCTCCCCTGACTCCCTCTGCTGTCTCCCTCTC
15 ACGAGAATAAAGTGTCAAGCAAGAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAA

Sequence ID 867

TTTTTTTTTTTTTTTTTTTTTTCAGAGTCACAGATATTGTATAGCTGAGGTAAGCATTT
20 TACAACCTTTTCAGACACAAGTAAGTACATAAATATTATTTTACAACCAACAATNTT
TAATATTTCCACATTGAANAATAGATGTGATAATTAAATCTTTTATAAGGTTTTAA
AAAGACATGAAACATAAACCTAATTATACATAAAAGAAAAGAATTTTAAACAAGAG
CTTATTGNGATGACATTACTCATAACTTTTACCTTTAAAACCTTTTCTTGGGTAGC
TATTCAAAAGTAAAGACCACAAGTTTTGTTGCCCANATTTCTTATGTTTNGTATAT
25 TTAAGCTCTTTATTTATTGAACAGATGNGTCATTAAATTCATTNGGAGCATTACTAT
TATCAGTAAATTTGATTTTTTTTTTCCCCTCAGTCATAGGTAAATCAGCTCCACCT
GGAATTTCTAAGGACCCAGTTTTAGTCAATATTTTCAAGTAATCATGACCTCAGAA
ATAGTCTTAATTAAGATAACAAATATTAGCCATCAAAATGGAACCAAGACAAGATT
CTAATGTTTGTAACAGTCAATCCATATTTATGAATATTAGCATATATTGGNGAAT
30 AGTTAAGGCAAAAGGGTCTAGCAG

Sequence ID - 869

nt: 667

TTGTGTTTTTTAGGACTCCTTATCTAAATTAAGGCAGAGAAGTTACAGTATTTATAT
CTGCATTAAATCTCAATTCCAGAAAAACCTTTTGAAAAATTATTTAATCCTCTGGA
35 AACTATTGATATGATACAGGAGAAATTTTCAGAAGTTTATTGAATAATTTAATATC
ATTTAATAGGACACTCTGGCTTGTATATAAGCAGATACGTTACTCAGACTTCTTGG
CTGTACTCTAAAATAATATATGTACTAGTCTCCTAAATATTACTAGCTCACCTTTC

- 225 -

AAAATGCATACTAATATTTCAATGTCTTTCTTCAATTTGAAAAGCTCTTGAATATC
TACTTGTGATAGCCCTAAGAGCTGAGATAATTATTTCCAGGAGGTTGAATCCCTGA
TTCTTAACTGTTTCAAGCAATGCATAAGCAAGAGAGAATATGACATAAGAGGACCATT
TCTACATTAGCCATTTTTTTTTTTCACAAGATACCTATGTGAATACAGGGCACCTGGGA
5 NGGTAAGTGGAGGACTATTTCTAACTATATTTATAAGCACATACTGATATTGNTGA
ATCAAAACCTACAGCAGTGCTTCTCAGATGGGAAGGGAGACAATGTGTAAGGAGAT
CAGGAATTCATTAGTCACCTTTTCAGATGGTTTAATGCATACAGCTGTACCG

Sequence ID 870

10 GGAGTTTGAGCAGATCCTTCAGGAGCGGAATGAACTCAAAGCCAAAGTGTTTCTGC
TCAAGGAGGAACTGGCCTACTTCCAGCGGGAGCTGCTCACAGACCACCGGGTCCCC
GGCCTTCTGCTCGAGGCCATGAAGGTGGCTGTCCGGAAGCAGCGGAAGAAGATCAA
GGCCAAGATGTTAGGGACACCAGAGGAAGCAGAGAGCAGTGAGGATGAGGCTGGCC
CATGGATCCTGCTCTCCGATGACAAGGGAGACCATCCCCACCCCCGGAGTCCAAA
15 ATACAGAGTTTCTTTGGCCTATGGTATCGGGGTAAAGCTGAATCCTCTGAGGATGA
GACCAGCAGCCCTGCACCCAGCAAGCTAGGGGGAGAAGAGGAGGCCCAACCACAGT
CTCCAGCTCCTGATCCGCCCTGTTCTGCCCTCCACGAACACCTTTGTCTGGGGGCC
TCAGCCGCCCCAGAGGCCTGACTTAGGGGTCTGGCTGTGGAAGGATGTGTGGCCTC
AAATGAGGACAGGGCTCCCGCCTTCACAGCCCTCGCCAGGGGTCTGCCCAATCCT
20 GGCCTGCATCAGGCAAGGACGGGGTCTCAGC

Sequence ID - 871

nt: 642

GCAAGTCTTCAGTATGTACATTTATCCCCTAGAAGAAGAAAAATTAGTTGTGCATG
AAAAAGAAACATTAAGTCAAGCTAAATGCTCACACTCTAAATCAGTGCTCTCCA
25 AAGTACAGCAGGCGGGAAAAGAAAATGGTAGATTTTTTTTCTTCCAATTACTTTAAC
TTATTCTTTTTTAATGGACACTTCATACATAAATATATTCACAATATATTAATATAT
ACATAATGTATAAGCATACATATTGAATGTGCAGTCAAAAAATGTACTAATGGAAT
GCTCTACCAAAACAAGTTCACGTTTCTGTAAAATGGGAATAATATTTTTTAAAG
GCATACAGTCTGAACATTTTTTAGATTATTCATAAAATCTATTCAGAAAGTTAAACT
30 AAAAAATTTAACGTATGCCTATAACAAATTTTGTACTTAATGTAATTGNTTTTCAT
CCTGAGATCTAATATCCTCGTTTTTAAGTAGAGCCACTTGTTTGCTACAGTTTAGT
CAAAACGTTAACATTAGATGGGTAAAGTAATATGAAATCTTTCTACTACTCCAAA
TAGAAAACAGAACATTAATAAGATAAAAATTCAAACATACTTACCAGTAGATTTTC
AACTGNGCAAAAGCTCATTGCATGGG
35

Sequence ID 873

GTTTTCCACCGTGAAGAGAACATTTCTCTGGGAATGACAAAGCCCTCAGGAACNG

- 226 -

CTTTTATTTCTATTGGAAGATGCCCATCATACTTCTGGCAGGATAAAATGATAAAT
TTATTTATTCAACAGATGATACTCAATTCCCTGCTGTTTTACTAAAGGTTCTTTAC
GTTTTATAGAAGCTAAATTTACTGTCATAGAAATTGCAATTGTAGATGTTACTGTA
ATCTAGTCAGAATATCCTTATCCTTCTAAAATAAACTAGTTAAAATTATTAACAT
5 ACGTACTGATATTAATTTTTAAGTTTAATGCTGCCACGTGCTTCTGCTAAGAACAT
TTATCACTACAAGTGGCAGAAAATTCCAACTCATCAAAACCAAAGTGTGCTTCT
TCCCTGCTTTTTTCAGAAAATGAGAAAGGATGACTTTATTCCAACATATTCTAAAAG
TATTCCAAGAACACTACCTTTATTCTAAATTCGTTATTTTCACAAAATAAAGGCTG
CAGATTGAAAGATAAAGGATTGCTATTAAAGAACAAGAAACAAAACCGAGAGA
10 GAAGGAGAGCTAGGGAAATCCCTGCANAANAACCGAATANGGTCCCTCTATTCTGG
GCCGGGGCCTGAAACTATGAAACAGGCCAACACAGAATCTTGGCA

Sequence ID 875

CCTCTGACTCGCTCAGCTCACCCACGCTGCTGGCCCTGTGAGGGGGCAGGGAAGGG
15 GAGGCAGCCGGCACCCACAAGTGCCACTGCCCGAGCTGGTGCAATTACAGAGAGGAG
AAACACATCTTCCCTAGAGGGTTCCCTGTANACCTAGGGAGGACCTTATCTGTGCGT
GAAACACACCAGGCTGTGGGCCTCAAGGACTTGAAAGCATCCATGTGTGGACTCAA
GTCCTTACCTCTTCCGGAGATGTAGCAAAACGCATGGAGTGTGTATTGTTCCAGT
GACACTTCANAGAGCTGGTAGTTAGTAGCATGTTGAGCCAGGCCCTGGGTCTGTGTC
20 TCTTTTCTCTTTCTCCTTAGTCTTCTCATAGCATTAACTAATCTATTGGGTTTCATT
ATTGGAATTAACCTGGTGCTGGATATTTTCAAATTGTATCTAGTGCAGCTGATTTT
AACATAACTACTGTGTTCTGGCAATAGTGTGTTCTGATTAGAAATGACCAATAT
TATACTAAGAAAAGATACGACTTTATTTTCTGGTAGATAGAAATAAATAGCTATAT
CCATGTACTGNAGTTTTTCTTCAACATCAATGGTCATTGNAATGTTACTGATCATG
25 CATTGGTGAGGNGGTCTGAATGTTCTGACATTAAACAATTTTCCAT

Sequence ID - 876

nt: 115

AAACTTTTGTGGCAACAGTGCCTAATTTGGATAATGTTTGTTCCTCAATAAATTAA
GAGCCAAATTGTAAA
30 AAA

Sequence ID - 878

nt: 634

GCCAGGCTTTGTGAATTACAGGACATTTGAGACAATCGTGAAACAGCAAATCAAGG
CACTGGAAGAGCCGGCTGTGGATATGCTACACACCGTGACGGATATGGTCCGGCTT
35 GCTTTCACAGATGTTTCGATAAAAAATTTTGAAGAGTTTTTTAACCTCCACAGAAC
CGCCAAGTCCAAAATTGAAGACATTAGAGCAGAACAGAGAGAGAAGGTGAGAAGC
TGATCCGCCTCCACTTCCAGATGGAACAGATTGTCTACTGCCAGGACCAGGTATAC

- 227 -

AGGGGTGCATTGCAGAAGGTCAGAGAGAAGGAGCTGGAAGAAGAAAAGAAGAAGAA
ATCCTGGGATTTTGGGGCTTTCCAATCCAGCTCGGCAACAGACTCTTCCATGGAGG
AGATCTTTCAGCACCTGATGGCCTATCACCAGGAGGCCAGCAAGCGCATCTCCAGC
CACATCCCTTTGATCATCCAGTTCTTCATGCTCCAGACGTACGGCCAGCAGCTTCA
5 AAAGGCCATGCTGCAGCTCCTGCAGGGACAAGGACACCTACAGCTGGCTCCTGAAG
GAGCGGAGCGACACCAGCGACAAGCGGAAGTTNCTGAAGGAGCGGCTTGCACGGCT
GACGCAGGCTCGGCGCCG

Sequence ID 879

10 GTTGCCGGGTCCTGTGATAACTCTGTTTAACATTTTGAGGAACTGTTGAATGGTTT
TTCACAGCAGCTGCCTCATTTTTTATTCCCATCAGCAGTACTTCTTGGTTCTAATA
CCTCCACGTTCTCGCCAACACTTGTTGTTGTCTGTAATTCGTTGTTAGCCATCCC
AGTGGGGATGAAGTAGTATCTTACTGTGGTTTTTCAGTTGCGTTTTCCCTGATAATTA
ATGATGGTGAACATCTTTTCATGTTCTTGTGGCCATTTGTATGTCTTCTTGGGAA
15 AAAAAAATGTCTGTTCAAATCCTTTACAAAGTATTTATTTTTTATGTCAACAATA
TAACCACTCAGTACACTGCTTTTTANACAATGATCTTTTAAAGGTTTGTTTACAAC
ATTTAGCACTTGAAATTTTAAGGTTATGCCCTCAAAAAAATTGCTGAGGGAGCTAA
GCTATGAAGATGCAAAGGCATAANAATTATACAATGGACTTTGGGGGAATCCAGGG
AAAGGGTGGGAGGGGGGTGANGGA

20

Sequence ID 881

TCGACTCTGATTTTTTTTTCTCCTTCCTCGCAGCCGCGCCAGGGAGCTCGCGGNGC
GCGGCCCCGTGCTCCTCCGGCCCCGAGATGAATCCTGCGGCAGAAGCCGAGTTCAACAT
CCTCCTGGCCACCGACTCCTACAAGGTTACTCACTATAAACAATATCCACCCAACA
25 CAAGCAAAGTTTATTCTACTTTGAATGCCGTGAAAAGAAGACAGAAAACCTCCAAA
TTAAGGAAGGTGAAATATGAGGAAACAGTATTTTATGGGTTGCAGTACATTCTTAA
TAAGTACTTAAAAGGTAAAGTAGTAACCAAAGAGAAAATCCAGGAAGCCAAAGATG
TCTACAAAGAACATTTCCAAGATGATGTCTTTAATGAAAAGGGATGGAACCTACATT
CTTGAGAAGTATGATGGGCATCTTCCAATANAAATAAAAGCTGTTCTGAGGGCTT
30 TGTCAATCCCAGAGGAAATGTTCTCTTCACGGTGGAAAACACAGATCCAGAGTGTT
ACTGGCTTACAAATTGGATTGAGACTATTCTTGTTTCAGTCCTGGTATCCAATCACA
GTGGCCACAAATT

Sequence ID 883

35 TCATTTACATTAAATACTCAAACTGCTCGATTAAGCAGGTGCTGTTCTTATCGCCA
TTTTGCATATGATGAGAAAGGGTAAGGTCAACCAGCTAGTATTTGGCTCACAGCAG
GCCTTAAGACTTGGTTTGTGTGACTCATCAGTCCACGCTCCTAAAACCACTAAGTT

- 228 -

5 GTTCTACCCTTTAATGTTGAATTAACATTGGATAGTGTTCAGTTTANATGGGTGG
GTGAGGGCCCAAGGACCTTTCAAACCTCAGATCTCTTATTTAATAACCTGGTCCCAG
ATCCATTCTCTGTCTGAAGAGGAAGTCATCCTTCAGTGGCTATTCATTGTGGGGTT
AAGAGCGCAGACTATGAATTCAGTCTTTTTTGGGTCCCAGTTTGCCAGACCTTGAGT
10 GAGTGCCCCGAGTTTACTTACTTGTAAAGGTAGGTGGAGGTAATATAATTAAATAA
ACTTAAAAAACTAATTAAAAACAAAACAAATGAACTAAGGTCTTAGGATATCTGGC
GTCTATTTTTCGCCAAATCACATAATGTCTATTGTTGTGTGTTGGACTATAGGATT
GTCCTTTAACAGGGAAGGGTTTATTTCTGTAATCAAGTCTGTCAATATTATGACCA
TGTTGATAATAGCTACCTTTAATTGAGGGCTTCCATGTGCCAA

Sequence ID 885

15 TCAGTGGAAAAGGGCAGGTTGAATCAAGGTGAATCAATCTGAAATTGAGCACACCT
GCCTGCCATCGCTGTTCCCTTCAACTGAGTGTCTGCACATCATGGGCTCTGTCTGTGA
GAGAAAAATCCCGGTGCTTGGTGTCTTGCATGACATGGAGTTTTCATGTAGATC
AATTTAAATGTACCTCTTGTTTACATAATTTGCATAATTTTAAAAGATAATGTTG
20 CCAAACCTTTGGAAATGTTAATGTTCANACTGAAAATCTCCACTACATGTAACCTTC
TTCCTCTGGATCAGTGGCATGGCTTATAATCCCAGCCAGTGGTTTGAACCTGTTCCA
GTGTCAACTGCCATGTGCTCTGCTTCAAGGGGGAAC TAGCCTTTTGTGAATTTTTT
GTACATAAGTATTTGTTACAAATATTTTAGCAAATGCTTTCTATTTCTCTTGCTTG
TGCATATCTTGGCTGGCGTTACAGAAAAATAGTGTAAACATTATTTCTTACCGGG
GAATGAGGGTTTT

Sequence ID 887

25 AGCACCTGGCACAGAGTAGTAGCTAACACAGATGTTAATTTTGCTGCGTCAAATGT
TTTCACTTTGAATCTCTCTTGAGTATTGTTCTCCTTATTGATTACATGATGACATC
CTGTTTTCTCTCCCTGACCTTTACTGTTTGTAGAAAAAAAAAAAAAAAAAAAA
AAAAAA

Sequence ID 889

30 CAGAGAGCTTGTTCCCTCCCTCCCTGTGCATGCAAACAAGAGGGCATGGGAGCACA
CAGAGAGATGGCAGCCACCTACAAGCCAAGAGGAGAAGCCTCACAATCAAACCTCTC
GCTGCTGGCGAGAGTCTTGGACTCTGTCTTGGACTTCCAGCCTCCAGACTGTGAGA
AACAAATTTCTGTTGTTTCAGCTTCTCAGTCTCTGGTGTTTTGTATTGTCAGCCTG
AGAACACAGCTGTACNATTATNAGGGAAACAGAAAACACTGATACTTAACAATGCT
35 AATGCAATTATTTATTTGCTTTTTCAGTCTCTACAAAACGTTCTAAAACACTAATCT
AAATATTAACAGTAAAATATTTGCATAACTAATGGAACTAAGAAATCATATGACC
AATATTTCACTTATTGGTAATCTTACTCTACTGATTTCCCCCAGACTGTGATTTT

- 229 -

TGAACTTCCTTGCCTTTCTCCTGTCTTTCTGNGTTTATTTCATGGAATTCCAGTTAT
CTGGGCTTGAAATTGCAGGCTCTCCTAACTTAAGCAAAATCTGACAGATCAGCAAA
ATGAGATAAATGTTTCTTTTTTCTTTCTGACTGCATTAAATCAGATACAACCTCAGC
ATTAAAAAGCTATCTTTGNAAAATGNTGGTACTAATAAATTAGTCTTA

5

Sequence ID 890

CCAGTTCACATTCAGTGAAGTCATGAACTTGAAATTGGCCATGATCAAAAAGTAT
TTAAATCACAGAAGTTGCAAATGCCACAAATCAAGGTCTTTTTCTCTTGAGAACC
TGTTAAACATTTACCAACTCACGACCGCCATGCACCCAATACTGCAATAGGTCTAT
10 AGATGCAGATACTGTCTCCATGAATCTTATAGGCTAGAAAGGAAATAGATAAGTAG
TCCTACCAGAAGACATGATGAAGGCATTTGTGGTAAACAGAATGATGGCCCCCA
AAGATGTCCACATCCTAATCCCTGAAGCCTATGAATATACTACTTTACTTGGA
AGGGACTTTGCCACAGGTTTTTAATTAAGGACCTTGAAATAGAGAGATTATCCTGG
ATAATCCAGATGGCCCCAGTGTAATCCCAAGGGTCCTCACAAAGGGTAGGAAGGAG
15 AGCCAGAGTCAGAGAAGGAGACGTAGCAATGGAGGCAGAGGTCANAGAGAGATCTG
CAGATGCTGCTGTGTTGGCTTTGAAAATGAGGAATGCAGGTGACCTCAANGNGCTA
GATGATGCAAGGAAACAAATAATCTCCTATGAACCCTAGGATGGGCATTATTATGA
GTCCTATTTTATAAACAAGGAACTGACNTCCAGAAAGATAAATGC

20

Sequence ID - 891

nt: 626

GGCAGAGGTTGCAGTGAAGTGAAGATCATGCCATTGCAATCCAGCCTGGGCAACANG
AGTGAGACTCCATCTCAAAAAAAAAAAAAAAAAAGACAAGAGTNTCCACTCTAAACA
CTTNTATTCAACATAGTCCTGAAAGTCGTAGCCACAGCAATTTAACAAGATAAAGC
AATAAAATGTATTCAAATAGAAAAAGAGGAAGTCAAATTATCTTCACTGGNGATAT
25 AATTCTCTACCTGGGAAACTTCACCGAAAAAGATTTACCAAAGATTTCTAAGCC
TAAATAATGACTTCAGCAAAGTCTCACCATACAAAATCAACATACACAAATGAGTA
GCATTTCTGTGCACCAATAATATTCAAGCTGAGAAAAAAGAACATGGTTCTATTT
ACAATAGCTACAAACAAAAAATATGTACCTAGTAATACATTAAATCAAGGNGGTA
AAATATCTNTACAACAAGAACTACAAACTGCTGAAAAAATAGAGACACGCAAA
30 TAAGTAAAAGGCACTCCATGCTCATGAATTTAAAGAATCAATATAATTAAATGT
CCGNGCTGCCTAAAGCAACTTACAGATTAAAGGCTATTTCTCTCAAACATATAAATG
CACCTTTTTTA

35

Sequence ID - 893

nt: 585

GTCATTGCTGGGTGGCGCCAGCCCTCAGACTTGCCCTCTTTCAGTAGGAAGAAGGC
CTCCCCACATACCTTCCCACACTCATCACCTTAAGCCAGACTCGGTGTCCAGTGAA
TATGACCATCTCTTGCCCATTTTCTAATGAGTGTTTTTATTAAATGAGTTATAAGAA

- 230 -

TGTGGTGGGTAAATCTATGGGCTTTGAACTAGTGAATCAACTTGGTTTCAGAATCT
GGCACTGCTACTTACTAGTGAATTTAAGCAAGTTATTTACCTTTCAGAGTGTGAG
TTCCCTCATGCATACAAGGAAGATAAAAAATAATGTNTACNAAAGTATTGGAGTAA
TTAATACATGGAGAAGTACATGTAAAGCGTTTAGCATGATGTCTGACATATTAAGC
5 ATCCAATATTAGTNGCTTGCAGAATTATTAGTAAAAGAGATTGCTTCTGAAAGCCA
TTCCAATTCTTAAATTTTATAATGCCACATTTGAGGTCACCTGAAGTCGTGTATAA
CATGTGTACATTTTTCGATTTATTTTTTCAATTCCCANATTAAAGGCATAGAGAT
ATCCTAGCNANGGACTCCAAGTGTG

10 Sequence ID 895

nt: 560

GTAATTGCAGCCTGGGCAACGGAGTGAGAGACTGTCTCAGGAAAAAAAAAAGAAAA
AAAAGTACTGAGGTAGTTGAATATATCCTCCATTCCTCCATTTGTGGATTAGTTAGT
AAATGGGGCATCTTAGGGTTTAAATATGTCCAGGGTCACTGAGGATCAGATCCTAG
GGTTCCTTTGACTCAAGGCTTTTGTCTCAGCAAAACGTCACCTTCAGCAGGAAGG
15 CTTTCTCAGGCAAGTAGCAGGGTGGCTACTATGTATCGCTTCTTTATTTTTTCTTT
TTTAAATAATGCAGGCACCGTGCGCATAATTTAAAAAATCAGTGCTAAAACCTT
AAAAAAAAAAGCTGTTCTCATCTCCTGTCTTTCTTTTTTTTTCTTTTTATTTTT
TTCTTTTATTATTATTATACTTTAAGTTTTAGGGTACATGTGCACAACGTGCAGGT
TTGTTACATATGTATACATGTGCCATGTNGGTGAGCTGCACCCATTAACTCGTCAT
20 TTAGCATTAGGTATATCTCCTAATGCTATCCCTCCCCCTCCCCCTTTTTTTTTTT

Sequence ID 896

GGGAATGTCTTAGGCACTGGGACTGTAAGTGCAAAGACCCTGTGGCACAAGGGAAT
GTTAATTATCTACCTTTCANAACTGGAANAAGGCCTAGCCTAGAGCATTGAAAAC
25 AATAAGGGAAAGGAGGAGTAAGGCTGGANAGATAGGAATGGTTTAAAGTCTTTGTT
AAAAATTTTTTTTAAAAAATCTTTATCACAAGAAGAGGATTGGCNTGATCAAATTT
GACTTTTAAAAANATTACTTGGGTTGGGCATGATCAAATACTACTTAGGGAGATTA
GTTTANATGATAATGGCATTCTGGACCANAGTGGAGTCAGAGGTGAAAAGAGGTAG
ATATTCCANAATTGAGGGATTTGTGAGGTGAAATCATTGTTACAGATATTAAAGG
30 ATAAGGAGCTTTGTCAAAGGGGATCTTAAGTTTCTGGTATGGTAACTGGGTTAGAG
AGCCCTGGAACATGACCAGCTTTAAGGGAAGAGAGCTTGAGCTCTGTTCTTGTTAA
GCTCAGTTTGAGATCTTTGTGGAATCAAGTGGAGAGGTCTAAGCAGGGAACCTGGCT
TGGCTAGGCTGTAAAGATGAATCTGAGAGTCCCAAGAATATGGTAATTATTAATAA
AAGCCTTAGGTANATGAAATTGTTTTGGG

35

Sequence ID - 897

nt: 509

GCAAATCTACACATTTGATTAAATGATAGGGAAGTATGCACACACATAATACATAT

- 231 -

5 AATGCTAGTTTCTTGGTTTTGATATTGTACCATAGTTATGTAAGATGTAACCATTG
GGGGAAACTGGGTGAAGGCTACATGAGACCTCTCTGTACTTAATCTTTGCAACTTA
TGTGAATCTATAATTATTCCAAAATAAAAAGTTTTAAAGAACCCTAAGTATCCTTAT
TACTGAGGGTCATCGTGCTAGACAGCAAGGTTGGGCCAGAGCTTCTAGTTATTTAA
AATACTAAATACCAGCCTGGGCAACATAGCAAGACCCTGCCTCTACAAAAGCAAA
AAAATTAGCTGGGCATGGTGGTACATGCCTGTGGTCCTAGTTACTCTTGGAGGAGT
CTGAGGTGGGGAGCTTGAGCCTAGGAGTTTGAGGCCGCAGTGAGCCTTGATTGTGT
CTCTGTACTCCAGTCTGGGCCACAGAGCAAGACCCGGTCTCTAAAAATAAATAAAT
AAATA

10

Sequence ID 898

15 ANTGCACTCCAGCTTGGTGACAGAGGGAGACTCCATNTTAAAAAAAAAAAAAAAAAA
AAAAAAGGGAGTAGCTTGAAGCCACATAGTAGTTAGTGGTAAAGGCCACCCCTTT
TCCCACAACTCACACCAGCACCAAGCTAGCCTTTNTAATTTCCAAGCCAGTGCC
CTTTCAACGCACACACCCCTGTGTGAGTTCCCTTTCTGCTGCAAGCTCTCTGGAGG
CAGATACTGTTGAGTCCCTGGCCTGCCTATGAGAACGGCTCATGATCTCTATTTCT
TCTGCTTAATGACCATCTCGAAGTAACAAGTTTAGCCTAAAAATAAAGTTGCTAAGT
TAGCAAAGGAAGTCCTTAGCAGCCACCATTTCTCGATTCCCTCCATCACCTCCCCTG
CCCCTCAACTCCCTCATTTCTCCCAAGATATGGGCTCCAGGCTGGGCGCGGTGGCT
20 CACGCCTATAATCCTAGCACTTTGGGAGGCTGAGGTGAGCAGATCACTGAGGTCAG
GAGTTCG

Sequence ID 899

TCNTTCGGAACGCGCC

25

Sequence ID 900

30 CTGGAGGGATGGGTAGGATTTTGACAAGAGTGGTTGAAGGTATTCTAATTCACCTA
GTACCTACATGTGCGAGGCAGCATGAAGGCAAAAAGCCTGGGGCATGTTTCAGAGA
ATAGCAAGTATTCTAGTTTGAGTGGCACCTGGTACGTATATAAGGGAATAGTAAAA
GATCTGGCTGGAAAGGAAAAGTAGGGGCAGGTTACGAAGGACCTCTGAAAGTCAGA
CTGTGGAACCTGGAACTTTTATCAGGAAGCAGTAGTTAGTTTTTTCAAGCAAAAGCT
AATTAGAGTTGATATTTAGGAGGATGAATCTAACAGTTGTGTGCAAGGATGCCTTC
AAACTGAGTGAGACTAGTACTGGAGACTGGTTAAGAGACTACAACAATAACCTGAG
TAAGAATTAATACAGGCCTGACCTAGTTTTGAGTGAGTAGGATTGGAAACAAGAGT
35 TTTAGGTATTATAGGATTTATGCATATAAAATGGACTTGACAGAACTTGAAGAAAG
AGAAAGTGTCAAAGGACACAGAAAGTGAGGCAGGATATCTTACAATGTTAAAGGA
AAGGAATAATAGAAGTTAC

- 232 -

Sequence ID 903

GGAAACATAAGCTTGTTCAGTACACTCACGCTGTAGATTAATTCTGATATTACAT
ATCTCCATCAGACTTTGTACCCTCTCTCTTCCATCCCTTACCCTTACCGATTAGGT
TGGTATTACCTAAAAATCCATAGAAAATGTCCAGGTGAATTGCCTTATGCTTTCTA
5 CCCATAAGGTATAATT

Sequence ID 904

CTCTGTGGTGTGAGAACACAGTGGGTGACCAAGGCTTTCCAGATGAACCCAAGGAA
AGTGAAAAAGCTGATGCTAATAACCAGACAACAGAACCTCAGCTTAAGAAAGGCAG
10 CCAAGTGGAGGCACTCTTCAGTTATGAGGCTACCCAACCAGAGGACCTGGAGTTTC
AGGAAGGGGATATAATCCTGGTGTTATCAAAGGTGAATGAAGAATGGCTGGAAGGG
GAGTGCAAAGGGAAGGTGGGCATTTTCCCCAAAGTTTTTGTGTAAGACTGCGCAAC
TACAGATTTGGAAAGCACTCGGAGAGAAGTCTAGGATGTTTCAAAACTACAAAGC
TGAAGAAAATGAAGCCCTATTACTTGTTTGTAAGATTTAGCACCCCTTCTGCTGTAT
15 ACTGTACTGAGACATTACAGTTTGGAAAGTGTTAACTATTTATTCCTGTAAAATT
TAACCTACTAGACAATGATGTGAGTACCCAGGATGATTTCCCTGGGGCACAGTGGGT
GAGGAGATGGGGACAGGTGAATGGAGGAGTTAGGGGAGAGGAAAAGTGGATGGAAG
TGTCTGGAAAGGGCACCAAAAAAGTCTTCCAGGTCTGATCCTGTTTCTTGCTCTGA
GTGCTAGCTACCACTGTGTACACTGTAAACATN
20

Sequence ID - 905

nt: 655

CTCAGCTCTTGCCTGGTCACCTTGTGGCTTTTACCATCCTCATCCCCTGTGCCACC
CACATCCTGCCACTTCTGCATGGAGTTGGGGTGGGGCCATTGGAGAAAAGAGGTTA
AACAAGCAGTAATTTACTTGAGTACAGTCTTTGAGCCAATGAAATGCCAGTCATCA
25 TTTCCCAGGGGTACTTGTCTATCTTGTCAACAACCCGCTGATAATGCTCCTTCAATG
TGAATAGCAAAAGTAGGGAGAGACGCTGAATGAAGAAGATGCCTACCCCTCAGGAA
GACTGCTGTCCGCCTCCAGGCCTGCATGCACACACCCATGCCACCTGCACCCCCA
GCACCACGCCCACACTCACTCGCACACACCCACATGCCAGTGTTTTGGGGTTGGCA
GCCTGGACACTGCTGAGGCAAACACAAGTCATCAAGCATAATTCTCATTTCTCTCCT
30 TCTGTCTCTGTTTTAGTTACAGGAATTTGGTCAGTTTAGAGGATTTAATAAGTCCG
TGGAATAATTTGTTTTCTGTCTCTTGCTACCCACGTGAAAAGTAAGTGCATGCTTCAT
GATGTGTTTTTCCCACTACCTTCCAGGCCAGCCGAGCCCACTGGCCANGGCCTGGCC
CGGTGACCTCGGTTGACACTGTCCTCANGCCACTCACTT

35 Sequence ID 906

CAGAATTTTCATGTTTATGCTGCACAAGGCCTGTATTTTATAATGGTGGCTCTTTTG
GACGATGACTTCCTCGATGGTGAAACTTCCAGTAATCTCCCTCATCATACTGAAAT

- 233 -

5 GATATCAGTATATCATCAGAACACCATGGAGCTTGTCAATTTGAGGGACACAGCTTG
CTTGTGTGCTTGGGAAAGAAGAGGTTTAGCATGGTTTCAGGTCAGTGATGAGTCCA
ATGATCTCTGCAAGTTCCCTTAGCTCTGANAATTCTGATGTCATATGCACTTCTGC
CGCCAGAGTTGCTGCTTACTGGATGCGTAAGAAGAAAAGAAAAAAAAAAAAAAAAA

Sequence ID - 907

nt: 582

10 CTTCCATTGGGGGTAAAGATCAAACTTTAGGCGAGCCAGGTCTGTATCTCCATTCC
TGTCTCTGACTGCTTCCCTGTAGGGATTGTCTGCAAGCGCACACCTGCATTTTCTT
GTCCACAAGTCTATGCTCTAACTCTGTCACCTGCATGGCTGCAAATTAGCTTCCTT
CTTCCTGCCCTCTTCTCTCTAGCTTGGATTTTGAATTTGAATGGCAGGCATGGGAT
GTCCGTGTGTGTGTACTGCTGATGTGTACAGCCGCTTGTTAGCGCTCTCATTGTCT
TCAAATGTAAGTCATTTTGGCTGGGTGCGGTGGCTCATGCGTATAATCCCACGCTT
TGGGAGGCTGAGGTGAGCTGATCATTTGAGGTTAGGAGTTCGAGACCAGCCTGGCC
AACATGGCAAAACTCCATCTCTACCAAAAATACAAAATTAGCTGGGTATGGTAGT
15 GCACGCCTGTAATCCCAGCTACTTGAATGCTGAAGCAGGAGAATTGCCTGAACCC
ANGAGGCGGAGGTTGCGGTGAGCCAAGATCACGCCACTGCACTCCAACCTGGGTGA
CAGAGCAAGGCTGTGTCTCAA

Sequence ID 908

20 ACCTGACTTCAAAC TATACTACGAGGCTACAGTAATCAAACAGCATGGTACTAGT
ACAAAAACAGACCAATGGAACAGAATAGAGATCTCAGAAATAAACTGCACATCTA
CAACCATCTGATCTTCAACAAACCTGACAAAACGAGCAATGGGGAAAGGATTCCCT
ATTTAATAAATGGTGCTGGGAGAACTGGCT-AGCCATGTGCAGAAAATTGAACTG
GACCCCTTCCCTTACACCTTATACAAAATTA ACTCAAGATGGATTAAAGACTTAAA
25 TG TAGAACC CAAAACGATAAAAACCCTAGAAGAAAATCTAGGCAATATCATTAAGG
ACATAGACATGGGCAAAAATTT CATGATGAAAACATCAAAGCAATGGCAACAAA
GCAGAACTGACAAATGGGCTTCTGCACAGCAAAAGAACTATCGTCAGAGTGAAC
AGACAACCTACAGAATGGGAGACAGTTTTTTGCAATCTATCCATCTGACAAAAGTCT
AATATCCAGAATCTACAAGGAATTTAA

30

Sequence ID 910

35 CAAAAACAAGAATTACCCGGGCTTGGTGGTGCATGTCTGTAGTCCTATCTACTCA
GGAGGCTGAGGCTGAAGGATCACTTGAGCCCAGGAGTTTGAGGCTGCAGTGAGTGA
GCCATGATCATGCCAGTGTACTCCAGCCTTGGCAGACTGAGCAAACTTGGTCCCT
CGCAAAATGTTGAAGCCCAGTTTTTCACTATTAACCTGTATTTTCACTTTCCCATGC
TAACTTTGAAACACTGGGGCTGGCCTGAGGGGTATAAAGGCTTATTCAAACCTCAGTA
ATTTAAACTTAAAATCCTAAGGAACCTCAAAGGTGTAATCTAGTCCAAATGGGGC

- 234 -

ATCAATTCTAAAGCATTGCTTGTTGAGCAGATTTTCTGTGTCTGAGGTATATAG
ATAACTTATCTTTTTATGACTAAATCCAAGTCCTTAGTTCCTGTTGGAATTCAAAA
TCATATTTAAAAATTGATGCTTTGTTCTATAATTAATGCTTTGATTGTATAAATAA
TAAGTATTCTTCCAAATCCCTTTTTACAGATGATGATTCTGATACCGAGACGTCAA
5 ATGACTTGCCAAAATTTGCAGATGGAATCAAGGCCNGAAACAGAAATCAGAACTAC
CTGGNTCCCAGTCCTGTNCTTAAAATTCTAACTCGAC

Sequence ID - 911

nt: 595

GAGGGTGTAGAAGAGAAGAAGAAGGAGGTTCCCTGCTGTGCCANAAACCCTTAAGAA
10 AAAGCGAAGGAATTTTCGCAGAGCTGAAGATCAAGCGCCTGAGAAAGAAGTTTGCCC
AAAAGATGCTTCGAAAGGCAAGGAGGAAGCTTATCTATGAAAAANCAAAGCACTAT
CACAAAGGAATATAGGCAGATGTACAAANCTGAAATTCGAATGGCGAGGATGGCAAG
AAAAGCTGGCAACTTCTATGTACCTGCAGAACCCAAATTGGCGTTTGTTCATCAGAA
TCAGAGGTATCAATGGAGTGAGCCCAAAGGTTTCGAAAGGTGTTGCAGCTTCTTCGC
15 CTTTCGTCAAATCTTCAATGGAACCTTTGTGAAGCTCAACAAGGCTTCGATTAACAT
GCTGAGGATTGTAGAGCCATATATTGCATGGGGGTACCCCAATCTGAAGTCAGTAA
ATGAACTAATCTACAAGCGTGTTATGGCAAAATCAATAAGAAGCGAATTGCTTTG
ACAGATAACGCTTTGATTGCTCGATCTCTTGGTAAATACNGCATCATCTGCATGGA
GGATTTGATTCATGAGATCTATACTGTTGGAAAAC
20

Sequence ID - 912

nt: 651

CATTTCCAGAGTTTATGTGAATTGAATTGAACTATGGTTTTATGTTACTGTCAGTA
GAATGAAGTACGAATATTTGAAAAATACACCTTCAACTTCAAAGTGATTCTTGACA
AAAATTATAAGGAATCATTTTGGACACATTTTCTGGTAGAGCCTTGTA AAAAATTAA
25 AACCAAGTGTTGTTTTCAAGAAGAACTGTAATACATAATCAGGAATTTGAGTAGGG
AGATTATTTTGTATTTTAAAAATTAAAGTGGCTGTGTAGTTTTAACTTTAGTATTGC
AGGTAGAGTAAGCTTACATGATAACAAAAATCTTGGTCTTAGTGACTTAATGATTC
TGATATTTATTGATTGATTGGTTATCATTCCAAATATTTTAAAAGATAATAGCTGG
CTGGGTGCGGTGGCTCATGCCTGTAATCCAGCACTTTGGGAGGCCAGGACGGGCG
30 GATCACGAGGTCAGGAGATCAAGACCATCCTGGCTAACACGGTGAAACCCCGTCTC
TACTAAAAATCAAAAAATTAGCCGGGTGTAGTGGCGGGCACCTGTAGTCCCAGCTA
CTCAGGAGGCTGAGGCAGGAGAATGGCATGAACCTGGGAGGCGGAGCTTGCAGTGA
GCTGAAATCGTGCCACTGCCTCCACCTGGCGACAA

35 Sequence ID 913

GTGAGGTGGGGACTTCATTCATTGTCCTATTTCTATCTCCACTTTGTGCCTGGAGA
GCTTTCAGGGGAGGTGGAGGAGGAGGGTCTGCCAAGCTACTGCAACATCTGTCACC

- 235 -

CACTATACCCAGTTACTTGGGGGAGGACAGACACTGTGGTGTCAATTAAAGTTGTTT
GAACCAAAGTGGCGGCTGCATCTTTGTCCCGATGCTAGCCGTGCCGGTCTCCCATC
ATCCGCTCGCCCTCCTTTNCCCTGGGCTGCGCCCACTTGTCTTCTGGATATTTGG
GGGTGACTCGCCATGCTTGGCACCCCTCTGCTTCTGGTGCTGCTCTGACTCGAAGA
5 CGGGACAGTCCCTGGTGCACATCCAGGGAAGAGGAGTGTCCGGTAGTTCTTGCAGTA
GGCACTTTATCAGGACCTGACCTGTTGCTGGGTGATTTTAGTCTCTACAAACAGAA
AGCGTTTCAAAGCGTCAGCTGTGGGAGCAGAGTGACCCTTTGCTGATGCTGGGGGG
AGGGGATCTAAATCCTCATTTATCTCT

10 Sequence ID 914

GGCGCCTGCTGGAGGAGGAGAGAGCTCTGCTGGCATGAGCCACAGTTTCTTGACTG
GAGGCCATCAACCCTCTTGTTGAGGCCTTGTTCTGAGCCCTGACATGTGCTTGGG
CACTGGTGGGCCTGGGCTTCTGAGGTGGCCTCCTGCCCTGATCAGGGACCCTCCCC
GCTTTCTGGGCCTCTCAGTTGAACAAAGCAGCAAAACAAAGGCAGTTTATATGA
15 AAGATTANAAGCCTGGAATAATCAGGCTTTTAAATGATGTAATTCCCACTGTAAT
AGCATAGGGATTTTGAAGCAGCTGCTGGTGGCTTGGGACATCAGTGGGGCCAAGG
GTTCTCTGTCCCTGGTTCAACTGTGATTTGGCTTTCCCGTGTCTTCTGGTGATG
CCTTGTTTGGGGTTCTGTGGGTTTGGGTGGGAAGAGGGCCATCTGCCTGAATGTAA
CCTGCTAGCTCTCCGAAGCCCTGCGGGCCTGCTTGTGTGAACCGTGTGGACAGTGG
20 TGGCCGCGCTGTGCCTGCTCGTGTGCTTACATGTCCCTGGCTGTTGAGGCGCTGC
TTTAACCTGCACCCCTNCCTTG-CTCATANATGCTCCTTTTGA

Sequence ID - 915

nt: 230

TTTGAGACCAGCCTAGCCAACATGGTGAAACCCCATCTCTACTAAAAATACAAAAA
25 TTAGCCGGGCGTGGCGGCACATGCCTATAATCCCACCTTACTTGGGAGGCTGANGTA
GGAGAATCGCTTGAACCCANANAGGCAGAGTTTGCAGTGAGCCGAGATTGTGCCAT
TGCACTCCAGCCTGGGCGACAGAGCGAGACTCCATCTAAANAAAATAAATGAATA
AAATAA

30 Sequence ID 917

NNCAGATTTTTTTTTTTTTTTCAGNGTTAGACCATCTTTCAATTCCTGGAACAAAC
TTAACTTTCCATGATATGTATTTTTTATACATTGCTGGATTTTATTTGCTAATATT
TTACTTAGGATTTAATTTTCTAAGTNGACCTATAATTNTCCTGTATAAAATTGCAT
TTGTCACATTTTAGTATCAAGGTTGTCCTANCNCCATGAAATGGATTTANAATGGT
35 TTATGTAANATAAAGTACATTTCTTCTAAAGGTTTGNGTGATTAACCTTCAAATC
TGCCANAGNGNGTTTTTTTTCTTTTTTTTTTTTTTTCATTTNAAGGGAGNGCAAGT
ANCTTTTCAAATNCTGATTTAATTTTAAATAATTTNCAAGTNTNTTANAGTTT

- 236 -

TATTTNTTNTNGAANGTTAACATTTTTATANAAAANGGTNTTATCTTTTTTAAATTC
TTTGACATCAGTTTCTTCANAATTCCTTCTTTTAA

Sequence ID 926

5 GTCATATCTCTTCCCAGGGAAAGCAGGAGCCCTTCTGGAGCCCTTCAGCAGGGTCA
GGGCCCCCTCGTCTTCCCCTCCTTTCCCAGAGCCATCTTCCCAGTCCACCATCCCCA
TCGTGGGCATTGTTGCTGGCCTGGCTGTCTAGCAGTTGTGGTCATCGGAGCTGTG
GTCGCTACTGTGATGTGTAGGAGGAAGAGCTCAGGTAGGGAAGGGGTGAGGGGTGG
GGTCTGGGTTTTCTTGTCCCCTGGGGGTTTCAAGCCCCAGGTAGAAGTGTTCCT
10 GCCTCATTAAGGGAAGCAGCATCCACACAGGGGCTAACGCAGCCTGGGACCCTGT
GTGCCAGCACTTACTCTTTTGTGCAGCACATGTGACAATGAAGGACGGATGTATCA
CCTTGATGGTTGTGGTGTGGGGTCTGATTTCAGCATTTCATGAGTCAGGGGAAGG
TCCCTGCTAAGGACAGACCTTAGGAGGGCAGTTGGTCCAGGACCCACACTTGCTTT
CCTCGTGTTCCTGATCCTGCCTTGGGTCTGTAG

15

Sequence ID 938

TGGCCATCCTTTTCCCCCAAACACACCCCCCTTAACCTATCTCTTGGGACTTAGCC
CGACCCTCCCTCTCATTTCCCATTAAGTCTGAGAGGCAAGAGCTAGGTTAGGCAAG
GAGGTGGTTGGCCAGAGATGGGGAACAGCCAGGTGCCCCAGTCTCTGATTTTTCC
20 TCCATCCTGCTTACCACCTCCCTGGGTACTTACAGCCTTCTCTTGGGAACAGCCGG
GGCCAGGACTGGGTACCTATGAGCTGAATCAGCATCTCCTCCTGAGTCCCAGGGC
CCCTGCAGTTCCCAGTCTCTTCTGTCTGCAGCCCTTGCCTCTTTCCCACAGGTTT
CACTTTATATCCACCTTTTCTTTTGTTCATTTTTTATTTTTTTTATTAT
TAAATGATGTGGTCTATGGAAAAAAAATAAAAATCTGACTTAGTTTT

25

Sequence ID - 939

nt: 513

GGAACCCAGTGTATTACCTGCTGGAACCAAGGAACTAACAATGTAGGTTACTAGT
GAATACCCCAATGGTTTCTCCAATTATGCCCATGCCACCAAAACAATAAAACAAAA
TTCTCTAACAAGAGTGAGCCATGCCTGTAACTGTAAAGAATGTAACA
30 TGTGGGGGACACACAGGGGCAGATGGGATGGTTTAGTTTAGGATTTTATTAGTGCA
TGCCCTACCCTCTGGGGGAACGTCCCATCTGAGGTTTTCTTCTCGGTGGGGGGATT
TAACTTCTGTCTAGGGAAAACAGTGTCTGATGAGGAGTGTTCCAACACAGGCTA
CATGAATTCCCCTATACCAGTGCGAAAGCAGCCAGGAGTCCCCGTTGGAAAAGAAC
AATGCCACTCTCTTTTATGTATCTTGGTCTGCAACTCATTTGTTGTAAGTAGGGT
35 TAATCGAGTATCAGGTTTACAGTATCCTGCCCTTATTATTTTATGATTCACTGACT
CAAGTTCCA

- 237 -

Sequence ID 947

5 GAGAGTGAAAAAATTCTGGTACAAATTGGGAAATTAGTATATAACAACATAGTGTT
AAATTCAATGGGAAAAGTTTAATAAGAGGATTTGGTATCAACTGGCTGTCCAAAGA
TAAAAATGGACCGTCCTATCACATACAAAATTGTTTTTTAGATAAAGATTTAAATA
CAGGCACTCCTTCATTTGCGTGGTGACCTTGAGGTGTTGCAGAAATGATGAGAGC
TGAAACTGCAAAGCAATTTTAATACTTTATCTGTTGGAAATCTTATAGTTTTCTG
TGACCGTTAAAAATTTTCATTAAACTATTAAAAACACCCATGACTGGTCACAAATGT
ATTGGGAAATGGAAAAGAATTAATACTAAAAATACAAAAATAGAAAATATTTA
AAATTATCTAAAAATTTGAAACATTAGAAAAATTGAGAACTAGGCAGGGCGTGGTG
10 GCTCACATCTGTAATTTTAGCCCTTTGGGAGGCTGANGCAGGTGGATCACCTGANG
TCAGGAGTTCGAGACCAGCCTGCCAACGTGGGGAAACCCCGTCTCTACTGAAAATA
CAAAAATTANCCGGGCATGGTGGCACAAGCCTGTAATNCTTGCTNACCAGGANGCT
GAGGCAGGAGAATCACTTGAACCCANGANG

15 Sequence ID 949

GTTCACATGAGAAGGTAGTATTATGTACAGTGACCTTGTTTAAAGTGTCNGTTTA
ATGTTACCACTAAGGCCCTGCCCCAGCTTTATCACCTGAGCACTAACAAAGTGCTGT
GTGGAGTTCAGTCCATGCTGGTAACTNTTGAGTATTCAGTGGGTCTTTTAAACAATT
ACCACCGTGGAGGANANAGCAAGGAAGAGAAATGCTGTGATCTTTTNCTGTTTTTA
20 ATTAGNGAAAGAGGGATTANATTAAACAAATGTTACAGAGNTGTGACTNTGATCCC
CCAGNGGTAAGCAATAATTGTANAGACTGGATTTNANAAGCCCTGAGAGTTTATTT
TCAACCTATNTATTATAGNNCAATCC

Sequence ID 1028

25 ACAAGGCTTGGGGGCTGGACTCCCTCTACTGCCTCTGGCCATACCCCTCCTGGAG
ATGGGGTCAAGGCACCAGGACTGA

Sequence ID - 1056

nt: 435

30 TCGCTTGTAAGCCTGAGACAGCTGCCTGTGTGGGACTGAGATGCAGGATTTCTTC
ACACCTCTCCTTTGTGACTTCAAGAGCCTCTGGCATCTCTTTCTGCAAAGGCATCT
GAATGTGTCTGCGTTCCTGTTAGCATAATGTGAGGAGGTGGAGAGACAGCCCACCC
CCGTGTCCACCGTGACCCCTGTCCCCACACTGACCTGTGTTCCCTCCCCGATCATC
TTTCCTGTTCCAGAGAAGTGGGCTGGATGTCTCCATCTCTGTCTCAACTTCATGGT
GCGCTGAGCTGCAACTTCTTACTTCCCTAATGAAGTTAAGAACCTGAATATAAATT
35 TGTTTTCTCAAATATTTGCTATGAAGGGTTGATGGATTAATTAAATAAGTCAATTC
CTGGAAGTTGAGAGAGCAAATAAAGACCTGAGAACCTTCCAGA

- 238 -

Sequence ID 1071

NGATATAGTNCCGCATGGGAAAGATGANCAGGTATAACCNAGCNTNATATAGCAAG
GACTAACCCCCCTGCCTTCTGCATAATGAATTAAGTAACTAGAAATAACTTNGCAAGGAG
AGCCAAAGCTAAGACCCCNAAACAGACGAGCTACCTAAGAACAGNTAAAAGAGC
5 ACACCCGTCTATGTAGCAAAATAGTGGGAAGATTTATAGGTAGAGGCGACAAACCT
ACCGAGCCTGGTGATAGCTGGTTGTCCAAGATAGAATCTTAGTTCAACTTTAAATT
NGCCACAGAACCCCTCTAAATCCCCTTGTAATTTAACTGTTAGTCCAAAGAGGAA
CAGCTCTTTGGACACTAGGAAAAAACCTTGTAAGAGAGAGTAAAAAATTTAACACCC
ATAGTAGGCCTAAAAGCAGCCACCAATTAAGAAAGCGTTCAAGCTCAACACCCACT
10 ACCTAAAAAATCCCAAACATATAACTGAACCTCTNACACCCAATTGGACCAATCTA
TCACCCTATAGAAGAACTAATGTTAGTATAAGTAACATGAAACATTCTCCTCCGC
ATAAGCCTGCN

Sequence ID - 1074

nt: 689

GGGAGGCGGAGGCTGCAGTGAGCTGAGATCGTGCCACTTCATTCCAGCCTGGGCAA
15 CAAAGCGAACTCTGTCTCAAAAAAAAAAAAAAAAAAATTTGTTGACTGTTGTAA
TTTAAAGCTTGTCAATTTTTTATTTAGTAATAAACTCATTAGTGTAGTATCTATGA
TGAACCAGGTTCTGCACAAAGTACCTTATGTTTCATGGCCTCATATCGTCTTCTCCA
AACTCTGCAAGATAGGATTCATCACCACCTTATAGGGAGAGATCTGAAAGTTTAA
20 ATTGTACCCAAGGTCACACAGCTGGTAAGTGCCAGAGCTGGGATTCCGTAGGGTGT
TCANAGTGCCTCTCCTGCCGTAGGCTTATCACAAAAGTCAAAGTTTGGTCATAAT
AAAGCCTGAAGTTTGGCAGGATTTAAAAATAGTCACCANACTTTTGAGTTGGAGCA
TCCCACCTCACTGCTGTTTACCTTCTGTGGCAGGGAGAGTCATCATTTCCATTTC
GCTTGTGGAATATCTTGTCAATTAACATTCTCATGCAAAAGCCATTTTATGGTGCCC
25 AATGAANATGGTTAAGCTACTGCCCCAAGCCTNTGGAAGCCTTCCTAATTTTGGAC
TTGCACTATGCAAATTGNATAATATTTTCTCTACCCTAAGCCAAATATTTTCTTCA
CTTTTCATTCAATTCTAC

Sequence ID 1081

CGCCGCCGCGCCGCGCTCTCCAACGCCAGCGCCGCTCTCGCTCGCCGAGCT
30 CCAGCCGAAGGAGAAGGGGGTAAGTAAGGAGGTCTCTGTACCATGGCTCGTACAA
AGCAGACTGCCCCGAAATCGACCGGTGGTAAAGCACCCAGGAAGCAACTGGCTACA
AAAGCCGCTCGCAAGAGTGCGCCCTCTACTGGAGGGGTGAAGAAACCTCATCGTTA
CAGGCCTGGTACTGTGGCGCTCCGTGAAATTAGACGTTATCAGAAGTCCACTGAAC
35 TTCTGATTTCGCAAACTTCCCTTCCAGCGTCTGGTGCGAGAAATTGCTCAGGACTTT
AAAACAGATCTGCGCTTCCAGAGCGCANCTATCGGTGCTTTGCAGGAGGCAAGTGA
GGCCTATCTGGTTGGCCTTTTTGAAGACACCAACCTGTGTGCTATCCATGCCAAAC

- 239 -

GTGTAACAATTATGCCAAAAGACATCCAGCTAGCACGCCGCATACGTGGAGAACGT
GCTTAAGAATCCACTATGATGGGAAACATTTTCATTCTC

Sequence ID - 1083

nt: 198

5 GCGCGTCGACTTTGTTTAGACATTGAATGACTTTGTTAAAGGCACAATTAATCACA
TTGGTTGTACTCTGNNGACAGCCTTCTTTAAAAAAAATAAACAATTTAAAACAA
AA
AAAAAAAAAAAAAAAAAAAAAAAAAANTTTTAACC

10 Sequence ID - 1084

nt: 198

GCGCGTCGACTTTGTTTAGACATTGAATGACTTTGTTAAAGGCACAATTAATCACA
TTGGTTGTACTCTGNNGACAGCCTTCTTTAAAAAAAATAAACAATTTAAAACAA
AA
AAAAAAAAAAAAAAAAAAAAAAAAAANTTTTAACC

15

Sequence ID - 1099

nt: 561

TGCATGCTTGTGGATTGGAAAACTTTGGAGACTGATTACTTTTCATTATATATGT
GTCACAGTGAAACAGCTTTTATGTGTCATGTAAGATTACTGCTTGCCTCTCTAAGG
AAGGTCGTGACTGTTTAAATAGACGGGCAAGGTGGAACCTTTTGAAAGATGAGCTT
20 TTGAATATAAGTTGTCTGCTAGATCATGGTTTGTATTGAACTAACAAGGTTTGCAG
ATCTGCTGACTTATATAAAGCTTTTGTATTCTCTACTAAGCTTTAAGATTTAAAAAA
TGTTCAATGTTGAAATTTCTGTGGGGCTCTATTTTGTCTTTGGCTTTCTGGTGAGA
GAGTGAGGAAGCATTCTTTCCTTCACTAAGTTTGTCTTTCTTGTCTTCTGGATAGA
TTGATTTTAAGAGACTAAGGGAATTTACAACTAAAGATTTTAGTCATCTGGTGGA
25 AAAGGAGACTTTAAGATTGTTTAGGGCTGGGCGGGGTGACTCACATCTGTAATCCC
AGCACTTTGGGAGGCCAAGGCAGGCAGAACACTTGAAGGAGTTCAAGACCAGCGTG
G

Sequence ID 1109

30 TTTGNCGGTNTTGGANNNNNANAANTTTCTTCCANNCNTNACNTNTTGGTGGNCTA
AATTAANATGGNTTTNGNGGGTTCNTTNCNTNNNTNNNNCATGGGANANAATTNATT
NTCNTNCNNNTTCTTNNCCCTNAANCTACCTTCCCCNATTTTCTCCCCTNTTCN
TNAATTANCATCCTCTCCNCNTANNTCNANACNTTAATGGCAANACTATCTAATAN
CNANNATAANANCTCCTGTNNNCCACATNTCTTATTNNNCGCNCNCANGTTNCANNC
35 CCNCAGAGTNAACTCATCCTCNCNNAANTTCATATCGTGNNCTNTNNNCNNNTNGC
GCGANATATTAANNANACCNGTANNTNNNANACANNANNTNNGNAANAANCCTTCT
NANNTTTTAGCNTCNGCNNTAACNNNNNTCTTNGTGNNNNCNCAGCTTTCNCNNC

- 240 -

ATNATNCTNCNNCGAANTNTCANNCNTCTCCNCTTNAATGNNTTCCCATGNATTAA
NTNCCTCGNNNANAGCACTATCGTNNNGAGNNNATTATNGNCNNTTTACNTCATG
TGGTCCANTNNCGTTNGNCGCNNNNAATNTTCGTNNNNCNN

5 Sequence ID 1118

GGATTTTAGAGGAAGGCGCTNGGTTACATTGGAGAACTGGAGTGGTCTGGAGTTCC
ACGGTGTTAGTGGACCAGAGGCCACCTCTCCTGGGCTTCTCAGTGTCTCGCCGGCGG
GGTTCGGCCTGAGCTGGATTGACATAGCCCTTGGCGGATTTAAACAACCTAAACAT
TAAGCAGTACAGCTGCCTCAAACCTTTGGGATTTTCAGAATGACTGACACTGCCGA
10 AGCTGTTCCAAAGTTTGAAGAGATGTTTGCTAGTAGATTACAGAAAATGACAAGG
AGTATCAGGAATACCTGAAACGCCCTCCTGAGTCTCCTCCAATTGTTGAGGAATGG
AATAGCANAGCTGGTGGGAACCAAAGAAACAGAGGCAATCGGTTGCAAGACAACAG
ACAGTTCAGAGGCAGGGACAACAGATGGGGGTGGCCAAGTGACAATCGATCCAATC
AGTGGCATGGACGATCCTGGGGTAACAACCTACCCGCAACACAGACAAGAACCTTAC
15 TATCCCCAGCAATATGGACATTATGGTTACAACCAGCGGCCTCCTTACGGTTACTA
CTGATAGAAATGTTGGCAGCTTTTAGTAAAAGCATTTACTCTGTTACCATGAGAAA

Sequence ID 1125

NGACTGGCTCCCGAAAAGAAGGGTGGCGAGAAANAAAAGGGCCGTTCTGCCATGGA
20 CGAAGTGGTAACCCGCGAATACACCATCAACATTNACAAGCGCATCCATGGAGTGG
GCTTCAAGAANCGTGCACCTCGGGCACTCAAAGAGATTCGGAAATTTGCCATGAAG
GAGATGGGAACTCCATATGTGCGCATTGACACCAGGCTCAACAAANCTGTCTGGGC
CAAAGGAATAAGGAATGTGCCATACCGAATCCGTGTGCGGCTGTCCANAAAACGTA
ATGAGGATGAAGATTCACCAAATAAGCTNTATACTTTGGTTACCTATGTACCTGTT
25 ACCACTTTCAAAAATCTACAGACAGTCAATGTGGATGANAAACNAATCGCTGATCGT
CAGATCAANAAANT

Sequence ID - 1139

nt: 503

CAGCACTGCCAGTGGAGATGGGCGTCACTACTGCTACCCTCATTTACCTGCGCTG
30 TGGACACTGAGAACATCCGCCGTGTGTTCAACGACTGCCGTGACATCATTCAGCGC
ATGCACCTTCGTCAGTACGAGCTGCTCTAAGAAGGGAACCCCCAAATTTAATTAAA
GCCTTAAGCACAATTAATTAAAAGTGAAACGTAATTGTACAAGCAGTTAATCACCC
ACCATAGGGCATGATTAACAAAGCAACCTTTCCCTTCCCCCGAGTGATTTTGCGAA
ACCCCTTTTCCCTTCAGCTTGCTTAGATGTTCCAAATTTAGAAAGCTTAAGGCGG
35 CCTACAGAAAAAGGAAAAAAGGCCACAAAAGTTCCCTCTCACTTTCAGTAAAAATA
AATAAAACAGCAGCAGCAAAACAAATAAAATGAAATAAAAGAAACAAATGAAATAAA
TATTGTGTTGTGCAGCATTAAAAAAATCAAAATAAAATTAATGTGAGCAAAG

- 241 -

Sequence ID - 1148

nt: 587

5 TGAAAAATAAAGTTTTTATGTATATTCTACATATGTATATGTTGGTAGAAAGCAAA
AACGCTAGGTAAAAATAAATGTAATACAATTTTAGCTATGAACCAAAAAACCATTT
GTGGTGTGGATGCAAGAAAGTCTGGATGGGTGCAGAGTTCTCCATGTTTCACCTTCT
GACATTTGAAAATACGCAGTTTGCATTTGATACGTCAAATGTTATTTTTAAGAAAA
CCAATAAAATCATTAACCGAAAAGGCAGTTTTGCTTGTTTTTACCTTAGTTGGA
GTTATCTGCAATTGCCGTATTAGTGTTTTAAGGAACTTGTAAGTAAGCTCCTTAGT
CCCCTTTAGAGCTACGAAACATGTCAATTTTACTTTTCTCCAGCTTTTTGGAATCT
TATCTAAATTACCATGTAGAGTTCTGCATAGCTTCAAATTTCTCTTAGCCAATGTGG
10 TCTGTAAGTGTCTATCGATGAATTTACCCGTTAATTGCCGTAGTATACTGTCCTGT
ACCGGATGTGAAGAGGAGCAACTCTGCACAGTGCACCTGGTTGCTCCCATGGTAGGA
ANGAATGGCTTATCAATGGTCGGATTT

Sequence ID - 1160

nt: 650

15 GGAGGATGGAGCAGTGAGCGGGTCTGGGCGGCTGCTGGCAGCGCCATGGAGACGGT
ACAGCTGAGGAACCCGCCGCGCCGGCAGCTGAAAAAGTTGGATGAAGATAGTTTAA
CCAAACAACCAGAAGAAGTATTTGATGTCTTAGAGAACTTGGAAGGGTGAGTG
TAAAGAACTATAGGTAGGTCATTGGGTCCCAGTCTTTTTCCTGCCCCAGAAGAAG
CAGAAGGATATGAACCTTTAGCATTTGTTCTAGGTGGGGTGGAAGGTAAATTTACA
20 GCTTGATGTCTTCTTCGCTTTACTCCAATCCCTATTATAGACAGATTTAGTGA
TTCCTGGTCTTTTAAACACGAAGAATATCTATTGTTTTCTCTTTTGTAGGATCTGT
ATGATTTTATCTACTTAACAGATAGCACTAATTAGATTAAAATTCTATAAGAACT
TTTAAATTTGCTGTTTATAATTTCTGATTGGTATGCAATAACTGTTTCAATGAAAA
TCAATGTAATTTAGTATTTTAAATTTGCACCTTTGTGAAATATAGTAAATAAATT
25 AAGCACTATCACCACCTTCACAGCTACTTAGGAGATCCACAATCCTGGGTTGGGAG
CCAGTGGATTTCTGAAACACAGATTTGTTAATG

Sequence ID - 1165

nt: 502

30 CTCAAGTGAATCCTGGCTTCTTGGAAGCGCTTGCCCTAGACGAGACACAGTGCATAA
AAACAACCTTTTGGGGGACAGGTATGTTTTCTTGCAGCTGCGGTGTGAAGGTCTTGG
CAAGACAAGCAGTGTGGCCAGAATTTTGAACCTTCTGATGAATGTGTAATGCAAAGG
ACCTTGATACATTTTTTTGTTTCAAGGTCCTCAAATGAGCACATGAAGAGGTTGCT
GTGAAACTTTAAGTGGCCCTACTGCGCAGAAGCATTCAGATGTCACTTGATGATCT
GTAAGGGAACCTTGCTGATTTGGGAATGTGCTTAGGGGAACACACATTCCTTTTGACA
35 GGGTCTGTCACTGGGTGGGTGATGAATTATACAGATGACATGTGCTTTTTTTTCTT
TTTTCAACCTCAATGGTATTCCTACAGGAAATGGATAACCATTTTAACTGTATTTT
TTGCAGCCCGTACCTTCTTGGGAATACAATTGTCTAACTTTTTATTTTTGGTCT

- 242 -

Sequence ID - 1172

nt: 648

CCACAATAATAAGAGAAAAACAGGAGCAAAAGGATATACAAAACCACCAGAAAACA
AATAACAAAGTGACAGGAGTAAGTCCTTAACTGGCAATAATAACCATGAATCTAAA
TGGATTCCATTTCCCACTTAAAAGATAAAGACATGCTGAATGGATAAAAAGCTGTC
5 ACCCAGTTATATGCTGCCTACAACAACTCACTTCACCTGTAAACATACATATGGA
TGGAAAGAGAAGGCATGGGAAAAGATACTCTACTCAAATGAAAACAAAACCAAAC
AAAGGTGGCTATTCTTATATGAGATAATACAGACATTAAATCAAAAACCTGGAAACA
AACACAAAGTCATTGTATAATGATGAATTCAATTATATCATGATGAATTCAATTAT
ATCCTCCTTCCTGATCAATTCAGAAAGGAGGATATAATCTTTTAAATATATATAC
10 ACCCAACACCAGAGCATATAAATATGTAAAGGAAGATAAAGGGAGTCCTGTGATCA
AGAATAAATATAACAATTATAAATATTTTATCTAAAGTGATAGATAGACTGTAATA
CAATAATAGGGTGGTGACATTAACACCCCCTCTCACATTGGACTGATCATCTAGAA
GGGAGAAAAAGCTTTATGATTGGAAAAGCCAT

15 Sequence ID 1178

ATTGTGTTGGCCACCCGGGAATTGCGGGCCGCGTCGACCTACGCACACGAGAACAT
GCCTCTCGCAAAGGATCTCCTTCATCCCTCTCCAGAAGAGGAGAAGAGGAAACACA
AGAAGAAACGCCTGGTGCAGAGCCCCAATTCTACTTCATGGATGTGAAATGCCCA
GGTGAGGAGACGGCTTGCTGTAGTGGGGAAAGCACTGGACCTCAACAGTTGGAAAA
20 TGTGTAGTGTTAGCTGTCTCGTATCCTTGAAGCTGTGCAGCAGCTTCAGTTTCTT
CGCCTGTGGAAAATATTTTCCCTGATACTCTTAAAATTTGAATGTATGAGACTGGC
AAAGTTTTGCATCTTAGGAGGAGTGATTCAATTTACCGTGATCTCTCATCACATTT
CACATACAACCCCTACGTTTTTTTGTGTTGGGAAACAATGTAATGGATGATGAGTT
GGGCATAAGTGCAGGAAAGACGGGTGTAATAGAGGAAAAAATGTTATCTGCTTTT
25 CTTTCAGGATGCTATAAAATCACACGGTCTTTAGCCATGCACAAACGGTAGTTTT
GTGTGTTGGCTGCTCCACTGTCTCTGCCAGCCTACAGGAGGAAAAGCAAGGCTTA
CAGAAGGATGTTCCCTTCAGGAGGAAGCAGCACTAAAAGCACTCTGAGTCAANATGA
GTGGGAAACCATCTCAATAAACACATTTTGGAT

30 Sequence ID - 1180

nt: 622

CTTTTCCTCCCGCTGTCCCCACGGGAGGGGACTGCTCTCCCCGCTGCATCCTTT
CTGTGAGGTACCTTACCCACCTCAGCACCTGAGAGGGTGAAATAGAATTCTAACCT
CGACATTCGGGAAGTGTTTTTGAGAAGTCTCGGTCCGTAAGGGAAGTCTTCCAAGT
CCGTGCAGCACTAACGTATTGGCACCTGCCTCCTCTTCGGCCACCCCCAGATGAG
35 GCAGCTGTGACTGTGTCAAGGGAAGCCACGACTCTGACCATAGTCTTCTCTCAGCT
TCCACTGCCGTCTCCACAGGAAACCCAGAAGTTCTGTGAACAAGTCCATGCTGCCA
TCAAGGCATTTATTGCAGTGTAATTTGCTTCAAAGGATCAGGCCCTGAGAACA

- 243 -

ATGACCTTATTTTCCTACAACAGTGTCTGGGTTGCGTGCCAGCAGATGCCTCAGATA
CCAAGAGATAACAAAGCTGCAGCTCTTTTGATGCTGACCAAGAATGTGGATTTTGT
GAAGGATGCACATGAAGAAATGGAGCAGGCTGTGGAAGAATGTGACCCTTACTCTG
GCCTCTTGAATGATACTGAGGAGAACAACCTCTGACANCCACAATCATGAGGATGAT
5 GTGTTG

Sequence ID - 1181 nt: 155
CGCCACTTATCCAGTGAACCACTATCACGAAAAAACTCTACCTCTCTATACTAAT
CTCCCTACAAATCTCCTTAATTATAACATTACAGCCACAGAATAATCATATTAA
10 AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Sequence ID 1182
CATTGTGTTGGCNCCTGGGAATTCGCGGCCGCGTCGACTTTTTGTGTTGTTTGGAG
CAGAAATACTAAAGAAGATTCGGGGCCGAGTATCCACAGAAGTGGACGCAAGGCTC
15 TCCTTTGATAAAGATGCGATGGTGGCCAGAGCCAGGCGGCTCATCGAGCTCTACAA
GGAAGCTGGGATCAGCAAGGACCGAATTCCTTATAAAGCTGTCAACCTGGGAAG
GAATTCAGGCTGGAAAGGAGCTCGAGGAGCAGCACGGCATCCACTGCAACATGACG
TTACTCTTCTCCTTCGCCCAGGCTGTGGCCTGTGCCGAGGCGGGTGTGACCCTCAT
CTCCCCATTTGTTGGGCGCATCCTTGATTGGCATGTGGCAAACACCGACAAGAAAT
20 CCTATGAGCCCCTGGAAGACCCTGGGGTAAAGAGTGTCACTAAAATCTACAACCTAC
TACAAGAAGTTTAGCTACAAAACCATTTGTTCATGGGCGCCTCCTTCCGCAACACGGG
CGAGATCAAAGCACTGGCCGGCTGTGACTTCCTCACCATCTCACCCAAGCTCCTGG
GAGAGCTGCTGCAGGACAACGCCAAGCTGGTGCCTGTGCTCTCAGCCAAGGCGGCC
CAAGCCAGTGACCTGGAAAAAATCCACCTGGATGAGAAGTCTTTCGTTGGTTGCA
25 CAACGAGGACCAGATGGCTGTGGAGAAG

Sequence ID - 1183 nt: 479
CGTGGCAGCCATCTCCTTCTCGGCATCATGGCCGCCCTCAGACCCCTTGTGAAGCC
CAAGATCGTCAAAAAGAGAACCAAGAAGTTCATCCGGCACCAGTCAGACCGATATG
30 TCAAAATTAAGCGTAACTGGCGGAAACCCAGAGGCATTGACAACAGGGTTTCGTAGA
AGATTCAAGGGCCAGATCTTGATGCCCAACATTGGTTATGGAAGCAACAAAAAAC
AAAGCACATGCTGCCCAGTGGCTTCCGGAAGTTCCTGGTCCACAACGTCAAGGAGC
TGGAAGTGCTGCTGATGTGCAACAAATCTTACTGTGCCGAGATCGCTCACAATGTT
TCCTCCAAGAACCGCAAAGCCATCGTGGAAGAGCTGCCCAACTGGCCATCAGAGT
35 CACCAACCCCAATGCCAGGCTGCGCAGTGAAGAAATGAGTAGGCAGCTCATGTGC
ACGTTTTCTGTTTAAATAAATGTAAAACTG

- 244 -

Sequence ID - 1185

nt: 628

CTTTGATTACCTTTGAGTATTAGGTTGAAAGCTTCTCTGTGCTTGATTGAACATTG
TGATGATGTTGATTGGGTCATGTCAGATTTAGACAGTGTTGTGTTTAAGATAAATG
TTTAATGGCTCTTAGCAGTGTTTCATGCCTCCCCTTTTCCCCTGATACTTTAAAAAC
5 AGAATATACAGAAAAGGGGAGTTGGGTGAAGAATCACCATATTCTCATTACCAGAG
TAGTGTCTACCAGCTGTTTTACATTTTTCTGTTTCCTTCTGTCCTTGGAATCCTT
TTTTTAGATCCTTGTAATACTAGTAAAGATATTCCACTCTGTGTTGTAAGCATTTT
TCCATTTTGCTCCATGGTCTTCATAATGCCCTGTGGTCCTTTATTAAGGGGATGCA
CCATGTAGAGGTGAAAGGCTTTCCTTGACTTGCCACCATTCTGTATTTTCCTTA
10 GAGGAGGAGGTTTCCAACATTTCTTTTTTAGAGACAGAGTCTCGTTCTGACACGCA
GGCAGGAGTGCAGTGGCATGATAACAGCTCACTGCAGCCTCGAACTCCTGGGCTCA
AGTTATCCTCCACCTCAGCTTCCTGAGTAGCTAGGACTGCAGGTGCCTGCCACCA
CACCCAGCTAAT

15 Sequence ID - 1186

nt: 494

CAGCCCTCCGTCACCTCTTCACCGCACCCCTCGGACTGCCCCAAGGCCCCCGCCGCC
GCCTCCAGCGCCGCGCAGCCACCGCCGCGCCGCGCCTCTCCTTAGTCGCCGCCA
TGACGACCGCGTCCACCTCGCAGGTGCGCCAGAACTACCACCAGGACTCAGAGGCC
GCCATCAACCGCCAGATCAACCTGGAGCTCTACGCCTCCTACGTTTACCTGTCCAT
20 GTCTTACTACTTTGACCGCGATGATGTGGCTTTGAAGAACTTTGCCAAATACTTTC
TTCACCAATCTCATGAGGAGAGGGGAACATGCTGAGAACTGATGAAGCTGCAGAA
CCAACGAGGGTGGCCGAATCTTCCTTCAGGATATCAAGAAACCAGACTGTGATGAC
TGGGAGAGCGGGCTGAATGCAATGGAGTGTGCATTACATTTGGAAAAAATGTGAA
TCAGTCACTACTGGAACGACAAACTGGCCACTGACAAAAATGAC

25

Sequence ID - 1188

nt: 599

GGGAGACAAGCCCAGCCTTTCGGCGAGNATACGTCTAACCTGTGCAACAGCCACT
ACATTACTTCAAACAGATCCTTCCTTTTGAGGGAGCAAGTCCTTCCCTTTCATT
TTTTCCAGTCTTCCCTCCCTGTGTATTCACTCTCATGATTATTATTTTAGTGGGGGC
30 GGGGTGGGAAAGATTACTTTTTCTTTATGTGTTTGACGGGAAACAAAACCTAGGTAA
AATCTACAGTACACCACAAGGGTCACAATACTGTTGTGCGCACATCGCGGTAGGGC
GTGGAAAGGGGCAGGCCANAGCTACCCGCAGAGTTCTCAGAATCATGCTGAGAGAG
CTGGAGGCACCCATGCCATCTCAACCTCTTCCCCGCCCCGTTTTACAAAGGGGGAGG
CTAAAGCCCAGAGACAGCTTGATCAAAGGCACACAGCAAGTCAGGGTTGGAGCAGT
35 AGCTGGAGGGACCTTGTCTCCAGCTCAGGGCTCTTTCCTCCACACCATTTCAGGTC
TTTCTTTCGAGGGCCCCTGTCTCAGGGTGAGGTGCTTGAGTCTCCAACGGCAAGGG
AACAAGTACTTCTTGATACCTGGGATACTGTGCCCAGAG

- 245 -

Sequence ID 1189

GGGAGACAAGCCCAGCCTTTCGGCGAGATACGTCTAACCCCTGTGCAACAGCCACTA
CATTACTTCAAACCTGAGATCCTTCCTTTTGAGGGAGCAAGTCCTTCCCTTTTCATTT
TTTCCAGTCTTCCTCCCTGTGTATTTCATTCTCATGATTATTATTTTAGTGGGGGCG
5 GGGTGGGAAAGATTACTTTTTCTTTATGTGTTTGACGGGAAACAAACTAGGTAAA
ATCTACAGTACACCACAAGGGTCACAATACTGTTGTGCGCACATCGCGGTAGGGCG
TGGAAAGGGGCAGGCCAGAGCTACCCGCAGAGTTCTCAGAATCATGCTGAGAGAGC
TGGAGGCACCCATGCCATCTCAACCTCTTCCCCGCCCGTTTACAAAGGGGGAGGC
TAAAGCCCAGAGACAGCTTGATCAAAGGCACACAGCAAGTCAGGGTTGGAGCAGTA
10 GCTGGAGGGACCTTGTCTCCAGCTCAGGGCTCTTTCCTCCACACCATTACAGGTCT
TTCTTTCCGAGGCCCTGTCTCAGGGTGAGGTGCTTGAGTCTCCAACGGCAAGGGA
ACAAGTACTTCTTGATACCTGGGATACTGTGCCAGAGCCTCGAGGAGGT

Sequence ID 1190

15 GTTTAAATTTGACAACTAAAGCTNATNACTGCTATAAGAGTAATAACTGCTCATT
TTCCATAACTCATTCTTAAAGTTTTAGTAATGTAAAAGTTATTTTTTTTGCAGTAAG
TTATAATGATAGAAGCTTACATGTTTTTTTCATGCCTCATCTGTTTCCCCTTAAAC
TATAATTATCAGTAAAGTCCTGTGGTATTTTTCAATTTGTAAGAACTAGGCTATA
TATACATTGGGAAAAACAGCCTTCATTTGTCAATGCACTAGTGTTCCAAAGGTTTC
20 TGGTAATTGTGTGCTATTGCTTTTTGTTGACTTGCAAAAAAAAAAAAAAAAAAATT
ACTATGACTTGNGGTAGCCCTGCAACCTTCGGAAGTGCTTAGCCCAGTCTGACCAT
ACATTTATATTTANAATGCTTAGGTAAATAAATAATATGCCTAAACCCAATGCTAT
AAGATACTATATAATATCTCATAATTTTAAAAATCACTGTTTTGTATAATAATAAA
ACAAGGCAGGCAAGCTGTTCTACAATGACTGTTGGTAAGGGTGCTGAGGAAGAAAA
25 ACAACAATCTTGATTGAGGGATAGTGAATAGACAAAAAATGTCCTAATCAATGAA
GCTGTGTGATGATTCTGATTGACAGAGA

Sequence ID 1191

GTGCAAAGTGTTATATCCACTTTCAACAAAGAGAGAAGCTGAAAAGCTAACCCAAT
30 GTTAATTTTGGATCACACACATTGAGTGAGACTTTAAGATTTTACTTCTGTTGGA
GTAGCTATATTATTTCTAGTTAAAAAATCTCTATATACATATTTATTTGTTTTTC
TACTTGTTTAATATTTTTCTCTTCCAATTAGGAACTCAATATGGAATAAAAAATAT
TTAAATGTATTTTACTCAAACGTGTGTGTATATATGTTTGTGTGCATGATAAGGAG
AGTGAGAGCAAGAGTAAGAGAGAGAGAGCACGCATAGATGGAAGCACACATTTAAT
35 GTCTATGAAATGAGAAAACATTAAGGCTAAGATATTTTTCTTCTGAACTAGCAGA
TTGTATCAATGGCTGGTCACTTAAATTAATCAGTTTGTAAAGATATTTAAAAGGTA
TGTCTACCTTCTTGCAATTAATTTGATTATGTTCTAATGGCATGGCAAGAGAAATG

- 246 -

AAAGAAGATAACTAAAAGTTAAAAGTCGTTGCATGTTTTTGTTCAGCATACCCCTT
CTTTCAGGCTACCGAATAACCTTGATTGACATTGGATTAGTAGTAGAATACCTCAT
TGGTAGAGCATATCGCAGCANCTACACTAGAAAACAT

5 Sequence ID 1192

GTCTGGAACCTCCAGACCTCAGGTGATACCCCTGCCTCAGCCTCCCAATGTGCTGGG
ATTACAGCTGTGAAGCCACCGCGCCCGGCTGCTGTGATAGTTGAGATGTAAACCAA
AAATAAAATTTCTAAGCCACCCAATCCGACTGAATGGACCCTTCCTGTTGAGCAAGG
ACATTCCAAAGTAAACTGAAAAGACCAGCTTAGGCCATGATGGGAAGGGGAGGTGT
10 CAACATGCCTCATTCTACCTTCCTCCCTCTGGAATCCAGACACAACCTGACCAGCAT
TAACATTAAAACAGAGATCTTAAGCTGGGCACGGTGGCTCATGCCTGTAAATCCCAG
CACTTTGGGAGGCCAAGGTGGGATCACCTGAGGTGCGAAGTTCAAGACCAGCCTGG
CCGGTATGGTGAAGCCATGTCTCTACTGAAAATGCAAATTTGGCCGGACATTGTGG
TGCA

15

Sequence ID 1193

TNCNTTTTTTTTCCCNCGGGAAAGCGCGCCATTGTGTTGGTCCCCGGGAATTCGCG
GCCGCGTCGACGAGAAATGGCTTGAACCCAGTAGGCAGAGGTTGTAGTGAGCCAG
AATNGGNCACCTGCACNTTANCCNTGGGTGACAAAANTGAAAACCTTGTCTNAAA
20 AAAAAAAAAAAAAAATTTTAANTNAAATNAAAAANCCTTTNCNTTNTTTTTNAAAN
NGGGGGGGGNNTTTTTNGGGNTTNGNNNTGGTAAAAANTNNNTTTTTTTTTTTTA
GGGGCCNANNCCCCNTTTTANAAAANCCNGNTTTTNAAAAAANTTTTTTNCCNCN
NTTNGGGGGGGGGGNTTTTNANCNNTNTTNGGGGGGGNNCCCCNTTANNACCNNC
AAANTTTTTANTTTTTTGNNAANNCCCCCTTTTTTNNTTTTTTTTTNGGGGGGGG
25 GGGNNGCCCCCNCCCTTTNGGGGGGGGGGNTTNGNAAAANNACTTTTNAAAANNA
AGGGNNGGGGNANATNNCCCCCCCCNGGNTTTTTTTTTTAAAAANTNAANNGGGGG
GGGNNNCTNANTNGGGGCNCCANNGGGGGNTTANAANNATTTTCTNCCCAAACCC
CCNGNTTTTATNNCCCCCCCCCCCCCNCNNNGAANGGGNGGNCCNTTTTTTTTTAT
TTTNNGGNGGGNAAAAAANTTTNAAAAANNANATNTTTTTTCCCCCCCCCCCCCNC
30 TTTTNGGNAAANCCNNGGGGGGNTCCTTTTTNAAANNNNCCCCCAAAAAAANTTT
TTTTNTTNTNTTTTTCTCTNGGGGNCCNNANTTNTANANTTTTNCNCCNAAAAAA
ANGGGNCCCCCTTTTTTNCNGGNNGGNNCCCAAAANNTTTTTTTTTNAAAAA
AAAA

35 Sequence ID 1195

GTTCGTGACNTTCGGAGCTACCTGACAGAGCAGAGTCAACCAGGNTCTGCCCAAAG
AGAGTGTTAGGCCTGAGCTTGAGAGCCCTGGAGAGACGTGTGCACAAAATGTGACC

- 247 -

TGAGGCCCTAGTCTAGCAAGAGGACATAGCACCCCTCATCTGGGAATAGGGAAGGCA
CCTTGCAGAAAATATGAGCAATTTGATATTAACATCTTCAATGTGCCATAGA
CCTTCCCACAAAGACTGTCCAATAATAAGAGATGCTTATCTATTTTA

5 Sequence ID - 1196

nt: 412

GTCTGACGCGGCCGCGGTCTGCTGGAGNCGATCAACTCTAGGCTCCAACCTCGTTATGA
AAAGTGGGAAGTACGTCTTGGGTACAAGCAGACTCTGAAGATGATCAGACAAGGC
AAAGCGAAATTGGTCATTCTCGCTAACAACCTGCCCAGCTTTGAGGAAATCTGAAAT
AGAGTACTATGCTATGTTGGCTAAAACCTGGTGTCCATCACTACAGTGGCAATAATA
10 TTGAACTGGGCAÇAGCATGCGGAAAATACTACAGAGTGTGCACACTGGCTATCATT
GATCCAGGTGACTCTGACATCATTAGAAGCATGCCAGAACAGACTGGTGAAAAGTA
AACCTTTTTCACCTACAAAATTTACCTGCAAACCTTAAACCTGCAAAATTTTCCTT
TAATAAAATTTGCTTGTTTT

15 Sequence ID 1197

CCGCCAACATGGGCGCGTTCGCACCAAAACCGTGAAGAAGGCGGCCCGGGTCATC
ATAGAAAAGTACTACACGCGCCTGGGCAACGACTTCCACACGAACAAGCGCGTGTG
CGAGGAGATCGCCATTATCCCCAGCAAAAAGCTCCGCAACAAGATAGCAGGTTATG
TCACGCATCTGATGAAGCGAATTCTAGAGAGGCCCAGTAAGAGGTATCTCCATCAAG
20 CTGCAGGAGGAGGAGAGAGAAAGGAGAGACAATTATGTTCTGAGGTCTCAGCCTT
GGATCAGGAGATTATTGAAGTAGATCCTGACACTAAGGAAATGCTGAAGCTTTTGG
ACTTCGGCAGTCTGTCCAACCTTCAGGTCACTCAGCCTACAGTTGGGATGAATTTTC
AAAACGCCTCGGGGACCTGTTTGAATTTTTTCTGTAGTGCTGTATTATTTTCAATA
AATCTGGGACAA

25

Sequence ID 1198

CAGAGGTGGGAGGATTGCTTCAGTTCAAGAGTTTGAGACCAGCCTGGGTAACATGG
CGAAACCCTGTCTTTACAAAAAATGCAAACCTTTGCCGCATGTGTTGGGGTGCGCC
TGTAAGTCCCAGCTTCTCGGGAGGCTGAGGTGGGGGGACCACCTGAGCCATGGAGGT
30 TGAGGCTGCAGTGAGCCGTGATACCACCACTGTACTCTAGCCTGGGCCATAGAGTG
AGACACCCTGCCTCAGAAATA

Sequence ID - 1199

nt: 439

CCCATCCCCTCGACCGCTCGCGTCGCATTTGGCCGCCTCCCTACCGCTCCAAGCCC
35 AGCCCTCAGCCATGGCATGCCCCCTGGATCAGGCCATTGGCCTCCTCGTGGCCATC
TTCCACAAGTACTCCGGCAGGGAGGGTGACAAGCACACCCTGAGCAAGAAGGAGCT
GAAGGAGCTGATCCAGAAGGAGCTCACCATTGGCTCGAAGCTGCAGGATGCTGAAA

- 248 -

TTGCAAGGCTGATGGAAGACTTGGACCGGAACAAGGACCAGGAGGTGAACTTCCAG
GAGTATGTCACCTTCCTGGGGGCCTTGGCTTTGATCTACAATGAAGCCCTCAAGGG
CTGAAAATAAATAGGGAAGATGGAGACACCCCTCTGGGGGTCTCTCTGAGTCAAAT
CCAGTGGTGGGTAATTGTACAATAAATTTTTTTTGGTCAAATTTAA

5

Sequence ID - 1200

nt: 526

CTGGAGACGACGTGCAGAAATGGCACCTCGAAAGGGGAAGGAAAAGAAGGAAGAAC
AGGTCATCAGCCTCGGACCTCAGGTGGCTGAAGGAGAGAATGTATTTGGTGTCTGC
CATATCTTTGCATCCTTCAATGACACTTTTGTCCATGTCACTGATCTTTCTGGCAA
10 GGAAACCATCTGCCGTGTGACTGGTGGGATGAAGGTAAAGGCAGACCGAGATGAAT
CCTCACCATATGCTGCTATGTTGGCTGCCCAGGATGTGGCCAGAGGTGCAAGGAG
CTGGGTATCACCGCCCTACACATCAAACCTCCGGGCCACAGGAGGAAATAGGACCAA
GACCCCTGGACCTGGGGCCCAGTCGGCCCTCANAGCCCTTGCCCGCTCGGGTATGA
AGATCGGGCGGATTGAGGATGTACCCCCATCCCCTCTGACAGCACTCGCAGGAAG
15 GGGGGTCCGCGTGGTCCGCGTCTGTGAACAAGATTCTCTAAAATATTTTCTGTTAA
TAAATTGCCTTCATGTAAACTG

Sequence ID - 1201

nt: 613

CTTAAGTATGCCCTGACAGGAGNATGAAGTAAAGAAGATTTGCATGCAGCGGTTCA
20 TTAAATCGATGGCAAGGTCCGAACTGATATAACCTACCCTGCTGGATTCATGGAT
GTCATCAGCATTGACAAGACGGGAGAGAATTTCCGTCTGATCTATGACACCAAGGG
TCGCTTTGCTGTACATCGTATTACACCTGAGGAGGCCAAGTACAAGTTGTGCAAAG
TGAGAAAGATCTTTGTGGGCACAAAAGGAATCCCTCATCTGGTGACTCATGATGCC
CGCACCATCCGCTACCCCGATCCCCTCATCAAGGTGAATGATACCATTGAGATTGA
25 TTTAGAGACTGGCAAGATTACTGATTTTCATCAAGTTTCGACACTGGTAACCTGTGTA
TGGTGA CTGGAGGTGCTAACCTAGGAAGAATTGGTGTGATCACCAACAGAGAGAGG
CACCTGGATCTTTTGGACGTGGTTCACGTGAAAGATGCCAATGGCAACAGCTTTGC
CACTCGACTTTCCAACATTTTGTATTATTGGCAAGGGCAACAAACCATGGATTTCTC
TTCCCCGAGGAAAGGGTATCCGCCTCACCATTGCTGAAGAGAGAGACAAAAGA

30

Sequence ID 1202

GGAATTTCGCGGCCGCGTCGACCTCTGCTCGAATTGACAGAAAAGGATTCTGTGAAG
AGTGATGAGATTTCCATCCATGCTGACTTTGAGAATACATGTTCCCGAATTGTGGT
CCCCAAAGCTGCCATTGTGGCCCCGCCACACTTACCTTGCCAATGGCCAGACCAAGG
35 TGCTGACTCAGAAGTTGTCATCAGTCAGAGGCAATCATATTATCTCAGGGACATGC
GCATCATGGCGTGGCAAGAGCCTTCGGGTTGAGAAGATCAGGCCTTCTATCCTGGG
CTGCAACATCCTTCGAGTTGAATATTCCTTACTGATCTATGTTAGCGTTCCTGGAT

nt: 692

Sequence ID 1204

Sequence ID 1205

[illegible]

- 250 -

TNNNNNTNNCNNCNANTAAANNCANNTCNANNNNNANNNNAATTACTTNNANGTNNNTC
ACN

Sequence ID - 1207

nt: 642

5 ACGAGAAGCCAGATACTAAAGAGAAGAANCCCGAAGCCAAGAAGGTTGATGCTGGT
GGCAAGGTGAAAAAGGGTAACCTCAAAGCTAAAAAGCCCAAGAAGGGGAAGCCCCA
TTGCAGCCGCAACCCTGTCCTTGTCAGAGGAATTGGCAGGTATTCCCGATCTGCCA
TGTATTCCANAAAGGCCATGTACAAGAGGAAGTACTCAGCCGCTAAATCCAAGGTT
GAAAAGAAAAAGAAGGAGAAGGTTCTCGCAACTGTTACAAAACCAGTTGGTGGTGA
10 CAAGAACGGCGGTACCCGGGTGGTTAAACTTCGCAAAATGCCTAGATATTATCCTA
CTGAAGATGTGCTCGAAAGCTGTTGAGCCACGGCAAAAAACCCTTCAGTCAGCAC
GTGAGAAAACTGCGAGCCAGCATTAACCCCGGGACCATTCTGATCATCCTCACTGG
ACGCCACAGGGGCAAGAGGGTGGTTTTCTGAAGCAGCTGGCTAGTGGCTTATTAC
TTGTGACTGGACCTCTGGTCCTCAATCGAGTTCCTCTACGAAGAACACACCAGAAA
15 TTTGTCAATTGCCACTTCAACCAAAATCGATATCAGCAATGTAAAAATCCCAAAACA
TCTTACTGATGCTTACTTCAAAAAGA

Sequence ID 1208

CCCTATACCTTCTGCATAATGAATTANCTAGAAATAACTTTGCAAGGGAGAGCCAA
20 AGCTAAGACCCCCGAAACCAGACGAGCTACCTAAGAACAGCTAAAAGAGCACACCC
GTCTATGTAGCAAAATAGTGGGAAGATTTATAGGTAGAGGCGACAAACCTACCGAG
CCTGGTGATAGCTGGTTGTCCAAGATAGAATCTTAGTTCAACTTTAAATTTGCCCA
CAGAACCCTCTAAATCCCCTTGTAATTTAACTGTTAGTCCAAAGAGGAACAGCTC
TTTGGACACTAGGAAAAAACCTTGTAGAGAGAGTAAAAAATTTAACACCCATAGTA
25 GGCCTAAAAGCAGCCACCAATTAAGAAAGCGTTCAAGCTCAACACCCACTACCTAA
AAAATCCCAAACATATACTGAACTCCTCACACCCAATTGGACCAATCTATCACCC
TATAGAAGAACTAATGTTAGTATAAGTAACATGAAAACATTCTCCTCCGCATAAG

Sequence ID - 1209

nt: 620

30 CTCTCCTGTCAACAGCGGCCAGCCTCCCAACTACGAGAATGCTCAAGGAGGAGCAG
GAAGTGGCTATGCTGGGGGCGCCCCACAACCTGCTCCCCCGACGTCCACCGTGAT
CCACATCCGCAGCGAGACCTCCGTGCCCCGACCATGTCGTCTGGTCCCTGTTCAACA
CCCTCTTCATGAACACCTGCTGCCTGGGCTTCATAGCATTTCGCCTACTCCGTGAAG
TCTAGGGACAGGAAGATGGTTGGCGACGTGACCGGGGGCCAGGCCTATGCCTCCAC
35 CGCCAAGTGCCTGAACATCTGGGCCCTGATTTTGGGCATCTTCATGACCATTTCTGC
TCGTCATCATCCCAGTGTTGGTTCGTCCAGGCCAGCGATAGATCAGGAGGCATCAT
TGAGGCCAGGAGCTCTGCCCCGTGACCTGTATCCCACGTACTCTATCTTCCATTCTCT

- 251 -

CGCCCTGCCCCCAGAGGCCAGGAGCTCTGCCCTTGACCTGTATTCCACTTACTCCA
CCTTCCATTTCCTCGCCCTGTCCCCACAGCCGAGTCCTGCATCAGCCCTTTATCCTC
ACACGCTTTTCTACAATGGCATTCAATAAAGTGTATATGTTTCTGGTGCTGCTGTG
ACTT

5

Sequence ID 1210

TTCGTAATTAGAATACTGTTTGGACTTGCTCAACAAGCACCTTATCTTAACAAAAA
GTAACCTTATAGAAAAGGGAGACATTCATTTAACTTCAAGCCCATATTATTCTTAAA
AGCTGACTCTTGAAATAGTATTTATTGAGTCATAGTGGAGTCATGGGACTTTTTTAA
10 GGGCCGGAAGGGACTATTTAGATCATCCAGTCCCACCCTGTCATTTTATGGAGGAG
GAAACTGAGGCCTAGATAAGATAACCAGTTAGTGGGTCCACTGACCTTTAGGACAG
TAGTCTATCCGTAAGAGACAACATGGAGAAAGAAATACAACGTTTTTATAGTGAAT
TATCATCTTACAAAGAATATTCTTCCCATATCGCACTTTTAAAAAGTGGGTACCTT
AGTCAAATAGGAGAAAAAACCACTTGAGTAGTTTCATCCTCAGGTTTTAGGTGAGG
15 AAAC TGATACTCAGATTAAATAACTTTAAGCACACAGAGCCTGAATGATAGTCTTA
TTTGAGCTCATCTGTGCTTTTAATGTGTACTACGTTAGGTGTTTTCACTTGCAATTT
CCTTTAGTCTTATTTGAGCTCATCTGTGCTTTTAATGTGTACTACGTTAGGTGTTT
TCACTTGCAATTCCTTGTTTGACGTTGACAATAAATCGTGAAGCTGCCTTATCTAA
GGAAGTCCTAAAGTAAATCATTGGAACACA

20

Sequence ID 1211

CCATTGTGTTGGNACCCGGGAATTCGCGGCCGCGTCGACGGAGTTTTACCTTATTA
CACTTTAATCTCTGGATTTACCCCATCTCATTTCTCTTTTAGGAAAACGTTTGTA
TGTGGTGGCTTTGATGGTTCTCATGCCATCAGTTGTGTGGAAATGTATGATCCAAC
25 TAGAAATGAATGGAAGATGATGGGAAATATGACTTCACCAAGGAGCAATGCTGGGA
TTGCAACTGTAGGGAACACCATTTATGCAGTGGGAGGATTTCGATGGCAATGAATTT
CTGAATACGGTGGAAGTCTATAACCTTGAGTCAAATGAATGGAGCCCCCTATACAAA
GATTTTCCAGTTTTTAACAAATTTAAGACCCTCTCAAATAACAGGCTTAGTGATGT
AATTATGGTTAGCAGAGGTACACTTGTGAATAAAGAGGGTGGGTGGGTATAGATGT
30 TGCTAACAGCAACACAAAGCTTTTGCATATTGCATACTATTAAACATGCTGTACAT
ACTTTTTGGGTTTATTTGGAAAGGAATGCAAAGATGAAGGTCTGTTTTGTGTACTT
TTAAGACTTTGGTTATTTTACTTTTGGAAAAGAATAAACCAAGAATTGATTGGGC
ACATCATTTCAAGAAG

35

Sequence ID - 1212

nt: 374

AGAGCAGCAGCCATGGCCCTACGCTACCCTATGGCCGTGGGCCTCAACAAGGGCCA

- 252 -

CAAAGTGACCAAGAACGTGAGCAAGCCCAGGCACAGCCGACGCCGCGGGCGTCTGA
CCAAACACACCAAGTTCGTGCGGGACATGATTCGGGAGGTGTGTGGCTTTGCCCCG
TACGAGCGGCGCGCCATGGAGTTACTGAAGGTCTCCAAGGACAAACGGGGCCCTCAA
ATTTATCAAGAAAAGGGTGGGGACGCACATCCGCGCCAAGAGGAAGCGGGAGGAGC
5 TGAGCAACGTACTGGCCGCCATGAGGAAAGCTGCTGCCAAGAAAGACTGAGCCCCT
CCCCTGCCCTCTCCCTGAAATAAAGAACAGCTTGACAG

Sequence ID - 1213

nt: 567

GAATTATTGACTTTGAATTGCATTTTCAGTACCATGAAGTCAAAGTCAGTGGTGTAT
10 TTGCTCATTTGTTTCATTCTTTCTTTTCCACCAACATTACTGCCTGCAGAGCCAGAG
GTGAGTGCAGAAATCCTGTCAATTCGTCACTTGTGGACAACCTGCAGCTTGCCACA
GCCTACAGTTCCACCACTGTGACCTCTGAAAACCTCCTGAACAAAAGGAAGGAGAC
TTGGAAATCCTGAATGGGCTTGGAGACATTAAGGGAGAAGTGCCTCCCTGGACCAA
GGCAGAATTCAATAGAACCAGCAAGAAATTTTCTATGAATGGGAAAGCAGGTGGC
15 AGGGGGCAGGGGTGGAAAAGCTTTGTACAGGAATTGTGGAAAAGCTTTTGCATTAT
CTCTAGTCTGAAAGTCACATTTCTCAGTTCTTTTCCACTCTCTTCTGTCAACTTGC
TGTGAGTAAATGACATCTGTCACCTGTGACACGGGCCAGGGACTATCACCATATGG
CCCCCACACATTATCTAGTACCAGCCTGCCTGGGCCATGCCTTTTCCAGTCACTGT
ACCAGCC

20

Sequence ID - 1214

nt: 620

CTCTCCTGTCAACAGCGGCCAGCCTCCCAACTACGAGAATGCTCAAGGAGGAGCAG
GAAGTGGCTATGCTGGGGGCGCCCCACAACCCTGCTCCCCCGACGTCCACCGTGAT
CCACATCCGCAGCGAGACCTCCGTGCCCCGACCATGTCGTCTGGTCCCTGTTCAACA
25 CCCTCTTCATGAACACCTGCTGCCTGGGCTTCATAGCATTGCGCTACTCCGTGAAG
TCTAGGGACAGGAAGATGGTTGGCGACGTGACCGGGGCCAGGCCTATGCCTCCAC
CGCCAAGTGCCTGAACATCTGGGCCCTGATTTTGGGCATCTTCATGACCATTCTGC
TCGTCATCATCCCAGTGTGGTTCGTCCAGGCCAGCGATAGATCAGGAGGCATCAT
TGAGGCCAGGAGCTCTGCCCCGTGACCTGTATCCCACGTACTCTATCTTCCATTCTC
30 CGCCCTGCCCCCAGAGGCCAGGAGCTCTGCCCTTGACCTGTATTCCACTTACTCCA
CCTTCCATTCCCTCGCCCTGTCCCCACAGCCGAGTCCTGCATCAGCCCTTTATCCTC
ACACGCTTTTCTACAATGGCATTCAATAAAGTGTATATGTTTCTGGTGCTGCTGTG
ACTT

35 Sequence ID 1215

CACAAGATAGAATGGTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAATTTTAA
GTGACAGTGCCATAGTTTGGACAGTACCTTTCAATGATTAATTTTAATAGCCTGTG

- 253 -

AGTCCAAGTAAATGATCACTTTATTTGCTAGGGAGGGAAGTCCTAGGGTGGTTTCA
GTTTCTCCCAGACATAACCTAAATTTTACATCAATCCTTTTAAAGAAAATCTGTAT
TTCAAAGAATCTTTCTCTGCAGTAAATCTCGCAGGGGAATTTGCACTATTACACTT
GAAAGTTGTTATTGTTAACCTTTTCGGCAGCTTTTAATAGGAAAGTTAAACGTTTT
5 AAACATGGTAGTACTGGAAATTTTACAAGACTTTTACCTAGCACTTAAATATGTAT
AAATGTACATAAAGACAAACTAGTAAGCATGACCTGGGGAAATGGTCAGACCTTGT
ATTGTGTTTTTTGGCCTTGAAAGTAGCAAGTGACCAGAATCTGCCATGGCAACAGGC
TTTAAAAAAGACCTTTAAAAAGACACTGTCTCAACTGTGGTGTAGCACCAGCCAG
CTCTCTGTACATTTGCTAGCTTGTAGTTTTCTAAGACTGAGTAAACTTCTTATTTT
10 TAGAAAGTGGAGGTCTGGTTTGTAACCTTCCTTGTACTTAATTGGGTAAAAGT

Sequence ID - 1216

nt: 484

CAACCTTAGCCAAACCATTACCCAAATAAAGTATAGGCGATAGAAATTGAAACCT
GGCGCAATAGATATAGTACCGCAAGGGAAAGATGAAAAATTATAACCAAGCATAAT
15 ATAGCAAGGACTAACCCCTATACCTTCTGCATAATGAATTAAGTAGAAATAACTTT
GCAAGGAGAGCCAAAGCTAAGACCCCCGAAACCAGACGAGCTACCTAAGAACAGCT
AAAAGAGCACACCCGTCTATGTAGCAAAATAGTGGGAAGATTTATAGGTAGAGGCG
ACAAACCTACCGAGCCTGGTGATAGCTGGTTGTCCAAGATAGAATCTTAGTTCAAC
TTTAAATTTGCCCACAGAACCTCTAAATCCCCTTGTAATTTAACTGTTAGTCCA
20 AAGAGGAACAGCTCTTTGGACACTAGGAAAAAACCTTGTAAGAGAGAGTAAAAAATT
TAACACCCATAGTAGGCCTAAAAGCAGCCACCAATT

Sequence ID 1217

GACAGGCGGGGGCCAGCGGCCGGGTGAAGGCCGGGTGGCTCTGTGAATCAAAGGA
25 GAGTCCCAGAAAACCTGTGACTGTTGAAGAAAATTCATCTGTGAATTTTTATATTC
AAGGAGTCAGTATTTATATTCATCTTTTAAACTGGGAAGATTTATATTTTACTTTA
AACTTCTTGATAATAATTTACAATGAATGGACACAGTGATGAAGAAAGTGTTAGA
AACAGTAGTGGAGAATCAAGGTAAGTAAGCACTTTGTTATCAATTGTTTACTATGA
AGAGAGTTGAAAACCTTGACTTTTTTCTTTATTGTTATTGTTGTTATTTAGTTTTCC
30 TCATAGGTAGCAGAGTTTTTCAGGTTTTCTCTTAGCTATCCAAATACTAAAAAAT
TCTGATATACGAACCTTTTTTTCATAATACAGGTTTAAATTATATTTTTTCATTGAGA
TACACAGTAGATCTTAAATATAGAAAGTTTTTGTTTACTTAAATCTATTTGGAAGT
TTATATTTGAGCTAATAATTAAGCTGGAGCATGTATAATAGATTTAAATTGTTTTG
ACTGTTAGTGAAATTT

35

Sequence ID 1218

CTCACTTGGTGGGTGAGCCTCCAATGACTACACCCAAGGAGGATTTAACACAGGGA

- 254 -

TTTTATGACTTGCAACAAGTCAGGAGGACATGGGGTTGGGGTAGTTTCAGCAGTGCC
TGTCTGAACAAAGGTGAAAATTGGGCTTTTATTGGGCTGATCAAGGGGGAGTAAAG
GCAGCCAGGAGCAGTCGCCTGTCATGCTTCTACCTATATTGCATGTATAGAAAAGG
GAAAATAAACTCCTTCCTGGGCAGGGTTTTAGTATGCTAAGGAGGGGAGTTATTCA
5 ACTTCAATCCAACCTCAAGCATCAGCATTGCTGCGTCCATCCCAGTTTTGTGTTTGCT
GGGGCTGAACTTCTTCCTATAACTTTTTGAAACAACAAGAACTCAAGGTGTGACAG
TTACAAGTGGGCCCTTTTTTCACAGTGTGTACCTAAACACGTGAGGACCCTGGATTA
CAGAATGACAGACTCGAAGTGAAGTCAAGTTCGGTTGTTTCATCTTTAGATGGTAAA
GATGGCTGTACGTACTATCCTTGCTTATTTCCAATCTATTGTTTAAACTCTTGAT
10 ATGTAATACCGCAGAGGCTAGAGATAAACCTTTGACCAAATGAGTGAATTCAAGT
AATCCATTACTAATGTGATCTGGAAACAAACATGGTGTGGAATGTGCATATGT

Sequence ID - 1219

nt: 559

CTTGGCAGCTCCGTTATGTGCCCAGCTCTTTGCAAGGGCATACTGGGAAATGAGTG
15 GAGATAAAGGACCCAATCATAAGCATTTTACAGTATGGATAACCCATTTTAAAAAG
GTAAACTGAGGCACAATGCAATTTTTTTTTTTTTTTAAGGAGTTTATTTGAGCAAA
CAGTGATTCATGAATCAGGCAGCACCAAACCAGAAGGAGGCTTTGCTGAANAAGGA
TGAGGGACAAGCATTTATAAAGTGAATGTAGATGTAATACAAAGAAAATATTTGAA
CCGGGTGCGGTGGCTTACACTTGTAATCCCAACACTTTGGGAGGCCAAGGCGGGCA
20 GATCACAAGATCAAGAGATCGAGACCATCCTGGTCAACATGGTGAAACCCCATCTN
TACTAAAAAATACAAAAATTANCTGGGCGTGGTGGTGCCTGTAGTCCCAGCT
ACTTGGGCGGCTGAGGCAGGANAATTGCTTGAACCCGGGAGGTGGAGGTTGCAGTA
AGCCGAGATTGCACCATTCGACTACTCCAGCCTGGTGACAGAGAGAGACTCCATC

25 Sequence ID 1220

GANNNGTGCGATANNATGNNTGTCTTTTTTTTAAAGTNTTTCNNATNGNAGNGAAN
CCCCNNANNTNNCATAANGAGAGATNACTACNGTACANATAGNGNCANACNGATA
GTAGTANCAANATTGTNTTAGCTANATNANTCAATAGATATCNAGATANANAANA
NCNNGGATATACAGCGATGTNTNANNGGNNNNNNNANGGAACGAACATCNACNTTA
30 ANNATAAGCTNGNGGAGAGAGACANGTANGTTATANANNAGAATNGNAGTAGGNGT
GATCATAATAGNNNNNANNTANTATATANGATNTTANTGNNCTNTNNTNNGTTTAT
CNNAATNTCTATNCTNGAGAGNAGCNNNATNNNNAGGCGANGANATTGGGNNNTN
CTCNTNATAGANANCTGGTGTCNNANAANTACNTCATCTATTNANCTCTCACNANA
TGGNANNATANAGNAGNGNNNTNNANAGGANTANGCATAGNGNNTNNCTNAAACAA
35 AANNNATAAGANNTCTCGNNAANANGGGCCTNTNNTNTAGCGAGGNNTTANTTTNT
ATANTTNTTCNCTCTTNNAATANNTANGATANATGANCTNGNNGTGATANATANNN
NNTACNGTNAANNTNTANTCNTATAATAGATANAAATATAGGATNTTNTCTGTCN

- 255 -

GGTNGAANANTTNNTNCNNTTTTNAATAATGNTGTTAGNGACNGNGNTNTNANANN
NNTTAGAAAGGTACTCTATATACTNNTATGNTNCGGCNNATAATANAACAGATGTT
TGTATNAATATNAAANAAGGTCNNTTTCGNCAGAGAANNNTGNCTGGTNATAGAA
TTAGCATAANTTANNTANTATGATNNANTNNTNCTACNANTNTTAGCNNTTNGCAG
5 NAGTCAT'TNNGNATNTATNNGNNTANTAGTNANTTGGGNC'TNNTNCAGANTATAT
TNTGNGAANATGAANNACGNANTCCTNNGNANTATNATNNTGANTANGANAANCN
ANANNTNTTNTANNANTGNCTATANATTGCCNNGATANATTNTNNAATGAANCGA
TAGCCCGCNCTAAGGANNTNNGTNANNTAAANNTCTCAGATAANNTACNTNTTNT
TATTAANCNANNATCACANTATANCNGNGACANNNGCGANANTATATGTATGNNAN
10 TATNACNGNTCCNNCCGNGAANNTANTCNTANNAGGCATT'CNNGNAGC'TNTTCT
NCTAGACNATTTNNANTGAAANNATGCNGNNAAAAACGACNNNCTTNAANTTNTGT
CTACANTCCGCNNTNTTNTACAGATNGCAGNTAAGNNNANTNANNGCTCTCANCT
NGCTNNNACT

15 Sequence ID - 1221 nt: 741
AAGCAGAANTNTCTCTAAAAACATTATCTCCTTAAATCTTGAGGTGCATATNAGA
GCCACAGGCAATCTCTGACATATAAAATTCAGTACAGGCCTTTCAAATTTGGCAT
TTCAGTGGTACAATAACAACAAGATATATAAATACTGTACAGTGCCTAGACAT
TCCAGTAAGAACCATTATTTTCTTTAATGTAGAATGATTAATACATATTCTACAAG
20 GGGCAGTAAGGTTAGTAATTCTATAGGGTATGTCCCGACATAATTTTCAAATTGTA
CAATAACACAAACAAC'TTTGTTAAGGCCATGTTTTATTTGCTGATTAATGGACAAA
AGGCAATGTAATTTATTTTCAAGTATTTTCTTGAAAGTCTGTGCTCATAAAAATCA
TGAAAAGTTGGAAAGACTGTTAAATCACTGAAACTTCAAATATATCTTACACAATC
TTGTTTGTACAAAAATACAAGTTAAATATAAACATAAAGCAATCATGGTAATTTTA
25 TGCAAATCTGTTTTATGTGATCATCAGTTATATATAAAAGTTTCTCAGTTCTGTTA
TTTGTGAAAAGATCAATACCAGATTGAATGACTACCTATTTGGCAAAGGGCCCTAAA
AAGCTTACTTTAGCACTCATCTTTTACATGGTTAAATGCATTTCCCTAATTTGAGAT
CACCTAAACACTGGAAAAGAAAAAAATGAAAGGGCAGTATGTCCATAAACCAACA
AATAATTTGGCTG

30 Sequence ID - 1224 nt: 485
CGAAATTTCTTGTGACACAGAGGAAGGGCAAAGGTCTGAGCCCAGAGTTGACGGA
GGGAGTATTTTCAGGGTTCAC'TTCAGGGGCTCCCAAAGCGACAAGATCGTTAGGGAG
AGAGGCCCAGGGTGGGGACTGGGAATTTAAGGAGAGCTGGGAACGGATCCCTTAGG
35 TTCAGGAAGCTTCTGTGCAAGCTGCGAGGATGGCTTGGGCCGAAGGGTTGCTCTGC
CCGCCGCGCTAGCTGTGAGCTGAGCAAAGCCCTGGGCTCACAGCACCCCAAAGCC
TGTGGCTTCAGTCCTGCGTCTGCACCACACAATCAAAGGATCGTTTTGTTTTGTT

Sequence ID 1228

nt: 741

30 AAGCAGAANTNTCTCTAAAAACATTATCTCCTTAAAATCTTGAGGTGCATATNAGA
GCCACAGGCAATCTCTGACATATAAAATTGCAGTACAGGCCTTTCAAATTTGGCAT
TTCCTGTTTACAATACAACAACCAAGATATATAAATACTGTACAGTGCCTAGACAT
TCCAGTAAGAACCATTTATTTTCTTTAATGTAGAATGATTAATACATATTCTACAAG
GGGCAGTAAGGTTAGTAATTTCTATAGGGTATGTCCCGACATAATTTTCAAATTGTA
35 CAATAACACAAACAACCTTTGTTAAGGCATGTTTTATTTGCTGATTAATGGACAAA
AGGCAATGTAATTTATTTTCAAGTATTTTCTTGAAAGTCTGTGCTCATAAAAATCA
TGAAAAGTTGGAAAGACTGTTAAATCACTGAAACTTCAAATATATCTTACACAATC

- 257 -

TTGTTTGTACAAAAATACAAGTTAAATATAAACATAAAGCAATCATGGTAATTTTA
TGCAAATCTGTTTTATGTGATCATCAGTTATATATAAAAGTTTCTCAGTTCTGTTA
TTTGTGAAAAGATCAATACCAGATTGAATGACTACCTATTGGCAAAGGGCCCTAAA
AAGCTTACTTTAGCACTCATCTTTTACATGGTTAAATGCATTTCCCTAATTTGAGAT
5 CACCTAAACACTGGAAAAGAAAAAAATGAAAGGGCAGTATGTCCATAAACCAACA
AATAATTTGGCTG

Sequence ID - 1231

nt: 203

TTGAGGAAGGGTCTACTGTCTTTTTAAATGGCACAATTTTAAGAGGTTTGAGAGGT
10 ACAGTCCCTTAACCTGCCACGGGAGAGGGGCCCCCAAACTTTCTTCCCCCACACT
TCTGGTTTTTCTGTGTGGAGGGGGAGCAGGGATATCTAAGCTGTGGTGTGAAAGGGT
AGGAGAGATGCTGGAGGTGGGGGTGCTGTGTTCTA

Sequence ID 1239

15 TTTCTCGGGAAGCGCGCCATTGTGTTGGTACCCGGAATTCGCGGCCGCGTTCGAC
ATTTTTTTTTTTTTTTTTTTTTTTTAGAATGATTAACAATTTATTGAGTTTTATTTATC
TACAAAAATATAGCAATACAGNGAATTACCAAATCCTAAATATTCAGTACCTGA
ACTGGCTACAACACCGNGTGCACACCCAGTTCCTGCAGAATCTCTTGCAGATATGG
GAGAGTCAGCCAGTGAAAAGATCCATTTCTTGGGAATCCTTGTCAACAAGACCAGT
20 TCAGAAATCCAGGATATATAGAAGCCTACTGTAATTTAAAAACAGTAACAAAAACC
CCAACAAAACCCAAATCAACAAAGACCAAGATAAAGNGTGATAAACATTAATTGT
AATGGTTTTTCCTTTACATGCAATACATGCATTTTAAATCACTAAGAAACACGAAA
TTTTGTAGAGCAAAGTTTGN GTTTCACGTAAGTGCAAATGAATATATATTTTATTT
TTTATACTATTAAATTATATATATTTTTTCCATACAAAAGCACACAGTGTTAATCT
25 ATAAATGACATCCAAGTGGATGATGATTGTTTTTGCATGTCCCCCTGCTTAGATT
TTTTTAAAATATATAGTCAAAAATTAACATCCTTCTTTAAAAATACAGAAGGGAAA
AANGGGCAAAAAAAAAAATCTAGACTCGAGCAAGCTTATGCATGCATGCGGCCGCA
ATTCGANCTCGNCGACTTGGCCAATTCGCCCTATAGNGAGTCGNATTACAATTCA
CTGGGCCGNCGN TTTACAACGTCGNGACTGGGAAAACCTGGCGTTACCCNNCTNA
30 TCGNCTTGNAACAATNCCN TTTNGCCAGNGGGG

Sequence ID 1255

TCACTTCGTATNGAANCTGTTTGGACTTGCTCAACAAGACCTTATCTTAACAAAAA
GTAAC TTATAGAAAAGGGAGACATTCAATTTAACTTCAAGCCCATATTATTCTTAAA
35 AGCTGACTCTTGAAATAGTATTTATTGAGTCATAGTGGAGTCATGGGACTTTTTAA
GGGCCGGAAGGGACTATTTAGATCATCCAGTCCCACCCTGTCAATTTTATGGAGGAG
GAAACTGAGGCCTAGATAAGATAACCAGTTAGTGGGTCCACTGACCTTTAGGACAG

- 258 -

5 TAGTCTATCCGTAAGAGACAACATGGAGAAAGAAATACAACGTTTTTATAGTGAAT
TATCATCTTACAAAGAATATTCTTCCCATATCGCACTTTTAAAAAGTGGGTACCTT
AGTCAAATAGGAGAAAAAACCACTTGAGTAGTTTCATCCTCAGGTTTTAGGTGAGG
AAACTGATACTCAGATTAAATAACTTTAAGCACACAGAGCCTGAATGATAGTCTTA
TTTGAGCTCATCTGTGCTTTTAATGTGTACTACGTTAGGTGTTTTCACTTGCAATTT
CCTTTAGTCTTATTTGAGCTCATCTGTGCTTTTAATGTGTACTACGTTAGGTGTTT
TCACTTGCAATTCCTTGTTTGACGTTGACAATAAATCGTGAAGCTGCCTTATCTAA
GNAGTCCTAAAGTAAATCATTGGAACACATGTANCCAGTTTGTGTGTTTTTATTTGC
CAGGTNTCAAATATAACTGAAAACCCATGCTAACTGACTNATTTTAAAAGNTGTNT
10 GGGGCATGAAANGATTGCTCTGCCTGGGCGGGNGGTTNANCCTGNGTCCCCCNTTT
NGGAGNCCACCCANGANGCGATATTTNAGGGNNGATTNAAACCCCTGGCACGNGN
NAACCCCNNTTTTAAANANAAAAANANCGGNG

Sequence 1256

15 TTGTTGTTGGTACCCGGGAATTCGCGGCCGCGTCGACGGAGTTTTACCTTATTACAC
TTTAATCTCTGGATTTACCCCATCTCATTCTCTTTTAGGAAAACGTGTTTGTATGT
GGTGGCTTTGATGGTCTCATGCCATCAGTTGTGTGGAAATGTATGATCCAAC TAG
AAATGAATGGAAGATGATGGGAAATATGACTTCACCAAGGAGCAATGCTGGGATTG
CAACTGTAGGGAACACCATTTATGCAGTGGGAGGATTTCGATGGCAATGAATTTCTG
20 AATACGGTGGAAGTCTATAACCTTGAGTCAAATGAATGGAGCCCCCTATACAAAGAT
TTTCCAGTTTTTAACAAATTTAAGACCTCTCAAATAACAGGCTTAGTGATGTAAT
TATGGTTAGCAGAGGTACACTTGTGAATAAAGAGGGTGGGTGGGTATAGATGTTGC
TAACAGCAACACAAAGCTTTTGCATATTGCATACTATTAAACATGCTGTACATACT
TTTTGGGTTTATTTGGAAAGGAATGCAAAGATGAAGGTCTGTTTTGTGTACTTTTA
25 AGACTTTGGTTATTTTTACTTTTTTGGAAAAGAATAAACCAAGAATTGATTGGGCACA
TCATTTCAAGAAGTCCCCCTCTCCTCCACATTTGTTTTGCCAATTTGCACATTAAAT
GACTCTTCCCTCAAATGTGTACTATGGGGTAAAAGGGGTAGGGNTTAAANATGTAA
ACAGTTGGGTTTTTTTAAGGGNCCTTTTTTCATAACTGGAACACTCTNTACAAGGNTN
CTTNTTAAATAAATAACTTGACTTTTTTTGTTTTNTAAANGNANCTTCNTGCTTCCA
30 TAAAAAATAAATTTAANTNGNCANCNTNGCTGCTGCGNCCANTTNGCTNGNCCNT
GGCATTCCCTAGGGANGNTNAATANTGGCENNNTTAACNNGGCNGNAACNNNNNCCA
NT

Sequence ID 1331

35 GGGCGATGCATGCTTTATTAAGGCTCTTGTTTCACCTGGCAGTGTACTGTATCAAC
GTATAATACAGAAAAAAATCTCTTTAAGGTCCTCCTTCACAAAGACATAGAGTGA
AACTCCCTTTACATGTCAGTATTTGTTCAACACTTTAGGCAACTTGACTGTCAGTG

- 259 -

TTAAAATGGAAAACAGGAAAATGGAAAATCTGACCAATTCTGCCACCTTGAGACT
TTCATATAGACCTTGCACAACAATTGTATAGATCACACACCGGCTGTATTTAATAT
GTAACATTTTTCACACATATTAAAGATACAGAAGTATTAAAAAACCCCAATGTTAA
TGTATTTGCTTAAAAGGCACAAGTTTCACATATCTGTCTAGCTATCTGTTGGTAAT
5 ACAGAAAGTATACTACTTTTTTTAAAAAAGTGGGCAGAATTCTTGTGTATGTATATT
TGTGTGTACAGTATGTGTATGTGTGTATATATATATATTATATATATAGATAATAT
ATAAATATTTTTTTTTTAAGGAGAACTAGAATGTTTAGCTAGAAAATTCCACAGCCT
GTGAAGAAATATTTCAAATGGCCATAAAGGAGGTAAAAATGAAAACCATAACCTA
ACTTTTATAGAGGCTTTATCTTTAATTTAACGATGTGCGGAGGACTTTCTTGCTTG
10 AATCTGTTCCGGGCTGTCTGCTCTGTCCATCAAATGGGCAGGTCTGGGAATGAGGC
ACCTTCGGCCGTTTCAAAGTGGCCTGAACAGAATGCTGGAACCCAGGCTGGACTCG
GAC

Sequence ID 1332

15 CAAACCTGCATGTTCTGCACATGTATCCAGGAACCTAAAAAAGATAGTT
TGTGTGTCTTAATTGAATAATAGTAGATTTATAGATTAAAGATCTATGGGTTTTTA
ATATGGATTAGAAATCTGTGGGTTTTTGATATGGATTAGAAATCTGTGGGTTTTTA
ATATGGATTGGAAATCTGTGGGTTTTTAATATGGATTAAAAACATCTGTGGGTTT
TTAATATGGATTAAACATCTGTGGGTTTTTAATATGGATTAAACATCTGGGTTTTT
20 AATATGGATTAAACATCTGTGGGTTTTTAATATGGGTTAAAAATCAAAGAAAATG
AACTATTTGCTCCAGTGCAGGAAAATACAGGCAATACTGGATACAATTAGATGGTC
AGGAGCGATAACCCGGTTGCCATTGTTTGAAGAAGAGAATAAGGTGCTAGCATTC
TATCCGTAGATAATTTGACAGCTAGGAAATAGGGGGAGTCTTCTATGTAGTTAGTG
AAGGCTAAATGAACTATTATATGCAGTTATCGTAGAAGAGTACTCAAAAAATCTG
25 TAAAAAATAAAGAAAGGCCGGGCGCGGTGGCTCACGCCTGTAATCCCAGCACTTTG
GGAGGCCGAGGCGGGTGGATCATGAGGTCAGGAGATCGAGACCATCCTGGCTACCA
NGGTGAAACCCCGTCT

Sequence ID 1335

30 CAAGACTCCATCTCAAAAAAATCTACAGTGCTGAGTATATAAAATTAT
TAACACATTTTACAACAATATGTGTTTGTGGAGTTAAATATTTTTTGTCTTTAAAA
CAGGTAATTTTAGTGCATACTTAATTTGATGATTAAATATGGTAGAATTAAGCATT
TTAAATGTAAATGTTTGTACATTGTTCAAGAAATAAGTAGAAATATATTCCTTTG
TTTTTTATTTAAATTTTTGTTTCTCTGTAACTAAAAGAACACGAAGTAATTGGTC
35 ACAATTAAGTGGTGTAACTGCCAAATATGGGTAAATAAGGGAAAATTTTGTTTAA
TATTTAGTCCTTCTGAGATGGCTTGAATATTTGAATTTTGTGTACGCTCTATACTG
GGTAGTCACAAGTCTTATAAACACTTTAGAGGAAAGATGGATTTTCACTCTGTATTT

- 260 -

TTAAACATCATTTATTTTAAATCTGGTGCTGAAAAATAAGAAAAAATTAAACTGC
ATTCTGCTGTTCTTCTTTAGAACATTCCTGCGTAAATACTGCTGTAATACTGTCA
TGCAAAGTGTATCCTTTCTTGTCGTATCCTTTTTGGGGCAGTGTTTTTTTGTTTTT
TTCTAGAAATGTTTGTCTTCCCCACCTGTTGATCCAGGTTAAGGAATACTTTT
5 TTACACTTTATTCAAA

Sequence ID 1336

CTTTTCCTCCCGCTGTCCCCACGGAGGGGACTGCTCTCCCCCGCTGCATCCTTTC
TGTGAGGTACCTTACCCACCTCAGCACCTGAGAGGGTGAAATAGAATTCTAACCTC
10 GACATTTCGGGAAGTGTTTTTGAGAAGTCTCGGTTCGGTAAGGGAAGTCTTCCAAGTC
CGTGCAGCACTAACGTATTGGCACCTGCCTCCTCTTCGGCCACCCCCCAGATGAGG
CAGCTGTGACTGTGTCAAGGGAAGCCACGACTCTGACCATAGTCTTCTCTCAGCTT
CCACTGCCGTCTCCACAGGAAACCCAGAAGTTCTGTGAACAAGTCCATGCTGCCAT
CAAGGCATTTATTGCAGTGTACTATTTGCTTCCAAAGGATCAGGCCCTGAGAACAA
15 TGACCTTATTTCTTACAACAGTGTCTGGGTTGCGTGCCAGCAGATGCCTCAGATAC
CAAGAGATAACAAAGCTGCAGCTCTTTTGATGCTGACCAAGAATGTGGATTTTGTG
AAGGATGCACATGAAGAAATGGAGCAGGCTGTGGAAGAATGTGACCCTTACTCTGG
CCTCTTGAATGATACTGAGGAGAACAACCTCTGACAACCACAATCATGAGGATGATG
TGTGTTGGGGTTTCCCAGCAATCAGGACTTGTATTGGTCAGAGGACGATCAAGAGCTC
20 ATAATCCCATGCCTTGCGCTGGTGAGAGCATCCAAAGCCTGCCTGAAGAAAA

Sequence ID 1337

CAAGAACTCTGGGACATTTGCAAAGGGTATGGCATATGTGTAATGGGAATACCAGA
GGAGAGGAAAGACAGGAAGTCAAAAAAGAATTTTTCCAAATTAATGATAGGTTCC
25 AAACCACAGATGCAGGAAGCTTAAACACCAACAGGATAAATAAAACAAAATCTACG
CTTAAGCATATCATACTTAACCTGCAGAAAATTACAGACAAAGAAAAACACCAGA
GGGGAAGCTGGCAGAAACATAACCACCTATAGCGGAAGAAGAATAAGAATTACATCA
GACTTCCCTTCAGAAATCTTGCAAACAAAAAGATGTAGCACAATATTTAAAGTATT
AAAGGAGGCCGGGCCCGGTGGCTCGGGCCTGTAATCCTAACACTTTGGGAGGCTGA
30 GGCAGGAGGACCATGAGGTGAGGATCGAGACCATCCTGGTGATGGTGATACCCC
ATCTCTACTAAAAATACAAAAAATTAACCGGGCATGGTGACACGCACCTGTAATCC
CAGCTACTTGGGAGGCTGAAGCAGGAGAATCGTTTGAGCCCAGGAGGTGGAGGTTG
CAGTGAGCCGAGATCACATCACTGCACGCCTGGGCAACAGAGCGAGACTCCATCTC
AAAAAA

35

Sequence ID 1338

CGACCCGTTTTAGTCAGGATGGTCTCGATCTCCTGACCTCGTGATCCGCCTGCCTC

- 261 -

GGCCTCCCAAAGTGCTGGGATTACAGGCGTGAGCCACCGCGCCCGGCGTAAATCAG
GTTTTTTTAAATGTTTGCCAAACCTTATCACTGACTTTTATAACAAAATTATTTACT
ATAATCATTAGGGAATATTTAAGTTCTGCTAATACTTAAAATTGCAGAGTGCTAAA
ACCAGCAGTGAGTTTAGAATCAAGCTAAGCTTTATTGTTGCTACTATTTGAGGCAT
5 ATTAGTTGACTGGTGTTTCATATGCAAGGCAGTCTACTGGGTGCAACAAGGGTTAGA
AGGATATTTTTTAAAAACTGACCCTATTCTCAGGATGAAAATAATACACTAGTAAT
AGTCTGCTCTGTTGGTTAACTCCTCGTAAGGAGGTACAATTAAAATGCTGTAGTGT
TGCAAGGGAAGGAGAGGAAGAATCATATTCCTTCACTAGCAGGATCAAGAAAGCTT
TTATAGAAATATACAAAATCTTCACTTCTTGAAGGATTGGTAAAATTTAATAGCCA
10 ACATTGGGCACTTATTCACTCTCTGAGTAAATATTTATTGCATGCTTATCTTGTAT
CAAGCATTTGTGATGAAAGCACAGAATGAAAGAGGAGGGAGAATGTTTAGAGAATA
AGGGCTGAAACACAGATTTTGTAGGGAGCGTAGGGGAGACTGANAAGACAGGTTCA
GGTTAGTAAGGGCGCTCATATTTTGACCCTGAATGTAACTATGTGCACATCATGC
TAGCTATTCTAAATCAGGCATTTTCAAATGGAAGCAGGCACTGACATTTT

15

Sequence ID 1344

CGTGAAGGGTCTTTATGTATTAGTATTAGAGTGATCTTTTGATTATTTTCCTCACT
ATAAGGAAATTATTTCTCAGGATGAGCTGCCATAACATTCCACTGTCTGATGGCA
ATTTTAAAGCCTGAAATTGAAGCCCATGGCTAGGCTATGAGAACCCTAGTTCGTAT
20 AGTAAAGTTGATATCTTCTGGATGTATACTAATTTTAGGCTTTATTTTAAACTGC
TGGAAACTGAACTTAGACAAAAGTATTTTCAGGACATCATTTACAATGTTTAGCC
CTAAAGAGTCAAGCTGTGGGATTCTGAGTCTTTCATATGTTACAGCAGAACTTAA
AAGCAAGAGGAAATTGGCTGGGCACAGTGGCTCTGTAATCCCAGCACTTTGGGAGG
CTGAGGTGGGTGGATCATGAGGTCAAGAGATTGAGACCATCCTAGCCAACATGGTG
25 AAACCCCATCTCTACTAAAAATACAAAATTAGCTGGGCGTGGTGGCACACGCCTG
TAATCCCAGCTAGTCAGGAGGCTGAGGCAGGAGAATATCTTGAAGTTGGGAGGCAG
AGGTTGCAGTGAGCCAAGATTACATCACTGCACTCCAGCCTGGTGACAGAGCGAGA
CTCCGACT

30

Sequence ID 1348

CTGAAACTGCACTGAACCCACAGGTAGGTTACATCACAGGACAGAAATCTGAGGAG
CTGGAGAAAGCAAAAGAATAAAGGATGGGCTGACACCAGAAGGAATTAAAGGAATT
TTTATACTGAACTTCAATTACTTGTTCAATTTGAAGTTTGTTTTTTTAATGAACGTT
TTTGCTGTTACTTAAATATAGTGTTTTGAAAGTGTTTCAAATGTATTCAAGTTGGG
35 ATTTTCCATATTTTACTACAGTTCTGTCTTAGTATGTTTACCATAAAACACTTATC
ATTAAAGCTCACAAAGTGCTTTTTTTGTAATATGAGGATAAAATGAAGCCATATAAG
AATTTTTTTTATATCTGTACATTTAACCACATTTGAGCTTTAGCCAAAATATATAG

- 262 -

CTTTTTTTTTTCTGACCTGGCCAACGTATTATCCAGCAAACATCAACTGAAGCAAT
ATGGAAACACTTCCAAATGTTTGCCAATAATGCTATTAAGTGACTGATGTCAACAT
TAGTTACATGGCAAACCTAAAGAGGCATTATACATTTTTTAAAACACACTAACATATA
ACTGTAGATAATGTAAGGTTTATTTATATGCATATTTTCATAGTATATTTAAATGTT
5 TAAATATAAAAAAGGGTTTTTAAACACTTTTAATTTTTATCTTTGATTTTTTTTTAT
TGATATCTCTTTCCAGGCTACTAATAAAATTGCCAGAACTAACTATCAGGTAAAG
GTTAAGGCATCAATTGACAAGTAAGTTTTCTAATTTTCGTTTTGAATTACAATTCCA
AATGTAAGACTTTTAAAAATGAATGGCCTTTATTTTATAGAATAATTTTGACCTTT
TAAATTTACTTATCTAACATTATATAATGAATGTACTTCAAATATTTGACTTTGAA
10 GTCAACATTAACAAATTCATGGATCCTAATTAAAATTTACTATAAACTGGAATCA
TTTATTACTTCCTT

Sequence ID 1351

TTTTTTTTTTTTTAAAAGAGATGGGTTCTCACTATGTTGCCCATAAATGTTTATGAG
15 ATTAAGTTCATCTTTTTTATCTGAGTAGTATTTTATTGTATGAATATAACCACCATT
TATTTATCTGTTGGTTATTTCCAGTTTTGGGCTATAATCCAAAATGCTTTTTTCAA
ACAATAGGCTATATATCATTAATGTCCGTTTATCAGCAGTATAAAATATCTTACCA
TAAATATTAATAAAAAGAAGCATTATATATAAAATATAGATATTTCAAACCTACA
GAGGGCCTTTTAATGATTAAATATTTTGTCCCTACAAAAGGTCCAGGTAATTACA
20 CCCATGAGGTTAACCTGCCTTAGTGCAGGACTTAAAATAAGGCTTCTCCTGCCATC
TCTCTCCATTTGTAGAATGTGAAATTCTTTAAAATGCATCCTATATTAGGAATACT
ATAGCTGTGCACTGGTGTGTTCTCTTCTTTAACTCGGGACCGTATATATCTGC
TCAAATTGCCCAAGTATACATATGCTGCACTCCATCAAGTGTCAGGCCACATTCTA
TCAGCACAGCGTGACTGCCTATCAGTGACAATATAAGTGAGCTCTATTTGGATCCC
25 TCTTACCCTACCTTTTATATTTATGACAGCATTATCATAAACTCCAATATTCTTC
AATAACTTACATGTTTGTGTAGGATAAAATTATTACCCTCAATGAACTACAT

Sequence ID 1352

ACCAGCTTCTTCACAGGTTCCACGAGTCATGTCAACACAGCGTGTGCTAACACAT
30 CAACACAGACAATGGGTCCACGTCTGCAGCTGCAGCCGCTGCAGCTACTCCTGCT
GTCCGCACCGTTCCACAGTATAAATATGCTGCAGGAGTTCGCAATCCTCAGCAACA
TCTTAATGCACAGCCACAAGTTACAATGCAACAGCCTGCTGTTTCATGTACAAGGTC
AGGAACCTTTGACTGCTTCCATGTTGGCATCTGCCCCCTCCTCAAGAGCAAAAGCAA
ATGTTGGGTGAACGGCTGTTTCTTATTCAAGCCATGCACCCTACTCTTGCTGG
35 TAAATCACTGGCATGTTGTTGGAGATTGATAATTCAGAACTTCTTCATATGCTCG
AGTCTCCAGAGTCACTCCGTTCTAAGGTTGATGAAGCTGTAGCTGTACTACAAGCC
CACCAAGCTAAAGAGGCTGCCCAGAAAGCAGTTAACAGTGCCACCGGTGTTCCAAC

- 263 -

TGTTTAAAATTGATCAGGGACCATGAAAAGAACTTGTGCTTCACCGAAGAAAAAT
ATCTAAACATCGAAAAACTTAAATATTATGGAAAAAAACATTGCAAAATATAAAA
TAAATAAAAAAAGGAAAGGAACTTTGAACCTTATGTACCGAGCAAATGCCAGGTC
TAGCAAACATAATGCTAGTCCTAGATTACTTATTGATTTAAAA

5

Sequence ID 1353

ACATTCTGGAAAAGGCCAAAAGGGAGGAAGAACTGATTAGTGGTTAGCCCAGGGTTA
GAGTTGGGGAGAGGATATAATGAGGGAACCTTTTGTGGATTCTGTACCATGATTATG
ATTACACAAACCTATGCATACATTGAAACACATAGAACTATACGTTGAAAAAAGTG
10 AATCTGCCTGTATGTAAATTTAAAAGAAAAATATTTTTTTAAAAAACAGATGCTT
CTTAACACATTATCATCTATGTCAGTTTAAACAGTTAGTAGACTTAGGCCAGGTGTC
ATGGCTCACTCCTGTAATCCCAGTGCTTTGGGAGTCTGAGGTGGGACGATCTCTTG
AGACTAGGAGGGAGTTTGAGACAAACCTAGGCAATGTAATGAGACTCTTCTCTAC
AAAAAATTTTAAAGTTATCTGGACATGGTGGTGCCTGCCTGTAGTCCCAGCTACTT
15 GGGAGGCTGAGGTGGGAGGATTCCTTGAGCCCAGAAGTTCAAGGCTACAGTGTGCT
ATGATAGAGCCACTGCACTCCAGCCTGGGCAACCAGGTGAGACCTTGTCTCTAAAA
TGAATAAATAAAT

Sequence ID 1355

20 TGGTCTTTCACCCAGCCAGGGAGAAGGTTCTTCGCTCAGTATGAAGAAAAGCAACC
CAAACCTCTCAATCTGATTTGTTTTTGTATTATGTCGATGCCCTGTAGTTTGAAAGT
GAAGTAAAGATTTAGAATTCACCTAAGTCCAAAGGAAAACACGTGGTTTTTTAAAGC
CATTAGGTAAAAAAAGTTCTCAATAAAGGCATTACAATTTTTTTAGGTTTAGAAAGA
TGGACTTTTCTGATAAATCTTGGCAGACATCTAAAAAACCATATTTTTTCACA
25 AGAAAATGCAAGTTACTTTTTTTGGAAATAATACTCACTGATTATGGATAAAATGG
AATATTTTCAGATACTATATTGGCTGTTTCAAAATAGTACTATTCTTTAACTTGT
AATTTTTTGCTAAGTTATTTGTCTTTGTTGTATCTATAAATATGTAAAAAATATTTA
AATAGATGTACCTGTTTTGCTTTCACACTTAATAAAAAATTTTTTTTTGT

30 Sequence ID 1359

CGGGATCCCTAGTATAACACATTCAAGTGTCCCTTTCAGTCTTACTACTTTGACC
GCGATGATGTGGCTTTGAAGAACTTTGCCAAATACTTTCTTCACCAATCTCATGAG
GAGAGGGAACATGCTGAGAACTGATGAAGCTGCAGAACCAACGAGGTGGCCGAAT
CTTCCTTCAGGATATCAAGAAACCAGACTGTGATGACTGGGAGAGCGGGCTGAATG
35 CAATGGAGTGTGCATTACATTTGGAAAAAATGTGAATCAGTCACTACTGGAAGTGC
ACAAACTGGCCACTGACAAAAATGACCCCCATGTGAGTATTGGAACCCCAGGAAAT
AAATGGAGGAAATCATTTGCCTTAGGGATTGGGAAAGCTGCCCACTAACTGTCTTC

- 264 -

CCCATTTGTTTTGCAGTTGTGTGACTTCATTGAGACACATTACCTGAATGAGCAGGT
GAAAGCCATCAAAGAATTGGGTGACCACGTGACCAACTTGCGCAAGATGGGAGCGC
CCGAATCTGGCTTGGCGGAATATCTCTTTGACAAGCACACCCTGGGAGACAGTGAT
AATGAAAGCTAAGCCTCGGGCTAATTTCCCCATAGCCGTGGGGTGACTTCCCTGGT
5 CACCAAGGCAGTGCATGCATGTTGGGGTTTCCTTTACCTTTTCTATAAGTTGTACC
AAAACATCCACTTAAGTTCTTTGATTTGTACCATTCCCTTCAAATAAAGAAATTTGG
TACC

Sequence ID 1360

10

TGCGCAGACCAGACTTCGCTCGTACTCGTGCGCCTCGCTTCGCTTTTCTCCGCAA
CCATGTCTGACAAACCCGATATGGCTGAGATCGAGAAATTCGATAAGTCGAAACTG
AAGAAGACAGAGACGCAAGAGAAAAATCCACTGCCTTCCAAAGAAACGATTGAACA
GGAGAAGCAAGCAGGCGAATCGTAATGAGGCGTGCGCCGCCAATATGCACTGTACA
15 TTCCACAAGCATTGCCTTCTTATTTTACTTCTTTTAGCTGTTTAACTTTGTAAGAT
GCAAAGAGGTTGGATCAAGTTTAAATGACTGTGCTGCCCCCTTTCACATCAAAGAAC
TACTGACAACGAAGGCCGCGCCTGCCTTTCCCATCTGTCTATCTATCTGGCTGGCA
GGGAAGGAAAGAACTTGCATGTTGGTGAAGGAAGAAGTGGGGTGGAAGAAGTGGGG
TGGGACGACAGTGAAAT

20

Sequence ID 1361

TATAAATACACTCCGGGATGATTTACCCCCGGAGGTCAGCTAGTAAAATACATGAG
TAGAATTCTTTAAAGTATGTGATAATTGCTCATCACTATCCAAGTGTGACATAAAT
CATAAAAAGAATTGACAAAATCAGGGTCGCAAAGAGAATTGAAAAAATCTGTAC
25 AACC AAAATTTAAATTGACCTCTGTCTAGAGTATGAGAGCCACACTGAACAGAAA
AACCAGATAAATCTTTTATAAAATATTCATTTGCAGCCCCATTAACGTTGCTTGTC
ACCCACCTCCCATGTCTTGGACAACTGAATGTATAGTAACATCATCCCAGGC
CAGGCGCGGTGGCTCATGCCTGTAATCCCAGCACTTTGTGAGGCTAAGGCAGGCAG
ATCAGGAGGTCAGGAGTTCAGGACCAGCCTGGCCAAAAGGTGAACTCCGTCTCT
30 ACTAACAATACAAAATTAGCTGGGTGCGGTAGTAGGCGCCTGTAATCCCAGCTAC
TCGGGAGGCTGAGGCAGGAGAATTGCTCAAACCCGGAAGGTGGAGGTTGCAGTGAG
CTGAGATCGTGCCACTGCACTCCAGCCTGGGTGACAGAGCAAGACTCTGTCTCGGG
GAGGGGGGTGGCGGAGATAAAGAAATAACATCATCTTATACTGTCAAGCTCAAGGT
GTCTGCAGCCTTATCTTCAGGGGAAGTTGTGTCTTTCTCAGGGGAAGATACAGATTT
35 CAATTTAGAGCAAGACAGAGAGAAGTTACATTCAGAGAGGAAAATGCAGTAGTCTA
ACTG

- 265 -

Sequence ID 1364

5 GCGGCCGCGCTCTTTTCAATTTTAAAAAGAAGTTTGTTCATTTTCAGTAATTT
CTGCTTTGATCTTCCTTATGTCCTCCTATTGAGTTGATCAGCTTTCTTTATTCTTG
CCTTTTCTCCTCTGTGTGCCCTTTCTATTAACGTATTTACCCTTAGGCTGGGCACA
ATGGCTGATGCCTGTAATCCCTGCACTTTGGGAGGCCGAGGCAGGTGGATCACCTA
AGGTCAGGAGTTCAAGACCAGCCTGGCCAACATGGTGAAACCTGGTCTCTACTAAA
AACACAAAAATTAGCCAGGCATGGTGGTGTGCACCTGTAATCCCAGCTACTCAGGA
GGCTGAGGCAGGAGAATTGCTTGAACCTGGGAGGCCGAGATTGTGCCAAAGCACTC
CAGCCTGGGCAACAAAATGAGACTTTGTGTC

10

Sequence ID 1365

15 CACCAGGCTGTCTTCAGATACTTCATACAGAAATGAGCCTCCCTGTGGGGTCCCTCT
TCCCTCCTTCAGCCTGTCCATCAACACAGCATTGCGGGATCCTTACCATGGCATCC
AGCCCTGGAGATGCTTCAGGAAAGTTGCAGGTCCATGCTGCAGGACAGGCTCAGAT
CAGCAGAGACGCATCTCACATCGGGCTGTGAAATTCAAGTTGAGCTGCAATTGGCA
ATGAGAA

Sequence ID 1366

20 GTTATTCACTGAGACCGTGCCCCGGTTATGAGGTTGTACCAGAAAGCAAGTATTCA
CTATGCACACTATTACCGCTCACCCCTAGCATTGAAGCCAGCCTGTAGCCTGAAAG
CCTTTGCTTTGAGGGCAGGTCTTTCCCCAAAATGCAGACACGAAGGTGCAAAGTGA
AGCTGCCAGTCTTGCAAAAGATGTAACCTGTCACGAAGGCCACGAGTGGCAGGGAG
AGCTGTCCCACATTTGCGGAAGTGGCTATGTGAGGACGGGGGAGGCGGGTCCCTTA
GAGATGAGACAATCATAAGGGGAGATATCAGAGAAAATCGTAAGGGGAGCAGATGG
25 TTGTCAAGAGAATAGGCTGACCATCGAAGGACTGGCAGAAGCTTTCAGAAAACCAC
TGGACGGCTGGGCACAGTGGCTTAGGCCTGTAATCCCAGCACTTTGGGAGGCTGAC
GCAGGTGAATCACTTGAGGTCAGGAGTTCAGACCAGCCTGGCCAACATGGTGAAA
CCCCATCTCTACAGAAAATATAAAAATTAGCCAGGCGTGGTGGCACAAGCCTAGAA
TCCCAGCTACTTGGGAGGCTGAGGCAGGCGAATGGCTTGAACCCAGGAGTCAGAGG
30 CTGCAGTGAGTCGAGATTGTTCCACTGCACTCCAGCCTGGGTGACAGTGCAAGACT
CCTTCCAAAAA

Sequence ID 1367

35 TTCGTGAGTGATGGCGTCCCGGGTTGCTTGCCGGTGCTGGCCGCCCGGGAGAGC
CCGGGGCAGAGCAGAGGTGCTCATCAGCACTGTAGGCCCGGAAGATTGTGTGGTCC
CGTTCCCTGACCCGGCCTAAGGTCCCTGTCTTGCACTGGATAGCGGCAACTACCTC
TTCTCCACTAGTGCAATCTGCCGATATTTTTTTTTTGTATCTGGCTGGGAGCAAGA

- 266 -

5 TGACCTCACTAACCAGTGGCTGGAATGGGAAGCGACAGAGCTGCAGCCAGCTTTGT
CTGCTGCCCTGTACTATTTAGTGGTCCAAGGCAAGAAGGGGGAAGATGTTCTTGGT
TCAGTGCGGAGAGCCCTGACTCACATTGACCACAGCTTGAGTCGTCAGAACTGTCC
TTTCCTGGCTGGGGAGACAGAATCTCTAGCCGACATTGTTTTGTGGGGAGCCCTAT
10 ACCCATTA CTGCAAGATCCCGCCTACCTCCCTGAGGAGCTGAGTGCCCTGCACAGC
TGGTTCAGACACTGAGTACCCAGGAACCATGTCAGCGAGCTGCAGAGACTGTACT
GAAACAGCAAGGTGTCCTGGCTCTCCGGCCTTACCTCCAAAAGCAGCCCCAGCCCCA
GCCCCGCTGAGGGAAGGGCTGTCACCAATGAGCCTGAGGAGGAGGAGCTGGCTACC
CTATCTGAGGAGGAGATTGCTATGGCTGTTACTGCTTGGGAGAANGGCCTAGAAAG
15 TTTTGCCCCCGCTGCGGCCCCAGCANAATCCAGTGTTGCCCTGTGGCTGGAGAAAGG
AATGTGCTCATCAGCAGTGCCCTCCNTTACGTCAACAATGTCCCCCACCTTGGGAA
CATCATTTGGTTGTGTGCTCAGTGCCCGATGTCTT

Sequence ID 1368

15 CAGTGAGCCAAGATCACACCACTGCACTCCAGCCTGGACAACAGAACGAGACTCCA
TATCAAAAAAATTAAATTAAAATATAATAAATTTCTTGCCGGGCGCAGTGGCTCAC
ACCTGTAATCCCAGCACTTTGGGAGGCCGAGGTGGGCGGATCACGAAGTCAGGAGA
TTGAGACCATCCTGGCTAATACAGTGAAACCCCGTCTCTACTATAAATACAAAAAA
TTAGCTGGGCATGGTGGCGGGCGTCTGTAGTCCCAGCTACTCAGGAGTCTGAGGCA
20 GGAGAATGGTGTGAACCCGGGAGGCGGAGCTTGCAAGTGAGCCGAGATCGTGCCACT
GCAATCCAGCCTGGGCAGCAGAACGAGACTCCATCTCAAATAAATAAATAAATAAA
ATGAATTTTCAGCTAGAAGAGCCTTATTCCATTTTCCTTTTTATTAAACATCTGGCA
TAAGTTGGTAAGTATGTGAAGTTTATCATATATTCTTATGCGAATTATTATTTTCG
CCTTTTTTTTTTATAATTCTGTCTGGGATTTGAATAGTAGAGTTTGAATTCAGGAAG
25 GACACCTGTGATAGGACAATAAAAT

Sequence ID 1369

CTGATTGCAAAAACATTACAACCTCAGTACTGCGGCTTTCATTCAAATAGGTTATAT
GTATAAACTGAGGTTCAACAATATTGTATTTGAGATGGGAAAGTTAAAGAAATGCA
30 ATAATGTAAATAATACTTAAGAAAATAAGATCTCAGGAAACTGTGTATACTCTGTA
CTTTTATGCAACTTTATCAGATCATTTTCAGTATATGCATCAAGGATATAGTGTATA
TGACATGAACTTTGAGTGCAAAAACGTACTATGTACCTTTTGTTTATTTTGCTGT
CAACATCTAAATAAAGGTTTTTTTG

35 Sequence ID 1370

CGAAAGGACTACAGAGCCCCGAATTAATACCAATAGAAGGGCAATGCTTTTAGATT
AAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTTTAAAAGTTGTAGGTGATTAA

- 267 -

AATAATTTGAAGGCGATCTTTTAAAAAGAGATTAAACCGAAGGTGATTAAAAGACC
TTGAAATCCATGACGCAGGGAGAATTGCGTCATTTAAAGCCTAGTTAACGCATTTA
CTAAACGCAGACGAAAATGGAAAGATTAATTGGGAGTGGTAGGATGAAACAATTTG
GAGAAGATAGAAGTTTGAAGTGGAAAACCTGGAAGACAGAAGTACGGGAAGGCGAAG
5 AAAAGAATAGATAAGATAGGGAAATTAGAAGATAAAAACATACTTTTAGAAGAAAA
AAGATAAATTTAAACCTGAAAAGTAGGAAG

Sequence ID 1371

GTCCAGNAGAAAGTTCAGTGACTTGTCCAGAGCTGCAGGTCTTAAGAGGCTGAAAT
10 CTCGCCTCTGCCCTCGAGGCTGCGGTTCCACTGACCCATACTACTTGCCTTCAGGAA
AGAGAAATGGTGTAGGAAGGCTGTGGATGAAGACGCTTACATTCATGAAGGATTTG
GATAGGCGAACATGAGCTTTTCCACCAAATTTTCAGAATTTTAAGAAATGCCTTAAA
TTATTTCTTAAAAATCAATTTGGGGCAGACGAGAAGTTCTGATAATAGTTTTTAGG
GAACATGATAAAATTTCTGACCTTAGAAGTGGTATACCAGTTTGAGAAGAAGAACAA
15 GCTATAAACGGTGTAGATAACATTCACGGCTATTTAAGAAAGAGTTACTAAGGGAA
ACCAGAATGACTTAAGAGTGTTACTCTTCTTTTTCTGAGAGAACAATAGCATCATC
TCAGAAAGCCTTTCATGCCATTAATAGGTAAGAATCTGGGCTTCTTGGACCATGGG
TTAGACTTTCTTACAAAACCATATATGCATTTCTTAGCAAAATTTATGCTATTAC
ATTTCTTATCTCAACAAAGACTGGTAAATTCAGTACTTATTCCTCAATTTTCTTA
20 CCCTTAAAATGGGGATATTCTGCCTCTCCAAGGAATGCTGGGAACAAGCAAGTCCT
CATGTTAGGGGTCTTTGAGTTTTTCATGGAAGTTTAGGTTATTTATATGATGACATA
GTTGTCAACTTACTTTCAGGATGGACTTTTCTTTTGTGAGTTTGTGACCTAAATAC
AATAGTTGTTATGCATGTCCAGTTTATGGAAGTACCACTGCAATANCAG

25 Sequence ID 1372

CAGTGCAGCCAAGTATCACACCACTGCACTCCAGTCCTGGACAACAGAAACGANTA
CTCCATATCAAAAAAATTAAATTAAANGATAATAAATTTCTTGCCGGGCGCAGTGG
CTCACACCTGTAATCCCAGCACTTTGGGAGGCCGAGGTGGGCGGATCACGAAGTCA
GGAGATTGAGACCATCCTGGCTAATACAGTGAAATCCCCGTCTCTACTATAAATAC
30 AAAAAATTAGCTGGGCATGGTGGCGGGCGTCTGTAGTCCCAGCTACTCAGGAGTCT
GAGGCAGGAGAATGGTGTGAACCCGGGAGGCGGAGCTTGCAGTGAGCCGAGATCGT
GCCACTGCAATCCAGCCTGGGCGAGCAGAACGAGACTCCATCTCAAATAAATAAATA
AATAAAATGAATTTTCAGCTAGAAGAGCCTTATTCCATTTTCCTTTTTATTAAACAT
CTGGCATAAGTTGGTAAGTATGTGAAGTTTATCATATATTCTTATGCGAATTATTA
35 TTTTCGCCTTTTTTTTTTATAATTCTGTCTGGGATTTGAATAGTAGAGTTTGAATTC
AGGAAGGACACCTGTGATAGGACAATAAAATCTA

- 268 -

Sequence ID 1374

GAAAGCACATATGATATACATGTGTGTCATATGTATTATTTTGTGTTGCCATCTGAG
TCTTCAAAATTTGTTACAGAATACCTGCATATTAATATTTCAAGGTATGGATTAAT

5 Sequence ID 1378

CTGAGTATTAATAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Sequence ID 1380

CCAAACCCAACTGGTCCAGTAGGATACTCACCTTACAGGGGGCGTCTCAAGAGTCT
10 CACAGTTCCCTTGCGTCTTAAGAGACTCACTGTTGGACCAGGCGTGGTGACTCACG
CCTGTAAAACCACTTTGGGAGGCCGAGGCGGGCGGATCAGTTGAGGTCAAGAG
TTCAAGACCAGCCTGACCAAGGTGCTGAAACCCCGTCTCTACTAAAAATACAAAA
TTAGCCAGGCATGGTGGTGTGCGCCTGTAATCCCAGCTACTCCAGAGGCTGAGGCA
GGAGAATCTCTTGAACCCAGGAGGTGGAGGTTGCAGTGAGTCGAGATCATGCCACT
15 GCACTCCAGCCTGGGTGACAGAGCGAGACTCCGTCTTAGAAAAAAAAAAAAAAAAAA
AAAAGAACCTCACAGTTCAGCAGGGTTCTAGCATGAGACAATGAGGACAAGGGTAG
GTGAGCAGGTGGAAAGAGTGAGAACAGGTCAATTGTGATGGAGAAAATAATAAAGA
CAGAAAAGGCAGAAGACTGCCTGGCAGAAGACCTGTCCCAGCAGATACAAAAATAC
AGACAACAGGAGCCAGCATAGACCCTTGACCTGTGTAAGTCTTTCTCAGGCCTTCT
20 TTTAAGTAGAAACATGCCTTTGAAAAAAGTTTAAATAAACAGGAAAATCATAAAT
CCCTATTTACATAAATAATATATCCTGGTCTTATTCTTAAAACCATTGATTTTTCA
CGGCTCATTAAANAAGCTGGGCGAGGTGGCTCACGCCCGTCATCCTAGCACTTTGG
GAGGCCGAGGCGGGCANATCACAAGGTGAGGAGTTGGGAGACCAGCCTGACCAACA
CGGTGAAACCCAGTCTCTACTAAAAATACAAAAATTANCTGGGGGTGGTGGTGTGT
25 GCCTGTAATCCAAGCTACTCGGGAGGCTGAGGCAGGA

Sequence ID 1382

CTTACTACCTCCAACATGAAACAAGCAGCCCCGCACTTCTCGAAGGTCTGAGTTAC
TTGGAATCGTTTTACCACATGATGGACAGAAGGAATATTTAGATATCTCTGAAAA
30 CCTCAAGCGTTACCTTCTTCAGTATTTTAAGCCAGTGATTGACAGGCAAAGCTGGA
GTGACAAGGGCTCAGTCTGGGACAGGATGCTCCGCTCGGCTCTCTTGAAGCTGGCC
TGTGACCTGAACCATGCTCCTTGATCCAGAAAGCTGCTGAACTCTTCTCCAGTG
GATGGAATCCAGTGGAATAATAATATACCAACAGATGTTTTAAAGATTGTGTATT
CTGTGGGTGCTCAGACAACAGCAGGATGGAATTACCTTTTAGAGCAATATGAACTG
35 TCAATGTCAAGTGCTGAACAAAACAAAATTCTGTATGCTTTGTCAACGAGCAAGCA
TCAGGAAAAGTTACTGAAGTTAATTGAACTAGGAATGGAAGGAAAGGTTATCAAGA
CACAGAACTTGGCAGCTCTCCTTCATGCGATTGCCAGACGTCCAAAGGGGCAGCAA

- 269 -

CTAGCATGGGATTTTGTAAAGAGAAAATTGGACCCATCTTCTGAAAAAATTGACTT
GGGCTCATATGACATAAGGATGATCATCTCTGGCACAACAGCTCACTTTTCTTCCA
AGGATAAGTTGCAAGAGGTGAACTATTTTTTTGAATCTCTTGAGGCTCAAGGATCA
CATCTGGATATTTTTCAAACGTCTGGAACGATAACCAAAAATATAAAATGGCT
5 GGAGAAGAATCTTCCGACTCTGAGGACTTGGCTAATGGTTAATACTTAAATGGTCA
ATAGAAAAAGTAGGCTGGGCGCGGTGGCTCACGCCTGTAATCCCAGCACTTTGGGA

Sequence ID 1387

AAAATT
10 TTTCAGT
GTTAAAGTAGGTTTGTTCGACGCGCCACGAATTTCCCGGGGACCAA

Sequence ID 1389

TTTTTTTTTTTTTTTGGGAGTCAGTTTTCTTTTCTTTTCTTTTCTTTTTTTTTTTTTT
15 GNTTTTCGGAAACGGAGTCTCGCTTCTCGCCCACTCTGGAGTGGNGCAGTGGGGN
GGTCTCAGCTCACCACAGCCTCCACCTCCTGGGCCCCAAGCGATCCTNTCACCTCAG
CCTCCTGCGTAGCTGGGACTACAGGCGTGCACCACCATTTCCAGGTAATTTTTGT
TTTTTGTANANACAGGGTTTCACTGTTGTTGCCCAGGCTGGTCTCGAACTCCTGC
TTCAGTCTGCCANAATGCTGGATTCTAGGCGTGAGCCACCGNGCCTGGCCCCAAAG
20 TTACTTTTCTTACAGAAGCAAAGCTTTAATGCATTTTACTGAATGCTTATAGCTTT
GTAGATACTGAAAAGAGTATGAGCGTCACATACAGACACATNTAACAGCACTGCCT
CCAACCAGCCCCTACCCACTGGTCAGGNGAGTAANAATCAAAATTTCTTTTCTGNGA
GTGGAACGGAAATTTTCATCTCTCCTCCTCAGGCAAGTAGTTAANAGGCTGGNGGGA
GTCATGGCCCCATTTTGTTCAAAATACAAGCTCCACAGGAACAAAAGGCTGAACTG
25 CTCACCTCCCAACTGATGAACCTCGTCTTTGTTCCATGTCAAAGGGGCCTTTGTGT
TACTGCAGCAGAACTCCAGCTATCAAACCATCAGGCACCAAAAGTAAAACCTCCTT
TCTCTAAAAGACCTCTCTTTACCTGAGCCTTTCAATGCATCTTTGCCCCCANATA
ATCCTGGATGAGATAATCCCCAGAGGAANACCAGCGCTTGCCTAGTGAAATTATAC
TATGAGACAAGGGTAAAAGACCTCAAANACCGGGTTGGCAGGTAAGGGAGTAGGGN
30

Sequence ID 1390

TCNGTGGCACCCGTTTCCGGCACCTTCAGACTCTGAAGAGCCACCTGCGAATCCAC
ACAGGAGAGAAACCTTACCATGTACGTAAGCCTCTTGAGGCCGCTCTCTGACCTGC
GGGGATGTGGAGGGCAGGGAAGGAGGTGGAGCGCAGGGAAGGAGGTGGAGCAGGGA
35 GGCAGTGGAACGTGTTTGCTCCCATCTCAAGCACACAGTGGGGCAACCACTACGCTA
ATGGTTGGAAGACCTAGATCTGGGCCCAATGGCCAGACACCCTGCTTGACCTTGGC
CCAAGCATTAGGGGACTCATCTTTAAATGAGGGTATGGGACTAGATGATCTGGGC

- 270 -

CTTAGGAGAGGAGT

Sequence ID 1391

CGGCTNCTACCCTGCGGAGATCACACTGACCTGGCAGTGGGATGGGGAGGACCAAA
5 CTCAGGACACCGAGCTTGTGGAGACCAGGCCAGCAGGAGATGGAACCTTCCAGAAG
TGGGCAGCTGTGGTGGTGCCTTCTGGAGAAGAGCAGAGATACACGTGCCATGTTCA
GCACGAGGGGCTGCCGGAGCCCCTCACCTTGAGATGGAAGCCGTCTTCCCAGCCCA
CCATCCCCATCGTGGGCATCGTTGCTGGCCTGGCTGTCCTGGCTGTCCTAGCTGTC
CTAGGAGCTATGGTGGCTGTTGTGATGTGTAGGAGGAAGAGCTCAGGTGGAAAAGG
10 AGGGAGCTGCTCTCAGGCTGCGTCCAGCAACAGTGCCAGGGCTCTGATGAGTCTC
TCATCGCTTGTAAGCCTGAGACAGCTGCCTGTGTGGGACTGAGATGCAGGATTTT
TTCACACCTCTCCTTTGTGACTTCAAGAGCCTCTGGCATCTCTTTCTGCAAAGGCA
TCTGAATGTGTCTGCGTTCCTGTTAGCATAATGTGAGGAGGTGGAGAGACAGCCCA
CCCCCGTGTCCACCGTGACCCCTGTCCCCACACTGACCTGTGTTCCCTCCCGATC
15 ATCTTTCCTGTTCCAGAGAAGTGGGCTGGATGTCTCCATCTCTGTCTCAACTTCAT
GGTGGCTGAGCTGCAACTTCTTACTTCCCTAATGAAGTTAAGAACCTGAATATAA
ATTTGTTTTCTCAAATATTTGCTATGAAGGGTTGATGGATTAATTAATAAGTCAA
TTCCTGGAAGTTGAGAGAGCAAATAAAGACCTGAGAACCTTCCANAATCCG

20 Sequence ID 1392

TGAAACAAAATGAATTTNTATGGGTAAGAGAGGGTAATATTTTAGAGTTGTGTTAC
AAACTACAAATTTTTATTAAATTAATAAATCAGAATACTAAATCCATGTGTTTTT
TTCTTTCTTAAAAAATATCTTTTGGCTGGGCACGGTAGCTCATGGCTGTAATCCCA
GCACTTTGGGAGGCTGAGGTGGGTGGATCGCCTGATGTCAGGAGTTCAAGACCAGC
25 CTGGTCAACATGTTGAAACCCCATCTCTACTAAAAATATAAAAATTAGCCGGTGTG
GTGGTGGGCGCCTGTAATCCCAGCTACTCAGGAGGCTAAGGCAGGAGAATTGCGTG
AACCAGGAGTTCAGTGATGTAGCGGGGAGCTGAGATTGTGCCACTACACTCCAGC
CTGGATGACAGAGTGAGACTCCATCTCAAAAAAAAAAAAAAAAAAAAA

30 Sequence ID 1394

GCATAATGTGAGGAGGTGGAGAGACAGCCCAACCCCGTGTCCACCGTGACCCCTGT
TCCCATGCTGACTTGTGTTTCTTCCCCAGTCATCTTTCCTGTTCCAGAGAGGTGGG
GCTGGATGTCTCCATCTCTGTCTCAACTTTATGTGCACTGAGCTGCAACTTCTTAC
TTCCCTACTGAAAATAAGAATCTGAATATAAATTTGTTTTCTCAAATATTTGCTAT
35 GAGAGGTTGATGGATTAATTAAATAAGTCAATTCTGGAATTTGAGAGAGCAAATA
AAGACCTGAGAACCTTCCAGAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AA

- 271 -

Sequence ID 1395

CTTACCATGTCAGTGCACAGAAATGCTGCTTGGGATGTAGGAAAAATAAATCCAC
AAAAGCTACCAAGTTTGAAGGGGACCATGAGTCTTCAGGCTGGAGCTTCCAAACCA
GATGAAAACCCCAACAATTAACCTGCAGTTTAAAGATCCAGCAGCTGGCCATTTCTGG
5 ACTCAAGGTGAATCGTCTGGATATGTATGGAGAAAAGTACAAACCCTTTAAGGGCA
TAAATAACATGACCAAAGCTGGGAAGTTCCAAGTTCGAACCTGAAGGGAGCATTG
CTGAGGGAATAGTCTTGCACATTTTTTTCATTTCTTACTTGTCTAAAAGTAAAAAAA
AATATCAGCCTGTCTCCTAGGTGAGTCCCCCTCCTGGACCCACCCGCTCCCTTTTTT
CCTTAGCCTTCAGTGCCATGGAACATAATCAAGGGAGGAAAAGGTCACCAGGGAGAA
10 CTGGACAGAACTGAAACACAGCAACACCAGTTCTCAAGGACAAGGTGTGTGATGGG
GGTAGGAAGCTTGGTGCTTATGTAACCATTTTAAACGTGGTTTCTATAGGAAAGAC
CAACATTTGTTTAGCTTGCTTGGCTTTAATTATCTAAAGCCAATGAAAGACTTCTT
TGTTGATTTTTTAAGATAGAAAGATT

15 Sequence ID 1396

CAAACACTATGTTATTTTATGAANAAGACTTGAACATCTATGGATTTTGGTATTTG
CAAGGGGTGAATGGGGTATTTGCAAGCAGTGAATGAGGAGGCCTGGAACCAATCTT
CTGCTGATATTGAGGCACAACCTGAAAAAGGTATATTACTTAAATCTCTTATTGTAT
TGTAAACTGTATAAGTAATGAAATTAAAAGGCAGAAATTGTCAGACTGAATAAAAT
20 GAAAAGACCAAACAATATGCTGCTTACAAGAAACACAATTCAAATATAAGGACACA
ATTAGTTTAAAGGAAAAGAACTGGAAAAGATATACCATGATAACACAAGTCAGAAG
AAAGCTGCTGTGGATATATTAATATGAGATGTAGATTTTCAGAGCAGTGAATATTGC
CAGGCATAAAGAAAGTTATTACATAATAATTAAGGTATCAGTTCATCAAGAAGATG
TAATAACCCTAAGTATTTTATACAACTAATATCAGAGCTTCAAATACATGAAGCAA
25 AAACCAGTGGAATTGATAGGAGAAACACACAATTACACAATTATAGTCAGAAATTT
CAACATATCTTTCTCAATGGAGAAAACAACTAGACAGGAAATCATTAAGGATATAG
ATGATTTAAATTATATGATCAACTACCTGGACGTAATTGGCATTATGGAACACTG
CACCACCAACAGCAGAGTACATATTATTTTCAAGTACACAGAAAACAGTTACCAAT
ATAGACCATTTTCTGGGTCATAAAACACATCTCAATAAATGTAAAACAATTAATGT
30 TATATAAAGTATGTGCTCTGACCNCAAAGGAATTAGAGATCAATAAAAGAACATCT
TTGAAAAATCTCACNTATTTAAAACTAATAACTCACTTCTAAATAACTCCTGTNT
CAAGAGAATNAAANGG

Sequence ID 1397

35 CCCAGCCTCACTGCGCCCCGTGAGGCCAGGCAGCTGCCCTCAGGGTCTGCCAAGGT
GGGGGTCAAGGGCCATGGGGGCAGGTAGCTCTGCCTGCAAAGCCCACAAGCATGTC

- 272 -

AGATCACCTGGGCTGCAGACAGACAAACACCTGAGCTGTTCTGAATACCTTCAGGT
TCCTGGCCTCGCTGAGCAAGTGCAGAAATTTTACCTTCAAGGATCAGGGTTTTTC
TGTTTGTTTGTTTTTTAACACACACATATGTGAACAAAGAGTATGCGTTTGTACTG
GCAGAAGAAGCGTCTGGTAAGACAACCAGCAAGTTAACAATGGTCACCTCCAGAAA
5 TGGGCTGGGTAAACCAAAGAATTTTTTTGTTTTTTGTTTTTTTTGAGTCAGGGTCTA
GCTCTGTCACCCAGGCTGGAACGCACTGGTGTGATCACGGCTCACTGCAGCCTTGA
CCTCCCTGGCTCAAGCAATCCTCCCAGCTCAGCCTCCTGAGTCGTTGGGACTACAG
GCACGTGCCACCACGCCTGACACATTTTTTAAATTTTTGTAGAGACAGTGTTTCAC
CATGTTGCCCAGGCAGGTCTCAAACCTCCTGGGCTCAAGTGGTCCTCCAGCTTCAGC
10 CTCCCAAAGTGCTAGGATTATAGGTGTGAGCCACAGTGCCAGCCCCGTAAGTGGAG
AATTTCTGTTGAATGAACCAAAGCAACTGCCAACCTCTCCATGCACCATGTGTTT
CAGAGGAGAAAGCACAGTGAAGAATGCAGTGTGTTCTGAGGTCCTGTCACCCCTGA
GGCTGTGTGTGTCCTTTGCCAAATTAAAGAGTCTTACTGAATGCGGTGCATCCAGG
AGACAGGCCNAGGTTTGGACTGGTAAAAAAA

15

Sequence ID 1399

CAGACACCTGGNAGAACGGGAAGGAGACGCTGCAGCGCGCGGACCCCCCAAAGACA
CATGTGACCCACCACCCCATCTNTGACCATGAGGCCACCCTGAGGTGCTGGGCCCT
GGGCTTCTACCCTGCGGAGATCACACTGACCTGGCAGCGGGATGGCGAGGACCAAA
20 CTCAGGACACCGAGCTTGTGGAGACCAGACCAGCAGGAGACAGAACCTTCCAGAAG
TGGGCAGCTGTGGTGGTGCCTTCTGGAGAAGAGCAGAGATACACATGCCATGTACA
GCATGAGGGGCTGCCGAAGCCCCTCACCTGAGATGGGAGCCATCTTCCAGTCCA
CCGTCCCCCATCGTGGGCATTGTTGCTGGCCTGGCTGTCCTAGCAGTTGTGGTCATC
GGAGCTGTGGTTCGCTGCTGTGATGTGTAGGAGGAAGAGTTCAGGTGGAAAAGGAGG
25 GAGCTACTCTCAGGCTGCGTCCAGCGACAGTGCCAGGGCTCTGATGTGTCTCTCA
CAGCTTGAAAAGCCTGAGACAGCTGTNTTGTGAGGGACTGAGATGCAGGATTTCTT
CACGCCTCCCCCTTTGTGACTTCAAGAGCCTCTGGCATCTCTTCTGCAAAGGCACC
TGAATGTGTCTGCGTCCTTGTAGCATAATGTGAGGAGGTGGAGAGACAGCCCACC
CTTGTGTCAACTGTGACCCCTGTTCCCATGCTGACCTGTGTTTCTCCCCAGTCA
30 TCTTTTTTGTTCNCAATAGGTGGGGCCTGGATGTCTCCATCTCTGTNTCA

Sequence ID 1440

TTATAAGGTACTTTTAAGGTATTTTAGTTGTCTTAGTCTATATTTCTGTACTCACC
35 TTTCTTTATCCACTCATCAGTTGATGGGCATGTAGGTTGGTTCCATATCTTTGCAA
TTCTGAATTGTGCTGTGATCAGGTGTCTTTTTAGTATAATGATTTACTCTCCTTTG
GGTAGATACCCAGTAGTGGGATTGCTGGATCGAATGGTTTTTATAATTTTCTATTT

- 273 -

- TACCACAGTTTCTCTCTGCATTTTTCCTCTTTGACCACTAACCATGTGAAATTCTC
ATATTGACCTTTATAATGATCATGAACCTCTTAGTATCATTGGGAAGGCCACATTTG
CCACTTATGATTGTAAACCTTATCCTCCATTTTTCTGTATTGTTGGTGCAAAAA
GCACCTATTATACCAGGACTTTAAAAATCAGTCTGATAAGTCTTTGATAAGTCTAA
5 TAATAATAACTGATAAGTCCATTGAATTTGCTTCTGATTACTTTTTCTTTAGTAGC
TAAACATGTATGTACTCCTATGATTACAATGAACACTCCTCTCCATTTAAATTAAT
TATTTACATTGATGAAATAGCAAAATGTTAATGACTAAATACTGTCTTGGTTTTTT
CGTTCCAGGTCAAGTCAATATTAACCTTCTTATAATTTCTTTTTTTCTTT
- 10 Sequence ID 1447
GCAAGGACTAACCCCTATACCTTCTGCATAATGAATTAAGTAGAAATAACTTTGCA
AGGAGAGCCAAAGCTAAGACCCCCGAAACCAGACGAGCTACCTAAGAACAGCTAAA
AGAGCACACCCGTCTATGTAGCAAAATAGTGGGAAGATTTATAGGTAGAGGCGACA
AACCTACCGAGCCTGGTGATAGCTGGTTGTCCAAGATAGAATCTTAGTTCAACTTT
15 AAATTTGCCACAGAACCCCTCTAAATCCCCTTGTAATTTAACTGTTAGTCCAAAG
AGGAACAGCTCTTTGGACACTAGGAAAAAACCTTGTAAGAGAGAGTAAAAAATTTAA
CACCCATAGTAGGCCTAAAAGCAGCCACCAATTAAGAAAGCGTTCAAGCTCAACAC
CCACTACCTAAAAAATCCCAAACATATAACTGAACTCCTCACACCCAATTGGACCA
ATCTATCACCCCTATAGAAGAACTAATGTTAGTATAAGTAACATGAAAACATTTCTCC
20 TCCGCATAAGCCTGCGTCAGATTAAACACTGAACTGACAATTAACAGCCCAATAT
CTACAATCAACCAACAAGTCATTATTACCTCACTGTCAACCCAACACAGGCATGC
TCATAAGGAAAGGT
- Sequence ID 1448
25 GGCCACCGGGTGCAAGGTCAGGGCTGGGGTGAGGCTGGGAAGCCCAGGGCTTGGC
CCACTGTGGCCGCCTTGTGTGGTCACTGCTTTCCTGGGCCTGCTGTGAGCTCCCTC
TAGGACCCCAGGCCTGTCTGGTGGGTCACTGTGACCACCACCTTGCACAGCACCTG
GCGCGTGGCAGGTGCTCAAACATTACTTGTTCGGAATGAACTTCATCTTGCTCTT
GGCTTTTTTGACTAATGCTGTGGAACATCTGACTAATTAGTGACTCTTTGGGGCCCC
30 CAGTTTCCCAGCTATAAAGTGGTAATATTAAGATAATAATTCTGGCCGGGCGCGGTG
GCTCACGCCTGTAATCCCAGCAGCACTTTGGGAGGCCGAGGTGGGCAGATCACGAG
GTCAGAAGATCGAGACCATCCTGGCTAACACGGTGAAACCCCATCTCTACTAAAAA
TACAAAAAATTANCCGGGCGTGGTGGCGGGCGCCTGTAGTCCCAGCTACTCANGAG
GCTGANGCAGGAGAATGGTGTGAACCCGGGAGGCAGAGGTTGCAGTGAACCAAGAT
35 CGNNCCACTGCACTCCAGCCTGGGCAACAGAGCGAGACTCCATCTTAAAAAA

Sequence ID 1449

- 274 -

AATCAGGGCCGCGAGTGTGTTCTGCGCCTGCCAGAGCTGACTCCTGATTTAACCGC
TGGCGTAACCGCGGGTTGCACGCATGCGTGCTGAAAAGCCTTTACCCCTCACGTGG
TTTCTTTTTTTAAACAGTCATCAAGCGAGGCTCGCGCGCAGGCCCCGCGTTGGAAAA
TGGCGGGGAAGCTGAAACCTCTGAATGTGGAGGCGCCAGAAGCTGCTGAGGAGGCT
5 GAAGGTAGTGAGGGCAAGTGGGCTGCACTCCTTTCTCTCCAACCAGGGCAGAAAGG
AGGGAGGATTTCGTCCCATTACAATAATGAAATAATGATATTCTAATTTTTTTTAAAT
AAAATGTTAAGCCTTTTGTTATTGAA

Sequence ID 1450

10 GGAAANCATGAGGCTTCGGGAGCCGCTCCTGAGCGGCAGCGCCGCGATGCCAGGCG
CGTCCCTACAGCGGGCCTGCCGCCTGCTCGTGGCCGTCTGCGCTCTGCACCTTGGC
GTCACCCCTCGTTTACTACCTGGCTGGCCGCGACCTGAGCCGCCTGCCCCAACTGGT
CGGAGTCTCCACACCGCTGCAGGGCGGCTCGAACAGTGCCGCCGCCATCGGGCAGT
CCTCCGGGGAGCTCCGGACCGGAGGGGCCCCGGCCGCCCTCCTNTAGGCGCCTCC
15 TCCCAGCCGCGCCCCGGGTGGCGACTCCAGCCCAGTCGTGGATTCTGGCCCTGGCCC
CGCTAGCAACTTGACCTCGGTCCCAGTGCCCCACACCACCGCACTGTGCTGCCCC
CCTGCCCTGAGGAGTCCCCGCTGCTTGTTAAGGACTCGGGTCCGCGCCAGTCCGAG
GATTGGGACCCCCCGGATTTCCCCGACAGGGTCCCCCANACATTCCTCAGGCTG
GCTCTTCTACGACAGCCAGCCTCCCTCTTCTGGATCAGAGTTTTAAATCCCANACA
20 GAGGCTTGGGACTGGATGGGAGAGAAGGTTTGCGAGGTGGGTCCCTGGGGAGTCCT
GTTGGAGGCGTGGGGCCGGGACCGCACAGGGAAGTCCCGAGGCCCTCTAGCCCCA
AAACCANAGAAGGCCTTGGAGACTTCCCTGCTGTGGCCCGAGGCTNAGGAAGTTTT
GGAGTTTTGGGTCTGCTTANGGCTTCNAGCAGCCTTGCACTGAGAACTTTGGTAGG
GACCTCGAGTAATCCACTCCNTTTTNGGGACTGACGTGAGGCTCCCGGTGGGGAAA
25 GANACTGACCTNTC

Sequence ID 1453

CCGACCTGTCTCGCTCCGTGGCCTTAGCTGTGCTCGCGCTACTCTCTCTTTCTGGC
CTGGAGGCTATCCAGCGTACTCCAAAGATTACAGGTTTACTCACGTCATCCAGCAGA
30 GAATGGAAAGTCAAATTTCTGAATTGCTATGTGTCTGGGTTTCATCCATCCGACA
TTGAAGTTGACTTACTGAAGAATGGAGAGAGAATTGAAAAAGTGGAGCATTTCAGAC
TTGTCTTTCAGCAAGGACTGGTCTTTCTATCTCTTGTACTACACTGAATTCACCCC
CACTGAAAAAGATGAGTATGCCTGCCGTGTGAACCATGTGACTTTGTCACAGCCCA
AGATAGTTAAGTGGGATCGAGACATGTAAGCAGCATCATGGAGGTTTGAAGATGCC
35 GCATTTGGATTGGATGAATTCCAAATTCTGCTTGCTTTTAAATATTGATATG
CTTATACACTTACACTTTATGCACAAAATGTAGGGTTATAATAATGTTAACATGGA
CATGATCTTCTTTATAATTCTACTTTGAGTGCTGTCTCCATGTTTGATGTATCTGA

- 275 -

GCAGGTTGCTCCACAGGTAGCTCTAGGAGGGCTGGCAACTTAGAGGTGGGGAGCAG
AGAATTCTCTTATCCAACATCAACATCTTGGTCAGATTTGAACTCTTCAATCTCTT
GCACTCAAAGCTTGTTAAGATAGTTAAGCGTGCATAAGTTAACTTCCAATTTACAT
ACTCTGCTTAGAATTTGGGGGAAAATTTAGAAATATAATTGACAGGATTATTGGAA
5 ATTTGTTATAATGAATGAAACATTTTTGTGCATATAAGATTCATATTTACTTCTTAT
ACA

Sequence ID 1454

TAAATAGGGAATCCTTTCCCATTTGCTTGTTTTCTCAGGTTTGTCAAAGATCAGA
10 TAGTTGTAGATATGCGACGTTATTTCTGAGGGCTCTGTTCTGTTCCATTGATCTAT
ATCTCTGTCACATGCACACGTATGTTTGTGTGGCACTATTCACAGTGGCAAAGAC
TTGGAACCAACCCAAATGTCCAACAATGATAGACCGGGTTAAGAAAATGCGGCACA
TATACACCATGGAATACTATGTAGCCATAAAAAATGATGAGTTCGTGTCCTTTGTA
GGGACATGGATGAAATTGGAAATCATCATTTCTCAGTAACTATCGCAGGAACAAAA
15 AACCAAACACTGCATATTCTCACTCATAGGTGGGAATTGAACAGTGGGAACACATG
GACACAGGAAGGGGAACATCACACTCTGAGGACTGTTGTGGGGTGGGGGGAGGGAG
GAGGGATAGCATTGGGAGATATACCTAGTGCTGGATGACGAGTTAGTGGGTGCAGC
GCACCAGCATGTCACATGTATACATATGTAACCTAACCTGCACATTGTGCACATGTA
CCCTAAAACTTAAGGTAT

20

Sequence ID 1456

CCGCAACAAACACGGGAGTGCAGATATCGCTGCGATGGGCTGATTTCTTTTATTTG
GGTATATACCCAGCAGTGGGATTGCTGGATTGTATGGTAGCTCTATTAGTTTTTTG
25 AGGAACCTCCAACTGTTCTNCATAGTGGTTGTACTCATTTACATTCCCAGTGTGA
ACCCTGAAAATTTGAGGCAGGTCTCAGTTAAATTAGAAAGTTGATTTTGCCAAGTT
GGGGACACGCACTCGTGACACAGCCTCAGGAGGAACTGATGACATGTGCCCAGGTG
GTCAGAGCACAGCTTGGTTTTATACATTTTAGGGAAACCTGAGCCATCAATCAACA
TACGTAAAATGGGCCGGGCACAGCAGCTCAAGCTGTAATCCCAGCACTCTGGGAGG
30 CCGAGGCGGGTGGATCACTTGAGGTGAGGAGTTTCGAGACCAGCCTGGCCAACATGG
TGAAACCCCGTCTCTATTAAAAATACAAAGCTTAGCTGGATGTGGTGGCGCATGCC
TGTAGTCCCAGCTGCTCTAGGAGGCTGAGGCATGAGAATTGCTTGAACCTGGGAGG
CAGAGGCTGCAGTGAGCCGAGATCGAGCCACTATACTCCAGCCTGGTCAACAGAGT
GAGACCCTGTCT

35

Sequence ID 1460

CCACAACCTGTGTTCACTAGCAACCTCAAACAGACACCATGGTGCACCTGACTCCTG

- 276 -

AGGAGAAGTCTGCCGTTACTGCCCTGTGGGGCAAGGTGAACGTGGATGAAGTTGGT
GGTGAGGCCCTGGGCAGGCTGCTGGTGGTCTACCCTTGGACCCAGAGGTTCTTTGA
GTCCTTTGGGGATCTGTCCACTCCTGATGCTGTTATGGGCAACCCTAAGGTGAAGG
CTCATGGCAAGAAAGTGCTCGGTGCCTTTAGTGATGGCCTGGCTCACCTGGACAAC
5 CTCAAGGGCACCTTTGCCCACTGAGTGAGCTGCACTGTGACAAGCTGCACGTGGA
TCCTGAGAACTTCAGGCTCCTGGGCAACGTGCTGGTCTGTGTGCTGGCCCATCACT
TTGGCAAAGAATTCACCCCAACAGTGCAGGCTGCCTATCAGAAAGTGGTGGCTGGT
GTGGCTAATGCCCTGGCCCAAGTATCACTAAGCTCGCTTTCTTGCTGTCCAATT
TCTATTAAAGGTTCTTTGTTCCCTAAGTCCAACTACTAACTGGGGGATATTATG
10 AAGGGCCTTGAGCATCTGGATTCTGCCTAATAAAAAACATTTATTTTCATTG

Sequence ID 1490

ATGGGCATCTCTCGGGACAACTGGCACAAGCGCCGCAAAACCGGGGGCAAGAGAAA
GCCCTACCACAAGAAGCGGAAGTATGAGTTGGGGCGCCAGCTGCCAACACCAAGA
15 TTGGCCCCCGCCGCATCCACACAGTCCGTGTGCGGGAGGTAACAAGAAATACCGT
GCCCTGAGGTTGGACGTGGGGAAATTTCTCCTGGGGCTCANAGTGTTGTACTCGTAA
AACAAGGATCATCGATGTTGTCTACAATGCATCTAATAACGAGCTGGTTCGTACCA
AGACCCTGGTGAAGAATTGCATCGTGCTCATCGACAGCACACCGTACCGACAGTGG
TACGAGTCCCACTATGCGCTGCCCTGGGCCGCAAGAAGGGAGCCAAGCTGACTCC
20 TGAGGAAGAAGAGATTTTAAACAAAAACGATCTAAAAAATTCAGAAGAAATATG
ATGAAAGGAAAAAGAATGCCAAAATCAGCAGTCTCCTGGAGGAGCAGTTCAGCAG
GGCAAGCTTCTTGCGTGCATCGCTTCAAGGCCGGGACAGTGTGGCCGAGCAGATGG
CTATGTGCTAGAGGGCAAAGAGTTGGAGTTCTATCTTAGGAAAATCAAGGCCCGCA
AAGGCAAATAAATCCTTGTTTTGTCTTCACCCATGTAATAAAGGTGTTTATTGTTT
25 TTGTT

Sequence ID 1491

CTTNCACATACTGATTGATGTCTCATGTCTCTCTAAAATGTGTAAAACCAAGCTGT
GCCCCAACCACCTTGGGNACATGTGGNGAGGACCTCCTGAGGCTGTGTCTATGGGCA
30 CACCTTAACCCTGGGAAAATAAACTTTCTAACTGACTTGAGAGCTGTCTCAGATA
TTCTGAGCTTACAGTTATTGTGAAATCATTTTAAATTATAAATTAAGTGGAGATTTA
CTTAAATCATGTGTAGTAAGTAGCCTGTGATATAGTCCTAGATACATACATTATCA
TCTTATGTATCTTCCCTCCCTCTTCCAGGTTCTGATAAAAACAGATGAAATCTGAA
AGACCATGACAGTAGTATTTTGAAAATGACAGTATTTGAAATTAATAAATTTGTAA
35 AGTGTTCTGTTCTATCACTGCCAAAGGATAAGTTACAAATTGGTTCTTGGAACGTA
ATATGTACTATGTGCTTGCTATTTAATAATTTACCAGTCTTAGTCTTTTTTATTCA
GACTAATTTTACCTTTTTTTAACCCTATGACTCTTTAGTTATAGTAGTACAAAAAAG

- 277 -

TAGTTTTAGTTATAGTTTTAGTTGTAGTACAAAAAGCATTTTCTGTAAGCTTAAT
TTCTTTCCCTTCCCGCTTTCCAGTCAGATGACTTTAGTGATTTGGAGTTGTGTG
CTTTATAAGTGCATTCCTCAGAGGACTTAATATTACTAAGATTTTAGCAACNCTGA
AATATGTT

5

Sequence ID 1492

TGTNCCTGTAGTCCTGTGTGGGAGGATTGCCTGAGCCTAGGAGCTCAAAGTTGCAG
TGAGCCCAGATCGNGNCATTGCAGTCCAGCCTGGGTGACAGAGTGAGACCCCATGT
CAAAAAAAAAAAAAACAAAAACAGGGGCCTGCCTCANCCAGCAGGTGAGGTCTGCC
10 ACTGAGAGCACTTCTAGCAGCAGGAACAGCCTCCACCCCCACACTGCAATCAAGTT
TTTTGGGTGAGCCTTAGGAGCTAANAAAGGGCCTAGTTTGNCTAAATAGCAGGAGT
TATATCCAGGGATCTTCAGGCCCAGGAATGCTAATGAGTAGGCATTCCATGGGCCC
TGGGAATGGCTTTGTGTGCCANAAATGATGGCCACAAAGGCCTTGCTGCCTTTTTT
CAAAATGGCTGCATCCAGCTGAGTGCTCTCTGCCAAAGGGGANAANAATAAGTC
15 TCCAGTGCATTTAGATTGGTCTCTCATCATCTCTCTCCTTTTTGTTTTTATTAGTC
TCCTTAACCAAACTGCCAAGAAAGGCTTGGAATTGAAACAAAACCTGATANAANA
GGTAAGAGGTTGTTCTTTT

Sequence ID 1493

TGTNTCAAAAAAAAAAAAAAGAACGGNAATGTACTGGAGATGTATTTGATAACCAA
20 GGNTTTAGGTAAATTTTCACCAGTATTAGTTNTATTTGCAAACCTGAAAAATGTTGT
AGGCTTAATATAAAATAACCACATTAGTGAACATTATATCTCTTAGAAGAAAGGCC
ATATTTTGCTCCTGCTTCTGTAAAAATATTATTTGTTTGAAGGGGAAATAATGGTA
GTGTGACCTTTCACCTAATTCCTACTCCCTTAATGTGAGAGAGACAAAATGAGCTG
25 AAGAAGGAAAATTCTGGAGTTACACTCCACAACCTTGAACATACTGACGGACATCT
CTGTTTTGACAACGATTTCTCCATGCCACCCATGCTNTAATGCCTTGTGGATCACG
GACAACCCTCTTTGCACAAGCTACAGCATCAGCGATGTTATCTTGACAGCAAAGCAC
TGCAGGATAAATGACAGGCATTAAGTCTCCTGGGGTTTTGCCATCATTACACCAG
TAGCGGCTATTGATCTGAAATATCCATAATCAGTGCTTCTGTCTCCAGCATTGTA
30 GTTTGTAGCTCGTGTGTTGTAACCACTCTCCCATTTGGCCAAACACATCCAGTTTG
CTAGGCTGATTCCTCTGTAGCCATCCATTCCCAATCTTTTCAGAGTTCTGGCCAA
TCACACCTTTCAAAGACCTTGCCCTGGACCGTAACAGAAAGGAGGACAAGCCCCAG
AACAAATGAGAGCCTTCATGTTGAC

35

Sequence ID 1494

TTGGTACCCGGGAAATTCCTTTGCCGCGTCGACGGCCGGTGAGGCAGATCACCTGAG
CCCAGGAGTTCAGGACCAGCCTGGGCAGCATACCGGGATTCCATCTNNACTAAAAA

- 278 -

5 CAGTAGGCTGGGTGTGGTGGCTCATGTCTGTAAGCTCAGGACTTTGGAAGGCCAAG
ATGGGAGGATCACTTGAGCCTGGGAGTTTGACACCAGCTTGAGCATCGTAGCCAGG
CCCTGACTCTACAAAAAAGTGAAATAATTAGCCGAGTGTGGTGGTTACACCTGTA
ATCCCAGCTGCTCAGGAGGCTGAGGTAGGAGAATCATTTGAACCCGGGAGGTGGAG
GTTGCAGTTAGCCGAGATCACGCCATTGCACTCCGGCCTGGGCGATAAAGCGAGAC
TCTGTCTCAAAAAAAAAAAAAA

Sequence ID 1495

10 ATTCGGGCCGAGATGTCTCGCTCCGTGGCCTTAGCTGTGCTCGCGCTACTCTCTCT
TTCTGGCCTGGAGGCTATCCAGCGTACTCCAAAGATTCAGGTTTACTCACGTCATC
CAGCAGAGAATGGAAGTCAAATTTCTGAATTGCTATGTGTCTGGGTTTCATCCA
TCCGACATTGAAGTTGACTTACTGAAGAATGGAGAGAGAATTGAAAAAGTGGAGCA
TTCAGACTTGTCTTTCAGCAAGGACTGGTCTTCTATCTCTTGTACTACACTGAAT
TCACCCCCACTGAAAAAGATGAGTATGCCTGCCGTGTGAACCATGTGACTTTGTCA
15 CAGCCCAAGATAGTTAAGTGGGATCGAGACATGTAAGCAGCATCATGGAGGTTTGA
AGATGCCGCATTTGGATTGGATGAATTCCAAATCTGCTTGCTTGTCTTTTAATAT
TGATATGCTTATACACTTACACTTTATGCACAAAATGTAGGGTTATAATAATGTTA
ACATGGACATGATCTTCTTTATAATTCTACTTTGAGTGCTGTCTCCATGTTTGATG
TATCTGAGCAGGTTGCTCCACAGGTAGCTCTAGGAGGGCTGGCACCTTAGAGGTGG
20 GGAGCAGAGAATTCTCTTATCCAACATCAACATCTTGGTCAGATTTGAACTCTT

Sequence ID G6

25 GGATTTTGGTCCGCACGCTCCTGCTCCTGACTCACCGCTGTTGCTCTCGCCGAG
GAACAAGTCGGTCAGGAAGCCCGCGCAACAGCCATGGCTTTTAAGGATACCGGA
AAAACACCCGTGGAGCCGGAGGTGGCAATTCACCGAATTCGAATCACCTAACAAG
CCGCAACGTAAAATCCTTGGAAGGTGTGTGCTGACTTGATAAGAGGCGCAAAAG
AAAAGAATCTCAAAGTGAAAGGACCAGTTCGAATGCCTACCAAGACTTTGAGAATC
ACTACAAGAAAACTCCTTGTGGTGAAGGTTCTAAGACGTGGGATCGTTTCCAGAT
GAGAATTCACAAGCGACTCATTGACTTGACACAGTCCTTCTGAGATTGTTAAGCAGA
30 TTACTTCCATCAGTATTGAGCCAGGAGTTGAGGTGGAAGTCACCATTGCAGATGCT
TAAGTCAACTATTTTAATAAATTGATGACCAGTTGTTAAAA

Sequence ID - 61 nt: 362

CTTATTGAAAATTTTACTAATTTCTTACTTTTTAGGTTTTAGGAGAATACTTTTGGA
TAATTGACTAGCCTCACATTATATTGATAGAGGFTCTTGAAAACTTTAATGCCAAT
TCATGTATCTTATGACTAAAATAGATAATCCATTTAGAAATTTAAGTCATTCTTGC
GTGCTTGATATGTGTTCAGCACTATCCAAGTTGCTAGGGGATAACAATGGTGAAGTG
AAAATATCAGCTAGGTGCCGGTGGCTCACACCTGTTATCCCAACAGTTTGGGAGG
CCAGGGTGGGAGGATCACTCAAGCACANGCGTTTCACACCAGCCTGGACAACAT
ACAAGACCCCATCTTTACCAAAAGTTAAG

Sequence ID - 490 nt: 382

TTTTCTTAGAACTTTATTTTTTCTGGCCAGGCGCAGTGGCTCACACCTGTAATCCC
AGCACTTTGGGAGGCCAAGGCAGGTCGATCACCTGAGGTGAGGAGCTCAAGACC
AGCCTGGCCAACATGGTGAAACCCTGTCTCTACTAAAAATACAAAAATTAGCTGG
GCGTGGTGGCGCATGCCTGTAATCCCANCTACTCAGGAGGCTGAGGCAGGAGAA
TTGTTTGAACCCGGGAGGCGGAGGTTGCANTGAGCCGAGATTGCGCCACTGCACT
CCAGCCTGGGCAACAGAGCGAAACTCCATCTCAAAAAAAAAAAAAAAAAACAAC
CTTTATTTTTTCTGATTTTAAAAGTAATAACTAGTTTGTAGAAACATTAAAAGT

Sequence ID - 892 nt: 559

TCTTTCGGAAGCGCGCCTTGTGTTGGTACCCGGGAATTCGCGGCCGCGTCGACGC
GGTCGTAAGGGGCTGAGGATTTTTGGTCCGCACGCTCCTGCTCCTGACTCACCGCT
GTTGCTCTCGCCGAGGAACAAGTCGGTCAGGAAGCCCGCGCGCAACAGCCATG
GCTTTTAAGGATACCGGAAAAACACCCGTGGAGCCGGAGGTGGCAATTCACCGA
ATTCGAATCACCCTAACAAGCCGCAACGTAAAAATCCTTGGAAGAGGTGTGTGCTG
ACTTGATAAGAGGCGCAAAAGAAAAGAATCTCAAAGTGAAAGGACCAGTTTCGAA
TGCCTACCAAGACTTTGAGAATCACTACAAGAAAACTCCTTGTGGTGAAGGTTT
TAAGACGTGGGATCGTTTCCAGATGAGAATTCACAAGCGACTCATTGACTTGCAC
AGTCCTTCTGAGATTGTTAAGCAGATTACTTCCATCAGTATTGAGCCAGGAGTTG
AGGTGGAAGTCACCATTGCAGATGCTTAAGTCAACTATTTTAATAAATTGATGAC
CAGTTGTTT

Sequence ID - 77 nt: 464

GCGGCTGCTGTTGGTTGGGGGCCGTCCCGCTCCTAAGGCAGGAAGATGGTGGCCG
CAAAGAAGACGAAAAAGTCGCTGGAGTCGATCAACTCTAGGCTCCAACCTCGTTAT
GAAAAGTGGGAAGTACGTCTGGGGTACAAGCAGACTCTGAAGATGATCAGACA
AGGCAAAGCGAAATTGGTCATTCTCGCTAACAACCTGCCCAGCTTTGAGGAAATCT
GAAATAGAGTACTATGCTATGTTGGCTAAAACCTGGTGTCCATCACTACAGTGGCA
ATAATATTGAACTGGGCACAGCATGCGGAAAAATACTACAGAGTGTGCACACTGG
CTATCATTGATCCAGGTGACTCTGACATCATTAGAAGCATGCCAGAACAGACTGG
TGAAAAGTAAACCTTTTCACCTACAAAATTTACCTOCAAACCTTAAACCTGCAA
AATTTTCCTTTAATAAAATTGCTTG

- 280 -

Claims:

1. A set of oligonucleotide probes, wherein said set comprises at least 10 oligonucleotides selected from:
5 an oligonucleotide as described in Table 1 or derived from a sequence described in Table 1, or an oligonucleotide with a complementary sequence, or a functionally equivalent oligonucleotide.
- 10 2. A set of oligonucleotide probes as claimed in claim 1 wherein said oligonucleotide probes are selected from: an oligonucleotide as described in Table 2 or derived from a sequence described in Table 2, or an
15 oligonucleotide with a complementary sequence, or a functionally equivalent oligonucleotide.
- 20 3. A set of oligonucleotide probes as claimed in claim 1 wherein said oligonucleotide probes are selected from: an oligonucleotide as described in Table 4 or derived from a sequence described in Table 4, or an
oligonucleotide with a complementary sequence, or a functionally equivalent oligonucleotide.
- 25 4. A set of oligonucleotide probes as claimed in any one of claims 1 to 3, wherein each probe in said set binds to a different transcript.
- 30 5. A set as claimed in any one of claims 1 to 4 consisting of from 10 to 500 oligonucleotide probes.
- 35 6. An oligonucleotide probe wherein said probe is selected from the oligonucleotides listed in Table 1, or derived from a sequence described in Table 1, or a complementary sequence thereof.
7. A set of oligonucleotide probes as claimed in any one of claims 1 to 5, or an oligonucleotide probe as

- 281 -

claimed in claim 6, wherein each of said oligonucleotide probes is from 15 to 200 bases in length.

- 5 8. A set of oligonucleotide probes as claimed in any one of claims 1 to 5 or 7 or an oligonucleotide probe as claimed in claim 6 or 7, wherein the transcript to which said probe binds is derived from a gene which is constitutively moderately or highly expressed.
- 10 9. A set of oligonucleotide probes as claimed in any one of claims 1 to 5, 7 or 8 or an oligonucleotide probe as claimed in any one of claims 6 to 8, wherein said probes are immobilized on one or more solid supports.
- 15 10. A set of oligonucleotide probes or an oligonucleotide probe as claimed in claim 9, wherein said solid support is a sheet, filter, membrane, plate or biochip.
- 20 11. A polypeptide encoded by the mRNA sequence to which an oligonucleotide as defined in claim 6 binds.
- 25 12. An antibody to a polypeptide as defined in claim 11.
- 30 13. A kit comprising a set of oligonucleotide probes immobilized on one or more solid supports as defined in claim 9 or 10.
- 35 14. A kit as claimed in claim 13 wherein said probes are immobilized on a single solid support and each unique probe is attached to different region of said solid support.
15. A kit as claimed in claim 13 or 14 further comprising standardizing materials.

- 282 -

16. The use of a set of probes as described in any one of claims 1 to 5 or 7 to 10 or a kit as described in any one of claims 13 to 15 to determine the gene expression pattern of a cell which pattern reflects the level of gene expression of genes to which said oligonucleotide probes bind, comprising at least the steps of:

- a) isolating mRNA from said cell, which may optionally be reverse transcribed to cDNA;
- b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotides or a kit as defined in any one of claims 1 to 5, 7 to 10 or 13 to 15; and
- c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce said pattern.

17. A method of preparing a standard gene transcript pattern characteristic of a disease or condition or stage thereof in an organism comprising at least the steps of:

- a) isolating mRNA from the cells of a sample of one or more organisms having the disease or condition or stage thereof, which may optionally be reverse transcribed to cDNA;

b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotides or a kit as defined in any one of claims 1 to 5, 7 to 10 or 13 to 15 specific for said disease or condition or stage thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and

c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce a characteristic pattern reflecting the level of gene expression of genes to which said oligonucleotides bind, in the sample with the disease, condition or stage thereof.

18. A method of preparing a test gene transcript pattern comprising at least the steps of:

- a) isolating mRNA from the cells of a sample of

- 283 -

said test organism, which may optionally be reverse transcribed to cDNA;

5 b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotides or a kit as defined in any one of claims 1 to 5, 7 to 10 or 13 to 15 specific for a disease or condition or stage thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and

10 c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce said pattern reflecting the level of gene expression of genes to which said oligonucleotides bind, in said test sample.

15 19. A method of diagnosing or identifying or monitoring a disease or condition or stage thereof in an organism, comprising the steps of:

- 20 a) isolating mRNA from the cells of a sample of said organism, which may optionally be reverse transcribed to cDNA;
- 25 b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotides or a kit as defined in any one of claims 1 to 5, 7 to 10 or 13 to 15 specific for said disease or condition thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation;
- 30 c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce a characteristic pattern reflecting the level of gene expression of genes to which said oligonucleotides bind in said sample; and
- 35 d) comparing said pattern to a standard diagnostic pattern prepared as described in claim 17 using a sample from an organism corresponding to the organism and sample under investigation to determine the degree of correlation indicative of the presence of said

- 284 -

disease or condition or a stage thereof in the organism under investigation.

20. A method as claimed in any one of claims 17 to 19
5 wherein said mRNA or cDNA is amplified prior to step b).

21. A method as claimed in any one of claims 17 to 20
10 wherein the oligonucleotides and/or the mRNA or cDNA are labelled.

22. A method as claimed in any one of claims 17 to 21
15 wherein said probes are as defined in claim 3 and said disease is Alzheimer's disease.

23. A method as claimed in any one of claims 17 to 21
20 wherein said probes are as defined in claim 2 and said disease is breast cancer.

24. A method as defined in any one of claims 17 to 23,
25 wherein said set of oligonucleotides as defined in any one of claims 1 to 5, 7 to 10 or 13 to 15 are replaced with a set of oligonucleotides which are randomly selected, preferably from a cDNA library.

25. A method of preparing a standard gene transcript
30 pattern characteristic of a disease or condition or stage thereof in an organism comprising at least the steps of:

a) releasing target polypeptides from a sample of
35 one or more organisms having the disease or condition or stage thereof;

b) contacting said target polypeptides with one or more binding partners, wherein each binding partner is specific to a marker polypeptide (or a fragment thereof)
40 encoded by the gene to which an oligonucleotide of Table 1 (or derived from a sequence described in Table 1) binds, to allow binding of said binding partners to said

- 285 -

target polypeptides, wherein said marker polypeptides are specific for said disease or condition thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and

- 5 c) assessing the target polypeptide binding to said binding partners to produce a characteristic pattern reflecting the level of gene expression of genes which express said marker polypeptides, in the sample with the disease, condition or stage thereof.

10

26. A method of preparing a test gene transcript pattern comprising at least the steps of:

- a) releasing target polypeptides from a sample of said test organism;
- 15 b) contacting said target polypeptides with one or more binding partners, wherein each binding partner is specific to a marker polypeptide (or a fragment thereof) encoded by the gene to which an oligonucleotide of Table 1 (or derived from a sequence described in Table 1)
- 20 binds, to allow binding of said binding partners to said target polypeptides, wherein said marker polypeptides are specific for said disease or condition thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and
- 25 c) assessing the target polypeptide binding to said binding partners to produce a characteristic pattern reflecting the level of gene expression of genes which express said marker polypeptides; in said test sample.

30

27. A method of diagnosing or identifying or monitoring a disease or condition or stage thereof in an organism comprising the steps of:

- a) releasing target polypeptides from a sample of said organism;
- 35 b) contacting said target polypeptides with one or more binding partners, wherein each binding partner is specific to a marker polypeptide (or a fragment thereof)

- 286 -

5 encoded by the gene to which an oligonucleotide of Table 1 (or derived from a sequence described in Table 1) binds, to allow binding of said binding partners to said target polypeptides, wherein said marker polypeptides are specific for said disease or condition thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and

10 c) assessing the target polypeptide binding to said binding partners to produce a characteristic pattern reflecting the level of gene expression of genes which express said marker polypeptides in said sample; and

15 d) comparing said pattern to a standard diagnostic pattern prepared as described in claim 25 using a sample from an organism corresponding to the organism and sample under investigation to determine the degree of correlation indicative of the presence of said disease or condition or a stage thereof in the organism under investigation.

20 28. A method as claimed in any one of claims 17 to 27 wherein said pattern is expressed as an array of numbers relating to the expression level associated with each probe.

25 29. A method as claimed in any one of claims 17 to 28 wherein said organism is a eukaryotic organism, preferably a mammal.

30 30. A method as claimed in claim 29 wherein said organism is a human.

35 31. A method as claimed in any one of claims 17 to 30 wherein the data making up said pattern is mathematically projected onto a classification model.

32. A method as claimed in any one of claims 17 to 31 wherein said disease is cancer or a degenerative brain

- 287 -

disorder.

33. A method as claimed in any one of claims 17 to 32 wherein said sample is tissue, body fluid or body waste.

5

34. A method as claimed in any one of claims 17 to 33 wherein said sample is peripheral blood.

10

35. A method as claimed in any one of claims 17 to 34 wherein the cells in the sample are not disease cells, have not been in contact with such cells and do not originate from the site of the disease or condition.

15

36. A method as claimed in any one of claims 19 to 35 for the diagnosis, identification or monitoring of two or more diseases, conditions or stages thereof in an organism, wherein said pattern produced in step c) is compared to at least two standard diagnostic patterns prepared as described in claim 17 or 25, wherein each standard diagnostic pattern is a pattern generated for a different disease or condition or stage thereof.

20

37. A method of identifying probes useful for diagnosing or identifying or monitoring a disease or condition or stage thereof in an organism, comprising the steps of:

25

- a) immobilizing a set of oligonucleotide probes, preferably as described hereinbefore, on a solid support;
- 30 b) isolating mRNA from a sample of a normal organism (normal sample), which may optionally be reverse transcribed to cDNA;
- c) isolating mRNA from a sample from an organism, corresponding to the sample and organism of step (b), which is known to have said disease or condition or a stage thereof (diseased sample), which may optionally be reverse

35



- 288 -

- transcribed to cDNA;
- 5 d) hybridizing the mRNA or cDNA of steps (b) and
 (c) to said set of immobilized oligonucleotide
 probes of step (a); and
- e) assessing the amount of mRNA or cDNA
 hybridizing to each of said oligonucleotide
 probes to determine the level of gene
 expression of genes to which said
10 oligonucleotide probes bind in said normal and
 diseased samples to generate a gene expression
 data set for each sample;
- f) normalizing and standardizing said data set of
 step (e);
- 15 g) constructing a calibration model for
 classification, preferably using the
 statistical techniques Partial Least Squares
 Discriminant Analysis (PLS-DA) and Linear
 Discriminant Analysis (LDA);
- 20 h) performing JackKnife analysis and identifying
 those oligonucleotide probes which are
 required for classification of said disease
 and normal samples into their respective
 groups.

1 / 12

FIG. 1

Effect of Direct Standardization (DS) on the Alzheimer
Data measured in two different series of experiments



2 / 12

FIG. 1_{CONT'D}





FIG. 2

Projection of normal (including benign) and breast cancer samples onto a
 classification model generated by PLSR-DA using the data of 44 informative genes

PC = Principal Components
 N = normal
 C = breast cancer patients




FIG. 3




FIG. 4

Projection of normal (women with suspected mammogram but no breast cancer) and breast cancer samples onto a classification model generated by PLSR-DA using the data of 719 cDNAs

PC = Principal Components

N = normal

B = breast cancer patients




FIG. 5

Prediction of normal (women with suspected mammogram but no breast cancer) and breast cancer samples based on 3 principal components using the data of 719 cDNAs




FIG. 6

Projection of normal (women with suspected mammogram but no breast cancer) and breast cancer samples onto a classification model generated by PLSR-DA using the data of 111 cDNAs already on file

PC = Principal Components
 N = normal
 B = breast cancer patients




FIG. 7

Prediction of normal (women with suspected mammogram but no breast cancer) and breast cancer samples based on 4 principal components using the data of 111 cDNAs already on file




FIG. 8

Projection of normal (women with suspected mammogram but no breast cancer) and breast cancer samples onto a classification model generated by PLSR-DA using the data of 345 sequenced cDNAs for breast cancer




FIG. 9

Prediction of normal (women with suspected mammogram but no breast cancer) and breast cancer samples based on 3 principal components using the data of 345 sequenced cDNAs




FIG. 10

Projection of non-alzheimer and alzheimer samples onto a classification model generated by PLSR-DA using the data of 520 sequenced cDNAs




FIG. 11

Prediction of non-alzheimer and alzheimer samples based on 4 principal components using the data of sequenced cDNAs

