

Product and Method

The present invention relates to oligonucleotide probes, for use in assessing gene transcript levels in a cell, which may be used in analytical techniques, particularly diagnostic techniques. Conveniently the probes are provided in kit form. Different sets of probes may be used in techniques to prepare gene expression patterns and identify, diagnose or monitor different states, such as diseases, conditions or stages thereof. Also provided are methods of identifying suitable probes and their use in methods of the invention.

The identification of quick and easy methods of sample analysis for, for example, diagnostic applications, remains the goal of many researchers. End users seek methods which are cost effective, produce statistically significant results and which may be implemented routinely without the need for highly skilled individuals.

The analysis of gene expression within cells has been used to provide information on the state of those cells and importantly the state of the individual from which the cells are derived. The relative expression of various genes in a cell has been identified as reflecting a particular state within a body. For example, cancer cells are known to exhibit altered expression of various proteins and the transcripts or the expressed proteins may therefore be used as markers of that disease state.

Thus biopsy tissue may be analysed for the presence of these markers and cells originating from the site of the disease may be identified in other tissues or fluids of the body by the presence of the markers. Furthermore, products of the altered expression may be

- 2 -

released into the blood stream and these products may be analysed. In addition cells which have contacted disease cells may be affected by their direct contact with those cells resulting in altered gene expression and their expression or products of expression may be similarly analysed.

However, there are some limitations with these methods. For example, the use of specific tumour markers for identifying cancer suffers from a variety of defects, such as lack of specificity or sensitivity, association of the marker with disease states besides the specific type of cancer, and difficulty of detection in asymptomatic individuals.

In addition to the analysis of one or two marker transcripts or proteins, more recently, gene expression patterns have been analysed. Most of the work involving large-scale gene expression analysis with implications in disease diagnosis has involved clinical samples originating from diseased tissues or cells. For example, several recent publications, which demonstrate that gene expression data can be used to distinguish between similar cancer types, have used clinical samples from diseased tissues or cells (Alon et al. 1999, PNAS, 96, p6745-6750; Golub et al. 1999, Science, 286, p531-537; Alizadeh et al, 2000, Nature, 403, p503-511; Bittner et al., 2000, Nature, 406, p536-540).

However, these methods have relied on analysis of a sample containing diseased cells or products of those cells or cells which have been contacted by disease cells. Analysis of such samples relies on knowledge of the presence of a disease and its location, which may be difficult in asymptomatic patients. Furthermore, samples can not always be taken from the disease site, e.g. in diseases of the brain.

In a finding of great significance, the present inventors identified the previously untapped potential of all cells within a body to provide information

- 3 -

relating to the state of the organism from which the cells were derived. WO98/49342 describes the analysis of the gene expression of cells distant from the site of disease, e.g. peripheral blood collected distant from a cancer site.

This finding is based on the premise that the different parts of an organism's body exist in dynamic interaction with each other. When a disease affects one part of the body, other parts of the body are also affected. The interaction results from a wide spectrum of biochemical signals that are released from the diseased area, affecting other areas in the body. Although, the nature of the biochemical and physiological changes induced by the released signals can vary in the different body parts, the changes can be measured at the level of gene expression and used for diagnostic purposes.

The physiological state of a cell in an organism is determined by the pattern with which genes are expressed in it. The pattern depends upon the internal and external biological stimuli to which said cell is exposed, and any change either in the extent or in the nature of these stimuli can lead to a change in the pattern with which the different genes are expressed in the cell. There is a growing understanding that by analysing the systemic changes in gene expression patterns in cells in biological samples, it is possible to provide information on the type and nature of the biological stimuli that are acting on them. Thus, for example, by monitoring the expression of a large number of genes in cells in a test sample, it is possible to determine whether their genes are expressed with a pattern characteristic for a particular disease, condition or stage thereof. Measuring changes in gene activities in cells, e.g. from tissue or body fluids is therefore emerging as a powerful tool for disease diagnosis.

- 4 -

Such methods have various advantages. Often, obtaining clinical samples from certain areas in the body that is diseased can be difficult and may involve undesirable invasions in the body, for example biopsy is often used to obtain samples for cancer. In some cases, such as in Alzheimer's disease the diseased brain specimen can only be obtained post-mortem. Furthermore, the tissue specimens which are obtained are often heterogeneous and may contain a mixture of both diseased and non-diseased cells, making the analysis of generated gene expression data both complex and difficult.

It has been suggested that a pool of tumour tissues that appear to be pathogenetically homogeneous with respect to morphological appearances of the tumour may well be highly heterogeneous at the molecular level (Alizadeh, 2000, supra), and in fact might contain tumours representing essentially different diseases (Alizadeh, 2000, supra; Golub, 1999, supra). For the purpose of identifying a disease, condition, or a stage thereof, any method that does not require clinical samples to originate directly from diseased tissues or cells is highly desirable since clinical samples representing a homogeneous mixture of cell types can be obtained from an easily accessible region in the body.

We have now identified a set of probes of surprising utility for identifying one or more diseases. Thus, we now describe probes and sets of probes derived from cells which are not disease cells and which have not contacted disease cells, which correspond to genes which exhibit altered expression in normal versus disease individuals, for use in methods of identifying, diagnosing or monitoring certain conditions, particularly diseases or stages thereof.

Thus the invention provides a set of oligonucleotide probes which correspond to genes in a cell whose expression is affected in a pattern characteristic of a particular disease, condition or

- 5 -

stage thereof, wherein said genes are systemically affected by said disease, condition or stage thereof. Preferably said genes are metabolic or house-keeping genes and preferably are constitutively moderately or highly expressed. Preferably the genes are moderately or highly expressed in the cells of the sample but not in cells from disease cells or in cells having contacted such disease cells.

Such probes, particularly when isolated from cells distant to the site of disease, do not rely on the development of disease to clinically recognizable levels and allow detection of a disease or condition or stage thereof very early after the onset of said disease or condition, even years before other subjective or objective symptoms appear.

As used herein "systemically" affected genes refers to genes whose expression is affected in the body without direct contact with a disease cell or disease site and the cells under investigation are not disease cells.

"Contact" as referred to herein refers to cells coming into close proximity with one another such that the direct effect of one cell on the other may be observed, e.g. an immune response, wherein these responses are not mediated by secondary molecules released from the first cell over a large distance to affect the second cell. Preferably contact refers to physical contact, or contact that is as close as is sterically possible, conveniently, cells which contact one another are found in the same unit volume, for example within 1cm^3 .

A "disease cell" is a cell manifesting phenotypic changes and is present at the disease site at some time during its life-span, e.g. a tumour cell at the tumour site or which has disseminated from the tumour, or a brain cell in the case of brain disorders such as Alzheimer's disease.

- 6 -

"Metabolic" or "house-keeping" genes refer to those genes responsible for expressing products involved in cell division and maintenance, e.g. non-immune function related genes.

5 "Moderately or highly" expressed genes refers to those present in resting cells in a copy number of more than 30-100 copies/cell (assuming an average 3×10^5 mRNA molecules in a cell).

10 Specific probes having the above described properties are provided herein.

Thus in one aspect, the present invention provides a set of oligonucleotide probes, wherein said set comprises at least 10 oligonucleotides selected from:

15 an oligonucleotide as described in Table 1 or derived from a sequence described in Table 1, or an oligonucleotide with a complementary sequence, or a functionally equivalent oligonucleotide.

20 "Table 1" as referred to herein refers to Table 1a and/or Table 1b. Table 1b contains reference to additional clones and sequences as disclosed herein. Similarly Tables 2 and 4 comprise 2 parts, a and b.

25 The invention also provides one or more oligonucleotide probes, wherein each oligonucleotide probe is selected from the oligonucleotides listed in Table 1, or derived from a sequence described in Table 1, or a complementary sequence thereof. The use of such probes in products and methods of the invention, form further aspects of the invention.

30 As referred to herein an "oligonucleotide" is a nucleic acid molecule having at least 6 monomers in the polymeric structure, ie. nucleotides or modified forms thereof. The nucleic acid molecule may be DNA, RNA or PNA (peptide nucleic acid) or hybrids thereof or modified versions thereof, e.g. chemically modified
35 forms, e.g. LNA (Locked Nucleic acid), by methylation or made up of modified or non-natural bases during synthesis, providing they retain their ability to bind

- 7 -

to complementary sequences. Such oligonucleotides are used in accordance with the invention to probe target sequences and are thus referred to herein also as oligonucleotide probes or simply as probes.

5 An "oligonucleotide derived from a sequence described in Table 1" (or any other table) refers to a part of a sequence disclosed in that Table (e.g. Table 1-4), which satisfies the requirements of the oligonucleotide probes as described herein, e.g. in
10 length and function. Preferably said parts have the size described hereinafter.

 Preferably the oligonucleotide probes forming said set are at least 15 bases in length to allow binding of target molecules. Especially preferably said
15 oligonucleotide probes are from 20 to 200 bases in length, e.g. from 30 to 150 bases, preferably 50-100 bases in length.

 As referred to herein the term "complementary sequences" refers to sequences with consecutive
20 complementary bases (ie. T:A, G:C) and which complementary sequences are therefore able to bind to one another through their complementarity.

 Reference to "10 oligonucleotides" refers to 10 different oligonucleotides. Whilst a Table 1
25 oligonucleotide, a Table 1 derived oligonucleotide and their functional equivalent are considered different oligonucleotides, complementary oligonucleotides are not considered different. Preferably however, the at least 10 oligonucleotides are 10 different Table 1
30 oligonucleotides (or Table 1 derived oligonucleotides or their functional equivalents). Thus said 10 different oligonucleotides are preferably able to bind to 10 different transcripts.

 Preferably said oligonucleotides are as described
35 in Table 1 or are derived from a sequence described in Table 1. Especially preferably said oligonucleotides are as described in Table 2 or Table 4 or are derived

- 8 -

from a sequence described in either of those tables. Especially preferably the oligonucleotide (or the oligonucleotide derived therefrom) has a high occurrence as defined in Table 3, especially preferably >40%, e.g. >80 or >90, e.g. 100%.

A "set" as described refers to a collection of unique oligonucleotide probes (ie. having a distinct sequence) and preferably consists of less than 1000 oligonucleotide probes, especially less than 500 probes, e.g. preferably from 10 to 500, e.g. 10 to 100, 200 or 300, especially preferably 20 to 100, e.g. 30 to 100 probes. In some cases less than 10 probes may be used, e.g. from 2 to 9 probes, e.g. 5 to 9 probes.

It will be appreciated that increasing the number of probes will prevent the possibility of poor analysis, e.g. misdiagnosis by comparison to other diseases which could similarly alter the expression of the particular genes in question. Other oligonucleotide probes not described herein may also be present, particularly if they aid the ultimate use of the set of oligonucleotide probes. However, preferably said set consists only of said Table 1 oligonucleotides, Table 1 derived oligonucleotides, complementary sequences or functionally equivalent oligonucleotides, or a sub-set thereof (e.g. of the size as described above), preferably a sub-set for which sequences are provided herein (see Table 1 and its footnote). Especially preferably said set consists only of said Table 1 oligonucleotides, Table 1 derived oligonucleotides, or complementary sequences thereof, or a sub-set thereof.

Multiple copies of each unique oligonucleotide probe, e.g. 10 or more copies, may be present in each set, but constitute only a single probe.

A set of oligonucleotide probes, which may preferably be immobilized on a solid support or have means for such immobilization, comprises the at least 10 oligonucleotide probes selected from those described

- 9 -

hereinbefore. Especially preferably said probes are selected from those having high occurrence as described in Table 3 and as mentioned above. As mentioned above, these 10 probes must be unique and have different sequences. Having said this however, two separate probes may be used which recognize the same gene but reflect different splicing events. However oligonucleotide probes which are complementary to, and bind to distinct genes are preferred.

As described herein a "functionally equivalent" oligonucleotide to those described in Table 1 or derived therefrom refers to an oligonucleotide which is capable of identifying the same gene as an oligonucleotide of Table 1 or derived therefrom, ie. it can bind to the same mRNA molecule (or DNA) transcribed from a gene (target nucleic acid molecule) as the Table 1 oligonucleotide or the Table 1 derived oligonucleotide (or its complementary sequence). Preferably said functionally equivalent oligonucleotide is capable of recognizing, ie. binding to the same splicing product as a Table 1 oligonucleotide or a Table 1 derived oligonucleotide. Preferably said mRNA molecule is the full length mRNA molecule which corresponds to the Table 1 oligonucleotide or the Table 1 derived oligonucleotide.

As referred to herein "capable of binding" or "binding" refers to the ability to hybridize under conditions described hereinafter.

Alternatively expressed, functionally equivalent oligonucleotides (or complementary sequences) have sequence identity or will hybridize, as described hereinafter, to a region of the target molecule to which molecule a Table 1 oligonucleotide or a Table 1 derived oligonucleotide or a complementary oligonucleotide binds. Preferably, functionally equivalent oligonucleotides (or their complementary sequences) hybridize to one of the mRNA sequences which corresponds

- 10 -

to a Table 1 oligonucleotide or a Table 1 derived oligonucleotide under the conditions described hereinafter or has sequence identity to a part of one of the mRNA sequences which corresponds to a Table 1
5 oligonucleotide or a Table 1 derived oligonucleotide. A "part" in this context refers to a stretch of at least 5, e.g. at least 10 or 20 bases, such as from 5 to 100, e.g. 10 to 50 or 15 to 30 bases.

In a particularly preferred aspect, the
10 functionally equivalent oligonucleotide binds to all or a part of the region of a target nucleic acid molecule (mRNA or cDNA) to which the Table 1 oligonucleotide or Table 1 derived oligonucleotide binds. A "target" nucleic acid molecule is the gene transcript or related
15 product e.g. mRNA, or cDNA, or amplified product thereof. Said "region" of said target molecule to which said Table 1 oligonucleotide or Table 1 derived oligonucleotide binds is the stretch over which complementarity exists. At its largest this region is
20 the whole length of the Table 1 oligonucleotide or Table 1 derived oligonucleotide, but may be shorter if the entire Table 1 sequence or Table 1 derived oligonucleotide is not complementary to a region of the target sequence.

25 Preferably said part of said region of said target molecule is a stretch of at least 5, e.g. at least 10 or 20 bases, such as from 5 to 100, e.g. 10 to 50 or 15 to 30 bases. This may for example be achieved by said functionally equivalent oligonucleotide having several
30 identical bases to the bases of the Table 1 oligonucleotide or the Table 1 derived oligonucleotide. These bases may be identical over consecutive stretches, e.g. in a part of the functionally equivalent oligonucleotide, or may be present non-consecutively,
35 but provide sufficient complementarity to allow binding to the target sequence.

Thus in a preferred feature, said functionally

- 11 -

equivalent oligonucleotide hybridizes under conditions of high stringency to a Table 1 oligonucleotide or a Table 1 derived oligonucleotide or the complementary sequence thereof. Alternatively expressed, said
5 functionally equivalent oligonucleotide exhibits high sequence identity to all or part of a Table 1 oligonucleotide. Preferably said functionally equivalent oligonucleotide has at least 70% sequence identity, preferably at least 80%, e.g. at least 90, 95,
10 98 or 99%, to all of a Table 1 oligonucleotide or a part thereof. As used in this context, a "part" refers to a stretch of at least 5, e.g. at least 10 or 20 bases, such as from 5 to 100, e.g. 10 to 50 or 15 to 30 bases, in said Table 1 oligonucleotide. Especially preferably
15 when sequence identity to only a part of said Table 1 oligonucleotide is present, the sequence identity is high, e.g. at least 80% as described above.

Functionally equivalent oligonucleotides which satisfy the above stated functional requirements include
20 those which are derived from the Table 1 oligonucleotides and also those which have been modified by single or multiple nucleotide base (or equivalent) substitution, addition and/or deletion, but which nonetheless retain functional activity, e.g. bind to the
25 same target molecule as the Table 1 oligonucleotide or the Table 1 derived oligonucleotide from which they are further derived or modified. Preferably said modification is of from 1 to 50, e.g. from 10 to 30, preferably from 1 to 5 bases. Especially preferably
30 only minor modifications are present, e.g. variations in less than 10 bases, e.g. less than 5 base changes.

Within the meaning of "addition" equivalents are included oligonucleotides containing additional sequences which are complementary to the consecutive
35 stretch of bases on the target molecule to which the Table 1 oligonucleotide or the Table 1 derived oligonucleotide binds. Alternatively the addition may

- 12 -

comprise a different, unrelated sequence, which may for example confer a further property, e.g. to provide a means for immobilization such as a linker to bind the oligonucleotide probe to a solid support.

5 Particularly preferred are naturally occurring equivalents such as biological variants, e.g. allelic, geographical or allotypic variants, e.g. oligonucleotides which correspond to a genetic variant, for example as present in a different species.

10 Functional equivalents include oligonucleotides with modified bases, e.g. using non-naturally occurring bases. Such derivatives may be prepared during synthesis or by post production modification.

 "Hybridizing" sequences which bind under conditions
15 of low stringency are those which bind under non-stringent conditions (for example, 6x SSC/50% formamide at room temperature) and remain bound when washed under conditions of low stringency (2 X SSC, room temperature, more preferably 2 X SSC, 42°C). Hybridizing under high
20 stringency refers to the above conditions in which washing is performed at 2 X SSC, 65°C (where SSC = 0.15M NaCl, 0.015M sodium citrate, pH 7.2).

 "Sequence identity" as referred to herein refers to the value obtained when assessed using ClustalW
25 (Thompson et al., 1994, Nucl. Acids Res., 22, p4673-4680) with the following parameters:
Pairwise alignment parameters - Method: accurate,
Matrix: IUB, Gap open penalty: 15.00, Gap extension
penalty: 6.66;
30 Multiple alignment parameters - Matrix: IUB, Gap open
penalty: 15.00, % identity for delay: 30, Negative
matrix: no, Gap extension penalty: 6.66, DNA transitions
weighting: 0.5.

 Sequence identity at a particular base is intended
35 to include identical bases which have simply been derivatized.

 The invention also extends to polypeptides encoded

- 13 -

by the mRNA sequence to which a Table 1 oligonucleotide or a Table 1 derived oligonucleotide binds. The invention further extends to antibodies which bind to any of said polypeptides.

5 As described above, conveniently said set of oligonucleotide probes may be immobilized on one or more solid supports. Single or preferably multiple copies of each unique probe are attached to said solid supports, e.g. 10 or more, e.g. at least 100 copies of each unique
10 probe are present.

 One or more unique oligonucleotide probes may be associated with separate solid supports which together form a set of probes immobilized on multiple solid support, e.g. one or more unique probes may be
15 immobilized on multiple beads, membranes, filters, biochips etc. which together form a set of probes, which together form modules of the kit described hereinafter. The solid support of the different modules are conveniently physically associated although the signals
20 associated with each probe (generated as described hereinafter) must be separately determinable.

 Alternatively, the probes may be immobilized on discrete portions of the same solid support, e.g. each unique oligonucleotide probe, e.g. in multiple copies,
25 may be immobilized to a distinct and discrete portion or region of a single filter or membrane, e.g. to generate an array.

 A combination of such techniques may also be used, e.g. several solid supports may be used which each
30 immobilize several unique probes.

 The expression "solid support" shall mean any solid material able to bind oligonucleotides by hydrophobic, ionic or covalent bridges.

 "Immobilization" as used herein refers to
35 reversible or irreversible association of the probes to said solid support by virtue of such binding. If reversible, the probes remain associated with the solid

- 14 -

support for a time sufficient for methods of the invention to be carried out.

Numerous solid supports suitable as immobilizing moieties according to the invention, are well known in the art and widely described in the literature and generally speaking, the solid support may be any of the well-known supports or matrices which are currently widely used or proposed for immobilization, separation etc. in chemical or biochemical procedures. Such materials include, but are not limited to, any synthetic organic polymer such as polystyrene, polyvinylchloride, polyethylene; or nitrocellulose and cellulose acetate; or tosyl activated surfaces; or glass or nylon or any surface carrying a group suited for covalent coupling of nucleic acids. The immobilizing moieties may take the form of particles, sheets, gels, filters, membranes, microfibre strips, tubes or plates, fibres or capillaries, made for example of a polymeric material e.g. agarose, cellulose, alginate, teflon, latex or polystyrene or magnetic beads. Solid supports allowing the presentation of an array, preferably in a single dimension are preferred, e.g. sheets, filters, membranes, plates or biochips.

Attachment of the nucleic acid molecules to the solid support may be performed directly or indirectly. For example if a filter is used, attachment may be performed by UV-induced crosslinking. Alternatively, attachment may be performed indirectly by the use of an attachment moiety carried on the oligonucleotide probes and/or solid support. Thus for example, a pair of affinity binding partners may be used, such as avidin, streptavidin or biotin, DNA or DNA binding protein (e.g. either the lac I repressor protein or the lac operator sequence to which it binds), antibodies (which may be mono- or polyclonal), antibody fragments or the epitopes or haptens of antibodies. In these cases, one partner of the binding pair is attached to (or is inherently

- 15 -

part of) the solid support and the other partner is attached to (or is inherently part of) the nucleic acid molecules.

As used herein an "affinity binding pair" refers to two components which recognize and bind to one another specifically (ie. in preference to binding to other molecules). Such binding pairs when bound together form a complex.

Attachment of appropriate functional groups to the solid support may be performed by methods well known in the art, which include for example, attachment through hydroxyl, carboxyl, aldehyde or amino groups which may be provided by treating the solid support to provide suitable surface coatings. Solid supports presenting appropriate moieties for attachment of the binding partner may be produced by routine methods known in the art.

Attachment of appropriate functional groups to the oligonucleotide probes of the invention may be performed by ligation or introduced during synthesis or amplification, for example using primers carrying an appropriate moiety, such as biotin or a particular sequence for capture.

Conveniently, the set of probes described hereinbefore is provided in kit form.

Thus viewed from a further aspect the present invention provides a kit comprising a set of oligonucleotide probes as described hereinbefore immobilized on one or more solid supports.

Preferably, said probes are immobilized on a single solid support and each unique probe is attached to a different region of said solid support. However, when attached to multiple solid supports, said multiple solid supports form the modules which make up the kit. Especially preferably said solid support is a sheet, filter, membrane, plate or biochip.

Optionally the kit may also contain information

- 16 -

relating to the signals generated by normal or diseased samples (as discussed in more detail hereinafter in relation to the use of the kits), standardizing materials, e.g. mRNA or cDNA from normal and/or diseased samples for comparative purposes, labels for incorporation into cDNA, adapters for introducing nucleic acid sequences for amplification purposes, primers for amplification and/or appropriate enzymes, buffers and solutions. Optionally said kit may also contain a package insert describing how the method of the invention should be performed, optionally providing standard graphs, data or software for interpretation of results obtained when performing the invention.

The use of such kits to prepare a standard diagnostic gene transcript pattern as described hereinafter forms a further aspect of the invention.

The set of probes as described herein have various uses. Principally however they are used to assess the gene expression state of a test cell to provide information relating to the organism from which said cell is derived. Thus the probes are useful in diagnosing, identifying or monitoring a disease or condition or stage thereof in an organism.

Thus in a further aspect the invention provides the use of a set of oligonucleotide probes or a kit as described hereinbefore to determine the gene expression pattern of a cell which pattern reflects the level of gene expression of genes to which said oligonucleotide probes bind, comprising at least the steps of:

- a) isolating mRNA from said cell, which may optionally be reverse transcribed to cDNA;
- b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotide probes or a kit as defined herein; and
- c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce said pattern.

The mRNA and cDNA as referred to in this method,

- 17 -

and the methods hereinafter, encompass derivatives or copies of said molecules, e.g. copies of such molecules such as those produced by amplification or the preparation of complementary strands, but which retain the identity of the mRNA sequence, ie. would hybridize to the direct transcript (or its complementary sequence) by virtue of precise complementarity, or sequence identity, over at least a region of said molecule. It will be appreciated that complementarity will not exist over the entire region where techniques have been used which may truncate the transcript or introduce new sequences, e.g. by primer amplification. For convenience, said mRNA or cDNA is preferably amplified prior to step b). As with the oligonucleotides described herein said molecules may be modified, e.g. by using non-natural bases during synthesis providing complementarity remains. Such molecules may also carry additional moieties such as signalling or immobilizing means.

The various steps involved in the method of preparing such a pattern are described in more detail hereinafter.

As used herein "gene expression" refers to transcription of a particular gene to produce a specific mRNA product (ie. a particular splicing product). The level of gene expression may be determined by assessing the level of transcribed mRNA molecules or cDNA molecules reverse transcribed from the mRNA molecules or products derived from those molecules, e.g. by amplification.

The "pattern" created by this technique refers to information which, for example, may be represented in tabular or graphical form and conveys information about the signal associated with two or more oligonucleotides. Preferably said pattern is expressed as an array of numbers relating to the expression level associated with each probe.

- 18 -

Preferably, said pattern is established using the following linear model:

$$y = Xb + f \quad \text{Equation 1}$$

wherein, X is the matrix of gene expression data and y is the response variable, b is the regression coefficient vector and f the estimated residual vector. Although many different methods can be used to establish the relationship provided in equation 1, especially preferably the partial Least Squares Regression (PLSR) method is used for establishing the relationship in equation 1.

The probes are thus used to generate a pattern which reflects the gene expression of a cell at the time of its isolation. The pattern of expression is characteristic of the circumstances under which that cells finds itself and depends on the influences to which the cell has been exposed. Thus, a characteristic gene transcript pattern standard or fingerprint (standard probe pattern) for cells from an individual with a particular disease or condition may be prepared and used for comparison to transcript patterns of test cells. This has clear applications in diagnosing, monitoring or identifying whether an organism is suffering from a particular disease, condition or stage thereof.

The standard pattern is prepared by determining the extent of binding of total mRNA (or cDNA or related product), from cells from a sample of one or more organisms with the disease or condition or stage thereof, to the probes. This reflects the level of transcripts which are present which correspond to each unique probe. The amount of nucleic acid material which binds to the different probes is assessed and this information together forms the gene transcript pattern standard of that disease or condition or stage thereof. Each such standard pattern is characteristic of the disease, condition or stage thereof.

- 19 -

In a further aspect therefore, the present invention provides a method of preparing a standard gene transcript pattern characteristic of a disease or condition or stage thereof in an organism comprising at least the steps of:

- a) isolating mRNA from the cells of a sample of one or more organisms having the disease or condition or stage thereof, which may optionally be reverse transcribed to cDNA;
- b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotides or a kit as described hereinbefore specific for said disease or condition or stage thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and
- c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce a characteristic pattern reflecting the level of gene expression of genes to which said oligonucleotides bind, in the sample with the disease, condition or stage thereof.

For convenience, said oligonucleotides are preferably immobilized on one or more solid supports.

The standard pattern for a great number of diseases or conditions and different stages thereof using particular probes may be accumulated in databases and be made available to laboratories on request.

"Disease" samples and organisms as referred to herein refer to organisms (or samples from the same) with an underlying pathological disturbance relative to a normal organism (or sample), in a symptomatic or asymptomatic organism, which may result, for example, from infection or an acquired or congenital genetic imperfection. Such organisms are known to have, or which exhibit, the disease or condition or stage thereof under study.

A "condition" refers to a state of the mind or body of an organism which has not occurred through disease,

- 20 -

e.g. the presence of an agent in the body such as a toxin, drug or pollutant, or pregnancy.

"Stages" thereof refer to different stages of the disease or condition which may or may not exhibit particular physiological or metabolic changes, but do exhibit changes at the genetic level which may be detected as altered gene expression. It will be appreciated that during the course of a disease or condition the expression of different transcripts may vary. Thus at different stages, altered expression may not be exhibited for particular transcripts compared to "normal" samples. However, combining information from several transcripts which exhibit altered expression at one or more stages through the course of the disease or condition can be used to provide a characteristic pattern which is indicative of a particular stage of the disease or condition. Thus for example different stages in cancer, e.g. pre-stage I, stage I, stage II, II or IV can be identified.

"Normal" as used herein refers to organisms or samples which are used for comparative purposes. Preferably, these are "normal" in the sense that they do not exhibit any indication of, or are not believed to have, any disease or condition that would affect gene expression, particularly in respect of the disease for which they are to be used as the normal standard. However, it will be appreciated that different stages of a disease or condition may be compared and in such cases, the "normal" sample may correspond to the earlier stage of the disease or condition.

As used herein a "sample" refers to any material obtained from the organism, e.g. human or non-human animal under investigation which contains cells and includes, tissues, body fluid or body waste or in the case of prokaryotic organisms, the organism itself. "Body fluids" include blood, saliva, spinal fluid, semen, lymph. "Body waste" includes urine, expectorated

- 21 -

matter (pulmonary patients), faeces etc. "Tissue samples" include tissue obtained by biopsy, by surgical interventions or by other means e.g. placenta. Preferably however, the samples which are examined are
5 from areas of the body not apparently affected by the disease or condition. The cells in such samples are not disease cells, e.g. cancer cells, have not been in contact with such disease cells and do not originate from the site of the disease or condition. The "site of
10 disease" is considered to be that area of the body which manifests the disease in a way which may be objectively determined, e.g. a tumour or area of inflammation. Thus for example peripheral blood may be used for the diagnosis of non-haematopoietic cancers, and the blood
15 does not require the presence of malignant or disseminated cells from the cancer in the blood. Similarly in diseases of the brain, in which no diseased cells are found in the blood due to the blood:brain barrier, peripheral blood may still be used in the
20 methods of the invention.

It will however be appreciated that the method of preparing the standard transcription pattern and other methods of the invention are also applicable for use on living parts of eukaryotic organisms such as cell lines
25 and organ cultures and explants.

As used herein, reference to "corresponding" sample etc. refers to cells preferably from the same tissue, body fluid or body waste, but also includes cells from tissue, body fluid or body waste which are sufficiently
30 similar for the purposes of preparing the standard or test pattern. When used in reference to genes "corresponding" to the probes, this refers to genes which are related by sequence (which may be complementary) to the probes although the probes may
35 reflect different splicing products of expression.

"Assessing" as used herein refers to both quantitative and qualitative assessment which may be

- 22 -

determined in absolute or relative terms.

The invention may be put into practice as follows. To prepare a standard transcript pattern for a particular disease, condition or stage thereof, sample
5 mRNA is extracted from the cells of tissues, body fluid or body waste according to known techniques (see for example Sambrook et. al. (1989), Molecular Cloning : A laboratory manual, 2nd Ed., Cold Spring Harbor Laboratory Press, Cold Spring Harbor, N.Y.) from a
10 diseased individual or organism.

Owing to the difficulties in working with RNA, the RNA is preferably reverse transcribed at this stage to form first strand cDNA. Cloning of the cDNA or selection from, or using, a cDNA library is not however
15 necessary in this or other methods of the invention. Preferably, the complementary strands of the first strand cDNAs are synthesized, ie. second strand cDNAs, but this will depend on which relative strands are present in the oligonucleotide probes. The RNA may
20 however alternatively be used directly without reverse transcription and may be labelled if so required.

Preferably the cDNA strands are amplified by known amplification techniques such as the polymerase chain reaction (PCR) by the use of appropriate primers.
25 Alternatively, the cDNA strands may be cloned with a vector, used to transform a bacteria such as E. coli which may then be grown to multiply the nucleic acid molecules. When the sequence of the cDNAs are not known, primers may be directed to regions of the nucleic
30 acid molecules which have been introduced. Thus for example, adapters may be ligated to the cDNA molecules and primers directed to these portions for amplification of the cDNA molecules. Alternatively, in the case of eukaryotic samples, advantage may be taken of the polyA
35 tail and cap of the RNA to prepare appropriate primers.

To produce the standard diagnostic gene transcript pattern or fingerprint for a particular disease or

- 23 -

condition or stage thereof, the above described oligonucleotide probes are used to probe mRNA or cDNA of the diseased sample to produce a signal for hybridization to each particular oligonucleotide probe species, ie. each unique probe. A standard control gene transcript pattern may also be prepared if desired using mRNA or cDNA from a normal sample. Thus, mRNA or cDNA is brought into contact with the oligonucleotide probe under appropriate conditions to allow hybridization.

When multiple samples are probed, this may be performed consecutively using the same probes, e.g. on one or more solid supports, ie. on probe kit modules, or by simultaneously hybridizing to corresponding probes, e.g. the modules of a corresponding probe kit.

To identify when hybridization occurs and obtain an indication of the number of transcripts/cDNA molecules which become bound to the oligonucleotide probes, it is necessary to identify a signal produced when the transcripts (or related molecules) hybridize (e.g. by detection of double stranded nucleic acid molecules or detection of the number of molecules which become bound, after removing unbound molecules, e.g. by washing).

In order to achieve a signal, either or both components which hybridize (ie. the probe and the transcript) carry or form a signalling means or a part thereof. This "signalling means" is any moiety capable of direct or indirect detection by the generation or presence of a signal. The signal may be any detectable physical characteristic such as conferred by radiation emission, scattering or absorption properties, magnetic properties, or other physical properties such as charge, size or binding properties of existing molecules (e.g. labels) or molecules which may be generated (e.g. gas emission etc.). Techniques are preferred which allow signal amplification, e.g. which produce multiple signal events from a single active binding site, e.g. by the catalytic action of enzymes to produce multiple

- 24 -

detectable products.

Conveniently the signalling means may be a label which itself provides a detectable signal. Conveniently this may be achieved by the use of a radioactive or
5 other label which may be incorporated during cDNA production, the preparation of complementary cDNA strands, during amplification of the target mRNA/cDNA or added directly to target nucleic acid molecules.

Appropriate labels are those which directly or
10 indirectly allow detection or measurement of the presence of the transcripts/cDNA. Such labels include for example radiolabels, chemical labels, for example chromophores or fluorophores (e.g. dyes such as fluorescein and rhodamine), or reagents of high electron
15 density such as ferritin, haemocyanin or colloidal gold. Alternatively, the label may be an enzyme, for example peroxidase or alkaline phosphatase, wherein the presence of the enzyme is visualized by its interaction with a suitable entity, for example a substrate. The label may
20 also form part of a signalling pair wherein the other member of the pair is found on, or in close proximity to, the oligonucleotide probe to which the transcript/cDNA binds, for example, a fluorescent compound and a quench fluorescent substrate may be used.
25 A label may also be provided on a different entity, such as an antibody, which recognizes a peptide moiety attached to the transcripts/cDNA, for example attached to a base used during synthesis or amplification.

A signal may be achieved by the introduction of a
30 label before, during or after the hybridization step. Alternatively, the presence of hybridizing transcripts may be identified by other physical properties, such as their absorbance, and in which case the signalling means is the complex itself.

35 The amount of signal associated with each oligonucleotide probe is then assessed. The assessment may be quantitative or qualitative and may be based on

- 25 -

binding of a single transcript species (or related cDNA or other products) to each probe, or binding of multiple transcript species to multiple copies of each unique probe. It will be appreciated that quantitative results
5 will provide further information for the transcript fingerprint of the disease which is compiled. This data may be expressed as absolute values (in the case of macroarrays) or may be determined relative to a particular standard or reference e.g. a normal control
10 sample.

Furthermore it will be appreciated that the standard diagnostic gene pattern transcript may be prepared using one or more disease samples (and normal samples if used) to perform the hybridization step to
15 obtain patterns not biased towards a particular individual's variations in gene expression.

The use of the probes to prepare standard patterns and the standard diagnostic gene transcript patterns thus produced for the purpose of identification or
20 diagnosis or monitoring of a particular disease or condition or stage thereof in a particular organism forms a further aspect of the invention.

Once a standard diagnostic fingerprint or pattern has been determined for a particular disease or
25 condition using the selected oligonucleotide probes, this information can be used to identify the presence, absence or extent or stage of that disease or condition in a different test organism or individual.

To examine the gene expression pattern of a test
30 sample, a test sample of tissue, body fluid or body waste containing cells, corresponding to the sample used for the preparation of the standard pattern, is obtained from a patient or the organism to be studied. A test gene transcript pattern is then prepared as described
35 hereinbefore as for the standard pattern.

In a further aspect therefore, the present invention provides a method of preparing a test gene

- 26 -

transcript pattern comprising at least the steps of:

a) isolating mRNA from the cells of a sample of said test organism, which may optionally be reverse transcribed to cDNA;

5 b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotides or a kit as described hereinbefore specific for a disease or condition or stage thereof in an organism and sample thereof corresponding to the organism and sample thereof under
10 investigation; and

c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce said pattern reflecting the level of gene expression of genes to which said oligonucleotides bind, in said test sample.

15 This test pattern may then be compared to one or more standard patterns to assess whether the sample contains cells having the disease, condition or stage thereof.

 Thus viewed from a further aspect the present
20 invention provides a method of diagnosing or identifying or monitoring a disease or condition or stage thereof in an organism, comprising the steps of:

a) isolating mRNA from the cells of a sample of said organism, which may optionally be reverse
25 transcribed to cDNA;

b) hybridizing the mRNA or cDNA of step (a) to a set of oligonucleotides or a kit as described hereinbefore specific for said disease or
30 condition or stage thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation;

c) assessing the amount of mRNA or cDNA hybridizing to each of said probes to produce a characteristic pattern reflecting the level
35 of gene expression of genes to which said oligonucleotides bind, in said sample; and

d) comparing said pattern to a standard

- 27 -

5 diagnostic pattern prepared according to the method of the invention using a sample from an organism corresponding to the organism and sample under investigation to determine the presence of said disease or condition or a stage thereof in the organism under investigation.

The method up to and including step c) is the preparation of a test pattern as described above.

10 As referred to herein, "diagnosis" refers to determination of the presence or existence of a disease or condition or stage thereof in an organism. "Monitoring" refers to establishing the extent of a disease or condition, particularly when an individual is
15 known to be suffering from a disease or condition, for example to monitor the effects of treatment or the development of a disease or condition, e.g. to determine the suitability of a treatment or provide a prognosis.

20 The presence of the disease or condition or stage thereof may be determined by determining the degree of correlation between the standard and test samples' patterns. This necessarily takes into account the range of values which are obtained for normal and diseased samples. Although this can be established by obtaining
25 standard deviations for several representative samples binding to the probes to develop the standard, it will be appreciated that single samples may be sufficient to generate the standard pattern to identify a disease if the test sample exhibits close enough correlation to
30 that standard. Conveniently, the presence, absence, or extent of a disease or condition or stage thereof in a test sample can be predicted by inserting the data relating to the expression level of informative probes in test sample into the standard diagnostic probe
35 pattern established according to equation 1.

Data generated using the above mentioned methods may be analysed using various techniques from the most

- 28 -

basic visual representation (e.g. relating to intensity) to more complex data manipulation to identify underlying patterns which reflect the interrelationship of the level of expression of each gene to which the various probes bind, which may be quantified and expressed mathematically. Conveniently, the raw data thus generated may be manipulated by the data processing and statistical methods described hereinafter, particularly normalizing and standardizing the data and fitting the data to a classification model to determine whether said test data reflects the pattern of a particular disease, condition or stage thereof.

The methods described herein may be used to identify, monitor or diagnose a disease, condition or ailment or its stage or progression, for which the oligonucleotide probes are informative. "Informative" probes as described herein, are those which reflect genes which have altered expression in the diseases or conditions in question, or particular stages thereof. Probes of the invention may not be sufficiently informative for diagnostic purposes when used alone, but are informative when used as one of several probes to provide a characteristic pattern, e.g. in a set as described hereinbefore.

Preferably said probes correspond to genes which are systemically affected by said disease, condition or stage thereof. Especially preferably said genes, from which transcripts are derived which bind to probes of the invention, are metabolic or house-keeping genes and preferably are moderately or highly expressed. The advantage of using probes directed to moderately or highly expressed genes is that smaller clinical samples are required for generating the necessary gene expression data set, e.g. less than 1ml blood samples.

Furthermore, it has been found that such genes which are already being actively transcribed tend to be more prone to being influenced, in a positive or

- 29 -

negative way, by new stimuli. In addition, since transcripts are already being produced at levels which are generally detectable, small changes in those levels are readily detectable as for example, a certain
5 detectable threshold does not need to be reached.

In preferred methods of the invention, the set of probes of the invention are informative for a variety of different diseases, conditions or stages thereof. A sub-set of the probes disclosed herein may be used for
10 diagnosis, identification or monitoring a particular disease, condition or stage thereof.

Thus the probes may be used to diagnose or identify or monitor any condition, ailment, disease or reaction that leads to the relative increase or decrease in the
15 activity of informative genes of any or all eukaryotic or prokaryotic organisms regardless of whether these changes have been caused by the influence of bacteria, virus, prions, parasites, fungi, radiation, natural or artificial toxins, drugs or allergens, including mental
20 conditions due to stress, neurosis, psychosis or deteriorations due to the ageing of the organism, and conditions or diseases of unknown cause, providing a sub-set of the probes as described herein are informative for said disease or condition or stage
25 thereof.

Such diseases include those which result in metabolic or physiological changes, such as fever-associated diseases such as influenza or malaria. Other diseases which may be detected include for example
30 yellow fever, sexually transmitted diseases such as gonorrhea, fibromyalgia, candida-related complex, cancer (for example of the stomach, lung, breast, prostate gland, bowel, skin, colon, ovary etc), Alzheimer's disease, disease caused by retroviruses such as HIV,
35 senile dementia, multiple sclerosis and Creutzfeldt-Jakob disease to mention a few.

The invention may also be used to identify patients

- 30 -

with psychiatric or psychosomatic diseases such as schizophrenia and eating disorders. Of particular importance is the use of this method to detect diseases, conditions, or stages thereof, which are not readily detectable by known diagnostic methods, such as HIV which is generally not detectable using known techniques 1 to 4 months following infection. Conditions which may be identified include for example drug abuse, such as the use of narcotics, alcohol, steroids or performance enhancing drugs.

Preferably said disease to be identified or monitored is a cancer or a degenerative brain disorder (such as Alzheimer's or Parkinson's disease).

In particular, a set of oligonucleotide probes, wherein said set comprises at least 10 oligonucleotides selected from:

an oligonucleotide as described in Table 4 or an oligonucleotide derived therefrom or an oligonucleotide with a complementary sequence, or a functionally equivalent oligonucleotide, may be used for diagnosis or identification or monitoring the progression of Alzheimer's disease. Similarly Table 2 probes and Table 2 derived probes and their functional equivalents may be used to diagnose, identify or monitor the progression of breast cancer. Especially preferably the probes used for breast cancer analysis are selected based on their occurrence as set forth in Table 3 and as described hereinbefore.

The diagnostic method may be used alone as an alternative to other diagnostic techniques or in addition to such techniques. For example, methods of the invention may be used as an alternative or additive diagnostic measure to diagnosis using imaging techniques such as Magnetic Resonance Image (MRI), ultrasound imaging, nuclear imaging or X-ray imaging, for example in the identification and/or diagnosis of tumours.

The methods of the invention may be performed on

- 31 -

cells from prokaryotic or eukaryotic organisms which may be any eukaryotic organisms such as human beings, other mammals and animals, birds, insects, fish and plants, and any prokaryotic organism such as a bacteria.

5 Preferred non-human animals on which the methods of the invention may be conducted include, but are not limited to mammals, particularly primates, domestic animals, livestock and laboratory animals. Thus preferred animals for diagnosis include mice, rats,
10 guinea pigs, cats, dogs, pigs, cows, goats, sheep, horses. Particularly preferably the disease state or condition of humans is diagnosed, identified or monitored.

As described above, the sample under study may be
15 any convenient sample which may be obtained from an organism. Preferably however, as mentioned above, the sample is obtained from a site distant to the site of disease and the cells in such samples are not disease cells, have not been in contact with such cells and do
20 not originate from the site of the disease or condition. In such cases, although preferably absent, the sample may contain cells which do not fulfil these criteria. However, since the probes of the invention are concerned with transcripts whose expression is altered in cells
25 which do satisfy these criteria, the probes are specifically directed to detecting changes in transcript levels in those cells even if in the presence of other, background cells.

It has been found that the cells from such samples
30 show significant and informative variations in the gene expression of a large number of genes. Thus, the same probe (or several probes) may be found to be informative in determinations regarding two or more diseases, conditions or stages thereof by virtue of the particular
35 level of transcripts binding to that probe or the interrelationship of the extent of binding to that probe relative to other probes. As a consequence, it is

- 32 -

possible to use a relatively small number of probes for screening for multiple disorders or diseases. This has consequences with regard to the selection of probes, discussed in relation to random identification of probes hereinafter, but also for the use of a single set of probes for more than one diagnosis. Table 9 which represents preferred probes of the invention discloses probes which are informative for both Alzheimer's and breast cancer.

Thus, the present invention also provides sets of probes for diagnosing, identifying or monitoring two or more diseases, conditions or stages thereof, wherein at least one of said probes is suitable for said diagnosing, identifying or monitoring at least two of said diseases, conditions or stages thereof, and kits and methods of using the same. Preferably at least 5 probes, e.g. from 5 to 15 probes, are used in at least two diagnoses.

Thus, in a further preferred aspect, the present invention provides a method of diagnosis or identification or monitoring as described hereinbefore for the diagnosis, identification or monitoring of two or more diseases, conditions or stages thereof in an organism, wherein said test pattern produced in step c) of the diagnostic method is compared in step d) to at least two standard diagnostic patterns prepared as described previously, wherein each standard diagnostic pattern is a pattern generated for a different disease or condition or stage thereof.

Whilst in a preferred aspect the methods of assessment concern the development of a gene transcript pattern from a test sample and comparison of the same to a standard pattern, the elevation or depression of expression of certain markers may also be examined by examining the products of expression and the level of those products. Thus a standard pattern in relation to the expressed product may be generated:

- 33 -

In such methods the levels of expression of a set of polypeptides encoded by the gene to which an oligonucleotide of Table 1 or a Table 1 derived oligonucleotide, binds, are analysed.

5 Various diagnostic methods may be used to assess the amount of polypeptides (or fragments thereof) which are present. The presence or concentration of polypeptides may be examined, for example by the use of a binding partner to said polypeptide (e.g. an
10 antibody), which may be immobilized, to separate said polypeptide from the sample and the amount of polypeptide may then be determined.

15 "Fragments" of the polypeptides refers to a domain or region of said polypeptide, e.g. an antigenic fragment, which is recognizable as being derived from said polypeptide to allow binding of a specific binding partner. Preferably such a fragment comprises a significant portion of said polypeptide and corresponds to a product of normal post-synthesis processing.

20 Thus in a further aspect the present invention provides a method of preparing a standard gene transcript pattern characteristic of a disease or condition or stage thereof in an organism comprising at least the steps of:

25 a) releasing target polypeptides from a sample of one or more organisms having the disease or condition or stage thereof;

30 b) contacting said target polypeptides with one or more binding partners, wherein each binding partner is specific to a marker polypeptide (or a fragment thereof) encoded by the gene to which an oligonucleotide of Table 1 (or derived from a sequence described in Table 1) binds, to allow binding of said binding partners to said target polypeptides, wherein said marker polypeptides
35 are specific for said disease or condition thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and

- 34 -

c) assessing the target polypeptide binding to said binding partners to produce a characteristic pattern reflecting the level of gene expression of genes which express said marker polypeptides, in the sample with the disease, condition or stage thereof.

As used herein "target polypeptides" refer to those polypeptides present in a sample which are to be detected and "marker polypeptides" are polypeptides which are encoded by the genes to which Table 1 oligonucleotides or Table 1 derived oligonucleotides bind. The target and marker polypeptides are identical or at least have areas of high similarity, e.g. epitopic regions to allow recognition and binding of the binding partner.

"Release" of the target polypeptides refers to appropriate treatment of a sample to provide the polypeptides in a form accessible for binding of the binding partners, e.g. by lysis of cells where these are present. The samples used in this case need not necessarily comprise cells as the target polypeptides may be released from cells into the surrounding tissue or fluid, and this tissue or fluid may be analysed, e.g. urine or blood. Preferably however the preferred samples as described herein are used. "Binding partners" comprise the separate entities which together make an affinity binding pair as described above, wherein one partner of the binding pair is the target or marker polypeptide and the other partner binds specifically to that polypeptide, e.g. an antibody.

Various arrangements may be envisaged for detecting the amount of binding pairs which form. In its simplest form, a sandwich type assay e.g. an immunoassay such as an ELISA, may be used in which an antibody specific to the polypeptide and carrying a label (as described elsewhere herein) may be bound to the binding pair (e.g. the first antibody:polypeptide pair) and the amount of label detected.

- 35 -

Other methods as described herein may be similarly modified for analysis of the protein product of expression rather than the gene transcript and related nucleic acid molecules.

5 Thus a further aspect of the invention provides a method of preparing a test gene transcript pattern comprising at least the steps of:

- a) releasing target polypeptides from a sample of said test organism;
- 10 b) contacting said target polypeptides with one or more binding partners, wherein each binding partner is specific to a marker polypeptide (or a fragment thereof) encoded by the gene to which an oligonucleotide of Table 1 (or derived from a sequence described in Table 1)
- 15 binds, to allow binding of said binding partners to said target polypeptides, wherein said marker polypeptides are specific for said disease or condition thereof in an organism and sample thereof corresponding to the organism and sample thereof under investigation; and
- 20 c) assessing the target polypeptide binding to said binding partners to produce a characteristic pattern reflecting the level of gene expression of genes which express said marker polypeptides, in said test sample.

25 A yet further aspect of the invention provides a method of diagnosing or identifying or monitoring a disease or condition or stage thereof in an organism comprising the steps of:

- a) releasing target polypeptides from a sample of said organism;
- 30 b) contacting said target polypeptides with one or more binding partners, wherein each binding partner is specific to a marker polypeptide (or a fragment thereof) encoded by the gene to which an oligonucleotide of Table 1 (or derived from a sequence described in Table 1)
- 35 binds, to allow binding of said binding partners to said target polypeptides, wherein said marker polypeptides are specific for said disease or condition thereof in an

- 36 -

organism and sample thereof corresponding to the organism and sample thereof under investigation; and

5 c) assessing the target polypeptide binding to said binding partners to produce a characteristic pattern reflecting the level of gene expression of genes which express said marker polypeptides in said sample; and

10 d) comparing said pattern to a standard diagnostic pattern prepared as described hereinbefore using a sample from an organism corresponding to the organism and sample under investigation to determine the degree of correlation indicative of the presence of said disease or condition or a stage thereof in the organism under investigation.

15 The methods of generating standard and test patterns and diagnostic techniques rely on the use of informative oligonucleotide probes to generate the gene expression data. In some cases it will be necessary to select these informative probes for a particular method, e.g. to diagnose a particular disease, from a selection of available probes, e.g. the probes described
20 hereinbefore (the Table 1 oligonucleotides, the Table 1 derived oligonucleotides, their complementary sequences and functionally equivalent oligonucleotides). The following methodology describes a convenient method for
25 identifying such informative probes, or more particularly how to select a suitable sub-set of probes from the probes described herein.

30 Probes for the analysis of a particular disease or condition or stage thereof, may be identified in a number of ways known in the prior art, including by differential expression or by library subtraction (see for example WO98/49342). As described hereinafter, in view of the high information content of most transcripts, as a starting point one may also simply
35 analyse a random sub-set of mRNA or cDNA species and pick the most informative probes from that sub-set. The following method describes the use of immobilized

- 37 -

oligonucleotide probes (e.g. the probes of the invention) to which mRNA (or related molecules) from different samples is bound to identify which probes are the most informative to identify a particular type of sample, e.g. a disease sample.

The immobilized probes can be derived from various unrelated or related organisms; the only requirement is that the immobilized probes should bind specifically to their homologous counterparts in test organisms. Probes can also be derived from commercially available or public databases and immobilized on solid supports or, as mentioned above, they can be randomly picked and isolated from a cDNA library and immobilized on a solid support.

The length of the probes immobilised on the solid support should be long enough to allow for specific binding to the target sequences. The immobilised probes can be in the form of DNA, RNA or their modified products or PNAs (peptide nucleic acids). Preferably, the probes immobilised should bind specifically to their homologous counterparts representing highly and moderately expressed genes in test organisms. Conveniently the probes which are used are the probes described herein.

The gene expression pattern of cells in biological samples can be generated using prior art techniques such as microarray or macroarray as described below or using methods described herein. Several technologies have now been developed for monitoring the expression level of a large number of genes simultaneously in biological samples, such as, high-density oligoarrays (Lockhart et al., 1996, Nat. Biotech., 14, p1675-1680), cDNA microarrays (Schena et al, 1995, Science, 270, p467-470) and cDNA macroarrays (Maier E et al., 1994, Nucl. Acids Res., 22, p3423-3424; Bernard et al., 1996, Nucl. Acids Res., 24, p1435-1442).

In high-density oligoarrays and cDNA microarrays,

- 38 -

hundreds and thousands of probe oligonucleotides or cDNAs, are spotted onto glass slides or nylon membranes, or synthesized on biochips. The mRNA isolated from the test and reference samples are labelled by reverse transcription with a red or green fluorescent dye, mixed, and hybridised to the microarray. After washing, the bound fluorescent dyes are detected by a laser, producing two images, one for each dye. The resulting ratio of the red and green spots on the two images provides the information about the changes in expression levels of genes in the test and reference samples. Alternatively, single channel or multiple channel microarray studies can also be performed.

In cDNA macroarray, different cDNAs are spotted on a solid support such as nylon membranes in excess in relation to the amount of test mRNA that can hybridise to each spot. mRNA isolated from test samples is radio-labelled by reverse transcription and hybridised to the immobilised probe cDNA. After washing, the signals associated with labels hybridising specifically to immobilised probe cDNA are detected and quantified. The data obtained in macroarray contains information about the relative levels of transcripts present in the test samples. Whilst macroarrays are only suitable to monitor the expression of a limited number of genes, microarrays can be used to monitor the expression of several thousand genes simultaneously and is, therefore, a preferred choice for large-scale gene expression studies.

A macroarray technique for generating the gene expression data set has been used to illustrate the probe identification method described herein. For this purpose, mRNA is isolated from samples of interest and used to prepare labelled target molecules, e.g. mRNA or cDNA as described above. The labelled target molecules are then hybridised to probes immobilised on the solid support. Various solid supports can be used for the

- 39 -

purpose, as described previously. Following hybridization, unbound target molecules are removed and signals from target molecules hybridizing to immobilised probes quantified. If radio labelling is performed, 5 PhosphoImager can be used to generate an image file that can be used to generate a raw data set. Depending on the nature of label chosen for labelling the target molecules, other instruments can also be used, for example, when fluorescence is used for labelling, a 10 FluoroImager can be used to generate an image file from the hybridised target molecules.

The raw data corresponding to mean intensity, median intensity, or volume of the signals in each spot can be acquired from the image file using commercially 15 available software for image analysis. However, the acquired data needs to be corrected for background signals and normalized prior to analysis, since, several factors can affect the quality and quantity of the hybridising signals. For example, variations in the 20 quality and quantity of mRNA isolated from sample to sample, subtle variations in the efficiency of labelling target molecules during each reaction, and variations in the amount of unspecific binding between different macroarrays can all contribute to noise in the acquired 25 data set that must be corrected for prior to analysis.

Background correction can be performed in several ways. The lowest pixel intensity within a spot can be used for background subtraction or the mean or median of the line of pixels around the spots' outline can be used 30 for the purpose. One can also define an area representing the background intensity based on the signals generated from negative controls and use the average intensity of this area for background subtraction.

35 The background corrected data can then be transformed for stabilizing the variance in the data structure and normalized for the differences in probe

- 40 -

intensity. Several transformation techniques have been described in the literature and a brief overview can be found in Cui, Kerr and Churchill
<http://www.jax.org/research/churchill/research/expression/Cui-Transform.pdf>). Normalization can be performed by dividing the intensity of each spot with the collective intensity, average intensity or median intensity of all the spots in a macroarray or a group of spots in a macroarray in order to obtain the relative intensity of signals hybridising to immobilised probes in a macroarray. Several methods have been described for normalizing gene expression data (Richmond and Somerville, 2000, Current Opin. Plant Biol., 3, p108-116; Finkelstein et al., 2001, In "Methods of Microarray Data Analysis. Papers from CAMDA, Eds. Lin & Johnsom, Kluwer Academic, p57-68; Yang et al., 2001, In "Optical Technologies and Informatics", Eds. Bittner, Chen, Dorsel & Dougherty, Proceedings of SPIE, 4266, p141-152; Dudoit et al, 2000, J. Am. Stat. Ass., 97, p77-87; Alter et al 2000, supra; Newton et al., 2001, J. Comp. Biol., 8, p37-52). Generally, a scaling factor or function is first calculated to correct the intensity effect and then used for normalising the intensities. The use of external controls has also been suggested for improved normalization.

One other major challenge encountered in large-scale gene expression analysis is that of standardization of data collected from experiments performed at different times. We have observed that gene expression data for samples acquired in the same experiment can be efficiently compared following background correction and normalization. However, the data from samples acquired in experiments performed at different times requires further standardization prior to analysis. This is because subtle differences in experimental parameters between different experiments, for example, differences in the quality and quantity of

- 41 -

mRNA extracted at different times, differences in time used for target molecule labelling, hybridization time or exposure time, can affect the measured values. Also, factors such as the nature of the sequence of
5 transcripts under investigation (their GC content) and their amount in relation to the each other determines how they are affected by subtle variations in the experimental processes. They determine, for example, how efficiently first strand cDNAs, corresponding to a
10 particular transcript, are transcribed and labelled during first strand synthesis, or how efficiently the corresponding labelled target molecules bind to their complementary sequences during hybridization. Batch to batch difference in the printing process is also a major
15 factor for variation in the generated expression data.

Failure to properly address and rectify for these influences leads to situations where the differences between the experimental series may overshadow the main information of interest contained in the gene expression
20 data set, i.e. the differences within the combined data from the different experimental series. Figure 1 provides one such example showing a classification based on Principal Component Analysis (PCA) of combined data from two experimental series where the main goal is to
25 distinguish between Alzheimer/non-Alzheimer patients.

PCA (also known as singular value decomposition) is a technique for studying interdependencies and underlying relationships of a set of variables. The data are modelled in terms of a few significant factors
30 or principal components (PC's), plus residuals. The PC's contain the main phenomena and define the systematic variability present in the data, while the residuals represent the variability interpreted as noise. Details on PCA can be found in Jolliffe (1986,
35 Principal Component Analysis, Springer-Verlag, NY), and Jackson (1991, A User's Guide to Principal Components, Wiley, NY). The results of Figure 1 show that two

- 42 -

clusters are formed representing the data from two experimental series rather than the Alzheimer/non-Alzheimer differentiation. There were eight samples in common between the two series of experiments, which ideally should have fallen on top of, or in near proximity to, each other if appropriately standardized.

We have now found that gene expression data between different experiments can be efficiently standardized by including a subset of samples from one experimental series in the next experimental series and using a direct standardization method (DS), originally described by Wang and Kowalski (Anal. Chem., 1991, 63, p2750 and J. Chemometrics, 1991, 5, p129-145). Although the method of DS is well known in the field of analytical chemistry, it remains undescribed and unused in the field of gene expression data analysis.

In DS, the secondary data representing for example experimental series 2 (secondary measurements, R_2) are corrected to match the data measured on the primary measurements representing data from series 1 (R_1), while the calibration model remains unchanged. In DS, response matrices for both experimental series are related to each other by a transformation matrix F , i.e.

$$R_1 = R_2 F \quad (1)$$

Where F is a square matrix dimensioned gene by gene. From (1), the transformation matrix is calculated as:

$$F = R_2^+ R_1 \quad (2)$$

The transformation matrix F in equation (2) is calculated using a relatively small subset of samples which are measured on both the master primary and the secondary series of data.

Finally, the response of the unknown sample measured on the secondary series $r_{2,un}^T$, is standardized

- 43 -

to the response vector $\hat{r}_{1,un}^T$ expected from the primary series

5

$$\hat{r}_{1,un}^T = r_{T_{21,un}}^T \hat{F} \quad (3)$$

From the preceding equation it can be seen that the column i of the transformation matrix contains the multiplication factors for a set of genes measured in the secondary series to obtain the intensity at spot i of the corrected series.

The number of samples that are repeated in the experimental series, R_1 and R_2 , should be equal to their ranks, which in this case is equal to the number of principal components retained for explaining the variation in the R_1 and R_2 . For example, if three principal components are retained for explaining the variation in the data set, a minimum of three samples should be repeated between R_1 and R_2 . The samples that should be repeated between different series should ideally be those that exhibit high leverages in the gene expression pattern. At times, two samples may suffice, while at other times, more than two samples should be ideally be included for good representativity. In some cases, the samples selected can be the same in all the experimental series to be compared (reference samples), while in other cases, representative samples can be selected sequentially by analyzing the expression pattern after each experiment. The selected samples with high leverages are then included in the next experimental series. The results of using Direct Standardization are shown in Figure 1.

Another approach for normalizing and standardizing the gene expression data set is to hybridize each DNA array with target molecules prepared from a test sample and an equal amount of labelled target molecules prepared from representative reference samples. In order to measure the intensity of labelled target

- 44 -

molecules hybridizing to the immobilized probes it is necessary that the labelled molecules are prepared from test and reference samples using different labels, for example, different fluorescent dyes can be used for
5 preparing the labelled material. The labelled molecules prepared from reference samples can be added to the hybridization solution together with the labelled material prepared from test samples. A data file from each array representing the expression pattern of
10 different genes in the test sample and reference samples can then be obtained, normalized and standardized by the direct standardization method as described above. An instant advantage of including the differentially labelled target molecules from reference samples during
15 hybridization is that it enables an efficient comparison of new test samples to the data sets already stored in a database.

Monitoring the expression of a large number of genes in several samples leads to the generation of a
20 large amount of data that is too complex to be easily interpreted. Several unsupervised and supervised multivariate data analysis techniques have already been shown to be useful in extracting meaningful biological information from these large data sets. Cluster
25 analysis is by far the most commonly used technique for gene expression analysis, and has been performed to identify genes that are regulated in a similar manner, and or identifying new/unknown tumour classes using gene expression profiles (Eisen et al., 1998, PNAS, 95,
30 p14863-14868, Alizadeh et al. 2000, supra, Perou et al. 2000, Nature, 406, p747-752; Ross et al, 2000, Nature Genetics, 24(3), p227-235; Herwig et al., 1999, Genome Res., 9, p1093-1105; Tamayo et al, 1999, Science, PNAS, 96, p2907-2912).

35 In the clustering method, genes are grouped into functional categories (clusters) based on their expression profile, satisfying two criteria: *homogeneity*

- 45 -

- the genes in the same cluster are highly similar in expression to each other; and *separation* - genes in different clusters have low similarity in expression to each other.

5 Examples of various clustering techniques that have been used for gene expression analysis include hierarchical clustering (Eisen et al., 1998, *supra*; Alizadeh et al. 2000, *supra*; Perou et al. 2000, *supra*; Ross et al, 2000, *supra*), K-means clustering (Herwig et
10 al., 1999, *supra*; Tavazoie et al, 1999, *Nature Genetics*, 22(3), p. 281-285), gene shaving (Hastie et al., 2000, *Genome Biology*, 1(2), research 0003.1-0003.21), block clustering (Tibshirani et al., 1999, Tech report Univ Stanford.) Plaid model (Lazzeroni, 2002, *Stat. Sinica*,
15 12, p61-86), and self-organizing maps (Tamayo et al. 1999, *supra*). Also, related methods of multivariate statistical analysis, such as those using the singular value decomposition (Alter et al., 2000, *PNAS*, 97(18), p10101-10106; Ross et al. 2000, *supra*) or
20 multidimensional scaling can be effective at reducing the dimensions of the objects under study.

 However, methods such as cluster analysis and singular value decomposition are purely exploratory and only provide a broad overview of the internal structure
25 present in the data. They are unsupervised approaches in which the available information concerning the nature of the class under investigation is not used in the analysis. Often, the nature of the biological perturbation to which a particular sample has been
30 subjected is known. For example, it is sometimes known whether the sample whose gene expression pattern is being analysed derives from a diseased or healthy individual. In such instances, discriminant analysis can be used for classifying samples into various groups
35 based on their gene expression data.

 In such an analysis one builds the classifier by training the data that is capable of discriminating

- 46 -

between member and non-members of a given class. The trained classifier can then be used to predict the class of unknown samples. Examples of discrimination methods that have been described in the literature include

5 Support Vector Machines (Brown et al, 2000, PNAS, 97, p262-267), Nearest Neighbour (Dudoit et al., 2000, supra), Classification trees (Dudoit et al., 2000, supra), Voted classification (Dudoit et al., 2000, supra), Weighted Gene voting (Golub et al. 1999, supra),

10 and Bayesian classification (Keller et al. 2000, Tec report Univ of Washington). Also a technique in which PLS (Partial Least Square) regression analysis is first used to reduce the dimensions in the gene expression data set followed by classification using logistic

15 discriminant analysis and quadratic discriminant analysis (LD and QDA) has recently been described (Nguyen & Rocke, 2002, Bioinformatics, 18, p39-50 and 1216-1226).

A challenge that gene expression data poses to

20 classical discriminatory methods is that the number of genes whose expression are being analysed is very large compared to the number of samples being analysed. However in most cases only a small fraction of these genes are informative in discriminant analysis problems.

25 Moreover, there is a danger that the noise from irrelevant genes can mask or distort the information from the informative genes. Several methods have been suggested in literature to identify and select genes that are informative in microarray studies, for example,

30 t-statistics (Dudoit et al, 2002, J. Am. Stat. Ass., 97, p77-87), analysis of variance (Kerr et al., 2000, PNAS, 98, p8961-8965), Neighbourhood analysis (Golub et al, 1999, supra), Ratio of between groups to within groups sum of squares (Dudoit et al., 2002, supra), Non

35 parametric scoring (Park et al., 2002, Pacific Symposium on Biocomputing, p52-63) and Likelihood selection (Keller et al., 2000, supra).

- 47 -

In the methods described herein the gene expression data that has been normalized and standardized is analysed by using Partial Least Squares Regression (PLSR). Although PLSR is primarily a method used for regression analysis of continuous data (see Appendix A), it can also be utilized as a method for model building and discriminant analysis using a dummy response matrix based on a binary coding. The class assignment is based on a simple dichotomous distinction such as breast cancer (class 1) / healthy (class 2), or a multiple distinction based on multiple disease diagnosis such as breast cancer (class 1) / Alzheimer (class 2) / healthy (class 3). The list of diseases for classification can be increased depending upon the samples available corresponding to other diseases or conditions or stages thereof.

PLSR applied as a classification method is referred to as PLS-DA (DA standing for Discriminant analysis). PLS-DA is an extension of the PLSR algorithm in which the Y-matrix is a dummy matrix containing n rows (corresponding to the number of samples) and K columns (corresponding to the number of classes). The Y-matrix is constructed by inserting 1 in the k th column and -1 in all the other columns if the corresponding i th object of X belongs to class k . By regressing Y onto X , classification of a new sample is achieved by selecting the group corresponding to the largest component of the fitted, $\hat{y}(x) = (\hat{y}_1(x), \hat{y}_2(x), \dots, \hat{y}_k(x))$. Thus, in a -1/1 response matrix, a prediction value below 0 means that the sample belongs to the class designated as -1, while a prediction value above 0 implies that the sample belongs to the class designated as 1.

An advantage of PLSR-DA is that the results obtained can be easily represented in the form of two different plots, the score and loading plots. Score plots represent a projection of the samples onto the principal components and shows the distribution of the

- 48 -

samples in the classification model and their relationship to one another. Loading plots display correlations between the variables present in the data set.

5 It is usually recommended to use PLS-DA as a starting point for the classification problem due to its ability to handle collinear data, and the property of PLSR as a dimension reduction technique. Once this purpose has been satisfied, it is possible to use other
10 methods such as Linear discriminant analysis, LDA, that has been shown to be effective in extracting further information, Indahl et al. (1999, Chem. and Intell. Lab. Syst., 49, p19-31). This approach is based on first decomposing the data using PLS-DA, and then using the
15 scores vectors (instead of the original variables) as input to LDA. Further details on LDA can be found in Duda and Hart (Classification and Scene Analysis, 1973, Wiley, USA).

 The next step following model building is of model
20 validation. This step is considered to be amongst the most important aspects of multivariate analysis, and tests the "goodness" of the calibration model which has been built. In this work, a cross validation approach has been used for validation. In this approach, one or
25 a few samples are kept out in each segment while the model is built using a full cross-validation on the basis of the remaining data. The samples left out are then used for prediction/classification. Repeating the simple cross-validation process several times holding
30 different samples out for each cross-validation leads to a so-called double cross-validation procedure. This approach has been shown to work well with a limited amount of data, as is the case in some of the Examples described here. Also, since the cross validation step
35 is repeated several times the dangers of model bias and overfitting are reduced.

 Once a calibration model has been built and

- 49 -

validated, genes exhibiting an expression pattern that is most relevant for describing the desired information in the model can be selected by techniques described in the prior art for variable selection, as mentioned
 5 elsewhere. Variable selection will help in reducing the final model complexity, provide a parsimonious model, and thus lead to a reliable model that can be used for prediction. Moreover, use of fewer genes for the purpose of providing diagnosis will reduce the cost of
 10 the diagnostic product. In this way informative probes which would bind to the genes of relevance may be identified.

We have found that after a calibration model has been built, statistical techniques like Jackknife
 15 (Effron, 1982, The Jackknife, the Bootstrap and other resampling plans. Society for Industrial and Applied mathematics, Philadelphia, USA), based on resampling methodology, can be efficiently used to select or confirm significant variables (informative probes).

20 The approximate uncertainty variance of the PLS regression coefficients B can be estimated by:

$$S^2B = \sum_{m=1}^M ((B-B_m)g)^2$$

25

where

S^2B = estimated uncertainty variance of B ;

30 B = the regression coefficient at the cross validated rank A using all the N objects;

B_m = the regression coefficient at the rank A using all objects except the object(s) left out in cross validation segment m ; and

g = scaling coefficient (here: $g=1$).

35 In our approach, Jackknife has been implemented together with cross-validation. For each variable the difference between the B -coefficients B_i in a

- 50 -

cross-validated sub-model and B_{tot} for the total model is first calculated. The sum of the squares of the differences is then calculated in all sub-models to obtain an expression of the variance of the B_i estimate for a variable. The significance of the estimate of B_i is calculated using the t-test. Thus, the resulting regression coefficients can be presented with uncertainty limits that correspond to 2 Standard Deviations, and from that significant variables are detected.

No further details as to the implementation or use of this step are provided here since this has been implemented in commercially available software, The Unscrambler, CAMO ASA, Norway. Also, details on variable selection using Jackknife can be found in Westad & Martens (2000, J. Near Inf. Spectr., 8, p117-124).

The following approach can be used to select informative probes from a gene expression data set:

- a) keep out one unique sample (including its repetitions if present in the data set) per cross validation segment;
- b) build a calibration model (cross validated segment) on the remaining samples using PLSR-DA;
- c) select the significant genes for the model in step b) using the Jackknife criterion;
- d) repeat the above 3 steps until all the unique samples in the data set are kept out once (as described in step a). For example, if 75 unique samples are present in the data set, 75 different calibration models are built resulting in a collection of 75 different sets of significant probes;
- e) select the most significant variables using the frequency of occurrence criterion in the generated sets of significant probes in step d). For example, a set of probes appearing in all sets (100%) are more informative than probes appearing in only 50% of the

- 51 -

generated sets in step d).

Once the informative probes for a disease have been selected, a final model is made and validated. The two most commonly used ways of validating the model are cross-validation (CV) and test set validation. In cross-validation, the data is divided into k subsets. The model is then trained k times, each time leaving out one of the subsets from training, but using only the omitted subset to compute error criterion, RMSEP (Root Mean Square Error of Prediction). If k equals the sample size, this is called "leave-one-out" cross-validation. The idea of leaving one or a few samples out per validation segment is valid only in cases where the covariance between the various experiments is zero. Thus, one sample at-a-time approach can not be justified in situations containing replicates since keeping only one of the replicates out will introduce a systematic bias in our analysis. The correct approach in this case will be to leave out all replicates of the same samples at a time since that would satisfy assumptions of zero covariance between the CV-segments.

The second approach for model validation is to use a separate test-set for validating the calibration model. This requires running a separate set of experiments to be used as a test set. This is the preferred approach given that real test data are available.

The final model is then used to identify a disease, condition or stage thereof in test samples. For this purpose, expression data of selected informative genes is generated from test samples and then the final model is used to determine whether a sample belongs to a diseased or non-diseased class or has a condition or stage thereof.

Thus viewed from a yet further aspect the present invention provides a method of identifying probes useful for diagnosing or identifying or monitoring a disease or

- 52 -

condition or stage thereof in an organism, comprising the steps of:

- 5 a) immobilizing a set of oligonucleotide probes, preferably as described hereinbefore, on a solid support;
- b) isolating mRNA from a sample of a normal organism (normal sample), which may optionally be reverse transcribed to cDNA;
- 10 c) isolating mRNA from a sample from an organism, corresponding to the sample and organism of step (b), which is known to have said disease or condition or a stage thereof (diseased sample), which may optionally be reverse transcribed to cDNA;
- 15 d) hybridizing the mRNA or cDNA of steps (b) and (c) to said set of immobilized oligonucleotide probes of step (a); and
- e) assessing the amount of mRNA or cDNA hybridizing to each of said oligonucleotide probes to determine the level of gene expression of genes to which said oligonucleotide probes bind in said normal and diseased samples to generate a gene expression data set for each sample;
- 20 f) normalizing and standardizing said data set of step (e);
- g) constructing a calibration model for classification, preferably using the statistical techniques Partial Least Squares Discriminant Analysis (PLS-DA) and Linear Discriminant Analysis (LDA);
- 30 h) performing JackKnife analysis and identifying those oligonucleotide probes which are required for classification of said disease and normal samples into their respective groups.
- 35

Preferably a model for classification purposes is

- 53 -

generated by using the data relating to the probes identified according to the above described method. Preferably the sample is as described previously. Preferably the oligonucleotides which are immobilized in
5 step (a) are randomly selected as described below or are the probes as described hereinbefore. Such oligonucleotides may be of considerable length, e.g. if using cDNA (which is encompassed within the scope of the term "oligonucleotide"). The identification of such
10 cDNA molecules as useful probes allows the development of shorter oligonucleotides which reflect the specificity of the cDNA molecules but are easier to manufacture and manipulate.

The above described model may then be used to
15 generate and analyse data of test samples and thus may be used for the diagnostic methods of the invention. In such methods the data generated from the test sample provides the gene expression data set and this is normalized and standardized as described above. This is
20 then fitted to the calibration model described above to provide classification.

The method described herein can also be used to simultaneously select informative probes for several related and unrelated diseases or conditions. Depending
25 upon which diseases or conditions have been included in the calibration or training set, informative probes can be selected for the said diseases or conditions. The informative probes selected for one disease or condition may or may not be similar to the informative probes
30 selected for another disease or condition of interest. It is the pattern with which the selected genes are expressed in relation to each other during a disease, condition, or stage thereof, that determines whether or not they are informative for the disease, condition or
35 stage thereof.

In other words, informative genes are selected based on how their expression correlates with the

- 54 -

expression of other selected informative genes under the influence of responses generated by the disease, condition or stage thereof under investigation. In examples 1 and 2 provided hereinafter, 139 informative probes were selected for breast cancer diagnosis and 182 probes were selected for Alzheimer's disease diagnosis by training the gene expression data set of genes representing 1435 or 758 randomly picked cDNA clones for breast cancer/non breast cancer samples, or Alzheimer/non-Alzheimer samples, respectively. Among the probes selected for breast cancer and Alzheimer, about 10 probes were informative both for breast cancer and Alzheimer disease diagnosis.

For the purpose of isolating informative probes or identifying several related and unrelated diseases, conditions and stages thereof simultaneously, the gene expression data set must contain the information on how genes are expressed when the subject has a particular disease, condition or stage thereof under investigation. The data set is generated from a set of healthy or diseased samples, where a particular sample may contain the information of only one disease, condition or stages thereof or may also contain information about multiple diseases, conditions or stages thereof. For example, if the isolation of informative probes for Alzheimer disease, breast cancer and diabetes is sought, whole blood samples can be obtained from an Alzheimer patient who has breast cancer and diabetes. Hence, the method also teaches an efficient experimental design to reduce the number of samples required for isolating informative probes by selecting samples representing more than one disease, condition or stage thereof.

As mentioned previously, in view of the high information content of most transcripts, the identification and selection of informative probes for use in diagnosing, monitoring or identifying a particular disease, condition or stage thereof may be

- 55 -

dramatically simplified. Thus the pool of genes from which a selection may be made to identify informative probes may be radically reduced.

5 Unlike, in prior art technologies where informative probes are selected from a population of thousands of genes that are being expressed in a cell, like in microarray, in the method described herein, the informative probes are selected from a limited number of randomly obtained genes. For example, from a population
10 of 1435 cDNA clones, randomly picked from a human whole blood cDNA library, we were able to select 139 informative probes for breast cancer diagnosis (see Example 1 and Table 2).

Thus in a preferred aspect of the above mentioned
15 method of identifying probes useful for diagnosing or identifying or monitoring a disease or condition or stage thereof in an organism, said set of oligonucleotides which are immobilized in step (a) are randomly selected from a larger set of oligonucleotides,
20 e.g. from a cDNA library or other oligonucleotide pool, which may be, but is preferably not selected from the set provided herein. Preferably said larger set comprises oligonucleotides which correspond to moderately or highly expressed genes. Thus preferably
25 in methods of the invention, the set of oligonucleotides according to the invention are replaced with a set of oligonucleotides which are randomly selected, e.g. from commercially available oligonucleotide or cDNA libraries.

30 As referred to herein "random" refers to selection which is not biased based on the extent of information carried by the transcripts in relation to the disease, condition or organism under study, ie. without bias towards their likely utility as informative probes.
35 Whilst a random selection may be made from a pool of transcripts (or related products) which have been biased, e.g. to highly or moderately expressed

- 56 -

transcripts, preferably random selection is made from a pool of transcripts not biased or selected by a sequence-based criterion. The larger set may therefore contain oligonucleotides corresponding to highly and moderately expressed genes, or alternatively, may be enriched for those corresponding to the highly and moderately expressed genes.

Random selection from highly and moderately expressed genes can be achieved in a wide variety of ways. A strategy used in this work, but not limiting in itself involves randomly picking a significant number of cDNA clones from a cDNA library constructed from a biological specimen under investigation. Since, in a cDNA library, the cDNA clones corresponding to transcripts present in high or moderate amount are more frequently present than transcripts corresponding to cDNA present in low amount, the former will tend to be picked up more frequently than the latter. A pool of cDNA enriched for those corresponding to highly and moderately expressed genes can be isolated by this approach.

To identify genes that are expressed in high or moderate amount among the isolated population for use in methods of the invention, the information about the relative level of their transcripts in samples of interest can be generated using several prior art techniques. Both non-sequence based methods, such as differential display or RNA fingerprinting, and sequence-based methods such as microarrays or macroarrays can be used for the purpose. Alternatively, specific primer sequences for highly and moderately expressed genes can be designed and methods such as quantitative RT-PCR can be used to determine the levels of highly and moderately expressed genes. Hence, a skilled practitioner may use a variety of techniques which are known in the art for determining the relative level of mRNA in a biological sample.

- 57 -

Especially preferably the sample for the isolation of mRNA in the above described method is as described previously and is preferably not from the site of disease and the cells in said sample are not disease cells and have not contacted disease cells.

The following examples are given by way of illustration only in which the Figures referred to are as follows:

Figure 1 shows the effect of Direct Standardization (DS) on the Alzheimer data measured in two different series of experiments in which AD denotes Alzheimer's samples and A,B are non-Alzheimer's samples. The samples in both series have been labelled systematically as (xx_7/xx_8), whereas the corrected samples from series 8 (in b,c,d) have been labelled as (xx_c), thus, for example, AD2-7 denotes Alzheimer disease sample number 2 in experiment series 7. The circled spots represent the samples chosen as the transfer samples. The connecting lines in figures b,c,d show the proximity of the replicated samples after applying DS. The dashed lines in figures a,c,d represent the decision boundary separating the classes. These lines have not been drawn on the basis of any statistical criteria, but serve the purpose of visually separating the classes. All the four figures show scores plot (PC1-PC2) from PCA analysis based on (a) non-standardized data, (b) scores plot after direct standardization using 3 transfer samples, (c) scores plot after direct standardization using 4 transfer sample, (d) scores plot after direct standardization using 8 transfer samples;

Figure 2 shows the projection of normal (including benign) and breast cancer samples onto a classification model generated by PLSR-DA using the data of 44 informative genes, in which PC is the principal components and N and C are normal and breast cancer samples, respectively;

Figure 3 shows the projection of individuals with

- 58 -

and without Alzheimer's disease onto a classification model generated by PLSR-DA using 182 informative genes;

Figures 4, 6 and 8 show projection plots as Figure 2 in which the classification model is generated using 719, 111 and 345 cDNAs, respectively, wherein PC is the principal components, N denotes normal and B denotes breast cancer samples;

Figures 5, 7 and 9 show prediction plots based on 3 principal components using the data of 719, 111 and 345 cDNAs, respectively;

Figure 10 shows a projection plot as Figure 3 in which the classification model is generated using 520 cDNAs; and

Figure 11 is the prediction plot corresponding to Figure 10.

Example 1: Diagnosis of Breast Cancer

Methods

Whole blood was obtained from the arms of breast cancer patients and patients with benign tumours (Ullevål and Haukland hospitals in Norway). All of the patients with breast cancer had a malignant tumour of the breast (disease samples). Healthy blood was collected from the above two hospitals, or collected at a Health station at Ås, Norway or at DiaGenic AS, Norway, from the arms of female donors with no reported signs of breast cancer. The blood from healthy individuals or with benign tumours comprise the normal samples. The blood was either collected in tubes containing EDTA and stored immediately at -80°C or was collected in PAXgene tubes and stored for 12-24 hours at room temperature before finally storing them at -80°C before use. Further details of the breast cancer and benign tumour patients from which blood was taken is provided in Table 5.

- 59 -

mRNA was isolated from the blood of the 29 breast cancer patients and 46 normal donors and used to prepare labelled probes by reverse transcribing in the presence of $\alpha^{33}\text{P}$ -dATP. The first strand cDNA of the normal and diseased samples was bound, separately to 1435 cDNA clones immobilized on a solid support (nylon membrane). These cDNA clones were randomly picked, without any prior knowledge of their gene sequences, from a cDNA library constructed using whole blood of 550 healthy individuals (Clontech, Palo Alto, USA). These methods were conducted as follows.

For amplification of inserts, bacterial clones were grown in microtiter plates containing 150 μl LB with 50 $\mu\text{g/ml}$ carbenicillin, and incubated overnight with agitation at 37°C. To lyse the cells, 5 μl of each culture were diluted with 50 μl H₂O and incubated for 12 min. at 95°C. Of this mixture, 2 μl were subjected to a PCR reaction using 20 pmoles of M13 forward and reverse primer in presence of 1.5 mM MgCl₂. PCR reactions were performed with the following cycling protocol: 4 min. at 95°C, followed by 25 cycles of 1 min. at 94°C, 1 min. at 60°C and 3 min. at 72°C either in a RoboCycler[®] Temperature Cycler (Stratagene, La Jolla, USA) or DNA Engine Dyad Peltier Thermal Cycler (MJ Research Inc., Waltham, USA). The amplified products were denatured by incubating with NaOH (0.2 M, final concentration) for 30 min. and spotted onto Hybond-N+ membranes (Amersham Pharmacia Biotech, Little Chalfont, UK), using MicroGrid II workstation according to the manufacturer's instructions (BioRobotics Ltd, Cambridge England). The immobilized cDNAs were fixed using a UV cross-linker (Hoefer Scientific Instruments, San Francisco, USA).

In addition to the 1435 cDNAs, the printed arrays also contained controls for assessing background level, consistency and sensitivity of the assay. These were

- 60 -

spotted at multiple positions and included controls such as PCR mix (without any insert); positive and negative controls of SpotReport™ 10 array validation system (Stratagene, La Jolla, USA) and cDNAs corresponding to constitutively expressed genes such as b-actin, g-actin, GAPDH, HOD and cyclophilin. Also, oligonucleotides corresponding to SIX1, b-tubulin, TRP-2, MDM2, Myosin Light C, CD44, Maspin, Laminin, and SRP 19 were included to detect disseminated cancer cells.

The total RNA from blood collected in EDTA tubes was purified using Trizol LS Reagent protocol (Invitrogen/Life Technologies). From blood contained in PAXgene tubes, the total RNA was purified according to the supplier's instructions (PreAnalytiX, Hombrechtikon, Switzerland). Contaminating DNA was removed from the isolated RNA by DNAase I treatment using DNA-free kit (Ambion, Inc. Austin, USA). RNA quality was determined visually by inspecting the integrity of 28S and 18S ribosomal bands following agarose gel electrophoresis. The concentration and purity of extracted RNA was determined by measuring the absorbance at 260 nm and 280 nm. mRNA was isolated from the total RNA using Dynabeads as per the supplier's instructions (DynaL AS, Oslo, Norway).

Labelling and hybridization experiments were performed in batches. The number of samples assayed in each batch varied from six to nine. In the case of samples that were assayed more than once (replicates), aliquots derived from the same mRNA pool were used for probe synthesis. For probe synthesis, aliquots of mRNA corresponding to 4-5 µg of total RNA were mixed together with oligodT_{25NV} (0.5 µg/ml) and mRNA spikes of SpotReport™ 10 array validation system (10 pg; Spike 2, 1 pg), heated to 70°C to remove secondary structures, and then chilled on ice. Probes were prepared in 35µl

- 61 -

reaction mixes by reverse transcription in the presence of 50 μ Ci [α^{33} P] dATP, 3.5 μ M dATP, 0.6 mM each of dCTP, dTTP, dGTP, 200 units of SuperScript reverse transcriptase (Invitrogen, LifeTechnologies) and 0.1 M DTT, labelling for 1.5 hr at 42°C. Following synthesis, the enzyme was deactivated for 10 min. at 70°C and mRNA removed by incubating the reaction mix for 20 min. at 37°C in 4 units of Ribo H (Promega, Madison USA). Unincorporated nucleotides were removed using ProbeQuant G 50 Columns (Amersham Biosciences, Piscataway, USA).

Prior to hybridization, the membranes were equilibrated in 4 x SSC for 2 hr at room temperature and prehybridized overnight at 65°C in 10 ml prehybridisation solution (4 x SSC, 0.1 M NaH_2PO_4 , 1 mM EDTA, 8% dextran sulphate, 10 x denhardt's solution, 1% SDS). Freshly prepared probes were added to 5 ml of the same prehybridisation solution, and hybridization continued overnight at 65°C. The membranes were washed at 65°C at increasing stringency (2 x 30 min. each in 2 x SSC, 0.1% SDS; 1 x SSC, 0.1% SDS; 0.1 x SSC, 0.1% SDS) to remove unspecific signals.

The amount of labelled first strand cDNA binding to each spot was assessed and quantified using a PhosphorImager to generate a gene expression data set. The data was generated using Phoretix software version 3 (Non Linear Dynamics, England). Background subtraction was performed on the generated data by subtracting the median of the line of pixels around each spot outline from the total intensity obtained from the respective spots.

The background-subtracted data was then normalized and transformed by selecting out 50 lowest and 50 maximum signals from each membrane. This step was to exclude genes that were expressed with a high degree of

- 62 -

variance. Since the genes varied from membrane to membrane, the expression data from 497 genes were removed from the data set. The values for the remaining 938 genes were then normalised by using different approaches such as external controls, dividing each spot by the median intensity of the observed signal in the respective membrane, range normalizing the data from each membrane, and then log transforming the data obtained.

10

The processed data obtained above was then used to isolate the informative probes by:

- a) keeping one unique sample (including all repetitions of the selected sample) out per cross validation segment;
- b) building a calibration model (cross validated) on the remaining samples using PLSR-DA;
- c) selecting the set of significant genes for the model in step b using the Jackknife criterion;
- d) repeating steps a), b) and c) until all the unique samples were kept out once (hence, in all 75 different calibration models were built (after repeating step b) 75 times), resulting in 75 different sets of significant probes (after repeating step c) 75 times));
- e) selecting significant variables using the frequency of occurrence criterion amongst the 75 different sets of significant probes.

25

20

30

35

The selected informative probes based on occurrence criterion were used to construct a classification model. The result of the classification model based on probes appearing in at least 90% of the generated sets after the step of isolating informative probes as described above is shown in Figure 2 in which it is seen that the expression pattern of these genes was able to classify most women with breast cancer and women with no breast cancer into distinct groups. In this figure PC1 and PC2

- 63 -

indicate the two principal components statistically derived from the data which best define the systemic variability present in the data. This allows each sample, and the data from each of the informative probes to which the sample's labelled first strand cDNA was bound, to be represented on the classification model as a single point which is a projection of the sample onto the principal components - the score plot.

The ability of the generated model, based on isolated informative probes, to predict future samples was determined by the double cross-validation approach. The performance of the diagnostic test for breast cancer based on the occurrence criterion is presented in Table 6.

Correct prediction of most breast cancer cells was achieved. These included all three samples obtained from women with ductal carcinoma in situ (DCIS), 11/15 samples obtained from women with stage I breast cancer, all five samples obtained from women with stage II breast cancer, and one of two samples obtained from women with stage III breast cancer. Interestingly, two correctly predicted stage I samples were obtained from women having a tumour size of <5 mm in diameter.

The model also correctly predicted the class of most non-cancer samples (41/46), including those that were obtained from women with non-cancerous breast abnormalities.

Confirmation that the gene transcripts are not from cells which are disseminated disease cells has been confirmed by several lines of evidences. Firstly, the informative genes were expressed constitutively at high or moderate levels in blood cells of women irrespective of whether they had cancer or not. Secondly, in the

- 64 -

assay described in this Example, in order to identify transcripts, at least 720 disseminated cells in blood samples would be required. Since, the average number of disseminated cells present in blood during different stages of breast cancer is much lower (organ confined breast cancer, 0.8 cells per ml; invasive breast cancer spread to lymph nodes only, 2.4 cells per ml; and metastatic breast cancer, 6 cells per ml; SD>100%) (29), we believe that the signals being detected originated from peripheral blood cells and could not have originated from disseminated cells. Thirdly, we were not able to detect any signal from the eight cancer markers known to have elevated expression in malignant cancer cells, including cancer cells that are disseminated in the blood.

Example 2: Diagnosis of Alzheimer's disease

Similar experiments were conducted with samples from Alzheimer's patients. In this method 7 patients diagnosed with Alzheimer's Disease at the Memory Clinic at Ullevål University Hospital were used in the trial. The patients were confirmed as having Alzheimer's disease based on the following criteria:

- * A standardized interview with a care-giver using IQCODE, an ADL scale and a scale measuring behaviour of the patient (Green scale).
- * Neuropsychological evaluation using MMSE, Clock drawing test, Trailmaking test A and B (TMT A and B), Kendrick object learning test (visual memory test), part of the Wechsler battery and Benton test.
- * A psychiatric evaluation using scales for detection of depression, MADRS for interviewing the patient and Cornell scale for interviewing the care-giver.
- * A physical examination.
- * Laboratory tests of blood samples to rule out other

- 65 -

diseases.

- * CT scan of the brain.
- * SPECT of the brain.

5 The mean age of the patients was 72.3 with an age range of 69-76. The mean MMSE score was 22.0 (the maximum score attainable being 30).

10 Six age-matched individuals without diagnosed Alzheimer's disease were used as a control. All had been tested with MMSE and had a minimum score of 28 (mean: 28.4). The mean age of the normal control group was 73.0 and the age range 66-81. A sample from a 16-year old individual, with a consequent minimal chance of
15 having Alzheimer's disease, was also included as an additional control.

Using the methods described above (except that hybridization to 758 rather than 1435 cDNA clones was
20 performed), informative probes were selected based on occurrence criterion and used to construct a classification model. The results of the classification model based on probes appearing at least once in the generated sets after the method to isolate informative
25 probes as described above is shown in Figure 3 in which it will be seen that the expression pattern of these genes was able to classify individuals with or without Alzheimer's disease into distinct groups. In this Figure PC1 and PC2 indicate the 2 principal components
30 statistically derived from the data which define the systematic variability present in the data. This allows each sample, and the data from each of the informative probes to which the samples' cDNA was bound, to be represented on the classification model as a single
35 point which is a projection of the sample onto the principal components - the score plot.

- 66 -

The ability of the generated model, based on isolated informative probes, to predict future samples was determined by the double cross-validation. The performance of the diagnostic test for Alzheimer's disease is presented in Table 7.

- 67 -

Appendix A

Partial Least Squares regression (PLSR)

5 Let a multivariate regression model be defined as:

$$Y = XB + F$$

where

10 X a $N \times P$ matrix with N predictor variables (genes);
 Y ($N \times J$) being the J predicted variables. In our case Y
 represents a matrix containing dummy variables;
 B is a matrix of regression coefficients; and
 F is a $N \times J$ matrix of residuals.

15

The structure of the PLSR model can be written as:

$$X = TP^T + E_A, \text{ and}$$

$$Y = TQ^T + F_A, \text{ where}$$

20

where

T ($N \times A$) is a matrix of score vectors which are linear
 combinations of the x-variables;

25 P ($P \times A$) is a matrix with the x-loading vectors p_a as
 columns;

Q ($J \times A$) is a matrix with the y-loading vectors q_a as
 columns;

E_a ($N \times P$) is the matrix for X after A factors; and

F_a ($N \times J$) is the matrix for Y after A factors.

30

The criterion in PLSR is to maximize the explained
 covariance of $[X, Y]$. This is achieved by the loading
 weights vector w_{a+1} , which is the first eigenvector of
 $E_a^T F_a F_a^T E_a$ (E_a and F_a are the deflated X and Y after a
 35 factors or PLS components).

The regression coefficients are given by:

- 68 -

$$B = W(P^T W)^{-1} Q^T$$

A PLSR model with full rank, i.e. maximum number of components, is equivalent to the MLR solutions. Further
5 details on PLSR can be found in Martens & Naes, 1989,
Multivariate Calibration, John Wiley & Sons, Inc., USA
and Kowalski & Seasholtz, 1991, supra.

- 69 -

Example 3: Validation of Example 1, diagnosis of breast cancer

The results in Example 1 were validated by using the
5 informative probes identified in Example 1 on new breast
cancer and control samples.

Methods

The methods, essentially as described in Example 1, were
10 used. Blood was taken from patients as described in
Table 8. However, blood was collected in PAXgene tubes
and the first strand labelled cDNAs were hybridized to
719 cDNAs spotted on nylon membranes along with other
15 controls as described in Example 1. After background
subtraction using control spots, the data of each
membrane was normalized using the inter quantile range.
The data was analysed as described in Example 1 and the
model validated by cross validation.

20 The 719 cDNAs which were spotted are a subset of the
cDNAs spotted in Example 1 and include 111 cDNAs
described in Table 2 and which were found to be
informative in Example 1.

Results

25 The results are shown in Figures 4 to 9. Figures 4, 6
and 8 are projection plots similar to Figure 2 and show
the projection of normal and breast cancer patients'
samples onto a classification model generated using all
30 719 cDNA. Figure 6 is similar but uses a classification
model generated with the 111 probes common to Example 1.
Figure 8 uses the 345 sequences of the 719 for which
sequence information is provided herein. In each case
35 classification of normal and breast cancer groups was
possible. Figures 5, 7 and 9 show prediction plots
which reflect the ability of the generated models to
correctly diagnose breast cancer. In the 3 prediction

- 70 -

plots shown, the disease samples appear on the x axis at +1 and the non-disease samples appear at -1. The y axis represents the predicted class membership. During prediction, if the prediction is correct, disease samples should fall above zero and non-disease samples should fall below zero. In each case almost all samples are correctly predicted.

Example 4: Validation of Example 2, diagnosis of
Alzheimers

The results in Example 2 were validated by using the informative probes identified in Example 2 on new Alzheimer's patient samples.

Methods

The methods, essentially as described in Example 2, were used. Twelve female patients diagnosed with Alzheimer's disease at the Memory Clinic at Ullevål University Hospital who were confirmed as having Alzheimer's disease based on the criteria of Example 2 were used in the trial. The mean age of the patients was 72.3 with an age range of 66-83. The mean MMSE score was 22.0 (the maximum score attainable being 30).

Sixteen age-matched female individuals without diagnosed Alzheimer's disease were used as the normal control group. All had been tested with MMSE and had a minimum score of 29. The mean age of the normal control group was 74.0 and the age range 66-86.

After transfer of the blood to PAXgene tubes, total mRNA was isolated from the blood of the Alzheimer's disease and from the control group donors according to the manufacturers's instructions (PreAnalytiX, Hombrechtikon, Switzerland). The isolated mRNA was labelled during reverse transcription in the presence of

- 71 -

$\alpha^{33}\text{P}$ -dATP, yielding a labelled first strand cDNA.

Hybridization was performed as described previously onto
730 cDNA clones picked from a cDNA library from whole
blood of 550 healthy individuals without knowledge of
5 the gene sequence of the random cDNA clones.

Results

The results are shown in Figures 10 and 11. Figure 10
is a projection plot generated using 520 probes which
10 have been sequenced. Figure 11 is a prediction plot and
shows correct prediction of almost all samples.

Table 1a

List of probes informative for disease diagnosis

	Clone ID	Sequence ID	No. of nucleotides
1	I-01	-	-
2	I-02	-	-
3	I-13	-	-
4	I-21	-	-
5	I-24	308	373
6	I-28	310	584
7	I-30	1180	622
8	I-34	313	554
9	I-37	-	-
10	I-42	-	-
11	I-52	-	-
12	I-54	1181	155
13	I-58	326	554
14	I-71	-	-
15	I-72	-	-
16	I-86	-	-
17	I-95	-	-
18	II-03	361	622
19	II-05	363	628
20	II-06	364	528
21	II-10	368	329
22	II-24	381	534
23	II-25	382	444
24	II-26	383	566
25	II-33	390	523
26	II-34	391	566
27	II-41	397	534
28	II-42	398	512
29	II-47	-	-
30	II-57	411	505
31	II-61	415	596
32	II-69	423	387
33	II-70	424	420
34	II-75	429	535
35	II-83	-	-
36	II-84	438	577
37	II-87	441	552
38	II-88	442	606
39	II-90	-	-
40	II-94	448	329
41	III-02	453	747
42	III-05	-	-
43	III-06	458	682

44	III-08	460	536
45	III-10	-	-
46	III-13	464	615
47	III-15	-	-
48	III-17	-	-
49	III-20	1183	479
50	III-23	473	694
51	III-26	476	476
52	III-35	485	551
53	III-39	487	224
54	III-40	488	349
55	III-43	490	382
56	III-44	491	382
57	III-53	500	390
58	III-56	503	109
59	III-57	504	374
60	III-60	-	-
61	III-60	-	-
62	III-61	507	521
63	III-63	509	575
64	III-68	-	-
65	III-74	518	502
66	III-80	523	585
67	III-82	-	-
68	III-85	526	516
69	III-89	530	660
70	III-92	-	-
71	III-96	-	-
72	IV-14	684	545
73	IV-15	1185	628
74	IV-23	-	-
76	IV-26	1186	494
75	IV-26	-	-
77	IV-29	-	-
78	IV-31	687	268
79	IV-32	688	569
80	IV-34	-	-
81	IV-35	-	-
82	IV-41	-	-
83	IV-45	-	-
84	IV-53	61	362
85	IV-62	-	-
86	IV-69	192	286
87	IV-80	701	579
88	IV-82	-	-
89	IV-93	-	-
90	IX-10	736	641
91	IX-12	-	-
92	IX-38	757	583
93	IX-39	758	424
94	IX-42	-	-
95	IX-48	764	626
96	IX-77	785	556
97	V-01	-	-
98	V-02	-	-
99	V-03	706	486

100	V-04	707	397
101	V-06	-	-
102	V-07	708	293
103	V-11	1188	599
104	V-12	711	498
105	V-15	-	-
106	V-17	-	-
107	V-21	-	-
108	V-25	-	-
109	V-32	-	-
110	V-35	-	-
111	V-39	-	-
112	V-42	-	-
113	V-43	-	-
114	V-47	-	-
115	V-49	-	-
116	V-52	-	-
117	V-54	-	-
118	V-55	77	412
119	V-58	-	-
120	V-59	-	-
121	V-65	-	-
122	V-68	-	-
123	V-71	-	-
124	V-75	-	-
125	V-79	-	-
126	V-80	726	260
127	V-90	-	-
128	V-91	-	-
129	V-92	-	-
130	V-94	-	-
131	VI-02	-	-
132	VI-04	865	122
133	VI-07	93	405
134	VI-09	-	-
135	VI-10	-	-
136	VI-12	869	667
137	VI-14	871	642
138	VI-17	-	-
139	VI-20	876	115
140	VI-21	-	-
141	VI-23	878	634
142	VI-34	-	-
143	VI-41	-	-
144	VI-42	-	-
145	VI-43	-	-
146	VI-44	-	-
147	VI-48	891	626
148	VI-49	-	-
149	VI-50	893	585
150	VI-52	-	-
151	VI-53	895	560
152	VI-55	897	509
153	VI-65	-	-
154	VI-70	108	550
155	VI-71	-	-

156	VI-72	-	-
157	VI-74	905	655
158	VI-76	907	582
159	VI-78	-	-
160	VI-79	-	-
161	VI-84	-	-
162	VI-87	911	595
163	VI-88	912	651
164	VI-90	-	-
165	VI-93	-	-
166	VI-95	915	230
167	VI-96	-	-
168	VII-02	-	-
169	VII-03	1196	412
170	VII-06	-	-
171	VII-10	-	-
172	VII-11	-	-
173	VII-15	1199	439
174	VII-19	562	580
175	VII-21	564	671
176	VII-25	-	-
177	VII-32	571	457
178	VII-36	575	209
179	VII-39	576	541
180	VII-42	579	502
181	VII-43	580	316
182	VII-46	583	631
183	VII-47	1200	526
184	VII-48	1201	613
185	VII-59	593	565
186	VII-60	-	-
187	VII-63	595	98
188	VII-66	598	362
189	VII-67	-	-
190	VII-72	600	595
191	VII-73	601	522
192	VII-75	-	-
193	VII-76	603	624
194	VII-77	1203	692
195	VII-80	605	338
196	VII-81	606	556
197	VII-83	-	-
198	VII-86	-	-
199	VII-88	-	-
200	VII-90	612	576
201	VII-91	613	341
202	VII-93	615	379
203	VIII-01	-	-
204	VIII-02	-	-
205	VIII-03	-	-
206	VIII-06	-	-
207	VIII-09	618	598
208	VIII-10	-	-
209	VIII-15	-	-
210	VIII-20	628	419
211	VIII-22	-	-

212	VIII-26	-	-
213	VIII-28	634	511
214	VIII-29	635	592
215	VIII-30	636	572
216	VIII-31	637	482
217	VIII-32	638	545
218	VIII-33	639	624
219	VIII-39	-	-
220	VIII-41	645	649
221	VIII-42	646	600
222	VIII-44	-	-
223	VIII-46	649	425
224	VIII-48	651	251
225	VIII-58	-	-
226	VIII-64	663	627
227	VIII-65	-	-
228	VIII-66	665	345
229	VIII-67	666	252
230	VIII-74	-	-
231	VIII-76	675	591
232	VIII-78	-	-
233	VIII-82	-	-
234	VIII-83	-	-
235	VIII-85	-	-
236	VIII-87	-	-
237	VIII-91	-	-
238	VIII-92	-	-
239	VIII-93	-	-
240	VIII-95	-	-
241	X-04	-	-
242	X-07	808	641
243	X-15	814	132
244	X-29	821	370
245	X-34	-	-
246	X-35	-	-
247	X-54	837	603
248	X-56	839	71
249	X-68	1207	642
250	X-72	849	622
251	X-94	860	501
252	XI-07	-	-
253	XI-13	1209	620
254	XI-50	-	-
255	XI-58	-	-
256	XI-81	1212	374
257	XII-07	1213	567
258	XII-17	-	-
259	XII-26	-	-
260	XII-27	-	-
261	XII-31	-	-
262	XII-32	-	-
263	XII-35	1214	620
264	XII-36	-	-
265	XII-52	-	-
266	XII-59	1216	484
267	XIII-19	1219	559

268	XIII-29	-	-
269	XIII-52	939	513
270	XIII-62	-	-
271	XIII-84	-	-
272	XIII-92	1221	741
273	XV-18	-	-
274	XV-22	1099	561
275	XV-24	-	-
276	XV-25	1224	485
277	XV-28	-	-
278	XV-34	-	-
279	XV-42	-	-
280	XV-68	-	-
281	XV-74	-	-
282	XV-93	-	-
283	XV-94	-	-
284	XV-96	-	-
285	XVI-36	1056	435
286	XVI-53	1230	741
287	XVI-59	-	-
288	XVI-66	1074	689
289	XVI-76	1083	198
290	XVI-77	1084	198
291	XVII-07	-	-
292	XVII-08	-	-
293	XVII-17	-	-
294	XVII-28	-	-
295	XVII-29	-	-
296	XVII-31	1139	503
297	XVII-36	-	-
298	XVII-39	-	-
299	XVII-40	1231	203
300	XVII-48	1148	587
301	XVII-55	-	-
302	XVII-58	-	-
303	XVII-67	-	-
304	XVII-72	-	-
305	XVII-76	1160	650
306	XVII-82	-	-
307	XVII-87	1165	502
308	XVII-95	1172	648

Table 1 b**List of sequences of probes informative for disease diagnosis**

Please see the note at the bottom

Clone ID	Sequence ID
I-09	298
I-10	299
I-13	1331
I-14	1178
I-15	300
I-16	301
I-17	302
I-19	304
I-20	305
I-22	306
I-23	307
I-24	308
I-25	309
I-28	310
I-30	1180
I-31	311
I-32	312
I-34	313
I-37	1440
I-38	314
I-39	315
I-40	316
I-42	1332
I-44	317
I-45	318
I-46	319
I-47	320
I-48	321
I-49	322
I-53	323
I-54	1181
I-56	324
I-57	325
I-58	326
I-60	327
I-64	328
I-67	330
I-69	331
I-71	332
I-72	333
I-73	334
I-77	335
I-79	336
I-80	337

I-81	338
I-82	339
I-86	1336
I-88	1182
I-95	1337
II-02	360
II-03	361
II-05	363
II-06	364
II-07	365
II-08	366
II-09	367
II-10	368
II-11	369
II-12	370
II-13	371
II-14	372
II-15	373
II-16	374
II-17	375
II-18	376
II-20	377
II-21	378
II-22	379
II-23	380
II-24	381
II-25	382
II-26	383
II-27	384
II-28	385
II-29	386
II-30	387
II-31	388
II-32	389
II-33	390
II-34	391
II-35	392
II-37	393
II-38	394
II-39	395
II-40	396
II-41	397
II-42	398
II-43	399
II-44	400
II-46	401
II-47	402
II-48	403
II-49	404
II-50	405
II-52	406

II-53	407
II-54	408
II-55	409
II-56	410
II-57	411
II-58	412
II-59	413
II-60	414
II-61	415
II-62	416
II-63	417
II-64	418
II-65	419
II-66	420
II-67	421
II-68	422
II-69	423
II-70	424
II-71	425
II-72	426
II-73	427
II-74	428
II-75	429
II-76	430
II-77	431
II-78	432
II-79	433
II-80	434
II-81	435
II-82	436
II-83	437
II-84	438
II-85	439
II-86	440
II-87	441
II-88	442
II-89	443
II-90	444
II-91	445
II-92	446
II-93	447
II-94	448
II-95	449
II-96	450
III-01	452
III-02	453
III-03	454
III-04	455
III-05	457
III-06	458
III-07	459

III-08	460
III-09	461
III-11	462
III-12	463
III-13	464
III-14	465
III-15	466
III-16	467
III-17	468
III-18	469
III-19	470
III-20	1183
III-21	471
III-22	472
III-23	473
III-24	474
III-25	475
III-26	476
III-27	477
III-28	478
III-29	479
III-31	481
III-32	482
III-33	483
III-34	484
III-35	485
III-37	486
III-39	487
III-40	488
III-42	489
III-43	490
III-44	491
III-45	492
III-46	493
III-47	494
III-48	495
III-49	496
III-50	497
III-51	498
III-52	499
III-53	500
III-54	501
III-55	502
III-56	503
III-57	504
III-58	505
III-59	506
III-61	507
III-62	508
III-63	509
III-64	510

III-65	511
III-66	512
III-67	513
III-69	514
III-70	515
III-71	516
III-73	517
III-74	518
III-75	519
III-77	520
III-78	521
III-79	522
III-80	523
III-81	524
III-82	1348
III-83	525
III-85	526
III-86	527
III-87	528
III-88	529
III-89	530
III-91	531
III-92	1351
III-93	532
III-94	533
III-95	534
III-96	535
IV-02	681
IV-04	682
IV-13	683
IV-14	684
IV-15	1185
IV-17	685
IV-23	1353
IV-26	1186
IV-28	686
IV-31	687
IV-32	688
IV-35	1355
IV-37	g6
IV-38	689
IV-40	690
IV-42	691
IV-43	1239
IV-44	692
IV-47	693
IV-53	61
IV-55	694
IV-56	695
IV-61	696
IV-64	697

IV-65	698
IV-69	192
IV-72	699
IV-73	700
IV-80	701
IV-82	196
IV-85	702
IV-93	703
IV-95	704
IV-96	705
IX-10	736
IX-12	738
IX-13	739
IX-24	747
IX-38	757
IX-39	758
IX-48	764
IX-50	766
IX-56	768
IX-62	773
IX-65	776
IX-72	782
IX-77	785
IX-91	796
IX-96	801
V-01	1361
V-03	706
V-04	707
V-07	708
V-08	709
V-09	710
V-11	1188
V1-16	873
V1-19	875
V-12	711
V-17	1364
V-18	712
V-20	713
V-24	714
V-25	1365
V-28	1189
V-35	1366
V-37	716
V-38	1190
V-39	1109
V-40	717
V-41	718
V-47	1368
V-48	719
V-49	1369
V-55	77

V-57	720
V-58	1370
V-61	721
V-64	722
V-65	723
V-68	1448
V-71	1495
V-74	724
V-75	1372
V-80	726
V-81	727
V-87	728
V-90	1374
VI-02	340
VI-03	341
VI-04	342
VI-06	343
VI-07	344
VI-08	345
VI-09	346
VI-11	347
VI-12	869
VI-13	870
VI-14	871
VI-16	873
VI-18	348
VI-19	349
VI-20	350
VI-21	351
VI-22	352
VI-23	878
VI-24	879
VI-25	353
VI-26	354
VI-27	355
VI-31	356
VI-32	885
VI-33	357
VI-35	358
VI-39	887
VI-43	1382
VI-44	1193
VI-45	889
VI-48	359
VI-49	892
VI-50	893
VI-53	895
VI-55	897
VI-58	899
VI-66	903
VI-67	904

VI-70	108
VI-71	1387
VI-74	905
VI-75	906
VI-76	907
VI-77	110
VI-79	1389
VI-80	908
VI-85	910
VI-87	911
VI-88	912
VI-90	1390
VI-93	1391
VI-95	915
VI-96	1392
VII-02	547
VII-03	548
VII-04	549
VII-05	550
VII-06	551
VII-07	552
VII-08	553
VII-09	554
VII-10	555
VII-11	556
VII-12	557
VII-14	558
VII-15	559
VII-17	560
VII-18	561
VII-19	562
VII-20	563
VII-21	564
VII-22	565
VII-23	566
VII-24	567
VII-25	1397
VII-26	250
VII-27	568
VII-28	569
VII-29	570
VII-32	571
VII-33	572
VII-34	573
VII-35	574
VII-36	575
VII-39	576
VII-40	577
VII-41	578
VII-42	579
VII-43	580

VII-44	581
VII-45	582
VII-46	583
VII-47	1200
VII-48	584
VII-49	585
VII-50	586
VII-52	587
VII-53	588
VII-54	589
VII-55	590
VII-57	591
VII-58	592
VII-59	593
VII-62	594
VII-63	595
VII-64	596
VII-65	597
VII-66	598
VII-67	1399
VII-71	599
VII-72	600
VII-73	601
VII-74	602
VII-76	603
VII-77	604
VII-80	605
VII-81	606
VII-82	607
VII-83	608
VII-84	609
VII-86	1453
VII-87	610
VII-89	611
VII-90	612
VII-91	613
VII-92	614
VII-93	615
VII-94	616
VII-96	617
VIII-09	618
VIII-10	619
VIII-11	620
VIII-12	621
VIII-13	622
VIII-15	623
VIII-16	624
VIII-17	625
VIII-18	626
VIII-19	627
VIII-20	628

VIII-21	629
VIII-22	1455
VIII-23	630
VIII-24	631
VIII-25	632
VIII-26	1456
VIII-27	633
VIII-28	634
VIII-29	635
VIII-30	636
VIII-31	637
VIII-32	638
VIII-33	639
VIII-34	640
VIII-36	641
VIII-37	642
VIII-38	643
VIII-40	644
VIII-41	645
VIII-42	646
VIII-43	647
VIII-45	648
VIII-46	649
VIII-47	650
VIII-48	651
VIII-50	652
VIII-51	653
VIII-53	654
VIII-54	655
VIII-55	656
VIII-56	657
VIII-57	658
VIII-58	659
VIII-59	660
VIII-60	661
VIII-61	662
VIII-64	663
VIII-65	664
VIII-66	665
VIII-67	666
VIII-68	667
VIII-69	668
VIII-70	669
VIII-71	670
VIII-72	671
VIII-73	672
VIII-74	673
VIII-75	674
VIII-76	675
VIII-77	676
VIII-78	677

VIII-79	678
VIII-80	679
X-07	808
X-15	814
X-20	817
X-29	821
X-34	825
X-46	833
X-54	837
X-56	839
X-68	1207
X-72	849
X-73	1208
X-94	860
XI-13	1209
XI-37	1460
XI-43	1210
XI-67	1211
XI-81	1212
XII-07	1213
XII-35	1214
XII-36	1215
XII-59	1216
XII-65	1028
XII-92	1217
XIII-03	917
XIII-04	1218
XIII-19	1219
XIII-24	926
XIII-51	938
XIII-52	939
XIII-67	947
XIII-69	949
XIII-88	1220
XIII-92	1221
XV-22	1099
XV-24	1101
XV-25	1224
XV-42	1108
XV-62	1226
XV-64	1118
XV-84	1125
XVI-19	1228
XVI-36	1056
XVI-53	1230
XVI-60	1071
XVI-66	1074
XVI-74	1081
XVI-76	1083
XVI-77	1084
XVII-31	1139

XVII-40	1231
XVII-48	1148
XVII-76	1160
XVII-87	1165
XVII-95	1172

Note

Sequences not available for sequence IDs in Table 1, and corresponding sequence IDs in Table 2 and 4.

298,301,305,307,312,317,318,319,320,332,333,334,336,340,341,342,343,344,345,346,347,348,349,350,351,352,353,354,355,356,357,358,359,367,372,375,376,377,379,385,392,393,404,437,439,440,443,444,445,449,455,457,465,466,467,468,470,486,498,501,511,514,516,517,520,522,528,531,535,547,548,549,550,551,552,553,554,555,556,557,558,559,573,584,604,608,616,620,623,640,659,662,664,667,668,673,677,678,679,681,695,702,712,716,825,886,894,902,909,916,1101,1108,1109,1177,1187,1193,1204,1220,1239,1255,1256,1342,1347,1354,1357,1362,1363,1364,1373,1375,1379,1403,1404,1405,1406,1413

Table 2a

List of informative probes for diagnosis of breast cancer

Clone ID	Sequence ID
I-24	308
I-28	310
I-30	1180
I-52	-
I-54	1181
II-41	397
II-70	424
II-87	441
III-06	458
III-20	1183
III-40	488
III-57	504
III-60	-
III-61	507
III-89	530
IV-14	684
IV-15	1185
IV-26	1186
IV-32	688
IV-41	-
IV-53	61
IV-62	-
IV-69	192
IV-80	701
IV-82	196
IX-10	736
IX-12	-
IX-38	757
IX-39	758
IX-42	-
IX-48	764
IX-77	785
V-11	1188
V-32	-
V-39	-
V-55	77
V-80	726
V-94	-
VI-07	93
VI-34	-
VI-41	-
VI-48	891
VI-49	-
VI-52	-
VI-55	897
VI-65	-
VI-70	108

Clone ID	Sequence ID
VI-72	-
VI-78	-
VI-84	-
VII-03	1196
VII-15	1199
VII-32	571
VII-39	576
VII-47	1200
VII-48	1201
VII-60	-
VII-73	601
VII-77	1203
VII-90	612
VIII-20	628
VIII-29	635
VIII-30	636
VIII-31	637
VIII-39	-
VIII-44	-
VIII-46	649
VIII-48	651
VIII-66	665
VIII-74	-
VIII-76	675
X-04	-
X-07	808
X-15	814
X-29	821
X-34	-
X-35	-
X-54	837
X-56	839
X-68	1207
X-72	849
X-94	860
XI-07	-
XI-13	1209
XI-50	-
XI-58	-
XI-81	1212
XII-07	1213
XII-17	-
XII-26	-
XII-27	-
XII-31	-
XII-32	-
XII-35	1214

Clone ID	Sequence ID
XII-36	-
XII-52	-
XII-59	1216
XIII-19	1219
XIII-29	-
XIII-52	939
XIII-62	-
XIII-84	-
XIII-92	1221
XV-18	-
XV-22	1099
XV-24	-
XV-25	1224
XV-28	-
XV-34	-
XV-42	-
XV-68	-
XV-74	-
XV-93	-
XV-94	-
XV-96	-
XVI-36	1056
XVI-53	1230
XVI-59	-
XVI-66	1074
XVI-76	1083
XVI-77	1084
XVII-07	-
XVII-08	-
XVII-17	-
XVII-28	-
XVII-29	-
XVII-31	1139
XVII-36	-
XVII-39	-
XVII-40	1231
XVII-48	1148
XVII-55	-
XVII-58	-
XVII-67	-
XVII-72	-
XVII-76	1160
XVII-82	-
XVII-87	1165
XVII-95	1172

Table 2b**List of sequences of probes informative for breast cancer**

Please see the note at the bottom of Table 1. Some sequences are missing.

Clone ID	Sequence ID
I-13	1331
I-14	1178
I-24	308
I-25	309
I-28	310
I-30	1180
I-37	1440
I-42	1332
I-48	321
I-54	1181
I-60	327
I-72	1335
I-81	338
I-82	339
I-86	1336
I-88	1182
I-95	1337
II-02	360
II-03	361
II-06	364
II-07	365
II-10	368
II-21	378
II-23	380
II-24	381
II-25	382
II-27	384
II-33	390
II-34	391
II-41	397
II-42	398
II-46	401
II-47	1338
II-48	403
II-52	406
II-57	411
II-58	412
II-59	413
II-60	414
II-61	415
II-62	416
II-64	418

II-67	421
II-69	423
II-70	424
II-74	428
II-80	434
II-82	436
II-84	438
II-87	441
II-88	442
II-96	450
III-01	452
III-02	453
III-06	458
III-08	460
III-12	463
III-13	464
III-17	1344
III-18	469
III-20	1183
III-21	471
III-23	473
III-24	474
III-25	475
III-26	476
III-27	477
III-28	478
III-29	479
III-32	482
III-33	483
III-35	485
III-39	487
III-40	488
III-42	489
III-45	492
III-46	493
III-47	494
III-48	495
III-56	503
III-57	504
III-58	505
III-59	506
III-61	507
III-62	508
III-63	509
III-64	510
III-66	512
III-67	513
III-70	515
III-74	518
III-75	519
III-78	521

III-80	523
III-81	524
III-82	1348
III-85	526
III-86	527
III-88	529
III-89	530
III-92	1351
III-93	532
III-95	534
III-96	1352
IV-04	682
IV-13	683
IV-14	684
IV-15	1185
IV-17	685
IV-23	1353
IV-26	1186
IV-31	687
IV-32	688
IV-35	1355
IV-37	96
IV-38	689
IV-42	691
IV-43	1239
IV-47	693
IV-53	61
IV-61	696
IV-64	697
IV-69	192
IV-72	699
IV-80	701
IV-82	196
IV-85	702
IV-93	1360
IV-96	705
IX-10	736
IX-12	738
IX-13	739
IX-24	747
IX-38	757
IX-39	758
IX-48	764
IX-50	766
IX-56	768
IX-62	773
IX-65	776
IX-72	782
IX-77	785
IX-91	796
IX-96	801

V-01	1361
V-03	706
V-04	707
V-07	708
V-08	709
V-11	1188
V-12	711
V-17	1364
V-24	714
V-25	1365
V-28	1189
V-35	1366
V-38	1190
V-39	1109
V-41	718
V-47	1368
V-49	1369
V-55	77
V-57	720
V-58	1370
V-61	721
V-64	722
V-65	1371
V-68	1448
V-71	1495
V-74	724
V-75	1372
V-80	726
V-90	1374
VI-03	864
VI-04	865
VI-07	93
VI-08	867
VI-09	1378
VI-12	869
VI-13	870
VI-14	871
VI-16	873
VI-19	875
VI-20	876
VI-21	1380
VI-23	878
VI-24	879
VI-25	1192
VI-26	881
VI-32	885
VI-39	887
VI-43	1382
VI-44	1193
VI-45	889
VI-48	891

VI-49	892
VI-50	893
VI-53	895
VI-55	897
VI-58	899
VI-66	903
VI-67	904
VI-70	108
VI-71	1387
VI-74	905
VI-75	906
VI-76	907
VI-77	110
VI-79	1389
VI-80	908
VI-85	910
VI-87	911
VI-88	912
VI-90	1390
VI-93	1391
VI-95	915
VI-96	1392
VII-02	1195
VII-03	1196
VII-06	1394
VII-08	1197
VII-09	1198
VII-10	1395
VII-11	1396
VII-15	1199
VII-17	560
VII-19	562
VII-21	564
VII-22	565
VII-23	566
VII-24	567
VII-25	1397
VII-26	250
VII-27	568
VII-29	570
VII-32	571
VII-33	572
VII-36	575
VII-39	576
VII-41	578
VII-42	579
VII-43	580
VII-46	583
VII-47	1200
VII-48	1201
VII-49	585

VII-54	589
VII-57	591
VII-58	592
VII-59	593
VII-62	594
VII-63	1202
VII-64	596
VII-66	598
VII-67	1399
VII-72	600
VII-73	601
VII-77	1203
VII-80	605
VII-82	607
VII-86	1453
VII-87	610
VII-90	612
VII-91	613
VII-92	614
VII-93	615
VII-96	617
VIII-09	618
VIII-10	619
VIII-13	622
VIII-16	624
VIII-20	628
VIII-21	629
VIII-22	1455
VIII-23	630
VIII-24	631
VIII-25	632
VIII-26	1456
VIII-27	633
VIII-28	634
VIII-29	635
VIII-30	636
VIII-31	637
VIII-32	638
VIII-33	639
VIII-34	1204
VIII-38	643
VIII-40	644
VIII-41	645
VIII-46	649
VIII-48	651
VIII-55	656
VIII-57	658
VIII-59	660
VIII-60	661
VIII-61	1205
VIII-64	663

VIII-66	665
VIII-73	672
VIII-74	673
VIII-76	675
VIII-80	679
X-07	808
X-15	814
X-20	817
X-29	821
X-34	825
X-46	833
X-54	837
X-56	839
X-68	1207
X-72	849
X-73	1208
X-94	860
XI-13	1209
XI-37	1460
XI-43	1210
XI-67	1211
XI-81	1212
XII-07	1213
XII-35	1214
XII-36	1215
XII-59	1216
XII-65	1028
XII-92	1217
XIII-03	917
XIII-04	1218
XIII-19	1219
XIII-24	926
XIII-51	938
XIII-52	939
XIII-67	947
XIII-69	949
XIII-88	1220
XIII-92	1221
XV-22	1099
XV-24	1101
XV-25	1224
XV-42	1108
XV-62	1226
XV-64	1118
XV-84	1125
XVI-19	1228
XVI-36	1056
XVI-53	1230
XVI-60	1071
XVI-66	1074
XVI-74	1081

XVI-76	1083
XVI-77	1084
XVII-31	1139
XVII-40	1231
XVII-48	1148
XVII-76	1160
XVII-87	1165
XVII-95	1172

Table 3

List of informative probes (Clone ID) selected for breast cancer diagnosis based on their occurrence criterion during variable selection.

Occurrence*	Clone ID
100%	XI-8, XVI-66, VIII-66, XVI-59, VII-03, XIII-19, XII-35, X-35, XI-50, XII-26, IV-53, XIII-29, XIII-62, I-30, III-06, XV-22, XV-94, VII-15, VII-39, IX-39, XVII-39, III-40, VII-32
90%	I-52, VI-65, VI-34, IV-62, XV-34, XVII-58, V-11, VI-78, XII-36, XIII-92, VIII-29, XVI-53, XVI-77, XI-13, XIII-84, IV-14, XII-31, V-80, VII-48, XVII-29, XVII-72
80%	III-60, VIII-74, IX-12, X-04, XIII-52, VIII-30, IX-38
70%	VI-49, X-29, VIII-48
60%	IV-82, IX-10, VI-52, X-68, VII-77
50%	IV-15
40%	XV-28, II-70, V-55
30%	XVII-17, XVII-67
20%	XI-58, XVI-36, VIII-39, VIII-44, III-61, IV-69, XV-68, X-72
10%	IX-42, IX-77, X-94, XV-96, XVII-55
5%	XII-59, XVI-76, I-54, XV-18, V-94, X-54, VI-07, VII-47, XVII-31, XVII-87, XVII-48
In at least one model	II-41, VI-41, III-57, III-89, VII-73, XV-25, IV-26, X-34, IV-41, VII-90, XV-42, XVII-82, XII-27, VIII-20, I-28, VII-60, VIII-76, III-20, VI-84, XI-07, XVII-28, XII-17, XVII-36, XII-52, XVII-76, VIII-46, VI-70, XV-74, XV-93, VIII-31, II-87, V-39, VI-55, X-07, X-15, XII-07, XVII-07, XVII-08, XVII-95, I-24, IV-32, V-32, VI-48, VI-72, IV-80, IX-48, X-56, XV-24, XII-32, XVII-40

*100% = Genes appearing in all the 75 cross validated models; 90% = Additional genes appearing in at least 68 out of 75 cross validated models; 5% = Additional genes appearing in at least 4 out of 75 cross validated models and so on.

Table 4a

List of informative probes for diagnosis of Alzheimer disease

Clone ID	Sequence ID
I-01	-
I-02	-
I-13	-
I-21	-
I-34	313
I-37	-
I-42	-
I-58	326
I-71	-
I-72	-
I-86	-
I-95	-
II-03	361
II-05	363
II-06	364
II-10	368
II-24	381
II-25	382
II-26	383
II-33	390
II-34	391
II-42	398
II-47	-
II-57	411
II-61	415
II-69	423
II-75	429
II-83	-
II-84	438
II-88	442
II-90	-
II-94	449
III-02	453
III-05	-
III-06	458
III-08	460
III-10	-
III-13	464
III-15	-
III-17	-
III-23	473
III-26	476
III-35	485
III-39	487
III-43	490
III-44	491
III-53	500
III-56	503

Clone ID	Sequence ID
III-60	-
III-63	509
III-68	-
III-74	518
III-80	523
III-82	-
III-85	526
III-92	-
III-96	-
IV-23	-
IV-26	-
IV-29	-
IV-31	687
IV-34	-
IV-35	-
IV-45	-
IV-80	701
IV-82	-
IV-93	-
V-01	-
V-02	-
V-03	706
V-04	707
V-06	-
V-07	708
V-12	711
V-15	-
V-17	-
V-21	-
V-25	-
V-35	-
V-42	-
V-43	-
V-47	-
V-49	-
V-52	-
V-54	-
V-58	-
V-59	-
V-65	-
V-68	-
V-71	-
V-75	-
V-79	-
V-80	726
V-90	-
V-91	-
V-92	-

Clone ID	Sequence ID
VI-02	-
VI-04	865
VI-09	-
VI-10	-
VI-12	869
VI-14	871
VI-17	-
VI-20	876
VI-21	-
VI-23	878
VI-41	-
VI-42	-
VI-43	-
VI-44	-
VI-48	891
VI-49	-
VI-50	893
VI-53	895
VI-71	-
VI-74	905
VI-76	907
VI-78	-
VI-79	-
VI-87	911
VI-88	912
VI-90	-
VI-93	-
VI-95	915
VI-96	-
VII-02	-
VII-03	-
VII-06	-
VII-10	-
VII-11	-
VII-19	562
VII-21	564
VII-25	-
VII-36	575
VII-42	579
VII-43	580
VII-46	583
VII-59	593
VII-63	595
VII-66	598
VII-67	-
VII-72	600
VII-73	601
VII-75	-
VI-02	-
VI-04	866
VI-09	-
VI-10	-
VI-12	873
VI-14	875
VI-17	-

Clone ID	Sequence ID
VII-91	613
VII-93	615
VIII-01	-
VIII-02	-
VIII-03	-
VIII-06	-
VIII-09	618
VIII-10	-
VIII-15	-
VIII-22	-
VIII-26	-
VIII-28	634
VIII-30	636
VIII-32	638
VIII-33	639
VIII-41	645
VIII-42	646
VIII-48	651
VIII-58	-
VIII-64	663
VIII-65	-
VIII-67	666
VIII-78	-
VIII-82	-
VIII-83	-
VIII-85	-
VIII-87	-
VIII-91	-
VIII-92	-
VIII-93	-
VIII-95	-

Table 4b**List of sequences of probes informative for Alzheimer disease**

Please see note to Table 1

Clone ID	Sequence ID
I-09	298
I-10	299
I-15	300
I-16	301
I-17	302
I-19	304
I-20	305
I-22	306
I-23	307
I-24	308
I-25	309
I-28	310
I-31	311
I-32	312
I-34	313
I-38	314
I-39	315
I-40	316
I-44	317
I-45	318
I-46	319
I-47	320
I-48	321
I-49	322
I-53	323
I-56	324
I-57	325
I-58	326
I-60	327
I-64	328
I-67	330
I-69	331
I-71	332
I-72	333
I-73	334
I-77	335
I-79	336
I-80	337
I-81	338
I-82	339
VI-02	340

VI-03	341
VI-04	342
VI-06	343
VI-07	344
VI-08	345
VI-09	346
VI-11	347
VI-18	348
VI-19	349
VI-20	350
VI-21	351
VI-22	352
VI-25	353
VI-26	354
VI-27	355
VI-31	356
VI-33	357
VI-35	358
VI-48	359
II-02	360
II-03	361
II-05	363
II-06	364
II-07	365
II-08	366
II-09	367
II-10	368
II-11	369
II-12	370
II-13	371
II-14	372
II-15	373
II-16	374
II-17	375
II-18	376
II-20	377
II-21	378
II-22	379
II-23	380
II-24	381
II-25	382
II-26	383
II-27	384
II-28	385
II-29	386
II-30	387
II-31	388
II-32	389
II-33	390
II-34	391
II-35	392

II-37	393
II-38	394
II-39	395
II-40	396
II-41	397
II-42	398
II-43	399
II-44	400
II-46	401
II-47	402
II-48	403
II-49	404
II-50	405
II-52	406
II-53	407
II-54	408
II-55	409
II-56	410
II-57	411
II-58	412
II-59	413
II-60	414
II-61	415
II-62	416
II-63	417
II-64	418
II-65	419
II-66	420
II-67	421
II-68	422
II-69	423
II-70	424
II-71	425
II-72	426
II-73	427
II-74	428
II-75	429
II-76	430
II-77	431
II-78	432
II-79	433
II-80	434
II-81	435
II-82	436
II-83	437
II-84	438
II-85	439
II-86	440
II-87	441
II-88	442
II-89	443

II-90	444
II-91	445
II-92	446
II-93	447
II-94	448
II-95	449
II-96	450
III-01	452
III-02	453
III-03	454
III-04	455
III-05	457
III-06	458
III-07	459
III-08	460
III-09	461
III-11	462
III-12	463
III-13	464
III-14	465
III-15	466
III-16	467
III-17	468
III-18	469
III-19	470
III-21	471
III-22	472
III-23	473
III-24	474
III-25	475
III-26	476
III-27	477
III-28	478
III-29	479
III-31	481
III-32	482
III-33	483
III-34	484
III-35	485
III-37	486
III-39	487
III-40	488
III-42	489
III-43	490
III-44	491
III-45	492
III-46	493
III-47	494
III-48	495
III-49	496
III-50	497

III-51	498
III-52	499
III-53	500
III-54	501
III-55	502
III-56	503
III-57	504
III-58	505
III-59	506
III-61	507
III-62	508
III-63	509
III-64	510
III-65	511
III-66	512
III-67	513
III-69	514
III-70	515
III-71	516
III-73	517
III-74	518
III-75	519
III-77	520
III-78	521
III-79	522
III-80	523
III-81	524
III-83	525
III-85	526
III-86	527
III-87	528
III-88	529
III-89	530
III-91	531
III-93	532
III-94	533
III-95	534
III-96	535
VII-02	547
VII-03	548
VII-04	549
VII-05	550
VII-06	551
VII-07	552
VII-08	553
VII-09	554
VII-10	555
VII-11	556
VII-12	557
VII-14	558
VII-15	559

VII-17	560
VII-18	561
VII-19	562
VII-20	563
VII-21	564
VII-22	565
VII-23	566
VII-24	567
VII-27	568
VII-28	569
VII-29	570
VII-32	571
VII-33	572
VII-34	573
VII-35	574
VII-36	575
VII-39	576
VII-40	577
VII-41	578
VII-42	579
VII-43	580
VII-44	581
VII-45	582
VII-46	583
VII-48	584
VII-49	585
VII-50	586
VII-52	587
VII-53	588
VII-54	589
VII-55	590
VII-57	591
VII-58	592
VII-59	593
VII-62	594
VII-63	595
VII-64	596
VII-65	597
VII-66	598
VII-71	599
VII-72	600
VII-73	601
VII-74	602
VII-76	603
VII-77	604
VII-80	605
VII-81	606
VII-82	607
VII-83	608
VII-84	609
VII-87	610

VII-89	611
VII-90	612
VII-91	613
VII-92	614
VII-93	615
VII-94	616
VII-96	617
VIII-09	618
VIII-10	619
VIII-11	620
VIII-12	621
VIII-13	622
VIII-15	623
VIII-16	624
VIII-17	625
VIII-18	626
VIII-19	627
VIII-20	628
VIII-21	629
VIII-23	630
VIII-24	631
VIII-25	632
VIII-28	634
VIII-29	635
VIII-30	636
VIII-31	637
VIII-32	638
VIII-33	639
VIII-34	640
VIII-36	641
VIII-37	642
VIII-38	643
VIII-40	644
VIII-41	645
VIII-42	646
VIII-43	647
VIII-45	648
VIII-46	649
VIII-47	650
VIII-48	651
VIII-50	652
VIII-51	653
VIII-53	654
VIII-54	655
VIII-55	656
VIII-56	657
VIII-57	658
VIII-58	659
VIII-59	660
VIII-60	661
VIII-61	662

VIII-64	663
VIII-65	664
VIII-66	665
VIII-67	666
VIII-68	667
VIII-69	668
VIII-70	669
VIII-71	670
VIII-72	671
VIII-73	672
VIII-74	673
VIII-75	674
VIII-76	675
VIII-77	676
VIII-78	677
VIII-79	678
VIII-80	679
IV-02	681
IV-04	682
IV-13	683
IV-14	684
IV-17	685
IV-28	686
IV-31	687
IV-32	688
IV-38	689
IV-40	690
IV-42	691
IV-44	692
IV-47	693
IV-55	694
IV-56	695
IV-61	696
IV-64	697
IV-65	698
IV-72	699
IV-73	700
IV-80	701
IV-85	702
IV-93	703
IV-95	704
IV-96	705
V-03	706
V-04	707
V-07	708
V-08	709
V-09	710
V-12	711
V-18	712
V-20	713
V-24	714

V-37	716
V-40	717
V-41	718
V-48	719
V-57	720
V-61	721
V-64	722
V-65	723
V-74	724
V-80	726
V-81	727
V-87	728
VI-13	870
VI-14	871
VI-16	873
VI-23	878
VI-24	879
VI-28	883
VI-32	885
VI-38	886
VI-39	887
VI-45	889
VI-46	890
VI-49	892
VI-50	893
VI-52	894
VI-53	895
VI-54	896
VI-55	897
VI-57	898
VI-58	899
VI-63	900
VI-65	902
VI-66	903
VI-67	904
VI-74	905
VI-75	906
VI-76	907
VI-80	908
VI-81	909
VI-85	910
VI-87	911
VI-88	912
VI-91	913
VI-94	914
VI-95	915
VI-96	916
I-13	1177
I-14	1178
I-30	1180

I-54	1181
I-88	1182
III-20	1183
IV-15	1185
IV-26	1186
IV-62	1187
V-11	1188
V-28	1189
V-38	1190
V-45	1191
VI-44	1193
VII-47	1200
I-42	1332
I-52	1333
I-86	1336
I-95	1337
III-10	1342
III-60	1347
III-82	1348
III-92	1351
IV-23	1353
IV-34	1354
IV-35	1355
IV-41	1356
IV-45	1357
IV-82	1359
V-01	1361
V-02	1362
V-06	1363
V-17	1364
V-25	1365
V-35	1366
V-42	1367
V-47	1368
V-49	1369
V-58	1370
V-75	1372
V-79	1373
V-90	1374
V-91	1375
V-94	1376
VI-10	1379
VI-41	1381
VI-43	1382
VI-71	1387
VI-72	1388
VI-79	1389
VI-90	1390
VI-93	1391
VII-25	1397
VII-60	1398

VII-67	1399
VIII-22	1403
VIII-26	1404
VIII-39	1405
VIII-44	1406
I-37	1440
V-32	1445
V-52	1447
V-68	1448
V-92	1449
VI-42	1450
VI-78	1452
VII-86	1453
VII-88	1454
IV-29	1490
V-15	1491
V-39	1492
V-54	1493
V-59	1494
V-71	1495

Table 5**Samples**

Diagnosis	No. of women
Normal /Benign	42*
DCIS	3
Invasive cancer	26

* From one woman, whole blood was collected at weeks 1,2,3,4,5 following menstruation. Hence, the number of unique normal/benign samples tested in the experiment is 75

Information about women with breast cancer

Sample	AGE	Stage	Cancer type	Size hist. (mm)	Nodes
1	51	II	IDC	20	1/7
2	84	II	IDC	22	2/2
3	50	I	DCIS+ 1 IDC	>50 DCIS; 5 x 14	0/7
4	47	I	IDC	15	0
5	69	III	ILC g.2 + tubular adenocarcinoma	50 + 3	1 av 12 + 1 av 7
6	50	II	IDC	24	0
7	65	I	IDC	15	0
8	63	II	IDC	23	0
9	55	I	IDC + DCIS	4	0 av 1
10	52	0	DCIS + small colloid carcinoma foci	50 + 3	0
11	60	II	IDC	24	0
12	54	I	IDC	11	0
13		0	DCIS	20	0
14	49	0	DCIS	9	0
15	48	I	IDC	4	0
16	56	I	IDC	4	0
17	68	I	IDC	14	0
18	68	I	IDC	7	0
19	63	I	IDC	10	0
20	45	I	IDC	19	1
21	57	III	IDC	60	8/20
22	55	II	IDC/DCIS	35 + 55	0
23	71	I	IDC/extensive DCIS	8	0
24	56	I	IDC	9	?

25	66	II	IDC	26	0
26	66	I	IDC	15	?
27	61	I	IDC	9	?
28	?	?	?	?	?
29	65	I	IDC	11	0

Other diseases /conditions present in the women tested

Other diseases /conditions present in the women tested

Disease/condition
Diabetes
Asthma
Ulcerous colitis
Hemochromatose
Crohn's disease
Fibromyalgia
Psoriasis
Atopic eczema
Rheumatism
Allergies

Prior history of cancer in the women tested

Cancer type	No. of women
Breast	3
Colon	2
Stomach	1
Skin	1

Table 6

Number of samples tested by double cross validation and success of the diagnostic test for breast cancer based on selected informative genes

Number of samples tested by double cross validation

Number of unique samples tested	76
Number of unique non cancer samples tested	48
Number of cancer samples tested	28

Success of the diagnostic test for breast cancer based on selected informative genes

Occurrence in percentage*	Number of informative probes	Specificity	Sensitivity	Accuracy	False Positive rate	False negative rate	Total error rate
100.00	23	84.78	75.86	81.33	15.22	24.14	18.67
90.00	44	91.30	79.31	86.67	8.70	20.69	13.33
80.00	51	86.96	79.31	84.00	13.04	20.69	16.00
70.00	54	88.13	75.86	84.00	10.87	24.14	16.00
60.00	58	89.13	75.86	84.00	10.87	24.14	16.00
50.00	59	89.13	75.86	84.00	10.87	24.14	16.00
40.00	63	88.13	75.86	84.00	10.87	24.14	16.00
30.00	66	86.86	75.86	82.67	13.04	24.14	17.33
20.00	74	89.13	75.86	84.00	10.87	24.14	16.00
10.00	79	89.13	75.86	84.00	10.87	24.14	16.00
5.00	90	86.96	79.31	84.00	13.04	20.69	16.00
1.33	139	84.78	72.41	80.00	15.22	27.59	20.00

*100% = Genes appearing in all the 75 cross validated models; 90% = Genes appearing in at least 68 out of 75 cross validated models; 5% = Genes appearing in at least 4 out of 75 cross validated models; and so on.

Table 7

Double cross-validation and details of the success of the diagnostic test for Alzheimer disease based on the expression 182 informative genes

Validation Result

Total number of samples tested	14
Number of Alzheimer's disease samples tested	7
Number of Alzheimer's disease samples incorrectly predicted	1
Number of non-Alzheimer's disease samples tested	7
Number of non-Alzheimer's disease samples incorrectly predicted	0

Success of diagnostic test

Performance	Description	%
Accuracy	Percentage of the total number of predictions that were correct	92.9
Sensitivity	Percentage of positive cases that were correctly identified	85.7
Specificity	Percentage of negatives cases that were correctly predicted	100
False positive rate	Percentage of negatives cases that were incorrectly classified as positive	0.0
False negative rate	Percentage of positives cases that were incorrectly classified as negative	14.3
Total error rate	Percentage of the total cases incorrectly predicted	7.1

Table 8

Some relevant features of the blood donors. B, Female donors with breast cancer; N, Female donors with suspected mammogram but no breast cancer; **IDC**, invasive ductal carcinoma; **DCIS**, ductal carcinoma in situ; na, not available nd, not determined; ++, no degradation of mRNA and no ribosomal contamination in the sample, +, no degradation of mRNA but ribosomal contamination in the sample.

		AGE	Cancer type /breast abnormality	Size Hist. (mm)	mRNA Quality
1	B1	na	IDC	5	++
2	B2	49	DCIS	8	nd
3	B3	54	IDC	18	++
4	B4	59	IDC	12	+
5	B5	61	DCIS+micro invasive cancer	15+1.5	++
6	B6	55	IDC	12+17	nd
7	B6		IDC	12+17	nd
8	N1	45	Fibroadenoma	-	nd
9	N2	52	na	-	+
10	N3	55	Cyst	-	++
11	N4	54	na	-	++
12	N5	51	Benign ductal epithelium	-	nd
13	N6	57	Benign	-	nd
14	N7	50	na	-	++
15	N8	52	na	-	+

Table 9

List of sequence of probes informative for both alzheimer and breast cancer disease

Clone ID	Sequence ID
I-24	308
I-25	309
I-28	310
I-48	321
I-60	327
I-72	333
I-81	338
I-82	339
II-02	360
II-03	361
II-06	364
II-07	365
II-10	368
II-21	378
II-23	380
II-24	381
II-25	382
II-27	384
II-33	390
II-34	391
II-41	397
II-42	398
II-46	401
II-47	402
II-48	403
II-52	406
II-57	411
II-58	412
II-59	413
II-60	414
II-61	415
II-62	416
II-64	418
II-67	421
II-69	423
II-70	424
II-74	428
II-80	434
II-82	436
II-84	438

II-87	441
II-88	442
II-96	450
III-01	452
III-02	453
III-06	458
III-08	460
III-12	463
III-13	464
III-17	468
III-18	469
III-21	471
III-23	473
III-24	474
III-25	475
III-26	476
III-27	477
III-28	478
III-29	479
III-32	482
III-33	483
III-35	485
III-39	487
III-40	488
III-42	489
III-45	492
III-46	493
III-47	494
III-48	495
III-56	503
III-57	504
III-58	505
III-59	506
III-61	507
III-62	508
III-63	509
III-64	510
III-66	512
III-67	513
III-70	515
III-74	518
III-75	519
III-78	521
III-80	523
III-81	524
III-85	526
III-86	527
III-88	529
III-89	530
III-93	532
III-95	534

III-96	535
IV-04	682
IV-13	683
IV-14	684
IV-17	685
IV-31	687
IV-32	688
IV-38	689
IV-42	691
IV-47	693
IV-61	696
IV-64	697
IV-72	699
IV-80	701
IV-85	702
IV-93	703
IV-96	705
V-03	706
V-04	707
V-07	708
V-08	709
V-12	711
V-24	714
V-41	718
V-57	720
V-61	721
V-64	722
V-65	723
V-74	724
V-80	726
VI-03	341
VI-04	342
VI-07	344
VI-08	345
VI-09	346
VI-12	869
VI-14	871
VI-19	349
VI-20	350
VI-21	351
VI-23	878
VI-25	353
VI-26	354
VI-48	359
VI-50	893
VI-53	895
VI-74	905
VI-76	907
VI-87	911
VI-88	912
VI-95	915

VII-02	547
VII-03	548
VII-06	551
VII-08	553
VII-09	554
VII-10	555
VII-11	556
VII-15	559
VII-17	560
VII-19	562
VII-21	564
VII-22	565
VII-23	566
VII-24	567
VII-27	568
VII-29	570
VII-32	571
VII-33	572
VII-36	575
VII-39	576
VII-41	578
VII-42	579
VII-43	580
VII-46	583
VII-48	584
VII-49	585
VII-54	589
VII-57	591
VII-58	592
VII-59	593
VII-62	594
VII-63	595
VII-64	596
VII-66	598
VII-72	600
VII-73	601
VII-77	604
VII-80	605
VII-82	607
VII-87	610
VII-90	612
VII-91	613
VII-92	614
VII-93	615
VII-96	617
VIII-09	618
VIII-10	619
VIII-13	622
VIII-16	624
VIII-20	628
VIII-21	629

VIII-23	630
VIII-24	631
VIII-25	632
VIII-28	634
VIII-29	635
VIII-30	636
VIII-31	637
VIII-32	638
VIII-33	639
VIII-34	640
VIII-38	643
VIII-40	644
VIII-41	645
VIII-46	649
VIII-48	651
VIII-55	656
VIII-57	658
VIII-59	660
VIII-60	661
VIII-61	662
VIII-64	663
VIII-66	665
VIII-73	672
VIII-74	673
VIII-76	675
VIII-80	679

- 124 -

Nucleotide sequences

Sequence ID - 93

nt: 405

5 GGATCCTGTGGCCACAGAGCTGCCCCAGCAGACGCTCCGCCCCACCCGGTGATGG
AGCCCCGGGGGGACAATCGTGCCCTGGGGAGGAGCAGGGTACAGCCCATTCCCCAG
CCCTGGCTGACCTGGCCTAGCAGTTTGGCCCTGCTGGCCTTAGCAGGGAGACAGGG
GAGCAAAGAACGCCAAGCCGGAGGCCCCGAGGCCAGCCGGCCTCTCGAGAGCCAGAG
CAGCAGTTGAATGTAATGCTGGGGACAGGCATGCTGCCGCCAGTAGGGCGGGGACC
CGGACAGCCAGGTGACTACCAGTCCTGGGGACACACTCACCATAAACACATCCCCA
10 GGCAGGACAGATCGGGGAAGGGGTGTGTACCAGGCTATGATTCTCTTGCATTAA
ATGTATTATTATT

Sequence ID - 108

nt: 550

15 GGCTTTGACAGAGTGCAAGACGATGACTTGCAAAATGTGCGCATCTGGAACGCAACA
TAGANACCATCATCAACACCTTCCACCAATACTCTGTGAAGCTGGGGCACCCAGAC
ACCTTGAACCAGGGGGAATTCAAAGAGCTGGTGCGAAAAGATCTGCAAAATTTTCT
CAAGAAGGAGAATAAGAATGAAAAGGTCATAGAACACATCATGGAGGACCTGGACA
CAAATGCAGACAAGCAGCTGAGCTTCGAGGAGTTCATCATGCTGATGGCGAGGCTA
ACCTGGGCCTCCCACGAGAAGATGCACGAGGGTGACGAGGGCCCTGGCCACCACCA
20 TAAGCCAGGCCTCGGGGAGGGCACCCCCTAAGACCACAGTGGCCAAGATCACAGTG
GCCACGGCCACGGCCACAGTCATGGTGGCCACGGCCACAGCCACTAATCAGGAGGC
CAGGCCACCCTGCCTNTACCCAACCAGGGCCCCGGGGCCTGTTATGTCAAACCTGTC
TTGGCTGTGGGGCTAGGGGCTGGGGCCAAATAAAGTCTCTTTCTCC

25 Sequence ID 110

ACGAAGACAGACATCTGTGGAATGATTCACATCCTCTCAAGTTAGGAGGATGGAGG
CCTGCTTCATTAAGAAGCTGGGGGTAGGGTGGGGGTGGGGAGAACACTTAACAACA
TGGGGACCAAGTCAGGGGAATCCCCTTATTTCTGTTTTGCATATGAGGAACCCTAGA
GCAGCCAGGTGAGGCTCTCTAGTTTAATAAAAATCATGGAAAGACTCTTAATGCAG
30 ACTCTTCTTAAGTGTTAATAGGGATTTTTTCAGCTTATTTTGGTTGCAGTTTCCAA
TTTTTAAAAATGTTGAGGTAATCTTTCCACCTTCCCAAACCTAATTCTTGTAGAT
GCATTAGTGTTGAACCAATGCTTTCTCATGTCTCAATTCTTTGTATATGCATTCTT
TTCAGATGTATTAAACAAACAAAAACCCTTC

35 Sequence ID - 192

nt: 286

CCGGTAATAGAATAGAAAAGGGAGAGTGTCTTCATGCAATGTGGCATCCTGGATTG
GGTCTCGNNACAAAAACAGGACATTAGTGGGAAAATTGGAAATCTGAAAAAGTCT

- 125 -

GAATTTTAGTTAATATACCAATTTCACTCTCTTGGTTTTGACAGATGTACCATGGT
GATGTAAGATGTTGACCTTGGGGTAGGCTGGGTGAAGGGTATACAGGAACTCTTTG
TACTATCTCTGCAACTTCTCTGTAAATCTAGTATCATTCCAAATAAAAGTTTATT
TAATTT

5

Sequence ID 250

GTGGAAGTGACATCGTCTTTAAACCCTGCGTGGCAATCCCTGACGCACCGCCGTGA
TGCCCAGGGAAGACAGGGCGACCTGGAAGTCCAACACTTCTTAAAGATCATCCAA
CTATTGGATGATTATCCGAAATGTTTCATTGTGGGAGCAGACAATGTGGGCTCCAA
10 GCAGATGCAGCAGATCCGCATGTCCCTTCGCGGGAAGGCTGTGGTGCTGATGGGCA
AGAACACCATGATGCGCAAGGCCATCCGAGGGCACCTGGAAAACAACCCAGCTCTG
GAGAACTGCTGCCTCATATCCGGGGGAATGTGGGCTTTGTGTTACCAAGGAGGA
CCTCACTGAGATCAGGGACATGTTGCTGGCCAATAAGGTGCCAGCTGCTGCCCGTG
CTGGTGCCATTGCCCCATGTGAAGTCACTGTGCCAGCCCAGAACACTGGTCTCGGG
15 CCCGAGAAGACCTCCTTTTTCCAGGCTTTAGGTATCACCCTAAAATCTCCAGGGG
CACCATTGAAATCCTGAGTGATGTGCACTGATCAAGACTGG

Sequence ID 299

CAGCGCAGGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCGCAGATAAG
20 TTTTTTCTCTTTGAAAGATAGAGATTGNTACAACTACTTAAAAAATATAGTCAAT
AGGTTACTAAGATATTGCTTAGCGTTAAGTTTTTAACGTAATTTTAATAGCTTAAG
ATTTTAAGAGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAAA
AGGTTTCTAAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAATG
TATTTAAAAGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACCAATAGA
25 AGGGCAATGCTTTTAGATTAAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTTT
AAAAGTTGTAGGTGATTAAAATAATTTGAAGGCGATCTTTTAAAAGAGATTAAAC
CGAAGGTGATTAAAAGACCTTGAAATCCATGACGCANGGAGAATTGCGCATTTAAA
GCCTAGTTACGCATTTACTAAACGCAGACGAAAATGGGAAGATTAATTGGGAGTGG
TAGGATGAAACAATTTTGGAGAAGATAGAAG

30

Sequence ID 300

CTCAAAGGAGAAAAAAACCTTGTAAAAAAAGCAAAAATGACAACAGAAAAACAAT
CTTATTCAGCATTCCAGTAACCTTTTTTGTGTATGTACTTAGCTGTACTATAAGT
AGTTGGTTTGTATGAGATGGTTAAAAGGCCAAAGATAAAAGGTTTCTTTTTTTTTT
35 CCTTTTTTGTCTATGAAGTTGCTGTTATTTTTTTTTTGGCCTGTTTGATGTATGTGT
GAAACAATGTTGTCCAACAATAAACAGGAATTTTATTTTGCTGAGTTGTTCTAAAA
AAAAAAAAAAAAAAAAAAAA

- 126 -

Sequence ID 302

AGTAGAGACGGGGTTTCACTGTGTTAGCCAGGATGGTCTCGATCTCCTGACCTCGT
GATCCGGCCACCTCGGCCTCCCGAAAAGTGCTGGGATTACAGGCGTGAGCCACGGCG
CCCAGCCCCAGCCTGTCACTTAAACTGATAAACGACAGATTAACAGTAGAAAAATT
5 TTATTTTGCATACATAATGAGGCTTCACAAAAGAGAAGTGAAAACCCAAGTAGGAG
TTTAGGGCTGGGGGCTTATATACCATTTAACAAGGGGTGATAAATTGTAAGAGAAT
AG

Sequence ID 304

10 TCCTTGGTTTTCGATTTGTGGCAACAATCCAGTCTTTTTGTTTTTTTTCAGGGATACC
ATATGTAACAGGTGCCATTGTTACTGTAACCTTTTCACACATGCCTTCAGTTTGATG
TCAAAGTCATCATTTAGTGTAACAGCAAGTTATCTGTTAGGCTGCACATCATGAA
CTTTACTTTTAGAAAAGTCTTATCTTTTATGCCACAGAAAATAGCATTGCGCTATTAG
TCATGGATGGCAAAGAAATTAATTTTGAGTTGTTTGGATAAAAATGTTTCAGTTGA
15 CTGTAGTGTGTATTGAGAGACACTGCCAGTAAACAAACTCTCTTGGTAGGTGGAAA
TCCCCTAGAAGTTACAGAAAATTGGGAGGAGGTGAACTTAATTAAATAACTTGAAT
TGTTTAGACATATTCAGAGCTTCTTATGACCTTGAAGAAATCACCCAACCTTCAAAA
GACCTCGGTTTCTTCATTTGTAAAATTAGGGAGTTTGACTAGATGTGTAAATCTAG
TTGTTAGTTAACTTCTAAGATGTAAAAACCCTCTTGTTTAAACAAAACCTACAAGA
20 TCAAGTTGCTTATCTGAAATCTTTATGAATCAACACTAGTCACTAAGTCTAGCTCG
ACC

Sequence ID 306

CTTTTCCTCCCGCTGTCCCCACGGAGGGGACTGCTCTCCCCGCTGCATCCTTTC
25 TGTGAGGTACCTTACCCACCTCAGCACCTGAGAGGGTGAAATAGAATTCTAACCTC
GACATTCGGGAAGTGTTTTTGTAGAAGTCTCGGTCGGTAAGGGAAGTCTTCCAAGTC
CGTGCAGCACTAACGTATTGGCACCTGCCTCCTCTTCGGCCACCCCCAGATGAGG
CAGCTGTGACTGTGTCAAGGGAAGCCACGACTCTGACCATAGTCTTCTCTCAGCTT
CCACTGCCGTCTCCACAGGAAACCCAGAAGTTCTGTGAACAAGTCCATGCTGCCAT
30 CAAGGCATTTATTGCAGTGTACTATTTGCTTCCAAAGGATCAGGCCCTGAGAACAA
TGACCTTATTTCTACAACAGTGTCTGGGTTGCGTGCCAGCAGATGCCTCAGATAC
CAAGAGATAACAAAGCTGCAGCTCTTTTGATGCTGACCAAGAATGTGGATTTTGTG
AAGGATGCNCATGAANAAATGGACNAGCTGTG

35 Sequence ID - 308

nt: 373

AAGTGGGTCTTGCCATCCCTGAACTGNAATCATCCCTAACATATTCATACCTGTTT
TCATTTTAAAAGTTGGGTCAGTTTTTTTATTAGTACATGTATTTCTATCCTACTGA

- 127 -

TTTATTTGCTATATCATCTAATTTAGTTTGAATATTCCATAATTTACTTAATTAGT
CCTGTATGGAGACCTAGCTCTTCTCAGTGTCTACTATTATAAACAATGCTACAGTG
AATATTGGTGNATAAATCCATACNCACCACGTACATATCTTAAGTTCTGGAAGAGA
TATTGCTAAACCAGAAGATAACCTGCATTTAAAATTTGACTGCTAGGGNCAGGGNC
5 ACATTTAATTAAATTAGAACAANGAATGCATAATGNC

Sequence ID 309

CCGGAATCGCGGCCGCGTCGACGAAAATATGTGCCCTGGCCAACTCCACAGGACTA
GTTCTAGGCAATCTGAAGGAAACCAGAAAATGTGAATTTCTCTTCCCTCAAAAAGC
10 TATACTGAAGTAGTATTTAATATTCAAGTACTTGTAATTTGCAGAACAGTACTTT
TTAATTTGACCCATGAATTTCTATTTAAATTTGTCACTTAATATTTAGCCAAGAAGC
AAACCATCTAAAAGATTTCTGGTTTATTTCTCCAACCTCCTAATAAATAGGGTCAC
ATATTTTTTAACTTTTTTCTAATTTGAAAAGTAATACAGGCATATGGTATTTTTAAA
AATGAAACAACACAAAGGGATATGTTTTGAAAAGTGGTCTTGCCATCCCTGAACTG
15 TAATCATCCCTAACATATTCATACCTGTTTTTCATTTTAAAAGTTGGGTGAGTTTTT
TTATTAGTACATGTATTTCTATCCTACTGATTTATTTGCTATATCATCTAATTTAG
TTTGAATATTTCCATAATTTACTTAATTTAGTCCTGTATGGAGACCTAGCTCTTCTCA
GTGTCTACTATTATAAACAATGCTACAGTGAATATTGGTGNATAAATCCTACACAC
CACGTAACATATCTTAAGTTCCTGGAAGAGATATTGCTAAACCAGAAGATAACCTG
20 CATTTAAAATTTGACTGCTAGGGTCAGGGTCACATTTAAATTAAATTAGAACAAGG
AATGCATAATGTCTTCGATAGCAATCTATTCAAGGTGCACCGTGGTCACAAAGGAA
AGCAAAACTGTC

Sequence ID - 310 nt:564

25 CCTGGNCAGAGGCCTCTATCCTGTANTGATAATTGCCATCAAAATTGTCAAAAANG
ATTTAATTTCTATGGGNAATAGTCCTTTTCTTAGCTTCTGCCNNTCACTTGCTTAT
TTTTTGTGTGGGAATGGGGTTGGATAAACCAATGAACTTTATTATAAACAATCCC
ACCTATATCTANCAATTTATATTTTCGGTGAAATACAGATATTTGCCTTTCTGGA
GTANTATAGAAGCTGTCAATATGTATCTACTGTACAGTACTAAATAGTATTCATTT
30 ATGAAATGAGTAGTGTTTGGGTGGCTGGGGTTAAGGAAAAATGAGACTTGGAATTG
TAGCTTTTATCCAAGTTTTGAGTATAAATAGGGTTTTGTGTTTTGTTTTTTTAACTT
AAAACTGAAATGCCATATAGAAAAACAGCATTGTTTTTACAGTTTGTAGTAAGTA
ACTTTTTTAAAGATTTTATCAAAAAGAATTTTGTCTATNGTGAGTAAAAGAAGTTCT
AATAATGGCCTAATCACTGCATTTTAAAAAACAAAGTTCAACACAAATGACATTT
35 GTTT

Sequence ID 311

- 128 -

CCTCTCCTCCATCTAAAGGCAACATTCCTTACCCATTAGTCTCAGAAATTGTCTTA
AGCAACAGCCCCAAATGCTGGCTGCCCCGGCCAAGCATTGGGGCCGCCATCCTGC
CTGGCACTGGCTGATGGGCACCTCTGTTGGTTCCATCAGCCAGAGCTCTGCCAAAG
GCCCCGCAGTCCCTCTCCCAGGAGGACCTTAGAGGCAATTAAATGATGTCCTGTTC
5 CATTGG

Sequence ID - 313

nt: 554

CCCGGAATCGCGGCCCGCGTCGACAACAAACCTGCATGTTCTGCACATGTATCCAG
GAACTTAAAAAAAAAAAAAGATAGTTTGTGTGTCTTAATTGAATAATAGTAGATTT
10 ATAGATTAAAGATCTATGGGTTTTTTAATATGGATTANAAATCTGTGGGTTTTTGAT
ATGGATTANAAATCTGTGGGTTTTTTAATATGGATTGGAAATCTGTGGGTTTTTAAT
ATGGATTAAAAACATCTGTGGGTTTTTTAATATGGATTAAACATCTGTGGGTTTTT
AATATGGATTAAACATCTGTGGGTTTTTTAATATGGATTAAACATCTGTGGGTTTTTAA
TATGGGTTAAAAATCAAAGAAAATGAACTATTTGCTCCAGTGCAGGAAAATACAG
15 GCAATACTGGATACAATTAGATGGTCAGGAGCGATAACCCGGTTGCCATTGTTTGA
AGAAGAGAATAAGGNGCTAGCATTCCTATCCGTAGATAATTTGACAGCTAGGAAAT
AGGGGGAGTCTTCTATGTAGTTAGTGAAGGCTAAATGAACTATTATATGC

Sequence ID 314

CTTTTCCTCCCGCTGTCCCCCACGGAGGGGACTGCTCTCCCCCGCTGCATCCTTTTC
20 TGTGAGGTACCTTACCCACCTCAGCACCTGAGAGGGTGAAATAGAATTCTAACCTC
GACATTCGGGAAGTGTTTTTTGAGAAGTCTCGGTCCGTAAGGGAAGTCTTCCAAGTC
CGTGCAGCACTAACGTATTGGCACCTGCCTCCTCTTCGGCCACCCCCCAGATGAGG
CAGCTGTGACTGTGTCAAGGGAAGCCACGACTCTGACCATAGTCTTCTCTCAGCTT
25 CCACTGCCGTCTCCACAGGAAACCCAGAAGTTCTGTGAACAAGTCCATGCTGCCAT
CAAGGCATTTATTGCAGTGTACTATTGCTTCCAAAGGATCAGGCCCTGAGAACAA
TGACCTTATTTCTTACAACAGTGTCTGGGTTGCGTGCCAGCAGATGCCTCAGATAC
CAAGAGATAACAAAGCTGCAGCTCTTTTGATGCTGACCAAGAATGTGGATTTTGTG
AAGGATGCACATGAAGAAATGGAGCAGGCTGTGGAAGAATGTGACCCTTACTCTGG
30 CCTCTTGAATGATACTGAGGAGAACAACCTCTGACAACCACAATCATGAGG

Sequence ID 315

TGGTACAGATACAACTGGACTCTCAGGACAAAACGACACCAGCCAAACCAGCAGC
CCCTCAGCATCCAGCAGCATGAGCGGAGGCATTTTCCTTTTCTTCGTGGCCAATGC
35 CATAATCCACCTCTTCTGCTTCAGTTGAGGTGACACGTCTCAGCCTTAGCCCTGTG
CCCCCTGAAACAGCTGCCACCATCACTCGCAAGAGAATCCCCCTCCATCTTTGGGAG
GGGTTGATGCCAGACATCACCAGGTTGTAGAAGTTGACAGGCAGTGCCATGGGGGC

- 129 -

AACAGCCAAAATAGGGGGGTAATGATGTACGGGCCAAGCACTGCCCAGCTGGGGGT
CAATAAAGTTACCCTTGTACTTG

Sequence ID 316

5 CGCCACTTATCCAGTGAACCACTATCACGAAAAAACTCTACCTCTCTATACTAAT
CTCCCTACAAATCTCCTTAATTATAACATTCACAGCCACAGAACTAATCATATTAA
AA

Sequence ID 321

10 CAGAACAGTACTTTTTAATTTGACCCATGAATTCTATTTAAATTTGTCACTTAATA
TTTAGCCAAGAAGCAAACCATCTAAAAAGATTTCTGGTTTATTTCTCCAACCTCCTA
ATAAATAGGGTCACATATTTTTTAACTTTTTTCTAATTTGAAAAGTAATACAGGCA
TATGGTATTTTAAAAATGAAACAACACAAGGGATATGTTTTGAAAAGTGGTTCTT
GCCATCCCTGAACTGTAATCATCCCTAACATATTCATACCTGTTTTCATTTTAAAA
15 GTTGGGTCAGTTTTTTTATTAGTACATGTATTTCTATCCTACTGATTTATTTGCTA
TATCATCTAATTTAGTTTGAATATTCCATAATTTACTTAATTAGTCCTGTATGGAG
ACCTAGCTCTTCTCAGTGTCTACTATTATAACAATGCTACAGTGAATATTGGTGT
ATAAATCCATACACACCACGTAACAATATCTTAAGTTCCTGGAAGAGATATTGCTAA
ACCAGAAGATAACCTGCATTTAAAATTTTGACTGCTAGGGTCAGGGTCACATTTAA
20 ATTAAATTAGAACAAGGAATGCATAATGTCTTCGATAGCAATCTATTCCAGGTGCA
CCGTGGTCACAAAGGAAAGCAAACTGTCAATAACTTTCTTCTCA

Sequence ID 322

25 TAGCATTGCGCCTTTTAAAACATTTGTTTTATTTTTTTCTGAGAATGGCTAACACA
CTTTATTGAGGTTGGAATTAATAAAGAAAATAAAAGAAATGTATCTTCATTTCATT
CTGTATGTTAGTGTTTTAATTACCCTTAGAATATATGGATAAAAAATACTATTCTT
TGTCTTGGAGAAGGTAAGAGTCTAGTTAGATGAATAAGGGTTATCTATGTAGAACA
ACTAGAGAATGAGAAGAGAGCTTATGAGATTGAGTACTACGTTATGCAGTAGAGTA
GCACGTCATCTGCTACTGAGTATGGTGTGATAACATTGTGTAACAGGAAAGTATGA
30 TCAATATCTACTTAAAATTAAGGACAATATTAGCACTACATTGCTTTATTTTAAAG
TAAAAATTAGAGAACTAAACACAAGCATTGTAAGTACAATAAAAGCTGATCTTTCT
AGTTAAGCAGAATAATACATGTTCAAGCATCTGCTAAATCATTAATATAAGAATA
TAGGGGTTTTCTATAATCTTATTTTCTTTGGAAGAGTACCTCATTTTCAAGANGAG
AAGTTTCTAATTGCCACTTCTTTAAAATAAAACAGGGTTTTAATGTTCCCAGCAC
35 AAAAATTAATATCTCTTCAAAAAGTCTCTTGTGATTAAAGTTTGAATCCCTTGTCTAT
ACTGCTTCTAATATTGACACTGACCTCCTTAGGTATTTTTTCAGGGGTTATAATCTT
TTCTTAAGGTATCTTTTTTCAAGAATTGGATACCTTGGGCTT

- 130 -

Sequence ID 323

CGCGTCGACTTTTAAAGTCATCTCTATAGGAAGGTGCTGGGCAGGGATCCCAGAGA
AAGAAAGGGTCCAAGACTCCATTAAGTGCCTTGGATGAAGGGCACTGCTACAGCAG
CTAGTACCAGAGACTCTCCTATCTCACGGTTGAGGCAGACCCAGGATAGAATAGAG
5 AATAAAAGGAATGCTTATAGGAAACAATTTTGTATGGAATGCTAGATGGCCAAGCC
TCAGCCTTTGGTCCAGTGCAACCCTTGCCTCGCTTGTCAACAGTGAAAAATTAGTT
TGGTTAGAAGAACCATCTGGAAACACACCAGCTTCTGCTACCTTCATGCTCATTGT
TAAAAAAGATTAACCAGTGTGAACATTCTGATCTGTTAATTCCAGGGACTGTTTT
CTTTCCAATGGACTGTTTGGTTGGTAGAATAACCCCCAAAAGCTCAAAGCTAAAATG
10 CATCATCAGTCCTAGTCGGCAGTTCCTTAAGAATGGACTGGCGGCGTGTTGAGCT
GATATGGAAAAGCTGCACCTTCCTGCAGAAGATCAACTGACCTGCTATCCCACCCC
AAATTCAACCTGAGGTATATTTTCAGTGAAGCAGGTAGCTGTGCTTCTCAAAGCAGA
GAAGCAGTTTTTAAGAACC AAAAAGGTAGAGGAAATCTA

15 Sequence ID 324

GTTTGTTACAGGCAGAATTGGATAGATACAGCCCTACAAATGTATATGCCCTCCCC
TGAAAAAATTGGATGAAAATCTGCACAGCAAAGTGAAACACACAGATAATAGGAA
CAAAATGTAGTTCCCATGTGCCAAACAAAATAAATGAAATCTCTGCATGTTTGCAG
CATATCTGCCTTTTGGGAATGTAATCAAGGNATAATCTTTGGCTAGTGTTATGTGC
20 CTGTATTTTTTTTAAATGGTACACCAGAAAAGGACTGGCAGTCTACTTCTACCATA
GTTAAACTTCACCCTCTTTAATTTCAACAATATTTCTTTGGAAGCAGGAAGAAATG
CTCATAAAGAGGATCAGACCTTCTTTCCCGTGAAACCAGTATTTGGCGCCATATAT
AAGCCTGGTTAAATTGGTCATCTAAAGCTGTCAAATAAGACATTCTGTGAAAGGTA
AACATCGAAACTGGTTATAAGTAAAACCATCAAGCCAACAACAGGGTCTTGAGATA
25 ACCTTTGAAGCTTATTGTCTGGCCTGCACCAGAAGATGTCTGCATTACTCATTGCT
AAAAATGTGTACACAGAACTGCACTAGGATTAATTGGTTCAAGAAGAAATTTAAAC
TTACGTTTGGGTTTCCATACAGCACTCTATTGAATACATGCATCTGAATTTAAGTT
GCAA

30 Sequence ID 325

GACCAGTAATGGCTTTTAAAGAGTCCATTTTGTCAATTGTCTCCCTAGTTAATTACAG
GTGGGGGATCTTTTGCCTCTATTCTCTTCATATTGAAATGAATCATACTCATGTTT
TGTGGAACCTCCTTAAAGTTGTAGCTGTCAATGATCAGATTTTTTTTATATTTCTCA
GCTTAACTCTGCTACTTGATTTACAGTGACCCATAACCTACTCATCCTTGTTTAT
35 AGTGACACATAATCTTATCTCTTTATAGAACCCTTAAATTTTATCATTATTTTCGCT
TAGAATACAGCATTCTTTGCTTCTGTTGCTGGTTTGACTTAAGAAATAAGGCAGT
AACTCTGATCAATCAATTATCCATAAGGAAGGGCTTTTCATGGGTTCTATTAATTT

- 131 -

GTTAGTACCCTAAGTATATCTGAAAAATATGTCTATTGAGAGAAGATTTTGGCATT
CCAGATGGTATAGTCTATATATATTTAAAGTTTGAATTTGCTTATATATACTCAG
CTTTCTTTTTCTAGCATTTTGCATTTACCTGTTAATTGAAGTATACCCCCACAT
ATAAAAGTTCCTCTTAAAGACACTGGACTCTTTCTGGGGGGCTAAAATA

5

Sequence ID - 326

nt: 554

CCCGGAATCGCGGCCCGCGTCGACAACAAACCTGCATGTTCTGCACATGTATCCAG
GAACTTAAAAAAAAAAAAAGATAGTTTGTGTGTCTTAATTGAATAATAGTAGATTT
ATAGATTAAAGATCTATGGGTTTTTAATATGGATTANAAATCTGTGGGTTTTTGAT
10 ATGGATTANAAATCTGTGGGTTTTTAATATGGATTGGAAATCTGTGGGTTTTTAAT
ATGGATTAAAAACATCTGTGGGTTTTTAATATGGATTAAACATCTGTGGGTTTTT
AATATGGATTAAACATCTGGGTTTTTAATATGGATTAAACATCTGTGGGTTTTTAA
TATGGGTTAAAAATCAAAGAAAATGAACTATTTGCTCCAGTGCAGGAAAATACAG
GCAATACTGGATACAATTAGATGGTCAGGAGCGATAACCCGGTTGCCATTGTTTGA
15 AGAAGAGAATAAGGNGCTAGCATTCCTATCCGTAGATAAATTTGACAGCTAGGAAAT
AGGGGGAGTCTTCTATGTAGTTAGTGAAGGCTAAATGAACTATTATATGC

Sequence ID 327

CGGCTACCGACAGAAGGACTATTTTCATCGCCACCCAGGGGCCACTGGCACACACGG
20 TTGAGGACTTCTGGAGGATGATCTGGGAGGGGAAGTCCCACACTATCGTGATGCTG
ACGGAGGTGCAGGAGAGAGAGCAGGATAAATGCTACCAGTATTGGCCAACCGAGGG
CTCAGTTACTCATGGAGAAATAACGATTGAGATAAAGAATGATAACCTTTTCAGAAG
CCATCAGTATACGAGACTTTCTGGTCACTCTCAATCAGCCCCAGGCCCGCCAGGAG
GAGCAGGTCCGAGTAGTGCGCCAGTTTCACTTCCACGGCTGGCCTGAGATCGGGAT
25 TCCCGCCGAGGGCAAAGGCATGATTGACCTCATCGCAGCCGTGCAGAAGCANCAGC
AGCAGACAGGCAACCACCCCATCACCGTGCACTGCAGTGCCGGAGCTGGGCGAACA
GGTACATTTCATAGCCCTCAGCAACATTTTGGAGCGAGTAAAAGCCGAGGGACTTTT
ANATGTATTTCAAGCTGTGAAGAGTTTACGACTTCAGAGACCACATATGGTGCAAC
CCTGGAACAGTATGAAATGTGCTACAAAGTGGTACAAGATTTATTGATATATTTCT
30 GATTATGCTAATTTCAATGAAGATCCTGCCTTAAATATTTTTTAATTTAATGGCAN
AT

Sequence ID 328

CAAGACTCCATCTCAAAAAAAAAAAAAATCTACAGTGCTGAGTATATAAAATTAT
35 TAACACATTTCAACAATATGTGTTTGTGGAGTTAAATATTTTTTGTCTTTAAAA
CAGGTAATTTTAGTGCATACTTAATTTGATGATTAAATATGGTAGAATTAAGCATT
TTAAATGTAAATGTTTGTACATTGTTCAAGAAATAAGTAGAAATATATTCCTTTG

- 132 -

TTTTTTATTTAAATTTTGTTCCTCTGTAAACTAAAAGAACACGAAGTAATTGGTC
ACAATTACTGGTGTTTAACTGCCAAATATGGGTAAATAAGGGAAAATTTTGTTTAA
TATTTAGTCCTTCTGAGATGGCTTGAATATTTGAATTTTGTGTACGTCTATACTG
GGTAGTCACAAGTCTTATAAACACTTTAGAGGAAAGATGGATTTTCAGTCTGTATTT
5 TTAACATCATTTATTTTAAATCTGGTGCTGAAAAATAAGAAAAAATTAACCTGC
ATTCTGCTGTTCTTCTTTANAAGCATTCTGCGTAAATACTGCTGTAATACTGTCA
TGCAAAGTGTATCCTTTCTTGTCTATCCTTTTTTGGGGCAGTGGTTTTT

Sequence ID 330

10 GCGGGAATCGCGGCCCGCTCGACCTCAAAGGAGAAAAAAACCTTGTAAAAAAG
CAAAAATGACAACAGAAAAACAATCTTATTCCGAGCATTCCAGTAACTTTTTTGTG
TATGTACTTAGCTGTACTATAAGTAGTTGGTTTGTATGAGATGGTTAAAAAGGCCA
AAGATAAAAGGTTTCTTTTTTTTTCTTTTTTGTCTATGAAGTTGCTGTTTATTTT
TTTTGGCCTGTTTGATGTATGTGTGAAACAATGTTGTCCAACAATAAACAGGAATT
15 TTATTTTGCTGAGTTGTTCTAAAAAATAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAATTTTAAAATTTTAAAATAAAACCCTTGGTTAT

Sequence ID 331

GCCGCGTCGACCTGCATGAGCCACAGTTTCTTGACTGGAGGCCATCAACCCTCTTG
20 GTTGAGGCCTTGTTCTGAGCCCTGACATGTGCTTGGGCACTGGTGGGCCTGGGCTT
CTGAGGTGGCCTCCTGCCCTGATCAGGGACCCTCCCCGCTTTCCTGGGCCTCTCAG
TTGAACAAAGCAGCAAAACAAAGGCAGTTTATATGAAAGATTANAAGCCTGGAAT
AATCAGGCCTTTTTAAATGATGTAATTCCCCTGTAATAGCATAGGGATTTTGGAA
CAGCTGCTGGTGGCTTGGGACATCANTGGGGCCAAGGGTTCTCTGTCCCTGGTTCA
25 ACTGTGATTTGGCTTTCCTGTCTTCTGCTGGTGATGCCTTGTTTGGGGTTCTGTG
GGTTTGGGTGGGAAGAGGGCCATCTGCCTGAATGTAACTGCTAGCTCTCCGAAGC
CCTGCGGGCCTGGCTTGTGTGAGCGTGTGGACAGTGGTGGCCGCGCTGTGCCTGCT
CGTGTTGCCTACATGTCCCTGGCTTGTTGAGGCGCTGCTTCAACCTGCACCCCTCC
TTGTCTCATAGATGCTCCTTTTGACCTTTTCAAAATTAATATGGATGGGAAAGCTC
30 CTATGCCTTTTGGCTTCTGCTAGTAAGGCGGGATGCCCAAGGGTCTGCCTGGGTGT
GGATTGGATGCTTGGGGTGTGGGGGTGGAACCTGTCTTGTGGCCCACTTGGGCCC
C

Sequence ID 335

35 CCCGCGTCGACTTTTAAAGTCATCTCTATAGGAAGGTGCTGGGCAGGGATCCCAGA
GAAAGAAAGGGTCCAAGACTCCATTAACCTGCCCTGGATGAAGGGCACTGCTACAGC
AGCTAGTACCAGAGACTCTCCTATCTCACGGTTGAGGCAGACCCAGGATAGAATAG

- 133 -

AGAATAAAAGGAATGCTTATAGGAAACAATTTTGTATGGAATGCTAGATGGCCAAG
CCTCAGCCTTTGGTCCAGTGCAACCCTTGCCTCGCTTGTCAACAGTGAAAAATTAG
TTTGGTTAGAAGAACCATCTGGAAACACACCAGCTTCTGCTACCTTCATGCTCATT
GTTAAAAAAGATTAAACCAGTGTGAACATTCTGATCTGTTAATTCCAGGGACTGTT
5 TTTCTTTCCAATGGACTGTTTGTGGTAGAATAACCCCAAAAGCTCAAAGCTAAAA
TGCATCATCAGTCCTAGTCGGCAGTTCCTTAAGAATGGACTGGCGGCGTGGGTGAG
CTGATTTGGAAAACCTGCCCTTCTGCAAAAAACACTGGCCTGCTTTCCA

Sequence ID 337

10 CAAGACTCCATCTCAAAAAAAAAAAAAAAAAATCTACAGTGCTGAGTATATAAAATTAT
TAACACATTTTACAACAATATGTGTTTGTGGAGTTAAATATTTTTTGTCTTTAAAA
CAGGTAATTTTAGTGCATACTTAATTTGATGATTAAATATGGTAGAATTAAGCATT
TTAAATGTTAATGTTTGTACATTGTTCAAGAAATAAGTAGAAATATATTCCTTTG
TTTTTTTATTTAAATTTTTGTTCCTCTGTAAACTAAAAGAACACGAAGTAATTGGTC
15 ACAATTACTGGTGTTTAACTGCCAAATATGGGTAAATAAGGGAAAATTTTGTTTAA
TATTTAGTCCTTCTGAGATGGCTTGAATATTTGAATTTTGTGTACGTCTATACTG
GGTAGTCACAAGTCTTATAAACACTTTAGAGGAAAGATGGATTTTCACTCTGTATTT
TTAAACATCATTTATTTTAAATCTGGTGCTGAAAAATAAGAAAAAAATTAACTGC
ATTCTGCTGTTCTTCTTTAGAAGCATTCCTGCGTAAATACTGCTGTAATACTGTCA
20 TGCAAAGTGATCCTTTCTTGTCTATCCTTTTTTGGGGCAGTGTT

Sequence ID 338

CTGGACTGCATGACCAGATCTGATGGGTGAGACTCAGGTGGCATGGAAGAGCCGAA
AGAGGATACCATATGTGGGTGCCGGGGGGGATAGGTGAGAAGTACTAGAAGGCGGA
25 ATGGAAGGACACTTCTGCTCAGCTCTGTGACACGGGCAGGGACCCTGCAGGGCTCA
GGTCCTTTAACACAGCAGCTTCATTCTAACACCAGCAGCGTTGGAACACACGTACA
AGTATGCAGACTAAGCTCTTGCTTGGCTGATACGGCTTTTTTGGGTTTTTAGAGAAC
ATGCATATATGTTCTCATTCATGGTACATGAACTCAGAAGCCTTACTGCCTATTTT
TGTTAATACTTCTGGGCAAACATTACCACTTACAACCTCACACCAGTTAGAAATCAT
30 TTGTAAAATGTTATTTAATAAAGCCAAAGAACTAAATCATATTTATTTTCCAAGGN
TTTCTAAGATCTCTGAACTAATGAGGTTTTTTTAAATCCCCATTAAGTACTCATCA
CTGCTAGTAAAAGCAGTTGTCTTTACCTTTAATTCCAGTGAGTCCCCTTAAATTTA
TTTTTTATTATCTTTGGCTACATTGCCTTAGACAAAAATGTGGTCACCCTAATTTAA
NGGATAAAATTCACATCCTCACAGATTTCTTATTAAGAGGGTCTAANCCTTGAATA
35 ATCANCAGTGGAATGGAAGTCTTCTTTACTGGNTTNNATCCTTTCCCTTTTTTAT
CCCATG

- 134 -

Sequence ID 339

TTTTTTTTTTAAATAAAGCTGTCGGCACTCAAGGGTAATTTTCATATCAGTGTGNTCT
ACAAGCTGGGGGAAAATGAGTTCTAATTGTCTCANAGCTACCAAATCCTTCACCTTTA
GCATAAAGGTTTAAAGATATCACAAAGATGCCAAGTGATTAATAATGTTTTAAACC
5 ACCCCTTTTTCTGTCTGAAAAACAATAAACAATATTACAACAGTATAGTTACA
GAAGGGTTCTATTTTCATATGTTTTATGCACACTGTGCCTCAAAGGTACTATTTAA
ATATATATACTTTTGAGGGGGTGGCTAATGCAGAAACACCCAAGACCTAAGGAAGA
TACAACCCCATTTCTAGGTGTGAGGTCTAAATGCTTCACACACCCACTTGTGACCT
TTTTTCATGAAGAATCATAACACTGTGCAGTGAGAAACAGTGGCAAAGCAATACTG
10 AAAGCATTTTAAATTATTTACTAGGTTAAAGGGTGAAGTACTTTAAATACAT
CAAATTTTCATCAT

Sequence ID 360

GCAAGTGAGAGCCGGACGGGCACTGGGCGACTCTGTGCCTCGCTGAGGAAAAATAA
15 CTAAACATGGGCAAAGGAGATCCTAAGAAGCCGAGAGGCAAATGTCATCATATGC
ATTTTTTGTGCAAACTTGTCTGGGAGGAGCATAAGAAGAAGCACCCAGATGCTTCAG
TCAACTTCTCAGAGTTTTCTAAGAAGTGCTCAGAGAGGTGGAAGACCATGTCTGCT
AAAGAGAAAGGAAAATTTGAAGATATGGCAAAGCGGACAAGGCCCGTTATGAAAG
AGAAATGAAAACCTATATCCCTCCCAAAGGGGAGACAAAAAGAAGTTCAAGGATC
20 CCAATGCACCCAAGAGGCCTCCTTCGGCCTTCTTCCTCTTCTGCTCTGAGTATCGC
CCAAAAATCAAAGGAGAACATCCTGGCCTGTCCATTGGTGATGTTGCGAAGAACT
GGGAGAGATGTGGAATAACACTGCTGCAGATGACAAGCAGCCTTATGAAAAGAAGG
CTGCGAAGCTGAAGGAAAAATACGAAAAGGTA

25 Sequence ID - 361

nt: 622

CTGTNATNGAATCTGCTTGTNACTNAAATGCTAAACTCAATTCTGTAATTCAATAG
GTGCACCTNTCTGAGAAACATANNAGACAATGAGGAAAAGGATTCANCATTCCGTG
GAATTTGTACCATGATCAGTGTGAATCCCAN TGCGTAATCCAAGTAAGATGTTCA
CAAAGATTTGTTTTTAATGTCTAATTAATAAAAATTTTAAAGGAAGAAACATTCTAA
30 TACTTTAATTATAAAAAGTTAACTATTTTCAAAGGTATCAAATACAGTTAAACCT
TTAAAATGTATATTTCTTAATATCTTGAAATTGTAATGCCTTTTTTTTTTCTTAAA
TTTTTTTTTGTGTCATGAAATGAGATAGTAACAGCAGATTGGGACAACAAGGTTATATT
CTTGTCTTGAATCAGGCCATGGCTTCTTTCATCCAAATTTTCAAGCCTCATTTATTT
ACTTTGTCCCTGCCTCCCATCCCTGGATATCANGTTTGTGGATATCTACAGTTAAT
35 AGAGTGACCAAATAGTAGGAATACTGTCTCTCTATTTCTGAATAAAATACTTTGAAT
CAGATTTAGAAATAATGAATAAAATACAAATCACCATTGAAATTGCTCTAATTTTG
AGAGCT

- 135 -

Sequence ID - 363

nt: 628

ATCACNTGAGGCAAGAGTTTGAGCCAGCCTAGCTAACATGGTGAAACCCCATCTCT
ACAAAAATATAAAAATTAGCCTGGGTGGTGATGGGCACCTGTAACCCAGCTACTC
GGGAGGCTGAGGTAGGAGAATCACTTGAACCCGGGAGATGGAGGTTGCAGTGAGCC
5 AAGATCGTGCCACTGCACTCCAGCCTGTGTGACAGAACAAGACTCTGTCTCAAAA
AAAATAATAATAATAATAATAATAAAAAGGAATAACATAGCTAGGAATAAATTTAA
TCAAAGAGGTGAAAGACTTATACACTTAAACTACAAAAAAAAAATCACTGAAGGA
ATTATAGACCCAAATAAAAATAATAAAAAGACATTCTGTGTTTTAGGGAAAGAAG
ACTTAATATTGTTAAGATGTCAATACTACCCAAAGTGATCTACAGATTCAACATAA
10 TCCCTATCAAAATTCCAACAGCCTACTTTGTAGAAATGGAAAAGCCAATTTTCAA
TTCAGATGGAATTGCGAGGGGTTCTGAATAACAAAAACAATCTTGGGGAAAAAAA
CAAAAAACAAAGTCAAAGAACTCACACTTCTCTATTTATAAATTTACTACAAAGTT
ATAGTAATCAAA

15 Sequence ID - 364

nt: 528

TGAACATCCAGCCATGTCATTTCTTCCATTCCCTGCCCTGGAGTAAAGTAGATTTAC
TGAGCTGATGACTTGTGTGCATTTGTACATTGCAACCTTAGCTTACCTCTTGAAGC
ATGTAGAGCATTCAATCACCCACCATTCATTCACTGCCTACTCCCACCACAGCTGTT
TCGTGGTCTGTCTGCTCCCTGTGCCACCCCCACCCCATCAGGTGGGCCTTTTGCAA
20 GTGATGAAGTCACCTGTGGGGGAAGAGCTTTCCTTTCCTCTCCTCAACTCAGAAGG
CCTCTTCTCTTGCTCAAGAGGGTGCTGCTGCTTCTGCCTCCTTCCCCGGCCGGC
CTCCATCCCAGTTCACCTTTTCAGAAATGGCCCTCAGTCAACTCTTCCCTTTTCT
CCTGGCTTTTTATTCTCCAGTCTCTTAAGAGTATCCTTAGCTTTAAAAACAATA
ACACAGAGGATGGGTGCAGTGGCTCATGCCTGTAATCCCAGCACTTGGAGCCTGG
25 GGCGGGCGGATCACTTGAGGNCA

Sequence ID 365

GTCCCGGAATCGCGGCCGCGTGCACCTTTTCTATGCCTGCTATATAAACAGTACCT
TGCAAGATGTCCTGTCTGATATCCACAAAGGGGTATTGTCAACCCCAAGTTCAGAC
30 AGCTTTGTATTCTTCTGTCCCTGGATACATGAATTACTGCCATCTTTACACAGCGC
CCTAAAATACCAACGCGAAGTTACCTGCTCAGCTTGAAGCTGCGCTGTACCCTGGA
ACCAGCACTTCTGCTGAATGACTCAGGATGAAGCCTCGACTTCTCCTTCCCATCCC
ATGCCAGACCCAGTGGCTCCTTTCCCAATCTGATCCAGTGACTTTAAGTCCAGC
TGTTGCAACCTGGGCATGAGGAGGAGTGCAAGATGGCTTTGTCCTACCTGGAAAGA
35 GGCTTTCTGGA

Sequence ID 366

- 136 -

CACCATTTACACACAGTGGGTCCTTGAATAGCATCGTTTTATTCAATGTCATTTTG
TTATAACATTGAGAAAAAATTGATTCCCGGCTGGGGCCACTGTCTGTGCACCGT

Sequence ID - 368

nt: 329

5 GAAAGATCTAAAATCGACACCCTAACATCACAATTAAAAGAACTAGAGAAGCAAGA
GCAAATTCAAAAGCTAGCAGAAGGCAAGAAATAACTAAGATCAGAGCAGAGCTGAA
AGAGATAGAGACACAAAAACCATTCAAAAAAACAATGAATCCAGGAGTTTTTTT
TTTTAAAAGATCAACAGAATTGACAGACTGCTAGCAAGACTAATAAAGAAGAGAG
AAGCATCAAATAGACTCAATAAAAAATGATAAAGGGGATATCACCACCAATCCCAC
10 AGAAATACAACTACCATCAGAGAACACTATAAACACCTCTATGCAAAT

Sequence ID 369

GAAAGATCTAAAATCGACACCCTAACATCACAATTAAAAGAACTAGAGAAGCAAGA
GCAAATTCAAAAGCTAGCAGAAGGCAAGAAATAACTAAGATCAGAGCAGAGCTGAA
15 AGAGATAGAGACACAAAAACCATTCAAAAAAACAATGAATCCAGGAGTTTTTTT
TTTTAAAAGATCAACAGAATTGACAGACTGCTAGCAAGACTAATAAAGAAGAGAG
AAGCATCAAATAGACTCAATAAAAAATGATAAAGGGGATATCACCACCAATCCCAC
AGAAATACAACTACCATCAGAGAACACTATAAACACCTCTATGCAAATAAAGTAG
AAAAT

20

Sequence ID 370

GAAAGATCTAAAATCGACACCCTAACATCACAATTAAAAGAACTAGAGAAGCAAGA
GCAAATTCAAAAGCTAGCAGAAGGCAAGAAATAACTAAGATCAGAGCAGAGCTGAA
AGAGATAGAGACACAAAAACCATTCAAAAAAACAATGAATCCAGGAGTTTTTTT
25 TTTTAAAAGATCAACA

Sequence ID 371

GCCCCGAATCGCGGCCGCGTCGACGTAAGCTCGGCTGAATCCACGGTTCAAGAACA
GGAAAGAAGGCCAAGGCATAGGGAGTGGGGCAGTTGGGTGAATATTAGTACCTTTC
30 CCTCAGNTNCATTAATTACCCCTGCCTACTCTGCACAAAAGGATNTAACAACAGTT
TCCTTTTTTAATGGCCAGGTACAGCTGCTTATATGGANGGGCATTNTNAATGATAT
CCTTNATCACTGTCTTAATCATCATNCTTAAACAATCACTTTATTGTGTTAAG
GAAGATAAAAATGGCTGGGTTCAATTTCCGTTCTGGAAGAAATCGANTNAAAAGGT
AACCATTTAATAATGCANAGGGCANTTCACTGCAGACCCTAATACTGGAAATTTT
35 TAAAAACAAATGAAAACTTCTACTTTTTCTTCTAAGCTTACTTAACCACCCAAAT
TTTCCAGCCACATATCTTCCTAGTCTACAACCTGCCTTTAAGAGATGCTCA
AAAAAATGTAAATTCTCAAATACATTCTTATTACAATTACTGCTAACCT

- 137 -

Sequence ID 373

CCAGTGTGCTGGGATTACAGGCATGAGCCCTGCACCCAGCCTCTTAAACTGATCAT
ATGATATTGGTTCTCAACCAAGGGTGACTTTGCCCCCAGAGGATACTTGGCAATGT
CTGGAGATACTCAGTTGTCTGACTTTGGACAGGTGCTACTGTCACCCAGTGGGTAG
5 AGGTCAGGGATGGTGCTAAACATAGGACAGCTGTCAAGAGAAAAGAATGTACCCAG
CCCCAAATGTCAGTAGGGCTGAGGTTGAGAAACCCAGCTGTAGCTGACGTGTGAAG
GACAGACTGGCCTGGAAGTGTGTTTTCTGCCCTTTCCACCCCTGCATATTAGTTA
AGGCCAAAGGAAAAAAGGAATGCAGGAAATGCCCGTTAAAAATCTTCAAAACAATA
TAAATGATCAATTCCACTAAAACCTTTACACATTTAAGTATAAAGGTATTGGTA
10 GGAAAATTTGTTATTCACTGCTTTTCTCAGTGT CATGAAATAATTATTTCTGCTGT
CAGTTT

Sequence ID 374

AAAAAAAAAATCACTGAAGGAATTATAGACCCAAATAAAAATAAATAAAAAGACAT
15 TCTGTGTTTTAGGGAAAGAAGACTTAATATTGTTAAGATGTCAATACTACCCAAAG
TGATCTACAGATTCAACATAATCCCTATCAAATTCCAACAGCCTACTTTGTAGAA
ATGGAAAAGCCAATTTTCAAATTCAGATGGAATTGCGAGGGGTTCTGAATAACAAA
AACAACTCTTGGGGAAAAAAAACAAAAACAAAGTCAAAGAACTCACACTTCTCTAT
TTATAATTTACTACAAAGTTATAGTAATCAAAGTCGACGCGGCCGCGATTCCGGG

20

Sequence ID 378

CGACTGCGGCTCTTCCTCGGGCAGCGGAAGCGGCGGGCGGTCTGGAGAAGTGGCCT
AAAACCTTCGGCGTTGGGTGAAAGAAAATGGCCCGAACCAAGCAGACTGCTCGTAAG
TCCACCGGTGGGAAAGCCCCCGCAAACAGCTGGCCACGAAAGCCGCCAGGAAAAG
25 CGCTCCCTCTACCGGCGGGGTGAAGAAGCCTCATCGCTACAGGCCCGGGACCGTGG
CGCTTCGAGAGATTCTGTCGTTATCAGAAGTCGACCGAGCTGCTCATCCGGAAGCTG
CCCTTCCAGAGGTTGGTGAGGGANATCGCCCAGG

30 Sequence ID 380

GCAATTTAATTTTTAATAACAAAGATACTGTATTTTAACATGGTGAAATATACTTG
GCTAAGTCCAGATTAAAAAAAAGTATCTAGCCCAACAGTACAATTATACAGCT
TTGTACAGAACATTCCATAGATCAACAGAAAAATACATTTGAGCGCAAAAATAAAAA
ATATTTAAGGAGAATCTCTAAGCAGCATTTTATTTCTGCAAAAGACATATCTTGTC
35 TGATTAAATATCTACAAGTGCTTTTCCTTTCAAAAATACATATATTCTTAATAGAC
TAAGTCATTAACAATGACCTGGTAATTCTTTCACCTTCAATTTGAATGATTTATAAG
CTAAATCTTCAACCACAAAAGGTTTTTATTTGTATTAAAGATGTTACCACTTTTGA

- 138 -

CAAAAAGCTTAAAATATTTTATATTTCAAAGGAAAATTAGCAACATAACTTTACAA
TATATTCTATGATATTTTGATTGTGAGGGCTACTCTATTTAAAACTGATGATCTCT
GTTGTGTTGCTCAGATGCAGGAAAGCAGCAAAA

Sequence ID - 381 nt: 534

5 GACTTANATCTAAATGGACCACATTCTCTACTTAAAAAATGCTATTAACCATGTG
ATCTTCTCAGTCATGAGGTAATCTGGTGACTACCCCTTCCCTCAAAGCCAGTTGGGAT
ATTCTTTGAATAGAGTAAACAGTGTTTCTAGGCTGGGAGACACCAGACATAGTTG
AGGACAGAGGTGCTAGAAAATAGGAAGTTTAAAGCATGTGCGGTGATGCTCAGAG
10 GAGGTAAACCCACCCCTCATGCTCATAGCTTCCAATCATTTTCTCTAGTTCTTAAC
TCTTAAATGTGAGAAATGCTTGAAGATTCTAGTCATCTGAAGAAAGTCTCTTTATT
AAAGATTTTCATAAAAGAGACCAAAGCAGACAAACAGAAAAAGACATCTTGGGGAA
AAAAACAAGGATAATGGGAAGAGAAGGAAAGTTTTAAAAATTATCAATATCCTCAG
GGGGACAAAATATTATATCCTATAAAGACAGATTTTTATTTTTTAAAAAATAGAA
15 AGCAAAACAAGCTCCTAAAAATAAAGTTTG

Sequence ID - 382 nt: 444

20 GTTAAGGAAGTCAGCACTTACATTAAGAAAATTGGCTACAACCCCGACACAGTAGC
ATTTGTGCCAATTTCTGGTTGGAATGGTGACAACATGCTGGAGCCAAGTGCTAACA
TGCCCTTGGTTC AAGGGATGGAAAGTCACCCGTAAGGATGGCAATGCCAGTGGAACC
ACGCTGCTTGAGGCTCTGGACTGCATCCTACCACCAACTCGTCCAACTGACAAGCC
CTTGCGCCTGCCTCTCCAGGATGTCTACAAAATTTGGTGGTATTGGTACTGTTCCCTG
TTGGCCGAGTGGAGACTGGTGTCTCAAACCCGGTATGGTGGTCACCTTTGCTCCA
GTCAACGTTACAACGGAAGTAAATCTGTGAAATGCACCATGAAGCTTTGAGTGA
25 AGCTTTTTCCTGGGGACAATGTGGGCTTCAATGTCAAGAATGTGTCTGTCAAG

Sequence ID - 383 nt: 566

CTTTGAAGAACTTTGCCAAATACTTTCTTACCAATCTCATGAGGAGAGGGGAACATG
CTGAGAAACTGATGAAGCTGCAGAACCAACGAGGTGGCCGAATCTTCCTTCAGGAT
30 ATCAAGAAACCAGACTGTGATGACTGGGAGAGCGGGCTGAATGCAATGGAGTGTGC
ATTACATTTGGAAAAAATGTGAATCAGTCACTACTGGAAGTGCACAACTGGCCA
CTGACAAAATGACCCCCATTTGTGTGACTTCATTGAGACACATTACCTGAATGAG
CAGGTGAAAGCCATCAAAGAATTGGGTGACCACGTGACCAACTTGCGCAAGATGGG
AGCGCCCGAATCTGGCTTGGCGGAATATCTCTTTGACAAGCACACCTGGGAGACA
35 GTGATAATGAAAGCTAAGCCTCGGGCTAATTTCCCCATAGCCGTGGGGTGACTTCC
CTGGTCACCAAGGCAGTGCATGCATGTTGGGGTTTCCTTTACCTTTTCTATAAGTT
GTACCAAAACATCCACTTAAGTTCTTTGATTTGTCCATTTCCTTCAAATAAAGAAAT

- 139 -

TTGGTA

Sequence ID 384

TTTTGGGGTTTATATATAAGCCTGGTTCTTGCTGAAACTGCTTATGTTGATAACCA
5 GTTAGTGAGTTCCTCTCTATTGACTTGCTGGGAAGTTTATAGAGACATTTTTTATG
CATTCAGAGATTTTCAGTACAAATCTTGAAAAAGGGACATTTAGGCCGGGCGCGGTG
GCTCACATCTGTAACCCTAGCACTCTGGGAGGCTGAGGTGGGTGGATCATGAAGTC
AAGAGATAGAGACCATCCTGGCAAAAATTAGCTGGGCGTGGTGGGGTGCGCCCGTA
GTCCCAGCTACTCGGGAGGCTGAGGCAGGAGAATTGCTTGAGCCCGGGAGGCGGAG
10 GTTTCATTGAGCCGAGATAGTGCCACTGCACTCCAGCCTGGACAACAGAGCGAGAC
TGTGTCTT

Sequence ID 386

CTAAGGGTTTAAAGATGGAAAGAGGCATTGATGAACAGCTGGGGAAGGAGTAGTTT
15 GAGGTAGATGTGCAGATGGAATGAAGAGAAGGTCTCAAGAAGAGGGTGGAGCCAAA
GAGGGCTGCAGATTTAGAAGGCTAAAGTCTTTAGATGGCTTTGGATAGCCTGTTGT
ATCTTGACCATGCAGGTTACAGTGAGCATGGAGTGGGGACAGAAGTGGAGGAAG
GAACCAGGGAACATGGAGTGAGAAGCTAAAGGAAAGTGATGCAGTAGATACATGGC
TCTAAAGTACTCAGGACTTTCAGAGGCTTAAACATAGGGTGACCAACTATCCCACT
20 ATGCCTGATACTAAGGGCATTCCCTGGATGTGGACCTTTCATTCCCCAAATTAGGA
AAGTCTTGGGCATACCAAGACAAGTTGGCCACCCTACTCAAAAGTATGTAAGCTAA
CATATCTGTTCTCTAAGAGGTAAAGCTGGATGGGGATACCAGATGTATGTACGTG
ATGCAGTTAAACAGCAATACAAGGGGGCAAGTCTACCTGATCGGCCAATTCAATGG
GA

25

Sequence ID 387

GAAGCCAAACCAAAGGAGCTTCTACTTCATGATGCCATTTATGTAAAGTTCAGGCA
GAGAAAATCAGTGGTTTAAGAAGTTAGAATAATGATTATCTTTGGAGGGATTGCAA
CTGGAAGAAGTCATGATTGGGATTTCTGGGTCTAATAGTGCTCTGTGTCTTGATC
30 TGAGTGCCGACTACATGAGTGGTTAGGTTTGCAAAATTCATTGAGTTATGCACTTA
ATGGTGTGTCTTATTAGAGCTGATGGAGGAGAGAGGGCTTCAATTTGCACAACTG
AGTAATCAGCTAGGCCCAGTCACTAGGTGAACAACTTACTGCTCCAATCAGCCTTA
GAGCAGGAATCAAACCTCATGTCTCAGAAAAGTTATTAATTCAGCTTGTCTTGGGAC
TTCCTTCAGAGTCACTCTTGAATAGCTGAAATAGTAAATGTTAAATCTGTGGATGC
35 AAGTGTGTAAATTATTTTAGTCATCAGCTCTAATAAGATGGCCTTTGGGGAAATGA
GTATAAGGTCACGAAAATGAAATGGCAAGAAGGAGGTCTACTATTTCTTCTGTAAT
ACTGATTTTTACCCCATCAGGGTCAGTCCCAGAGGTGTAAATGTGAAGCTTG-T

- 140 -

CTTTTCTTTAATAA

Sequence ID 388

CTTTGGACACTAGGAAAAACCTTGTAGAGAGAGTAAAAATTTAACACCCATAGT
5 AGGCCTAAAAGCAGCCACCAATTAAGAAAGCGTTCAAGCTCAACACCCACTACCTA
AAAAATCCCAAACATATAACTGAACTCCTCACACCCAATTGGACCAATCTATCACC
CTATAGAAGAACTAATGTTAGTATAAGTAACATGAAAACATTCCTCCTCCGCATAAG
CCTGCGTCAGATTAAACACTGAACTGACAATTAACAGCCCAATATCTACAATCAA
CCAACAAGTCATTATTACCCTCACTGTCAACCCAACACAGGCATGCTCATAAGGAA
10 AGGTTAAAAAAAGTAAAAGGAACTCGGCAAATCTTACCCCGCCTGTTTACCAAAA
CATCACCTCTAGCATCACCAAGTATTAGAGGCACCGCCTGCCAGTGACACATGTTT
AACGGCCGCGGTACCCTAACCGTGCAAAGGTAGCATAATCACTTGTTCTTAAATTA
GGGACCTGTATGAATGGCTCCACGAGGGTTCAGCTGTCTTTACTTTTAACCAGTG
AAATTGACCTGCCCGTGAAGAGGGCGGCATAACACAGCAAGACGAGAAGACCCTAT
15 GGAGCTTTAATTTATTAATGCAAACAGTCCTAACAAACCCAGGTCCTAAACTCCA
AACCTGCATTAAA

Sequence ID 389

CGACCCGGAATTCGCGGCCGCGTCACTGAGTTCTTGACAAGAGTGTTTTTCCCTT
20 CCCGTACAGAGTGGGCCCAACGACCTACGGCACTTTGACCCCGAGTTTACCGAAG
AGCCTGTCCCAACTCCATTGGCAAGTCCCCTGACAGCGTCCTCGTCACAGCCAGC
GTCAAGGAAGCTGCCGAGGCTTTTCTAGGCTTTTCTCTATGCGCCTCCCACGGACTC
TTTCTCTGAACCCTGTTAGGGCTTGGTTTTAAAGGATTTTATGTGTGTTTCCGAA
TGTTTTAGTTAGCCTTTTGGTGGAGCCGCCAGCTGACAGGACATCTTACAAGAGAA
25 TTTGCACATCTCTGGAAGCTTAGCAATCTTATTGCACACTGTTTCGCTGGAAGCTTT
TTGAAGAGCACATTCTCCTCAGTGAGCTCATGAGGTTTTCATTTTTATTCTTCTCT
CCAACGTGGTGCTATCTCTGAAACGAGCGTTAGAGTGCCGCTTAGACGGAGGCAG
GAGTTTCGTTAGAAAGCGGACGCTGTTCT

30 Sequence ID - 390

nt: 523

GAATCCCTAGAAAAAGAGAATTCCCAACTTGATGAGGAAAACCTTAGAACTGCGAAG
GAATGTAGAATCTTTGAAGTGTGCAAGCATGAAAATGGCTCAGCTACAGCTAGAAA
ACAAAGAACTGGAAAGTGAAAAAGAGCAACTTAAGAAGGGTTTGGAGCTCCTGAAA
GCATCTTTCAAGAAAACAGAACGCTTAGAAGTTAGCTACCAGGGTTTAGATATAGA
35 AAATCAAAGACTGCAAAAACCTTTAGAGAACAGCAATAAAAAAATCCAGCAATTAG
AGAGTGAAGTACAAGACTTAGAGATGGAAAATCAAACATTGCAGAAAAACCTAGAA
GAACTAAAAATATCTAGCAAAAGACTAGAACAGCTGGAAAAAGAAAATAAATCATT

- 141 -

AGAGCAAGAGACTTCTCAACTGGAAAAGGATAAGAAACAATTGGAGAAGGAAAATA
AGAGACTCCGACANCAAGCAGAAATTAAAGATCCACATTTGAAGAAAATAATGTGA
AGATTGGAAATTTGGAAAA

5 Sequence ID - 391

nt: 566

CTTTGAAGAACTTTGCCAAATACTTTCTTACCAATCTCATGAGGAGAGGGAACATG
CTGAGAAACTGATGAAGCTGCAGAACCAACGAGGTGGCCGAATCTTCCTTCAGGAT
ATCAAGAAACCAGACTGTGATGACTGGGAGAGCGGGCTGAATGCAATGGAGTGTGC
ATTACATTTGGAAAAAATGTGAATCAGTCACTACTGGAACCTGCACAACTGGCCA
10 CTGACAAAAATGACCCCCATTTGTGTGACTTCATTGAGACACATTACCTGAATGAG
CAGGTGAAAGCCATCAAAGAATTGGGTGACCACGTGACCAACTTGCGCAAGATGGG
AGCGCCCGAATCTGGCTTGGCGGAATATCTCTTTGACAAGCACACCCTGGGAGACA
GTGATAATGAAAGCTAAGCCTCGGGCTAATTTCCCCATAGCCGTGGGGTGACTTCC
CTGGTCACCAAGGCAGTGCATGCATGTTGGGGTTTCCTTTACCTTTTCTATAAGTT
15 GTACCAAAACATCCACTTAAGTTCTTTGATTGTCCATTCCTTCAAATAAAGAAAT
TTGGTA

Sequence ID 394

GACCCGGAATCGCGGCCGCGTCGACCATTTTAGCCAAGGTGCCTCTATAGGGGTCA
20 AGACATCATGTGCCCAGACCTAAGGTGAGGAATGTCATATTTTTCTGTAAAATCA
TTTTATTTCTGTGTATCTTACCTTTAAATCATTGTGGTTTACTCTGAGATTCTGTA
GTCCTAATATTGTATCATTGTGCTGTCTGCAAAACAACCTTGAATCTATTTTGTTTG
CATCTTTTGTTACATGTAACGCAGCTGTACTTTATGTTCTTTGCAACTGTTTCCAT
TATGAGAACGCTGTGCTATTTACAAGGTTACATTTTTCTTGGCCAGGCGAGGTGGT
25 CATGCCTGTGATCCCAGCACTTTGGGAGGCCAAGGTGGGCGGATCACTTGAGGTAA
AGAGTTGAGACCAGCCTGGCTAGCATGGCGAAGCCCAGTCTCTACTAAAATACAA
AAATTGGCCGGGTGAAATTAGCCGGGCGTGTTGGTGTGTGCTTGTAATCCCAGCTA
CTCGGGAGGCTGAGGCAGGAGAATCGCTTGAATCCGGGAGGCAGAGGTTGCAGTGA
GCCAAGATCANGCCACTGCACTCCACCTCGGGGTCAAGAGCGAAACTCTGTCTCAA

30

Sequence ID 395

CCGTTTTAGTCAGGATGGTCTCGATCTCCTGACCTCGTGATCCGCCTGCCTCGGCC
TCCCAAAGTGCTGGGATTACAGGCGTGAGCCACCGCGCCCGGCGTAAATCAGGTTT
TTTAAATGTTTGCCAAACCTTATCACTGACTTTTATAACAAAATTATTTACTATAA
35 TCATTAGGGAATATTTAAGTTCTGCTAATACTTAAATATGCAGAGTGCTAAAACCA
GCAGTGAGTTTAGAATCAAGCTAAGCTTTATTGTTGCTACTATTTGAGGCATATTA
GTTGACTGGTGTTCATATGCAAGGCAGTCTACTGGGTGCAACAAGGGTTAGAAGGA

- 142 -

TATTTTAAAAAACTGACCCTATTCTCAGGATGAAAATAATACACTAGTAATAGTC
TGCTCTGTTGGTTAACTCCTCGTAAGGAGGTCAATTAAAATGCTGTAGTGTTGCAA
GGGAAGGAGAGGAAGAATCATATTCCTTCACTAGCAGGATCAAGAAAGCTTTTATA
GAAATATACAAAATCTTCACTTCTTGAAGGATTGGTAAAATTTAATAGCCAACATT
5 GGGCACTTATTCATTCTCTGAGTAAATATTTATTGCAT

Sequence ID 396

CTTAAATCTAAATGGACCACATTCTCTACTTAAAAAATGCTATTAACCATGTGAT
CTTCTCAGTCATGAGGTAATCTGGTGACTACCCTTCCTCAAAGCCAGTTGGGATAT
10 TCTTTGAATAGAGTAAACAGTGTTTCTAGGCTGGGAGACACCAGACATAGTTGAG
GACAGAGGTGCTAGAAAATAGGAAGTTTAAAAGCATGTGCGGTGATGCTCAGAGGA
GGTAAACCCCAACCCTCATGCTCATAGCTTCCAATCATTTTCTCTAGTTCTTAACTC
TTAAATGTGAGAAATGCTTGAAGATTACTAGTCATCTGAAGAAAGTCTCTTTATTA
AAGATTTTTCATAAAAGAGACCAAAGCAGACAAACAGAAAAAGACATCTTGGGGAAA
15 AAAACAAGGATAATGGGAAGAGAAGGAAAGTTTAAAAATTATCAATATCCTCAGG
GGGACAAAATATTATATCCTATAAAGACAGATTTTATTTTTTAAAAAATAGAAA
GCAAAACAAGCTCCTAAAAA

Sequence ID - 397

nt: 534

20 GACCCGGAATCGCGGCCGCGTCGACGGAAGCTCCTGCCCTCCTAAAGCTGAAGCC
AAAGCGAAGGCTTTAAAGGCCAAGAAGGCAGTGTTGAAAGGTGTCCACAGCCACAA
AAAGAAGGAGATCCGCACGTCACCCACCTTCCGGCGGCCGAAGACACTGCGACTCC
GGAGACAGCCCAAATATCCTCGGAAGAGCGCTCCCAGGAGAAACAAGCTTGACCAC
TATGCTATCATCAAGTTTCCGCTGACCACTGAGTCTGCCATGAAGAAGATAGAAGA
25 CAACAACACACTTGTGTTTATTGTGGATGTTAAAGCCAACAAGCACCAGATTAAAC
AGGCTGTGAAGAAGCTGTATGACATTGATGTGGCCAAGGTCAACACCCTGATTTCGG
CCTGATGGAGAGAAGAAGGCATATGTTTCGACTGGCTCCTGATTACGATGCTTTGGA
TGTTGCCAACAAAATTGGGATCATTTAAACTGAGTCCAGCTGCCTAATTCTGAATA
TATATATATATATATATCTTTTCACCATAA

30

Sequence ID - 398

nt: 512

GGGGAGCCCCCTCTTCCCTCAGTTGTTCCCTACTCAGACTGTTGCACTCTAAACCTA
GGGAGGTTGAAGAATGAGACCCTTAGGTTTAAACACGAATCCTGACACCACCATCT
ATAGGGTCCCAACTTGGTTATTGTAGGCAACCTTCCCTCTCTCCTTGGTGAAGAAC
35 ATCCCAAGCCAGAAAGAAGTTAACTACAGTGTTTTCCTTTGCACCGATCCCCACCC
CAATTCAATCCCGGAAGGGACTTACTTAGGAAACCCTTCTTTACTAGATATCCTGG
CCCCCTGGGCTTGTGAACACCTCCTAGCCACATCACTACAGTACAGTGAGTGACCC

- 143 -

CAGCCTCCTGCCTACCCCAAGATGCCCCCTCCCCACCCTGACCGTGCTAACTGTGTG
TACATATATATTCTACATATATGTATATTAAACTGCACTGCCATGTCTGCCCTTT
TTTGTGGTGTCTAGCATTAACCTATTGTCTAGGCCAAAGCGGGGGTGGGAGGGGAA
TGCCACAG

5

Sequence ID 399

TTTTGGCATTACTTAATCCAATTATAAAACTGAATTTTTTAAAAACAGCACTTGT
TTTTTCTTCCAAGATTAAATTTGAATTTTTTTATGGACATTAGAAAACATTGCAGTT
TAGTCATAATCAAAAATAAATCTTGAGGCTGGTAGAGCAGCTTTGTTGCTGTTTAT
10 ATTTTTATTGCTTACTGGATTTCAAGTGTACCTAGTGCCATCAGTTTGGTATTTTG
CCACCTTGCACATTCAAGTGTATTTGATTTTTCTTTTTCTTTTTTTCATATTACT
TTTAAATCCTGAATAGTTTGTGGCAGCTGGAGATCACCTAGTCCACCACTGTCCAA
CATGGCAATGGTAAGTAATATTGAGTAAAGAATAGAAAATTAGTAAAATGCATGGC
TTCAGAATTATAGCAATTTGCAAAATAGGTTAATGGATGAAAATTAGAATGACCAG
15 TTTAACTTTCCCCCAGCAGATTCTTCTGTAAACAATGCCCTTCAAAATAAAGG
AAGAACAAGTGGGTGTTATACCTATGTTATTTGGCTATGTTAGCACAATATGATGG
ACTAATTTGAGAAAAAGCATTACTTCCTTTACTATTACTTCTTTTCTTTATAGGG
CTAAGTCTGCCTTCTGGGTCTTTGAA

20

Sequence ID 400

GAAGAAGCGCGAAGAGCCGTTAGTCATGCCGGTGTGGTGGCGGCGGCGGAGACTGC
GGGCCCGTAGCTGGGCTCTGCGAGGTGCAAGAAAGCCTTTGAGGTGAAGGTGTATG
AAAGTCATCATAACAGATGTTTTCCAAAACTTGTAGAAGGTGTGAAAAAACTAC
TAGGATCACGCGGCATGTATTGAGCATATAGGTTGCTGTAGATGAATGTTCTTAGC
25 TGTATGTTTAAAAATACTTCTGCTTCGTTACCTCAAGTGTGGCATGCAGCATTTT
GGAAGGAAAATTGAAGACGTGTTCAAGAAAACATGAACAGAAGCAAATGATGAAAA
TGAGCATTTTACTTGATGTTGATAACATCACAATAAATTATGGAGAAAAATACATA
TTTGGCTAACTTTTAATTGCTGAACAATAAAGTGTTCCTTTTAAATCNAAAAA

30

Sequence ID 401

GAAGCCAAACCAAAGGGAGCTTCTACTTCATGATGCCATTTATGTAAAGTTCAGGC
AGAGAAAATCAGTGGTTTTAAGAAGTTAGAATAATGATTATCTTTGGAGGGATTGCA
ACTGGAAGAAGTCATGATTGGGATTTCTGGGTCTTAATAGTGCTCTGTGTCTTGAT
CTGAGTGCCGACTACATGAGTGGTTAGGTTTGCAAAATTCATTGAGTTATGCACTT
35 AATGGTGTGTCTTATTAGAGCTGATGGAGGAGAGAGGGCTTCAATTTGCACAACT
GAGTAATCAGCTAGGCCAGTCACTAGGTGAACAACCTTACTGCTACCAATCAGCCT
TAGAGCAGGAATCAAACCTCATGTCTCAGAAAAGTTATTAATTCAGCTTGTCTTGGG

- 144 -

ACTTCCTTCAGAGTCACTCTTGAATAGCTGAAATAGTAAATGTTAAATCTGTGGAT
GCAAGTGTGTAAATTATTTTAGTCATCAGCTCTAATAAGATGGCCTTTGGGGAAAT
GAGTATAAGGTCACGAAAATGAAATGGCAAGAAGGAGGTCTACTATTTCTTCTGTA
ATACTGATTTTTTACCCCATCAGGGTCAGTCCCCAAAGGTTGTAAATGTGAAGCTTG
5 GTCTTTTTCTTTA

Sequence ID 402

GACCCCTATTCTCAGGATGAAAATAATACACTAGTAATAGTCTGCTCTGTTGGTTAA
CTCCTCGTAAGGAGGTACAATTAAAATGCTGTAGTGTGCAAGGGAAGGAGAGGAA
10 GAATCATATTCCTTCACTAGCAGGATCAAGAAAGCTTTTATAGAAATATACAAAAT
CTTCACTTCTTGAAGGATTGGTAAAATTTAATAGCCAACATTGGGGCACTTATTCAT
TCTCTGAGTAAATATTTATTCATGCTTATCTTGTATCAACATTGNGATGAAAGCN
CAAGAATGAAAGAGGAGGGAGAATGTTTANAGAATAAGGCTGAAACACAGATTTTG
TAGGGAGCGTAGGGGAGACTGANAAAACAG

15

Sequence ID 403

AAGACACCTGATAGATTGTCTTGTATTATTTTTCTTTGCCTTCTTACAATCTCAG
TGATTAGAATTGGGCTGAAAACAATACATCAAATTCTCAGCAAATCCTTATGGGT
TGCTGGATACCGAGGGTTTTTAAGATCTTTAGACTTCACTATATAGAACAAATGTT
20 GAATGGGAATTTTCTTTATTTCTATANCGTTTNG

Sequence ID 405

CCCGGAATCGCGGCCGCGTCGACGATGAGCATTTTTTTCATGTGTCTTTTGGCTGCA
TAAATGTCTTCTTTTGAGAAGTGTGCGTTCATATCCTTTGCCCACTTTTGGATGGG
25 GTTGTTTTTTTCTTGTAATTTGTTTGAGTTCATTGTAGATTCTGGATATTAGCCC
TTTGTGATGAGTAGGTTGCGAAAATTTTCTCCCATTTTGTAGGTTGCCTGTTCA
CTCTGATGGTAGTTTCATTTGCTGTGCAGAAGCTCTTTAGTTTAATTAGATCCCAT
TTGTCAATTTTGGCTTTTGTGTCATTGCTTTTGGTGTTTGTAGACTTGAAGTCCTT
GCCCATGCCTATGTCCTGAATGGTAATGCCTAGGTTTTCTTCTAGGGTTTTGATGG
30 TTTTAGGTCTAACGTTTCAGTCTTTAATCCATCTTTTAAAAGTCTTTCACAGTAC
ATGAGTAGTAGTGACACCAATAATGTCAGAGCAGGGAACTCCCAGGTTCTGCCCAT
CCACAAAACAACAAATAAGCTGGCAAAAACCTTTAAGAATCAACTTTTGCAGATCT
CTGAAATCTAGTCAAACTTAAACAGAGGAAAGATTAATAAAGACNGGCTGCCTGA
GATAACACTAACACACAC

35

Sequence ID 406

CATCAAATAAATAAATAAATAAATTTTAAAAGTCACAGCATTGAATTTTTAAATGT

- 145 -

TTGGGATGATAAAGCACCTGCTTATCATGAAGCTANAGAAATTCAATGACACGTTT
GCCAGGGTCTTTGCTAGTGATGTTGGAACAAGTCTGTAATGCTGATGAAACATCAC
TGTTCCGGGCATTATTGCCCCAGAAAGACACTGACTGCAGCTGATGAAACAGCCCTT
CCAAGAATTAAGGATGCCAAAGACCAAATAACTGTGCTGAGATATACTTACGCAGC
5 AGGCATGCATAAGTGTAACCTTGCTGTTATAAGCAAAAGCTTGCGTTCTCACTGTT
TTCAAGGAGTGAATTTCATACCAATCCATTATTATGCTAATAAAAAGGCATGGATC
ACCAGGGACATCTTTTCAGATTGGTTTCACAAACATTTTGTACCAGCAGCTTGTGC
TTACTGCAGGGAAGCTGACTGGATGATGACTGCAAGATTTTGTATATCTTAACAA
CTGTTGTGCTCATCCTCCAGCTGAAATTCTCATCAAAAATAATGTTTATGGCTCAC
10 ACCTGTAATCTCAACACTTTGGGAGGATTGCCTGACCCAGGAGTTCAAGCCCACCC
TGGGCAACACAGCAAGACCCAACCTNTC

Sequence ID 407

TTTTAAAAATCATAAAACGTTTCTTACAAAAGAGCATTACATTNTGCACACTGCTC
15 TGAACAGATGCCAGGGACATGTGGACTATTGTTACTTTTCCTCCCTGTCCCACCCC
CCAAATGTTACAGTGACCACAAAGCAAGGTGTTCAACAATAATTACATGGGGGGAAT
TTTTTAAACCACCAACAATAACGAAAAATAAAATCCACTCACTCTGCTGCTGTTTC
AAAATTTCAATGTTAGTTTTTGCACGCCCTTCCCCCCCCCAACCCTGTTTGTAAAGG
AACTAAACATTACATCTGGTGAACAGCAAAGATTTCACTACACCTCAAATGCAGA
20 ACACCTATGAAGCAGAGGAATGTTGGCTTTTTTAAACAGAAGCAGATAAAAAAAAAA
GATGCAGGACTCCTTCAGTTCCTCACTAGTCTTAGAAAACTTTCCAGAATACTGC
TTCACACTATAAAAAAGAAAAAATATCTTGCATTAGAATCCTTCAACATCTGCATA
CTGCTTCACACTGTTTCGTTTCTAGGAGCACTTTGTGCACAGGACACTTCTGCTTATA
TTTCTTTAATCAGAACTTAGTTGGATGGGCCGGGCATGGTGGCTCACGCCTGTAAT
25 CCCAGCACTTTGGGAGGCCGAGG-GGGTGGATCACC

Sequence ID 408

CCATCTCCAAATTTAGTATTCATTCTGTTTAGCATATTATCAGTTGCCATCTATTT
GTTTTAACTGATTACTTGAATCTGATTAAACATCACAGAAATGGGCTTTGATAAGA
30 ACAATATTGAATAAGAAATTTTAAATAACAAAACAGCTTATAGAAAAATTCAGCAT
AACTTTTCCATCACCTTCACCACCCCTTGCCTTTTATTATCCTGTCTGTATCACTG
CTTTCTGTTAGCAGTGTTGTGTGAGTTAGGATTTGGGCAGGAAAGCAAAAGCAACC
ACCCGTCATTTTCCCAGAATGAAGGGTTTGACGTAGGATGTAGACTTTGTATAGTA
GTTGGGAGAGCTGTGGGAGTGAAGGTCAGGGATGTCACCTACAGAAGTCAGGGAAT
35 CTGCCACCAGAGATCCTGCATCAGAAACAGCCAACAGCGTGCTTCTGAAGAACTAG
TGGGGAAGTGGCTATAATTCTTAGGAATCCCAGCAAGTCCGCACCACTGTCTCAGT
CTACAGCAGTGGAGAAAGGGTTTCCAGGAGCTCTCTGGAAAGTTCCTGCCACAC

- 146 -

TTTGCAACAATCTTCAGAGGATAATGGGCTTCTCTTCCAGCTTCCACACCCAACAA
GAGTGCCTTTTCATCGGCCAACTCTAACCTGGAACCTATGGCAGAGGGGATTTAGG
AGACAGTTTGTNATGTCTGTGGAATGCAAATGAANANGTANCAATGCTTANTTGAC
AGCGGNCATACACAAATNTNGAAA

5

Sequence ID 409

GATCCGTNGACT

Sequence ID 410

10 CTCTTCCCAGCCCCTGAGCCCAGCCCCTTCCCAAGTGGTGCCAGACAAAAAACTAC
ATGGCCCTTTTCGTGTCTTGGGGGTGGAAAGGGAGGGATGAATTGGGGTGATAGAAC
CCTGGTGAATTCAGAGTAATCTTTCTTTAGAAAACTGGTGTCTTCTAAAGAAACAG
GATAGGAGTTTAGAGAAGGCACCAAAGCTTTCACCTTGGTTTGGCACCAGTTTCTA
ACCATCTGTTTTTTCTACCCTAGCTATCTTTTATTTGGTAAAATATAAATGTATAAT
15 TATGTTTGTAGAGCTTTACCAAGGAGTTTCCCTCCTTTTTTGTGTTGATTAGCA
AATTTTTGATTCTCCATTTTCCAAAAGTAAGAGACTCCAGCATGGCCTTCTGTTTG
CCCCGCAGTAAAGTAACTTCCATATAAAATGGTATTTGAAAGTGAGAGTTCATGAC
AACAGACCGTTTTCCATTTTCATCTGTATTTTATCTCCGTGACTCCACTTGTGGGTT
T

20

Sequence ID - 411

nt: 505

TGGAGCTGAAAAATTCCTATTACCTAGGGGCATCACAACGCATTGCATTTGCCCCG
TGTTTGGGATGATGCTGGTGTAAACCTACTATGCTGCCAGTCATGTAAAAGTATAG
CACACACAATTAGTAGGTAATGCTTGCAAATAATAATGAAAGACTCTGCTACTGGT
25 TTATGTATTTACTATGCTATACTTTTTTGTCTATTACTTTAGAGTGTACTCCTACTTT
TTTTTTTTTTTTTTTTTGGAGATGGAGTTTCACTCTTGTCTGTAGGCTGGAGCGAAN
TGGCGCGATCTCGGCTTACTGCAACCTCCACCTCCTGGGTTCAAGCGATTCTCCTG
CCTCANCTTCCCAGAGTAGCTGAGATTACAGGCATGCACCGCCACGCACGGGTAAT
TTTGTATTTTTTGGTAGAGACAGGGTTTACCATGTTGGCCAGGCTGGTCACCAACT
30 CCTGACCTCAGGTGACCCGCCTCCTCACCTCCAGAGTGTTGGGATTACAGGNGTGA
G

Sequence ID 412

ATAAAAAATTAGCTGGGGGTGATGGGCCCTGTACCCAGCTACTCGGGAGGTGAGGT
35 AGGAGAATCACTTGAACCCGGGAGATGGAGGTTGCAGTGAGCCAAGATCGTGCCAC
TGCACTCCAGCCTGTGTGACAGAACAAGACTCTGTCTCAAAAAAATAATAATAA
TAATAATAATAAAAAGGAATAACATAGCTAGGAATAAATTAATCAAAGAGGTGAA

- 147 -

AGACTTATACACTTAAACTACAAAAAAAATCACTGAAGGAATTATAGACCCAA
ATAAAAATAAATAAAAAGACATTCTGTGTTTTAGGGAAAGAAGACTTAATATTGTT
AAGATGTCAATACTACCCAAAGTGATCTACAGATTCAACATAATCCCTATCAAAT
TCCAACAGCCTACTTTGTAGAAATGGAAAAGCCAATTTTCAAATTCAGATGGAATT
5 GCGAGGGGTCTGAATAACAAAACACAATCTTGGGGAAAAAACAACAAAG
TCAAAGAACTCACACTTCTCTATTTATAATTTACTACAAAGTTATAGNATCAAAGT
CGACGCGCCGCGATCCGGGC

Sequence ID 413

10 CACAGTACTCCATTTTGGGGTCCAACTGTAATGCTCAAAATAATAAATGCTTACA
CGAAAATTATTTATTGAGAATATTCATATAAAAATTACCTAAAGCAAAGTAAAAAA
AGTAAAATCAAGGTGGTATATTTGAAGTGAATGGTGATTGGAAATTTTGTAGCTGTA
ACAAAAAGAAAGAAAACAACTTTTTTTAAAGCCTCATTCTCTTTTCTTTCAAATG
TACCTTATTTCCACACACTCTTGGGCTGACCTTTATTTTATCAATAAGCTCAATAT
15 TACTTTGTTTAAAATAAGATGCTTCAGCAAAAGTCATTCTCTCTTTAACCATATAA
TTTAAAAACTCCTCTTCACGATTGATAGCAAAATCAGAAACGTTAGGGCACCAGTG
AGTTGAAAAAACTGGTCTTAAGTTGAAAAACTATTATTAATAATATTATCCTATC
CATCCATATCTATTGAAATTGTGAGGTCCATAATTTTCATTTTAATTAATTATAGGA
AAGAAGAAAAGATAATACCCATTTGTTCTAT

20

Sequence ID 414

CTCAGACTCTTTCTGCCCTAATGGCCATTACTATCCAGTCTGTATTGCTACAAGGG
ACCCACTGGTACCCCTTTTAGATTCTATCAAAAGGAACAGGGTTTTCTAGAGGCA
GGCAGCCTGGTGGTATGGCACAGCAGAAGCTTACTGCTAATGAAATGGGAACCTCC
25 CCCTCCCTTGTGGTTTCAGCACAGAACCTGAATGCCAGGAAAAATTCCTGGGCCAA
GAAGCTAAAGCTAAAGAAACCTTCCTTTTTTTCAACGTTTTTTTTCTTTCAAACCTG
TAGGGTCACTTTTGATTGAGGCAAAGGGGTCTACTGTAAGTGGAAGAACTCACT
CCCCTAACATAAGTTTTTCACTGTGGTGGGATGGTGCCGCCGATATGCTTGATATG
CTTTTCCTTCCACATGTTAAGCTAGGAAACCTAACAGGATGTCAGCAGGGCAGTTA
30 ACTCTGGACTCANAGCCCTCAAGGGCATGTGGCANAACCTCATGGCATNCAAGACC
A

Sequence ID - 415

nt: 596

GTATAATTGATTCTTTTGAACCTAAAGTATAAGACTTCACGATTAGAAAAAATTA
35 TCCAAAGACTAATGTAATTAAGTGAGGAAAAGGTGCTGGAGGAACTGGATAACCAC
ATGGAAATGTATGAACCATGACCTCTATGTACATACTATATATAAACTTAATTT
GAGGTGTATCACAGAGCTAACTGTGGGGGCTAAACGTTGAAGCCTTTGGATGGCC

- 148 -

GCACAAGAGATGTCTGCATTTCATAACCTTG GGGAGGGTATGAACATTTCTTG GTAA
CATGCAAAAAGCACTAACTGTAAAAGAGAACAGTTGGTCAGTTGAATTCATGAAA
CATTGTAAACTTCTGCTAAACAACTGACACCATTAAAGAATGTGGAAAAAGGCTGGG
CACAGTGGCTCATGCCCTATAATCCCAGCATTTTGGGAGGCCGGGGCGGGAGAATCA
5 CTTGAGGCCAGGAGTTTGAAACCAGCCTGGGCAACATGGCAAGACCCCGACTCTAC
AAAAATATTTTTTAAAAATTAGTTGGGTGTGGTGATGCACTCCTGTAGTCCTAGCTG
CCAGGANGCTAAGGNGGAAGGATCACTTAACCCTGG

Sequence ID 416

10 CTGGTGGCGGCGGTCGTGCGGACGCAAACATGCAGATCTTTGTGAAGACCCTCACT
GGCAAAACCATCACCCCTTGAGGTCGAGCCCAGTGACACCATTGAGAATGTCAAAGC
CAAAATTCAAGACAAGGAGGGTATCCCACCTGACCAGCAGCGTCTGATATTTGCCG
GCAAACAGCTGGAGGATGGCCGCACTCTCTCAGACTACAACATCCAGAAAGAGTCC
ACCCTGCACCTGGTGTTGCGCCTGCGAGGTGGCATTATTGAGCCTTCTCTCCGCCA
15 GCTTGCCCAGAAATACAACCTGCGACAAGATGATCTGCCGCAAGTGCTATGCTCGCC
TTCACCCTCGTGCTGTCAACTGCCGCAAGAAGAAGTGTGGTCACACCAACAACCTG
CGTCCCAAGAAGAAGGTCAAATAAGGTTGTTCTTTCCTTGAAGGGCAGCCTCCTGC
CCAGGCCCCGTGGCCCTGGAGCCTCAATAAAGTGTCCCTTTCATTGACTGGAGCAG

20 Sequence ID 417

GCAGGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCCGCAGATAAGTTT
TTTTCTCTTTGAAAGATAGAGATTAATACAACCTACTTAAAAAATATAGTCAATAGG
TTACTAAGATATTGCTTAGCGTTAAGTTTTTAAACGTAATTTTAATAGCTTAAGATT
TTAAGAGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAAAAGG
25 TTTCTAAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAATGTAT
TTAAAAGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACCAATAGAAGG
GCAATGCTTTTAGATTAAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTTTAAA
AGTTGTAGGTGATTAAAATAATTTGAAGGCGATCTTTTAAAAAGAGATTAAACCGA
AGGTGATTAAAAGACCTTGAAATCCATGACGCAGGGAGAATTGCGTCATTTAAAGC
30 CTAGTTAACGCATTTACTAAACGCAGACCAAATGGAAAGATTAATTGGGAGTGGT
AGGA

Sequence ID 418

35 CCCGGAATCGCGGCCGCGTCGACGGGAGGTGATAGCATTGCTTTCGTGTAAATTAT
GTAATGCAAAATTTTTTTAATCTTCGCCTTAATACTTTTTTTATTTTGTTTTATTTT
GAATGATGAGCCTTCGTGCCCCCCTTCCCCCTTTTTTGTCCCCCAACTTGAGATG
TATGAAGGCTTTTGGTCTCCCTGGGAGTGGGTGGAGGCAGCCAGGGCTTACCTGTA

- 149 -

CACTGACTTGAGACCAGTTGAATAAAAGTGCACACCTTATAAAAAA

Sequence ID 419

5 CCCGGAATCGCGGCCGCGTCGACGGGAGGTGATAGCATTGCTTTCGTGTAAATTAT
GTAATGCAAAATTTTTTTAATCTTCGCCTTAATACTTTTTTATTTTGTTTTATTTT
GAATGATGAGCCTTCGTGCCCCCCTTCCCCCTTTTTTGTCCCCCAACTTGAGATG
TATGAAGGCTTTTGGTCTCCCTGGGAGTGGGTGGAGGCAGCCAGGGCTTACCTGTA
CACTGACTTGAGACCAGTTGAATAAAAGTGCACACCTTATAAAA

10 Sequence ID 420

CTTCATTTGAAATGGTTGAATCTGCTGTGTAATAAAGTGGTTCAACCATGATTAGG
AACTGAAATTTAGTAGAAGAGGGAAAAGGAGTTAATGTAACAAATTATTTTAGCTA
CAAACCCCGGTAATAGAGCACTTGGGGGATGGGATGGGGTGGGTGGTGAGACAAT
CAGAATGGTAAATTGATTAAATGCTCCTAACCTGTAATTTTGTGCATAGAGCACC
15 CTATGCTGTGGAAATAACTGTTCTTAGATTTTATTGTAACCTGGACTGTTTCAAGTTG
CCCAGAGGGAAAGAACATTCCTAATTCTAATAAAATAAACTTTTTATTTTGTTTA

Sequence ID 421

TGTCATTGAATCTGCTTGTACTTAAATGCTAAACTCAATTCTGTAATTCAATAGG
20 TGCACCTCTCTGAGAAACATAAGAGACAATGAGGAAAAGGATTGAGCATTCCGTGG
AATTTGTACCATGATCAGTGTGAATCCCAGTGGCGTAATCCAAGTAAGATGTTTAC
AAAGATTTGTTTTTAATGTCTAATTAATAAAATTTTAAAGGAAGAAACATTCTAAT
ACTTTAATTATAAAAAGTTAACTATTTTCAAAGGTATCAAATACAGTTAAACCTT
TAAAATGTATATTTCTTAATATCTTGAAATTGTAATGCCTTTTTTTTTTCTTAAAT
25 TTTTTTTGTTCATGAAATGAGATAGTAACAGCAGATTGGGACAACAAGGTTATATTC
TTGTCTTGAATCAGGCCATGGCTTCTTTTCATCCAAATTCAGACCTCATTTATTTA
CTTTGTCCCTGCCTCCCATCCCTGGATATCAGTTTGTGGATATCTACAGTTAATAG
AGTGACCAAATAGTAGGAATACTGTCTCTCTATTCTGAATAAAATCTTTGAATCAG
ATTTAGAAATAATGAATAAAATACAAATCAGCCATTGAAATTGCTCTAATTTTGAG
30 AGCTTATGATTTATTTCATCTTTGGTTTCCAAGTTCAAGTTATATGTAGACATTTTA
ATT

Sequence ID 422

GCTTCCTAGGTGAGGTCACGAGGAAACCTGCTGGCCAAGTGACCTGGCAGGGTGTG
35 GCCAGTGTGGCCAGGGCCGCCGAGCCTGCTTTCCTTCCCTGCAGCAGGAACCCTTC
TGGGGCTGTGATCCTGCGATGGTGCCTGGGTGGGAGTGGGGGTGGGGGGCGGGATG
GTCTCCCTACCTGCCAGCTTCTTGGTTTGAGGTGAGGACAGCCCCGGAAGCTCANA

- 150 -

CTTGGCTCCTGTCCATGTACTTGGGGCCATGAGCTCTGCAGGGACCTTGGAAAGAN
AGAGACGGGTGGTGTANGGCANGGGAAGGCATTGTCTTCAAACAGGAAAAAGCTGA
NAAATGGAAACAGGCGAAACTTACCAAGTGTAAACATCACCTGGAAGTGAAGGAGGGT
GGGAAAGGTTTTAATTATTTTAAAAATAGAGATGGGGTCTCACTATGTTGCCCAGGC
5 TGGTCTCAAACACTACTGGGCTCAAGTGAACCTCCTTCT

Sequence ID - 423

nt: 387

TGTTTTCTCNAGGGCGAGAGGCTGTCTTANAGCACCATTCTCTGGCCCTNGTCCCAT
GAGAAGGAACCGCACTCAGGAGCCACACTCTCCCACTNCCCTTGCCCANAAAGACTC
10 ACAGAGGGCACGGAGCTGGCTGTGGTGAGAGGAGGTCCANCAAATTCCTGTCTGCA
NAAGGGTTCTGAACACCACCGCCTGGCAGCGTGCTGGAGGAGGGATTCTCTTTTC
CTCACAGCAATTCTGACCAGAAACCTGTCAAATCAGGAATGGCTAAAATAAGACCA
GGGTATGAATGACCATCAGCCACAGTAAAACCAAGGCACAGCTCTCCTGAGCCAC
CCAAGCTGCTGTGGCCAGACTGGTGACATCACCTCAGGGCAAACAAAAA

15

Sequence ID - 424

nt: 420

CGCAGAATGGCTCCCGCAAAGAAGGGTGGCGAGAAGAAAAAGGGCCGTTCTGCCAT
CAACGAAGTGTAACCCGAGAATACACCATCAACATTCAAGCGCATCCATGGAG
TGGGCTTCAAGAAGCGTGACCTCGGGCACTCAAAGAGATTTCGAAATTTGCCATG
20 AAGGAGATGGGAACTCCAGATGTGCGCATTGACACCAGGCTCAACAAAGCTGTCTG
GGCCAAAGGAATAAGGAATGTGCCATACCGAATCCGTGTGCGGCTGTCCAGAAAC
GTAATGAGGATGAAGATTCACCAAATAAGCTATATACTTTGGTTACCTATGTACCT
GTTACCACCTTTCAAAAATCTACAGACAGTCAATGTGGATGAGAACTAATCGCTGAT
CGTCAGATCAAATAAAGTTATAAAATTG

25

Sequence ID 425

GGAAACTGATGCCAGTCAGAACTCAGATCAAATGAAGGGGTGAAGAGAACCAGAA
TTGATCTCTCTGTAGGAGAATATAAATGACTTTTTTAAAGTACATATTTTCTGTGA
AAGACAGTTTTTTGTTTAATGCAAAAATGTTAACAATGTTTATATCATGTAGAAGT
30 AAAAGATCGTGAAACAGCACAGAGAACAGTAGTAAGACAGATTGAATTGCACTGTT
GTAAGATGATGAACTTACAATATTAAGTGAAGGTAGACTGTGATAGATTAAGGATA
TATATTGTAATCCCTAGAGCAATTGTCAAAGTGGTACAGGTAAAAAGCCAATAGAG
GTGATAAAATGGAATACTAAAAAATATCAGATGAATAATAAAGAAGACAGGAAATG
AGGAACAGTGGAACAGAATGAATAAAAAACAAGACCATTAACTTAATCATTAATAA
35 TTACTTTAAATGGGTAAACATTATGGTTATAAGGCAGAGATTTTCAGACTAGATA
AAAGAGCAAGCTCCACTATATACTGTCTACAAGAGATATACTTTAAAGTGTATATT
ATATTTAAATATAAAGATTTGGAATAAATAAACCTAAGAATAAGCTTACTAGGGAA

- 151 -

GTGAAAGATCTGTACAACAAGAATTACAAAACACTGCTGAACGAAATCATAGGTGA
CCA

Sequence ID 426

5 GTCCCGGAATCGCGGCCGCGTCGACGTTTCCTCAAAATTTATCTTCCTGTTAATGT
CAGGCATGTATCTCCTTAGCTTGCCACAAATAACTATATATACCACAGACCTTCCT
TTGTAGGGCTAACAGTGTTCATTGTAAGTGGAGGCCTCATAGATACCTGGCCTTT
TCCTACCTTATTCCAAAGATGGTTGCATCTTATAAATAATGTCATTCTTCAGCAAA
TGGTATGGAAATGAGATTGTAATGTCATTATTTTCCTCTTTAAATAATCAGGACAAC
10 TCATGATACAAAGAGCTCTTCTCTATAAAAGGTGGGACTTTTTTTTTTTAGTAATAG
CAAAAATAAAATTGTACCTCCTTAATCTTCTACAGAAAGATGGATTTCATTTTCAA
CATTAAAGAGGTAGTTTTAAGAAGCAGTAGAAGTCAGCCTGGGCAGCATGGTGAAAC
CCCGTCTCTACAAAAAAGTTAGCTGGGCTTAGTAGTTGCAATCCCAGCTACTCTGG
AGGCTGAGGTTGGAGATCATCTGANCTGGGGAGGTCNAGGCTGCAATGATACANT
15 GAGCCCTGATTGTGCCACTCCACCTGGTTGCAGA

Sequence ID 427

TTCCAATCTTCGTGTTCACTTTAAGAACACTCGTGAAACTGCTCAGGCCATCAAGG
GTATGCATATACGAAAAGCCACGAAGTATCTGAAAGATGTCACCTTTACAGAAACAG
20 TGTGTACCATTCCGACGTTACAATGGTGGAGTTGGCAGGTGTGCGCAGGCCAAGCA
ATGGGGCTGGACACAAGGTCGGTGGCCCCAAAAAGAGTGCTGAATTTTTGCTGCACA
TGCTTAAAAACGCAGAGAGTAATGCTGAACTTAAGGGTTTAGATGTAGATTCTCTG
GTCATTGAGCATATCCAAGTGAACAAAGCACCTAAGATGCGCCGCCGGACCTACAG
AGCTCATGGTTCGGATTAACCCATACATGAGCTCTCCCTGCCACATTGAGATGATCC
25 TTACGGAAAAGGAACAGATTGTTCTTAAACCAGAAGAGGAGGTTGCCCAGAAGAAA
AAGATATCCAGAAGAACTGAAGAAACAAAACTTATGGCACGGGAGTAAATTCA
GCATTAAAATAAATGTAATTAAAAGG

30 Sequence ID 428

TGCAGGATCCGTCGACTCTAGATAACATGGCTAGAAAAGAGAATGAAAAAGTTGGA
ATTTTTAATTGCCATGGTATGGGGGGTAATCAGGTTTTCTCTTATACTGCCAACAA
AGAAATTAGAACAGATGACCTTTGCTTGGATGTTTCCAACTTAATGGCCAGTTA
CAATGCTCAAATGCCACCACCTAAAAGGCAACCACTCTGGGAGTATGACCCAGTG
35 AAATTAACCCTGCAGCATGTGAACAGTAATCAGTGCTGGATAAAGCCACAGAAGA
GGATAGCCAGGTGCCCAGCATTAGAGACTGCAATGGAAGTCGGTCCCAGCAGTGGC
TTCTTCGAAACGTCACCCTGCCAGAAATATTCTGAGACCAAATTT

- 152 -

Sequence ID - 429

nt: 535

5 CACAGTACTCCATTTTGGGGTCCAACTGTAATGCTCAAAATAATAAATGCTTACA
CGAAAATTATTTATTGAGAATATTCATATAAAAATTACCTAAAGCAAAGTAAAAAA
AGTAAAATCAAGGTGGTATATTTGAAGTGAATGGTGATTGGAAATTTTGTAGCTGTA
10 ACAAAAAGAAAGAAAACAACCTTTTTTTAAAGCCTCATTCTCTTTTCTTTCAAAATG
TACCTTATTCCCACACACTCTTGGGCTGACCTTTATTTTATCAATAAGCTCAATAT
TACTTTGTTTAAAATAAGATGCTTCAGCAAAGTCATTCTCTCTTTAACCATATAA
TTTAAAAACTCCTCTTCACGATTGATAGCAAAATCAGAAACGTTAGGGCACCAGTG
AGTTGAAAAAACTGGTCTTAAGTTGGAAAACTATTATTAATAATATTATCCTATC
15 CATCCATATCTATTGAAATTGTCAGGTCCATAATTTTCAATTTTAATTAATTATAGGA
AAGAAGAAAAGATAATACCCATTTGTTCTAT

Sequence ID 430

15 CAGGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCGCAGATAAGTTTTT
TTCTCTTTGAAAGATAGAGATTAATACAACCTACTTAAAAAATATAGTCAATAGGTT
ACTAAGATATTGCTTAGCGTTAAGTTTTTAACGTAATTTTAATAGCTTAAGATTTT
AAGAGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAAAAGGTT
TCTAAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAATGTATTT
AAAAGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACCAATAGAAGGGC
20 AATGCTTTTAGATTAAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTTTAAAAG
TTGTAGGTGATTAAAATAATTTGAAGGCGATCTTTTAAAAGAGATTAAACCGAAG
GTGATTAAAAGACCTTGAAATCCATGACGCAGGGAGAATTGCGTCATTTAAAGCCT
AGTTAACGCATTTACTAAACGCAGACGAAAATGGAAAGATTAATTGGGAGTGGTAG
GATGAAACAATTTGGAGAAGATAGAAGTTTGAAGTGGAAAACCTGGAAGACAGAAGT
25 ACC

Sequence ID 431

30 CGCTGGGTGCCTGCAGCGCCTCCCTTGTCTCATATGGTGTGTCCAGCACTCTATTG
TTGTAAACTGTTGNTTTGNCTGACCTAAATTNTCTTTACTAAACANATTTAATAGT
TNAAAAAANANCA

Sequence ID 432

35 TTTTAAAGTCATCTCTATAGGAAGGTGCTGGGCAGGGATCCCAGAGAAAGAAAGGG
TCCAAGACTCCATTAACTGCCCTGGATGAAGGGCACTGCTACAGCAGCTAGTACCA
GAGACTCTCCTATCTCACGGTTGAGGCAGACCCAGGATAGAATAGAGAATAAAAGG
AATGCTTATAGGAAACAATTTTGTATGGAATGCTAGATGGCCAAGCCTCAGCCTTT
GGTCCAGTGCAACCCTTGCCTCGCTTGTCAACAGTGAAAAATTAGTTTGGT'TAGAA

- 153 -

GAACCATCTGGAAACACACCAGCTTCTGCTACCTTCATGCTCATTGTTAAAAAAG
ATTAACCAGTGTGAACATTCTGATCTGTTAATTCCAGGGACTGTTTTCTTTCCAAT
GGACTGTTTGTGGTAGAATAACCCCCAAAAGCTCAAAGCTAAAATGCATCATCAG
TCCTAGTCGGCAGTTCCTTAAGAATGGACTGGCGGCGTGTTGAGCTGATATGGAA
5 AAGCTGCACCTTCCTGCAGAAGATCAACTGACCTGCTATCCCACCCCAAATTTCAA
CCTGAGGTATATTTCAATGAAGGCAGGTAGCTGTGCTTCTCAGAGCA

Sequence ID 433

TCCCGGAATCGCGGCCGCGTCGACCCGCCGCCGAGGATTCAGCAGCCTCCCCCTTG
10 AGCCCCCTCGCTTCCCGACGTTCCGTTCCCCCCTGCCCGCCTTCTCCCGCCACCGC
CGCCGCCGCCCTTCCGCAGGCCGTTTCCACCGAGGAAAAGGAATCGTATCGTATGTC
CGCTATCCAGAACCTCCACTCTTTCGACCCCTTTGCTGATGCAAGTAAGGGTGATG
ACCTGCTTCTGCTGGCACTGAGGATTATATCCATATAAGAATTCAACAGAGAAAC
GGCAGGAAGACCCTTACTACTGTCCAAGGGATCGCTGATGATTACGATAAAAAGAA
15 ACTAGTGAAGGCGTTTAAGAAAAAGTTTGCCTGCAATGGTACTGTAATTGAGCATC
CGGAATATGGAGAAGTAATTCAGCTACAGGGTGACCAACGCAAGAACATATGCCAG
TTCTCGTAGAGATTGGACTGGCTAAGGACGATCAGCTGAAGGTTTCATGGGTTTTTA
AGTGCTTGTGGCTCACTGAAGCTTAAGTGAGGATTTCTTGCAATGAGTAGAATTT
CCCTTCTCCCTTGTCACAGGTTTAAAAACCTCACAGCTTGTATAATGTAACCATT
20 TGGGGTCCGCTTTTAACTTGGACTAGTGTAACCTNCTTCATGCAATAAACTGAAAAG
ACCATGCTGCTANTC

Sequence ID 434

TTCGGACGCAAGAAGACAGCGACAGCTGTGGCGCACTGCAAACGCGGCAATGGTCT
25 CATCAAGGTGAACGGCGGCCCTGGAGATGATTGAGCCGCGCACGCTACAGTACA
AGCTGCTGGAGCCAGTTCTGCTTCTCGGCAAGGAGCGATTTGCTGGTGTAGACATC
CGTGTCCTGTAAAGGGTGGTGGTCACGTGGCCCANATTTATGCTATCCGTCAGTC
CATCTCCAAAGCCCTGGTGGCCTATTACCANAAATATGTGGATGAGGCTTCCAAGA
AGGAGATCAAAGACATCCTCATCCAGTATGACCGGACCCTGCTGGTAGCTGACCCT
30 CGTCGCTGCGAGTCCAAAAAGTTTGGAGGCCCTGGTGCCCGCGCTCGCTACCAGAA
ATCCTACCGATAAGCCCATCGTGACTCAAACTCACTTGTATAATAAACAGTTTTT
GAGGGATTTTAAAA

Sequence ID 435

CTGCAATGTGCAATAGTTGCACCACTGCACTCCAGCCTGGGTGACAGAGTGAGAAC
35 CTATCTCTTAAAAAAGGGAAGAGACATGAGAGGGGCCCAAGT
CACTTGCTCACTCACTTTCCGTGTACATGTACCAAGAAAAGGCCATGTGGGAAAGA

- 154 -

GCAAGAAGGCAGCCGCCTTCAAGACAGGAAGAGAGCCCTCACCAGAAACTGAGCCA
GAACCTTGGAATTCCAGCCTCCANAAGTGTGAGAAAAGAATTTTCTGTTGTTTCAG
TCCCCACACTATGGCATTGTTTACGGCAGCCTGAGCTAATACTCCTACTTTGTC
CTGCATTTACTTGGTCTTCCAGTTAGTTTTTTAGACTTTGGGAATCAGAGCAGTCA
5 GTTGTGTCAGATTTTAGCTTACAGTTGTCCTACCTGTGCAACTGAAATTTCTTCCATT
TTAAACCAGAGCAGAGTTTTAGAGTCAAAGAAACCAGATCTTTTAGTGCAGAAGC
TTTCCACTGTATTANAAGTGAGGAAGTTGGT

Sequence ID 436

10 AAAAAAACTCCAGAGAAGTTTATAGAAAGAGATGACATGTAAACCCTGCTGAAAAA
TAGTTTCATTTGTTAGAATATAATTGTCCTTCCACTAAAAAAGAAAAAAGCA
TTTAAGGCTCTAAGATCTCTTGAAGTACCACTTTTCCCTGAATCCAGAGTTTTTAT
GTGCATTATTTTTATGCGTTTGTAGTTTGATATGTTGTATTTATAAGTAGTTTTAG
CTTTCCATTATGAATTCTTCTTTGACCCATGAGTTATTTAGGTAAGTGTTAAAAA
15 TTTACAATAGTTTATATATGCAAATATTATGTTGTTAGAGTTGGTTTTTCATGTCAT
TTTTACATATACAGGGGCAGTTTCCCCAACTAAATTGTATATTCCTTAAAGCAGCA
CTCTTAAATTTTTATTTCTGTGTCAATTTCTTGNCTGTGTTTCCCTGGCATGGAATAC
ATGGCATAAAATTTGTTATGTAATTAAATGAAATATTATTATACCTTCTATTTTTT
AGAAAAAA

20

Sequence ID - 438

nt: 577

GTGACAGGGATGACATAACTATTAGTGGCAGGTTAGTTGTTGGTCACTTTCAACT
CTGGGTTCAAGCGATTCTCCTACCTCAGCCTCCCGAGTAGCTGGGATTACAGGCAT
25 GCACCGCCACACCTAATTTTCTATTCTTAGTAGAGACGGGGTTTCTCCCTGTTGGT
CAGGCTGGTCTCGAACTCCCGACCTCAGGTGATCTGCCTGCCTCAGTCTCCCAAAG
TCCTGGAACCACAGACATGAGCCACCACGCCTGGCCCCCTTTTAAATATTTCTGCT
CATTGATGATGCACCCAGTCACCCAAGTGCTCTGATGGAGATGTATAAGGAGATGA
ATGCTGTTTTTCATGGCTGCTAATACAACATTCATTCTGCAACCCCCAAATCAAGAA
30 GTAATTTTGAATTTCAAGTCTTATTATTTAAGAAATATATTTTGCAAGACTATAGC
TGCCATAGACCGTGATTCCCTCTGATGGATCAGACAACTAAAATGAAAACCTCCTG
CAACGTATTCATCATTCTAGATCCCTGAGGAATCGCCACACTGACTTNCACAATGG
GTGAACTGGGTTACAGT

35 Sequence ID - 441

nt: 552

AAACAAAATTATTCTCTGAGAGGGGAAAGGACATTTGAGGGAAACATCAAATTTCCC
CATAAATAAATGAATGGAGTTTGCAGGAAGGTGAGGGTGAGCAGAGATGTGTGTGG

- 155 -

ACATCTCTGACCATCCATCGCTGTATTCAAATGGATTGTTTTATTCCATTCTGGTC
TCAGGCATGACCACGTCCAGTGAAGACATTTGAGGCAGCACATCTCAGGACCCAGG
CAATAGACTGGCCCCAACTCAGGCTGGACTAAGGTGTGATTAATTCTTTGTTTTTT
GTGTGGAACAGCTCACCTTGT CAGACAGCCTCAGGGCATCTCTGAGACACAGGGGC
5 AGAAAATGACATTCATCTTTTGAGTCCTCATCCATGGAGTGCTGTGTTTGGGGGGC
TGCATCTGCTGAAGCGAGAACCCCATCTGCCACCCACCAGGATGCCCATTCCTCC
AGGACTTCTCCAACCTTACTATTAGACTAAACCAGAACAACTGTATTTA
TGCAAGCAAAATTGATGAGAAAATTATATTCAAATAAAGCAAAAATTA

```

10      Sequence ID - 442                               nt:       606
TCGTGCCACTGCACTCCAGCCTGGACGACAGAGTGAGACTCCATCTCAAATAAAT
AAATAAATAAATAAATAAATAAATAAATAAAAAATAAAAAATACTTCTGCTATGA
AAAACCTAGTTGGTATTTTTGCTTATTTTAATACTATAGAAATATGGTGATCTCATC
TTTAATAGAGTGCTTTTAAGGTCCCCAGTGATAATCTCCTAAAATCATGAACTTTA
15     AGAATTTATAATGTTAATATGAGGAAATGAAATCTGGATTATCTCACCATATTA
TATAATTCATTAGTGACAGAGCAAGAACTCCAGGTCACCTGTCTATTCCATGTTTT
TCCTATCTGCCTTTAAATGTTGAGATACTACCCTTATCTCATGTGAATGGAGAAC
TGCCTAAAATGCTAAAACTGACTCAGAGGCACCCAGACATAAGTGAAGTGTGATTA
GAAAATCCTGGTCAGTTGAGTCTTAGCCAAATGTGTACCTACTGTGTCTGCCTCTA
20     TCAAGTCAATGAAAACATGATCTGAGAACTGTAAGTCCATTTATGGAAAGGGTTGA
TTTANAGATATTTTGAAC TTNCAGTGATGAGCCCCTTCTCAAATAG

```

Sequence ID 446

CGGACTCCTGTGCTAATTGTCTAGCTTACATATCATTGTATAGAGACTGTTTATTCT
25 GTACCAAACCTGATTTCAAAAGTACTACATNGAAAATAAACCGGTGACTGTTTTTCT
TCATAAAGTTCTGCGTTTGGCATCTTCACTCTTCCAAAATGTATCTGTACATCAN
AAATGTCACTATTCCAAGTGTCTTTTTTAGTGTGGCTTTAGTATGGCTTCCTTTTAA
TATTGNACATACATTGNATCTTTGTTTTATGGNAATAAGTAATAAAAATGTAGACT
TCATATTTTGTACAAAATGTCTTATGTACAGAATAAAAAAGTTCATAGAAACAGCC
30 NANAA

Sequence ID 447

35 AGGCCGAGGCAGGCAGATCNCNTGAGGTCAAGAGTTTGAGACCAGCNTAGCTAACA
TGGTGAAACCCCATCTCTACAAAAATATA-AAAATTAGCCTGG-GTGGTGATGGGC
ACCTGTAACCCCAGCTACTCGGGAGGCTGAGGTAGGAGAATCACTTGAACCCGGA
GATGGAGGTTGCAGTGAGCCAAGATCGTGCCACTGCACTCCAGCCTGTGTGACAGA
ACAAGACTCTGTCTCAAAAAAAAAATAATAATAATAATAATAAAAAAGGAATAAC

- 156 -

ATAGCTAGGAATAAATTTAATCAAAGAGGTGAAAGACTTATACACTTAAAACTACA
AAAAAAAATCACTGAAGGAATTATAGACCCAAATAAAAAATAAATAAAAAGACATT
CTGTGTTTTAGGGAAAGAAGACTTAATATTGTTAAGATGTCAATACTACCCAAAGT
GATCTACAGATTCAACATAATCCCTATCAAAATTCCAACAGCCTACTTTGTAGAAA
5 TGGAAAAGCCAATTTTCAAATTCAGATGGAATTGCGAGGGGTTNTGAATAACAAAA
CACNATCTTGGGGAAAAAAAACAAAAACAAAGTCAAAGAACTCACACTTCTNTAT
TTATAAATTTACTACAAAGTTATAGTAATCNAA

Sequence ID - 448 nt: 329
10 TACGCACACGAGAACATGCCTCTCGCAAAGGATCTCCTTCATCCCTCTCCAGAAGA
GGAGAAGAGGAAACACAAGAAGAAACGCCTGGTGCAGAGCCCCAATTCCTACTTCA
TGGATGTGAAATGCCCAGGATGCTATAAAATCACCACGGTCTTTAGCCATGCACAA
ACGGTAGTTTTGTGTGTTGGCTGCTCCA₂CTGTCCTCTGCCAGCCTACAGGAGGAAA
AGCAAGGCTTACAGAAGGATGTTCTTCAGGAGGAAGCAGCACTAAAAGCACTCTG
15 AGTCAAGATGAGTGGGAAACCATCTCAATAAACACATTTTGGGTTAAAA

Sequence ID 450
GAGCAGTGGCATGATCACACCTTACTGCGGCCTCCAACCCCTGAGCTTAAGTGATT
CTCCCGCATTATCCTCCTGAGTAGCTGAGACTACAGGTGCATGCCACCATACTACTA
20 CTAAATTTGGGTGCGGTGGTGGTGGTGAATTTTTTAATATTTTGTAGAGACAGGGT
CTCACTGTGATGCCCAGGCTGGTCTTGA₂ACTCCTGGGCTCAAGCAGTCACCCACCT
CAGCCTCCCAAAGCACTGGGATTACAGGTGTGAGCCACCACACTGGCCAGCTTTGT
TTTGT₂TTTGTATGACTAAGCTGCTCTTGCTAAAAGGGCTTCTCTCTGA₂ACTTCCCTA
CCTTTCTTCTGTTTCCCTGGGCTAGGGCTCCATGTTGGCAGTCCTACTCCCAATTA
25 ACCTGGGGCTGTCTGGTTAACCTTTATAAGATCTGCAGTCATTGGGAGACCCGGGG
ACCAGGAATATTGTTGTTGAGGGAGCTACCCTGGAAAGTGGATGGGTGGCCAAAGG

Sequence ID 452
TTTGGCTTTGCCTCTAGGCATTAGATGTTATCTTTGGAGGCATCCTTCTATGAGCA
30 TTCATTTTTGGACCAAGCCTGGATTTACAATTCTATTACTGGCCCAGACTTCATTT
CTATCCAATTTCA₂TTCCACTGTGCTATAGTTTACAACATATAATTTGACTTATAAA
TAATTCCTGACTATGGGTTTAAAGACTGAAAATGGATCAATAGAAACTTTGAAAAT
GTTAACATCTTGATTGCTTTTCTCAGTGTAGAAATGGACAATGTTTAGCTTAAAAA
CTGCATGTTTTTAATGAGATACGGGGTTGAAAGACTTATTCTTGGAAATTTATTGTT
35 CTGGAGAAAGCCTGTTGCTATCTGCCATACCTTGGTTTACTTTGTGCAAAATGAGC
TTCTTTTAAAGTAATGAGCTCTTCCATGTTTCTAGCTTAAATTGCTGTCTTAGACAC
TTCATCAGGGTTCCCTGCTCTGCCTCATTCCCCCTTTTGTCTCACTTGCAGCCTTTG

- 157 -

ACATAATCCTGGGAGGCAATTGGCATCATACATATTTTGCTTTGTAATCTCCTGCT
TTGATTCTGACTGGGACCCAGC

Sequence ID - 453

nt: 747

5 GGATCTAAGACCAGCCTGGCAGCCACCAGATGGTGATTCTAGTCCTGGCTCAGTCA
GTAATAGGTCACTGACCCCAGAGAAATCAATTCAGCCTCCCCAGGTCCTTGGATT
CTTTCTGTGAAAATGAAAGCATAGGTAGGAATTTCCCATGGAACAGCTAGCAGAGG
AGAAATATTAAAAGTCAGGAGACTCATGCTATAGTTTTTCATACTTCATTACAACAA
TGTTGTTTAGGACAAGTGAGTTAACCTGTTAGCTTCCTCTATATAAAATGGAAAGT
10 CATTAAAAACCTACATAGCAGGGTTCTTGTGAAGATCAAGTGATAATGTAGGAAGC
ATGTACAAATGTCACATTCTGCCGTCACGTAATGGTCCTCACAGCTTGAGGTAGCA
TTTAGCATGTGTCATGATTTAGTACAAGGGTTGGCAAACCTGTTGCTCTTGGATTAA
GTCTGGCTCATTGCCTGTTTTTCAAAGAAAAAAATTGTATATGTGTGTATATATGT
TATATATAGGTACACACACATATGTGCTATATATAGCATATATACACACATAATAT
15 ATAAACATGTACATATATAGCATTATATATATACCGTGTATAATATCTCCAGTCCT
CATGACCAGCCATGCTTGTTTACATTTGCATACTCTATGATTGCTTTTCATGC
AACAAATGGCAGAGTTGAGTGATTGTTTTGCACAGANACTGTATGGCCCACTAAACC
TAAAATATTAATCTCTGCC

20 Sequence ID 454

CTCCTGCCGGGCTCGTGGCGGCTTCTGTCCGCTCCGCGGAGGGAAGCGCCTTCCCC
ACAGGACATCAATGCAAGCTTGAATAAGAAAAACAAATTCTTCCTCCTAAGCCATG
GCATATCAGTTATACAGAAATACTACTTTGGGAAACAGTCTTCAGGAGAGCCTAGA
TGAGCTCATAAGTCTCAACAGATCACCCCCCACTTGCCCTTCAAGTTCTACTTC
25 AGTTTGATAAGGCTATAAATGCAGCACTGGCTCAGAGGGTCAGGAACAGAGTCAAT
TTCAGGGGCTCTCTAAATACGTACAGATTCTGCGATAATGTGTGGACTTTTGTACT
GAATGATGTTGAATTCAGAGAGGTGACAGAACTTATTAAAGTGGATAAAGTGAAAA
TTGTAGCCTGTGATGGTAAAAATACTGGCTCCAATACTACAGAATGAATAGAAAA
ATATGACTTTTTTACACCATCTTCTGTTATTCAATTGCTTTTGAAGAGAAGCATAGA
30 AGAGACTTTTTTATTTATT

Sequence ID - 458

nt: 682

TGCCACTGAAGATCCTGGTGTGCGCCATGGGCGCGCCCGCCCGTTGTTACCGGT
ATTGTAAGAACAAGCCGTACCCAAAGTCTCGCTTCTGCCGAGGTGTCCCTGATGCC
35 AAGATTGCGATTTTTTGACCTGGGGCGGAAAAAGGCAAAAGTGGATGAGTTTCCGCT
TTGTGGCCACATGGTGTGAGATGAATATGAGCAGCTGTCCTCTGAAGCCCTGGAGG
CTGCCCCGAATTTGTGCCAATAAGTACATGGTAAAAAGTTGTGGCAAAGATGGCTTC

- 158 -

CATATCCGGGTGCGGCTCCACCCCTTCCACGTCATCCGCATCAACAAGATGTTGTC
CTGTGCTGGGGCTGACAGGCTCCAAACAGGCATGCGAGGTGCCTTTGGAAAGCCCC
AGGGCACTGTGGCCAGGGTTCACATTGGCCAAGTTATCATGTCCATCCGCACCAAG
CTGCAGAACAAGGAGCATGTGATTGAGGCCCTGCGCAGGGCCAAGTTCAAGTTTCT
5 GGCCGCAGAAGATCCACATCTCAAAGAAGTGGGGCTTCACCAAGTTCAATGCTGAT
GAATTTGAAGACATGGTGGCTGAAAAGCGGCTCATCCCANATGGCTGTGGGGTCAA
GTACATCCCCAATCGTGGCCCTCTGGACAAGTGGCGGCCCTGCACTCATGAAGGCT
TTCAATGTGC

10 Sequence ID 459

TCCCGGAATCGCGGCCGCGTCGACCTTGTCCTTGAGCGTCAACCTTCTTTCCCTGA
AGTGGCTGGGGTTCTGTTCCTTCTTTGATTGACAACTTGTGTTAACCTTCGCAC
ATCTCTGGGCCAATTTTGTCTTGTAAGTCTTTCCGGAGACCCCTGGAATTTAAATC
ATTAGCACCGCGCCCTTCCCCGAAGAGTCTTCGAAGGGTTGCCGCTTTTCGGTGGC
15 GCAGTTCTCGCGAGAAGGTGACTTTCTTTCTCGGTATTTCTTGTTTCCAGAATCC
TTAGCGCGAGGCGGAAAAATATTTCTCCAGCTTGTTGTTGATGCCGCGATTTTGA
CTGAGACTTCTTCCCACGATTTCTGTTTTTGTCTTCTCCAAGGAAAATGGCAGCTCC
CGAGCAGCCGCTTGCGATATCAAGGGGATGCACGAGCTCCTCCTCGCTTTCCCCGC
CTCGGGGCGACCGAACCCTTCTGGTCAGGCACCTGCCGGCTGAGCTTACTGCTGAG
20 GAGAAAGAGGACTTGCTGAAGTACTTCGGGGCTCAGTCTGTGCGGGTCTGTGAGA
TAAGGGGCGACTGAAACATACAGCTTTTGCCACATTCCCTAATGAAAAAGCAGCTN
TAAAGGCATTGACAACTNCATCAACTGAAACTTTTAGTCATACTTTAATCG

Sequence ID - 460

nt: 536

25 CAGAGATCAAAATAGGCCTTACACAGTGCACGCGAATTTAAAAGATTACCCCAT
CAGGTGTATGGATTTTGCAGTATTAAAGATGCTGCCTGGAATAGGTCATTATCTTC
TCCAAGTACTCTGTAAAGTCAATGAGTCACATAGAGTATAAGGTTTATTATCTGCT
TTTCTTTTCATTAAATAAATCTTTATTGAATTTCTACTACATTAAAAAACCAACCA
AAACAAAACAAACAAAAAAACACTTCCCTGAGCCATAAAGGAGAAGGTAGTTTGT
30 ACTGGAACCTTGAAGGATGGGTAACTTTTCCAGCAGATAAAGATTGAGAGAAGACCT
TCCAGGTAGAGAAAGCAGTGTGGGCACAGGCAAAGATGGAAGAACACACGTGGCTG
TGGGAAACACAGCTAGAAGCCAGTGCAGGATAGAGAGTAGGCTATGATGTGCAAAGG
TTANACACTGGGAGAGACAGGTCCATGAGAGTAGCTTGGACTAACACAGGGAGGGT
TTGGAATCCCAACTGGGGAACCTANAAATCAA

35

Sequence ID 461

TAGGAGGCTTATTCAGTATTTCCCCTATTCTCAGGCTACACCCTAGACCAAACCT

- 159 -

ACGCCAAAATCCATTTCACTATCATATTCATCGGCGTAAATCTAACTTTCTTCCCA
CAACACTTTCTCGGCCTATCCGGAATGCCCCGACGTTACTCGGACTACCCCGATGC
ATACACCACATGAAACATCCTATCATCTGTAGGCTCATTCATTTCTCTAACAGCAG
TAATATTAATAATTTTCATGATTTGAGAAGCCTTCGCTTCGAAGCGAAAAGTCCTA
5 ATAGTAGAAGAACCCTCCATAAACCTGGAGTGACTATATGGATGCCCCCACCCTA
CCACACATTCGAAGAACCCGTATACATAAAAT

Sequence ID 462

TCTTTATCAAGTTGAGAAAGTTCCTCCCCTCTATTCCCTAGTTTGCTAAGAGTCCTT
10 CTATCCTATTTCTTAATGGTTTAGTAGATGACTCTGTGGTACTTTGAAGGTTGTTT
GCAGAATTTCCATGCCATAGGCAATTTACCTTTCTTGACATTTGAAGGATTGATG
TTGGTGCCAAGTATAGAATCTTCACAGAGTCCTCCTGTAGCTTCTAAAGGTTTAGC
TTGAAAATGTTAATTGCTTAACGCTAGTAAGTGAGTGAAAAAGCTGGGGATAAATT
TTGTATCTTGCTTATATTTTCAGTTCCACCTCTGTCCNGACNAAACCCCATATAT
15 AA

Sequence ID 463

TAGTTTACATATCCCAACCTTTAAAAATATTCCTCTTATTAGCTTTATATTCACTT
TATAGAAGTTGAGTTTAAATTAATAATTCTTGGCATCCTGAAGTATGTCACATAGCA
20 TGTGCTCCTTATAAATATGTTGATATCTCAGAAGACAGCATCCCGGTTTTTCATTTT
ATAAAGTACCATACTTAAGAATGCTGTAATACTTATCTTTTATAACATGTTTCCTT
CGCTTTGCTTGNCCTTTTATGNCATCAGTTTTAACTGTTTACTTCATTTAACAGNTT
ACATCATNCAACAGTTTACTTCATTAAACAGTAGGTGGAAAAATAGATGCCAGTCT
ATGAAAATCTTCCCATCTATATCAAAATACTTTCAAGGATATACTTT

25

Sequence ID - 464

nt: 615

CGACTTTCAACCATCAAGTGAGGAATACCTTCACATAACTGAGCCTCCCTCTTTAT
CTCCTGACACAAAATTAGAACCTTCAGAAGATGATGGTAAACCTGAGTTATTAGAA
GAAATGGAAGCTTCTCCACAGAACTTATTGCTGTGGAAGGAAGTGAATTTCTCCA
30 AGATTTCCAAAACAAAACCTATGGTCAAGTTTCTGGAGAAGCAATCAAGATGTTTC
CCACCATTAACACCTGAGGCTGGAAGTGTATTACAAGTCCCGATGAAATTGAA
TTAGAAGGTGCTACACAGTGGCCACACTCTACTTCTGCTTCTGCCACCTATGGGGT
CGAGGCAGGTGTGGTGCCTTGGCTAAGTCCACAGACTTCTGAGAGGCCACGCTTT
CTTCTTCTCCAGAAATAAACCTGAACTCAAGCAGCTTTAATCAGAGGGCAGGAT
35 TCCACGATAGCAGCATCAGAACAGCAAGTGGCAGCGAGAATTCTTGATTCCAATGA
TCAGGCAACAGTAAACCTGTGGAATTTAATACTGAGGGTGCAACACCCCATTTTC
CCTTCTGGGAGACTTCTAATGAAACANATTTCTGATTGGCATTAATGAANAGTCA

- 160 -

Sequence ID 469

GATTTTAAAAATACATATAGCAAAAATATTACAGGGTCAGGGGAGACAATTAGAA
TGATATAATTCAAAGTGGATTAAAAAAAACCTGTCACCCAGAATAACAATACCCAG
CAAAGTTGTCCTTCATAAATGAAAGAAAATNAAATCTTTNCCNAACNA

5

Sequence ID 471

TCCCGGGAATCTGCAGGATCCGTCGACT

Sequence ID 472

10 GACAGTGCCCAGGGCTCTGATATGTCTNTCACANCTTGNAAGTGTGAGACAGCTG
CCTTGTGTGGGACTGAAAGGCAAGATTTGTTCTTGCCCTTCCCTTTGTGACTTGAA
GAACCCTGACTTTGTTTCTGCAAAGGCACCTGCATGTGTCTGTGTTCTTGTAGGCA
TAATGTGAGGAGGTGGGGANACCACCCACCCCCATGTCCACCATGACCCTCTTNC
CACNCTNACCTGTGCTCCCTCCCCAATCATNTTT

15

Sequence ID - 473

nt: 694

TGGGCTTTGGGCTGGCTGCAGTCTGTCTGAGGGCGGCCGAAGTGGCTGGCTCATTT
AAGATGAGGCTTCTGCTGCTTCTCCTAGNGGCGGCGTCTGCGATGGTCCGGAGCGA
GGCCTCGGCCAATCTGGGCGGCGTGCCAGCAAGAGATTAAAGATGCAGTACGCCA
20 CGGGGCCGCTGCTCAAGTTCCAGATTTGTGTTTCTTGAGGTTATAGGCGGGTGT
GAGGAGTACATGCGGGTTATTAGCCAGCGGTACCCAGACATCCGCATTGAAGGAGA
GAATTACCTCCCTCAACCAATATATAGACACATAGCATCTTTCCTGTGCTCTTCA
AACTAGTATTAATAGGCTTAATAATTGTTGGCAAGGATCCTTTTGCTTTCTTTGGC
ATGCAAGCTCCTAGCATCTGGCAGTGGGGCCAAGAAAATAAGGTTTATGCATGTAT
25 GATGGTTTTCTTCTTGAGCAACATGATTGAGAACCAGTGTATGTCAACAGGTGCAT
TTGAGATAACTTTAAATGATGTACCTGTGTGGTCTAAGCTGGAATCTGGTCACCTT
CCATCCATGCAACAACTTGTTCAAATCTTGACAATGAAATGAACTCAATGTGCA
TATGGGATTCAATCCCCACCATCGATCATAGCACCCCCTATCAGCACTGNAACTC
TTTTGCATTAAGGGATCATTGC

30

Sequence ID 474

GGCAGCGCGGGGAGCCCGTCCGCGCCGGCGGGCGGGCCGGTTTCGAAGTTGATGCA
ATCGGTTTTAAACATGGCTGAACGCGTGTGTACACGGGACTGACGCAACCCACGTGT
AACTGTCAGCCGGGCCCTGAGTAATCGCTTAAAGATGTTCTACGGGCTTGTGCT
35 GTTGATGTTTTGTTTTGTTTTGTTTTTGGTCTTTTTTTGTATTATAAAAAATAAT
CTATTTCTATGAGAAAAGAGGCGTCTGTATATTTTGGGAATCTTTTCCGTTTCAAG
CATTAAAGAACACTTTTAATAAACTTTTTTTTGATAATGGTTAAAAA

- 161 -

A

Sequence ID 475

5 CATAATAAAAAACAATCAACAAACAGGGAATGGAAAGAACTTCCTCAGCATGGTG
AAGGCCACATATGAAAATCCCACAGCTAACATCATACTCAATGATGAAAGACTGAA
AGCTTTTCTCCTGAGATCAGGAACAAGACAAAGATGTCACCTTTTGTCACTTCTAT
TCAACTCATATTATTGGAAGTTTTTGCCAGAGCAATTAGGTAAG

10 Sequence ID - 476

nt: 476

GAGAACTCTTTTCATAGGCTGAATGTTGCTCCACAATGTGTCCTTTGACTATCTCTG
GCTAATTATTATTTTAAATCTCTTCTCAGCTTTTCCAAGAACATAACGTTAACCAA
GATCTTAGGCCATTCACAACTCTTTTGTAAAAATTAATGTGGATGTGAAACGAGGC
AACAAATCCTGAAGTAGAAAGTTATTCCTGGCCAGGCACGGTGGCTCACGCCTGTA
15 ATCCTGGCACTTTGGGAGGCCGAGGTGGGTGGATCATGAGGACAGGAGATCGAGAC
CATCCTGGCCAACATGATGAAACCCCATCTCTACTAAAATACAAAAAATTAGCTGG
GCATGGTGACGCGTGCCTGTAGTCCCAGTTACTCGGGAGGCTGAGGCAGGGGAATT
GCTTGAACCTCGGAGGTGGGAGGTTGCAGTGTGCCGAGATCACGCTACTGCACTCC
AGCCTGGCAACAGAGCAAGACTCCATCT

20

Sequence ID 477

AAACAGAAAGTTTCTTCTAAAGGCATGATTCAAGTTAAGTCATTCTTAAGTGTTAAA
AAATTGTGAAAAATGTGCCTGTAATCCCAACACTTTGGGAGGCCGAGGCAGGCAGA
TCACGAGGTCAGGAGATCAAGACCATCCTGGCTAACAAGGTGAAACCCCGTCTCTA
25 CGAAAAATACCAAAAACATTAGCCGGGCGTGGTTGTGGGCGCCTGTAGTCCCAGCT
ACTTGAGAGGCTGAGGCAGGAGAATG

Sequence ID 478

TTCTTGGGATATTGATGACTACTGTCTGAGAGGTGCTGTGGGGAGATTTTCAGGAT
30 TGTGTGGTCTTTGAGGGGGGTGTTTTTTTAAAGACAACATTGACCACTGTCCACTGT
CCACATGATCATTGTAAAATTGCAATGCCGCATGCTAGTTGGTTACATAAGACATA
ATTCCAGTGATTGAAGGTGGTTACACTGTATGGTGGTGTGTTCAAGATGGCACTGG
CATCTTTGAGCAGAGCCTGGCTATGCAGCATCATTGAGTTTTTTTAAACACCCTAN
AGGTCTGGTTGTTGTTGCTGTTGTCTTTTCTGTGAAAGTCACAAANANAAGTTACA
35 GTCCAGGTGAACCTGGAGTTTATAGGTTGGTTTTGTTTCTGNTATATATATATATA
TATATATTTTTTTTTTTTTTTTAAACATTTACCTGTAGTGCTGTAGCTGTTGATACTA
TCACCTGCATGCTATTTCTAGTGAGTGCTAAATACAGTATGGTCCAATGACAATAA

- 162 -

CAGCCCATGGTACTGCCAG

Sequence ID 479

5 CATCAGTCTGTTATCCATGCTGACTTTCGGAAGACTTGCAGCTACTGCATTGATAT
CTTTCCTGCCAATAAGCAAAGTGTTGAACACTTTCACAAAATATTTTACTGAGGCAG
GCTTGAAAGAGCTTTCAGAATATGTTTCGGAATCAGCAAACCATCGGAGCTCGTAAG
GAGCTCCAGAAAGAACTTCAAGAACAGATGTCCCGTGGTGATCCATTTAAGGATAT
AATTTTATATGTCAAGGAGGAGATGAAAAAAAAACAACATCCCAGAGCCAGTTGTCA
TCGGAATAGTCTGGTCAAGTGTAATGAGCACTGTGGAATGGAACAAAAAGAGGAG
10 CTTGTAGCAGAGCAAGCCATCAAGCACTTGAAGCAATACAGCCCTCTACTTGCTGC
CTTTACTACTCAAGGTCAGTCTGAGCTGACTCTGTTACTGAAGATTAGGGAGTATT
GCTATGACAACATTCATTTTCATGAAAGCCTTCCANAAAA

Sequence ID 481

15 CACACTTTCATGATAAAAACAGAACCTAGGAATGAAAAGAAATTATAGCAACATAA
TAAAGACCATATATGAGAAGCCCACAGCTAACATACTGTATGGTGAAAACTGAAA
GCTCTTCTCTAAGATCAGGAACAAGGCAAGGATGCCCATTTCTTGCCACTTCTATC
GAACGTAGTACTGGAAGCCCTAGCCAGAACAACCTAGGCAATAGAAAGAAATTAAAG
GCATCCATNTCAGAAAGGAAGAANCAAAATGCTGTCTGTTTAAANATGACA

20

Sequence ID 482

TTTCTATANAAAAAAATTTTTTAAATAATTGTAAAGTTAGATTTTAAATTTGTAA
ATATAAAATCACAAAGGAATGTACCCAATAAAATGTAAATGCNCCATAAAAAAA
AAAAAAAAAAAAAAAAAAAA

25

Sequence ID 483

CGNTAACGTGCAATCCGCCGCACGCCAGCAAACCTGGACAAACTCCGGGATCTCATC
GAAGCGATTGAGCACCAGTACCAGAGTAATACCGGACTGATGTAACGAGGCGAGTC
GCTCATCCAGCTTGCTGACGTGAGGCAACATCCAGGCCATCGAACGGNTCATCAAG
30 AATCAACAAGTCAGGCTCCGACATCAGCGCCTGACACAGCAGGGTTTTTCGCGTCT
CGCCAGTGGAAGGTATTTAAAGCGTCNGTCGAGGAGGGCGGTAATACCGAACTGC
TGCGCCAGTTGCATGCAACGCGGTGCATCCTTTACTTTCATCCTGAATGATCTCAGC
CGTAGTGCGTCCGGTGCCATCTTCGCCAGGGCCGAGCATATCGGTGTTATTCGCT
GCCATTCGTCGCTGACGAGTTTTTGCAATTGCTCGAAGGAGAGACGAGTGATGTGG
35 GAAACTGGCTTTGCCGTTACCTTTCAAAGCGGGAAGTTCCCCCGCCAGCGCGC
GGGCCAGGGCCCGAT

- 163 -

Sequence ID 484

TTTTTTTTTTTATTCTATTAAAAAATGTTNNTGAAAAAAGATACTTAAATTTTAA
AGATAACTNAATTCCTAANGATTTAAAATAATCCAAGCAGAGATGAAAGANCAAAT
GCAAATGCNTAAAAAGACCCCANAGCATTGTTAGCAAAAAGCAAATATAGTTAGCC
5 AAGCATATATATNTCATAAAAGCAATAANAAGGCNTAAAGCAAGTTTGGGGAGAGC
TTATTTAAAACTTGTA AAAATCATTG AATTTTAAAGTTTTC AAAC

Sequence ID - 485

nt: 551

10 TTTGGAACACAAAGTTCCCTTTT TAGAAGAATAGGTATTGAGCCCTTGAGCGTGGG
TAGAAAGATAGAGACAGAGTGATTTGCAAAATAATGGAGGATCATATTTATATATG
AATTTTCACTTATTTGAACTTTCAGATATCANCTTNAAAANCTTTGGTTTAAGTAA
AGTNTNTTAATGAGACTCCTTGGATGAAAGTAACCAAACCAGTAAAAATAAGGTA
ATAAGGATGTAATAGTTTCTTATGGACACTCAACAGCTAGAATGCAGTTAGTCTCA
GAAAAGAATTAGAACAAATAACTGGAAGGCCATCAGGAGTCCAAAACCATCACTCT
15 TTTATATTTTATATTTTATTTTCTCTCTTCANATGAGCATTCTCTTTCTATGTCC
ATATGGTANAAGGCGGCAGCTCCATAGATTATGGCTTCAGATGTTACAGTTCCGCT
NAATGCAGGGACAGACTTGCTATCTTTCAGTCCCCTTACATATCCTGGGGAGAGAG
CAAATGATTGACTGGCTTGAGTCAGGTGCCCGTTCCCTTTCCAATCT

20 Sequence ID - 487 nt:224

GTTTGNTTGTGACCATCTGTACTTGTAATTTCTTTACNTTCATTGGTATGAAAAAT
ATGTTCTTAGAAGCANGAAAAAGAATTCAGNTTTGCTTTGTATACTAAATTAAATG
CTGTAATTTTGATAAAATGAAAAATCTGCTTTATTTGCAACAATTGGTTTCTTCCT
TGACGTCAGCCTCACTCTTGGA CTTTGGTATT CAGCCNGNCACCCCTGGGAATTCC

25

Sequence ID - 488

nt: 349

GTGCCTCCCTGTGTGAGTAGCCTAAGGTGCATTGAAAAAGACTGGGATGTGTTTTA
TTTTTTTGTATTAGATAGCATTAAACCTTACTGTTGAAGTATTTTTGGTGAGTATT
AGTGACAAGCCATTGAGTCTTAAGCCTTACGGCTTCCTATAAAATCACTAATTTTCG
30 TGTGTGTTTGTGTGTAGGTTACGTTATATATAGGATTCGTGTTCCGCCGTGGTGGCC
GAAAACGCCCAGTTCCTAAGGGTGCAACTTACGGCAAGCCTGTCCATCATGGTGTT
AACCAGCTAAAGTTTGCTCGAAGCCTTCAGTCCGTTGCAGAGGANCGAGCTGGACN
CCCTGGGGGGCTC

35 Sequence ID 489

TTAACAGCTGCATAGAGTTTTAAAAGTACATTATATTTTGTGAGACAAGTAAAATA
TCTGTTTTTTCACGCAAAAAAGCCATGAAATACGTAATTTTTTAAAGACAAAAAAT

- 164 -

CATCTTTTGAGTTTGCTCTTTGGTTTTCTTCATTCCTTTGAGGATTGGGAAAAC
AGAAAGATTCTTTGATTGTTGGGTAATGAAGAGGTAATTTGGGACAGTGTGGTGGTAC
CAGGAAGAAAGAGGATTGGAAAGGCCAGTACTGTTTTAGTTGCTCGGCACTGTTGG
TTTTGTTTTAATGTGGTTGCCCTGTCCACTACATGGTTCTATCAGTAGTGTAATCC
5 ATTTTCAATGTAAAGCTCTTTTAGTTTTTGTTCATAGACATAAATTAATATTTTGAG
AGGCATCCCTCACCTGTTTCATTTCTTCTGTGTTGAAATGAAGTACTTAAAATTACC
GTTATACATGAACTTTGTGGACTGTAAGATTGTTTATATATGTTCAAATGCCTTTT
AGCTGGCTTTTTAATTAATATGCCTGTTTTGAGTGCTTAATACAATGTAATGNGGA
TTGTAAATCATACCTATTTTAAATCATTCCTTCCTGTATATTTGNACTCAGAGAGC
10 CTTATTTTATTCTTCCAGC

Sequence ID - 491 nt: 382
TTTTCTTAGAACTTTATTTTTCTGGCCAGGCGCAGTGGCTCACACCTGTAATCCC
AGCACTTTGGGAGGCCAAGGCAGGTCGATCACCTGAGGTGAGGAGCTCAAGACCAG
15 CCTGGCCAACATGGTGAAACCCTGTCTCTACTAAAAATACAAAATTAGCTGGGCG
TGGTGGCGCATGCCTGTAATCCCANCTACTCAGGAGGCTGAGGCAGGAGAATTGTT
TGAACCCGGGAGGCGGAGGTTGCANTGAGCCGAGATTGCGCCACTGCACTCCAGCC
TGGGCAACAGAGCGAACTCCATCTCAAAAAAAAAAAAAAAAAACAACCTTTATTT
TTTCTGATTTTAAAAGTAATAACTAGTTTGTAGAAACATTAAAAGT
20

Sequence ID 492
ACCCTAACATAACTTAAAATTTGTTNGGAATTTGAAAGTACAGAATTTTCCTGTA
ATTGAGACTNTTTAACTTTTGTGGTTGGAGAAGGTATTCTATTTTTTGAAAATAT
CTGTAAGTTTTATCTAAATAGTAACTCTAAGTATTCTTCCCCTTTACTTACAGCC
25 ACCCTGGGAATCTGAGACTAGAGAAAATAAAGTTTGTCTCTTGTCTAAGGAGGGT
CTGGTTTAGAAATCTGATTTAGACATAGAAAAATTGCAAGAAGCTTGAGGTGATTG
GAAGATACGATTTTGTATCAAAGNATGTTTCTGTTTTATAGATTTTATTATCTA
CAACTCCTTATTAATATATTTAAGAAGTCATTAACCCACCATTGATTACTTGATAT
AAAAGGAGAANC GGTTAAAGGTGAAATANAATTTTTAATTTTTTTTTTTTAA
30 GTTTAGGATTTTTTTTTTAAATCTAAGAGTTTCTGTCATTTGGGGACAATCAGAA

Sequence ID 493
TGGAATCATAATTNGTTAACTGAAGCTNATAAGATGAGAGCATTCANAGAGAAAA
GAACGGAAAGATTGAATATCAGTTTCCCTTCTTTAAAAAAATTGTGGATATGTGAT
35 CTAGCTTCTTGAGCATCACAGTGAAGTCTGTTGCTCGTGGTAATTGATCGCTATGCT
GACAATCTTATCTCCACCTATGTCATTCAATTTTCTAAGAGGCAAAATCCTTAATC
AGGAGGAGAGTTTAGCTCTAGCTAAATTTCCCTTGTCCAGCATGCTCCTGCTCCCC

- 165 -

CAACTTGTGGAAACAGCTAAAGGATTGGACTAGGAGCANAAGTTTGGAATGGTTAA
AATGTAGCAACATGTGTTTCCTGAAACAAAATTCCACTATAATAAAAAAGCATTT
GAATGCTCCCTTGTAATTCTGTTGGAGCTTGTGTCCTTTTTTATGACACAACCATA
ATCAGTGATAGACAGTAGCATAAAGAAGCAAGAGCAAAGCAATTAAGTAATAATAG
5 CACTACAAAATGTGTGCTGTACTTACCAAACACGACATTTATGAATTATTANATA
GGAATAAGGGGATGGT

Sequence ID 494

10 GACCCAGCCATCTAAATAAGTTRTACATGTTGCGTATTTTTTTGTTAGGGACTTAT
CTTCCGAAGAGGAAAGGTTTATGAAACCTAAAGTAACAATGATAGCTTGAATCAA
AATGATAGCATTGTTGGCACAGCTGTGAATGATCATGTCCTCAAAGTGTGGAATTC
TTACACTGGACAACCTGCTTCATAACTTAATGGGACATGCTGATGAAGTATTTGTTC
TGGAGACACATCCCTTTGATTCCAGAATTATGTTATCTGCAGGACATGATGGCAGC
15 ATATTTATATGGGATATTACAAAAGGTACCAAGATGAAACATTATTTTAATATGGT
AAGTGAAGTGAGATGTACCTTGATACATGCTTGATAATTTGTTTAGAGTATTTGGG
TTATGCGGCTTACCCAGAAATTGATCTGCTTGTTTTGGCAGTTTGTTTTTACAAAT
CAACATATTCAAAGCCTGCTAAATATTAGACAGCTACATGTATATACGTACATACA
TGAA

20

Sequence ID 495

TTTC

Sequence ID 496

25 CTCGCTGGCGGGAGGCCACGGGCTTTCCACAGCGCGGGGGAACGGGAGGCTGCAGG
ATGGTCAAGCTGACGGCGGAGCTGATCGAGCAGGCGGCGCAGTACACCAACGCGGT
GCGCGACCGGGAGCTGGACCTCCGGGGGTGATCTGGACCCTCTGGCATCTCTCAA
TCGCTGACTTACCTAAGTATCCTAAGAAATCCGGTAACCAATAAGAAGCATTACAG
ATTGTATGTGATTTATAAAGTTCCGCAAGTCATAGTACTGGATTTCCAGAAAGTGA
30 AACTAAAATTTTAATCCAGGTGCTGGTTTGCCAACTGACAAAAAGAAAGGTGGGCC
ATCTCCAGGGGATGTAAAAGCAATCAAGAATGCCATAGCAAATGCTTNAACTCTGG
CTGAAGTGGANAGGCTGAANGGGTTGCTGCAGTCTGGTC

Sequence ID 497

35 GAAGACCTCACATCTGAGAGCTCATCTGCGTTGGCATTCTGGAGAACGCCCTTTTG
TTTGTAAGTGGATGTACTGTGGTAAAAGATTTACTCGAAGTGATGAATTACAGAGG
CACAGAAGAACACATACAGGTGAGAAGAAATTTGTTTGTCCAGAATGTTCAAACG

- 166 -

CTTTATGANAAGTGACCACCTTGCCAAACATATTTAAAACACACCAGAATAAAAAAG
GTATTCACCTCTANCAGTACAGTGCTGGCATCTGTGGAAGCTGCGCGAGATGATACT
TTGATTACTGCAGGAGGAACAACGCTTATCCTTGCAAATATTCAACAAGGTTCTGT
TTCAGGGATAGGAACTGTTAATACTTCCGCCACCAGCAATCAAGATATCCTTACCA
5 AACTGAAATACCTTTACAGCTTGTACAGTTTCTGGAAATGAGACAATGGGAGTA
AATATTACACAAATACTTATTCATTGNGGTTATTTTATACAGTAGTGAGAAGAAT
ATTGTTCTTAAGTTCTTAGATATCTTTTTTTGGATGTGCAAAAATTTTGGATTGA
CAGTAACTTGGGTATACATGACACTGAAATGCCTTACTTTGGATGA

10 Sequence ID 499

TGCCTGCGGGCCAGGACCTCGCCCAGCCCATGTTTCATCCAGTCAGCCAACCAGCCC
TCCGANGGGCAGGCCCCCAGGTGACCGGCGACTGAGGGCCTGAGCTGGCAAGGCC
AAGGACACCCAACAACAATTTTGGCATAACAGCCCCAGGCAATGGGCACAGCCTTCC
TCCCANAGGACCCGGCCGACCTCAGCGCCTCCTGCAGGCTAGGACACTGGTGCAC
15 TACACCCCATGCCTGGGGGCGGAGATTCTCCAGCAGAAAGATGCAATATTTTTTGT
TTCCTTTTTTTCCATTTTTTTCTCTAAGGAATCAATATTTCAATATGTTGAGTGTG
TGTCCAATGCTATGAAATTAAATATTTAAATAACATATTTATGGCATTTTTCTTGAA
GAGTGTGGTTGAAGAAATATTTCTCCTTTTGTTTTTCTTTTTTTTTGNTTGNTAC
TGCCACTTCTTTTTTAGGAGCAAATCTCCCCAGGGGTGTACGGNATTTCTTGACTCT
20 GGGAACAGCTGCTACCCCCAAGACTTGCCACGTTGTTCTGCCCTCAAATGGAATTA
AGTG

Sequence ID - 500

nt: 390

GGAATATGGTCAGGATCTTCTCCATACTGTCTTCAAGAATGGCAAGGTGACAAAAA
25 GCTATTCATTTGATGAAATAAGAAAAAATGCACAGCTGAATATTGAACTGGAAGCA
GCACATCATTAGGCTTTATGACTGGGTGTGTGTTGTGTGTATGTAATACATAATGT
TTATTGTACANATGTGTGGGTTTGTGTTTTATGATACATTACAGCCAAATTATTT
GTTGGTTNATGGACATACTGCCCTTTCATTTTTTTCTTTTCCAGTGTTTAGGTGAT
CTCAAATTAAGAAATGCATTTAACCATGTAAANATGANTGCTAAAGTCAGCTTTT
30 TAGGGCCCTTTGCCAATAGGTANTCATTCATCTGGTATTGATCTTTTCACAAA

Sequence ID 502

ACCCGCCATCTTCCAGTAATTCGCCAAAATGACGAACACAAAGGGAAAGAGGAGAG
GCACCCGATATATGTTCTCTAGGCCTTTTANAAAACATGGAGTTGTTCCCTTTGGCC
35 ACATATATGCGAATCTATAAGAAAGGTGATATTGTAGACATCAAGGGAATGGGTAC
TGTTCAAAAAGGAATGCCCCACAAGTGTTACCATGGCAAACTGGAAGAGTCTACA
ATGTTACCCAGCATGCTGTTGGCATTGTTGTAAACAAACAAGTTAAGGGCAAGATT

- 167 -

CTTGCCAAGAGAATTAATGTGCGTATTGAGCACATTAAGCACTCTAAGAGCCGAGA
TAGCTTCCTGAAACGTGTGAAGGAAAATGATCAGAAAAAGAAAGAAGCCAAAGAGA
AAGGTACCTGGGTTCAACTAAAGCGCCAGCCTGCTCCACCCAGAGAAGCACACTTT
GTGAGAACCAATGGGAAGGAGCCTGAGCTGCTGGAACCTATTCCTATGAATTCAT
5 GGCATAATAGGTGTTAAAAAAAAAAAAATAAAGGACCTCTGGG

Sequence ID - 503

nt: 109

ACATTTTCCGGNCCTTTTGCCATACACAGTTACAGAGATCAGTCAAATCCATACCA
CCACTGAGATCTCATTTATTGCCACAGATGCACAAAATAAATAACCCAAAATC

10

Sequence ID - 504

nt: 374

CCAGCAACGACCCATACCTCAGACCCGACGGCCCGGAGCGGAGCGCGCCCTGCCCT
GGCGCAGCCAGAGCCGCCGGGTGCCCGCTGCAGTTTCTTGGGACATAGGAGCGCAA
AGAAGCTACAGCCTGGACTTACCACCCTAACTGCGAGAGAAGCTAAACGTGTTT
15 ATTTTCCCTTAAATTATTTTTGTAAATGGTAGCTTTTCTACATCTTACTCCTGTTG
ATGCAGCTAAGGTACATTTGTAAAAAGAAAAAAACCAGACTTTTCANACAAACCC
TTTGTATTGTANATAAGAGGAAAAGACTGAGCATGCTCCTTTTATATTAATTT
TTACAGTATTTGTAAGAATAAAGCANCATTTGAAATCG

20

Sequence ID 505

GTACAGGAGGTAAATTGGATACCCCATCTAAGGGGATCTGTGAGACCAGGTAGTTA
TTTGGAATGAAAGAGTAAGATATTAAACCAGCCAGCATGTCAACAGGTGGGTGATA
GTCTTGTTCTCACAGACAACAGATGGCCATCATCTTAAAACAACATTTATGTTAAC
25 CAGCAGATAAGGGACTCCTGCATTGTCAGTGGACTTTGAGCCTGAGTTTTTCTACT
TGCATAGGTGAAAGTGGACTGCAATGCTAGTATAAATGCCGTATGATGACTAGTAC
CCCTTAGGGAGCTCCAGTTTGCCTTCCTGGGGAACCACAGACCCCAAGTGTAATTT
CCTGAGGACAGCCCGACTTCT

30

Sequence ID 506

GTTACTGTGAGCCTGTCAGTAGTGGGTACCAATCTTTTGTGACATATTGTCATGCT
GAGGTGNGACACCTGCTGCACTCATCTGATGTAAACCATCCCANAGCTGGCGAGA
GGATGGAGCTGGGTGGAACTGCTTTGCACTATCGTTTGCTTGGTGTGTTTGT
ACGCACAACCTTGCTTGTAAGTAACTGTCTTCTGTACTATTTAACTGTAAATGG
35 AATTTTGAAGTATTTGTTACAATAATAACTCTGAGATGTGTGAAAAAAAAAAAA
AAAAAAAAAAAAA

- 168 -

Sequence ID - 507

nt: 521

CTGCGGTGGAGCCGCCACCAAAATGCAGATTTTTCGTGAAAACCCCTTACGGGGAAGA
CCATCACCCCTCGAGGTTGAACCCCTCGGATACGATAGAAAATGTAAAGGCCAAGATC
CAGGATAAGGAAGGAATTCCTCCTGATCAGCAGAGACTGATCTTTGCTGGCAAGCA
5 GCTGGAAGATGGACGTACTTTGTCTGACTACAATATTCAAAGGAGTCTACTCTTC
ATCTTGTGTTGAGACTTCGTGGTGGTGCTAAGAAAAGGAAGAAGAAGTCTTACACC
ACTCCCAAGAAGAATAAGCACAGAGAAAGAAGGTTAAGCTGGCTGTCTGAAATA
TTATAAGGTGGATGAGAATGGCAAAATTAGTCGCCTTCGTGAGAGTGCCCTTCTG
ATGAATGTGGTGTCTGGGGTGTATTATGGCAAGTCACTTTGACAGACATTATTGTGGC
10 AAATGTTGTCTGACTTACTGTTTCAACAAACCAGAAGACAAGTAACTGTATGAGTT
AATAAAAGACATGAACT

Sequence ID 508

AAGCTCATGATTTTAAATGTATTTTCTAATAAACTATACTCCCATTTAAAAATCA
15 CCAATACCTTAATGTTTCAATTATATAAGCTAATTAAAAATAAAGGCTGGGCGTGG
TGGCTCACTTTTGAAGACCGAGGCAGGCAGATCACCTGAGGTCAGGAGTTCGAGAC
CAGCCTGCCCAACATGGAGAAACCCCATCTCTACTAAAAATACAAAATTAGCCAGG
CATGGTGGCACATGCCCATAATCCCAGCTACTGGGGAAGCTGAGGCAGGAGAATCA
CTTGAACCTGGGAGGCAGGGGCTGCAGTGAGCCGAGATCATGCCATTGCACTCCAG
20 TCTGGGCAACAATAGTGGAACCTCCATCTCAAAAATAAATAAAAAAATAAAATAAAA
ATAAAATTCAAACCTAAATAGATGCTCTACTTCAGGAGTGGGCAAATTAATCACC
TGCATCCTTTTTTTTGGGCTTTC

Sequence ID - 509

nt: 575

TTTTTTTCTAAATGGNGATTACTAATATATGTGGAGACTATTAATCTCTTTTCTGT
25 TGCCATTAGTTCATTTTTCCCCAAAAGCCAATACATGTTTCATTACAAAATGAATT
ATAAAATATAAGTTAAAAGAAAAACATAAAACCCTACAATCTTACCCACCCAGACA
ACTACTATTAATACCTTAGTATTAACATATACACATCATGTATATGTATAAATTTA
TCTTAAACAAAAATAAAATTATTCTTTACATATTGTTTTAAACCTATTTATCTGG
30 CCAGGTGCCGTGGCTCACGCTTGTAATCCCAGCACTTTGGGAGGCTGAGGCACGTG
GATCACCTGAGGTCAGGAATTCGAGACCAGCCCAGCCAACATGGTGAAACCCTGTC
TCTAATGGTTTAAATACCAAAAATTAGCTGGGCATGGTGGCACATGCCTGTAATA
TCAGCTAACATGGGAGGCTGAGGCAGGAGAATCACTTGAACCANGGAGGGGGAGGT
TGCAGTGAGCCGAAATCACACCACTTCACTGCAGCCTGGGCAACAAAGCAAGACTG
35 TCTCAAAAAGAAAAA

Sequence ID 510

- 169 -

CACTGTCATTCCCAGGAGGCTTTGGAGTCAGAACTGGATTCAAATTCTGACTNTAT
GTTGTGTGACTTGGGCCAATAGCTTCTTTNTGTGCCTCAGTTTCTTTAGCTGTAAA
TANACGGGTAGGTCACCCCTTACCCCATAGGTTATGGGGAAAGTTACAGAAAATGG
TCAGCTGGGCNCAGTGGCTCAAGCCTGTGGTCCAGCNCCTTGGGAGGCCAAGGTG
5 AGCAGATTGCTTGAGCCCAGGAGTTTGACACCAGTNTGGCAACGTGACGAAACCCT
ATCNCCTGTGAAAAATACAAAAAATTAGCCAGGCATGGTGGTGTGTGTCTGTGGTTC
CAGCTGCTTGAGAGTTTGAAGTGGGAGGATCACCTGAGCCCAGAAGGTCGAGGCTG
CAGTGAGCTGTGATCGCGTCACTGCACTCCAGCCTGGC-GACAGAGTGAGA-CCCC
T-TTTGAAAAAAAAAAAAAAAAAAT

10

Sequence ID 512

GTGAGCGGTGGTGGTTTATTCTTCCGTGGAGTTAAGGGCTCCGTGGACATCTCAGG
TCTTCAGGGTCTTCCATCTGGAATATATAAAGTTTCAAGAAACATGTCTCGAAGAT
ATGACTCCAGGACCACTATATTTTCTCCAGAAGGTCGCTTATACCAAGTTGAATAT
15 GCCATGGAAGCTATTGGACATGCAGGCACCTGTTTGGGAATTTTAGCAAATGATGG
TGTTTTGCTTGCAGCAGAGAGACNCAACATCCACAAGCTTCTTGATGAAGTCTTTT
TTTCTGAAAAAATTTATAAACTCAATGAGGACATGGCTTGCAGTGTGGCAGGCATA
ACTTCTGATGCTAATGTTCTGACTAATGAACTAAGGCTCATTGCTCAAAGGTATTT
ATTACAGTATCAGGAGCCAATACCTTGTGAGCAGTTGGTTACAGCGCTGTGTGATA
20 TCAAACAAGCTTATACACAATTTGGAGGAAAACGTCCCTTTGGTGTTCATTGCTG
TACATTGGCTGGGATAAGCACTATGGCTTTCAGCTCTATCAGAGTGACCCTAGTGG
AAATTCGGGGGATGGGAAGGCCACATGCATTGGAAATAATANCGCTGCAGCTGTGT
CAATGTTGAAACAAG

25

Sequence ID 513

TTTTTTTTTTTATAAACTCCAATCATTTCCAGAGCTACTTAGCTCAGCATCTTTTTT
TTCCACGCTCTTAAGTTGTGTTTATACATTTTTTGATACAGTTAGATTGTTTTTGTG
ACATTCTTCATTCTATCCTGGGATCCCCCAACCACCTAAGTGGATTTTTTGATAAT
TTGCATGCTTTAAGGATAACTCTTCATTCTGNAAAGGGCTATGGGTTTTGGCAAAT
30 GCAGAGTCATGTATCCAAGATTACAATATCGCACAGAAGAGTTTCATCACTATATA
AAACTCACCAGTCTTCCTCCTATTCAACCATCTCCATGCCTTCTTCCAGCCCTAA
CTCCTTAAACCACCTCATATCTTTACTATTGCTATAGTATTGCCTCTTCCACCATG
TCATATAAATGGAAACATACAGTATTAGTCTTCTCAAAGTAGTTTCTTTTACCTAA
CAACATGCATTTAAGATTCATAGTGTCTTTTAATGACTTGATAGATTATTTCTTTG
35 TAGCTGAATAATATTGCATCTTATAGATGTAACCGTTTGTATATCCATATTTTCTC
ACAGCCTATGACTTGNCTTTTGATTCTCTGAACAGGCCATTACAAAGCAGAAGTT
TTAATTTTTTATAAAGCTAATGNATCAACTT

- 170 -

Sequence ID 515

CCTGGATGACAGCATATCTGTTTATAGCTCAGTTTACTGAATACTTTAAGCCCACT
GTTGAAACCTGCT

5 Sequence ID - 518

nt: 502

GATGCATGTCCAGCATAGGCAGGATTGCTCGGTGGTGAGAAGGTTAGGTCCGGCTC
AGACTGAATAAGAAGAGATAAAATTTGCCTTAAAACTTACCTGGCAGTGGCTTTGC
TGCACGGTCTGAAACCACCTGTTCCCACCCTCTTGACCGAAATTTCTTGTGACAC
AGAGAAGGGCAAAGGTCTGAGCCCAGAGTTGACGGAGGGAGTATTTAGGGGTTTAC
10 TTCAGGGGCTCCCAAAGCGACAAGATCGTTAGGGAGAGAGGCCAGGGTGGGGACT
GGGAATTTAAGGAGAGCTGGGAACGGATCCCTTAGGTTTACAGGAAGCTTCTGTGCAA
GCTGCGAGGATGGCTTGGGCCGAAGGGTTGCTCTGCCCGCCGCGCTAGCTGTGAGC
TGAGCAAAGCCCTGGGCTCACAGCACCCCAAAGCCTGTGGCTTCAGTCCTGCGTC
TGCACCACACATTCAAAGGATCGTTTGTGTTTTGTTTTTAAAGAAAGGTGANAT

15

Sequence ID 519

CTGCGATNGAGTTTTGAGAGGAAGGANTAAAGTNCTCATCTCNGACGGTGAGAAAG
ATCATNACTAAGGAAACGCAGGGTTGGAAGCAGTGCTGANTGTCCAGTTGAGTTTC
ATGANCAAACATTTGCTGTGGGACCAGTTTTTCATGGNGGTTTGTCATTTTGTCCAG
20 CTGCCTGGAGCTGCTTGGTTGAAGGCACAGAATAATCAGGATTAATTGTTNAACTT
GTATGAATTTCTTTATTTTAAATAGGAATAATATCTGCCTTGGGAGCAAGTTGTA
AGAGTTAACTGAAAGCTTNAGGAAAACTTCCCTTGCTATTTAAGTAGGGCTTTA
CAAGTTACAATTCTATCACAGTTTTTAAGATTATAAAC

25 Sequence ID 521

GCGGCGCANCTGCGGATCCANAAGGNCATAAACGANCNAACTGCCCAANNCGTG
TGATATCACCTTCTNAGATCCAGACNACCTCCTCAACTTCAAGCTGGTCATCTGTC
CTGATGAGGGCTTCNACAAGAGTGGGAAGTTTGTCTCAAAAAA

30 Sequence ID - 523

nt: 585

GATTTACTGTGGGAATTTGCTCATGCAATTATGGAAACCTAGAAGTCCCATATAT
GCCATCTTCAAGCTGGAATCCCAGGAAAGCAGGTGGTGTAAATCTGAGATTGAAGT
CTTGAGAACCGGGGGAGTCAATGGTGTAACTCCCAATCTAGGGCTTAAGGCCCAAG
GACCAGGGCTGCTGGTGTGACAGATGCAATCCTGGAGTTCAAAGGATTGAGAACCA
35 GGAGCTCTGGTGTCTGAGGGCAGTAGAAGATGGATGTTCCAGCTCAAGAAGGGAAA
GTAAGAATCCGTCCTTCCCTCCACTTTTTTGTCTATTTCAGATGAGCCCTCAATGGA
CTGAACGATGCTCACCCACACTGTGAGGGCTGGTCTTCTTTATTCAATCCACTGAC

- 171 -

TTAAGTGCTGATCTCTTCTGGAAACACCTTCACAGACACACCCAGAAATAATGTTCT
TACCAGCCATGGGCCTGTTACTTAGCCCAGTCAAGTTGACACAGAAAATTAGCTAT
CACAAATCTGTGTGTGTATATACATATGTATTTGCATGTGTGTGTATATATGGNG
TATATATATTTCATGTGTGTGTATAT

5

Sequence ID 524

CTTTTGCCAGTAGGCCCCCTGAGTAGGTTCCCTCTATCTTTTGGCATGACCCCAGAA
GTCTTTGATAACTTCCTTGCTTTCTGATGTGACAAGACATCCAGGGCCAGATTGTC
CATATCCTGCCCCGGATGCACGATGCACTGTTTCTCCAAGAATCCCTGTGTCCTTT
10 GCTGATGATGCCATGATTTTAAGTTCTCTAATATAGTTTTATCTCTTTGTTTCAGA
TAATGCTTTTGTGTTCTCACATGTCCTGCTCTCTCTCTCTCATTTTTGGTGTT
GATCAGTCTTTCCATAAGATTGTTTATTTCACTAGTCCTTCATTCTTCTTTTTTCT
AAATTTACTCTTCTTGACTAGTATCCTGTCACTTCTGAGGACTCATATTTTTGCAA
CTTGAAAATTATTCTTATTTATTTAAGTATATGTTNCTGAAACTCTCATTAGACAC
15 ATTTTG

Sequence ID 525

[illegible]

25

Sequence ID - 526

nt: 516

CTTTTCATGGTCTCTTGTTCAATTAATCATCTAAAATCCAAGCNCAGAGAATTCAAT
TTTAGATGGTCTCCAGAGCAGAATTTGATGTATAATCTTAATTACAAATCATAGAT
AATTAATATTGNTTACAAAATCANAATACGATTAGAGGTAGGGATCCTGCACACAC
30 CCTATTTTCTCCCCAGTGTTCTGACCGAGAGACTAATTAATAATTCAAGGAACTT
ACAGTGAATGANAAACCCATGGTTTTGCTTAATTATCAGAACAGCTAGATCTGAGAA
CAGCTGTCTCCACATGGATAGACACTTATTCCACCCATTTGCAGGTAGAATAGCT
GGCAATAATAAGTCCTTCCCATTGGATATGTTGAAAGGTGCCTGCCATGGCATAGT
TGCCACAAGAGAGGAAGAAATGGACACAAATGTAGGCTGTTTTCAGGGGCANAGGGA
35 AGGTGGGAGGAACCAANTTGCTGGTTTTACACACCCTCTGGGGAACACCCATGC
ACCTATGANATG

- 172 -

Sequence ID 527

GACAAAAGCTGAGAGAATTTTTTTCTTGAATATTTGCACTAAAAGATAGGTTAAAA
TTCTTCAGGCTGAAGAGAGCATAACCAGGTGGAGATTTGGATCTACAAAAGGAAGG
AAGATTTGGAAATGGATTTGGCACCATTGACTCAATTTCCAGAACAAGAAAGCAGG
5 GACAGTTTTGGGAAGCTCAAGACACACTGCCCATGAGCAGCAATTTGGACCTCCTG
CTGCATCCACTGTGCATCAAACACACACTGTACAGACAAAGACTCCCAGGAAAAGA
AGTATAAACATGGACTAACACAGAGATGGGCAAACCTACAGCCTGTGACCCAGCCAC
CTGTTTATGTAGAATCCAAAGTAAGAATCTTTAACTTACACATAAACTT

10 Sequence 529; 660nt

GACAGCAGAGCACACAAGCTTNTAGGACAAGAGCCAGGAAGAAACCACCGGAAGGA
ACCATCTCACTGTGTGTAAACATGACTTCCAAGCTGGCCGTGGCTCTCTTGGCAGC
CTTCCTGATTTCTGCAGCTCTGTGTGAAGGTGCAGTTTTTGCCAAGGAGTGCTAAAG
AACTTAGATGTCAGTGCATAAAGACATACTCCAAACCTTTCCACCCCAAATTTATC
15 AAAGAACTGAGAGTGATTGAGAGTGGACCACACTGCGCCAACACAGAAATTATTGT
AAAGCTTTCTGATGGAAGANAGCTCTGTCTGGACCCCAAGGAAAACCTGGGTGCANA
GGGTGTGGANAAGTTTTTTGAAGAGGGCTGAGAATTCATAAAAAAATTCATTCTCT
GTGGTATCCAAGAATCAGTGAAGATGCCAGTGAACTTCAAGCAAATCTACTTCAA
CACTTCATGTATTGTGTGGGTCTGTTGTAGGGTTGCCAGATGCAATACAAGATTCC
20 TGGTTAAATTTGAATTTTCAGTAAACAATGAATAGTTTTTTCATTGTACCATGAAATA
TCCAGAACATACTTATATGTAAAGTATTATTTATTGTAATCTACAAAAACAACAA
ATAATTTTTAGATATAAGGATTTTCTGGATATTGCACGGGAGA

Sequence ID 529

25 GACAGCAGAGCACACAAGCTTNTAGGACAAGAGCCAGGAAGAAACCACCGGAAGGA
ACCATCTCACTGTGTGTAAACATGACTTCCAAGCTGGCCGTGGCTCTCTTGGCAGC
CTTCCTGATTTCTGCAGCTCTGTGTGAAGGTGCAGTTTTTGCCAAGGAGTGCTAAAG
AACTTAGATGTCAGTGCATAAAGACATACTCCAAACCTTTCCACCCCAAATTTATC
AAAGAACTGAGAGTGATTGAGAGTGGACCACACTGCGCCAACACAGAAATTATTGT
30 AAAGCTTTCTGATGGAAGANAGCTCTGTCTGGACCCCAAGGAAAACCTGGGTGCANA
GGGTGTGGANAAGTTTTTTGAAGAGGGCTGAGAATTCATAAAAAAATTCATTCTCT
GTGGTATCCAAGAATCAGTGAAGATGCCAGTGAACTTCAAGCAAATCTACTTCAA
CACTTCATGTATTGTGTGGGTCTGTTGTAGGGTTGCCAGATGCAATACAAGATTCC
TGGTTAAATTTGAATTTTCAGTAAACAATGAATAGTTTTTTCATTGT

35

Sequence ID - 530

nt: 660

GACAGCAGAGCACACAAGCTTNTAGGACAAGAGCCAGGAAGAAACCACCGGAAGGA

- 173 -

ACCATCTCACTGTGTGTAAACATGACTTCCAAGCTGGCCGTGGCTCTCTTGGCAGC
CTTCCTGATTTCTGCAGCTCTGTGTGAAGGTGCAGTTTGGCCAAGGAGTGCTAAAG
AACTTAGATGTCAGTGCATAAAGACATACTCCAAACCTTTCCACCCCCAAATTTATC
AAAGAACTGAGAGTGATTGAGAGTGGACACACTGCGCCAACACAGAAATTATTGT
5 AAAGCTTTCTGATGGAAGANAGCTCTGTCTGGACCCCAAGGAAAACCTGGGTGCANA
GGGTTGTGGANAAGTTTTTGAAGAGGGCTGAGAATTCATAAAAAAATTCATTCTCT
GTGGTATCCAAGAATCAGTGAAGATGCCAGTGAAACTTCAAGCAAATCTACTTCAA
CACTTCATGTATTGTGTGGGTCTGTTGTAGGGTTGCCAGATGCAATACAAGATTCC
TGGTTAAATTTGAATTTTCAGTAAACAATGAATAGTTTTTTCATTGTACCATGAAATA
10 TCCAGAACATACTTATATGTAAAGTATTATTTATTTGAATCTACAAAAACAACAA
ATAATTTTTAGATATAAGGATTTTCCTGGATATTGCACGGGAGA

Sequence ID 532

GAATTGTGATAGTTTCAGCTTGAATGTCTCTTAGAGGGTGGGCTTTTGTGATGAGG
15 GAGGGGAAACTTTTTTTTTTTTCTATAGACTTTTTTTCANATAACATCTTCTGAGTCA
TAACCAGCCTGGCAGTATGATGGCCTANATGCAGAGAAAACAGCTCCTTGGTGAAT
TGATAAGTAAAGGCAGAAAAGATTATATGTCATACCTCCATTTGGGGAATAAGCATA
ACCCTGAGATTCTTACTACTGATGAGAACATTATCTGCATATGCCAAAAAATTTTA
AGCAAATGAAAGCTACCAATTTAAAGTTACGGAATCTACCATTTTAAAGTTAATTG
20 CTTGTCAAGCTATAACCACAAAAATAATGAATTGATGAGAAATACAATGAAGAGGC
AATGTCCATCTCAAATACTGCTTTTACAAAAGCAGAATAAAAGCGAAAAGAAATG
AAAATGTTACACTACATTAATCCTGGAATAAAAGAAGCCGAAATAAATGAGAGATG
AGTTGGGATCAAGTGGGATTGANGANGCTGTGCTGTGT

25 Sequence ID 533

CTTGAACCTCGGAGGCAGAGGTTGCAGTGAGCCGAGATCACGCCACTGCACTCCAG
CCTCGGGGACAGAGCAAGACTCCATCTCAAACACACACACACACACACACACA
CACACACACACACAAAACAGATATACACTGAACACAGCACAAGTGGGACATAAGAG
ATTTAAAAGGGTTAGAGATGTAAAATGGATCTAGGAATGGAAACCATAAGGNGGGA
30 TTTATCAACTGGATTCTGCANAATGCTGTTAAGGCCAGATGTTAGCAGGTGTTACA
TAAAAAAGGGATACCATGAGCAAAAGTATTTGAACATGGGCAATGGTTGAAACAAG
TTTAAACAGATTATNTTTATTACCAAATCTCTCAAACCTTTAATATGCTATAAACA
TTGTGAAACAATAAAAAAACTTTCCAAAA

35 Sequence ID 534

GGGAAGGGAGCTATGAGTGTGTGTGTGTGTATGGACTCACTCCCAGGTTTACCTG
GCCACAGGTGCACCCCTTCCCACACCCTTTACATTCCCAGAGCCAAGGGAGTTTAA

- 174 -

GTTTGCAGTTACAGGCCAGTTCTCCAGCTCTCCATCTTANAGAGACAGGTCACCTT
GCAGGCCTGCTTGCAGGAAATGAATCCAGCAGCCAACTCGAATCCCCCTAGGGCTC
AGGCACTGAGGGCCTGGGGACAGTGGAGCATATGGGTGGGAGACAGATGGAGGGTA
CCCTATTTTACAACTGAGTCAGCCAAGCCACTGATGGGAATATACAGATTTAGGTGC
5 TAAACCGTTTATTTTCCACGGATGAGTCACAATCTGAAGAATCAAACCTCCATCCT
GAAAATCTATATGTTTCAAACCACTTGCCATCCTGTTAGATTGCCAGTTCCTGGG
ACCAGGCCTCANACTGTGAAAGTA

Sequence ID 560

10 GGCGGAGGTTGCAGTGAGCTGAGATGGCGCCATTGCTCTCCAGCCTGGGTGACAA
GAGCAAAACTCCGTCTCAAAAAAAAAAAAAAAAAAAAAAGCAATTTACTTAAAAAC
ATACAAACACAGAGACAAGTATTTTTTGAGAAACAAATACCTTTTTTCATTTTTTATA
CCAATGTAACAATAATCCATTAAACACACCTTTACTAACTGTTTTCTAGGAGTCTG
ATATGATGAGGAAATAGGTAAACCTTTAATAGCCAGTACTAAATTAGAGTGGCACA
15 ACTTTCCTGAGGAAAAAGATGGGTATTTTACTTTTCTGTTTTAGAAAAGTGGCTT
GACAACAGTATGCTTATGTCTTAGAGTTTGAAATTCAAGTTCTTGAACATTATTAA
TGGCTACAATCATTCATACCCACATTGGGCTGTATTCTTGATGAATCCAAAGTGAT
TTTCACCTCAACTCTGAATTTTATTCTCCTCTTTTGAATATAATAACAACCATCTCA
CTAGAGGAAGCATTTTCTGATTGGAGATTCATTATTGTTTTAGATAAT
20 GTTTTCATTTGCTTATGGGTATATAAAAAATTTTATCTTAAAAATATTTCTCTCA
TTTAGCTAGCAACATTGTTTTT

Sequence ID 561

CTCAGGGTGATCTCTGAACCCAAACTTGCCCCAAAGAAGGTTGCTCTGTCTCTCC
25 ACATCCCCATCTCCTCCCTAGGGCCTTGTTGGGGAGAGGCTCCTCCATCTTTCCCA
AGTCACACCATCGTTTCTACGTGGTCTGGACAAGAGCAAGAGCACACCTTGTCCC
CACCTTCTCCAGAGCAGCCAGAACCCACCTCAGGTGCCTTCCCCATCCGGTGCAGT
TAAGGCACTTCTGCCAGCACCATGGTATGAGCACTAGACTTGGAGTTAAGATTTGA
GAGCCCCCTCTGTCACTGTGGAAGCTTGAGCATGTTGCTTGATCTCTCTGAACCTT
30 GTGTTTCTCATCTGTGAAAGGTGATAATGTGGGGCTGCTGTGAGATTTAAAGGACA
TAATGCACCTACGGTCCAAGCACTGCCTGGAATACAGCANAAGCTCAACAGATACT
GGACAACCCATCCCCTTAGTAGAGGCCTAACCATGTGACCCAAGGCAAAAGTGCT
TAAAAAAA

35 Sequence ID - 562

nt: 580

ATTGCATGCAAGTTTGCTGAGCTGAAGGAAAAGATTGATCGCCGTTCTGGTAAAAA
GCTGGAAGATGGCCCTAAATTCTTGAAGTCTGGTGATGCTGCCATTGTTGATATGG

- 175 -

5 TTCCTGGCAAGCCCATGTGTGTTGAGAGCTTCTCAGACTATCCACCTTTGGGTGCG
TTTGCTGTTTCGTGATATGAGACAGACAGTTGCGGTGGGTGTCATCAAAGCACTGGA
CAAGAAGGCTGCTGGAGCTGGCAAGGTCACCAAGTCTGCCCAGAAAGCTCAGAAGG
CTAAATGAATATTATCCCTAATACCTGCCACCCCACTCTTAATCAGTGGTGGGAAGA
ACGGTCTCAGAACTGTTTGTTCATTTGGCCATTTAAGTTTAGTAGTAAAAGACTG
GTTAATGATAACAATGCATCGTAAACCTTCAGAAGGAAAGGAGAATGTTTTGTGG
ACCACTTTGGTTTTCTTTTTTGCCTGTGGCAGTTTTAAGTTATTAGTTTTTAAAAAT
CAGTACTTTTTAATGGAAACAACCTTGACCAAAAATTTGTCACAGAATTTTGAGACC
CATTA AAAAAGTTAAATGAG

10

Sequence ID 563

15 GCAACCTGCACAACCCCGCCCTGTTTCGAGGGCCGGAGCCCTGCCGTGTGGGAGCTG
GCCGAGGAGTATCTGGACATCGTGCGGGAGCACCCCTGCCCCCTGTCCTACGTCCG
GGCCACCTCTTCAAGCTGTGGCACCACACGCTGCAGGTGCACCAGGAGCTGCGAG
AGGAGCTGGCCAAGGTGAANACCTGGAGGGCATCGCTGCTGTGAGCCAGGAGCTG
AAGCTGCGGTGTCAGGAGGAGATATCCAGGCAGGAGGGAGCGAAGCCCACCGGCGA
CTTGCCCTTCCACTGGATCTGCCAGCCCTACATCCGGCCGGGGCCAGGGAGGGGA
GCAAGGAGAAGGCAGGTGCGCGCAGCAAGCGGGCCCTGGAGGAAGAGGAGGGTGGC
ACGGAGGTCCTGTCCAAGAACAAGCAAAAGAAGCAGCTGAGGAACCCCAACAAGAC
20 CTTGACCCCTCTCTGAACCAAAATATGCAAAGTGTGACCAGTGTGGAAACCCAAA
GGGCAACAGATGTGTGTTTCAGCCTGTGCCGCGGNTTG

Sequence ID - 564

nt: 671

25 GGAATAGAATTTTAAATAGTAATAACTGCTTGTTTTTTTTTGTGCAAGTACTTTTAT
ACATAAGATAAAACAAAACCTTACCACCAAACATACCAAAATGCACCTCTTTCATA
AGTGAGTTACTAAGATTTCTATACCTGGAATATCATGTATGTTTCATTTACTGGAT
GTTTACATTTTAGGAAGGAAAATAGTTTTTGTTTATTTTAAACAACCTGAATACTTATA
AACTGTTGTTCCCTGGAAGTTATTTATTCATAAAAAATTTGTTCTTTTGTGTCATGAA
TTTATAATTCCTAAATGAAGACCAGAAAGTACAAATTGCTGGGAGGAAGAATAGGC
30 TTTATTAATCAACTGATGTCTTGATTTTTCTAAATGGGAAGATTGCTTTATTTTAA
ACACTAATTATGGGAGCAGATTCTTAGCAAACCTCTTTGGAAAAGTTAATGTTATG
ATGTGCATTAGGCTGCCCCATCGTGTATATAAATGAAGCAGATTTGATTTTTGTAT
TCTTACGTTTCTCTGCTTTGTAGTTGTGGCTGTACTTAAAGAAATACAGAATTTCA
TATATTTAAAAATGTTTAAATGTGACCCACAGACATTGTAAATGGATTNAAAAC
35 AACATGAAAAATATTCAACCTAAAAGAATTCTTAACTTCACAAGTGTTTTACTTC

Sequence ID 565

- 176 -

CTTGGTTCCGCGTTCCCTGCACAAAATGCCCGGCGAAGCCACAGAAACCGTCCCTG
CTACAGAGCAGGAGTTGCCGCAGCCCCAGGCTGAGACAGGGTCTGGAACAGAATCT
GACAGTGATGAATCAGTACCAGAGCTTGAAGAACAGGATTCCACCCAGGCAACCAC
ACAACAAGCCCAGCTGGCGGCAGCAGCTGAAATCGATGAAGAACCAGTCAGTAAAG
5 CAAAACAGAGTCGGAGTGAAAAGAAGGCACGGAAGGCTATGTCCAAACTGGGTCTT
CGGCAGGTTACAGGAGTTACTAGAGTCACTATCCGGAAATCTAAGAATATCCTCTT
TGTCATCACAAAACCAGATGTCTACAAGAGCCCTGCTTCAGATACTTACATAGTTT
TTGGGGAAGCCAAGATCGAAGATTTATCCCAGCAAGCACAACCTAGCAGCTGCTGAG
AAATTCAAAGTTCAAGGTGAAGCTGTCTCAAACATTCAAGAAAACACACAGACTCC
10 AACTGTACAAGAGGAGAGTGAAAGAGGAAGAGGTCGATGAAACAGGTGTAGAAGTTA
AGGACATAGAATTTGGTCATTGTCCAAAGCAAATGTGTGCGAGAGCA

Sequence ID 566

GTCACCAAGAGCTTGTTGTCAGGTTTTCACTTGCTATTCGCAGAGATTTTTTTTAA
15 AGGCACTATTTGTAGTGTTAAAAGGGTGAATTTATCANAAGGCATAATAATCATAA
ATGTGTATATGCCTAATAATAGAACTTTAAAAGGCATGAAGCAACACTCAAAGGA
TTAAAGGGAGATCATCTCACCCCTTCTTACCAATTGATAGAATGATCTGATGAAA
ACAGTAAAATAACAACAGATCTGAACACTGTCAACCATCTTGACAAATACTTATGC
CTAGTGTTCCATTATTGGAACACTAAACATGTGGAATGATTTATATCCTACTGCTC
20 AAGGTCATCACCAAGGTCTAATTGTAAAATTTCAAAAAATTGCAACCTCAGGCATA
AATGGGTTAATCGACATTTATAGCACACACATGCAACATGTACCAGAGATTCCTTC
TTTTCTATGAACATGGTACTTCCACCAAGATAGACCACATTGTGAACTATAAAACA
AATCTAAAAACATTTGAAATGAAGGAAATTATATAAAATATGTTCTCTTGATCTCA
ATGAAATTAAATTAATACTATAT

25

Sequence ID 567

CTCATGGCGGCCAATGTAGGCCCAAACCTTCCTCAAGTCAAACCTCTCCAGGCCAC
CTTCTGCTTCCCGGTGGCATCAACAGGCCAGCTTTGACTTGAGAACAGCCTCTGC
AGGCCCTGCTCTTGCCCTCCAGGGGCTTTTTCCAGGCCAGCTCTTGCCCTCATGGC
30 AGCTGCCCCAGGCCAAATTTCTGCCTGCCTGCCAGCAGCCTCAACAGGCACAGCTC
CTCCCTCACAGTGGCCCATTTAGGCCCAAACCTCATGACTGTGAGGCCATTTCCAGGC
CTAGTGCCCTGCCTCGTGGCTGACTCTTGAAGCCCAAACCTTCCTCAAATCAGCCTT
TTGCCCAAACCTTCTGTCTACTGTGCGACTCTACAGGTCAGCCTCTGCCTCACAGTGG
ACCTCCAGACCCAGATGGTGTCTNCTGTGGCATCCTCAGGCGAAGCTCCTGCCTT
35 TCGGCAGCCTCTCCAGGCCAGCTCCTCCTGCTCCAGCCTTCTCTCCAGGCTCTGA
ACTTTCTCAGGTCTCCCTCTGTTGTCCAAGGCTGGAGTGTAGTAG

- 177 -

Sequence ID 568

TATATATGTAATGCCCTTAACCTAGTGTGTTGGCATGATCGTTGCTGAAAGGGAAGC
TTGTGGGTACAGTGTCCCCCTCAGAAGCCAAAGCCCAGGGAAGGTCGCCTGCCCAGG
TCAGGCTCCCAGCGAGTTTGTCTGGGGAGGGG¹CCATTCATACCTCCAGGTCAGGAC
5 AGAGGCTCGGGCTGAGGGAACCCTACACAGGTCTTGGGAAGCAGATCCTTCCTGCCT
AAGCCAGCAGGACAGCTCAACAGGAAGCATCTTCAGCCACGGGAGGAGAGGCAGC
ACCTTTTTTTGGAACCATAACAGAGCTAAGAATGGTGGTACAAGTAATAGATTCTGTA
CTGGCAACCCCACTTGGTGGAGCAAGTTCTAGGAAAAGGGGGCTGTCCTTGAGTCA
GCCATGGGGTTCAGCCACACAGTCACCGCAGCTGCTCTTTGGCACCGGGCGCTGGAA
10 AGACCTAGGATGACACAGCCTGGAAAGAGCTTGGGAAAAGCTCATCTTCACAGAA
CTACCTGCTATACCAGCCAGGGCAGGTGCTTATTCCCACAACAGCCCTCTGTTGTA
GGCGGCAGTGCCATCCTGAANGTGCCGTGGTACCTTCTGAANACCCAGCTGAGGGC
CTGTAATGGCACTTGCATGCCACATGGNACACCCTTTCCCGGTAA

15 Sequence ID 570

ACCGCGGCCGCGTNAANAAAAAAAAAAAAAAAAAGAATTCCACTTGATCAACTTAATTCC
TTNTCTTTATCTTCCCTCCCTCACTTCCCTTTTCTCCCACCCTCTTTTCCAAGCTG
TTTCGCTTTGCAATATATTACTGGTAATGAGTTGCAGGATAATGCAGTCATAACTT
GTTTTCTCCTAAGTATTTGAGTTCAAACTCCTGTATCTAAAGAAATACGGTTGGG
20 GTCATTAATAAAGAAAATCTTTCTATCTTAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AA
AAA

Sequence ID - 571

nt: 457

25 TTAGAGAGGTGAGGATCTGGTATTTCTGGAATAATCCCCCTGGGGAAGACGAA
GGGATGCTGCAGTTCCAAAAGAGAAGGACTCTTCAGAGTCATCTACCTGAGTCCC
AAAGCTCCCTGTCTGAAAGCCACAGACAATATGGTCCCAATGACTGACTGCACC
TTCTGTGCCTCAGCCGTTCTTGACATCAAGAATCTTCTGTTCCACATCCACACAGC
CAATACAATTAGTCAAACCACTGTTATTAACAGATGTAGCAACATGAGAAACGCTT
30 ATGTTACAGGTTACATGAGAGCAATCATGTAAGTCTATATGACTTCAGAAATGTTA
AAATAGACTAACCTCTAACAACAAATTAAAAGTGATTGTTTCAAGGTGATGCAATT
ATTGATGACCTATTTTATTTTCTATAATGATCATATATTACCTTTGTAATAAAAC
ATTTTTCCC

35 Sequence ID 572

CGTCTATTTGNGTTTCTTCTCACAATTGGTAAGTTCTCTGTATTGATTGATGGCTA

- 178 -

AGTTTGATTAGTGTTTTTCTCTAGTTGGTAATTATATTCTAGTATTTTATCATCTT
ATTGTTTACTCAACTNAAAGTGNCACAGAAGAGTTGCCAGGTTTCTCTTTGATATG
AGATCTCTNNTTGATTGGAATGCAAATCANAAGTGTCATGTTTTGAATAAAGGGA
CCAGATGACTTATAGGTATTCTTTCTCTAAATATAACTAAGGTAAGATTTTGT
5 TGAGGTACTTAATCTATATAAGTGGTAAAGAATTTACTTGAATTTCTCCAAATTCT
CATGTCTAAAGTCTGATTGATTAAATTCATTCTTGGTATTTTCATTTTGAAAAGAAT
GTAGCTTTAGCAAACCTCTTTGTATAAATGCAGTGGGATTAAGGTCATTTAAAAA
TTGTTATATCATTTGTATTTTTTAAATTTACCAGTTTTATTTTTCTTTTTACCCTT
AGCCCGGCCTCAGAAAGTGTGTTTGTGTCCATTTCTCCCAGCGCACCTCTGCATA
10 TCTCTACCCACTTGTCTATAATTCAGCATCCAGCAGAGGAAAACAAAGTGTTCGTA
CAGTTCCTCTACTAGCAGCATGCCTCCCCCAGGACAAGTGTA

Sequence ID 574

TTATTGCTGACATAAAAATGGTGCACATCGGCCAGGGCCAGGATGAATCAGCCAA
15 TCTGCACCATTTATACATGGAAGTGGAGAACATTGTGCCAATAATCATTTAATATA
TGCCAAATCTTACACGTCTACTCTAAACTGCTCTAATGAAGTTTCAGTGACCTTGA
GGGCTAAAGATTGTTCTTCTGGGTAAGAGCTCTTGGGCTGGTTTTTTCANAGCAGAG
TTCTTGTTGTGGGTAGACTGTGACTAGGTTTCACAGCCTTTGTGGAACATTCCGTAT
AACGGCATTGTGGAAGCAATAACTAGTTTCTATGAAAGAACCAGAGCTGGGAAGAT
20 GGCTGGGAAGCCAGGCCAAAGTGGGGGCAACAGCTTGCTTCTCTTTCTCTTCTCAC
CCTCAGTTTGTATGGGAAAATGGAGATGTCCTCTCCACTTTATCCACGATATCTA
AATG

Sequence ID - 575

nt: 209

25 CAGGATATCGAGACCATCCAGACAGCATGGTGAAGTCCGTCTCTACTGGAATAC
AAAAAGTTAGCCGTGTGTGGTGGCACGCGCTCTAATCCCAGCTATTCGGGAGGCT
TAGGCAGGAGAATTACTTGAACCCGGGAGGCGAAGGTTGCAGTGAGCTGAGATCGC
ACCATTGCACCTCCACCCTGG-CGACAGAGCAAGACTCCGTCT

30 Sequence ID - 576

nt: 541

CAGCCAACCCAGAAGGAGCCAGTCTACAACTATGCCTGATCCTCCTCATGGCAGGC
CACGAAGCATTGCTGCCATGTGTTGAATTATAAAACCCACATTGCTTTTTTGAACCC
TGTTGCGGGTAAAAATAACCAAATTATCAGTCCTTGGAACCCAGGCAATCAAGTG
AGTACAAGGTAAAGATAAGTATGGTTTAGAGGAGAAATTATGTTTCTGAACTGGTG
35 TCCTTTGATGGCAGCGTCAGCCTTGCTAAGTCAGAGTAGAGGGAGCAGTGACCTTA
ATAAGCTTTGGTGAGCATCATGTGCACGCGTGGGTGGGAGTCCCTTTTCACTGATGC
TTTTAAAGTGCTTTTGCAGACCCTGGAAGGGATCCTCCACACATATGAGGTGTGG

- 179 -

GACAGGTAGGCCAGAGAGGATTAGCCCTGCTTTCGAGACTAGAAATCTACAGTCCT
GAAGGAGCAGTAATTAATTGGTACACCTGTCAGGGCCAGCCCCAGGTCTCCTGGC
TTTTTCCAGGTTTTCTGTCTCACATGATTTTGCTTTT

5 Sequence ID 577

CTTTAATTTTTCAAGTGTTTAAAAACAATTTTATACTTAAGCCAGCCTTGAAGAT
AAGCACAAAATTTACCAGTTTACATTTAAAAAACAAACAAAAACGACAACAACCTC
AAGCACCCGCTCTGTGCATAGCACTATTCTAGGTGCAATAAAAGGGAATCTTAACC
TTAGAAATATGAGTTCACTTTCTGGAATTGTATTATCTCCTTTTCCAGAGAGTAAA
10 AATAAATAAAATCACCATTGTTTACTACAGATCTGCCCCAAACCACATCTGGTTCA
CAGAAAGGCTAATTTCTGCCAAATTAAAGATGTAATGAACTCAGTTCCTGCTTTCC
CAAAAACACGAAAGCAGAATTCCTTTTCACTGAAAAAATAAACAGTTTTTCCATGC
AAGGGCAGTTTGCTTCTAATAAGTATTTTTTAAAAAATTTTTTTTTCTCTAGCTT
TTCTTTTAAATTTTCTTCTCTAATATTGCCTTTTCTTGTACAAGGCAGACCAGGTA
15 TCTTTTTATGCTGTTTTTCTTTTACTAAGAAAAGTATTGCATCTTGAAGACAAACC
ATTTCCAGAGTAGTGATAAAAAATAACACTAAAAAACTTTAAAGGTGAGTCACT
TCATCACCTTGATGAAGTAAAAAA

Sequence ID 578

20 GGAAAAAATATTTCCACTTAGATATTTTACATGGTTTTGTTTAAAATTACCATTAC
TTGTTTTTTTAAAAACACATGACCACATATGTATATGTATATCTACCTAAACATTGT
ATCATGGTTTTCAGTATGTTATTTCATGTATTACTGGGAGATGCTACCAAGAAACCAA
CCCAAAGAAAATCTGGAAAATACATTTCTATTTATAGAATAAATGTTTCATTTAT
ATAAAAGCAAAAGAACTTAGAGTTCTAATAAATGGGATGTCTAATAAATTATGAAG
25 TTAGTGATTTGAATATATTATATTTTTTATAACTTCCTTGCCAAAGTCCTGATTTAG
TACATTAGAGAACCTGTGTTTCTCTCTCTCTACCATTCATCTCTCTTCCATACA
GTCATTTGGGCTTTTTTACTCAAAGAGAATCAAGAAATAATAAGGTATAACAAGCTT
GGCAAAGTGTTGGCTTTTTTAAAAAAAATTTTTTTAATCTCTAGCAGTTTGGTAAAT
TTAGCAGCATCATTTATTTGGGATTCTTTTATCTGATTTCAACAGTGAAAAACATC
30 CCTATGATAAAGCCTAATGACCCATTTCCAAAAGATGGAATTGCCCTTCTAGAAA
ATATGACGGAGAAAAGT

Sequence ID - 579

nt: 502

35 CGAATAGCCAAGTGGTCTGACAAGATCGAGAGTAATGAGGCCCATACTTTAGTACA
GTCTTGAATGGCCAGATGGTGCTGGGCATACCCCAACCAGAGATATGTAAGTCTTT
ATGTTGTCAAAATTTCCCAGAAACATGAATTTCCCACTAAGATTCATTAAGGAAAA
CTAGAATGAAAACAAAACGTTCTTGTATAATATTCATTANAAAGAAATGAAGAA

- 180 -

GGCCGGGCATGGTGGCTCACGCCTGTAATCCCAGCACTTTGAGAGGCCAAGGTAGG
CAGATCATGAGGTCAGGAGTTTGAGACCAGCCTGGCCAACATAGTGAAATCCCGTC
TCTACCAAAAATACAAAAAATTAGCCGGGCATGGTGGCACACACCTGTCATCCCA
GCTACTCAGGAGGCTGAGGCAGGAGAATTGCTTGAACCTGGGAGGTGGAGGTTGCA
5 GTGAGCTGAGATTGCACCACTGTACTACAGCCTAGGTGACAGTGCAAGACTCTG

Sequence ID - 580

nt: 316

CCTATGCCAAACTAAAGAAAGCTTGCCTGGCCTACAGGCCTAAAGGTTCAAATGNG
GATTAAAAAACACAGTAGTCACATAAAATGTCTGCTGGCTGGCTGGAATTCCATC
10 ACCTACAATTTACCTGCTTTCAAAACTGTGTTCAACATTGAGAAAACAGAAAACC
ACTTATCTTGAGCTTAATATGGGCTTCTTTTTCCTTAACTGTAGAACACTTACTGA
AATATCAAATCAATGGTTAGGATATGTATCCTAGGCAGGCCTAAACCATTAACT
TGGTTTAAGCAACTTTGTATAATTNACCTCCTAAAT

15 Sequence ID 581

CTTCATGAGTGCCCGGTTGCCCAAGTCAAAAACCTGGGAGTGATATAAACTCCCCA
CACATCCAGTCAGTCACTCATCAACTCTATTGATTCTG-CTGCTAAATATATCTCA
ATTGTATTAACTTAAACATATGCATAATACATCTTCTTCTTCACTGCATTTTTGTG
GGCTGCACTTACCTTTCAGGTAACAACAACACTGGCCCCCTCTTGCCCTTCTAGTCA
20 GAAGTGCCAAAATGATGAGAGCTAGCCATGACAAACCCACAGCCAACATTACACTG
AATGTGCAAACTGGAAGGGCATCCAAACAGAGGAGG

Sequence ID 582

TAGAATTCTCGCCTGCCTTGGCTTCTCCCTCTAGTTGTTTCCTTCTCTGTCTTCTGT
25 GGGCTTCTTATTGTCTGCTCACTCCTTCTTCAGTGTCTCTCATGGGCTTCCCTCC
CTTCTCAGCTGATGCCATCACCTGGGGAATCACAGTTACTCAGCAGCACTGGGGCC
TCTCTATCTCTATGCTGGTCATGCCTATGTGTGAGCTGCAGACCCAGTGGAATTTT
CATTTGTGCATCCCATGCCCAGCCCACCCTCCACCAGCCTCGAATGCAGCTGTTCA
GCCCTACCCAGTCCTCAGAAAAGTTCTCTCCCTGGATCCTCTTTTTCCTTCATG
30 AGTGCCCGGTTGCCCAAGTCAAAAACCTGGGAGTGATATAAACTCCCCACACATCC
AGTCAGTCACTCATCAACTCTATTGATTCTGTCTGCTAAATATATCTCAATTGTAT
TAACTTAAACATATGCATAATACATCTTCTTCTTCACTGCATTTTTGTGGGCTGCA
CTTACCTTTCAGGTAACAACAACACTGGCCCCCTCTTGCCCTTCTAGTCAGAAGTGC
CAAAATGATGAGAGCTAGCCATGACAAACCCACAGCCAACATTACACTGAATGTGC
35 AAAACTGGAAGGGCATCCAAACAGAGGA

Sequence ID - 583

nt: 631

- 181 -

CTGAGGTGGGAGGATTCCACTCTCACCCATTTCTTCTTTTCATTTTCAGTTTCTCCA
GTTAGTAACTGAAGATGTTCTTTGAGTAATTAAGTGAGTGAGAAAATTTTAAAGTG
AGAAATCTATAAAAAGAACCATGTTAACATAAATATTTTCAGTCCTTACAAGTTGGT
ATTGACTTTTCTCATTGGTAATCTGACTGATTTAATACTGCTCATTCCAATATCTG
5 GTGATGTAATTCTGGTTATGAATCCTTGTATTAATAACACCTCCTGGGAGGTTTTT
TTTCCCAACATTACATTCAGAATATTAGAGCTGAAAATACCTTTTTTAAAGGTTAT
CAGGAGGAGGGAGCTTATGTTTAATGTGGTGGATAAACTTAACTGCTGGTTAATA
CAATTGTTATTTCAGGTGAAATTCCTTAACTTTTCACGTGCAAAGTTTTGTATGTA
TACAGACATTTGGGGAAAAGTTTTATCATCCCTAAAACCGTTACTGTCCAGAAAA
10 TGATAAGAATCCCTGGGTTCCAAATCCTTCATAAGGTATTTATTCAATTTATTTATT
CAACACATTTACTCAATGCCTCCGCTCTGCTGCAACTACACTGACATTCTGCTTCT
AATCTAACCGAAAAT

Sequence ID 585

15 TTTCAAATTGTACAATAACACAAACAACTTTGTAAAGGCCATGTTTTATTTGCTGA
TTAATGGACAAAAGGCAATGTAATTTATTTTCAAGTATTTTCTTGAAAGTCTGTGC
TCATAAAAATCATGAAAAGTTGGAAAGACTGTTAAATCACTGAAACTTCAAATATA
TCTTACACAATCTTGTTTGTACAAAAATACAAGTTAAATATAAACATAAAGCAATC
ATGGTAATTTTATGCAAATCTGTTTTATGTGATCATCAGTTATATATAAAAGTTTC
20 TCAGTTCTGTTATTTGTGAAAAGATCAATACCAGATTGAATGACTACCTATTGGCA
AAGGGCCCTAAAAGCTTACTTTAGCACTCATCTTTTACATGGTTAAATGCATTTTC
CTAATTTGAGATCACCTAAACACTGGAAAAGAAAAAAATGAAAGGGCAGTATGTC
CATAAACCAACAAATAATTTGGCTGTAATGTATCATAAAACACAAACCCACACAT
CTGTACAATAAACATTATGTATTACATACACACACACACCCAGTCATAAAGCC
25 TAATGATGTGCTGCTTCCAGTTCAATATTCAGCTGTGCATTTTTTCTTATTTTCATC
AAATGAATAGCTTTTTTGTACC

Sequence ID 586

GTAAACTGTTCTCTCCGAGGGAAAAAATGGAAGTTATCCTCACAGTTCACTGCCGT
30 GGTATTTCTTCTGTCCCATGCTTTGCATGACTGCCATGGTACAGCCTTGTTTTCAA
CTGTTCACTGTGATCTGTGGGTCTTTGAGTTTCAGTGAGTTTGCTGAAATGTCGAA
GAAGTAGTTCCAACTTCAATGTTCAATGAAATTTTTGTTCAAGTTTGAAATGGAG
AGAGCAGCTTTAAAAGGTACTAAGCCTTTTACAAATTGGTGAGTACTGGCACATGA
GAT

35

Sequence ID 587

TTTTTTTTTTTCCCTAAAAGGTAACCCCTAAACACAGCTAAAACATGCCATCAGC

- 182 -

TGACTCCAAGGNACACACAGTCCTGTATCTGGA ACTACTGAGTGGCAGGCATCTTT
CTCTGCCTCTGACAGTGGAGTCCCCATCACTGCAGAGCATAGCCAAAGGAGTCAAA
GGTCTCAGCGGGTCACTGCCTTATCAACCCTCACCAGTCCCTTATGTTTTTTAATA
TTTTATAATCTTGACATGACACCAAGATGCTTTAATAAAAAAGCACCTCTAACTCG
5 GTCTTGATTCACTTACCTTGAGCCTGGGACTTCTCTAGGCTCCTGAGGCAAAAAC
AGGTAGAGGGGAGATGGTGGAACATAAAACACAATTTTGCTTGGCACCCACCTTGG
CGTCTGTCCCATGACCAGGTCTTTC AATTCGATGATTTTGTCA TTGATGGAGGAG
CGATATCGTTTCTCAATGATATTATGGGTGTCCGCCTTTCTCCTTCTTTGGGGGG
CTCAAGCTGCTTGACTCCCCCAGGTACCTGCTTAATGGGGCACTTTCTCTTGCCCC
10 ATCATTACAGGCATTGTGGTCAGAATGGTCCCAC TGTGCCCACCAGGGTCTA

Sequence ID 588

CTAGTCTTTTTCATAGTCTGCATAGAGTCTGGCCATTACCATCAGTTTTTAAAGATGT
CCATATTGTGGCCGGGCGCGGTGGCTCACGCCTGGTAGTCCCAGCACTTTGGGAGG
15 CTGAGGCAGGTGGATCATGAGGT CAGGAGATCGAGACCATCCTGGCTAACACGGTG
AAACCCGTCTCTACTAAAAAAATATTAAAAAATTGGCCAGGCCTGGTGGTGGGCG
CCTGTGGTCCCGGTGCTTGGGAGGCTGAGGCAGGANAATGGTGTGAACCCGGAAG
TCGGAGGTTGCAGTGAGCCAAGATTGCACCTGGGCAACACAGCGAGACTCCGTCTC
AAAAAAAAAAAAAA

20

Sequence ID 589

CAATTATTTATTACCTTTCCATTTGTTTCGCCTGATGATGTGACAATGCATGGTCTT
TGTGCATGCTGCTAGACACTTTTCTTTCCAGCCGAAAAGTCTATTATGTAATTTT
TACATTCATAATTTTAATGTGGATGATCAGGATTAAATCAAGATATATATCTGGAA
25 CCTCTTATAAATGGAGCACTTAGAAATTTGTTGTTCTGCACTTAACCTAGAGAGAG
AAAAAATGCTTTTCTTTGTGAAAAATCTGAATTCCTGTCCTGACCTTCTGTGATGT
GGAAACCCTAGGCTCTGAGACACACTCTCTGGTGTCTGAGACAGAACCAAAGCAAT
AACGTTGTGATGCCCACAGGCCTGGAGCCAGCTAGCGACCTTGTGCCGCCAGCTG
TCCATGGCCCGTGCAGAGCAGAGGACAGTGAGTGTCTGCACTGAGAACCTTAAACC
30 ACAGTTGAACATACCACACCTGTTTGTCTTAAGCTATAGTGTA AAAACAAAGTTT
GGGCTCTGAAAATTTAACTGAAAAAGATTTCCCTTGTT

Sequence ID 590

GTGGCAGCAGGCGCAGCCCAGCCTCGAAATGCAGAACGACGCCGGCGAGTTCGTGG
35 ACCTGTACGTGCCGCGGAAATGCTCCGCTAGCAATCGCATCATCGGTGCCAAGGAC
CACGCATCCATCCAGATGAACGTGGCCGAGGTTGACAAGGTCACAGGCAGGTTTAA
TGGCCAGTTTAAACTTATGCTATCTGCGGGGCCATTTCGTAGGATGGGTGAGTCAG

- 183 -

ATGATTCCATTCTCCGATTGGCCAAGGCCGATGGCATCGTCTCAAAGAACTTTTGA
CTGGAGAGAATCACAGATGTGGAATATTTGTCATAAATAAATAATGAAAACCTAAA

Sequence ID 591

5 CAGCAGCAGAAATGTTTGCAAGATAGGCCAAAATGAGTACAAAAGGTCTGTCTTCC
ATCAGACCCAGTGATGCTGCGACTCACACGCTTCAATTCAAGACCTGACCGCTAGT
AGGGAGGTTTATTTCANATCGCTGGCAGCCTCGGCTGAGCAGATGCACAGAGGGGAT
CACTGTGCAGTGGGACCACCTCACTGGCCTTCTGCAGCAGGGTTCTGGGATGTTT
TCAGTGGTCAAATACTCTGTTTtagagcaagggtcagaaaacagaaatactgtca
10 TGGAGGTGCTGAACACAGGGAAGGTCTGGTACATATTGGAAATTATGAGCAGAAACA
AATACTCAACTAAATGCACAAAGTATAAAGTGTAGCCATGT

Sequence ID 592

TACTCAATGAAAAACCATGATAATTCTTTGTATATAAAATAAAcatttgAAAAAA
15 AAAAAA

Sequence ID - 593

nt: 565

CAGGATCAAGGTGAAAAGGAGAACCCCATGCGGGAACCTTCGCATCCGCAAACCTCTG
TCTCAACATCTGTGTTGGGGAGAGTGGAGACAGACTGACGCGAGCAGCCAAGGTGT
20 TGGAGCAGCTCACAGGGCAGACCCCTGTGTTTTCCAAAGCTAGATACACTGTCAGA
TCCTTTGGCATCCGGAGAAATGAAAAGATTGCTGTCCACTGCACAGTTCGAGGGGC
CAAGGCAGAAGAAATCTTGGAGAAGGGTCTAAAGGTGCGGGAGTATGAGTTAAGAA
AAAACAACCTTCTCAGATACTGGAAACTTTGGTTTTTGGGATCCAGGAACACATCGAT
CTGGGTATCAAATATGACCCAAGCATTGGTATCTACGGCCTGGACTTCTATGTGGT
25 GCTGGGTAGGCCAGGTTTCAGCATCGCAGACAAGAAGCGCAGGACAGGCTGCATTG
GGGCCAAACACAGAATCAGCAAAGAGGAGGCCATGCGCTGGTTCCAGCAGAAGTAT
GATGGGATCATCCTTCCTGGCAAATAAATTCCTGTTTCTATCCAAAAGAGCAATAA
AAAGT

30 Sequence ID 594

CAGAAGAGTAAGCAAATCTCAAAGCAGCGAAAGGGAAGAACTAAAAAAGGTAGAG
CAGAAATAAGAGAAAATAGAGAAGAGAACAATTGAGAAAAATAATTGAAACCAAAA
GGTGGTTCTTTGAAAAGCCTAACAAAATGGACACATCTTTAGTTAGAGTGACCAAG
AAAAAAGGGCAGTGA CTGACTCAGATTACTTCATTCAAGAGTGAAAGAGGGCACATCACT
35 ACCAATTTACAGAAATAAAAAGGATTATGAGGAAATACTACAGATAATTGATGACA
TTAACTTAGAAGAATATATTTCAAGAAAGACACAACTACTGAAACCGACTCAAGA
AGAAACAGAAAATCTGAACAGACCTATAAAAAATAGAGATTTAATTGATATTCAGA

- 184 -

AAGTTTCCCAAAAAGAAAAGCACTGGCCAAGATGACTTCACTGGTGAATTCTATCA
AGTGTCAAAGATGAATTACTGACATTCAATTCACACTCCTTTAAGAAATAGAAGAGG
GGACATCACTTTTCAAAGCATCGACATTCTAATCATTAGTCCCTTGGTTTCCTGCT
CCCAAAGCCAGGTGATGTATCACAAAAAACCCCTACAGACCCACTGGGCACAATG
5 GCTTTATGCCTAT

Sequence ID - 595

nt: 98

CTTTGCTCGAATNGTCAGATAAGGATTCTGTGAANGGAGATGAGATTTCCATCCAT
GCTGACTTTGANAATACATGTTCCCGAATTGGGGNCCCCAAA

10

Sequence ID 596

CTCAAGTGTTCCCTCAGCTTAGGCTTTGTTTAAATGATCCCACCCAGGGGCGATGG
TAGGGAACAACAGGGTCACTAAACTATTTGGCTGGCTACAATCTGGGAAATGGTA
AGACAGGGAAAGGCCATGTTGTTTCATTCCCTTGTGCAGATCTAGGGAGAACCGCAG
15 AGAGAACAGTTAGCATTTCTTGTTCAATGAATTATCCTATTAAGAACACTGGATGT

Sequence ID 597

CGGNCGCGGTCGACGCTACTCCTACCTATCTCCCCTTTTATACTAATAATCTTATA
AAAAAAAAAAAAAAAAAAAAAAAAAAAA

20

Sequence ID - 598

nt: 362

GGCATGTGCCTGTAGTCCTAGTTGCTGAGGTAAGAGGATTGCTTGAGCCCAAGAGT
TCAAGGCTGCAACAAGCTTTGATTGCGCCACTGCACTCCANCCTTGCGGACAGACT
AAAACGCTGTCTCAAAAAAAAAAACAAAAACGACNAAAAAAAAACAAAACAGAAAAA
25 ATTAACCTTAGGCAATGACAGTCCCTGGCAAATGCTGGGAGGGAGGCAACANTGGTC
AAGGAAGGTAACCCTGAANCAGGACTTGTAAGCAAATAANATTGGGAGGCCAAGG
TGGGTGGATCACNAGGTCAGGAGTTCGAGACCAACCTGGCCAACATAGTGAAACCC
CGTCTTTCTAAAAATACAAAAAATT

30

Sequence ID 599

GACAAAAGAACCATTTGGATACATAGGTATGGTCTGAGCTATGATATCAATTGGCT
TCCTAGGGTTTATCGTGTGAGCACACCATATATTTACAGTAGGAATAGACGTAGAC
ACACGAGCATATTTACCTCCGCTACCATAATCATCGCTATCCCCACCGGCGTCAA
AGTATTTAGCTGACTCGCCACACTCCACGGAAGCAATATGAAATGATCTGCTGCAG
35 TGCTCTGAGCCCTAGGATTATCTTTCTTTTCACCGTAGGTGGCCTGACTGGCATT
GTATTAGCAAACCTCATCACTAGACATCGTACTACACGACACGTACTACGTTGTAGC
TCACTTCCACTATGTCCTATCAATAGGAGCTGTATTTGCCATCATAGGAGGCTTCA

- 185 -

TTCACTGATTTCCCCTATTCTCAGGCTACACCCTAGACCAAACCTACGCCAAAATC
CATTTCACTATCATATTCATCGGCGTAAATCTAACTTTCTTCCCACAACACTTTCT
CGGCCTGTCCGGAATGCCCCGACGTTACTCGGACTACCCCGATGCATACACCACAT
GAAACATCCTATCATCTGGAG'

5

Sequence ID - 600

nt: 595

TTCAAATTCTTGNTAANAGTCTTTGTTCTGAATTTTACTTTGTCTGTTATTCCTAT
AGCCTTTCCAATTTTCTTTGCTTGGATTTTACGTGATAAGTTTTTCCCCCATTT
TACTTTTANCAACTCTATATTTTTTAGTTGAGGTTGGGTTTCTTGTAACAGCATA
10 TAATTTGGGTTTTTTAATCCAATCTGAAAATTAATGTCCTTAATTTTGTGTTTATA
CCATTTACACATAATGTACTCATATATAAGGTTTAACTGAAACCTACTATCTTGCT
AGTTGTGCTCTACTTGAATTTTTTTTTTAGTATTCTGTTTTAATTGACCAACATTTG
ACTGTATCTCTTTGTGTAATTCTTTTACAGGTTGCTGTAGGCATGACAATATATAC
ACTTAACTTTTCTCAGTACACTGAGAGTTGAAATTGTAGTACTTCGAGGAAAACAT
15 AGAAACTTGAATGATATCGGTTACATTTTACCACCTCCATATGTTGCAATTATT
AAATGTATTAGATCTGCCTACCTCGAAAACCCATCAGTCTTTTAACTTTGCTCTCA
ATGGTGATTCATATTTTTTAAAAAACTTGAGGCAA

20 Sequence ID - 601

nt: 522

TCGACCGGGTTTGGAGCAGTGCCTTGTTTGCTGTGCAGCGGATACTCTACAGGTAC
ATTTCCTTTTTGGAAACCAAAGGGAGGGATTGACAATATTGATGGTAGATCTTTT
TTCTTTAGCAAGAATTAAGGATTTTGGTGGGTGGGGGGAGGCTTCTGTGGGGACCA
AGACAATGTACTGTCAGTCAGGATTTAAGTCGAACTACCTCATCCCTTGCCCCAGA
GAACAGTTGATCGTGTTTTAAACCAAAGGTGCGGAATGGAGAGAGGGAGGCGGTG
25 CATTGCAGCTTCCGATAGAGCTTTTTATTTTTTGGATATCAGGAACCAATTTTGAAG
ATTTCTTAAGAAAGTCATTTACATCAGGGACATGAAGAGCAAAGTAGGTATTTTTG
GTCAGTACTTGAATTTGATAGGCTTTATGCAAACAACTCTCCCTCTGCTGGAGTCT
GGCAAGTTTGCTTTTCACTGGACGCTAATTCAAGTGCCATACAAAACCTAAAATAAN
AGTTTTACTTATAACACA

30

Sequence ID 602

CAGAAATCGCAATTGAAGACCAGATTTGTCAAGGTTTGAAACTGACATTTGATACT
ACCTTCTCACCAAACACAGGAAAGAAAAGTGGTAAAATCAAGTCTTCTTACAAGAG
GGAGTGTATAAACCTTGGTTGTGATGTTGACTTTGATTTTGCTGGACCTGCAATCC
35 ATGGTTCAGCTGTCTTTGGTTATGAGGGCTGGCTTGCTGGCTACCAGATGACCTTT
GACAGTGCCAAATCAAAGCTGACAAGGAATAACTTTGCAGTGGGCTACAGGACTGG
GGACTTCCAGCTACACACTAATGTCAATGATGGGACAGAATTTGGAGGATCAATTT

- 186 -

ATCAGAAAGTTTGTGAAGATCTTGACACTTCAGTAAACCTTGCTTGGACATCAGGT
ACCAACTGCACTCGTTTTGGCATTGCAGCTAAATATCAGTTGGATCCCACTGCTTC
CATTTCTGCAAAAGTCAACAACCTCTAGCTTAATTGGAGTAGGCTATACTCAGACTC
TGAGGCCTGGTGTGAAGCTTACACTCTCTGCTCTGGTAGATGGGAAGAGCATTAAAT
5 GCTGGAGGCCACAAGGTTGGGCTCG

Sequence ID - 603 nt: 624
GACACACGAGCATATTTACCTCCGCTACCATAATCATCGCTATCCCCACCGGCGT
CAAAGTATTTAGCTGACTCGCCACACTCCACGGAAGCAATATGAAATGATCTGCTG
10 CAGTGCTCTGAGCCCTAGGATTCATCTTTCTTTTACCGTAGGTGGCCTGACTGGC
ATTGTATTAGCAAACCTCATCACTAGACATCGTACTACACGACACGTACTACGTTGT
AGCCCACTTCCACTATGTCCTATCAATAGGAGCTGTATTTGCCATCATAGGAGGCT
TCATTCACTGATTTCCCTATTCTCAGGCTACACCCTAGACCAAACCTACGCCAAA
ATCCATTTCACTATCATATTCATCGGCGTAAATCTAACTTTCTTCCACAACACTT
15 TCTCGGCCTATCCGGAATGCCCGACGTTACTCGGACTACCCCGATGCATACACCA
CATGAAACATCCTATCATCTGTAGGCTCATTCAATTTCTTAACAGCAGTAATATTA
ATAATTTTCATGATTTGAGAAGCCTTCGCTTCGAAGCGAAAAGTCCTAATAGTAGA
AGAACCCTCCATAAACCTGGAGTGACTATATGGATGCCCCCACCTACCACACAT
TCGAAGAA

20
Sequence ID - 605 nt: 338
ACCTGAGGCCTCGGTGGGGCCAGTGCGACGCTGGCTTAAGGAGCTGGAGGGGTTC
TAATACACATTTAATTCACTTTCTCTTCCCTAAGAGGCTGCCGGAGTTGGGGCCTC
CTCCAGCAGAGACCCTCGGACCCCTGCAGGGCCTGGACTTGGGGTGAACAGGGCTT
25 CAGTCAGCGCAAGTATTCCATTTGCATTTGGTAATTTTTTCATGCCACCTATTTATG
AATATATAAATCTTTATACCAAATCTATTTTTTAAACATGGAAAAGTTGCCTTTA
TGGAACCTTGGCAGAGCCAGAGTGACACATTCCTAAACCATTAACAGATTTCTA
TA

30 Sequence ID - 606 nt: 556
GGATAATGATACCTCTGACCTTTCTTCCTTTTGGGAAGTACTTGAGTGTGCAGCTG
CATGAGGCCTCAGCAGGAGAGAGATTTTAGGTCCAAGAAGCTATACCAGTAGGACA
AGGCAGGAAAATACTACACTTTCAAGATCAAGCCCCTCTGACTCTCATTTGGAAAC
TGGATGTTTGCTAAGCACCTGCTTCTTAAGGATGCCGAGGGATTTAATGATACTCC
35 CAGAAACCTGGAGAGATTAATGGGGCCTATGGAGAAGTGCTCTGAACTCAGTGTTG
GGACTTGAATAAAATTAACCATTTGTATGTTTTTCAGAACAACTAAGCTGTTTTATA
TTTCATGTGCATGAAAGCCCTAGAACTAAGTTGTGTTATTTCCAGAAATGAAATAG

- 187 -

ATCCACAGTTAGATGATGTGGCCATTAGGAAGTACCAAATTTATAAAAATCACTG
GAGGTCTGTCTGAGCAGTACCTAATAAAATATAGTATACTGAAAGTGAACAGATCT
TTGTCTCTTTCTTTGGCTGCTTGATACTTTATCTGTGTCTGCCGGACAGTGC

5 Sequence ID 607

CAATAAAAGCAGGTTAACCTCAATGATAGCAGTTAAAATGTTCTATCTTATGTATT
TCTTTTAAGTATTACCATTATGGTGCTACTGAGCGTTTCTTTTGGTAAAAAGAAA
AATGCCATGGGCTGCAGTCTTCTTCCATCACTTTTCCCTACCAGGTCCATTAATAT
GCTTATAACACTAGTGCCAGTTATTTTATTTGATAATGCTTATGGTATTTGTATAT
10 TTGTTTGCATTCCAATTTTGTTTAATAATGAGTGTGTAAACTGCATACGTAAATA
AATGTAAATACTAATGTACTGCTGC

Sequence ID 609

TTTATTACCCAAGTTTAACTCTGTCTGGTGATTTGTTGTTGTTGTTGTTGTNG
15 TTGTTGTTGAAGTTCAGGCTGCATGTGGGATAGGTTTGCTCAGGCATACTTCTTAG
GAAGTAGTCACCTGCATGACTGTTTTTGGGATAACTCTTTGAGTATTTGGAGAGGT
CTATTGTAACCTTCTGAAAGGCATTGTTTTTACGTATGAATGTTCTAAAATTCATTC
TAAATGGTCATGAAAAGAAAAGGATTACATTTTAGAATGGCAATAGTCCCTGAGG
ACTATTATGTCTTTTAGATTTCTGTGGGTTTCTAGGAATGTTAGTGTAACCTTANA
20 TTTCCACCTACCTGATTTCTGGATGTGCCTATTGGAACCTTGCTGAGATCTTTTTTTT
TTCCTTAACATGTTGTCCCCTTGACCCGTACTTCGAAACTAAACATATTATTTTAT
TTGCTTACACTTCAGGAGGCAATTGGCAGACACCAGGCCAACAGTCT

Sequence ID 610

GCTCTGACCCCAGTTGGAAATGTATCTGTACTTTGTCCGGCTTCCACTCAAGGACC
25 ATTTATGACATTGCTTGGTGTGTCAGCTGACAGGGGCTCTGGCCACAGCTTGTGGGGA
TGACGCGATCCGCGTGTTTCAGGAGGATCCCAACTCGGATCCACAGCAGCCCACCT
TCTCCCTGACAGCCCCTTGATCAGGCCCATTCCAGGATGTCAACTGTGTGGCC
TGGAACCCCAAGGAGCCAGGGCTACTGGCCTCCTGCAGTGATGATGGGGAGGTGGC
30 CTTCTGGAAGTATCAGCGGCCTGAAGGCCTCTGAGCTACCTCGACTTTGGACAGAG
TAATGACTCCCCAGAAAACGTCATATAAGACTTTACCAGCCCCTGAGAGGACCAGG
AGGAGCATCCTTGACCTTCATTTAACTTGGCTCACTTCTCTTCANACTTGGGTTAGA
AGTGCAGAGCCACAAAATTGCTTTCCTTCCCCGCTTTGACATGAGGCCTTCAGTA
AAG

35

Sequence ID 611

TGCAGGATCCGTCGACT

- 188 -

Sequence ID - 612

nt: 576

5 GAGAAATATAAGATTATGTATAGATCAAATCTACCTCTATTTGGTGTCTGAAAGA
GATGAGGAGAATGGGACAACTTGGAAAGCTTATTTCAAGATAACATTCCTGAGAA
CTTCCCCAATCTTGCTAGAGAGGGCCAACATTAAAATTCAGTAAATGCTGAAAACCTC
10 CAGTAAGATATTTCTTAAGAAAATTATTCCCAAGATATATACTCATCAAATTATCT
AAGGTCAAATGAAGGAAAAAATTTTATAGGCAGCTAGAGAGAAATGTCAGGTCACC
TACAAAGAGAATGGCATAAGACAAAAAGTAGAACTCCCAGCAGAACTCTAAAAGC
CAGAAGAGATTAGGGGCCAATATTTAACAATTCTGAAAGAAATTCCAACAAGGAATT
TCATATCCAGCCAACTAAGCTTCATAATTGAAGGAGAAATAAGATATTTCCAGA
15 CAAGCAAATGCTGATGAAATCCATCACCACCAGACCTGCCTTATAAGAGCTCCTGA
GGGAAGCACTAAATATTGAAAGGGAAGAACTTTATGAACCATTTCAAAAACACATT
TAAGTNCACAAAGCAG

Sequence ID - 613

nt: 341

15 CCTTATTTTACAGGTGAAAAACCACGAATCAGATAGATTTTTATTTGCCCAAGTCA
CATAATATTAAGAACAGGCCAAGTGTGGTGGCTCATGTCTGTAATCTGAGCACTTT
GGGAGGCTAAGGCGGGTGGATTTCTGAGCCTAGGAGTTTGAGATCAGCCTGGGCA
ACATGGCGAAACCTCATCTCTACAAAACATACAAAATTAGTCAGTGTGGTGGTGA
GAGCCTGTAGTCCTGGCTACTCGTGAGGCTGAGGTGGGAGCATCACCTGAGCCTGG
20 GAAGTCGAGGCTGCAGTGGCAACAGAATGGGTAACCTGGACATCAGAGTGAGACCC
TGTCT

Sequence ID 614

25 CTCACACCTGTAATTCCATTACTTTGGAAGGCTGAGAGAGGAGGATCAGTGGAGCC
CAGGAGTTTGAGACCAGCCTGGGCAATATAGGGAGACCTGTCTCTACAAAATGA
AATAGCCAGGCGAGGTGGCATGTGCCTGTGGTCCCAGCTACTTGGGAGACTGAGGT
GGAAGGCTGCCTTGAGCCCAGGAGTTCCAGGCTGCAGTGAGCCATCATTATGCCAC
TGCACTCCAACCTGGGAGACAGAGTGAGAGAGACCCTGTCTCAAACAAACAAACCC
AAAATAGGCCAGGCACAGTGACTCATGCCTGTAATCCCAGCACTTTGGGAGGCTGA
30 AATAGGCGGATCATTTGAGGTGAGGAGTTCAAATTCAGACCAGCCCGGCAACAT
GGCAAAACCACATCTCTACTACAAATAAAAAATTAGTTGGGTGTGGNGGAGCATTC
CTGTAATCACAGCTATTTCAGGAGGCTGAGGCATGANAACCGCTTCA

Sequence ID - 615

nt: 379

35 TAAATTTAAAACATTTTAAATTAGCTGGCATGATGGCATGCACCTGTAGTCCTACCT
ACTTGGGAGGCCAAGGCAGGAAGATTGCTTGAGCCCAGGAGTTTGAGCTTACTGTG
AGCTGTGATCACACCACTGCACTCCAGCCTGGGTGACAAAGGAAGACCGTATTTCT

- 189 -

AAAAAATAAAAAATACAAATACAACTACAACTAGCACTAGACCAACAGTGACTAT
GTACCATGAACTGAGGAATATTATTAATTCCACCATTTCATCTGAGGTTAACAAT
ATGTCAATGACTTAAATAACATCATATCTCTGAGAGTAATTTCTCCTATATTTCCA
TGACAAATGTTAGATAATTTTCCATTTTTTCCATTCAACAAAA

5

Sequence ID 617

TTTTCAGGCATGTCAGAGAAGGGAGGACTCACTAGAATTAGCAAACAAAACCACCC
TGACATCCTCCTTCAGGAACACGGGGAGCAGAGGCCAAAGCACTAAGGGGAGGGCG
CATACCCGAGACGATTGTATGAAGAAAATATGGAGGAACTGTTACATGTTCCGTAC
10 TAAGTCATTTTCAGGGGATTGAAAGACTATTGCTGGATTTTCATGATGCTGACTGGC
GTTAGCTGATTAACCCATGTAAATAGGCACTTAAATAGAAGCAGGAAAGGGAGACA
AAGACTGGCTTCTGGACTTCCTCCCTGATCCCCACTCTTACTCATCACCTGCAGTG
GCCAGAATTAGGGACTCAGAATCAAACCACTGTAAGGCAGTGCTGGCTGCCATTGC
CTGGTCACATTGAAATTGGTGGCTTCATT

15

Sequence ID - 618

nt: 598

GATTAACTTTTCATTTTAAGCTCTTCTCTACTAATTCTGTTCGTATGTTTATTCATT
TTGCGTTGATCATATTTTGTACACCAGGCCTCTTCTCAGTTTTATATGTGTGTTA
ATTTACTCCTTTCAAGAGCCCTATGATACATGAATTTATCTCCATTTTATAGATGA
20 GGAAATTAAGACCTAGAGTTACTGAACTTGCCCAAGGTTATACAGCTGATGGGTAG
GGCCAGAACTTTGCCTCAGAGAATCTGAATTTCCAAAAAATAACCTAAAAGAGAAA
TTTAAGTACTAATTAGTAAGCAAAGAAATGCACATTTAAGGAAGACAGTGCACATT
TAAGGAAGACAGTAACCTTTTATCTATTAGAGAAAAACACACATTCTGTCTTTAAC
ACACACATAAATCTTATATTGGCAGGGATTTTCTTTATTTCAGCAATTATTTATTGG
25 TTGTCTGCTTTGTGGTACACATAAATGCTGGGGATAAACACTTAATAAAATATACT
TCCTTCTCTTGAATATCTTGCACTTTAAGTGGGAAGGTAAGTCAACAGAGTAGAGG
TGATATATCCAAGTGATAGACTGTTTCATTGCCAGTAG

Sequence ID 619

GTTGCCTGAGAGTGACCTTTGCATCTGCCTGTCCAGCCAGCATGGAACCAAAGCGG
ATCAGAGAGGGCTACCTTGTGAAGAAGGGGAGCGTGTTCAATACGTGGAAACCCAT
GTGGGTGTGATTGTTAGAAAGATGGAATTGAATTCATAAGAAGAAAAGTGACAACA
GCCCCAAAGGAATGATCCCGCTGAAAGGGAGCACTCTGACTAGCCCTTGTCAAGAC
TTTGGCAAAGGATGTTTGTGTTAAGATCACTATGACCAAACAGCAGGACCACTT
35 CTTCCAGGCAGCCTTCCTGGAGGAGAGAGATGCCTGGGTTCGGGATATCAATAAGG
CCATTAAATGCATTGAAGGAGGCCAGAAATTTGCCAGGAAATCTACCAGGAGGTCC
ATTGCACTGCCAGAAACCATTGACTTAGGTGCCTTATATTTGTCCATGAAAGACAC

- 190 -

TGAAAAAGGAATAAAAGAACTGAAT

Sequence ID 621

5 TGGTACTGAACCTACGAGTACACCGACTACGGCGGACTAATCTTCAACTCCTACAT
ACTTCCCCCATTATTCTTAGAACCAGGCGACCTGCGACTCCTTGACGTTGACAATC
GAGTAGTACTCCCGATTGAAGCCCCCATTCGTATAATAATTACATCACAAGACGTC
TTGCACTCATGAGCTGTCCCCACATTAGGCTTAAAAACAGATGCAATTCCCGGACG
TCTAAACCAAACCACTTTCACCGCTACACGACCGGGGGTATACTACGGTCAATGCT
CTGAAATCTGTGGAGCAAACCACAGTTTCATGCCCATCGTCCTAGAATTAATTCCC
10 CTA AAAATCTTTGAAATAGGGCCCGTATTTACCCTATAGCACCCCTCTACCCCT

Sequence ID 622

TTTTTCTTGTTTTTGTGTGTCTACCTTGGCATATACTAAAGGAAGGTGTGTATTCA
TTTATTACATGATATCTCTGGGTATAATTATTTACATATATGAATTTGAAAGAAA
15 GATTGAGAGGGATATGTGTGACCTTTGTTTCATTATGATCATTTACATGACTAAAG
ATAAGATCATATGTCTGATTTTTCAGTTTAATGGCAAGTTACTTAAAATAAATGAA
ATATGTTTTTATTGTTTTTCGTGGGTTTGATGCTTTGTGTTTTATTTCAAGTAACTT
GAGAATGCATTGTGTTTTGGTACTGTTTTTTATGAATATCATTA AAAATTTATTTAA
GGAGAGAGTAATTTTGCAATAATATTTTTGATTTATTTGAAAATAAAATTCAAGAT
20 AAATGAAATAATTGAAATTTTCTAAAGAAGGAATTGAATATATTTTTTACATTTGAA
TGA ACTAAGGATTAACTGAACCATTTATATATAGTACTTTCAGAACTGAATGTCTT
AAATGATAAAGCTCTAATTGGTTAAAGTGACTTTCTTTCAAGTCAAAGAACCAGA
AACTGAATAGATGATCTAACTACTGCCACTGAGGTTTTGGATTAGTGAGTATAAAT
TT

25

Sequence ID 624

TGCAGGATCCGTCGACT

Sequence ID 625

30 GACAATCAGAGCAGATCTTGGGCTTCTGTGGCTCATCTCAGCCCTTTATAACTGGC
CTGAGAAGAGGGTTTATCTACTTGTGCAAGTGGCCCAGAAATCTCACTCGTACATG
AGGCTTTGGAACATCCTTGCAAAGGTACGCTGAAAGCAAATTGCTGTTTTCTGGT
GGTTCTGCACGTTTCCTAACTTTTATCATAGTTTGATTTTCATTATTTAAGAAAAA
ATAAAAAATCCAAAGACCATAAGATGGCATTAGATTTTTTACCATTAAATTATTAA
35 TGCCTATTTGGTGCTCATAAAGATTAATCATGTACGCATGTTTCCAATCTTTCTT
TTGCAGTATATTATTTTCTAAAAATTGTTACATGCAAATTTAAACCAAGATTTATC
AGTA

- 191 -

Sequence ID 626

TTGGAAGAAATAAACCAAGGCAGAAAAATTTTAAATGGCCAAAATAAATTGTATTG
CTAACTTAGATGGCCACAGATGGGGGCAGGGGTGGAGAGAGGAGAAATTGAAAACN
CCACAAAGACCCCGCAATGGCTAGAACTTGAAATCTCTGGATATTGCAACAATAGC
5 AGCCTCCTTAAGTCAGCAAAAAGATAAAGATTGATCCAATGTTCTATATTACAGAA
CAGAGCAGATTGTCAATATAGCAAATAAAGTTACCGTTGAGTGGACTGCGCTGTNT
AAGCTGCTTGGTTGGCCTTAAGTGCCGACAATTAAGAGATGAAGGCAATGAGAACT
GAAACAAACATTTAAGTTCAAGACCCAGTTTACTGACACTGGGACTATTACTATAT
CTCTTTGGGCCTCAGTTTACTTATCTGTAACATTAAGAGGTTGGATTACATGATGT
10 CTCACGATTCTTTTTTTTTTATTTAGAGATGGGGTTTTTGCTCTGTTGCCCAGGCTGG
AGTGCAGTGGCATGATCATAGCTCACAGCAG

Sequence ID 627

CCAGCCTGTCACTGGCCTGGCCAAGGAGGAGAGACAGGCCAGGGATTCTGGTCCTA
15 ACTCTACTGGCCACACTGTGTGGCCTGAGACCCCCCTTCCCTCCCAAGCCCCTGC
CTCCGCATCTGCGTGGTGAAGGCCATTGGCCCTCATCGGTGGATCTGCGTTTCCTC
GGGCCTACACTGTCTAGGATTGTGCGGGGCTGGTGAGAGAACAAAGATCTCTTCCGT
GTTCAAGGCAGACTTCCTGCCCCCTGCACCCTGCTCTCTCCAGGCCCTTGAGGTCA
GTGTGAGCCCCAAGGGCAAGAACACTTCTGGAAGGGAGAGTGGATTTGGCTGGGCC
20 ATCTGGATGGAAGGTAAAAAAAAGAAAATCCCTTGAAAGGAGATTGAGGGAAGTTT

Sequence ID - 628

nt: 419

AAGAGAAAGGACTCAGTGTGTGATCCGGTTTCTTTTTGCTCGCCCCTGTTTTTTGT
AGAATCTCTTCATGCTTGACATACCTACCAGTATTATTCCCGACGACACATATACA
25 TATGAGAATATACCTTATTTATTTTTGTGTAGGTGTCTGCCTTCACAAATGTCATT
GTCTACTCCTAGAAGAACCAAATACCTCAATTTTTGTTTTTGAGTACTGTACTATC
CTGTAAATATATCTTAAGCAGGTTTGTTTTCAGCACTGATGGAAAATACCAGTGT
GGGTTTTTTTTTTAGTTGCCAACAGTTGTATGTTTGCTGATTATTTATGACCTGAAA
TAATATATTTCTTCTTCTAAGAAGACATTTTGTTACATAAGGATGACTTTTTTTATA
30 CAATGGGAATAAATTATGGCATTTTTTT

Sequence ID 629

CTGAGAGTCACTGTGTTTTTAGCCAAATCTAAGGGAGAAAATGAATATTGATAGCA
GCATGCTGTAGCCAGCTCCTTAAAGGAAGGATGGTGCCTGGTACAGAGTTAGAGTT
35 AGTGCTTCAGTAAATAATGAATGTGTGCTAGGTAGGTTCTGCTGGGTAGGCTGCAT
GCATTGACCAATTTATTCCTCCTTGTTTCAAAACAGGATTTAAGGGCACTTATATA
TATATATTTTTTAGTTTTTTTAAATGTAAATGAGAGAATAAAGATATATATATATGT

- 192 -

CTATATATGTATATATGTATATATATGTCTATATGTCTATATGTATATATGTCTAT
ATGTATATATGTGTGTGTGTATATATATATATATATATAAGTTTTCTGTTGCTA
GCATAACAACTACCAGAACTTAGCAACTGAAACAACATGAATTTATCTTACGGT
TCTATAGTTCAGAAGTCTAACGTGTCACTGGGATGAAATCCAGGTTTCAACAGGAC
5 TGGGTTCCCTTCTAGCTCATTACAGCTACCTGGCTCATTACAGGTTGTNNGCAGAATA
TACTTCCATGAACTGTAGGGCTGAGACCCCGTTCCTTCCTGGCTATCATCTGAAA
ACTTTC

Sequence ID 630

10 AGGCGCAGCCCAGCCTCGAAATGCAGAACGACGCCGCGAGTTCGTGGACCTGTAC
GTGCCGCGGAAATGCTCCGCTAGCAATCGCATCATCGGTGCCAAGGACCACGCATC
CATCCAGATGAACGTGGCCGAGGTTGACAAGGTCACAGGCAGGTTTAATGGCCAGT
TTAAAACTTATGCTATCTGCGGGGCCATTTCGTAGGATGGGTGAGTCAGATGATTCC
ATTCTCCGATTGGCCAAGGCGGATGGCATCGTCTCAAAGAACTTTTGACTGGAGAG
15 AATCACAGATGTGGAATATTTGTCATAAATAAATAATGAAAACCTAAAAAAAAAAAA
AAAAAAAAAAAAAAAA

Sequence ID 631

TNCACTCACACACTCCCAAACCTTAACAAACACATACATGTGCAGCCAACCCAATG
20 GGCCAGCCTCTTTTATGCTCCTCACATGTTTCCTTTAACTGGAATACCCATGACAG
CTCCCTACATAGTTACTTGTAAACTCCTCCTCTCTGTATAAGTTTTCTGAATTTT
TTTGATAAAATTAAGTTGTGCCACCCCTTTATGCTCTCTTANAACCTTGTTCTGTT
CTCATGGCTGTTCTGCAACGAATCTCATTGTGTTCTCCTACTCAATTACATTCCTG
CGTCTCCCACTAGATGGCAGACTCTTTGAGAGTAGGAGATTCCCTTGTTATCTCTG
25 GATCCCTGGCACTTGCGAGAAAGCCTGTTACGTAATAATTGCTCAACAATTAGTTTT
TAAATAAATGAATTATTTTAAACGCCAAAATTACAATGATTGTGCATTAAGTGA
AAGATGACCATCTAAAAACATAAAGCCATGCTTCATGACATTGGC

Sequence ID 632

30 GACCATTACAGGGAATTTTATAAAAAATGCAGATACTGTCTTGAGCAGATCGAAAT
GCCGATGAGGTGGATGCAATTTCTTTTGTGCAAGCAGTGACGGTGCCCCCCCCCT
CGGGTGTCCGTGCTGTGCCTTAGCTTCCCCAGGTGCCGGGACTCACACCTGCTAGG
GGCTGGGCAAGGCCCGGCTCTGCTTTCTCTGAAGGGCTTGTCCAAGTTCAATTGCC
CTGTTACAGGTGGTCAAGACGTCCGGCCGCCTTGACCCAGGCTACCCTTAGCCAAT
35 ATCCTCTGCCCCTGGGTGGTTGGTGGCTGGGCCTCAGGGTGGGCAACGTTAGGGGT
TTGGCGAAAGCCCCGCCCATGGGATTGAGGGACGGGGCTGCACTCCAACCGTCTGC
ACCTGCTCTTCCCCCACCCTGTGGGACCTCATCTTCACGTGCCATGTGTGCTGAA

- 193 -

GGCCCAGGGCCCAGCAGGGGGCAGTGGCACCTGTTGACGGAAAAGCCGAGGTGCTT
ACCAATGGACCTTCTGGCCCGCCCTCCCCTGTACTTGTGCGGGCATTACAGGGCCCCG
ACCTGTGCCTACCCGCA

5 Sequence ID 633

CAACTGTGTTCACTAGCAACCTCAAACAGACACCATGGTGCACCTGACTCCTGAGG
AGAAGTCTGCCGTTACTGCCCTGTGGGGCAAGGTGAACGTGGATGAAGTTGGTGGT
GAGGCCCTGGGCAGGCTGCTGGTGGTCTACCTTGGACCCAGAGGTTCTTTGAGTC
CTTTGGGGATCTGTCCACTCCTGATGCTGTTATGGGCAACCCTAAGGTGAAGGCTC
10 ATGGCAAGAAAGTGCTCGGTGCCTTTAGTGATGGCCTGGCTCACCTGGACAACCTC
AAGGGCACCTTTGCCCACTGAGTGAGCTGCACTGTGACAAGCTGCACGTGGATCC
TGAGAACTTCAAGCTCCTGGGCAACGTGCTGGTCTGTGTGCTGGCCCATCACTTTG
GCAAAGAATTACCCCCACCAAGTGCAGGCTGCCTATCANAAAGTGGTGGCTGGTGTG
GGCTAATGCCTGGCCCCACAAGTATCACTAAGCTCGCTTTCTTGCTGTCCAATTTT
15 TATTAAAGGTTCTTTGTTCCCTAAGTCCAACCTACTAACTGGGGGATATTATGAA
GGGCCTTG

Sequence ID - 634

nt: 511

TTTTTTAATTTACCAAAATTTGTTGACGTCCCTTGATTTGCTGATAGGGACAATA
20 ATTAAATATTTTCCACTTGTTTTTATAAAAACCTGTAATGGTGATTTGTTTAACAGA
TGTTGACTTAGCACCTTCTCTCTTTTTTTTTTTTTTTTTTTTGGAGTTGGAGTCTTGC
TCTGTCACCCAGCTGGAGTGCAGTGGCACGATTTCCGGCTCACTGCAACCTCCGCCT
CCCAGGTTCCGGGCGCTTCTCCTGCCTCAGCCTCCCANATAGTTGGGATTACAGGTG
CATGCCGCCACNCCTAGCTAATGTTTTTTGTATCTTGGTANANATGGNGTTTCCACC
25 TTGTTGCCCATGCCGCTCTTGAACCTCCTTGGCCTCCCAAAGTGTTAGGATTACAGG
CGTGAGCCACTGTGCCTGGCCCCAATTTANCACCTTACTGGGTGCTGAGGCTGTGA
GCCATAGTAGAATGCATGTGATCCAGGGCCTTGCTGAATTCATGGGCTAATAGGGA
GCCTGAC

30 Sequence ID - 635

nt: 592

TGAGCGTTGGGCTGTAGGTGCTGTGCTGTGTGATCCCCAGAGCCATGCCCGAGA
TAGTGGATACCTGTTTCGTTGGCCTCTCCGGCTTCCGTCTGCCGGACCAAGCACCTG
CACCTGCGCTGCAGCGTCGACTTTACTCGCCGGACGCTGACCGGGACTGCTGCTCT
CACGGTCCAGTCTCAGGAGGACAATCTGCGCAGCCTGGTTTTGGATACAAAGGACC
35 TTACAATAGAAAAAGTAGTGATCAATGGACAAGAAGTCAAATATGCTCTTGAGAA
AGACAAAGTTACAAGGGATCGCCAATGGAAATCTCTCTTCCCTATCGCTTTGAGCAA
AAATCAAGAAATTGTTATAGAAATTTCTTTTGAGACCTCTCCAAAATCTTCTGCTC

- 194 -

TCCAGTGGCTCACTCCTGAACAGACTTCTGGGAAGGAACACCCATATCTCTTTAGT
CAGTGCCAGGCCATCCACTGCAGAGCAATCCTTCCTTGTGAGGACACTCCTTCTGN
GAAATTAACCTATACTGCAGAGGTGTCTGTCCCTAAAGAACTGGTGGCACTTATGA
GTGCTATTTCGTGATGGAGAAACACCTGACCCA

5

Sequence ID - 636

nt: 572

CTTANAAGAGTTGCTCATTACACCCACGCCCTTGCCCAAGGCTGGCCCACTCAGA
GCGAAACTTAACTTTTGTCTGGATGGGAAGAGAAGTAAGTCTACCCCGAGGTTGCC
ATGTTGAAGAGTGAGAGGTCCAAGTGATTCTGTGCATTGAAACCAAGACACCCAC
10 CCAGAACACTTCTTCCCTCCCTCAGCCCAAACCAAAGGCTGGGGTTCTCATCTCCA
AGTGGCTGTTCTCCAACCTTCCCAAGCCGCTTGCACTCCCGAGACTGGACTACTGT
GGCGGTTAGGTTAGATTTGAAGACGGGGCCAGGCTGGGTATGAACGGGTGCAGCC
CTCTTCTCCTCTTCCCCCCCCACATCTCTCATGAGAGAGGTAGTGGCATTTCCTTCT
CAGGGAGCTTCAATGGGAAAGGTCTCGAAAGCTTCAGGAGGAGCAGAATACCAACG
15 CAGGGGGATGGCTGTAACGATCTCACCGTCTCCTAACCTCAGTCCCTTTTTTGAGA
GTGAATGGTGGAGGGTGGGAAAGGGACCCAAATTTGTAGATCTCTTTGTCTGGGGG
AGGGGAANGATG

20 Sequence ID - 637

nt: 482

TTAAAACAGGCGCAGGGGTAAAAATGAGAATGAATCTGAAAAAGAGAGTTGGTGT
TTAAAGAGGATGGACAAGAGTATGCTCAGGTAATCAAAATGTTGGGAAATGGACGA
TTGGAAGCATTGTGTTTTGATGGTGTAAAGAGGTTATGCCATATCAGAGGGAAAT
GAGAAAAAAGGTTTGGATAAATACATCAGACATTATATTGGTTGGTCTACGGGACT
ATCAGGATAACAAAGCTGATGTAATTTTAAAGTACAATGCAGATGAAGCTAGAAGC
25 CTGAAGGCATATGGCGAGCTTCCAGAACATGCTAAAATCAATGAAACAGACACATT
TGGTCCTGGAGATGATGATGAAATCCAGTTTGACGATATTGGAGATGATGATGAAG
ACATTGATGATATCTAAATTGAACCAAGTGTTTTACATGACAAGTTCTCTGAGGA
TGGTTCTACAGTTGGGATTTTGGCCATCATCAAC

30 Sequence ID - 638

nt: 545

TTTGAAGGCAAAGAGGGATTAACTGTGCTGGCATCATGTAAGGAGACTTGATAGA
TAAGAAAAAGCTTTACCTAAGTTTTGAAGAATAGGTTTTTCATAATGGAAAATTTA
AGGGAAAAATCTCCAAAAAGTGCTACTCAAGTTTTATCCATTTGTATTTCCAACA
CAGCCTAGGACAGTACCTGCACATAGTAGGTGATTAATAAAAATTTAGAAAGCATT
35 AATACTAAAGAGGAAAAATAGCAATGGCAAGAAAACACATGTAGGGAACACATGTA
GCCAAAAAATAATATATAATCAGAGAAATAATAGGACTTCTGGAAAAAAAAGATGA
GATCAGATTGGTTAGGATCTTTACTAACATGACAAGAGCATGAATTTTTTTTCTGT

- 195 -

AGATAATAAGTATGAAAGAATTTTAGCTTAAAAATTAGCATAATTTGGATCCACAT
ATGCAAATCAATGAATGTAATTCATAATATAAACAGAACTAAACACAAAAACCACG
TGATTATCTCAATAGACACAGAAAAGGCCTTCAAAAAAATT

5 Sequence ID - 639

nt: 624

GACACACGAGCATATTTACCTCCGCTACCATAATCATCGCTATCCCCACCGGCGT
CAAAGTATTTAGCTGACTCGCCACACTCCACGGAAGCAATATGAAATGATCTGCTG
CAGTGCTCTGAGCCCTAGGATTCATCTTTCTTTTCACCGTAGGTGGCCTGACTGGC
ATTGTATTAGCAAACCTCATCACTAGACATCGTACTACACGACACGTACTACGTTGT
10 AGCCCACTTCCACTATGTCCTATCAATAGGAGCTGTATTTGCCATCATAGGAGGCT
TCATTCACTGATTTCCCCTATTCTCAGGCTACACCCTAGACCAAACCTACGCCAAA
ATCCATTTCACTATCATATTCATCGGCGTAAATCTAACTTTCTTCCCACAACACTT
TCTCGGCCCTATCCGGAATGCCCCGACGTTACTCGGACTACCCCGATGCATACACCA
CATGAAACATCCTATCATCTGTAGGCTCATTCATTTCTCTAACAGCAGTAATATTA
15 ATAATTTTCATGATTTGAGAAGCCTTCGCTTCGAAGCGAAAAGTCCTAATAGTAGA
AGAACCCTCCATAAACCTGGAGTGACTATATGGATGCCCCCACCCTACCACACAT
TCGAAGAA

Sequence ID 641

20 CAAGATGACAAAGAAAAGAAGGAACAATGGTCGTGCCAAAAGGGCCGCGGCCACG
TGCAGCCTATTCGCTGCACTAACTGTGCCCGATGCGTGCCCAAGGACAAGGCCATT
AAGAAATTCGTCATTCGAAACATAGTGGAGGCCGCAGCAGTCAGGGACATTTCTGA
AGCGAGCGTCTTCGATGCCTATGTGCTTCCCAAGCTGTATGTGAAGCTACATTACT
GTGTGAGTTGTGCAATTCACAGCAAAGTAGTCAGGAATCGATCTCGTGAAGCCCGC
25 AAGGACCGAACACCCCCACCCCGATTTAGACCTGCGGGTGCTGCCCCACGTCCCCC
ACCAAAGCCCATGTAAGGAGCTGAGTTCTTAAAGACTGAAGACAGGCTATTCTCTG
GAGAAAAATAAAATGGAATTGTACTTAA

Sequence ID 642

30 TGCTTGGCCCTCTACCTCCTGCCCTCTTCTGTTCATCTCCCAACCACTGCACTCT
TGATTTTTATACCACACAGAAGGTAAGAAAATTCTAGGAACCCTAAGGATCAATCC
TCTCCATTTTCACTCAAATGCCTGGGGCCCAGCTCTGCAATGACTGACTCCAGGGC
CTCTTTCTCACTGCCAGCATAGAAGTCAGGGGAGCCAGCTGGGCCCTGCGGTCAG
GAAGGTTCTCATTTTTGGAGCATTCCCTGAGCCCAGATCATAGGAGCAGCTGTCCC
35 TGGTGGGACACAGGAGTCATGACTCCTACCCTCCACCCTCCACACCCACCAGGCAT
TTAGCAGTCTGTCTATGCAAGACAGATGAATTCTCAGCCAGGATACCTCAAGGCA
GGCAAAGGTGAGTGGAGGGAAAATTCACAAACATTACAGGGTGTGTGGTGCTGGCAT

- 196 -

CACCATGGCCAAATCCAAGAGGTCTTCCTGGAAGAGGGCCCAAACCTGGAACCAAAA
GAATGCTGTCAGCAGTTGGAATAGAGCTGTGAATT

Sequence ID 643

5 CTTTCCAAGAGGAATCCTCGGCAGATAAACTGGACTGTCCTCTACAGAAGGAAGCA
CAAAAAGGGACAGTCGGAAGAAATTCAAAAGAAAAGAACCCGCCGAGCAGTCAAAT
TCCAGAGGGGCCATTACTGGTGCATCTCTTGCTGATATAATGGCCAAGAGGAATCAG
AAACCTGAAGTTAGAAAGGCTCAACGAGAACAAGCTATCAGGGCTGCTAAGGAAGC
10 AAAAAAGGCTAAGCAAGCATCTAAAAAGACTGCAATGGCTGCTGCTAAGGCACCTA
CAAAGGCAGCACCTAAGCAAAAGATTGTGAAGCCTGTGAAAGTTTCAGCTCCCCGA
GTTGGTGGAAAACGCTAAACTGGCAGATTAGATTTTAAATAAAGATTGGATTATA
ACTCT

15 Sequence ID 644

CTTTGATAGAGAAGAAAATTCTCCTAGGATACAAGAGCCTCAACATTTTAAAGATT
TTCTGCATCTCAAAAGCGTAGGCTCCTTGCTGGGCAAGGTGAGCCTCTGTGAGTCC
TCATAGGACCGAGCAAATCTGATTCACCCAGAAAATCCAATATCGAAGCTGAGCT
TTGGCCTGAGCGGGTTCCATTTCCTCCCAGATCCTATTTAGGAAGTGTCTCCTGA
20 CAACCTCCAAAAGGTGCTAACATGCAACGTTCTGAAGGGTTATTGCTCAAAAACAA
GATTTTCTTGTGGTCAAGACTCTGCGAGCCTCGAACACGATGAATCCGCTCGAAT
GGGCTTGGGCTTTGCCCCGGGTGGCGCACGCTCACACGCTGGAAGCACAGCTTTGAC
GATCTCCACACACGCACAGGCACACAGCCACAGATGATGCCGGCTCATTCTCAGG
GGGTGTCTAAGTTCTGCTTTAAATATTTACCCCTAATTGTACAAACAATAGGGGC
25 ATGAGCCTGGTACTCGATAAATGGGGACTTNCTTAAAA

Sequence ID - 645

nt: 649

CTACAGCCTGGGCAGCGCGCTGCGCCCCAGCACCAGCCGAGCCTCTACGCCTCGT
CCCCGGGCGGCGTGTATGCCACGCGCTCCTCTGCCGTGCGCCTGCGGAGCAGCGTG
30 CCCGGGGTGCGGCTCCTGCAGGACTCGGTGGACTTCTCGCTGGCCGACGCCATCAA
CACCGAGTTCAAGAACACCCGCACCAACGAGAAGGTGGAGCTGCAGGAGCTGAATG
ACCGCTTCGCCAACTACATCGACAAGGTGCGCTTCCTGGAGCAGCAGAATAAGATC
CTGCTGGCCGAGCTCGAGCAGCTCAAGGGCCAAGGCAAGTCGCGCCTGGGGGACCT
CTACGAGGAGGAGATGCGGGAGCTGCGCCGGCAGGTGGACCAGCTAACCAACGACA
35 AAGCCCGCGTCGAGGTGGAGCGCGACAACCTGGCCGAGGACATCATGCGCCTCCGG
GAGAAATTGCAGGAGGAGATGCTTCAGAGAGAGGAAGCCGAAAACACCCTGCAATC
TTTCAGACAGGAAATCCAGGAGCTGCAGGCTCAGATTCAGGAACAGCATGTCCAAA

- 197 -

TCGATGTGGATGTTTCCAAGCCTGACCTCACGGCTGCCTTGCGTGACGTACGTANC
AATATGAAAGTGTGGCTGCCAAAAACCTTGCGAG

Sequence ID - 646 nt: 600

5 GAGATGTCTCGCTCCGTGGCCTTAGCTGTGCTCGCGCTACTCTCTCTTTCTGGCCT
GGAGGCTATCCAGCGTACTCCAAAGATTGAGGTTTACTCACGTCCATCCAGCAGAGA
ATGGAAAGTCAAATTTCTGAATTGCTATGTGTCTGGGTTTCATCCATCCGACATT
GAAGTTGACTTACTGAAGAATGGAGAGAGAATTGAAAAAGTGGAGCATTGAGACTT
GTCTTTTCTGCAAGGACTGGTCTTTCTATCTCTTGTACTACACTGAATTCACCCCCA
10 CTGAAAAAGATGAGTATGCCTGCCGTGTGAACCATGTGACTTTGTACAGCCCAAG
ATAGTTAAGTGGGATCGAGACATGTAAGCAGCATCATGGAGGTTTGAAGATGCCGC
ATTTGGATTGGATGAATTCCAAATCTGCTTGCTTGCTTTTTAATATTGATATGCT
TATACACTTACACTTTATGCACAAAATGTAGGGTTATAATAATGTTAACATGGACA
TGATCTTCTTTATAATTCTACTTTGAGTGCTGTCTCCATGTTTGATGTATCTGAGC
15 AGGGTGCTCCACAGGTAGCTCTAGGAGGGCTGGCAACTTA

Sequence ID 647

CGAATGTGCAGGTTTGTTACATAGGTATATATATGCCATGATGGAAATATTTATTT
TTTTAAGCGTAATTTTGCCAAATAATAAAAACAGAAGGAAATTGAGATTAGAGGGA
20 GGTGTTTAAAGAGAGGTTATAGAGTAGAAGATTTGATGCTGGAGAGGTTAAGGTGC
AATAAGAATTTAGGGAGAAATGTTGTTTATTATTGGAGGGTAAATGATGTGGTGCC
TGAGGTCTGTACGTTACCTCTTAACAATTTCTGTCCTTCAGATGGAAACTCTTTAA
CTTCTCGTAAAAGTCATATACCTATATAATAAAGCTACTGATTTCCAAAAA

25 Sequence ID 648

AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Sequence ID - 649 nt: 425

CAAAAAACGAAGAAAAGTGACGACAGTCTGAGGGACTTATGGGAGATCATCAAGT
30 GAACCACTATATGTGTAATGTAAGTCTTGGAATGAGAAGAGAGAAGGAGAAGGAGG
AGAGAGCTTATTTGTAGAAATAATGGCTGAAAACATCCCAAACCTTTCCTTTTTTTG
AGGAAAGAAATAGGCATACAAGTTCAAGAACTCAAGGAACTCCAGAGAGGACAAT
TCTAAAGACACCCCTCTAACATACATTATAATCAAATTGTCAAAAGTAAATACA
AAGAGAATCTTTTAAATTGACAAGAGAAAAGCAGCTGGTCACGTTCAAGGGAGTTC
35 TATAAGAATTTGAGCAGATTTCTCAGCAGAAACCTTGCAGGCCAACAGGCAGTGGG
ATGATACATTCAAAGTGCAAAAAAAAAAAAAAAAAA

- 198 -

Sequence ID 650

CGAGAGTTTACCAGTNGCCTAATAATGCAATAAAAAATGCTTTGAGATAGCTAACN
GCCCATAAACAACTCAAATTGCTTATAAAGTTTCTTCCCATGTTCCCATTTGAT
5 GAAAAGTCTTACATCACATATAACTGGGAAGCAGGGGTCCCTCCTCAATTTTCAGA
CATTTTGAAAGGATGACAGTTCTGTTTGTAGATGAGTAAACCTCTATATTCATAA
GTTCTAAAATCCTTCATTATGAGGGATTCAAAGTATTTATAAAAAACACTGCCCTCT
AAAAATTTCTCAGATCTGAAGTATGGNCTTGGNCCTGAATATACAGTGTTATCCT
ATGTTTAAAGGGTGATCCAGACATGAGACGCAACTAGTTGGTGCATAAGAAGGCC
CCACTTGGCTATTTTCATATCTACCTACAATTGACCAAAAAAATTTTTTAGGCCAG
10 CAATTATTATTTAGCTTCGCTCTTTCTAGTGCAAGAACTGCAGGCTGGATCAGTA
GTTCAACAGCTAAACAGTCATAAAATAGTCATTGGCATGTTAAATTTCTTTCAATG
CTTCAAAGATAAATTCCAATTCTATTACTTATTATTGNGACNGNATTACTAAAC
AGGTAAGGATGGGAATA

15 Sequence ID - 651

nt: 251

CTTTGGGAGGCCGAGGCGGGCGGATCACTTGAGGTCAGGGGTTCGAGACCAGTCTG
GCCAACATGGTGAAACCCCAACTCTACTAAAAATACAAAAGTTAGCCAAGTGTGGT
GGCAAGTGCCTGTAATCCCAGCTACTCGGGAGGCTGAGACAGGAGAATCACTTTGA
ACCTGGGAGGCGGAGGTTGCAGTGAGCCAAGATCGTGCCACTGCACTTCAGCCTGG
20 GCAACAGAGCAAGATTCCGTCCATCTC

Sequence ID 652

CTTTCCTTCAGCCTTGCAGACACCTAAACATCATGTAATTACCTAAGGAATTCCCAA
GTGCCTCTTCCAGGTTATACGTGTAAATAGCTGTTTTTATGCAAGATTAGTTAGAT
25 ACTGCTCTTTACAGGATGAGTGGTGTGTCTTTGGCTGGGGGGGNCTTAAATGTGT
TTCTAATGTGTGTGTCAAATAATTACCTGTAAACAGACTGCCAATCTGGCTGAAG
CCAATGCTTCTGAAGAAGATAAAATTAAAGCAATGATGTCGCAATCTGGCCATGAA
TACGACCCAATCAATTACATGAAGAAACCTCTAGGTCCACCACCTCCATCTTACAC
GTGTTTCCGTTGTGGTAAACCTGGACATTATATTAAGAATTGCCCAACAAATGGGG
30 ATAAAAAATTTGAATCTGGTCCTAGGATTAAAAAGAGCACTGGAATTCACAGAAGT
TTCATGATGGAAGTGAAAGATCCTAATATGAAAGGTGCAATGCTTACCAACACTGG
AAAATATGCAATCCAATATAGATGCAGAAGCATATGCAATTGGGAAGAAAGAGAA
ACCTCCTTNTTACCAGAGAGCCATCTTNTTTCT

35 Sequence ID 653

GTTGTGACTCGTTGGCATGTGATCTGAAGTTCCTGCCCTGCAGCTGACGAGCCAGT
GTTTCAATAATTAAAAACAACCTCAACTCACTGTCCTCCTGCCTTGAATTTGATCAT

- 199 -

TGCGCTTTGCATGTATGTATCACAATACCACATGTACCCCATAAATATGTACAAAG
ATTATGTGTCAATAAAAAACAAAAATTAAAAATCCCAATTTTTA

Sequence ID 654

5 GTTGCTAGTAGCGGCAGGAAGATGTCAGGCTCACTTTCCTCTGATTCCCGAAATGG
GGGGAACCTCTAACCATAAAGGAATGGTAGAACAGTCCATTCTCGGATCAGAGAA
AAATGCAGACATGGTGTACCTGGATTTTTTTCTGCCCATGAATGTTGCCAGTCAG
TACCTGTCCTCCTTGTTTCTCTATTTTTGGTTATGAATGTTGGGGTTACCACCTGC
ATTTAGGGGAAAATTGTGTTCTG

10

Sequence ID 655

GTCCCCGGGAATCGCGGCCGCGTCGACGGTTTATTTTCAGTGCTTGAAGATACATT
CACAAATACTTGGTTTGGAAGACACCGTTTAATTTTAAGTTAACTTGCATGTTGT
AAATGCGTTTATGTTTAAATAAAGAGGAAAATTTTTTGAAAAAAAAAAAAAAAAAA
15 AA
AAAAAAAAAAAAAAAAAATTTTT

Sequence ID 656

TAGAGGCCTGAATAGGTAGACAATGGCAGCAGCGTTTTTAATCACAGTCCTATTCA
20 TGCCCTAATTCGGGAGTGATGATTAAAGGACATTAGAGGGAGCACTTTGACATCTG
ATCCTTTGAACTGACGTCTGTGCAGGCTGCACTCCATAGAGCTCACTTGGCCAAAC
TGATTTCTTTAAATAAAGTGCTGTGATTTCCAATGTAGGAAATATTACATTAGAGC
CTATTGAAATGATTAGGAATTGAGGAGCTTTTCTTTAGGTGGGAATGTGGTGTATG
CTGTATACTCACAAAAGTGAGATCATTAATATTGCATGTACTACTTTGAATATCAG
25 GGACCACAGAGAAATAGCATGAGAAACGCCTTCCTGCAGTCATGCACTTAAAATGA
ATATGAACAAAAATGTGGAACCTCTGCTGTCATAGCTCTCCG

Sequence ID 657

GGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCGCAGATAAGTTTTTTT
30 CTCTTTGAAAGATAGAGATTAATACAACCTTTAAAAAATATAGTCAATAGGTACT
AAGATATTGCTTAGCGTTAAGTTTTTAACGTAATTTTAATAGCTTAAGATTTTAAG
AGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAAAAGGTTTCT
AAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAATGTATTTAAA
AGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACCAATAGAAGGGCAAT
35 GCTTTTAGATTAAAATGAAGGTGACTTAAACAGCTTAAAGTTTA

Sequence ID 658

- 200 -

GACCTTTGAGAAAATTAATTTAAATCCTAGAACTTTGGGTGAACCGAAGAAATTTA
TAATATTTGTTTAGTTAATAACAGATAAAAAGGAAAGATTCAAGCCTATTGGATGA
GAATTTGTACATTATTTTAGAGCTAATAATAATGGTTTTTCAGTTTAGTGAGGATTT
AAAAAATGTTTTTGAATCAAACTTTTTTTCTTTATAATCCTTTTTTAACCTCAG
5 GAAATAAGGTATTATGAAATCCACACACTGTTACCTCCTTAAAGTATGAGGATACT
TCCCACTGTTTGGTCCACTAGTGGCTGATTATTTTGTGTTGTGGATTATTTGTAATT
TTCTTTTTAATTCTTCCCTTAAAGAGCATGGCATTGAGTCACAGACCTATATTTG
AATCCTGTCATTTACTAGCGTTTTGACCTTGAACAATTATGCTCAGAGTCTCAGTT
TTTTCTTGTAAGTGATGATGATACTACTTAACTCACAGGGTTGTAGTGAAGATCA
10 AATGAGATCATGTCTGTANAACACCCTGCCCCGGCACTCAATAAGTATTAATAGGAA
CCCATATACCTC

Sequence ID 660

TGTTTTTATTTTTTTAAAGGTATAAACACCAAAAAAAAAAATTAAACATTGTATGAAG
15 ATGGAAAATAAGAAGATGCACTTTCTGTAACTTTGTCTAAGGATTTAAATTACTAA
CTTATGAACTCCAATTTGAATTGAACTTAACTATCGGCTTTCTTACTGGTAAAATT
ATATGGTTTATTTTTAAATGCGTACATATTGACCAATGGCCTCTGAAAAAGCACATT
TTAGATACTGAAATTGAAGGAAAGAAAATGCATCTTCAAACATTTTTTGGAAATCTC
ACCACATATACTTTGTTANATTTGTGTATTGTAGGGTGTGTTTTGTATTTTTGT
20 ATTGTATATGAACTTTTTTTAAATGTGACAGTTAAACACATCTTTAAAAGCATAGT
CACAGACAAAAGCATACAGTATAAAAATTTCTTGAAAACCTCACAATATTATAT
TTGGAGGCAGCTTCAGACTGTTTTATTGG

Sequeunce ID 661

CTCTGGCACACATTAGTTCCTCTTATATTACATTGATATAAGCAAGTCATATGGAT
25 TTATCTGAGTGTAAGGAGAGCTGGAAAAAATAGTTTCTAGCAGGTCAGCCACCTCC
CAGTGAGGGCTGCATACCATAGAAGGGGAGAATGAATTTTGGGAAAACAGGTAATT
ATCTCTGTCACAGAAGGGGATGAAAAGTATGGTAGTTACNCAAGTTANACATCTGT
ATGGAAAATACCACTTGGTTCTACAAATGNNG

30

Sequence ID - 663

nt: 627

GCCTCCCGGGTTTCAGGGATTTCTCCTGCCTCAGCCTCCTGAGTGGCTGCATTGCAG
GCACCTGCCACCACGCCTTGCAAATTTTTGTGTTTTTAGTGAGATGGGGTTTTGC
CATGTTGGCCAGGCTGGTCTCGGACTCCTGACCTCAGGTGATCCGCCCCGCTCAGC
35 CTCCCAGAGGGCTGGGATTACAGGCGTGAGCCACTGTGCCTGGCCCCAAGTTTTGC
ATCTTTTAATGCCCTCTGAACAAATACATAGAGAAAACCTCAGAACAAATTAAAC
CTGCAGAGCAACAGTGTCTCCATGTCTTAGGTTTCAAGTTTGCCTCTAAAATTCT

- 201 -

AATCCATATTTTTCTACTTCTCAGATAATTTATGTGTGTGTACTCTTCCTAGACGT
ACAAGAGACTTTTAAATGCTAAATATTTGTCAGTGCTTAACAAAACTCAATTTCA
CATTACTCATATTGTTTTTGTTTTAATTGAATGTGAATTAAATTTTTATTAGTTAT
TTGATTTGGAATGTTATGTATGCCATTAACACTATTAGGGGAATCTCTAGCATTTTC
5 TGTATTTTTTAAAGAATTTGATTCTTTTGTANATTCTGCCTGTGTGGCATTTTAAAC
ATGTGTGACAT

Sequence ID - 665 nt: 345
ACCGGCGACATGGCCAAACGTACCAAGAAAGTCGGGATCGTCGGTAAATACGGGAC
10 CCGCTATGGGGCCTCCCTCCGGAATAATGGTGAAGAAAATTGAAATCAGCCAGCACG
CCAAGTACACTTGCTCTTTCTGTGGCAAAACCAAGATGAAGAGACGAGCTGTGGGG
ATCTGGCACTGTGGTTCCTGCATGAAGACAGTGGCTGGCGGTGCCTGGACGTACAA
TACCACTTCCGCTGTCACGGTAAAGTCCGCCATCAGAAGACTGAAGGAGTTGAAAG
ACCAGTAGACGCTCCTCTACTCTTTGAGACATCACTGGCCTATAATAAATGGGTTA
15 ATTTATGTA

Sequence ID - 666 nt: 252
ATAATTCAGAACTTCTTCATATGCTCGAGTCTCCAGAGTCACTCCGTTCTAAGGTT
GATGAAGCTGTAGCTGTACTACAAGCCCACCAAGCTAAAGAGGCTGCCCAGAAAGC
20 AGTTAACAGTGCCACCGGTGTTCCAAGTGTTTAAAATTGATCAGGGACCATGAAAA
GAACTTGTGCTTCACCGAAGAAAAATATCTAAACATCGAAAACTTAAATATTAT
GGAAAAAAACATTGCAAAATATAAAAT

Sequence ID 669
25 TTAATTTTAACCAGNGAAATTGACCTGCCCCGTGAANAGGCGGGCNTGACACAGCAA
GACGAGAAGACCCTATGGAGCTTTAATTTATTAATGCAAACGGTACCTAACAAACC
CACAGGTCCTAAACTACCAAACCTGCATTAAAAATTTCCGGTTGGGGCGACCTCGGA
GCAGAACCCAACCTCCGAGCAGTACATGCTAAGACTTCACCAGTCAAAGCGAACTA
CTATACTCAATTGATCCAATAACTTGACCAACGGAACAAGTTACCCTAGGGATAAC
30 AGCGCAATCCTATT

Sequence ID 670
GGCTGATTCTTGAGCTATAAAAGCATAATTGCTTTATATTTTGGATCATTTTTTAC
TGGGGGCGGACTTGGGGGGGGTTGCATACAAAGATAACATATATATCCAACCTTTCT
35 GAAATGAAATGTTTTTAGATTACTTTTTCAACTGTAAATAATGTACATTTAATGTC
ACAAGAAAAAATGTCTTCTGCAAATTTTCTAGTATAACAGAAATTTTTGTAGATG
AAAAAATCATTATGTTTAGAGGTCTAATGCTATGTTTTCATATTACAGAGTGAAT

- 202 -

TTGTATTTTAAACAAAAATTTAAATTTTGGGAATCCTCTAAACATTTTTGTATCTTTA
ATTGGTTTATTATTAAATAAATCATATAAAAAATT

Sequence ID 671

5 CAGGAAGTCACCTGGGATTGGCTGCCTCACCCACTCACAGTGCCATCCCTGCCCCA
GGCCTCCCAGTGGCAATTCCAAACCTGGGTCCCTCCCTGAGCTCTCTGCCTTCTGC
TCTGTCTTTAATGCTACCAATGGGTATTGGGGATCGAGGGGTGATGTGTGGGTAC
CTGAAAGAACTACACCCTACCTCCACCACCTTACCCTCACCTGGAGAGCAGTTAT
TTCAGAACCATTCTACCTGGCATTTTATCTTATTTAGCTGACAGACCACCTCCACA
10 GTACATCCACCCTAACTCTATAAATGTTGATGGTAATACAGCATTATCTATCACCA
ATAACCCTTCAGCACTA

Sequence ID 672

CAGGAAGTCACCTGGGATTGGCTGCCTCACCCACTCACAGTGCCATCCCTGCCCCA
15 GGCCTCCCAGTGGCAATTCCAAACCTGGGTCCCTCCCTGAGCTCTCTGCCTTCTGC
TCTGTCTTTAATGCTACCAATGGGTATTGGGGATCGAGGGGTGATGTGTGGGTAC
CTGAAAGAACTACACCCTACCTCCACCACCTTACCCTCACCTGGAGAGCAGTTAT
TTCANAACCATTCTACCTGGCATTTTATCTTATTTAGCTGACAGACCACCTCCACA
GTACATCCACCCTAACTCTATAAATGTTGATGGTAATACAGCATTATCTATCACCA
20 ATAACCCTTCAGCACTAGATCCCTATCAGTCCAATGGAAATGTTGGATTANAACCA
GGCATTGTTTCAATANACTCTCGCTCTGTGAACACACATGG

Sequence ID 673

GGGTTTTCTTTTCGGAAGCGCGCCTTGTGTTGGTACCCGGGAATTCGCGGCCGCGTC
25 GACTGCTAAACAGAATACTGCTATTTTGAGAGAGTCAAGACTCTTTCTTAAGGGCC
AAGAAAGCCACNTGNNCCCTNGGNCTAATCTGGCTGAGTAGTCAGTTATAAAAGCC
NTAATNGCTTNNTNTTTGGNNTCNTTTTTNNCNGGGGNCGGNCTTGGGGGGGGTTG
CNTCCAAAGATANCATNTNTTCCAACCTTNTNAANNNAANNNGTTTTAAATCCCT
TTTCCNCCNGAAAANANNGCCCTTTAAGNGCCNCAAAAAAANNGTNTTCTGCAN
30 NTTTTCTANTATNACAAANNTTTTNGTAGAANAAAAATTTTTTTTTTAGNGGCTACC
CTTTNTTTNTTANNCANNGGAGTTTNTTTTTTACAAAAAANATTTGGGNCCCCT
CCACAACCTTGGGTCTNTAATNGGGGGGTTTTTAAATAAANCNTNTNTAAATCCCC
CNNNNNNNNNCNNNNNNNNNCCNNNNNNNNNNNNNNNCCNNNNAAAAATTTTTNC
TCCCCCNCCCTTTTTCTTCTGCGGCCCCAATTTAAGCCCGGCGCTTGGGGCAA
35 ATCCCCCTTTAGNGGGGGGGTTTTANAAAAACCNGGGGCGGGGNTTTAAACCNCGG
GGNNNGGGGAA

- 203 -

Sequence ID 674

ACCTCTAGCATCACCAGTATTAGAGGCACCGCCTGCCAGTGACACATG-TTTAAC
GGCCGCGGTACCCTAACCGTGCAAAGGTAGCATAATCACTTGTTCCTTAATTAGGG
ACCTGTNTGAATGGCTCCACNAGGGTTCACCTTGTCTCTTACTTTTAACCAGTGAAA
TTGACCTGCC

Sequence ID - 675

nt: 591

GTATAGAAAATAATGTCCCCAGNGCATAGAAAAAATGAGTCTCTGGGCCAGTGAAT
ACAAAACATCATGTTCGAGAATCATTTGGAAGATATACAGAGTTCGTATTTTACAGCTTT
10 GTTTATCCTTCCTGTTAAGAGCCTCTGAGTTTTTTAGTTTTTAAAAGGATGAAAAGCT
TATGCAACATGCTCAGCAGGAGCTTCATCAACGATATATGTCAGATCTAAAGGTAT
ATTTTCATTCTGTAATTATGTTACATAAAAGCAATGTAAATCAGAATAAATATGTT
AGACCAGAATAAAATTAATTATATTCTGGTCTTCAAAGGACACACAGAACAGATAT
CAGCAGAATCACTTAATACTTCATAGAACAAAAATCACTCAAAACCTGTTTATAAC
15 CAAAGAATTTCATGAAAAGAAAGCCTTTGCCATTTGTCTTAGAAAGTTATTTTTTTA
AAAAAAAATCATACTTACTATTAGTATCTATGGAAGTATATGTAACAATTTTTTATG
TAAAGGTCATCTTTCTGTGATAGTGAAAAAATATGTCTTTACTAAGTTGAAATGAA
TACTTTCTGNCTTTGCTAATGGATAGTTATT

20 Sequence ID 676

CTCAATTCTACTAAAAAGCCCCCAAGAAAAGCGAATGAGAAAACAGAGTCATCCT
CTGCACAGCAAGTAGCAGTGTACGCCTTAGCGCTTCAGCTCCAGCTCAGATTCC
AGCTCCTCCTCTTCCTCGTCGTCTTCAGACACCAGTGATTAGACTCAGGCTA
AGGGGT CAGGCCAGATGGGGCAGGAAGGCTNCGCAGGACCGGACCCCTAGACCACC
25 CTGCCCCACCTGCCCCCTTCCCCCTTTGCTGTGACACTTCTTCATCTCACCCCCCCC
TGCCCCCCTCTAGGAGAGCTGGCTCTGCAGTGGGGGAGGGATGCAGGGA

Sequence ID 679

GNANCNTTTCCTNTCGNAAANCGCGCCTTGTGTTGGTACCCGGGAATTCGCGGCCG
30 CGTCGACAAAAAAAAAAAAAAAAAAAAAAAAAANTNTAGACTCGANCAAGCTT
ATGCANGCNTGCGGCCGCAATTCGAGCTCGGCCGACTTGGCCAATTCGCCCTATAG
NGAGTCGTATTACAATTCACTGGCCGTCGTTTTACAACGTCGNGACTGGGAAAACC
CTGGCGTTACCCAACTTAATCGCCTTGCAGCACATCCCCCTTTCGCCAGCTGGCGT
AATANCGAANAGGCCCGCACC GATCGCCCTTCCCAACAGTTGCGCAGCCTGAATGG
35 CGAANGGAAATTGTAAGCGTTAATATTTTGT TAAAATTCGCGTTAAATTTTGTTA
AATCAGCTCATTTTTTAAACCAATAGGCCGAAATCGGCAAAATCCCTTATAAATCAA
AAGAATAGACCGAGATAGGGTTGAGNGTTGTTCCAGTTTGGAACAANAGTCCACTN

- 204 -

TTAAAGAACGNGGACTCCAACGTCAAAGGGCGAAAAACCGTCTATCAGGGCGATGG
CCCCTACGTGAACCATCNCCCTAATCAAGTTTTTTGGGGTCGAGGNGCCGTAAAG
CACTAAATCGGAACCTTAAAGGGAGCCCCGATTTAAAGCTTGACGGGGAAAGCCC
GGCGAACGTGGCGAAA

5

Sequence ID 682

CACCTGCAGTCCAAGTACATCGGCACGGGCCACGCCGACACCACCAAGTGGGAGTG
GCTGGTGAACCAACACCGCGACTCGTACTGCTCCTACATGGGCCACTTCGACCTTC
TCAACTACTTCGCCATTGCGGAGAATGAGAGCAAAGCGCGAGTCCGCTTCAACTTG
10 ATGGAAAAGATGCTTCAGCCTTGTGGACCGCCAGCCGACAAGCCCCGAGGAAAAC TG
AAACTTTGCTTAACNACCGAATGGNGGGGANCTTTTCCAACGNTTTT

Sequence ID 683

TTGGTTTCATACTGNTGGGGNTTGAATGNTCCCTNCAACACTNATGTTGANACTTA
15 ATCCCTAATGNGGCAATACTGAAAGGTGGGGCCTTTGAGATGTGATTGGATCGTAA
GGCTGTGCCTTCATTCATGGGTTAATGGATTAATGGGTTATCACAGGAATGGGACT
GGTGGCTTTATAAGAAGAGGAAAAGAGAACTGAGCTTGCATGCCC

Sequence ID - 684

nt: 545

20 GTGGAAGNGACATCGTCTTTAAACCCTGCGTGGCAATCCCTGACGCACCGCCGTGA
TGCCCANNGAAGACAGGGCGACCTGGAAGTCCAACACTTCCCTTAAGATCATCCAA
CTATTGGATGATTATCCGAAATGTTTCATTGTGGGAGCAGACAATGTGGGCTCCAA
GCAGATGCAGCAGATCCGCATGTCCCTTCNCGGGAAGGCTGTGGTGCTGATGGGCA
AGAACACCATGATGCGCAAGGCCATCCGAGGGCACCTGGAAAACAACCCAGCTCTG
25 GAGAACTGCTGCCTCATATCCGGGGGAATGTGGGCTTTGTGTTACCAAGGAGGA
CCTCACTGANATCAGGGACATGTTGCTGGCCAATAAGGTGCCAGCTGCTGCCCCGTG
CTGGTGCCATTGCCCCATGTGAAGTCACTGTGCCAGCCCAGAACAACACTGGTCTCGGG
CCCGATAAGACCTCCTTTTTCCAGGCTTTAGGTATCACCACTAAAATCTCCAGGGG
CACCATTGAAATCCTGAGTGATGTGCACTGATCAAGACTGG

30

Sequence ID 685

GGAAAGGGCCATTTTATTGCCTAAAACCACCTGGNTTTTNAGGTAACAGTTCCAAC
ATGTCCTTTTTTTGAATAGCTGTTCTAATTATTATATATTAGCTGATTAATAGGAG
TACTTGATAGGTGGACTGTGTCAGGTAGCCTCAGGCAATCCTACTTCAACAAGCTG
35 TCAGGGAGCCATGCCATGCTTCTTTATGACATAGGTGAATTTGATAGGCTCACTAG
CAGAACATGGGATCACAAGGTGGAACCNNTCCNNTT

- 205 -

Sequence ID 686

GACCCCTTCCTTACACCTTATACAAAAAACTGAAACTGGACCCCTTCCTTACACC
TTATACAAAAATTAACCTCAATTTTATTATGTTGTATTAAATTAAGTTGGGTTTAAT
TAAGATGGATTAAAGACTTAATTATAAGACCTAAAACCATAAAAACCCTAGAAGAA
5 AACCTAGGCCATACCATTTCAGGACACGGGTATGGGCAAAGACTTCATAACTAAAAC
ACCAAAGCAATGGCAACGAAGTCCAATAGACAAATTGGACCTGATTAAACTAAA
GAGCTTCAGCACAGCAGAAGAGACTATCGTCAGAGTGAACAGGCAACCCACAGAAT
GGAAGAAAATTCTTGCAATCTATCCATCTGACAAGGGGCTAATATCCAAAATCTAC
AAAGAACTTAAACAAATTTACAAGGAAAAACACAAACAACCCCATCAAAAAGTGGG
10 CTAAGGATGTGAACAGACACTTCTCAAAAGAAAAACATTTATGCAGCCAACAAACAT
GAAAAAAAGTTTCATCATCACTGCTCATTAGAGACATGCAAATCAAAACCACAATGA
GATCCCATCCACACCAGTTAGAATGGCAATCATTA AAAATGT

Sequence ID - 687

nt: 268

15 TTTATGTGTTTTTGGCTTGGGGGGCGCTGGGCCTAGCCCAGAGTAGTGCTTGCTCCC
CCTGCCTTGTCCCACCAGGGAGGCAGCAGACTCAGGCCCTCCATGGTCTCTTTGT
CATTTTGTGACATGCATTCCTCCTTTTGTTCATCTTGTGTTGGGGGGAGGGGATTAAAC
CAAAGGCCACCCTGACTTTGTTTTTTGTGGACACACAATAAAAGCCCCGTTTATTGT
TAAA

20

Sequence ID - 688

nt: 569

CTTTAGCCAGCCTGATCAGAAAAAAACAAAAGAAGAGGAAAGACGTAGATTACCAA
CATCAAGAATGTGAGTTATGATATCACTACAGACTCTCCAGGTATTAAAAGCATAA
TTAGAGAATGATATGAGCAGCTATATGCAAATAAGTTCAACATTGGACAAATGGAC
25 AAATTTCTTGAAAGATAAATTATGAAATTTCAATTCTGAAAGAACTACATGACCTTA
ATTGTCTTACATCTATTAAATAAGTGGAATTTGTAGTTTAGAACTTTCCACAAA
GAAAACTCTAGGCCCAGATGGCATCAAATAATATTCAGATGAATGAAATGGAGAA
AGGATAGCCTTTTCAACAAATGGTGGTGGAAACAATTGGATTTCCATATGCAAAAAA
ATAGAGATGGACGCAGAGGTGTGTGCTTAGGAGGCTGAGGTGAGAGGATTGTTTGA
30 GGCCAGCCTGGGCAACATAGCAAGACCCCATTTCAAAAACAAAAATAAAGAACTTG
TAGCCTTACCTTGTGCCATATTATGAAAATGTATCATAGGCTTAAATGTGAAACGT
AAAACAAA

Sequence ID 689

35 CGCAGGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCGCAGATAAGTTT
TTTTCTCTTTGAAAGATAGAGATTAATACAACTACTTAAAAAATATAGTCAATAGG
TTACTAAGATATTGCTTAGCGTTAAGTTTTTAACGTAATTTTAATAGCTTAAGATT

- 206 -

TTAAGAGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAAAAGG
TTTCTAAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAATGTAT
TTAAAAGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACTAATAGAAGG
GCAATGCTTTTAGATTAAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTTTAAA
5 AGTTGTAGGTGATTAAAATAATTTGAAGGCGATCTTTTAAAAAGAGATTAAACCGA
AGGTGATTAAAAGACCTTGAAATCCATGACGCAGGGAGAATTGC

Sequence ID 690

CGAAAAGCAAATATAACTTGCCACTAACCAAGATCACCTCTGCAAAAAGAAATGAA
10 AACAACTTTTGGCAGGATTCTGTTTCATCTGACAGAATTCAGAAGCAGGAAAAAA
GCCTTTTAAAAATACCGAGAACATTAAAAATTTCGATTTGAAGAAATCAGCATTTTC
TAACTGAAGTGAGCCAAAAGGAAAATTATGCTGGGGCAAAGTTTAGTGATCCACCT
TCTCCTAGTGTTCTTCCAAAGCCTCCTAGTCACTGGATGGGAAGCACTGTTGAAAA
TTCCAACCAAAACAGGGAGCTGATGGCAGTACACTTAAAAACGCTCCTCAAAGTTC
15 AAACTTAGATTTTCAGATTT

Sequence ID 691

CCGGTCTCTACACAATATATAGAAATCTGGGCATGGTGGTGCCTGGCTGTAGTCTC
AGCTACCTAGTTGGGTGAGGTGGGAGAGTCGCTTGAGTCCTGGAGGTTGAGGCTGT
20 AGTGAGCCAGGGCTGCACCACTGCATTCCAGCCTGGGTAAACAGAGTGAGACCCTGT
CTCAAAAAGAAAAAAAAAAATTGCTAATTTTAACAAATCACAAAACCTGACTCAGGC
AAGTTGTCTGACTCAAAAGCCCTTGAAAAACCATCAAAGACAGTAGAATGTTAACT
GGTCATTTACGTAAAATAGTGTTTCATTAAATTTTTGGTTCATTTAGGATAATCATT
TTAAATGAGACTGTATTTGAGACTGTATACACATACATATACATGTTTACACACAT
25 ATACGTACAATATATGTACATTCTATCTAAAAGATCATACATGTGTGTACATATAT
GTTTTTAAAAGTCAAACCTGACATATTAATGGAAACAGTGCTTACATCTCTGGTAGT
GATTTTCTATTAGCAGCAGCCCTACATATGCTGCGTCTCTGAACAGCATGTCAGTG
CCATGACTGTCTAAACATGCAAATATGACTGACAGACTCTTGAGACAGCTTTCACC
TTG

30

Sequence ID 692

AATTCGNGGCCGCGTCNNCCTANGAGGCACCAGGAAATCCCGCGGGGTGGCCCATG
CAGACCAGGCGCACGTGGCTCATGGGGCANAATTGCCAAGGACAGCTCACGACAGT
GCCACCTTCTCACCATTCCAGCCAAGGAGAGATGTGACGTTGGAACCTGCTCTGGCA
35 CTTCTGTCAAGCCTCCCCGCCCCAATTGCCTTGAGATCTCTGCTCTTTGTCAGAG
ATTTGCAAAGACTCACGTTTTTGTGTTGTTTTCTCATCATTCCATTGTGATACTAAGA
AACTAAGAAGCTTAATGAAAAGAAATAAAATGCCTATGTTGTTGTTCT

- 207 -

Sequence ID 693

CTAGAACCCATGACTCCTAGGTCTTATACTGCAACCACAGTATCAGCAAATAATCT
TTCATAAGGGGATTATTCTCTGATTAAACAGGAAATACAGGAATTTAATTTGTGAAC
ACGCTAGGTAGAAGCAGAAACCCAAATCCAAATCCAAATTTAAACATTTAAAATTC
5 ATTCTATAACTAAGATCTAACAGTCATTTTCTTCCCAGTAAGAAATAACCAAAGCA
TGCTAAAAATCACTGGACTAAATTGGTGTCAAACTGCCACATTGCCAGGCATGGG
GGGGTCATACTTGTAATCCCAGCACTTTGGGAGGCCGAGGTGGGAAAATTGCTTGA
GGCCAGGAGTTGAAACCAGCCTGGGCAACACAGTGAGACCCCATCTCCACAAAAA
AAAAAAATTAAAAAACAAACAAAACATTAGCTGGGCATGGTGGTACACGCCTGTA
10 GTCCCAGCTACTCAGGAGCCTGAAGTGAGAGGATCACTGAAGCCCAGGAGGTAGAG
CTATGACTGTAGTGAGCTATGACTGTGCCACTACACTCCACCTGGGTGACAGGGGA
CTC

Sequence ID 694

15 CGACTTCCATTTGTATTAATGGAATACTAAGTCCCTCTGTGATTTCTGAACCAAGC
TATTCCTAGGCCTGAGTTTTATTTTGTGACACAGAAATAAATTANAAGGCCAAGC
GTGGTGGCATGTGCCTGTAGTCCTAGTTGCTGAGGTAAGAGGATTGCTTGAGCCCA
GGAGTTCAAGGCTGCAGCAAGCTTTGATTGCGCCACTGCACTCCAGCCTTGCGCAC
AGACTAAGACGCTGTCTCAAAAAAAAAACAAAAA

20

Sequence ID 696

GGTTATCAATGAGATTAAGAGACAACCTAGAGTAAAAACAAAAGAAAAGAAAAGAAA
NGAAAACAACAGAAGCTCTATTAAGTACCTCTAACCAATACAAACAGGTAACTGA
TGTTCTCCATTCTGTATATAAAAAATCCAGTGGACACCCACAACACAGGCTTCAGG
25 CTTGTAGGACACTTTCTAGTTCATCTGAGCACTTTTGTCTCAGCAGTTGAGCTGT
ATACTTAGCAACATTTGGTGCCTTCCAAACCCATTTGTGCCTGTAGCACTTACTATT
GAAATACATAATTTAATTAAATATTATATAAAGGAATGGAATACGAGTTGGACAAG
AAAAAGAGTTAAATCTGAAGGTTAGGTAAAAAGAGCAACTTCTTTTCTCTGTTTTG
CAGGTTGGCAAAATCATTTAAAAACAATTGGAAGTATTATATGTTCTGCATTAAGT
30 TGTCATTTTACTTAAAACTAGGCATCAAAGATGATGCATAATAAATTTAGTGTAT
GCAAGAATGACTGCTTGGGACCTCAATATATGAATTCTTAATCCAAGGAAAGTCCT
TGGCCTTACATTTAAAAGTCGGCAAATAAGTGTACGTTTCATT

Sequence ID 697

35 GAACATTTAAAAATAATGCAAATAAGGCTGGGCGTGGGGGCTCACACCTGTAATCC
CAGCACTTTGGGAGGCCGAGGCAGGCAGATCACGAGGTCAGGAGATTGAGACCATC
CTGGCTAACACAGTGAAACCCTGTCTCTACTTAAAAAATAAAAAAATTAGCCAGGC

- 208 -

GTGGTGGTGGGCGCCTGTAGTCCCAGCTACTCAGGAGGCTGAGGCAGGAGAATGGT
GTGAACCCGGGAGGCGGAGCTTGCANTGAGCTGAGATCGTGCCACTGCACTCCAGC
CTGAGCGACAGAGCGAGACTCTGTCT

5 Sequence ID 698

TCATTAGAATCCAAGCTTTGAAAATTTCTGATTAATGCTCATGTATTTCTTTATCT
TTGTTTTTCTTGTGAAGAAAGACTTTCACCACTGTCTGAGTGATGATGCTGTTGA
TAAGGATGATGTCGATGACTACTATATTGCATCTCTCAGGAACAGCTGATGGGAAG
GGAGGGGCTGCTGAGTTCCTTGTCTAGCTAGCAGCACGCTCCTCANAGAGGGGG
10 CCGAGTTACAGACAGCAGCCGCATTCTCATGCAAAATTAGTTTTAACTGCTAGTG
TGGGCATCGGTACCTTTTGCTGGGTGATACCGAAGAATTGTTGAGGATTTAGTAT
GCTCCGTAGAGACAGTTCAGCCAGTCATTTCTGCATTGGAGAGACTTCTCATACTT
TCTTTGAAGACTCATAGAAAGCTGGAT

15 Sequence ID 699

ATTAAGGTTTGTNCCCAACAAGAATAGATGTAATTAGAAAAAANTGNCTTCCTTAC
CTATTGCCTCTGATNTTTACTTGCTTAAATTTTTTTTATTGNAAATCCAGAAAAAG
NGGATTTAGAGAACAACTAACTCCACCTAATCTATGACAGANATGTACAANAN
AGTACCTGTGAAAAATGTGAAAGNATNTGAAAAATGTAACCTTTGGCAGCCTGAGC
20 ATAGTCAACCAGAAAACTATCTGAATTAAAATAATTGGTCCATAGGTACTATTTT
ATTTGGTCCATAAGGATTATTTTTTTCAACTTTTTTTTCAAGTGTATTATTATGTCA
TTTCCCACGTAGGTTACTGATACCTGAAGACTTTTTTNCACCTTTAACCTTNCTCGT
TGAGGAGCTTTGTANTCTAATAAAAGAGAAATATAAGTAAATGTTAGATATATGGG
NGGATAATGGTAACTATGTGCTTAAAGAGGTATAAAAGAAGGGTAGGGAGCAGATA
25 AGACAAAGGAAGGGCTATATTATAANGAAGAATATTCCAAGTAGGGAAGAGAAAAA
GATATGTTATCCATATAATATTTTATGTGCAGTAGAGAACATGTTCTATAGAANAG
ACAGAAGATG

Sequence ID 700

CTTGAGCCCAGGAATTCCAGCCTGGGCAATATAGTAAGACTCCGTCTCTACAAAAG
30 ATACAAAATTAGCCAGATGTGGTGGTGCGTGCCGTGTAGTTCCAGATACTGGAAAG
ACTGAGGCAGGAGGATTGCTTGAGCATGGGAAGTTGAGGCTGCAATGAGCTGTGAT
TACGCCACTACACTCCAGCCTGGGCAACAGAGTAAGATCTTGTCTCAAAAAAAAAA
TTGAATTCAGCTAAAAATAATAAAATTTTAAAAATAATTTTAAAAAGCCCTCAACAG
35 CTTTGTCTTTCTCTCCTTGCCAGCTTCTCTGCAGCCTATAGCCTGCAGGCTGGCTG
CTGCGAGCCAGGACAAGCGGTGGGAAATGCAATCACAGCGTGAAATCTCTGTGTTC
AGAGACACGCAGGAAGCAGGTGAACCATGAAGGGCCAACACATGCCCCCAGTTAGC

- 209 -

AGGGTGTAGAGACCGGGGCAGGGCTTTCTTCTTCCTTCTGGGTTATAAATATCCAT
GTCCTGCCATTTGAAGCTGCAAGTGGCACACATGGATGCTGGACAGGCGCTCGCAC
TTTCTGGGCAGGGCANGGGGCTCAAAGGCAGGACAGCTGGGCAAAGCACCTTGCG
TGGGCCC

5

Sequence ID - 701

nt: 579

CTTTGGAGCTTCTGTCTGTGCTGTGGACCTCAATGCAGATGGCTTCTCAGATCTGC
TCGTGGGAGCACCCATGCAGAGCACCATCAGAGAGGAAGGAAGAGTGTTTGTGTAC
ATCAACTCTGGCTCGGGAGCAGTAATGAATGCAATGGAAACAAACCTCGTTGGAAG
10 TGACAAATATGCTGCAAGATTTGGGGAATCTATAGTTAATCTTGGCGACATTGACA
ATGATGGCTTTGAAGGTAATTAAAATTATCAAATTGGTGCTTGATTTCTGCTTTTA
AAATGGTTTATGGAAGAAAATATGATTAAAGTTTTGTATTGTTTTCTTCCTTCTATAG
AAGATGGAGCCAGAATGGCATGCTAAGTTTTTTCTTTTCTTTAGTGTTATATATGA
CTTCTCCTCAATTGTCACCCATTGATCTTTACCACTGTTAATAATGGATGATATTC
15 AAAATACCTTATTTTCAAGTATTCTAAGGCACCATTGATTAGAAACTGCATTATTAT
TTATGTGTCCCTAAAAGCTACCTATTAAGCTGTTACACCCACCATTTTTCTGTAA
GAAAATCCTGATTTTCAGAA

Sequence ID 702

20 GTNNTCCTCTCGGAACGCGCCTTNTGTAGCCAGGTGCTACCAGACCNAATACACGG
TTGTTCCAGCTTGCGCATTACCGATGGCGTAGATATCCGGATCGGAAGTCTGGCA
GGAATCATTAATGACAATACCCCCACGCGGAGCAACGTCCAGACCACACTGGGTTG
CCAGCTTATCGCGCGGACGGATAACGGTAGAGAAGACGATAAAGTCGACTTCCAGT
TCGCTGCCGTCGGCAAAACGCATGGTTTTACGCGCTTCAACACCTTCTGCACAAT
25 CTCAAGGGTGTTTTTGTGCTGGTGTGAACGCGCACGCCCATACTTTTGATTTTGCAGC
GCAGCTGCTCGCCGCCCATCTGATCAAGCTGTTCTGCCATCAGCATAGGGGCAAAT
TCGATAAC
GTGGGTTTCAATACCTAAGTTTTTTCAGCGCGCCTGCGGCTTCCAGACCTAACAGGC
CGCAATTCGAGCTCGGCCGACTTGGCCAATTGCGCCTATAGTGAGTCGTATTACAA
30 TTCACTGGCCGTCGTTTTACAACGTCGTGACTGGGAAAACCTTGGCGTTACCCAAC
TTAATCGCCTTGCAGCACATCCCCCTTTCGCCAGCTGGCGTAATAGCGAAAGAGGC
CCGCACCCGATCGCCCTTTCCAACAGTTGCGCACCTGAATGGCGAATGGAAATTGT
AAGCGTTAATATTTTGTAAATTCGCGT

35

Sequence ID 703

CTGCGCAGACCAGACTTCGCTCGTACTCGTGCGCCTCGCTTCGCTTTTCTCCGCA
ACCATGTCTGACAAACCCGATATGGCTGAGATCGAGAAATTCGATAAGTCGAAACT

- 210 -

GAAGAAGACAGAGACGCAAGAGAAAAATCCACTGCCTTCCAAAGAAACGATTGAAC
AGGAGAAGCAAGCAGGCGAATCGTAATGAGGCGTGCGCCGCCAATATGCACTGTAC
ATTCCACAAGCATTGCCTTCTTATTTTACTTCTTTTAGCTGTTTAACTTTGTAAGA
TGCAAAGAGGTTGGATCAAGTTTAAATGACTGTGCTGCCCTTTTACATCAAAGAA
5 CTTACTGACAACGAAGGCCGCGCCTGCCTTTCCCATCTGTCTATCTATCTGGCTGGC
AGGGAAGGAAAGAACTTGCATGTTGGTGAAGGAAGAAGTGGGGTGAAGAAGTGGG
GGTGGGACGACAGTGAAATCTAA

Sequence ID 704

10 CTTGTATTCAAGAACTACTGTAATGCATTAGTGGTCTGGCTTCATTTTGTATGATG
CCAGATCCTTAATTTACCCAGCACAATCATTTTCAAGTAGTTTCCCTATGGCTCCTGCA
AAAATGCAAACAGAAACCACCACAGGAACAGCCCCCTTGCTGCCTCCTGTTGCTGAG
GTAGTAGTCGCTAAAGAAAATTGAAGGCTCCTTACAATCTATATTTGAAAACCTAGA
ACTTCTGTAGAAACACACAGATCCCGATCTTAGAAGTTGTACAGGACAATCTGGTA
15 AAAGTACATAATTGTGATTTTATTAACATGAATTAAAATGCCCAACCAGTGCTTCA
GTGTGACAGTATATTTAAAATAAAAAAGAAATTAAAGGTCATATACTGTACTACTT
TCACAAAGATCCACAGTTTTTGCAAAGACTTGTTCATATGTACAATGCTATATATCA
AATGAGAAAAGCTGTAAGCAATTATATACGCAAAGAAATGGCAGTA

20 Sequence ID 705

TTCCAGTCCTTTTCATTTAGTATAAAAGAAATACTGAACAAGCCAGTGGGATGGAAT
TGAAAGAACTAATCATGAGGACTCTGTCTTGACACAGGTCCTCAAAGCTAGCAGAG
ATACGCAGACATTGTGGCATCTGGGTAGAAGAATACTGTATTGTGTGTGCAGTGCA
CAGTGTGTGGTGTGTGCACACTCATTCTTCTGCTCTTGGGCACAGGCAGTGGGTG
25 TAGAGGTAACCAGTAGCTTTGAGAAGCTACATGTAGCTCACCAGTGGTTTTTCTCTA
AGGAATCACAAAGGTAACTACCCAACCACATGCCACGTAATATTTTACGCCATTCA
GAGGAAACTGTTTTTCTCTTTATTTGCTTATATGTTAATATGGTTTTTTAAATTGGTA
ACTTTTATATAGTATGGTAACAGTATGTTAATACACACATACATATGCACACATGC
TTTGGGTCTTCCATAATACTTTTATATTTGTAAATCAATGTTTTTTGGAGCAATCC
30 CAAGTTTAAGGGAAATATTTTTGTAAA

Sequence ID - 706

nt: 496

CAACCCTCTCTCCTCAGCGCTTCTTCTTTCTTGGTTTGATCCTGACTGCTGTCATG
GCGTGCCCTCTGGAGAAGGCCCTGGATGTGATGGTGTCCACCTTCCACAAGTACTC
35 GGGCAAAGAGGGGTGACAAGTTCAAGCTCAACAAGTCAGAACTAAAGGAGCTGCTGA
CCCGGGAGCTGCCCAGCTTCTTGGGGAAAAGGACAGATGAAGCTGCTTTCCAGAAG
CTGATGAGCAACTTGGACAGCAACAGGGACAACGAGGTGGACTTCCAAGAGTACTG

- 211 -

TGTCTTCCTGTCCTGCATCGCCATGATGTGTAACGAATTCTTTGAAGGCTTCCCAG
ATAAGCAGCCCAGGAAGAAATGAAAACCTCCTCTGATGTGGTTGGGGGGTCTGCCAG
CTGGGGCCCTCCCTGTCGCCAGTGGGCACTTTTTTTTTTCCACCCTGGCTCCTTCA
ACACGTGCTTGATGCTGAGCAAAGTTCAATAAAGATTTTGGGAAGTTT

5

Sequence ID - 707

nt: 397

CGGATGTGGTGGCAGGCGCCTCTAGTCCCAGCTACTCGGCAGGCTGAGGTAGGAGA
ATGGCTTGAACCCAGGAGGTGGAGCTGACAGTGAGCCGAGATCGCGCCACTGCACT
CCAGCCTGGGCGGCAGAGCGAGACTCCATCTCAAAAAAAAAAAAAAAAAAATAGA
10 CTTTGAGACCAGCCTGACCAACATAGTGAAACCCGTCACTACTAAAAATACAAAAA
TTACCCGGGCGTGGTGACGGGCGCCTGTAATCCCAGCTACTTGGGAGGCTGAGACA
GGAGAATCACTTGAACCAGGGAGGCGGAGGTTGTAGTGAAGTGAATCGTGCCCCCT
GCACTCCAGCCTGGGTAACAAGAGCGAAACTCCGTCTCAAAAATAAATAAATAAAT
AAAAT

15

Sequence ID - 708

nt: 293

CCAGCTTTTTATGGTGTTTAATCTAATACACTTAAGCTGCAGTCCCAAAATTAGGG
GTCCTTCAGTCTTGGAGACTATAAGGGAGCCTCTGCACCCAGGGAAAATGTTACCC
TTTACAGGGGGGAAGGGTAAACCAGTAGGGAATACAGTACAATCCCAACCCTACTG
20 GGAGGGGCGGGAGGGAGGTGTTGCCGTCACTGTATTAAGTCGATGTTGGGAAACGT
TTTAACATCTGGAGCCTTTGTGGGTGGAAATATGTCTCCAGTTACAACCTCCGCAGT
GGATGTGAAGAAG

Sequence ID 709

25 GGAAGCTACAATGATTTTTGGGAATTACAACAATCAGTCTTCAAATTTTGGACCCAT
GAAGGGAGGAAATTTTGGAGGCAGAAGCTCTGGCCCCCTATGGCGGTGGAGGCCAAT
ACTTTGCAAAACCACGAAACCAAGGTGGCTATGGCGGTTCCAGCAGCAGCAGTAGC
TATGGCAGTGGCAGAAGATTTTAATTAGGAAACAAAGCTTANCAGGAGAGGAGAGC
CAGAGAAGTGACAGGGAAGCTACAGGTTACAACAGATTTGTGAACTCAGCCAAGCA
30 CAGTGGTGGCAGGGCCTAGCTGCTACAAAGAAGACATGTTTTAGACAAATACTCAT
GTGTATGGGCAAAAAACTCGAGGACTGTATTTGTGACTAATTGTATAACAGGTTAT
TTTAGTTTCTGTTCTGTGGAAAGTGTAAGCATTCCAACAAAGGGGTTTTAATGTA
NATT

Sequence ID 710

35 TGGATTCCCGTCGTAACTTAAAGGGAAACTTTTCACAATGTCCGGAGCCCTTGATGT
CCTGCAATGAAGGAGGAGGATGTCCTTAAGTTCCTTGACAGCAGGAACCCACTTAG

- 212 -

GTGGCACCAATCTTGACTTCCAGATGGAACAGTACATCTATAAAAGGAAAAGTGAT
GGCATCTATATCATAAATCTCAAGAGGACCTGGGAGAAGCTTCTGCTGGCAGCTCG
TGCAATTGTTGCCATTGAAAACCTGCTGATGTCAGTGTTATATCCTCCAGGAATA
CTGGCCAGAGGGCTGTGCTGAAGTTTGCTGCTGCCACTGGAGCCACTCCAATTGCT
5 GGCCGCTTCACTCCTGGAACCTTCACTAACCAGATCCAGGCAGCCTTCCGGGAGCC
ACGGCTTCTTGTGGTTACTGACCCAGGGCTGACCACCAGCCTCTCACGGAGGCAT
CTTATGTAAACCTACCTACCATTGCCCTGTGT

Sequence ID - 711 nt: 498
10 GTGGTACATATACACAAAGGAAAACCTATGTAGCCATTAAAAGAAAAGGAACTCCTA
TCATTTGTAACAACATAAATAAATCTGGAGGAGATTAGGCTAAGGTGAAATAAGCC
AGGCACAAAAGACAACCTACCATATGATCTTACTTATACGTGTGTGGAATCTAAAA
AGGTGGAATTTACAGAAGCAGAGAGTAGAATGGTGATTACCAGAGGCTGGGGAGTG
AGGGCAGGAGGTTGGAGAAATGTTGGTCAAAGGATACAAAGTTTCAGTTATACAGG
15 ATGAATAAGTTCAAGAGATCTATTGTACAACGTGGTGGCTATAGTTGATAACAATG
TATTGTGTTCTTGAAAAATGCTGAGAGAGTAGATTTTAAGTGTTCTCACCACAAAA
CATAAGTATGTGAGGTAATGCATGTGTTAATTANCTTAATTTAGACATTTTATAAT
GTATTATACATATTTCAAACCCACGTTGTACATGAGAAAGATACACAATT

20 Sequence ID 713
GCCCAGTCGACCCATGTTCTCCTTTCTACACCAGCATTAGACGCTGTCTTCACAGA
TTTGGAATCCTGGCTGCCATTTTTCAGCTGCCATCCATGACGTTGATCATCCTG
GAGTCTCCAATCAGTTTCTCATCAACACAAATTCAGAACTTGCTTTGATGTATAAT
GATGAATCTGTGTTGGAAAATCATCACCTTGCTGTGGGTTTCAAACCTGCTGCAAGA
25 AGAACACTGTGACATCTTCATGAATCTCACCAAGAAGCAGCGTCAGACACTCAGGA
AGATGGTTATTGACATGGTGTAGCAACTGATATGTCTAAACATATGAGCCTGCTG
GCAGACCTGAAGACAATGGTAGAAACGAAGAAAGTTACAAGTTCAGGCGTTCTTCT
CCTAGACAACTATACCCGATCGCATTACAGGTCCTTCGCAACATGGTCACTGTGCAG
ACCTGAGCAACCCACCAAGTCCTTG

30
Sequence ID 714
CTGTAACAGAGATTCCTTTTTTCAATAATCTTAATTCAAAAGCATTATTAGACTTG
AAAGGGTTTGATAATCTCCCAGTCCTTAGTAAAGATTGAGAGAGGCTGGAGCAGTT
TTCAGTTTTAAATGAGTCTGCAGTTAATATCAAATGTGAGTTTGGGACTGCCTGGC
35 AACATTTATATTTCTTATTCAGAACCTTGATGAGACTATTTTTAAACATACTAGT
CTGCTGATAGAAAGCACTATACATCCTATTGTTTCTTTCTTTCCAAAATCAGCCTT
CTGTCTGTAAACAAAATGTACTTTATAGAGATGGAGGAAAAGGTCTAATACTACAT

- 213 -

AGCCTTAAGTGTTCCTGTCATTGTTCAAGTGTATTTTCTGTAACAGAAACATATTT
GGAATGTTTTTCTTTTCCCCTTATAAATTGTAATTCCTGAAATACTGCTGCTTTAA
AAAGTCCCACCTGTCAGATTATATTATCTAACAATTGAATATTGNAAATATACTTGG
CTTACCTCTCAATAAAAGGGTCTTTTCTATT

5

Sequence ID 717

TCCACCCACCTTGACCTCCCAAAGTGCTGGGATTATAGGCGTGAGCCACCTCGCCC
AGCCCGATACTAGGACTTATGCAGAAAAAACCTTGACATGGAGGAAAGTAAGATCT
AAATAAATACTGTATTCATAGATTAAGAACTCAGCATAATAATATACCATTTCT
10 CCCCAGATTGATGTACAGATTTAACACAATTCCTATCAAGATCCCAGCAAGATTTT
TGTAGATATGTAAAAGATTATTCAAAAATGTAAAAGGAAGGACAAAGGACTAGAAT
AGATAAAACAAAATGGAGAAAGATTTAATAGGAATCACTGTAACCTGATTTTAAGAC
ATACAGAACAATAATAGAACTGCTTGTATTAGTCCATTTTCACGCTGCTGATAAA
GACATACCTGAGATTGGCAATTACAAAGGAAAGANGTTTATTGGCTTACAGTTCCC
15 ATGGCTGGGGAGGCCT

Sequence ID 718

CTCCTCTGGGTTGAAACCCGGGCGCCGCAAGATGCCGGCTTACCACTCTTCTCTC
ATGGATCCTGATACCAAACCTCATCGGAAACATGGCACTGTTGCCTATCAGAAGTCA
20 ATTCAAAGGACCTGCCCCCAGAGAGACAAAAGATACAGATATTGTGGATGAAGCCA
TCTATTACTTCAAGGCCAATGTCTTCTTCAAAAACCTATGAAATTAAGAATGAAGCT
GATAGGACCTTGATATATATAACTCTCTACATTTCTGAATGTCTGAAGAACTGCA
AAAGTGCAATTCCAAAAGCCAAGGTGAGAAAGAAATGTATACGCTGGGAATCACTA
ATTTTCCCATTCTTGAGAGCCTGGTTTTTCCACTTAACGCAATTTATGCCAAACCT
25 GCAAACAAACAGGAAGATGAAGTGATGAGAGCCTATTTACAACAGCTAAGGCAAGA
GACTGGACTGAGACTTTGTGAGAAAAGTTTTCGACCCTCAGAATGATAAACCCAGC
AAGTGGNGGGCTTGCTTTGTGAAGAGACAGTTTATGAACAANAGTCTTTCAGGACC
TGGACAGTGAAGGGAGCCCGGGCAGCCA

30 Sequence ID 719

CGNGGCCGCGTNAACTTTTGATCGTCAGCTGGGGCTGGCAGGCACCTAAATGGGAA
GGGTGATAGCAGTGTGTTGGGGGGAGTTTAGGGAACGGTCCTCTACCGATAGAGGC
AGCANCTCATTGGAATTTCTCCTGAAGTTGTCTTGCCCCCTTGAATCCTGCAGGAA
GGCTGGCAAATGGCCATTTCCCTTCCACTTGAATAGAGACCCATAACTCAAGTATC
35 TGCCCTTAAGACACCACAGGACTGTTCTTCGCGGGGCCCTGCCCTGGATTTGGGAG
AGGCAGTCCANCTCACCCAACTAGGCTCTGCANGGGGACCANGAGGGATGGGTTGT
GTCCACAGGACCAGCCAGACTGATGAGGGATGCGGCAAGCATATTCTCACCACCTT

- 214 -

CTTTCACGTTTACAACANACCAGCNTTCCCTGTGTGGCAGGGGTACATTGGTCAC
CGAGGACCTANAATCATGGAGTGCTCTGGGGATCCGGGCTTGA

Sequence ID 720

5 TCAGTGTGTAATTTTGTTCAGACACTTTCTCTGCATCAATTGGTATGACCATGTGAT
TTTTTTTCTGTAGCCTGTTAATATGGTTAATTTTCAAATATTGAGCTGATTAATTT
TCAAATATTGAGCTCTCCTTGCATCTCTGGAATAAGTACCACTTGGTCGTGGTATA
TATTTCTTTTAATATATTGCTGAATTCTGTTTGATCATGTTTTCTTAAAGACTTTC
GTGTCTGTTTTTCATGATAGATACTGGTCTATAGTTTTGTGTTGAATATCTTGGTTTG
10 ATTTTGATATCAGGATAATGCTACCTTAATAGAATGAATTGGAGCCAAGTATGGTG
GCAAATGCCTATAGTCCTAGCTACTCAGGAGGCTGAGGTGGTGGGGACTGCTTGAC
CCANGAGTTCAAATCTAGCTTGGGCAATGTAGCAAGAC

Sequence ID 721

15 TAGAAGGAATGACTATTCATGTCCAAAGTGAATGGTTTTGTGCAGTGAACAACACA
TGGCGAGGTACTAACTGAGAACTTTTTTCATGCTTTATGCCTACCTCTTGTAGTTG
TTGCAGAGCAAATATAAATTGTAATAAGATAGCTAGGCCTTGCAGAAACAAACAGA
AAAACCTAAAAAAAATGATATAAGAGCTGGAGTCTAGTATTTATATGAATCTGTG
AGAGATAATTTTTTTTGGTCTCACTGCAATGAACCAAAGCGGCTGAGTTTGGTTTT
20 TAATTGTAGCCATGTATTGAAGGCATCTTTTTTGACCAACTCTTGTTGGTTCTGTCT
TGAACCATTTGTTAATCACTGTGCTGTAATTAGTATAGCTAAATCTTTTCCTTCCTT
GCTCCTCCCCCAGCCACCCCGTCTTCCCTTAACATTTTTTTCAGGGGGGGTTGGGA
GTGGTTTCATTTTAATGTGAGTGGATGTTTTTGATAGTTGTAAGGAAAAAATGCATT
TCAGACACATTTACACATGAGCTATTTTCTTACACAGTATGTCTTATTGGTAATA
25 AGAATGTAATTCAT

Sequence ID 722

CNTTCCNTAAGAATACAAAAAATTAGCTGGGCGTGGTGGCAGGCGCCTGTAATCCC
ATCTACTCAGGAAGCTGAGGCTGGAGAATCGCTTGAACCCGGGAGGCGGAGGTTGC
30 AGTGAGCAGAGATCACGCCACTGCAGTCCAGCCTGGGCAACAGTGCGAGACTCTGT
CTCAAAAAAAAATAAATAAATTACCTGGGTGTGGCAGCGCGTGCCTGTAATCCCA
GCTACCCAGGAGGCTGAGGCAAGAGAAGTCTTGAACCCAGGAGGCAGAGGTTGCA
TGGAGCTGAGATGGCGCCACTGCACTCCAGTCTGGTGACAGAGTGAG

35 Sequence ID 724

CTCTCTACTAAAAATACAAAAATTAGCTGGGCACGGNGGTGCATGCCTGTAAACCC
AGCTACCAGGTACTCGGGAGGCTGAGGCAGGAGAAATCGCTTGAACCCAGGGAGTCTGG

- 215 -

AGGTTGCGGCGAGCTGAGATCATGCCACTGCACTGCGGCCTGGAGACAAGAGCAAG
ACTCCGTCTCAAAAAAAAAAAAAAAAAAAAAAAAAAAGACNTCACCTAATTGCAG
NGNGNGGACCTTATTTGGCTNTTAATTCAAACCTATTAAAAATGTGAACN

5 Sequence ID - 726

nt: 260

CGGGGTCTGTACCGGGCTGGCCTGTGCCTATCACCTCTTATGCACACCTCCCACCC
CCTGTATTCCCACCCCTGGACTGGTGGCCCCCTGCCTTGGGGAAGGTCTCCCCATGT
GCCTGCACCAGGAGACAGACAGAGAAGGCAGCAGGCGGCCTTTGTTGCTCAGCAAG
GGGCTCTGCCCTCCCTCCTTCTTCTTGTCTCATAGCCCCGGTGTGCGGTGCAT
10 ACACCCCCACCTCCTGCAATAAAATAGTAGCATCGG

Sequence ID 727

CTGAGTNTAGAAATGATGCCATTAACTGATTGCAAAAACATTACAACCTCAGTAC
TGCAGCTTTCATTCAAATAGGTTATATGTATAAACTGAGTTCAACAATATTGTATT
15 TGAGATGGTAAAGTTAAAGAAATGCAATAATGTAAATAATACTTAAGAAAATAAGA
TCTCAGGAAACTGTATATACTCTGTACTTTTATGCAACTTTATCAGATCATTTTCAG
TATATGCATCAAGGATATAGTGTATATGACATGAACCTTGAGTGCAAAAACCTGTAC
TATGTACCTTTTGTATTATTTTGTCTGTCAACATCTAAATAAAGGTTTTTTTGTGTTGT
TTTTTGTTTTTTTAATTGTTTTGTTTTTAAAGATTGTTTTAATTAATTAATAAATA
20 ATTGTTTTAATTAACAATTGTTTTAATTGTTTTTAAAGTCGCCAGGCTGAGGCAGGT
GAATCACAAGCTTAGGAGTTGGAGGCTAGCCTGCCAACATGGTGAAACCCCGTCTC
TACTAAAAATACAAAAAATTAACCTGGGTGTGGG

Sequence ID 728

CCCATCTGCACCAGTACACAGGCAGGCATTATCATTCTTCACCTACTTTTTTAAATA
GTGGCAACTTGGGATTCATTCTGGTGATTCTGAACCTTGCCTCATAGCTTAAAGTA
TAAAAAAGATTCAAGAGCAGTGAGGTTTGTCTTTCCAGTGAATGGTGGACTGAGT
GGTGCGAGGTGGAGGGCTAACAAGAGGAAAGAACTACATTCTTCAGAATACAGTGA
TGAAAATTCATTTTGAACTCAAATATTTTCATTTTGGATATTCTCCTGTTTTTAT
30 TAAACCAGTGATTACACCTGGCCATCCCTCTAAATGTTCTAGGAAGGCATGTCTAT
TGTGATTTTGATGAAGACAGAATTATTTTTCTCTGTAGAAACACAGATAACCACTTT
ATCAGGGGAAGTTAGTCAAATGAAATGGAAATTGGTAAATGGACAAAAGCTAGCTA
GTAAAAAGGACGACCCAGCAACATGCTTTAACCCCATTTGTATGTTTGTGGAAAGAG
CATAGTTTAAACATCTTGAGAAATTTGGGACATAAAAGTTTTTCATNGGTAGACAGTT
35 CATGGCAGTATATGAATTGACATAATGGAAATAATCTGATTTTATTTTTTACAAC
ACATCCTTTCCCC

- 216 -

Sequence ID - 736

nt: 641

5 GGAATTCCAAGTGCTTGGGGATAATGATACCTCTGACCTTTCTTCCTTTTGGGAAG
TACTTGAGTGTGCAGCTGCATGAGGCCTCAGCAGGAGAGAGATTTTAGGTCCAAGA
AGCTATACCAGTAGGACAAGGCAGGAAAATACTACACTTTCAGGATCAAGCCCCTC
TGACTCTCATTTGGAAACTGGATGTTTGCTAAGCACCTGCTTCTTAAGGATGCCGA
GGGATTTAATGATACTCCCAGAAACCTGGAGAGATTAATGGGGCCTATGGAGAAGT
GCTCTGAACTCAGTGTGGGACTTGAATAAAATTAACCATTGTCATGTTTTTCAGAA
CAACTAAGCTGTTTTATATTTTCATGTGCATGAAAGCCCTAGAACTAAGTTGTGTTA
TTTCCAGAAATGAAATAGATCCCACAGTTAGATGATGTGGCCATTAGGAAGTACCA
10 AATTTATAAAAATCACTGGAGGTCTGTCTGAGCAGTACCTAATAAAATATAGTATA
CTGAAAGTGAACAGATACTTTGTCTCTTTCTTTGGCTGCTTGATCTTTATCTGTGT
CTGCCGTACAGTGCACCCTTAAAGTATTCTACACCAGTGCTTCTCAAACCTGGAAAT
GTGCATGTAAGTCACCCANGGGTCT

15 Sequence ID 739

TGCATGCCCATAGTCCCAGCTATTTGGGAGGCTGAGGCAGGAAAATCGCTTGAACC
CGGGAGCCAGAGGTTGCAGTGAGCCGAGATCGCACTCCAGCTTGGCGACAGAACAA
GACTCTGTCTCAAAAAAAAAAAAAAAAAAGAAATCTTGGGATCCTGAACCCCTTACTC
GAAGGGCTAAGGTAGCATCTCAGCATGTCTTATTTCGAGACTTCGTANAACCAGACC
20 TGCTGTTTGTAGATGTTAATTAATCAAACCTTTCTCTACTCATTTCTGGACCAGTTA
AGGTTTTCTCCTTCTCCGTATGAGTTTTGATTTTCGTCCTCCTTGGTTGGAGATCA
CACTTTGGTCTGCTGCTAAGTTGGATGCCTCCCACTGTCTTTCCCTAAGTCTAGGG
CTTCANACCCCAAGTGTGGGGAGAGGGACTTTCGTTTCCTGCCCCCTCACCACATCAG
ACACAGGCAGGCAAGAATAAGATGGCCAAAAGGCCGATGAACTTCTTGACCTAGCC
25 TGGGACATTACCTGTTACTAGGTGGACTTCACTGCCTGTGAATGGAAGCTGAAGGG
CTGTTTTTTTTGGTTTGTATTTGGACAGGCCAGGCTTANAGAGGGAGAGAACTGGGC
TACTCTTCAGCAGTGATCTTTAAAATGCC

Sequence ID 747

30

CAGAGTGCAAGACGATGACTTGCAAAATGTCGCAGCTGGAACGCAACATAGAGACC
ATCATCAACACCTTCCACCAATACTCTGTGAAGCTGGGGCACCAGACACCCTGAA
CCAGGGGGAATTCAAAGAGCTGGTGCGAAAAGATCTGCAAAATTTTCTCAAGAAGG
AGAATAAGAATGAAAAGGTCATAGAACACATCATGGAGGACCTGGACACAAATGCA
35 GACAAGCAGCTGAGCTTCGAGGAGTTCATCATGCTGATGGCGAGGCTAACCTGGGC
CTCCACGAGAAGATGCACGAGGGTGACGAGGGCCCTGGCCACCACCATAAGCCAG
GCCTCGGGGAGGGCACCCCCTAAGACCACAGTGGCCAAGATCACAGTGGCCACGGC

- 217 -

CACGGCCACAGTCATGGTGGCCACGGCCACAGCCACTAATCAGGAGGCCAGGCCAC
CCTGCCTCTACCCAACCAGGGCCCCGGGGCCTGTTATGTCAAAC TGTCTTGCTGT
GGGGCTAGGGGCTGGGGCCAAATAAAGTCTCTTTCCTC

```

5      Sequence ID - 757                               nt:      583
      GAACCCTGCGGAGGGACTTCAATCACATCAATGTAGAACTCAGCCTTCTTGGAAG
      AAAAAAAAAAGAGGCTCCGGGTTGACAAATGGTGGGGTAACAGAAAGGAACTGGCTAC
      CGTTCGGACTATTTGTAGTCATGTACAGAACATGATCAAGGGTGTTACACTGGGCT
      TCCGTTACAAGATGAGGTCTGTGTATGCTCACTTCCCCATCAACGTTGTTATCCAG
10     GAGAATGGGTCTCTTGTTGAAATCCGAAATTTCTTGGGTGAAAAATACATCCGCAG
      GGTTCGGATGAGACCAGGTGTTGCTTGTTCAGTATCTCAAGCCCAGAAAGATGAAT
      TAATCCTTGAAGGAAATGACATTGAGCTTGTTTCAAATTCAGCGGCTTTGATTGAG
      CAAGCCACAACAGTTAAAAACAAGGATATCAGGAAATTTTTGGATGGTATCTATGT
      CTCTGAAAAAGGAACTGTTGAGCAGGCTGATGAATAAGATCTAAGAGTTACCTGGC
15     TACAGAAAGAAGATGCCAGATGACACTTAAGACCTACTTGTGATATTTAAATGATG
      CAATAAAAGACCTATTGATTTGG

```

```

Sequence ID - 758                               nt:      424
CTTGGCTCCTGTGGAGGCCTGCTGGGAACGGGACTTCTAAAAGGAAGCTATGTCTGG
20 AAGGCTGTGGTCCAAGGCCATTTTTGCTGGCTATAAGCGGGGTCTCCGGAACCAA
GGGAGCACACAGCTCTTCTTAAAATTGAAGGTGTTTACGCCCAGATGAAACAGAA
TTCTATTTGGGCAAGAGATGCGCTTATGTATATAAAGCAAAGAACAACACAGTCAC
TCCTGGCGGCAAACCAAACAAAACCAGAGTCATCTGGGGAAAAGTAACTCGGGCCC
ATGGAAACAGTGGCATGGTTCGTGCCAAATTCCGAAGCAATCTTCCTGCTAAGGCC
25 ATTGGACACAGAATCCGAGTGATGCTGTACCCCTCAAGGATT'TAAACTAACGAAAA
ATCAATAAATAAATGTGGATTTGTGCTCTTGT

```

```

Sequence ID - 764                               nt:      626
GATTTTTTTTTTTTTTTTTTGGAGATGGAGTCTTTCTCTGTCGCCCAGGCTGGAGTGCA
30 GTGGTGAAATCTCGACTCACTGCAACCTCCGTCTCCTGGGTTCAGCAATTCTCCT
GCCTCAGCCTCCTGAGTAGCTGGGATTACAGGCACCAGCCACCACGCCCGGCTAAT
TTTTGTATTTTAGTAGAGACAGGTTTTACCATGTTGGCTAGGCTGATTTTGAAC
TCATGACCCCAAGTGATCTGCCCGCCTCGGCCCTCCCAAAGTGCTGGAATTACAGGT
GTGAGCTACCACTCCAGCCAATGATTACATTTATAAGGTAAAATAACTTGTGCCA
35 ATCTGTACAAGTGAATTCAGATTTAAAATTTTAATTGTAAAAAGATATCCAGGTGA
TATTTCTCCCTGAATAATTTAGTTTTCTTTTCTATTTCTTGATATAAAAGTACTCA
GCATTGAAGTAAATTGCTATCTTCACATTTCTTCCTATTTGAGCTGCTCTAAATAAGT

```

- 218 -

AGTCCTACATATTTTCCCCCAACACAAAAAACCAGAAAAGAATTATTTTATACT
GGATTTTTTTGGTTGTAGCAGGAACCTAAAGGNGCCAATTGTAACATGCATGTTCT
TTTTGGCAAA

5 Sequence ID 766

GTCCATCCTGCAGGCCACAAGCTCTGGATGAGGAACTTGAGGCAAGTCACCAGCCC
CTGATCATTTTCGCCTAAAAGAGCAAGGACTAGAGTTCCTGACCTCCAGGCCAGTCC
CTGATCCCTGACCTAATGTTATCGCGGAATGATGATATATGTATCTACGGGGGCT
GGGGCTGGGCGGGCTCCTGCTTCTGGCAGTGGTCTTCTGTCCGCCTGCCTGTGTT
10 GGCTGCATCGAAGAGTAAAGAGGCTGGAGAGGAGCTGGGCCAGGGCTCCTCAGAG
CAGGAACTCCACTATGCATCTCTGCAGAGGCTGCCAGTGCCCAGCAGTGAGGGACC
TGACCTCAGGGGCAGAGACAAGAGAGGCACCAAGGAGGATCCAAGAGCTGACTATG
CCTGCATTGCTGAGAACAACCCACCTGAGCACCCAGACACCTTCTCAACCCAG
GCGGGTGGACAGGGTCCCCCTGTGGTCCAGCCAGTAAAAACCATGGTCCCCCACT
15 TCTGTGTCTCAGTCCTCTCAGTCATCTCGAGCCTCCGTTCAAATGATCATCATCA
AAACTTATGTGGCTTTTTTGACCTTTGAATAGGGAATTTTTTAAAATTTTTTAAAAA
TT

Sequence ID 768

20 CCAGCGCAGGGGCTTCTGCTGAGGGGGCAGGCGGAGCTTGAGGAAACCGCAGATAA
GTTTTTTTCTCTTTGAAAGATAGAGATTAATACAAC TACTTAAAAAATATAGTCAA
TAGGTTACTAAGATATTGCTTAGCGTTAAGTTTTTTAACGTAATTTTAAATAGCTTAA
GATTTTAAAGAGAAAATATGAAGACTTAGAAGAGTAGCATGAGGAAGGAAAAGATAA
AAGGTTTCTAAACATGACGGAGGTTGAGATGAAGCTTCTTCATGGAGTAAAAAAT
25 GTATTTAAAGAAAATTGAGAGAAAGGACTACAGAGCCCCGAATTAATACCAATAG
AAGGGCAATGCTTTTAGATTAAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTT
TAAAGTTGTAGGTGATTAAAATAATTTGAAGGCGATCTTTTAAAAGAGATTAAA
CCGAAGGTGATTAAAAGACCTTGAAATCCATGACGCAGGGAGAATTGCGTCATTTA
AAGCCTAGTTAACGCATTTACTAAACGCAGACGAAAATGGAAAGATTAAATTGGGAG
30 TGGTAGGATGAAACAATTTGGAGAAGATAGAAGTTT

Sequence ID 773

GAGGAAAGGGGAGTTAATATTTAGTGGACAGAATTTTCACTTTTACAGATGAAAAGA
GTTCTGGAGATAGACGGTGTGTGATAGTTGCACAGCAGTGTGAATGTGCTCATTTGTT
35 ACCGAACCTTAAAAATGTTTAAACATAGTATTATGTGATTTTTTATTTTGGCACTTAAA
AAAAAGAATGAAGTACTGATACATGCTACAACATGGGTGAGCTTTAAATACATTC
TGCTCAGTGAAATAAGCCAGATGCAAAAGATCACATATTATATAATCCACTTATAC

- 219 -

GAGATACCTAGAATAGGCCAAATTCATAGAGACAGAAAGTAGAATAGTGGTTCCCAG
GGGCTGGGGACAAGGGGGCAGTGAGAGATTGAGAGTTATTATTAATGCGTACAGAG
TTTCAGTTTGGGCTGATAAAAAAGTTCTGAAGATGGATGGTGATGATGGTTGTACA
TCAATGTGAGTGTAATTACCGCCACTGAACTGCCCTTAAAAACGTTTAAAAGAGTA
5 AATTTTATGTTGNGTATATTTTACCATAAT

Sequence ID 776

TTTTTTTTTTCATAAGAGGCAAGTACAAGAAAAAGCTTAATTACTTTAACTTCTAAG
TAGTTTGGAACTCTAAATAAATAGGAGTTACCAAATATATGCGCTTCTGTGAATAGT
10 TTTCCCCCACATGTTTATTTATATTTTTTGCATCTCATCAAACCTAACAGATTCTAA
AGTCTCTGGTGATAATGACAATATCTGCTACGGAGAGACTAGCCTGGGGGAAGAGG
ATCTCCCTGAACAAGGATAGCGGAGTTGCTGCAGCTTTCAAATGAAGCTGGACATT
TAGCTGCGGGGGTAGCACCCCTTTGATCAAGGCAGCCCAAAGATGAGTTTCAGGGAT
GGGACTGACAGAAGAGAAAAAGTTCTTCCCAGCCCTTTCTACTTTTTCTCTTTGTTT
15 CTCAGGCTTCTGGCCGTCTTCAGTTTTTCACAAGTTTCACTCTCAACCCTAAACAGT
ACTTCTGTGAAGTACCCTTTGGCCCCCTCGTTTTTCAGCTCCTAAACTCACCTGGAAA
TAGATGTCAATCTAATTTTGGGTCTGACTAGTGCAGTAGGCATTTTTTGGTGA

Sequence ID 782

CTCACACAGAACAATAATGAATGAGTGTGGCTGTGTGCCACTATCACTGTGTCTAC
20 AAAACAGCCAGTGGGCCTGATTTGGCCCTTGGCTGCAGTGCGCCCGTCTCTGTTT
TTGAGGAATAAATCGCATCATTTTCATATGGCTAATGCAATTTTTTTTCCCATCTGG
AAGCAACATCTGATTGGACTCATCTTGATGGTGCTTGTTACAGTCTCTGTAAATG
GGAGAGGGTCCGAGAATAGCTCTTCCTGTTTTTCATCAGGACTGTTTTTGGGGATGG
25 CAAAGAAGTCAGTGTGTCCAGCCTGTGTCTCCTCCTCACCACGTGGCTGATTCCTGAA
TCTGCATGTGCANCACNTGCCGTTGTCTGGGGCATGATCTGTGTGA

Sequence ID - 785

nt: 556

CTTTTCTCTGGGTATAGATTTACCCTAGCACCTATCTCATTATATTGAATTTTCCA
30 GCATATTTAAATAAACTATTAATTAGTCACACTATTTCTTAAAAGTCACACTATCA
ACTAATCGTGACCGCAATTATCTAGGGGTGATAATCTGCTGAGTCTACTCTTTAAA
TACACTGGGACCCAGCATATTGAGTTATATTGGCACAGAACTTCACTCTGGGTAT
AGATTTACCCTAGTACCTTGCCGGCAGGATCCTATTATTATCATGGTTGTACAAGCAA
GGTTCAGGGAAGAGGCTGGCACAGAGAAGGTACCTGGTAACTGTTGTTTGAGGCTG
35 AATTCAGCTCAACTCAGCTCCAGTAGAGATGGTGTCCCCTTCTCTACCGTGTTGAG
ATAGTGTGCAGTCCCTTCCTAAGGGCTGTTACCCACCGCAATAGGACTTGTCAGCT
TCAACTTTTAAATTTCTCTGCTCCCGCTGGGACCCACCCGCTTCAAAAATCATCAT

- 220 -

GGNGGNTTTAGCACCAATTTAGTAAACACAACTGTCTGAAATATTTTGGAT

Sequence ID 796

GAACATTCAAGATAGTGAGAGGAAGAAAAAGATATGGCTGTACGGGACCGAGGTCT
5 CTTCTATTATCGCCTCCTCTTAGTTGGCATTGATGAAGTTAAGCGGATTCTGTGTA
GCCCTAAATCTGACCCTACTCTTGGAATTTTGGAGGATCCGGCAGAAAGACCTGTG
AATAGCTGGGCCTCAGACTTCAACACACTGGTGCCAGTGTATGGCAAAGCCCCTG
GGCAACTATCTCTAAATGCCAGGGGGCAGAGCGTTGTGACCCAGAGCTTCCTAAAA
CTTCATCCTTTGCCGCATCAGGACCCCTTGATTCCCTGAAGAGAACAAGGAGAGGGTA
10 CAAGAACTCCCTGATTCTGGAGCCCTCATGCTAGTCCCCAATCGCCAGCTTACTGC
TGATTATTTTGGAGAAAATTTGGCTTAGCCTTAAAGTTGCTCATGAGCAAGTGTTGC
CTTGCGGGGAGAATTCCATCCTGACACCCTCCAGATGGCTCTTCAAGTAGTGAAC
ATCCAGACCATCGCAATGAGTAGGGCTGGGTCTCGGCCATGGAAAGCATACTCAG
TGCTCANGATGATACTGGCTGTCTGTTCTTAACAGAACTGCTATTGGAGCCTGGAA
15 ACTCAGAATGCAGATCTTTTGTGAACAAAATGAAGCAAGAACCGGAGACNCTGAAT
AGTTTTATTTCTGTATTAAAACTGNGATTGGAACAATTGAAGA

Sequence ID 801

CCACTCCACCTTACTACCAGACAACCTTAGCCAAACCATTTACCCAAATAAAGTAT
20 AGGCGATAGAAATTGAAACCTGGCGCAATAGATATAGTACCGCAAGGGAAAGATGA
AAAATTATAACCAAGCATAATATAGCAAGGACTAACCCCTATACCTTCTGCATAAT
GAATTAAGTAGAAATGAGGATTCTGACCTTGACTTTGATATCAGCAAATTGGAACA
GCAGAGCAAGGTGCAAAACACAGGACATGGAAAACCAAGAGAAAAGTCCATAATAG
ACGAGAAATTCTTCCAACCTCTCTGAAATGGAGGCTTATTTAGAAAACAGAGAAAAA
25 GAAGAGGAACGAAAAGATGATAATGATGATGAGTCAGGTAAAAGTTCCAGAAATGT
GAACAACAAAGATTTTTTTTGATCCAGTTGAAAGTGATGAAGACATAGCAAGTGATC
ATGATGATGAGCTGGGTTCAAACAAGATGATGAAATTGCTGAAGAAGAAGCAGAAG
AAGGAAGCATTTCTGAAATATGAATGAAAAAATTACATCTTTAGAAAAAGAGTTA
TTAGAAAAAAGCCTTGGCAGCCGTCNGGGGGAAGTGACGCACAGAAGAGACCAGAG
30 AATAGCTTCCTGGANGAGACCCTGCACTTTACCCATGCTGCTGGATGG

Sequence ID - 808

nt: 641

CCGGGTTTTAGTATTTAACCAAGAGCCTTTTAAATATTGAAAACCCATAGTTCAGA
AAATGTTAGTATTGCTGCCCTTCTTCACATAAATTTTTTTTAAATTATACTATTA
35 TTTTGCTTAATTTTATATTGGGTTAAAACAACCTTCAAGAAGGTTAACTAGGAAAG
AAGACCTTTTTGTTTTATTTTTACTATTTATATATAGAAGACAAATCAGCATTTGG
TGATAGTTTTACATGACCAGTTATCAAACGGTCATAGTATGAAGTGTGCAGTTGTT

- 221 -

CATTATTAGTAAATTATGTTTGATTTTTTAAACTATTTAGTACTAATAGTTGAGATG
AAAACCTGAAGAAAAATGCCAATGTGACGTTTGTGTATAGCTAGCCTTAAAAAACTT
CCCATGTTTTTAGGTGACTTTTTTCCCCCTCTTAGTACTCTGGAGAAACAATGAAG
ATGGGCCATCTCAATTCCAGATGTAAACAAAAAGTAATTTTTATTTC AACATTTAA
5 TGTAACCTGCTATTATTGNGGATTCTTGNCCTGNGTATTTTCTTTCCCTTATTCAAG
TAATATAGAATAACTTTCCTTAAAATGATTTGATCCAAGATACGTCATTTCTGTAT
TGGCAAAATGCCNCTATTAAAGTGT

Sequence ID - 814 nt: 132
10 GTTAAAGTGATACATTTTTTATACCAAATGTGTTTATTTTTTTTGTGCAAGTAATCCT
TAAAATTGCAATTGTATTAGGTGTTAAAAATAAAGTTTTTAAAAAAATAAAAAAA
AAAAAAAAAAAAAAAAAAAA

Sequence ID 817
15 GACAACCTTAGCCAAACCATTTACCCAAATAAAGTATAGGCGATAGAAATTGAAAC
CTGGCGCAATAGATATAGTACCGTAAGGGAAAGATGAAAAATTATAACCAAGCATA
ATATAGCAAGGACTAACCCTATACCTTCTGCATAATGAATTAAC TAGAAATAACT
TTGCAAGGAGAGCCAAAGCTAAGACCCCCGAAACCAGACGAGCTACCTAAGAACAG
CTAAAGAGGCACACCCGTCTATGTAGCAAAATAGTGGAAGATTTATAGGTAGAGG
20 CGACAAACCTACCGAGCCTGGTGATAGCTGGTTGTCCAAGATAGAATCTTAGTTCA
ACTTTAAATTTGCCACAGAACCCCTCTAAATCCCCTTGTAATTTAACTGTTAGTC
CAAAGAGGAACAGCTCTTTGGACACTAGGAAAAAACCTTG TAGAGAGAGTAAAAAA
TTTAACACCCATAGTAGGCCTAAAAGCAGCCACCAATTAAGAAAGCGTTCAAGCTC
AACACCCACTACCTAAAAAAATCCCAAACATATAACTGAACTCCTCACACCCAATT
25 GGACCAATCTATCACCTATAGAAGACTAATGTTAGTATAAGTAACATGAAAACAT
TCTTCTNCGCATAAGCCTGCGTCAGATTAAACACTGAACTGACAATTAA

Sequence ID - 821 nt: 370
AAAGAGCTCCCAAATGCTATATCTATTTCAGGGGCTCTCAAGAACAATGGAATATCA
30 TCCTGATTTANAAAATTTGGATGAAGATGGATATACTCAATTACACTTCGACTCTC
AAAGCAATACCAGGATAGCTGTTGTTTCANAGAAAGGATCGTGTGCTGCATCTCCT
CCTTGCGCCTCATTGCTGTAATTTTGGGAATCCTATGCTTGGAATACTGGTGAT
AGCTGTGGTCCTGGGTACCATGGCTGGTTTCAAAGCTGTGGAATTCAAAGGATAAA
TTAATGAAGAAAACAAGCGGAGCTGAAGAAGAAAGTACAATATGGTGCTGTCTTCC
35 TAATGAAATAAATTCACTAAATGGACATTAAAAA

Sequence ID 825

- 222 -

AGACTCGAGCAAGCTTATGCATGCATGCGGCCGCAATTCGAGCTCGGCCACTTGGC
CAATTCGCCCTATAGTGAGTCGTATTACAATTCAGTGGCCGTCGTTTTACAACGTC
GTGACTGGGAAAACCTGGCGTTACCCAACCTAATCGCCTTGCAGCACATCCCCCT
TTCGCCAGCTGGCGTAATAGCGAAGAGGCCCGCACCGATCGCCCTTCCCAACAGTT
5 GCGCAGCCTGAATGGCGAATGGAAATTGTAAGCGTTAATATTTTGTAAAATTTCGC
GTTAAATTTTTGTAAATCAGCTCATTTTTTAACCAATAGGCCGAAATCGGCAAAA
TCCCTTATAAATCAAAGAATAGACCGAGATAGGGTTGAGTGTGTTCAGTTTGG
AACAAAGAGTCCACTATTAAAGAACGTGGACTCCAACGTCAAAGGGCGAAAAACCGT
CTATCAGGGCGATGGCCCACTACGTGAACCATCACCTAATCAAGTTTTTTGGGGT
10 CGAGGTGCCGTAAAGCACTAAATCGGAACCTAAAGGGAGCCCCCGATTTAAAGCT
TGACGGGGAAAGCCGCGCAACGTGGCGAGAAAGGAAGGGAAAAAGCCAAANGGAG
CCGGCGCTAGGGCCTGGCAAGTGTACGGGCAACGCTGCGCGTAACCACCCACACCCC
GCCGNGCTTAATGCCCCNTTCAGGGCGCGTNCTGATGCCGNATTTTNTCTTACNCA
TNTGTGCNGGNTT

15

Sequence ID 833

TAAAATAATGGCAAAAAACAAACAAAAACAAGTTCTCTAAACAGAAAGGAAATTA
CTAAAGAAGGAATCTTGAAATAACAGGAAAGAGGAAATACCACAGTAGGCAACATT
ATGGGTAAATAAAACAGACTTTCCTTCTTTAGTTTCCTAAATATGTTTGATGATT
20 AATGCAAAAATTACAATATTTTCTTATGTAGCACTAAAGGTATGTAGAGAAAATAT
TTAAGATAATTGTACTGTAAGCGGGAGATGACAGTGACATAAAGGCAACGTTTTTA
TACTTCACTCAAACCTTATGTATTAATGTAATCCATAAAGCAACCAAAAAAGCTAT
ACTAAGTACATTCAAAAACACAATAGATAAACCAACAAAAATTCTAAAGGATGTAC
AAGTAACCCACTGGAAGCTGCAAAAAATGTAAACAGAACTAAAAACAGAGAATAA
25 ATGAAAAATTAAAAACGAAATGGCAGACTTAGGCCCTAATATACAAATTATCACAT
TAAATATAAATGGTCTAAATACACCAACTGTAAGACAGAGATTAGCAAAGTCGATT
TAAAAACATGACTCAACTACGTGCTGTCTACAAGAACTCACTTCAAATATACCAA
GATAGGAAGGTTGAAAGTAAAACGATGGAAAAAGATGTATCATGTGAACATTAATC
AAAGGAAAGCAGGGGTGGCTATATTAACATCAGGTAAAATAAACTTT

30

Sequence ID - 837

nt: 603

TGAGGNTGGTCATGATGCANAAGCTACTCAAATGCAGTCGGCTTGTCTGGCTCTT
GCCCTCATCCTGGTTCTGGAATCCTCAGTTCAAGGTTATCCTACGCGGAGAGCCAG
GTACCAATGGGTGCGCTGCAATCCAGACAGTAATTCTGCAAACCTGCCTTGAAGAAA
35 AAGGACCAATGTTTCGAACTACTTCCAGGTGAATCCAACAAGATCCCCCGTCTGAGG
ACTGACCTTTTTCCAAAGACGAGAATCCAGGACTTGAATCGTATCTTCCCACCTTTC
TGAGGACTACTCTGGATCAGGCTTCGGCTCCGGCTCCGGCTCTGGATCAGGATCTG

- 223 -

GGAGTGGCTTCCTAACGGAAATGGAACAGGATTACCAACTAGTAGACGAAAGTGAT
GCTTTCCATGACAACCTTAGGTCTCTTGACAGGAATCTGCCCTCAGACAGCCAGGA
CTTGGGTCAACATGGATTAGAAGAGGATTTTATGTTATAAAAGAGGATTTTCCCAC
CTTGACACCAGGCAATGTAGTTAGCATATTTTATGTACCATGGNTATATGATTAAT
5 CTTGGGACAAAGAATTTTATAGAAATTTTAAACATCTGAAAA

Sequence ID - 839

nt: 71

ATTTATCTAATATTTGGTTTAATAAAATGTGAATAATGAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAA

10

Sequence 849

nt: 622

15

TTTTTTTTTATTTTTTTGAGAATGGAGTCTTGCTCTGCCGTCCAGGCTAGAGTTCAG
TGGTGCGATCTCAGCTCACTGCCACCTCACCTCCTAGGTTCCAGAGATTCTTGTC
TTCAGCCTCCTCAGTAGTTGAGAATACAGGAACACGCCACCACGCCTAGCTAATTT
20 TTGTATTTTATAGTAGAGATGGGGTTTCACCATGTTGGCCAGGCTGGTCTCAAATCTC
CTGGCCTAAGTGACCCACCTGCCTCAGCCTCCCAAAGTGCTGGGATTATAGGCGTG
AGTCATTGTCCCCAGCCGGATGTTTTTCATCTTGATTTGCCTTAGTTTCTAAATCTC
ATCCTCTCCATTTTCTCCTGTTAGTAGTCACAGAGAACCAAATTCTGTCAAGTTAT
GAAACTAAAGTCTCTCTTCCACAAGTCTTCCTGTGTTCTGCCTCAAGTGAACCTGA
25 AAGAACATCAGTTTGTGGGAAGGTTGAAGACCGAATGATCTGCTGGGAAATCACTG
AGGCATTGCCATTCTCTTGAGGAATTTCAATTTTCATCGAAGTTTCGGTTTATATCC
CTTTCTTGGTGAGTACTATTGCTGTATGTAAATTAAATGAGTCGTCATCCTTCTT
NTGAGC

25

Sequence ID - 860

nt: 501

30

GTGAAATCACTTTCATGGATTATTAATGGATTTAAGAGGGCATCAATCAGCTCAAC
TCAAGATTTTATAATCATTTTTTAGTATTTAGATTGTGCCTCAAAGTTGTAGTACCT
CACAATACCTCCACTGGTTTCCTGTTGTAAAAACCTTCAGTGAGTTTGACCATTGT
GCTCTTGGCTCTTGGGCTGGAGTACCGTGGTGAGGGAGTAAACACTAGAAGTCTTT
35 AGTACAAAACCTGCTCTAGGGACACCTGGTGATTCCTACACAAGTGATGTTTATATT
TCTCATAAAGAGTCTTCCCTATCCCAAGGTCTTCATGATGCCAGTAGCCATATATG
ATAAATTATGTTCAGTGATAACTTAGTTATCAGAAATCAGCTCAGTGGTCTTCCCC
GCCATGATTCACATTTGATGAGTTTTTAAAAATCAAAGTGATTTTGAAAATCTCTA
ATGGCTCAGAAAATAAAAACATCCAGTTTGTGGATGACTATATTTAGATTTCT

35

Sequence ID 864

TTGTGTTTTTAGGACTCCTTATCTAAATTAAGGCAGAGAAGTTACAGTATTTATAT

- 224 -

CTGCATTAAATCTCAATTCCAGAAAAACCTTTTGAAAAATTATTTAATCCTCTGGA
AACTATTGATATGATACAGGAGAAATTTTCAGAAGTTTATTGAATAATTTAATATC
ATTTAATAGGACACTCTGGCTTGTATATAAGCAGATACGTTACTCAGACTTCTTGG
CTGTACTCTAAAATAATATATGTACTAGTCTCCTAAATATTACTAGCTCACCTTTC
5 AAAATGCATACTAATATTTCAATGTCTTTCTTCAATTTGAAAAGCTCTTGAATATC
TACTTGTGATAGCCCTAAGAGCTGAGATAATTATTTCCAGGAGGTTGAATCCCTGA
TTCTTAACTGTTTCAGCAATGCATAAGCAAGAGAGAATATGACATAAGAGGACCATT
TCTACATTAGCCATTTTTTTTTTTCACAAGATACCTATGTGAATACAGGGCACCTGGGA
GGGTAAGTGGAGGACTATTTCTAACTATATTTATAAGCACATACTGATATTGGTGA
10 ATCAAAACCTACAGCAGTGCTTCTCAGATGGGAAGGGAGACAATGTGTAAGGAGAT
CAGGAATTCATTAG

Sequence ID - 865

nt: 122

CCANAATCCACTCTCCAGTCTCCCTCCCCTGACTCCCTCTGCTGTCTCCTCCCCTCTC
15 ACGAGAATAAAGTGTCAAGCAAGAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAA

Sequence ID 867

TTTTTTTTTTTTTTTTTTTTTTCAGAGTCACAGATATTGTATAGCTGAGGTAAGCATTT
20 TACAACTTTTCAGACACAAGTAAGTACATAAAATATTATTTTACAACCAACAATNTT
TAATATTTCCACATTGAANAATAGATGTGATAATTAAATCTTTTATAAGGTTTTAA
AAAGACATGAAACATAAACCTAATTATACATAAAAGAAAAGAATTTTAAACAAGAG
CTTATTGNGATGACATTACTCATAACTTTTACCTTTAAAACCTTTTCTTGGGTAGC
TATTCAAAAGTAAAGACCACAAGTTTTGTGTTGCCCANATTTCTTATGTTTNGTATAT
25 TTAAGCTCTTTATTTATTGAACAGATGNGTCATTAATTCATTNGGAGCATTACTAT
TATCAGTAAAATTTGATTTTTTTTTTCCCCTCAGTCATAGGTAAATCAGCTCCACCT
GGAATTTCTAAGGACCCAGTTTTAGTCAATATTTTCAAGTAATCATGACCTCAGAA
ATAGTCTTAATTAAGATAACAAATATTAGCCATCAAAATGGAACCAAGACAAGATT
CTAATGTTTGTAACAGTCAATCCATATTTATGAATATTAGCATATATTGGNGAAT
30 AGTTAAGGCAAAGGGTCTAGCAG

Sequence ID - 869

nt: 667

TTGTGTTTTTAGGACTCCTTATCTAAATTAAGGCAGAGAAGTTACAGTATTTATAT
CTGCATTAAATCTCAATTCAGAAAAACCTTTTGAAAAATTATTTAATCCTCTGGA
35 AACTATTGATATGATACAGGAGAAATTTTCAGAAGTTTATTGAATAATTTAATATC
ATTTAATAGGACACTCTGGCTTGTATATAAGCAGATACGTTACTCAGACTTCTTGG
CTGTACTCTAAAATAATATATGTACTAGTCTCCTAAATATTACTAGCTCACCTTTC

- 225 -

AAAATGCATACTAATATTTCAATGTCTTTCTTCAATTTGAAAAGCTCTTGAATATC
TACTTGTGATAGCCCTAAGAGCTGAGATAATTATTTCCAGGAGGTTGAATCCCTGA
TTCTTAACTGTTTCAGCAATGCATAAGCAAGAGAGAATATGACATAAGAGGACCATT
TCTACATTAGCCATTTTTCACAAGATACCTATGTGAATACAGGGCACCTGGGA
5 NGGTAAGTGGAGGACTATTTCTAACTATATTTATAAGCACATACTGATATTGNTGA
ATCAAAACCTACAGCAGTGCTTCTCAGATGGGAAGGGAGACAATGTGTAAGGAGAT
CAGGAATTCATTAGTCACCTTTCAGATGGTTTAATGCATACAGCTGTACCG

Sequence ID 870

10 GGAGTTTGAGCAGATCCTTCAGGAGCGGAATGAACTCAAAGCCAAAGTGTTCCCTGC
TCAAGGAGGAACTGGCCTACTTCCAGCGGGAGCTGCTCACAGACCACCGGGTCCCC
GGCCTTCTGCTCGAGGCCATGAAGGTGGCTGTCCGGAAGCAGCGGAAGAAGATCAA
GGCCAAGATGTTAGGGACACCAGAGGAAGCAGAGAGCAGTGAGGATGAGGCTGGCC
CATGGATCCTGCTCTCCGATGACAAGGGAGACCATCCCCACCCCCGGAGTCCAAA
15 ATACAGAGTTTCTTTGGCCTATGGTATCGGGGTAAAGCTGAATCCTCTGAGGATGA
GACCAGCAGCCCTGCACCCAGCAAGCTAGGGGGAGAAGAGGAGGCCCAACCACAGT
CTCCAGCTCCTGATCCGCCCTGTTCTGCCCTCCACGAACACCTTTGTCTGGGGGCC
TCAGCCGCCCCAGAGGCCTGACTTAGGGGTCTGGCTGTGGAAGGATGTGTGGCCTC
AAATGAGGACAGGGCTCCCGCCTTCACAGCCCTCGCCAGGGGTCTGCCCAATCCT
20 GGCCTGCATCAGGCAAGGACGGGGTCTCAGC

Sequence ID - 871

nt: 642

GCAAGTCTTCAGTATGTACATTTATCCCCTAGAAGAAGAAAAATTAGTTGTGCATG
AAAAAGAAACATTAAGTCAAAGCTAAATGCTCACACTCTAAATCAGTGCTCTCCA
25 AAGTACAGCAGGCGGGAAAGAAATGGTAGATTTTTTTCTTCCAATTACTTTAAC
TTATTCTTTTAAATGGACACTTCATACATAAATATATTACAAATATATTAATATAT
ACATAATGTATAAGCATACATATTGAATGTGCAGTCAAAAAATGTACTAATGGAAT
GCTCTACCAAAACAAGTTCACGTTTCTGTAAAATGGGAATAATATTTTTTAAAAG
GCATACAGTCTGAACATTTTTAGATTATTCATAAAATCTATTCAGAAAGTTAAACT
30 AAAAAATTTAACGTATGCCTATAACAAATTTTGTACTTAATGTAATTGNTTTTCAT
CCTGAGATCTAATATCCTCGTTTTTAAAGTAGAGCCACTTGTTTGCTACAGTTTAGT
CAAAACGTTAACATTAGATGGGTAAAGTAATATGAAATCTTTCTACTACTCCAAA
TAGAAAACAGAACATTAAAAAGATAAAAATTCAAACATACTTACCAGTAGATTTTC
AACTGNGCAAAGCTCATTTGCATGGG

35

Sequence ID 873

GTTTTCCACCGTGAAGAGAACATTTCTCTGGGAATGACAAAGCCCTCAGGAACNG

- 226 -

CTTTTATTTCTATTGGAAGATGCCCATCATACTTCTGGCAGGATAAAATGATAAAT
TTATTTATTCAACAGATGATACTCAATTCCCTGCTGTTTTACTAAAGGTTCTTTAC
GTTTTATAGAAGCTAAATTTACTGTCATAGAAATTGCAATTGTAGATGTTACTGTA
ATCTAGTCAGAATATCCTTATCCTTCTAAAATAAAACTAGTTAAAATTATTAACAT
5 ACGTACTGATATTAATTTTTAAGTTTAATGCTGCCACGTGCTTCTGCTAAGAACAT
TTATCACTACAAGTGGCAGAAAATTCCAACTCATCAAAACCAAACCTGTTGCTTCT
TCCCTGCTTTTTTCAGAAAATGAGAAAGGATGACTTTATTCCAACATATTCTAAAAG
TATTCCAAGAACAACCTACCTTTATTCTAAATTCGTTATTTTCACAAAATAAAGGCTG
CAGATTGAAAGATAAAGGATTGCTATTAAAGAACAAAAGAAAACAAAACCGAGAGA
10 GAAGGAGAGCTAGGGAAATCCCTGCANAANAACCGAATANGGTCCCTCTATTCTGG
GCCGGGGCCTGAAACTATGAAACAGGCCAACACAGAATCTTGCCA

Sequence ID 875

CCTCTGACTCGCTCAGCTCACCCACGCTGCTGGCCCTGTGAGGGGGCAGGGAAGGG
15 GAGGCAGCCGGCACCACAAAGTGCCACTGCCCGAGCTGGTGCATTACAGAGAGGAG
AAACACATCTTCCCTAGAGGGTTCCCTGTANACCTAGGGAGGACCTTATCTGTGCGT
GAAACACACCAGGCTGTGGGCCTCAAGGACTTGAAAGCATCCATGTGTGGACTCAA
GTCCTTACCTCTTCCGGAGATGTAGCAAAACGCATGGAGTGTGTATTGTTCCAGT
GACACTTCANAGAGCTGGTAGTTAGTAGCATGTTGAGCCAGGCCTGGGTCTGTGTC
20 TCTTTTCTCTTTCTCCTTAGTCTTCTCATAGCATTAACTAATCTATTGGGTTTCAAT
ATTGGAATTAACCTGGTGCTGGATATTTTCAAATTGTATCTAGTGCAGCTGATTTT
AACATAACTACTGTGTTCTTGGCAATAGTGTGTTCTGATTAGAAATGACCAATAT
TATACTAAGAAAAGATACGACTTTATTTTCTGGTAGATAGAAATAAATAGCTATAT
CCATGTACTGNAGTTTTTCTTCAACATCAATGGTCATTGNAATGTTACTGATCATG
25 CATTGGTGAGGNGGTCTGAATGTTCTGACATTAACAATTTTCCAT

Sequence ID - 876

nt: 115

AAACTTTTGTGGCAACAGTGCACATAATTTGGATAATGTTTGTTCCCAATAAATTAA
GAGCCAAATTGTAAA
30 AAA

Sequence ID - 878

nt: 634

GCCAGGCTTTGTGAATTACAGGACATTTGAGACAATCGTGAAACAGCAAATCAAGG
CACTGGAAGAGCCGGCTGTGGATATGCTACACACCGTGACGGATATGGTCCGGCTT
35 GCTTTCACAGATGTTTCGATAAAAAATTTGAAGAGTTTTTTAACCTCCACAGAAC
CGCCAAGTCCAAAATTGAAGACATTAGAGCAGAACAAAGAGAGAGAAGGTGAGAAGC
TGATCCGCCTCCACTTCCAGATGGAACAGATTGTCTACTGCCAGGACCAGGTATAC

- 227 -

AGGGGTGCATTGCAGAAGGTCAGAGAGAAGGAGCTGGAAGAAGAAAAGAAGAAGAA
ATCCTGGGATTTTGGGGCTTTCCAATCCAGCTCGGCAACAGACTCTTCCATGGAGG
AGATCTTTCAGCACCTGATGGCCTATCACCAGGAGGCCAGCAAGCGCATCTCCAGC
CACATCCCTTTGATCATCCAGTTCTTCATGCTCCAGACGTACGGCCAGCAGCTTCA
5 AAAGGCCATGCTGCAGCTCCTGCAGGGACAAGGACACCTACAGCTGGCTCCTGAAG
GAGCGGAGCGACACCAGCGACAAGCGGAAGTTNCTGAAGGAGCGGCTTGCACGGCT
GACGCAGGCTCGGCGCCG

Sequence ID 879

10 GTTGCCGGGTCCTGTGATAACTCTGTTTAACATTTTGAGGAACTGTTGAATGGTTT
TTCACAGCAGCTGCCTCATTFTTTTATTTCCCATCAGCAGTACTTCTTGTTCTAATA
CCTCCACGTTCTCGCCAACACTTGTGTGTCTGTAATTTCTGTTGTAGCCATCCC
AGTGGGGATGAAGTAGTATCTTACTGTGGTTTTTCAGTTGCGTTTCCCTGATAATTA
ATGATGGTGAACATCTTTTCATGTTCTTGTGTTGGCCATTTGTATGTCTTCTTGGGAA
15 AAAAAAATGTCTGTTCAAATCCTTTACAAAGTATTTATTTTTTATGTCAACAATA
TAACCACTCAGTACACTGCTTTTTANACAATGATCTTTTAAAGGTTTGTTTACAAC
ATTTAGCACTTGAAATTTTAAGGTTATGCCCTCAAAAAAATTGCTGAGGGAGCTAA
GCTATGAAGATGCAAAGGCATAANAATTATACAATGGACTTTGGGGGAATCCAGGG
AAAGGGTGGGAGGGGGGTGANGGA

20

Sequence ID 881

TCGACTCTGATTTTTTTTTTCTCCTTCCTCGCAGCCGCGCCAGGGAGCTCGCGGNGC
GCGGCCCCCTGTCCTCCGGCCCCGAGATGAATCCTGCGGCAGAAGCCGAGTTCAACAT
CCTCCTGGCCACCGACTCCTACAAGGTTACTCACTATAAACAATATCCACCCAACA
25 CAAGCAAAGTTTATTCCTACTTTGAATGCCGTGAAAAGAAGACAGAAAATCCAAA
TTAAGGAAGGTGAAATATGAGGAAACAGTATTTTATGGGTTGCAGTACATTCTTAA
TAAGTACTTAAAAGGTAAAGTAGTAACCAAAGAGAAAATCCAGGAAGCCAAAGATG
TCTACAAAGAACATTTCCAAGATGATGTCTTTAATGAAAAGGGATGGAACACTACATT
CTTGAGAAGTATGATGGGCATCTTCCAATANAAATAAAAGCTGTTCTGAGGGCTT
30 TGTCATTCCCAGAGGAAATGTTCTCTTCACGGTGGAAAACACAGATCCAGAGTGTT
ACTGGCTTACAAATTGGATTGAGACTATTCTTGTTTCAGTCCTGGTATCCAATCACA
GTGGCCACAAATT

Sequence ID 883

35 TCATTTACATTAATACTCAAAACTGCTCGATTAAGCAGGTGCTGTTCTTATCGCCA
TTTTGCATATGATGAGAAAGGGTAAGGTCACCCAGCTAGTATTTGGCTCACAGCAG
GCCTTAAGACTTGGTTTGTGTGACTCATCAGTCCACGCTCCTAAAACCACTAAGTT

- 228 -

5 GTTCTACCCTTTAATGTTGAATTAACATTGGATAGTGTTCAGTTTANATGGGTGG
GTGAGGGCCCAAGGACCTTTCAAACCTCAGATCTCTTATTTAATAACCTGGTCCCAG
ATCCATTCTCTGTCTGAAGAGGAAGTCATCCTTCAGTGGCTATTTCATTGTGGGGTT
AAGAGCGCAGACTATGAATTCAGTCTTTTTGGGTCCCAGTTTGCCAGACCTTGAGT
GAGTGCCCCGAGTTTACTTACTTGTAAAGGTAGGTGGAGGTAATATAATTAAATAA
ACTTAAAAAACTAATTAAAAACAAAACAAATGAACTAAGGTCTTAGGATATCTGGC
GTCTATTTTGGCCCAAATCACATAATGTCTATTGTTGTGTGTTGGACTATAGGATT
GTCCTTTAACAGGGAAGGGTTTATTTCTGTAATCAAGTCTGTCAATATTATGACCA
TGTTGATAATAGCTACCTTTAATTGAGGGCTTCCATGTGCCAA

10

Sequence ID 885

15 TCAGTGGAAAAGGGCAGGTTGAATCAAGGTGAATCAATCTGAAATTGAGCACACCT
GCCTGCCATCGCTGTTCTTCAACTGAGTGCTGCACATCATGGGCTCTGTCTGTGA
GAGAAAAATCCCGGTGCTTGGTGTCTTGCATGACATGGAGTTTTCATGTAGATC
AATTTAAAATGTACCTCTTGTTTACATAATTTGCATAATTTTAAAAGATAATGTTG
CCAACTTTGGAAATGTTAATGTTCANACTGAAAATCTCCACTACATGTAACTTTC
TTCTCTGGATCAGTGGCATGGCTTATAATCCCAGCCAGTGGTTTGAACGTTCCTCA
GTGTCAACTGCCATGTGCTCTGCTTCAAGGGGGAAGTACCTTTTGTGAATTTTTT
GTACATAAGTATTTGTTACAAATATTTTAGCAAATGCTTCTATTTCTCTTGCTTG
20 TGCATATCTTGGCTGGCGTTACAGAAAAATAGTGTAACATTATTTCTTACCGGG
GAATGAGGGTTTT

Sequence ID 887

25 AGCACCTGGCACAGAGTAGTAGCTAACACAGATGTTAATTTTGCTGCGTCAAATGT
TTTCACTTTGAATCTCTCTTGAGTATTGTTCTCTTATTGATTACATGATGACATC
CTGTTTTCTCTCCCTGACCTTTACTGTTTGTGTTAGAAAAAAAAAAAAAAAAAAAA
AAAAAA

Sequence ID 889

30 CAGAGAGCTTGTTCCCTCCCTCCCTGTGCATGCAAACAAGAGGGCATGGGAGCACA
CAGAGAGATGGCAGCCACCTACAAGCCAAGAGGAGAAGCCTCACAAATCAAACCTCTC
GCTGCTGGCGAGAGTCTTGGACTCTGTCTTGGACTTCCAGCCTCCAGACTGTGAGA
AACAAATTTCTGTTGTTTCAGCTTCTCAGTCTCTGGTGTTTTGTATATGTCAGCCTG
AGAACACAGCTGTACNATTATNAGGGAAACAGAAAACACTGATACTTAACAATGCT
35 AATGCAATTATTTATTTGCTTTTTCAGTCTCTACAAAACGTTCTAAACACTAATCT
AAATATTAACAGTAAAATATTTGCATAACTAATGGAAACTAAGAAATCATATGACC
AATATTTCACTTATTGGTAATCTTACTCTACTGATTTCCCCCCAGACTGTGATTTT

- 229 -

TGAACTTCCTTGCCTTTCTCCTGTCTTTCTGNGTTTATTTCATGGAATTCCAGTTAT
CTGGGCTTGAAATTGCAGGCTCTCCTAACTTAAGCAAAATCTGACAGATCAGCAAA
ATGAGATAAATGTTTCTTTTCTTTCTGACTGCATTAAATCAGATACAACTCAGC
ATTAAAAAGCTATCTTTGNAAAATGNTGGTACTAATAAATTAGTCTTA

5

Sequence ID 890

CCAGTTCCACATTCAGTGAAGTCATGAACTTGAAATTGGCCATGATCAAAAAGTAT
TTAAATCACAGAAGTTGCAAATGCCACAAATCAAGGTCTTTTTCTCTTGGAGAACC
TGTTAAACATTTACCAACTCACGACCGCCATGCACCCAATACTGCAATAGGTCTAT
10 AGATGCAGATACTGTCTCCATGAATCTTATAGGCTAGAAAGGAAATAGATAAGTAG
TCCTACCAGAAGAACATGATGAAGGCATTTGTGGTAAACAGAATGATGGCCCCCA
AAGATGTCCACATCCTAATCCCTGAAGCCTATGAATATACTACTTTACTTGGCAAA
AGGGACTTTGCCACAGGTTTTTAATTAAGGACCTTGAAATAGAGAGATTATCCTGG
ATAATCCAGATGGCCCCAGTGTAATCCCAAGGGTCCTCACAAAGGGTAGGAAGGAG
15 AGCCAGAGTCAGAGAAGGAGACGTAGCAATGGAGGCAGAGGTANAGAGAGATCTG
CAGATGCTGCTGTGTTGGCTTTGAAAATGAGGAATGCAGGTGACCTCAANGNGCTA
GATGATGCAAGGAAACAAATAATCTCCTATGAACCCTAGGATGGGCATTATTATGA
GTCCTATTTTATAACAAGGAACCTGACNTCCAGAAAGATAAATGC

20

Sequence ID - 891

nt: 626

GGCAGAGGTTGCAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGT
AGTGAGACTCCATCTCAAAAAAAAAAAAAAAAAAGACAAGAGTNTCCACTCTAAACA
CTTNTATTCAACATAGTCCTGAAAGTCGTAGCCACAGCAATTTAACAAGATAAAGC
AATAAAATGTATTCAAATAGAAAAAGAGGAAGTCAAATTATCTTCACTGGNGATAT
25 AATTCTCTACCTGGGAACTTCACCGAAAAAGATTTACCAAAAGATTTCTAAGCC
TAAATAATGACTTCAGCAAAGTCTCACCATACAAAATCAACATACACAAATGAGTA
GCATTTCTGTGCACCAATAATATTCAAGCTGAGAAAAAAGAACATGGTTCTATTT
ACAATAGCTACAAACAAAAAATATGTACCTAGTAATACATTAAATCAAGGNGGTA
AAATATCTNTACAACAAGAACTACAAAACCTGCTGAAAAAATAGAGACACGCAAA
30 TAAGTAAAAAGGCACTCCATGCTCATGAATTTAAAGAATCAATATAATTAAATGT
CCGNGCTGCCTAAAGCAACTTACAGATTAAAGGCTATTTCTCTCAAACTATAAATG
CACCTTTTTTA

35

Sequence ID - 893

nt: 585

GTCATTGCTGGGTGGCGCCAGCCCTCAGACTTGCCTCTTTGCAGTAGGAAGAAGGC
CTCCCCACATACCTTTCCACACTCATCACCTTAAGCCAGACTCGGTGTCCAGTGAA
TATGACCATCTCTTGCCCATTTTCTAATGAGTGTTTTTCATTAATGAGTTATAAGAA

- 230 -

TGTGGTGGGTAAATCTATGGGCTTTGAACTAGTGAATCAACTTGGTTTCAGAATCT
GGCACTGCTACTTACTAGTGAATTTAAGCAAGTTATTTACCTTTCAGAGTGTCAG
TTCCCTCATGCATACAAGGAAGATAAAAAATAATGTTNTACNAAAGTATTGGAGTAA
TTAATACATGGAGAACTACATGTAAAGCGTTTAGCATGATGTCTGACATATTAAGC
5 ATCCAATATTAGTNGCTTGCAGAATTATTAGTAAAAGAGATTGCTTCTGAAAGCCA
TTCCAATTCTTAAATTTTATAATGCCACATTTGAGGTCACCTGAAGTCGTGTATAA
CATGTGTACATTTTTCGATTATTTTTTTCAATTCCCANATTAAAGGCATAGAGAT
ATCCTAGCNANGGACTCCAAGTGTG

10 Sequence ID - 895 nt: 560
GTAATTGCAGCCTGGGCAACGGAGTGAGAGACTGTCTCAGGAAAAAAAAAAGAAAA
AAAAC TACTGAGGTAGTTGAATATATCCTCCATTCCCCATTTGTGGATTAGTTAGT
AAATGGGGCATCTTAGGGTTTAAATATGTCCAGGGTCACTGAGGATCAGATCCTAG
GGTTCCTTTGACTCAAGGCTTTTGTCTCAGCAAAACGTCACCTTCCAGCAGGAAGG
15 CTTTCTCAGGCAAGTAGCAGGGTGGCTACTATGTATCGCTTCTTTATTTTTTCTTT
TTTAAAATAATGCAGGCACCGTGCGCATAATTTAAAAAATCAGTGCTAAAACCCTT
AAAAAAAAAAAGCTGTTCTCATCTCCTGTCTTTCTTTTTTTTTTCTTTTTATTTTT
TTCTTTTATTATTATTATACTTTAAGTTTTAGGGTACATGTGCACAACGTGCAGGT
TTGTTACATATGTATACATGTGCCATGTNGGTGAGCTGCACCCATTAACTCGTCAT
20 TTAGCATTAGGTATATCTCCTAATGCTATCCCTCCCCCTCCCCCTTTTTTTTTTT

Sequence ID 896
GGGAATGTCTTAGGCACTGGGACTGTAAGTGCAAAGACCCTGTGGCACAAGGGAAT
GTTAATTATCTACCTTTCANAACTGGAANAAGGCCTAGCCTAGAGCATTGAAAAC
25 AATAAGGGAAAGGAGGAGTAAGGCTGGANAGATAGGAATGGTTTAAAGTCTTTGTT
AAAAATTTTTTTTAAAAAATCTTTATCACAAGAAGAGGATTGGCNTGATCAAATTT
GACTTTTTAAAAANATTACTTGGGTTGGGCATGATCAAATACTACTTAGGGAGATTA
GTTTANATGATAATGGCATTCTGGACCANAGTGGAGTCAGAGGTGAAAAGAGGTAG
ATATTCCANAATTGAGGGATTTGTGAGGTGAAATCATTGTTACAGATATTAAAGG
30 ATAAGGAGCTTTGTCAAAGGGGATCTTAAGTTTCTGGTATGGTAACTGGGTTAGAG
AGCCCTGGAACATGACCAGCTTTAAGGGAAGAGAGCTTGAGCTCTGTTCTTGTTAA
GCTCAGTTTGAGATCTTTGTGGAATCAAGTGGAGAGGTCTAAGCAGGGAAGTGGCT
TGGCTAGGCTGTAAAGATGAATCTGAGAGTCCCAAGAATATGGTAATTATTAATAA
AAGCCTTAGGTANATGAAATTGTTTTGGG

35 Sequence ID - 897 nt: 509
GCAAATCTACACATTTGATTAAATGATAGGGAAC TATGCACACACATAATACATAT

- 231 -

AATGCTAGTTTCTTGGTTTTGATATTGTACCATAGTTATGTAAGATGTAACCATTG
GGGGAAACTGGGTGAAGGCTACATGAGACCTCTCTGTACTTAATCTTTGCAACTTA
TGTGAATCTATAATTATTCCAAAATAAAAAGTTTTAAAGAACCTAAGTATCCTTAT
TACTGAGGGTCATCGTGTCTAGACAGCAAGGTTGGGCCAGAGCTTCTAGTTATTTAA
5 AATACTAAATACCAGCCTGGGCAACATAGCAAGACCCTGCCTCTACAAAAGCAAA
AAAATTAGCTGGGCATGGTGGTACATGCCTGTGGTCCTAGTTACTCTTGGAGGAGT
CTGAGGTGGGGAGCTTGAGCCTAGGAGTTTGAGGCCGAGTGAGCCTTGATTGTGT
CTCTGTACTCCAGTCTGGGCCACAGAGCAAGACCCGGTCTCTAAAAATAAATAAAT
AAATA

10

Sequence ID 898

ANTGCACTCCAGCTTGGTGACAGAGGGAGACTCCATNTTAAAAAAAAAAAAAAAAAAAA
AAAAAAGGGAGTAGCTTGAAGCCACATAGTAGTTAGTGGTAAAGGCCACCCCTTT
TCCCACAACCTCACACCAGCACCACAAGCTAGCCTTTNTAATTTCCAAGCCAGTGCC
15 CTTTCAACGCACACACCCCTGTGTCAAGTTCCCTTTCTGCTGCAAGCTCTCTGGAGG
CAGATACTGTTGAGTCCCTGGCCTGCCTATGAGAACGGCTCATGATCTCTATTTCT
TCTGCTTAATGACCATCTCGAAGTAACAAGTTTAGCCTAAAAATAAAGTTGCTAAGT
TAGCAAAGGAAGTCCTTAGCAGCCACCATTTCTCGATTCCCTCCATCACCTCCCCTG
CCCCTCAACTCCCTCATTTCTCCCAAGATATGGGCTCCAGGCTGGGCGCGGTGGCT
20 CACGCCTATAATCCTAGCACTTTGGGAGGCTGAGGTGAGCAGATCACTGAGGTCAG
GAGTTTCG

Sequence ID 899
TCNTTCGGAACGCGCC

25

Sequence ID 900

CTGGAGGGATGGGTAGGATTTTGACAAGAGTGGTTGAAGGTATTCTAATTCACCTTA
GTACCTACATGTGCGAGGCAGCATGAAGGCAAAAAGCCTGGGGCATGTTTCAGAGA
ATAGCAAGTATTCTAGTTTGAGTGGCACCTGGTACGTATATAAGGGAATAGTAAAA
30 GATCTGGCTGGAAAGGAAAAGTAGGGGCAGGTTACGAAGGACCTCTGAAAGTCAGA
CTGTGGAACCTGGAACCTTTTATCAGGAAGCAGTAGTTAGTTTTTTCAAGCAAAAGCT
AATTAGAGTTGATATTTAGGAGGATGAATCTAACAGTTGTGTGCAAGGATGCCTTC
AAACTGAGTGAGACTAGTACTGGAGACTGGTTAAGAGACTACAACAATAACCTGAG
TAAGAATTAATACAGGCCTGACCTAGTTTTGAGTGAGTAGGATTGGAACAAGAGT
35 TTTAGGTATTATAGGATTTATGCATATAAAATGGACTTGACAGAACTTGAAGAAAG
AGAAAGTGTCAAAGGACACAGAAAGTGAGGCAGGATATCTTACAATGTAAAGGA
AAGGAATAATAGAAGTTAC

- 232 -

Sequence ID 903

GGAAACATAAGCTTGTTCAGTACACTCACGCTGTAGATTAATTCTGATATTACAT
ATCTCCATCAGACTTTGTACCCTCTCTCTTCCATCCCTTACCCTTACCGATTAGGT
TGGTATTACCTAAAAATCCATAGAAAATGTCCAGGTGAATTGCCTTATGCTTTCTA
5 CCCCATAAGGTATAATT

Sequence ID 904

CTCTGTGGTGTGAGAACACAGTGGGTGACCAAGGCTTTCAGATGAACCCAAGGAA
AGTGAAAAAGCTGATGCTAATAACCAGACAACAGAACCTCAGCTTAAGAAAGGCAG
10 CCAAGTGGAGGCACTCTTCAGTTATGAGGCTACCCAACCAGAGGACCTGGAGTTTC
AGGAAGGGGATATAATCCTGGTGTATCAAAGGTGAATGAAGAATGGCTGGAAGGG
GAGTGCAAAGGGAAGGTGGGCATTTTCCCCAAAGTTTGTGTTGAAGACTGCGCAAC
TACAGATTTGGAAAGCACTCGGAGAGAAGTCTAGGATGTTTCACAACTACAAAGC
TGAAGAAAATGAAGCCCTATTACTTGTGTTGTAAGATTTAGCACCCCTTCTGCTGTAT
15 ACTGTACTGAGACATTACAGTTTGGAAAGTGTAACTATTTATTCCTGTAAATTT
TAACCTACTAGACAATGATGTGAGTACCCAGGATGATTTCTGGGGCACAGTGGGT
GAGGAGATGGGGACAGGTGAATGGAGGAGTTAGGGGAGAGGAAAAGTGGATGGAAG
TGTCTGGAAAGGGCACCAAAAAGTCTTCCAGGTCTGATCCTGTTTCTTGCTCTGA
GTGCTAGCTACCACTGTGTACACTGTAACATN

20

Sequence ID - 905

nt: 655

CTCAGCTCTTGCCTGGTCACCTTGTGGCTTTTACCATCCTCATCCCCTGTGCCACC
CACATCCTGCCACTTCTGCATGGAGTTGGGGTGGGGCCATTGGAGAAAAGAGGTTA
AACAAGCAGTAATTTACTTGAGTACAGTCTTTGAGCCAATGAAATGCCAGTCATCA
25 TTTCCCAGGGGTACTTGTCTCTTGTCAACAACCCGCTGATAATGCTCCTTCAATG
TGAATAGCAAAAGTAGGGAGAGACGCTGAATGAAGAAGATGCCTACCCCTCAGGAA
GACTGCTGTCCGCCTCCAGGCCTGCATGCACACACCCATGCCACCTGCACCCCCA
GCACCACGCCCACACTCACTCGCACACACCCACATGCCAGTGTTTTGGGGTTGGCA
GCCTGGACACTGCTGAGGCAAACACAAGTCATCAAGCATAATTCTCATTCTCTCCT
30 TCTGTCTCTGTTTTAGTTACAGGAATTTGGTCAGTTTAGAGGATTTAATAAGTCCG
TGGAATTTTGTCTCTCTTGCTACCCACGTGAAAAGTAAGTGCATGCTTCAT
GATGTGTTTTCCCACTACCTTCCAGGCCAGCCGAGCCCACTGGCCANGGCCTGGCC
CGGTGACCTCGGTTGACACTGTCCTCANGCCACTCACTT

35

Sequence ID 906

CAGAATTTTATGCTGCACAAGGCCTGTATTTTATAATGGTGGCTCTTTTG
GACGATGACTTCCTCGATGGTGAACTTCCAGTAATCTCCCTCATCATACTGAAAT

- 233 -

GATATCAGTATATCATCAGAACACCATGGAGCTTGTCATTTGAGGGACACAGCTTG
CTTGTGTGCTTGGGAAAGAAGAGGTTTAGCATGGTTTCAGGTCAGTGATGAGTCCA
ATGATCTCTGCAAGTTCCTTAGCTCTGANAATTCTGATGTCATATGCACTTCTGC
CGCCAGAGTTGCTGCTTACTGGATGCGTAAGAAGAAAAGAAAAAAAAAAAAAA

5

Sequence ID - 907

nt: 582

CTTCCATTGGGGGTAAAGATCAAACCTTTAGGCGAGCCAGGTCTGTATCTCCATTCC
TGTCTCTGACTGCTTCCCTGTAGGGATTGTCTGCAAGCGCACACCTGCATTTTCTT
GTCCACAAGTCTATGCTCTAACTCTGTCACCTGCATGGCTGCAAATTAGCTTCCTT
10 CTTCTGCCCCCTTCTCTCTAGCTTGGATTTTGAATTTGAATGGCAGGCATGGGAT
GTCCGTGTGTGTGTACTGCTGATGTGTACAGCCGCTTGTTAGCGCTCTCATTGTCT
TCAAATGTAAGTCATTTTGGCTGGGTGCGGTGGCTCATGCGTATAATCCCACGCTT
TGGGAGGCTGAGGTGAGCTGATCATTTGAGGTTAGGAGTTCGAGACCAGCCTGGCC
AACATGGCAAAACTCCATCTCTACCAAAAATACAAAATTAGCTGGGTATGGTAGT
15 GCACGCCTGTAATCCCAGCTACTTGAATGCTGAAGCAGGAGAATTGCCTGAAACC
ANGAGGCGGAGGTTGCGGTGAGCCAAGATCACGCCACTGCACTCCAACCTGGGTGA
CAGAGCAAGGCTGTGTCTCAA

Sequence ID 908

20 ACCTGACTTCAAACCTATACTACGAGGCTACAGTAATCAAAACAGCATGGTACTAGT
ACAAAAACAGACCAATGGAACAGAATAGAGATCTCAGAAATAAACTGCACATCTA
CAACCATCTGATCTTCAACAAACCTGACAAAACGAGCAATGGGGAAAGGATTCCCT
ATTTAATAAATGGTGCTGGGAGAACTGGCT-AGCCATGTGCAGAAAATTGAAACTG
GACCCCTTCCTTACACCTTATACAAAAATTAACCTCAAGATGGATTAAAGACTTAAA
25 TGTAGAACCCAAAACGATAAAAACCTAGAAGAAAATCTAGGCAATATCATTAAGG
ACATAGACATGGGCAAAAATTTTCATGATGAAAACATCAAAGCAATGGCAACAAAA
GCAGAACTGACAAATGGGCTTCTGCACAGCAAAAGAACTATCGTCAGAGTGAAC
AGACAACCTACAGAATGGGAGACAGTTTTTGCATCTATCCATCTGACAAAAGTCT
AATATCCAGAATCTACAAGGAATTTAA

30

Sequence ID 910

CAAAAAACAAGAATTACCCGGGCTTGGTGGTGTCATGTCTGTAGTCCTATCTACTCA
GGAGGCTGAGGCTGAAGGATCACTTGAGCCCAGGAGTTTGAGGCTGCAGTGAGTGA
GCCATGATCATGCCAGTGTACTCCAGCCTTGGCAGACTGAGCAAACTTGGTCCCT
35 CGCAAAATGTTGAAGCCCAGTTTTTCACTATTAACCTGTATTTTCAGTTTCCCCATGC
TAACCTTTGAAACACTGGGGCTGGCCTGAGGGTATAAAGGCTTATTCAAACCTCAGTA
ATTTAAACTTAAAATCCTAAGGAACCTCAAAAAGTGTAATCTAGTCCAAATGGGGC

- 234 -

ATCAATTCTAAAGCATTTGCTTGTTTGAGCAGATTTTCTGTGTCTGAGGTATATAG
ATAACTTATCTTTTTATGACTAAATCCAAGTCCTTAGTTCTGTGGAATTCAAAA
TCATATTTAAAAATTGATGCTTTGTTCTATAATTAATGCTTTGATTGTATAAATAA
TAAGTATTCTTCCAAATCCCTTTTTACAGATGATGATTCTGATACCGAGACGTCAA
5 ATGACTTGCCAAAATTTGCAGATGGAATCAAGGCCNGAAACAGAAATCAGAACTAC
CTGGNTCCCAGTCCTGTNCTTAAAATTCTAACTCGAC

Sequence ID - 911 nt: 595
GAGGGTGTAGAAGAGAAGAAGAAGGAGGTTCTGTGCTGTCANAAACCCCTTAAGAA
10 AAAGCGAAGGAATTTTCGCAGAGCTGAAGATCAAGCGCCTGAGAAAGAAGTTTGCCC
AAAAGATGCTTCGAAAGGCAAGGAGGAAGCTTATCTATGAAAAANCAAAGCACTAT
CACAAAGGAATATAGGCAGATGTACAAANCTGAAATTGGAATGGCGAGGATGGCAAG
AAAAGCTGGCAACTTCTATGTACCTGCAGAACCCAAATTGGCGTTTGTTCATCAGAA
TCAGAGGTATCAATGGAGTGAGCCCAAAGGTTGAAAGGTGTTGCAGCTTCTTCGC
15 CTTTCGTCAAATCTTCAATGGAACCTTTGTGAAGCTCAACAAGGCTTCGATTAACAT
GCTGAGGATTGTAGAGCCATATATTGCATGGGGGTACCCCAATCTGAAGTCAGTAA
ATGAACTAATCTACAAGCGTGGTTATGGCAAATCAATAAGAAGCGAATTGCTTTG
ACAGATAACGCTTTGATTGCTCGATCTCTTGGTAAATACNGCATCATCTGCATGGA
GGATTTGATTCATGAGATCTATACTGTTGGAAAAC

20
Sequence ID - 912 nt: 651
CATTTCAGAGTTTATGTGAATTGAATTGAACTATGGTTTTATGTTACTGTCAGTA
GAATGAAGTACGAATATTTGAAAAATACACCTTCAACTTCAAAGTGATTCTTGACA
AAAATTATAAGGAATCATTTTGGACACATTTTCTGGTAGAGCCTTGTAATAATTAA
25 AACCAAGTGTTGTTTTCAAGAAGAACTGTAATACATAATCAGGAATTTGAGTAGGG
AGATTATTTTGTATTTTAAATTAAGTGGCTGTGTAGTTTTAACTTTAGTATTGC
AGGTAGAGTAAGCTTACATGATAACAAAAATCTTGGTCTTAGTGACTTAATGATT
TGATATTTATTGATTGATTGGTTATCATTCCAAATATTTTAAAAGATAATAGCTGG
CTGGGTGCGGTGGCTCATGCCTGTAATCCCAGCACTTTGGGAGGCCAGGACGGGCG
30 GATCACGAGGTCAGGAGATCAAGACCATCCTGGCTAACACGGTGAAACCCCGTCTC
TACTAAAAATCAAAAAATTAGCCGGGTGTAGTGGCGGGCACCTGTAGTCCCAGCTA
CTCAGGAGGCTGAGGCAGGAGAATGGCATGAACCTGGGAGGCGGAGCTTGCAGTGA
GCTGAAATCGTGCCACTGCCTCCACCTGGCGACAA

35 Sequence ID 913
GTGAGGTGGGGACTTCATTCATTGTCCTATTTCTATCTCCACTTTGTGCCTGGAGA
GCTTTCAGGGGAGGTGGAGGAGGAGGTCTGCCAAGCTACTGCAACATCTGTCACC

- 235 -

CACTATACCCAGTTACTTGGGGGAGGACAGACACTGTGGTGTCAATAAGTTGTTT
GAACCAAAGTGGCGGCTGCATCTTTGTCCCGATGCTAGCCGTGCCGGTCTCCCATC
ATCCGCTCGCCCTCCTTTNCCCTGGGCTGCGCCCACTTGTCTTCCTGGATATTTGG
GGGTGACTCGCCATGCTTGGCACCCCTCTGCTTCCTGGTGCTGCTCTGACTCGAAGA
5 CGGGACAGTCCCTGGTGCACATCCAGGGAAGAGGAGTGTGGTAGTTCTTGCACTA
GGCACTTTATCAGGACCTGACCTGTTGCTGGGTGATTTTAGTCTCTACAAACAGAA
AGCGTTTCAAAGCGTCAGCTGTGGGAGCAGAGTGACCCTTTGCTGATGCTGGGGGG
AGGGGATCTAAATCCTCATTTATCTCT

10 Sequence ID 914

GGCGCCTGCTGGAGGAGGAGAGCTCTGCTGGCATGAGCCACAGTTTCTTGACTG
GAGGCCATCAACCCTCTTGGTTGAGGCCTTGTTCTGAGCCCTGACATGTGCTTGGG
CACTGGTGGGCCTGGGCTTCTGAGGTGGCCTCCTGCCCTGATCAGGGACCCTCCCC
GCTTTCCTGGGCCTCTCAGTTGAACAAAGCAGCAAAACAAAGGCAGTTTATATGA
15 AAGATTANAAGCCTGGAATAATCAGGCTTTTTTAAATGATGTAATTCCCACTGTAAT
AGCATAGGGATTTTGAAGCAGCTGCTGGTGGCTTGGGACATCAGTGGGGCCAAGG
GTTCTCTGTCCCTGGTTCAACTGTGATTTGGCTTTCCCGTGTCTTTCCTGGTGATG
CCTTGTTTGGGGTTCTGTGGGTTTGGGTGGGAAGAGGGCCATCTGCCTGAATGTAA
CCTGCTAGCTCTCCGAAGCCCTGCGGGCCTGCTTGTGTGAACCGTGTGGACAGTGG
20 TGGCCGCGCTGTGCCTGCTCGTGTTCCTACATGTCCCTGGCTGTTGAGGCGCTGC
TTTAACCTGCACCCCTNCCTTG-CTCATANATGCTCCTTTTGA

Sequence ID - 915

nt: 230

TTTGAGACCAGCCTAGCCAACATGGTGAAACCCCATCTCTACTAAAAATACAAAAA
25 TTAGCCGGGCGTGGCGGCACATGCCTATAATCCCACTTACTTGGGAGGCTGANGTA
GGAGAATCGCTTGAACCCANANAGGCAGAGTTTGCAGTGAGCCGAGATTGTGCCAT
TGCACTCCAGCCTGGGCGACAGAGCGAGACTCCATCTAAAANAAAATAAATGAATA
AAATAA

30 Sequence ID 917

NNCAGATTTTTTTTTTTTTTTTTCAGNGTTAGACCATCTTCAATTCTTGGAACAAAC
TTAACTTTCATGATATGTATTTTTTATACATTTGCTGGATTTTATTTGCTAATATT
TTACTTAGGATTTAATTTTCTAAGTNGACCTATAATTNTCCTGTATAAAATTGCAT
TTGTCACATTTTAGTATCAAGGTTGTCCTANCNCCATGAAATGGATTTANAATGGT
35 TTATGTAANATAAAGTACATTTCTTCTAAAGGTTTGNGTGGATTAACTTTCAAATC
TGCCANAGNGNGTTTTTTTTTCTTTTTTTTTTTTTTTTTCATTTNAAGGGAGNGCAAGT
ANCTTTTCAAATNCTGATTTAATTTTAAAATATTTNCAAGTNTNTTTANAGTTTTT

- 236 -

TATTNTNTNTNGAANGTTAACATTTTTATANAAAANGGTNTTATCTTTTTAAATTC
TTTGACATCAGTTTCTTCANAATTCCTTCTTTTAA

Sequence ID 926

5 GTCATATCTCTTCCCAGGGAAAGCAGGAGCCCTTCTGGAGCCCTTCAGCAGGGTCA
GGGCCCCCTCGTCTTCCCCTCCTTTCCCAGAGCCATCTTCCCAGTCCACCATCCCCA
TCGTGGGCATTGTTGCTGGCTGGCTGTCCTAGCAGTTGTGGTCATCGGAGCTGTG
GTCGCTACTGTGATGTGTAGGAGGAAGAGCTCAGGTAGGGAAGGGGTGAGGGGTGG
GGTCTGGGTTTTCTTGTCCCCTGGGGGTTTTCAAGCCCCAGGTAGAAGTGTTCCTT
10 GCCTCATTACTGGGAAGCAGCATCCACACAGGGGCTAAGCAGCCTGGGACCCTGT
GTGCCAGCACTTACTCTTTTGTGCAGCACATGTGACAATGAAGGACGGATGTATCA
CCTTGATGGTTGTGGTGTGGGGTCTGATTTTACGCATTTCATGAGTCAGGGGAAGG
TCCCTGCTAAGGACAGACCTTAGGAGGGCAGTTGGTCCAGGACCCACACTTGCTTT
CCTCGTGTTTCCTGATCCTGCCTTGGGTCTGTAG

(15

Sequence ID 938

TGGCCATCCTTTTTCCCCCAAACACACCCCCCTTAACCTATCTCTTGGGACTTAGCC
CGACCCTCCCTCTCATTTCCCATTAAGTCTGAGAGGCAAGAGCTAGGTTAGGCAAG
GAGGTGGTTGGCCAGAGATGGGGAACAGCCAGGTGCCCCAGTCCTCTGATTTTTCC
20 TCCATCCTGCTTACCACCTCCCTGGGTACTTACAGCCTTCTCTTGGGAACAGCCGG
GGCCAGGACTGGGTACCTATGAGCTGAATCAGCATCTCCTCCTGAGTCCCAGGGC
CCCTGCAGTTCCAGTCTCTTCTGTCCTGCAGCCCTTGCCTCTTTCCCACAGGTTC
CACTTTATATCCACCTTTTCCTTTTGTTCATTTTTTATTTTTTATTATTTAT
TAAATGATGTGGTCTATGGAAAAAAAAATAAAATCTGACTTAGTTTT

25

Sequence ID - 939

nt: 513

GGAACCCAGTGTATTACCTGCTGGAACCAAGGAACTAACAATGTAGGTACTAGT
GAATACCCCAATGGTTTCTCCAATTATGCCCATGCCACCAAACAATAAAACAAAA
TTCTCTAACACTGCAAAGAGTGAGCCATGCCTGTAAACACTGTAAAGAATGTAACA
30 TGTGGGGGACACACAGGGGCAGATGGGATGGTTTAGTTTAGGATTTTATTAGTGCA
TGCCCTACCCTCTGGGGGAACGTCCCATCTGAGGTTTTCTTCTCGGTGGGGGGATT
TAACTTCTGTCCTAGGGAAAACAGTGTCTGATGAGGAGTGTTC AACACAGGCTA
CATGAATTTCCCTATAACAGTGCGAAAGCAGCCAGGAGTCCCCGTTGGAAAAGAAC
AATGCCACTCTCTTTTATGTATCTTGGTTCTGCAACTCATTTGTTGTAAGTAGGGT
35 TAATCGAGTATCAGGTTACAGTATCCTGCCCTTATTATTTTATGATTCACTGACT
CAAGTTCCA

- 237 -

Sequence ID 947

GAGAGTGAAAAAATTCTGGTACAAATTGGGAAATTAGTATATAACAACATAGTGTT
AAATTCAATGGGAAAAGTTTAATAAGAGGATTTGGTATCAACTGGCTGTCCAAAGA
TAAAAATGGACCGTCCTATCACATACAAAATTGTTTTTTAGATAAAGATTTAAATA
5 CAGGCACTCCTTCATTTGCGTGGTGCACCTTGAGGTGTTGCAGAAATGATGAGAGC
TGAAACTGCAAAGCAATTTTAATACTTTATCTGTTGGAAATCTTATAGTTTTCTCTG
TGACCGTTAAAATTTTCATTAAACTATTAAAAACACCCATGACTGGTCACAAATGT
ATTGGGAAATGGAAAAGAATTAATACACTAAAAATACAAAAATAGAAAATATTTA
AAATTATCTAAAAATTTGAAACATTAGAAAAATTGAGAACTAGGCAGGGCGTGGTG
10 GCTCACATCTGTAATTTTAGCCCTTTGGGAGGCTGANGCAGGTGGATCACCTGANG
TCAGGAGTTGAGACCAGCCTGCCAACGTGGGGAAACCCCGTCTCTACTGAAAATA
CAAAAATTANCCGGGCATGGTGGCACAAGCCTGTAATNCTTGCTNACCAGGANGCT
GAGGCAGGAGAATCACTTGAACCCANGANG

15 Sequence ID 949

GTTTCACATGAGAAGGTAGTATTATGTACAGTGACCTTGTTTAAAGTGTCNGTTTA
ATGTTACCACTAAGGCCCTGCCCCAGCTTTATCACCTGAGCACTAACAAGTGCTGT
GTGGAGTTCAGTCCATGCTGGTAACTNTTGAGTATTCAGTGGGTCTTTTAAACAATT
ACCACCGTGGAGGANANAGCAAGGAAGAGAAATGCTGTGATCTTTTNCTGTTTTTA
20 ATTAGNGAAAGAGGGATTANATTAAACAAATGTTACAGAGNTGTGACTNTGATCCC
CCAGNGGTAAGCAATAATTGTANAGACTGGATTTNANAAGCCCTGAGAGTTTATTT
TCAACCTATNTATTATAGNNCAATCC

Sequence ID 1028

25 ACAAGGCTTGGGGGCTGGACTCCCTCTACTGCCTCTGGCCATACCCCCTCCTGGAG
ATGGGGTCAAGGCACCAGGACTGA

Sequence ID - 1056

nt: 435

TCGCTTGTAAGCCTGAGACAGCTGCCTGTGTGGGACTGAGATGCAGGATTTCTTC
30 ACACCTCTCCTTTGTGACTTCAAGAGCCTCTGGCATCTCTTTCTGCAAAGGCATCT
GAATGTGTCTGCGTTCCTGTTAGCATAATGTGAGGAGGTGGAGAGACAGCCCACCC
CCGTGTCCACCGTGACCCCTGTCCCCACACTGACCTGTGTTCCCTCCCCGATCATC
TTTCCTGTTCCAGAGAAGTGGGCTGGATGTCTCCATCTCTGTCTCAACTTCATGGT
GCGCTGAGCTGCAACTTCTTACTTCCCTAATGAAGTTAAGAACCTGAATATAAATT
35 TGTTTTCTCAAATATTTGCTATGAAGGGTTGATGGATTAATTAAATAAGTCAATTC
CTGGAAGTTGAGAGAGCAAATAAAGACCTGAGAACCTTCCAGA

- 238 -

Sequence ID 1071

NGATATAGTNCCGCATGGGAAAGATGANCAGGTATAACCNAGCNTNATATAGCAAG
GACTAACCCCCCTGCCTTCTGCATAATGAATTAAGTAGAAATAACTTNGCAAGGAG
AGCCAAAGCTAAGACCCCNAAACAGACGAGCTACCTAAGAACAGNTAAAAGAGC
5 ACACCCGTCTATGTAGCAAAATAGTGGGAAGATTTATAGGTAGAGGCGACAAACCT
ACCGAGCCTGGTGATAGCTGGTTGTCCAAGATAGAATCTTAGTTCAACTTTAAATT
NGCCACAGAACCCCTCTAAATCCCCTTGTAATTTAACTGTTAGTCCAAAGAGGAA
CAGCTCTTTGGACACTAGGAAAAAACCTTGTAGAGAGAGTAAAAAATTTAACACCC
ATAGTAGGCCTAAAAGCAGCCACCAATTAAGAAAGCGTTCAAGCTCAACACCCACT
10 ACCTAAAAAATCCCAAACATATAACTGAACTCCTNACACCCAATTGGACCAATCTA
TCACCCTATAGAAGAACTAATGTTAGTATAAGTAACATGAAAACATTCTCCTCCGC
ATAAGCCTGCN

Sequence ID - 1074

nt: 689

15 GGGAGGCGGAGGCTGCAGTGAGCTGAGATCGTGCCACTTCATTCCAGCCTGGGCAA
CAAAGCGAACTCTGTCTCAAAAAAAAAAAAAAAAAAATTTGTTGACTGTTGTAA
TTTAAAGCTTGTCAATTTTTTATTTAGTAATAACACTCATTAGTGTAGTATCTATGA
TGAACCAGGTTCTGCACAAAGTACCTTATGTTTCATGGCCTCATATCGTCTTCTCCA
AAACTCTGCAAGATAGGATTCATCACCCTTATAGGGAGAGATCTGAAAGTTTAAA
20 ATTGTACCCAAGGTCACACAGCTGGTAAGTGCCAGAGCTGGGATTCCGTAGGGTGT
TCANAGTGCCTCTCCTGCCGTAGGCTTATCACAAAAGTCAAAGTTTGGTCATAAT
AAAGCCTGAAGTTTGGCAGGATTTAAAAATAGTCACCANACTTTTGAGTTGGAGCA
TCCCACCTCACTGCTGTTTACCTTCTGTGGCAGGGAGAGTCATCATTTCCATTTCA
GCTTGTGGAATATCTTGTCAATTAACATTCTCATGCAAAAGCCATTTTATGGTGCCC
25 AATGAANATGGTTAAGCTACTGCCCCAAGCCTNTGGAAGCCTTCCTAATTTTGGAC
TTGCACTATGCAAATTGNATAATATTTTCTCTACCCTAAGCCAAATATTTTCTTCA
CTTTTCATTCATTCTAC

Sequence ID 1081

30 CGCCGCCGCGCCGCCGTGCTCTCCAACGCCAGCGCCGCTCTCGCTCGCCGAGCT
CCAGCCGAAGGAGAAGGGGGGTAAGTAAGGAGGTCTCTGTACCATGGCTCGTACAA
AGCAGACTGCCCGCAAATCGACCGGTGGTAAAGCACCCAGGAAGCAACTGGCTACA
AAAGCCGCTCGCAAGAGTGCGCCCTCTACTGGAGGGGTGAAGAAACCTCATCGTTA
CAGGCCTGGTACTGTGGCGCTCCGTGAAATTAGACGTTATCAGAAGTCCACTGAAC
35 TTCTGATTGCAAACTTCCCTTCCAGCGTCTGGTGCGAGAAATTGCTCAGGACTTT
AAAACAGATCTGCGCTTCCAGAGCGCANCTATCGGTGCTTTGTCAGGAGGCAAGTGA
GGCCTATCTGGTTGGCCTTTTTGAAGACACCAACCTGTGTGCTATCCATGCCAAAC

- 239 -

GTGTAACAATTATGCCAAAAGACATCCAGCTAGCACGCCGCATACGTGGAGAACGT
GCTTAAGAATCCACTATGATGGGAAACATTTTCATTCTC

Sequence ID - 1083 nt: 198

5 GCGCGTCGACTTTGTTTAGACATTGAATGACTTTGTTAAAGGCACAATTAATCACA
TTGGTTGTACTCTGNNGACAGCCTTCTTTAAAAAAAATAAACAATTTAAAACAA
AA
AAAAAAAAAAAAAAAAAAAAAAAAAANTTTTAACC

10 Sequence ID - 1084 nt: 198

GCGCGTCGACTTTGTTTAGACATTGAATGACTTTGTTAAAGGCACAATTAATCACA
TTGGTTGTACTCTGNNGACAGCCTTCTTTAAAAAAAATAAACAATTTAAAACAA
AA
AAAAAAAAAAAAAAAAAAAAAAAAAANTTTTAACC

15 Sequence ID - 1099 nt: 561

TGCATGCTTGTGGATTGGAAAACTTTGGAGACTGATTACTTTTCATTATATATGT
GTCACAGTGAAACAGCTTTTATGTGTCATGTAAGATTACTGCTTGCCTCTCTAAGG
AAGGTCGTGACTGTTTAAATAGACGGGCAAGGTGGAACCTTTTGAAAGATGAGCTT
20 TTGAATATAAGTTGTCTGCTAGATCATGGTTTGTATTGAACATAACAAGGTTTGCAG
ATCTGCTGACTTATATAAAGCTTTTTGATTCCCTACTAAGCTTTAAGATTTAAAAAA
TGTTCAATGTTGAAATTTCTGTGGGCTCTATTTTGTCTTTGGCTTTCTGGTGAGA
GAGTGAGGAAGCATTCTTTCCTTCACTAAGTTTGTCTTTCTTGTCTTCTGGATAGA
TTGATTTTAAGAGACTAAGGGAATTTACAACTAAAGATTTTAGTCATCTGGTGGA
25 AAAGGAGACTTTAAGATTGTTTAGGGCTGGGCGGGGTGACTCACATCTGTAATCCC
AGCACTTTGGGAGGCCAAGGCAGGCAGAACACTTGAAGGAGTTCAAGACCAGCGTG
G

Sequence ID 1109

30 TTTGNCGGTNTTGGANNNNNANAANTTTCTTCCANNCNTNACNTNTTGGTGGNCTA
AATTAANATGGNTTTNNGGGGTTCTTNTNCTNNNTNNNNCATGGGANANAATTNATT
NTCNTNCNNNTTCTTNNCCCTNAANCTACCTTCCCCCNATTTTCTCCCCTNTTCN
TNAATTANCATCCTCTCCNCNTANNTCNANACNTTAATGGCAANACTATCTAATAN
CNANNATAANANCTCCTGTNNNCCACATNTCTTATTNNNCGCNNCANGTTNCANNC
35 CCNCAGAGTNAACTCATCCTCNNCNAANTTCATATCGTGNNCTNTNNNCNNNTNGC
GCGANATATTAANNANACCNGTANNTNNNANACANNANNTNNGNAANAANCCTTCT
NANNTTTTAGCNTCNNGCNNTAACNNNNNTCTTNGTGNNNNCNCAGCTTTCNCNNC

- 240 -

ATNATNCTNCNNCGAANTNTCANNCNTCTCCNCTTNAATGNNTTCCCATGNATTAA
NTNCCTCGNNNANAGCACTATCGTNNNGAGNNNATTATNGNCNNTTTACNTCATG
TGGTCCANTNNCGTTNGNCGCNNNNAATNTTCGTNNNNCNN

5 Sequence ID 1118

GGATTTTAGAGGAAGGCGCTNGGTTACATTGGAGAAGTGGAGTGGTCTGGAGTTCC
ACGGTGTAGTGGACCAGAGGCCACCTCTCCTGGGCTTCTCAGTGTCTCGCCGGCGG
GGTTCGGCCTGAGCTGGATTGACATAGCCCTTGGCGGATTTAAACAACCTAAACAT
TAAGCAGTACAGCTGCCTCAAACCTTTGGGATTTTCAGAATGACTGACACTGCCGA
10 AGCTGTTCCAAAGTTTGAAGAGATGTTTGCTAGTAGATTACAGAAAATGACAAGG
AGTATCAGGAATACCTGAAACGCCCTCCTGAGTCTCCTCCAATTGTTGAGGAATGG
AATAGCANAGCTGGTGGGAACCAAAGAAACAGAGGCAATCGGTGCAAGACAACAG
ACAGTTCAGAGGCAGGGACAACAGATGGGGGTGGCCAAGTGACAATCGATCCAATC
AGTGGCATGGACGATCCTGGGGTAACAACTACCCGCAACACAGACAAGAACCTTAC
15 TATCCCCAGCAATATGGACATATATGGTTACAACCAGCGGCCTCCTTACGGTTACTA
CTGATAGAAATGTTGGCAGCTTTTAGTAAAAGCATTTACTCTGTTACCATGAGAAA

Sequence ID 1125

NGACTGGCTCCCGAAAAGAAGGGTGGCGAGAANAAAAGGGCCGTTCTGCCATGGA
20 CGAAGTGGTAACCCGCGAATACACCATCAACATTNACAAGCGCATCCATGGAGTGG
GCTTCAAGAANCGTGCACCTCGGGCACTCAAAGAGATTCGGAAATTTGCCATGAAG
GAGATGGGAACCTCCATATGTGCGCATTGACACCAGGCTCAACAAANCTGTCTGGGC
CAAAGGAATAAGGAATGTGCCATACCGAATCCGTGTGCGGCTGTCCANAAAACGTA
ATGAGGATGAAGATTCACCAAATAAGCTNTATACTTTGGTTACCTATGTACCTGTT
25 ACCACTTTCAAAAATCTACAGACAGTCAATGTGGATGANAAACNAATCGCTGATCGT
CAGATCAAANAAANT

Sequence ID - 1139

nt: 503

CAGCACTGCCAGTGGAGATGGGCGTCACTACTGCTACCCCTCATTTACCTGCGCTG
30 TGGACACTGAGAACATCCGCCGTGTGTTCAACGACTGCCGTGACATCATTACGCGC
ATGCACCTTCGTCAGTACGAGCTGCTCTAAGAAGGGAACCCCCAAATTTAATTAAA
GCCTTAAGCACAATTAATTAAGTGAACGTAATTGTACAAGCAGTTAATCACCC
ACCATAGGGCATGATTAACAAAGCAACCTTTCCTTCCCCCGAGTGATTTTGCGAA
ACCCCTTTTCCCTTCAGCTTGCTTAGATGTTCCAAATTTAGAAAGCTTAAGGCGG
35 CCTACAGAAAAAGGAAAAAGGCCACAAAGTTCCCTCTCACTTTCAGTAAAAATA
AATAAACAGCAGCAGCAACAAATAAAATGAAATAAAAGAAACAAATGAAATAAA
TATTGTGTTGTGCAGCATTAAAAAAATCAAAATAAAATTAATGTGAGCAAAG

- 241 -

Sequence ID - 1148

nt: 587

5 TGAAAAATAAAGTTTTTATGTATATTCTACATATGTATATGTTGGTAGAAAGCAAA
AACGCTAGGTAAAAATAAATGTAATACAATTTTAGCTATGAACCAAAAACCATTT
GTGGTGTGGATGCAAGAAAGTCTGGATGGGTGCAGAGTTCTCCATGTTTCACCTTCT
10 GACATTTGAAAATACGCAGTTTGCATTTGATACGTCAAATGTTATTTTTAAGAAAA
CCAATAAAATCATTAAAACCGAAAAGGCAGTTTGTCTTGTTTTTACCTTAGTTGGA
GTTATCTGCAATTGCCGTATTAGTGTTTTAAGGAACTTGTAAGTAAGCTCCTTAGT
CCCCTTTAGAGCTACGAAACATGTCAATTTTACTTTTCTCCAGCTTTTTTGGAATCT
TATCTAAATTACCATGTAGAGTTCTGCATAGCTTCAAATTCTCTTAGCCAATGTGG
15 TCTGTAAGTGTCTATCGATGAATTTACCGTTAATTGCCGTAGTATACTGTCCTGT
ACCGGATGTGAAGAGGAGCAACTCTGCACAGTGCCTGGTTGCTCCCATGGTAGGA
ANGAATGGCTTATCAATGGTTCGGATTT

Sequence ID - 1160

nt: 650

15 GGAGGATGGAGCAGTGAGCGGGTCTGGGCGGCTGCTGGCAGCGCCATGGAGACGGT
ACAGCTGAGGAACCCGCCGCGCCGGCAGCTGAAAAAGTTGGATGAAGATAGTTTAA
CCAAACAACCAGAAGAAGTATTTGATGTCTTAGAGAAACTTGGAGAAGGGTGAGTG
TAAAGAACTATAGGTAGGTCATTGGGTCCCAGTCTTTTTCTGCCCCAGAAGAAG
CAGAAGGATATGAACCTTTCAGCATTGTTCTAGGTGGGGTGAAGGTAAATTTACA
20 GCTTGTGATGTCCTTCTTCGCTTTACTCCAATCCCTATTATAGACAGATTTAGTGA
TTCCTGGTCTTTTTTAACACGAAGAATATCTATTGTTTTCTCTTTTGTAGGATCTGT
ATGATTTTATCTACTTAACAGATAGCACTAATTAGATTAAATTTCTATAAGAACT
TTTTAATTTGCTGTTTCAATAATTTCTGATTGGTATGCAATACTGTTTCAATGAAAA
TCAATGTAATTTAGTATTTTAATATTTGCACCTTTGTGAAATATAGTAAATAAATT
25 AAGCACTATCACACCTTCACAGCTACTTAGGAGATCCACAATCCTGGGTGGGAG
CCAGTGGATTTCTGAAACACAGATTTGTTAATG

Sequence ID - 1165

nt: 502

30 CTCAAGTGAATCCTGGCTTCTTGGAAGCGCTTGCCTAGACGAGACACAGTGCATAA
AAACAACTTTTGGGGGACAGGTATGTTTTCTTGCAGCTGCGGTTGTAAGGTCTTGG
CAAGACAAGCAGTGTGGCCAGAATTTTGAACCTTCTGATGAATGTGTAATGCAAAGG
ACCTTGTACATTTTTTTTGTTC AAGGTCTCAAATGAGCACATGAAGAGGTTGCT
GTGAACTTTAAGTGGCCCTACTGCGCAGAAGCATTTCAGATGTCACTTGATGATCT
GTAAGGGAACCTTGCTGATTTGGGAATGTGCTTAGGGAACACACATTCCTTTTGACA
35 GGGTCTGTCACTGGGTGGGTGATGAATTATACAGATGACATGTGCTTTTTTTTCTT
TTTTCAACCTCAATGGTATTCCTACAGGAAATGGATAACCATTTTAACTGTATTTT
TTGCAGCCCGTACCTTCTTGGGAATACAATTGTCTAACTTTTTAATTTTGGTCT

- 242 -

Sequence ID - 1172

nt: 648

CCACAATAATAAGAGAAAAACAGGAGCAAAAGGATATACAAAACCACCAGAAAACA
AATAACAAAGTGACAGGAGTAAGTCCTTAACTGGCAATAATAACCATGAATCTAAA
TGGATTCCATTTCCCACTTAAAAGATAAAGACATGCTGAATGGATAAAAAGCTGTC
5 ACCCAGTTATATGCTGCCTACAACAACTCACTTCACCTGTAAACATACATATGGA
TGGAAAGAGAAGGCATGGGAAAAGATACTCTACTCAAATGAAAACAAAAACCAAAC
AAAGGTGGCTATTCTTATATGAGATAATACAGACATTAAATCAAAAACCTGGAAACA
AACACAAAGTCATTGTATAATGATGAATTCAATTATATCATGATGAATTCAATTAT
ATCCTCCTTCCTGATCAATTCAGAAAGGAGGATATAATCTTTTTTAAATATATATAC
10 ACCCAACACCAGAGCATATAAATATGTAAAGGAAGATAAAGGGAGTCCTGTGATCA
AGAATAAATATAACAATTATAAATATTTTATCTAAAGTGATAGATAGACTGTAATA
CAATAATAGGGTGGTGACATTAACACCCCTCTCAGATTGGACTGATCATCTAGAA
GGGAGAAAAAGCTTTATGATTGGAAAAGCCAT

15 Sequence ID 1178

ATTGTGTTGGCCACCCGGAATTCGCGGCCGCGTCGACCTACGCACACGAGAACAT
GCCTCTCGCAAAGGATCTCCTTCATCCCTCTCCAGAAGAGGAGAAGAGGAAACACA
AGAAGAAACGCCTGGTGCAGAGCCCCAATTCCTACTTCATGGATGTGAAATGCCCA
GGTGAGGAGACGGCTTGCTGTAGTGGGGAAAGCACTGGACCTCAACAGTTGGAAAA
20 TGTGTAGTGTTAGCTGTCTCGTATCCTTGAAGCTGTGCAGCAGCTTCAGTTTCTT
CGCCTGTGGAAAATATTTTCCCTGATACTCTTAAAATTTGAATGTATGAGACTGGC
AAAGTTTTCATCTTAGGAGGAGTGATTCAATTCACCGTGATCTCTCATCACATTT
CACATACAACCCCTACGTTTTTTTTGTGTTGGGAAACAATGTAATGGATGATGAGTT
GGGCATAAGTGCAGGAAAGACGGGTGTAATAGAGGAAAAAAATGTTATCTGCTTTT
25 CTTTCAGGATGCTATAAAATCACACGGTCTTTAGCCATGCACAAACGGTAGTTT
GTGTGTTGGCTGCTCCACTGTCCTCTGCCAGCCTACAGGAGGAAAAGCAAGGCTTA
CAGAAGGATGTTCTTCAGGAGGAAGCAGCACTAAAAGCACTCTGAGTCAANATGA
GTGGGAAACCATCTCAATAAACACATTTTGGAT

30 Sequence ID - 1180

nt: 622

CTTTTCCTCCCGCTGTCCCCACGGGAGGGGACTGCTCTCCCCGCTGCATCCTTT
CTGTGAGGTACCTTACCCACCTCAGCACCTGAGAGGGTGAAATAGAATTCTAACCT
CGACATTCGGGAAGTGTTTTTGAGAAGTCTCGGTGCGTAAGGGAAGTCTTCCAAGT
CCGTGCAGCACTAACGTATTGGCACCTGCCTCCTCTTCGGCCACCCCCCAGATGAG
35 GCAGCTGTGACTGTGTCAAGGGAAGCCACGACTCTGACCATAGTCTTCTCTCAGCT
TCCACTGCCGTCTCCACAGGAAACCCAGAAGTTCTGTGAACAAGTCCATGCTGCCA
TCAAGGCATTTATTGCAGTGTACTATTTGCTTCCAAAGGATCAGGCCCTGAGAACA

- 243 -

ATGACCTTATTTCTACAACAGTGTCTGGGTTGCGTGCCAGCAGATGCCTCAGATA
CCAAGAGATAACAAAGCTGCAGCTCTTTTGATGCTGACCAAGAATGTGGATTTTGT
GAAGGATGCACATGAAGAAATGGAGCAGGCTGTGGAAGAATGTGACCCTTACTCTG
GCCTCTTGAATGATACTGAGGAGAACAACTCTGACANCCACAATCATGAGGATGAT
5 GTGTTG

Sequence ID - 1181 nt: 155
CGCCACTTATCCAGTGAACCACTATCACGAAAAAACTCTACCTCTCTATACTAAT
CTCCCTACAAATCTCCTTAATTATAACATTCACAGCCACAGAACTAATCATATTA
10 AA

Sequence ID 1182
CATTGTGTTGGCNCCCGGAATTCGCGGCCGCGTCGACTTTTTGTGTTGTTTGGAG
CAGAAATACTAAAGAAGATTCGGGGCCGAGTATCCACAGAAGTGGACGCAAGGCTC
15 TCCTTTGATAAAGATGCGATGGTGGCCAGAGCCAGGCGGCTCATCGAGCTCTACAA
GGAAGCTGGGATCAGCAAGGACCGAATTCTTATAAAGCTGTCATCAACCTGGGAAG
GAATTCAGGCTGGAAAGGAGCTCGAGGAGCAGCACGGCATCCACTGCAACATGACG
TTACTCTTCTCCTTCGCCCAGGCTGTGGCCTGTGCCGAGGCGGGTGTGACCCTCAT
CTCCCCATTTGTTGGGCGCATCCTTGATTGGCATGTGGCAAACACCGACAAGAAAT
20 CCTATGAGCCCCTGGAAGACCCTGGGGTAAAGAGTGTCACTAAAATCTACAACCTAC
TACAAGAAGTTTAGCTACAAAACCATTTGTCATGGGCGCCTCCTTCGCAACACGGG
CGAGATCAAAGCACTGGCCGGCTGTGACTTCCTCACCATCTCACCCAAGCTCCTGG
GAGAGCTGCTGCAGGACAACGCCAAGCTGGTGCCTGTGCTCTCAGCCAAGGCGGCC
CAAGCCAGTGACCTGGAAAAAATCCACCTGGATGAGAAGTCTTTCCGTTGGTTGCA
25 CAACGAGGACCAGATGGCTGTGGAGAAG

Sequence ID - 1183 nt: 479
CGTGGCAGCCATCTCCTTCTCGGCATCATGGCCGCCCTCAGACCCCTTGTGAAGCC
CAAGATCGTCAAAAAGAGAACCAAGAAGTTCATCCGGCACCAGTCAGACCGATATG
30 TCAAAATTAAGCGTAACTGGCGGAAACCCAGAGGCATTGACAACAGGGTTCTGTAGA
AGATTCAAGGGCCAGATCTTGATGCCCAACATTGGTTATGGAAGCAACAAAAAAC
AAAGCACATGCTGCCCAGTGGCTTCCGGAAGTTCCTGGTCCACAACGTCAAGGAGC
TGGAAGTGCTGCTGATGTGCAACAAATCTTACTGTGCCGAGATCGCTCACAATGTT
TCCTCCAAGAACCGCAAAGCCATCGTGGAAAGAGCTGCCCAACTGGCCATCAGAGT
35 CACCAACCCCAATGCCAGGCTGCGCAGTGAAGAAAATGAGTAGGCAGCTCATGTGC
ACGTTTTCTGTTTAAATAAATGTAAAACTG

- 244 -

Sequence ID - 1185

nt: 628

CTTTGATTACCTTTGAGTATTAGGTTGAAAGCTTCTCTGTGCTTGATTGAACATTG
TGATGATGTTGATTGGGTCATGTCAGATTTAGACAGTGTTGTGTTTAAGATAAATG
TTTAATGGCTCTTAGCAGTGTTTCATGCCTCCCCTTTTCCCCTGATACTTTAAAAAC
5 AGAATATACAGAAAAGGGGAGTTGGGTGAAGAATCACCATATTCTCATTACCAGAG
TAGTGTCTACCAGCTGTTTTCACATTTTTCTGTTTCCTTCTGTCCTTGGAATCCTT
TTTTTAGATCCTTGTAATACTAGTAAAGATATTCCACTCTGTGTTGTAAGCATT
TCCATTTTGCTCCATGGTCTTCATAATGCCCTGTGGTCCTTTATTAAGGGGATGCA
CCATGTAGAGGTGAAAGGCTTTCCTTGACTTGGCCACCATTTCTGTATTTTCCTTA
10 GAGGAGGAGGTTTCCAACATTTCTTTTTTAGAGACAGAGTCTCGTTCTGACACGCA
GGCAGGAGTGCAGTGGCATGATAACAGCTCACTGCAGCCTCGAACTCCTGGGCTCA
AGTTATCCTCCACCTCAGCTTCTGAGTAGCTAGGACTGCAGGTGCCTGCCACCA
CACCCAGCTAAT

15 Sequence ID - 1186

nt: 494

CAGCCCTCCGTACCTCTTCACCGCACCTTCGGACTGCCCCAAGGCCCGCCGCC
GCCTCCAGCGCCGCGCAGCCACCGCCGCCGCCGCTCTCCTTAGTCGCCGCCA
TGACGACCGCGTCCACCTCGCAGGTGCGCCAGAACTACCACCAGGACTCAGAGGCC
GCCATCAACCGCCAGATCAACCTGGAGCTCTACGCCTCCTACGTTTACCTGTCCAT
20 GTCTTACTACTTTGACCGCGATGATGTGGCTTTGAAGAACTTTGCCAAATACTTTC
TTCACCAATCTCATGAGGAGAGGGGAACATGCTGAGAACTGATGAAGCTGCAGAA
CCAACGAGGGTGGCCGAATCTTCCTTCAGGATATCAAGAAACCAGACTGTGATGAC
TGGGAGAGCGGGCTGAATGCAATGGAGTGTGCATTACATTTGGAAAAAATGTGAA
TCAGTCACTACTGGAAGTGCACAACTGGCCACTGACAAAAATGAC

25

Sequence ID - 1188

nt: 599

GGGAGACAAGCCCAGCCTTTCGGCGAGNATACGTCTAACCTGTGCAACAGCCACT
ACATTACTTCAAACCTGAGATCCTTCCTTTTGAGGGAGCAAGTCCTTCCCTTTTATT
TTTTCCAGTCTTCCTCCCTGTGTATTATTCTCATGATTATTATTTTAGTGGGGGC
30 GGGGTGGGAAAGATTACTTTTTCTTTATGTGTTTGACGGGAAACAAAAGTAGGTAA
AATCTACAGTACACCACAAGGGTCACAATACTGTTGTGCGCACATCGCGGTAGGGC
GTGGAAAGGGGCAGGCCANAGCTACCCGCAGAGTTCTCAGAATCATGCTGAGAGAG
CTGGAGGCACCCATGCCATCTCAACCTCTTCCCCGCCCGTTTACAAAGGGGGAGG
CTAAAGCCCAGAGACAGCTTGATCAAAGGCACACAGCAAGTCAGGGTTGGAGCAGT
35 AGCTGGAGGGACCTTGTCTCCAGCTCAGGGCTCTTTCCTCCACACCATTCAGGTC
TTTCTTTCCGAGGCCCTGTCTCAGGGTGAGGTGCTTGAGTCTCCAACGGCAAGGG
AACAAGTACTTCTTGATACCTGGGATACTGTGCCCAGAG

- 245 -

Sequence ID 1189

GGGAGACAAGCCCAGCCTTTCGGCGAGATACGTCTAACCCTGTGCAACAGCCACTA
CATTACTTCAAAC TGAGATCCTTCCTTTTGAGGGAGCAAGTCCTTCCCTTTCATTT
TTTCCAGTCTTCCTCCCTGTGTATTCAATCTCATGATTATTATTTTAGTGGGGGCG
5 GGGTGGGAAAGATTACTTTTTCTTTATGTGTTTGACGGGAAACAAAAGTAA
ATCTACAGTACACCACAAGGGTCACAATACTGTTGTGCGCACATCGCGGTAGGGCG
TGGAAAGGGGCGAGGCCAGAGCTACCCGAGAGTTCTCAGAATCATGCTGAGAGAGC
TGGAGGCACCCATGCCATCTCAACCTCTTCCCCGCCCGTTTACAAAGGGGGAGGC
TAAAGCCCAGAGACAGCTTGATCAAAGGCACACAGCAAGTCAGGGTTGGAGCAGTA
10 GCTGGAGGGACCTTGTCTCCAGCTCAGGGCTCTTTCCTCCACACCATTGAGGTCT
TTCTTTCCGAGGCCCTGTCTCAGGGTGAGGTGCTTGAGTCTCCAACGGCAAGGGA
ACAAGTACTTCTTGATACCTGGGATACTGTGCCAGAGCCTCGAGGAGGT

Sequence ID 1190

GTTTAAATTTGACAAACTAAAGCTNATNACTGCTATAAGAGTAATAACTGCTCATT
15 TTCCATAACTCATTCTTAAAGTTTTAGTAATGTAAAGTTATTTTTTTGTCAGTAAG
TTATAATGATAGAAGCTTACATGTTTTTTCATGCCTCATCTGTTTCCCCTTAAAC
TATAATTATCAGTAAAGTCCTGTGGTATTTTTCAATTTGTAAGAACTAGGCTATA
TATACATTGGGAAAAACAGCCTTCATTTGTCAATGCACTAGTGTTCAAAGGTTTC
20 TGGTAATTGTGTGCTATTGCTTTTTTGTGACTTGCAAAAAAAAAAAAAAAAAAATT
ACTATGACTTGNGGTAGCCCTGCAACCTTCGGAAGTGCTTAGCCCAGTCTGACCAT
ACATTTATATTTANAATGCTTAGGTAAATAAATAATATGCCTAAACCAATGCTAT
AAGATACTATATAATATCTCATAATTTTAAAAATCACTGTTTTGTATAATAATAA
ACAAGGCAGGCAAGCTGTTCTACAATGACTGTTGGTAAGGGTGCTGAGGAAGAAAA
25 ACAACAATCTTGATTGAGGGATAGTGAATAGACAAAAAATGTCCTAATCAATGAA
GCTGTGTGATGATTCTGATTGACAGAGA

Sequence ID 1191

GTGCAAAGTGTTATATCCACTTTCAACAAAGAGAGAAGCTGAAAAGCTAACCCAAT
30 GTTAATTTTGGATCACACACATTCAGTGTAGACTTTAAGATTTTACTTCTGTTGGA
GTAGCTATATTATTTCTAGTTAAAAAACTCTCTATATACATATTTATTTGTTTTTC
TACTTGTTTAAATATTTTTCTCTTCCAATTAGGAACTCAATATGGAATAAAAAATAT
TTAAATGTATTTTACTCAAACGTGTGTGTATATATGTTTGTGTGCATGATAAGGAG
AGTGAGAGCAAGAGTAAGAGAGAGAGAGCAGCATAGATGGAAGCACACATTTAAT
35 GTCTATGAAATGAGAAAACATTAAGGCTAAGATATTTTTCCTTCTGAACTAGCAGA
TTGTATCAATGGCTGGTCACTTAAATTAATCAGTTTGTAAGATATTTAAAAGGTA
TGTCTACCTTCTTGCAATTAATTTGATTATGTTCTAATGGCATGGCAAGAGAAATG

- 246 -

AAAGAAGATAACTAAAAGTTAAAAGTCGTTGCATGTTTTTGTGTCAGCATACCCTT
CTTTCAGGCTACCGAATAACCTTGATTGACATTGGATTAGTAGTAGAATACCTCAT
TGGTAGAGCATATCGCAGCANCTACACTAGAAAACAT

5 Sequence ID 1192

GTCTGGAAGTCCAGACCTCAGGTGATACCCCTGCCTCAGCCTCCCAATGTGCTGGG
ATTACAGCTGTGAAGCCACCGCGCCCGGCTGCTGTGATAGTTGAGATGTAAACCAA
AAATAAAATTCTAAGCCACCCAATCCGACTGAATGGACCCTTCCTGTTGAGCAAGG
ACATTCCAAAGTAAACTGAAAAGACCAGCTTAGGCCATGATGGGAAGGGGAGGTGT
10 CAACATGCCTCATTCTACCTTCCTCCCTCTGGAATCCAGACACAACCTGACCAGCAT
TAACATTAAACAGAGATCTTAAGCTGGGCACGGTGGCTCATGCCTGTAATCCCAG
CACTTTGGGAGGCCAAGGTGGGATCACCTGAGGTCGGAAGTTCAAGACCAGCCTGG
CCGGTATGGTGAAGCCATGTCTCTACTGAAAATGCAAAATTGGCCGGACATTGTGG
TGCA

15

Sequence ID 1193

TNCNTTTTTTTTCCCNCGGGAAAGCGCGCCATTGTGTTGGTCCCCGGGAATTGCGG
GCCGCGTCGACGAGAAATGGCTTGAACCCAGTAGGCAGAGGTTGTAGTGAGCCAG
AATNGGNCACCTGCACNTTANCCNTGGGTGACAAAANTGAAAACCTTGTCTNAAA
20 AAAAAAAAAAAAAAATTTTAANTNAAATNAAAAANCCTTTNCNTTNTTTTTTNAAAN
NGGGGGGGGNNTTTTNGGGNTTNGNNNTGGTAAAAANTNNNTTTTTTTTTTTTA
GGGGCCNANNCCCCNTTTTANAAAAANCCNGNTTTTNAAAAAANTTTTTTNCCNCN
NTTNGGGGGGGGGGNTTTTNANCNNTNTTNGGGGGGGGNNCCCCNTTANNACCNNC
AAANTTTTTANTTTTTTGNNNAANNNCCCCCTTTTTTNNTTTTTTTGNGGGGGGG
25 GGGNNGCCCCCNCCCTTTNGGGGGGGGGGNTTNGNAAAANNACTTTTNAANAANNA
AGGGNNGGGGGNANATNNCCCCCNGGNTTTTTTTTTTAAAAANTNAANNGGGGG
GGGNNNCTNANTNGGGGCNCCCANNGGGGNTTANAANNATTTCTNCCCAAACCC
CCNGNTTTTATNNCCCCCCCCCCCCCNNNNGAANGGGNGGNCCNTTTTTTTTATT
TTTNGGNGGGNAAAAAANTTTNAAAAANNANNAATNTTTTTTCCCCCCCCCCCCNC
30 TTTTNGGNAAANCCNNGGGGGGNTCCTTTTTNAAANNNNCCCCCAAAAAAANTTT
TTTTNTTNTNTTTTTCTCTNGGGGNCNNANTTNTANANTTTTNCNCNAAAAAA
ANGGGNCCCCTTTTTTNCNGGNGGNNCCCCAAANNTTTTTTTNAAAAA
AAAA

35 Sequence ID 1195

GTTCTGTACNTTCGGAGCTACCTGACAGAGCAGAGTCAACCAGGNTCTGCCCAAAG
AGAGTGTTAGGCCTGAGCTTGAGAGCCCTGGAGAGACGTGTGCACAAAATGTGACC

- 247 -

TGAGGCCCTAGTCTAGCAAGAGGACATAGCACCCCTCATCTGGGAATAGGGAAGGCA
CCTTGCAGAAAATATGAGCAATTTGATATTAACTAACATCTTCAATGTGCCATAGA
CCTTCCCACAAAGACTGTCCAATAATAAGAGATGCTTATCTATTTTA

5 Sequence ID - 1196

nt: 412

GTGACGCGGCCGCGGTCGCTGGAGNCGATCAACTCTAGGCTCCAACTCGTTATGA
AAAGTGGGAAGTACGTCCTGGGGTACAAGCAGACTCTGAAGATGATCAGACAAGGC
AAAGCGAAATTGGTCATTCTCGCTAACAACTGCCAGCTTTGAGGAAATCTGAAAT
AGAGTACTATGCTATGTTGGCTAAAACCTGGTGTCCATCACTACAGTGGCAATAATA
10 TTGAACTGGGCACAGCATGCGGAAATACTACAGAGTGTGCACACTGGCTATCATTT
GATCCAGGTGACTCTGACATCATTAGAAGCATGCCAGAACAGACTGGTGAAAAGTA
AACCTTTTACCTACAAAATTTACCTGCAAACCTTAAACCTGCAAAATTTTCCTT
TAATAAAATTTGCTTGTTTT

15 Sequence ID 1197

CCGCCAACATGGGCGCGTTTCGCACCAAACCGTGAAGAAGGCGGCCCGGGTCATC
ATAGAAAAGTACTACACGCGCCTGGGCAACGACTTCCACACGAACAAGCGCGTGTG
CGAGGAGATCGCCATTATCCCCAGCAAAAAGCTCCGCAACAAGATAGCAGGTTATG
TCACGCATCTGATGAAGCGAATTCAGAGAGGGCCAGTAAGAGGTATCTCCATCAAG
20 CTGCAGGAGGAGGAGAGAGAAAGGAGAGACAATTATGTTCTGAGGTCTCAGCCTT
GGATCAGGAGATTATTGAAGTAGATCCTGACACTAAGGAAATGCTGAAGCTTTTGG
ACTTCGGCAGTCTGTCCAACCTTCAGGTCACTCAGCCTACAGTTGGGATGAATTTT
AAAACGCCTCGGGGACCTGTTTGAATTTTTTCTGTAGTGCTGTATTATTTTCAATA
AATCTGGGACAA

25

Sequence ID 1198

CAGAGGTGGGAGGATTGCTTCAGTTCAAGAGTTTGAGACCAGCCTGGGTAACATGG
CGAAACCCTGTCTTTACAAAAAATGCAAACCTTTGCCGCATGTGTTGGGGTGCGCC
TGTAGTCCCAGCTTCTCGGGAGGCTGAGGTGGGGGGACCACCTGAGCCATGGAGGT
30 TGAGGCTGCAGTGAGCCGTGATACCACCACTGTACTCTAGCCTGGGCCATAGAGTG
AGACACCCTGCCTCAGAAATA

Sequence ID - 1199

nt: 439

CCCATCCCCTCGACCGCTCGCGTCGCATTTGGCCGCTCCCTACCGCTCCAAGCCC
35 AGCCCTCAGCCATGGCATGCCCCCTGGATCAGGCCATTGGCCTCCTCGTGGCCATC
TTCCACAAGTACTCCGGCAGGGAGGGTGACAAGCACACCCTGAGCAAGAAGGAGCT
GAAGGAGCTGATCCAGAAGGAGCTCACCATTGGCTCGAAGCTGCAGGATGCTGAAA

- 248 -

TTGCAAGGCTGATGGAAGACTTGGACCGGAACAAGGACCAGGAGGTGAACTTCCAG
GAGTATGTCACCTTCCTGGGGGCCTTGGCTTTGATCTACAATGAAGCCCTCAAGGG
CTGAAAATAAATAGGGAAGATGGAGACACCCTCTGGGGGTCCTCTCTGAGTCAAAT
CCAGTGGTGGGTAATTGTACAATAAATTTTTTTTGGTCAAATTTAA

5

Sequence ID - 1200

nt: 526

CTGGAGACGACGTGCAGAAATGGCACCTCGAAAGGGGAAGGAAAAGAAGGAAGAAC
AGGTCATCAGCCTCGGACCTCAGGTGGCTGAAGGAGAGAATGTATTTGGTGTCTGC
CATATCTTTGCATCCTTCAATGACACTTTTGTCCATGTCACTGATCTTTCTGGCAA
10 GGAAACCATCTGCCGTGTGACTGGTGGGATGAAGGTAAAGGCAGACCGAGATGAAT
CCTCACCATATGCTGCTATGTTGGCTGCCCAGGATGTGGCCCAGAGGTGCAAGGAG
CTGGGTATCACCGCCCTACACATCAAACCTCCGGGCCACAGGAGGAAATAGGACCAA
GACCCCTGGACCTGGGGCCAGTCGGCCCTCANAGCCCTTGCCCGCTCGGGTATGA
AGATCGGGCGGATTGAGGATGTCACCCCATCCCCTCTGACAGCACTCGCAGGAAG
15 GGGGGTCGCCGTGGTGCCTGCTGTGAACAAGATTCTCTCAAATATTTTCTGTAA
TAAATTGCCTTCATGTAAACTG

Sequence ID - 1201

nt: 613

CTTAAGTATGCCCTGACAGGAGNATGAAGTAAAGAAGATTTGCATGCAGCGGTTCA
20 TTAAATCGATGGCAAGGTCCGAACCTGATATAACCTACCCTGCTGGATTTCATGGAT
GTCATCAGCATTGACAAGACGGGAGAGAATTTCCGTCTGATCTATGACACCAAGGG
TCGCTTTGCTGTACATCGTATTACACCTGAGGAGGCCAAGTACAAGTTGTGCAAAG
TGAGAAAGATCTTTGTGGGCACAAAAGGAATCCCTCATCTGGTGAATCATGATGCC
CGCACCATCCGCTACCCCGATCCCCTCATCAAGGTGAATGATACCATTTCAGATTGA
25 TTTAGAGACTGGCAAGATTACTGATTTTCATCAAGTTTCGACACTGGTAACCTGTGTA
TGGTGAATGGAGGTGCTAACCTAGGAAGAATTGGTGTGATCACCACAGAGAGAGG
CACCCCTGGATCTTTTGACGTGGTTCACGTGAAAGATGCCAATGGCAACAGCTTTGC
CACTCGACTTTCCAACATTTTTGTATTGGCAAGGGCAACAAACCATGGATTTCTC
TTCCCCGAGGAAAGGGTATCCGCCTCACCATTGCTGAAGAGAGAGACAAAAGA

30

Sequence ID 1202

GGAATTCGCGGCCGCGTCGACCTCTGCTCGAATTGACAGAAAAGGATTCTGTGAAG
AGTGATGAGATTTCCATCCATGCTGACTTTGAGAATACATGTTCCCGAATTGTGGT
CCCCAAAGCTGCCATTGTGGCCCGCCACACTTACCTTGCCAATGGCCAGACCAAGG
35 TGCTGACTCAGAAGTTGTCATCAGTCAGAGGCAATCATATTATCTCAGGGACATGC
GCATCATGGCGTGGCAAGAGCCTTCGGGTTCAGAAGATCAGGCCTTCTATCCTGGG
CTGCAACATCCTTCGAGTTGAATATTCCTTACTGATCTATGTTAGCGTTCCTGGAT

- 249 -

CCAAGAAGGTCATCCTTGACCTGCCCCCTGGTAATTGGCAGCAGATCAGGTCTAAGC
AGCAGAACATCCAGCATGGCCAGCCGAACCAGCTCTGAGATGAGTTGGGTAGATCT
GAACATCCCTGATACCCAGAGCTCCTCCCTGCTATATGGATGTCATTCTGAAG
ATCACCGATTGGAGAGCCCAACCACTCCTCTGCTAGATGACATGGATGGCTCTCAA
5 GACAGCCCTATCTTTATGTATGCCCCTGAGTTCAAGTTCATGCCACCACCGACTTA
TACTGAGGTGGATCCCTGCATCCTCAACAACAATGTGCAGTGAGCAT

Sequence ID - 1203

nt: 692

10 TGCAGAGGGGTCCATACGGCGTTGTTCTGGATTCCCGTCGTAACCTTAAAGGGAAAC
TTTCACAATGTCCGAGCCCTTGATGTCCTGCAAATGAAGGAGGAGGATGTCCTTA
AGTTCCTTGCAGCAGGAACCCACTTAGGTGGCACCAATCTTGACTTCCAGATGGAA
CAGTACATCTATAAAAGGAAAAGTGATGGCATCTATATCATAAATCTCAAGAGGAC
CTGGGAGAAGCTTCTGCTGGCAGCTCGTGCAATTGTTGCCATTGAAAACCTGCTG
ATGTCAGTGTTATATCCTCCAGGAATACTGGCCAGAGGGCTGTGCTGAAGTTTGCT
15 GCTGCCACTGGAGCCACTCCAATTGCTGGCCGCTTCACTCCTGGAACCTTCACTAA*
CCAGATCCAGGCAGCCTTCCGGGAGCCACGGCTTCTTGTGGTTACTGACCCAGGG
CTGACCACCAGCCTCTCACGGAGGCATCTTATGTTAACCTACCTACCATTGCGCTG
TGTAACACAGATTCTCCTCTGCGCTATGTGGACATTGCCATCCCATGCAACAACAA
GGGAGCTCACTCAGTGGGTTTAATGTGGTGGATGCTGGCTCGGGAAGTTCTGCGCA
20 TGCGTGGCACCATTTCCTGTAACACCCATGGGAGGTCATGCCTGATCTGTACTTC
TACAGAGATCCTGAAGAGAT

Sequence ID 1204

25 TTTTTTTTTTTTTTCTGCGGGAAGCGCGCCATTGTGTTGGTACCCGGGAAATTCTG
CGGCCGCGTCGACACAGGCCCCAGCATCAAGATCTGGGATTAGAGAGGAAAGATC
ATTGTAGATGAACTGAAGCAAGAAGTTATCAGTACCAGCAGCAAGGCAGAACCACC
CCAGTGCACCTCCCTGGCCTGGTCTGCTGATGACACAGGTGTTGGGCNNGNNCNNGG
GGNNGNNGNNGNNGNNGCNGNNGGNNCNGNNGNNGNNGNNGCNGNNGNNGNNTNNGNNGN
GNNCNGNNGNNGNNGNNGNNGNTCNGNNGCNGGGGCCNNGNCGNCGCGGNCGCG
30 NNTNNGGTTNCNNGCNCNNGGCGCGC

Sequence ID 1205

35 CAGACTCTGACCCAGCCTCAGTCCTAACTCCTGGGGCTGGGCTGAGGGGAACAAGC
ATTTGCTGAACTTGAAAAACAAGCAAATCAAAAACAGGAAAAAATTGTACCTG
GTACTTTTTTTTAGAAAAAAGATTAAAAAAGAAAGATAAATTCTTGTTTGGAAA
CTTGAAAAA
AAAAAATTTTAACTC

- 250 -

TNNNNNTNNCNNCNANTAANNCANNTCNANNNNANNNNAATTACTTNNANGTNNNTC
ACN

Sequence ID - 1207

nt: 642

5 ACGAGAAGCCAGATACTAAAGAGAAGAANCCCGAAGCCAAGAAGGTTGATGCTGGT
GGCAAGGTGAAAAAGGGTAACCTCAAAGCTAAAAAGCCCAAGAAGGGGAAGCCCCA
TTGCAGCCGCAACCCTGTCCTTGTGAGAGGAATTGGCAGGTATTCCTGATCTGCCA
TGTATTCCANAAAGGCCATGTACAAGAGGAAGTACTCAGCCGCTAAATCCAAGGTT
GAAAAGAAAAAGAAGGAGAAGGTTCTCGCAACTGTTACAAAACCAGTTGGTGGTGA
10 CAAGAACGGCGGTACCCGGGTGGTTAAACTTCGCAAAATGCCTAGATATTATCCTA
CTGAAGATGTGCCTCGAAAGCTGTTGAGCCACGGCAAAAACCCCTTCAGTCAGCAC
GTGAGAAAACCTGCGAGCCAGCATTACCCCGGGACCATTCTGATCATCCTCACTGG
ACGCCACAGGGGCAAGAGGGTGGTTTTCTGAAGCAGCTGGCTAGTGGCTTATTAC
TTGTGACTGGACCTCTGGTCTCAATCGAGTTCCTCTACGAAGAACACACCAGAAA
15 TTTGTTCATTGCCACTTCAACCAAAATCGATATCAGCAATGTAAAAATCCCAAAACA
TCTTACTGATGCTTACTTCAAAAAGA

Sequence ID 1208

CCCTATACCTTCTGCATAATGAATTANCTAGAAATAACTTTGCAAGGGAGAGCCAA
20 AGCTAAGACCCCCGAAACAGACGAGCTACCTAAGAACAGCTAAAAGAGCACACCC
GTCTATGTAGCAAAATAGTGGGAAGATTTATAGGTAGAGGCGACAAACCTACCGAG
CCTGGTGATAGCTGGTTGTCCAAGATAGAATCTTAGTTCAACTTTAAATTTGCCCA
CAGAACCCCTCTAAATCCCCTTGTAATTTAACTGTTAGTCCAAAGAGGAACAGCTC
TTTGGACACTAGGAAAAACCTTGTAGAGAGAGTAAAAAATTTAACACCCATAGTA
25 GGCCTAAAAGCAGCCACCAATTAAGAAAGCGTTCAAGCTCAACACCCACTACCTAA
AAAATCCCAAACATATAACTGAACTCCTCACACCCAATTGGACCAATCTATCACCC
TATAGAAGAACTAATGTTAGTATAAGTAACATGAAAACATTCTCCTCCGCATAAG

Sequence ID - 1209

nt: 620

30 CTCTCCTGTCAACAGCGGCCAGCCTCCCAACTACGAGAATGCTCAAGGAGGAGCAG
GAAGTGGCTATGCTGGGGGCGCCCCACAACCCTGCTCCCCCGACGTCCACCGTGAT
CCACATCCGCAGCGAGACCTCCGTGCCCCGACCATGTCGTCTGGTCCCTGTTCAACA
CCCTCTTCATGAACACCTGCTGCCTGGGCTTCATAGCATTCGCCTACTCCGTGAAG
TCTAGGGACAGGAAGATGGTTGGCGACGTGACCGGGGCCAGGCCTATGCCTCCAC
35 CGCCAAGTGCCTGAACATCTGGGCCCTGATTTTGGGCATCTTCATGACCATTCCTGC
TCGTTCATCATCCAGTGTTGGTTCGTCCAGGCCAGCGATAGATCAGGAGGCATCAT
TGAGGCCAGGAGCTCTGCCCGTGACCTGTATCCACGTAATCTATCTTCCATTCCCT

- 251 -

CGCCCTGCCCCCAGAGGCCAGGAGCTCTGCCCTTGACCTGTATTCCACTTACTCCA
CCTTCCATTCCCTCGCCCTGTCCCCACAGCCGAGTCCTGCATCAGCCCTTTATCCTC
ACACGCTTTTCTACAATGGCATTCAATAAAGTGTATATGTTTCTGGTGCTGCTGTG
ACTT

5

Sequence ID 1210

TTCGTAATTAGAATACTGTTTGGACTTGCTCAACAAGCACCTTATCTTAACAAAAA
GTAACCTTATAGAAAAGGGAGACATTCAATTTAACTTCAAGCCCATATTATTCTTAAA
AGCTGACTCTTGAAATAGTATTTATTGAGTCATAGTGGAGTCATGGGACTTTTAA
10 GGGCCGGAAGGGACTATTTAGATCATCCAGTCCCACCCTGTCAATTTTATGGAGGAG
GAAACTGAGGCCTAGATAAGATAACCAGTTAGTGGGTCCACTGACCTTTAGGACAG
TAGTCTATCCGTAAGAGACAACATGGAGAAAGAAATACAACGTTTTTATAGTGAAT
TATCATCTTACAAAGAATATTCTTCCCATATCGCACTTTTTAAAAAGTGGGTACCTT
AGTCAAATAGGAGAAAAAACCACTTGAGTAGTTTCATCCTCAGGTTTTAGGTGAGG
15 AAAGTACTCAGATTAAATAACTTTAAGCACACAGAGCCTGAATGATAGTCTTA
TTTGAGCTCATCTGTGCTTTTAAATGTGTACTACGTTAGGTGTTTTCACTTGCAATT
CCTTTAGTCTTATTTGAGCTCATCTGTGCTTTTAAATGTGTACTACGTTAGGTGTTT
TCACTTGCAATTTCTTGTGTTGACGTTGACAATAAATCGTGAAGCTGCCTTATCTAA
GGAAGTCCTAAAGTAAATCATTGGAACACA

20

Sequence ID 1211

CCATTGTGTTGGNACCCGGGAATTGCGCGCCGCGTCGACGGAGTTTTACCTTATTA
CACTTTAATCTCTGGATTTACCCCATCTCATTTCTCTTTTAGGAAAAGTGGTTGTA
TGTGGTGGCTTTGATGGTTCTCATGCCATCAGTTGTGTGGAAATGTATGATCCAAC
25 TAGAAATGAATGGAAGATGATGGGAAATATGACTTCACCAAGGAGCAATGCTGGGA
TTGCAACTGTAGGGAACACCATTTATGCAGTGGGAGGATTGATGGCAATGAATTT
CTGAATACGGTGGAGTCTATAACCTTGAGTCAAATGAATGGAGCCCCATACAAA
GATTTTCCAGTTTTAACAAATTTAAGACCCTCTCAAATAACAGGCTTAGTGATGT
AATTATGGTTAGCAGAGGTACACTTGTGAATAAAGAGGGTGGGTGGGTATAGATGT
30 TGCTAACAGCAACACAAAGCTTTTGCATATTGCATACTATTAAACATGCTGTACAT
ACTTTTTGGGTTTTATTTGGAAAGGAATGCAAAGATGAAGGTCTGTTTTGTGTACTT
TTAAGACTTTGGTTATTTTACTTTTTGGAAAAGAATAAACCAAGAATTGATTGGGC
ACATCATTTCAAGAAG

35

Sequence ID - 1212

nt: 374

AGAGCAGCAGCCATGGCCCTACGCTACCCTATGGCCGTGGGCCTCAACAAGGGCCA

- 252 -

CAAAGTGACCAAGAACGTGAGCAAGCCCAGGCACAGCCGACGCCGCGGGCGTCTGA
CCAAACACACCAAGTTCGTGCGGGACATGATTCTGGGAGGTGTGTGGCTTTGCCCCG
TACGAGCGGCGCGCCATGGAGTTACTGAAGGTCTCCAAGGACAAACGGGCCCTCAA
ATTTATCAAGAAAAGGGTGGGGACGCACATCCGCGCCAAGAGGAAGCGGGAGGAGC
5 TGAGCAACGTACTGGCCGCCATGAGGAAAGCTGCTGCCAAGAAAGACTGAGCCCCT
CCCCTGCCCTCTCCCTGAAATAAAGAACAGCTTGACAG

Sequence ID - 1213

nt: 567

GAATTATTGACTTTGAATTGCATTTTCAGTACCATGAAGTCAAAGTCAGTGGTGTAT
10 TTGCTCATTTGTTTCATTCTTTCTTTTCCACCAACATTACTGCCTGCAGAGCCAGAG
GTGAGTGCAGAAATCCTGTCAATTTCGTCACTTGTGGACAACCTGCAGCTTGCCACA
GCCTACAGTTCACCACTGTGACCTCTGAAAACCTCCTGAACAAAAGGAAGGAGAC
TTGGAAATCCTGAATGGGCTTGGAGACATTAAGGGAGAACTGCCTCCCTGGACCAA
GGCAGAATTCAATAGAACCAGCAAGAAATTTTCTATGAATGGGAAAGCAGGTGGC
15 AGGGGGCAGGGGTGGAAAAGCTTTGTACAGGAATTGTGGAAAAGCTTTTGCATTAT
CTCTAGTCTGAAAGTCACATTTCTCAGTTTCTTCCACTCTCTTCTGTCAACTTGC
TGTGAGTAAATGACATCTGTACCTGTGACACGGGCCAGGGACTATCACCATATGG
CCCCACACATTATCTAGTACCAGCCTGCCTGGGCCATGCCTTTTCCAGTCACTGT
ACCAGCC

20

Sequence ID - 1214

nt: 620

CTCTCCTGTCAACAGCGGCCAGCCTCCCAACTACGAGAATGCTCAAGGAGGAGCAG
GAAGTGGCTATGCTGGGGGCGCCCCACAACCCTGCTCCCCCGACGTCCACCGTGAT
CCACATCCGCAGCGAGACCTCCGTGCCCCGACCATGTCGTCTGGTCCCTGTTCAACA
25 CCTCTTTCATGAACACCTGCTGCCTGGGCTTCATAGCATTTCGCTACTCCGTGAAG
TCTAGGGACAGGAAGATGGTTGGCGACGTGACCGGGGCCAGGCCTATGCCTCCAC
CGCCAAGTGCCTGAACATCTGGGCCCTGATTTTGGGCATCTTCATGACCATTCTGC
TCGTCATCATCCCAGTGTGGTTCGTCCAGGCCAGCGATAGATCAGGAGGCATCAT
TGAGGCCAGGAGCTCTGCCCCGTGACCTGTATCCCACGTAATCTATCTTCCATTCCCT
30 CGCCCTGCCCCCAGAGGCCAGGAGCTCTGCCCTTGACCTGTATTCCACTTACTCCA
CCTTCCATTCCCTCGCCCTGTCCCCACAGCCGAGTCCTGCATCAGCCCTTTATCCTC
ACACGCTTTTCTACAATGGCATTCAATAAAGTGTATATGTTTCTGGTGCTGCTGTG
ACTT

35

Sequence ID 1215

CACAAGATAGAATGGTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAATTTTAA
GTGACAGTGCCATAGTTTGGACAGTACCTTTCAATGATTAATTTTAATAGCCTGTG

- 253 -

AGTCCAAGTAAATGATCACTTTATTTGCTAGGGAGGGAAGTCCTAGGGTGGTTTTCA
GTTTCTCCAGACATACCTAAATTTTACATCAATCCTTTTAAAGAAAATCTGTAT
TTCAAAGAATCTTTCTCTGCAGTAAATCTCGCAGGGGAATTTGCACTATTACACTT
GAAAGTTGTTATTGTTAACCTTTTCGGCAGCTTTTAATAGGAAAGTTAAACGTTTT
5 AAACATGGTAGTACTGGAAATTTTACAAGACTTTTACCTAGCACTTAAATATGTAT
AAATGTACATAAAGACAACTAGTAAGCATGACCTGGGGAAATGGTCAGACCTTGT
ATTGTGTTTTTGGCCTTGAAAGTAGCAAGTGACCAGAATCTGCCATGGCAACAGGC
TTTAAAAAAGACCCTTAAAAAGACACTGTCTCAACTGTGGTGTAGCACCAGCCAG
CTCTCTGTACATTTGCTAGCTTGTAGTTTTCTAAGACTGAGTAAACTTCTTATTTT
10 TAGAAAGTGGAGGTCTGGTTTGTAACTTTCCTTGTACTTAATTGGGTAAAAGT

Sequence ID - 1216 nt: 484
CAACCTTAGCCAAACCATTACCCAAATAAAGTATAGGCGATAGAAATTGAAACCT
GGCGCAATAGATATAGTACCGCAAGGGAAAGATGAAAAATTATAACCAAGCATAAT
15 ATAGCAAGGACTAACCCCTATACCTTCTGCATAATGAATTAAGTAACTAGAAATAACTTT
GCAAGGAGAGCCAAAGCTAAGACCCCCGAAACCAGACGAGCTACCTAAGAACAGCT
AAAAGAGCACACCCGTCTATGTAGCAAAATAGTGGGAAGATTTATAGGTAGAGGCG
ACAAACCTACCGAGCCTGGTGATAGCTGGTTGTCCAAGATAGAATCTTAGTTCAAC
TTTAAATTTGCCACAGAACCCCTCTAAATCCCCTTGTAATTTAACTGTTAGTCCA
20 AAGAGGAACAGCTCTTTGGACACTAGGAAAAAACCTTGTAGAGAGAGTAAAAAATT
TAACACCCATAGTAGGCCTAAAAGCAGCCACCAATT

Sequence ID 1217
GACAGGCGGGGGCCCAGCGGCCGGGTGAAGGCCGGGTGGCTCTGTGAATCAAAGGA
25 GAGTCCCAGAAAACCTGTGACTGTTGAAGAAAATTCATCTGTGAATTTTTATATTC
AAGGAGTCAGTATTTATATTCATCTTTTAACTGGGAAGATTTATATTTTACTTTA
AACTTCTTGATAATAATTTACAATGAATGGACACAGTGATGAAGAAAGTGTTAGA
AACAGTAGTGGAGAATCAAGGTAAGTAAGCACTTTGTTATCAATTGTTTACTATGA
AGAGAGTTGAAAACCTTGACTTTTTTCTTTATTGTTATTGTTGTTATTTAGTTTTCC
30 TCATAGGTAGCAGAGTTTTTCAGGTTTTCTTCTTAGCTATCCAAATACTAAAAAAT
TCTGATATACGAACCTTTTTTCATAATACAGGTTTTAATTATATTTTTCATTCAGA
TACACAGTAGATCTTAAATATAGAAAGTTTTTGTTTACTTAAATCTATTTGGAAGT
TTATATTTGAGCTAATAATTAAGCTGGAGCATGTATAATAGATTTAAATTGTTTTG
ACTGTTAGTGAAATTT

35

Sequence ID 1218
CTCACTTGGTGGGTGAGCCTCCAATGACTACACCCAAGGAGGATTTAACACAGGGA

- 254 -

TTTTATGACTTGCAACAAGTCAGGAGGACATGGGGTTGGGGTAGTTCAGCAGTGCC
TGTCTGAACAAAGGTGAAAATTGGGCTTTTATTGGGCTGATCAAGGGGGAGTAAAG
GCAGCCAGGAGCAGTCGCCTGTATGCTTCTACCTATATTGCATGTATAGAAAAGG
GAAAATAAACTCCTTCCTGGGCAGGGTTTTAGTATGCTAAGGAGGGGAGTTATTCA
5 ACTTCAATCCAACCTCAAGCATCAGCATTGCTGCGTCCATCCCAGTTTTTGTTTTGCT
GGGGCTGAACCTTCTTCCTATAACTTTTTGAAACAACAAGAACTCAAGGTGTGACAG
TTACAAGTGGGCCCTTTTTACAGTGTGTACCTAAACACGTGAGGACCCTGGATTA
CAGAATGACAGACTCGAAGTGACTCAAGTTCGGGTGTTTCATCTTTAGATGGTAAA
GATGGCTGTACGTACTATCCTTGCTTATTTCCAATCTATTGTTTAAACTCTTGAT
10 ATGTAATACCGCAGAGGCTAGAGATACAACCTTTGACCAAATGAGTGAATTCAAGT
AATCCATTACTAATGTGATCTGGAAACAACATGGTGTGAATGTGCATATGT

Sequence ID - 1219

nt: 559

CTTGGCAGCTCCGTATGTGCCAGCTCTTTGCAAGGGCATACTGGGAAATGAGTG
15 GAGATAAAGGACCCAATCATAAGCATTTTACAGTATGGATACCCCATTTTAAAAAG
GTAAACTGAGGCACAATGCAATTTTTTTTTTTTTTTTAAAGGAGTTTATTTGAGCAA
CAGTGATTCATGAATCAGGCAGCACCAAACCAGAAGGAGGCTTTGCTGAANAAGGA
TGAGGGACAAGCATTTATAAAGTGAATGTAGATGTAATACAAAGAAAATATTTGAA
CCGGGTGCGGTGGCTTACACTTGTAATCCCAACACTTTGGGAGGCCAAGGCGGGCA
20 GATCACAAGATCAAGAGATCGAGACCATCCTGGTCAACATGGTGAACCCCATCTN
TACTAAAAAATACAAAAATTANCTGGGCGTGGTGGTGCCTGTAGTCCCAGCT
ACTTGGGCGGCTGAGGCAGGANAATTGCTTGAACCCGGGAGGTGGAGGTTGCAGTA
AGCCGAGATTGCACCATTGCACTACTCCAGCCTGGTGACAGAGAGAGACTCCATC

25 Sequence ID 1220

GANNNGTGCGATANNATGNNTGTCTTTTTTTTAAAGTNTTTCNNATNGNAGNGAAN
CCCCNNANNTNNCATAANGAGAGATNACTACNGTACANATAGNGNCANACNGATA
GTAGTANCAANATTGTNTTAGCTANATNANTCAATAGATATCNAGATANAANAANA
NCNNGGATATACAGCGATGTNTNANNGGNNNNNNANGGAACGAACATCNACNTTA
30 ANNATAAGCTNGNGGAGAGAGACANGTANGTTATANANNAGAATNGNAGTAGGNGT
GATCATAATAGNNNNNANNTANTATATANGATNTTANTGNNCTNTNNTNNGTTTAT
CNNNAATNTCTATNCTNGAGAGNAGCNNNATNNNNAGGCGANGANATTGGGNNNTN
CTCNTNATAGANANCTGGTGTCTNNANAANTACNTCATCTATTNANCTCTCACNANA
TGGNANNATANAGNAGNGNNNTNNANAGGANTANGCATAGNGNNTNNCTNAAACAA
35 AANNNATAAGANNTCTCGNNAANANGGGCCTNTNNTNTAGCGAGGNNTTANTTTNT
ATANTTNTTCNCTCTTNNAATANNTANGATANATGANCTNGNNGTGATANATANN
NNTACNGTNAANNTNTANTCNTATAATAGATANAAATATAGGATNTTNTCTGCGCN

- 255 -

GGTNGAANANTTNNTNCNNTTTNAATAATGNTGTTAGNGACNGNGNTNTNANANN
NNTTAGAAAGGTACTCTATATACTNNTATGNTNCGGCNNATAATANAACAGATGTT
TGTATNAATATNAAANAAGGTCNNTTTCGNCAAGAGAANNNTGNCTGGTNATAGAA
TTAGCATAANTTANNTANTATGATNNANTNNTNCTACNANTNTTAGCNNTTNGCAG
5 NAGTCATTNNGNATNTATNNGNNTANTAGTNANTTGGGNCTNNTNCAGANTATAT
TNTGNGAANATGAANNACGNANTCCTNNGNANTATNATNNTGANTANGANAANCN
ANANNTNTTNTANNANTGNCTATANATTGCCNNGATANATTNTNNAATGAANCGA
TAGCCCGCNCTAAGGANNTNNGTNANNTAAANNTCTCAGATAANNTACNTNTTNT
TATTAANCNANNATCACANTATANCNGNGACANNNGCGANANTATATGTATGNAN
10 TATNACNGNTCCNNNCCGNGAANNANTCNTANNAGGCATTGNGNNGAGCTNTTCT
NCTAGACNATTTNNANTGAAANNATGCNGNNAACGACNNNCTTNAANTTNTGT
CTACANTCCGCNNTNTTNTACAGATNGCAGNTAAGNNNANTNANNGCTCTCANCT
NGCTNNNACT

15 Sequence ID - 1221 nt: 741
AAGCAGAANTNTCTCTAAAAACATTATCTCCTTAAATCTTGAGGTGCATATNAGA
GCCACAGGCAATCTCTGACATATAAAATTGCAGTACAGGCCTTTCAAATTTGGCAT
TTCAGTGGTACAATAACAACAACCAAGATATATAATAACTGTACAGTGCCTAGACAT
TCCAGTAAGAACCATTATTTTCTTTAATGTAGAATGATTAATACATATTCTACAAG
20 GGGCAGTAAGGTAGTAATTCTATAGGGTATGTCCCGACATAATTTCAAATTGTA
CAATAACACAAACAACCTTTGTTAAGGCCATGTTTTATTGCTGATTAATGGACAAA
AGGCAATGTAATTTATTTTCAAGTATTTTCTTGAAAGTCTGTGCTCATAAAAATCA
TGAAAAGTTGGAAAGACTGTTAAATCACTGAACTTCAAATATATCTTACACAATC
TTGTTTGTACAAAAATACAAGTTAAATATAAACATAAAGCAATCATGGTAATTTTA
25 TGCAAATCTGTTTTATGTGATCATCAGTTATATATAAAAGTTTCTCAGTTCTGTTA
TTTGTGAAAAGATCAATACCAGATTGAATGACTACCTATTGGCAAAGGGCCCTAAA
AAGCTTACTTTAGCACTCATCTTTTACATGGTTAAATGCATTTCTTAATTTGAGAT
CACCTAAACACTGGAAAAGAAAAAAATGAAAGGGCAGTATGTCCATAAACCAACA
AATAATTTGGCTG

30 Sequence ID - 1224 nt: 485
CGAAATTTCTTGTGACACAGAGGAAGGGCAAAGGTCTGAGCCCAGAGTTGACGGA
GGGAGTATTTAGGGTTCACCTCAGGGGCTCCCAAAGCGACAAGATCGTTAGGGAG
AGAGGCCAGGGTGGGGACTGGGAATTTAAGGAGAGCTGGGAACGGATCCCTTAGG
35 TTCAGGAAGCTTCTGTGCAAGCTGCGAGGATGGCTTGGGCCGAAGGGTTGCTCTGC
CCGCCGCGCTAGCTGTGAGCTGAGCAAAGCCCTGGGCTCACAGCACCCCAAAGCC
TGTGGCTTCAGTCCTGCGTCTGCACCACACAATCAAAGGATCGTTTTGTTTTGTT

- 256 -

TTTAAAGAAAGGTGAGATTGGCTTGGTTCTTCATGAGCACATTTGATATAGCTCTT
TTTCTGTTTTTCCTTGCTCATTTTCGTTTTGGGGAAGAAATCTGTACTGTATTGGGA
TTGTAAAGAACATCTCTGCACTCAGACAGTTTACAGA

5 Sequence ID 1226

GGTTTTTATACTTGCCATGAACTGTTCTTTGGGATATTATTTTGTTCAGGTTCCC
CACTTGGACAGCAGAGGGGGTGACTCTGCCCATCCCTGCCACTGGTAGCCAGGCGG
GCAATGTCTGCTAGCAGTCTGCTTCTGTCTGAACTCAGCCAGCAGAGGCAAACCTCC
CGGTTCCCCGAGAAACACTCTGAAGGCAGGGTGGGTGACTCCACCCACCACCGCCT
10 CTCCTAGCCATGCAGGCCATGTCTGCTAGAGCTTCCAGCGCAGTGGTCCTAATTCT
GTCTGAATCCGGCTGAGGGGTGCAGCCTCCTGTTACTGCCCAGGGAAACACCCAGA
TGGCAGGGTGGGTGACTCCAACCACCTCTGCCTGTGGTAGCCAGATGGGCCACACC
TGCTAGAGCTTCCAGCCCAGCAGTCCCGCTACTCTGTGGGTGGGTGCCATCCCCTG
TTCCTCTGGGAAGCACCCAGACAGCTGATTACGTGACCCCAACCACTTCTGCAGAT
15 CCTAGCTGAGCAGGACTTGCTGGTTTGGACAATGCCCAAGCAGGGAAGAGCCCTCA
TTCTCTTATCACTGACAGAGGTGAGATGTCCGANTTTGTANGCTGGTGGAGGAGTG
AGGTGGAGGAGGTATGCCTCT

Sequence ID 1228

20 GTTATTCAGGTATCCATCAAAATTTTATAAGAGGGCCGGAAACATCGGCTCACACC
TGTAATCCCAGCACTTTGGGAGGCTGAGGCAGGTGGTTCACTTGAGGTCAGGAGTT
CGAGACCAGCCTGGCCAACATGGCAAAACCCCGTCACTATTAAAAATACAAAACAT
TAGCTGGGTGTAGTGGCAGGTGCCTGTAATCCCAGCTATTTCGGGAGGCCTAGGAAG
GAAAATGGCTTGAACCTGGGGGTGGAGGTGGAGTGAGGCAAGATCACACCACTGC
25 ACTCCAGCCTGGGCGACAGAGCGAGACTCCATCTCAAAGAAGAAAAAAAAAACAA
CAAAAAACCTTTATCAGATTATCAGAGGTTATCACTACAGAGGGAGGTAAAATTG
GAGGGAAAAGGGTACAAATTTATTTAC

Sequence ID - 1230

nt: 741

30 AAGCAGAANTNTCTCTAAAAACATTATCTCCTTAAAATCTTGAGGTGCATATNAGA
GCCACAGGCAATCTCTGACATATAAAATTGCAGTACAGGCCTTTCAAATTTGGCAT
TTCAGTGGTACAATAACAACAACCAAGATATATAATAACTGTACAGTGCCTAGACAT
TCCAGTAAGAACCATTATTTTCTTTAATGTAGAATGATTAATACATATTCTACAAG
GGGCAGTAAGGTTAGTAATTCTATAGGGTATGTCCCGACATAATTTTCAAATTGTA
35 CAATAACACAAACAACTTTGTTAAGGCCATGTTTTATTTGCTGATTAATGGACAAA
AGGCAATGTAATTTATTTTCAAGTATTTTCTTGAAAGTCTGTGCTCATAAAAATCA
TGAAAAGTTGGAAAGACTGTTAAATCACTGAACTTCAAATATATCTTACACAATC

- 257 -

TTGTTTGTACAAAAATACAAGTTAAATATAAACATAAAGCAATCATGGTAATTTTA
TGCAAATCTGTTTTATGTGATCATCAGTTATATATAAAAGTTTCTCAGTTCTGTTA
TTTGTGAAAAGATCAATACCAGATTGAATGACTACCTATTGGCAAAGGGCCCTAAA
AAGCTTACTTTAGCACTCATCTTTTACATGGTTAAATGCATTTCCTAATTTGAGAT
5 CACCTAAACACTGGAAAAGAAAAAAATGAAAGGGCAGTATGTCCATAAACCAACA
AATAATTTGGCTG

Sequence ID - 1231 nt: 203
TTGAGGAAGGGTCTACTGTCTTTTTAAATGGCACAAATTTAAGAGGTTTGAGAGGT
10 ACAGTCCCTTAACCTGCCACGGGAGAGGGGCCCCCAAACTTTCTTCCCCCACACT
TCTGGTTTTCTGTGTGGAGGGGGAGCAGGGATATCTAAGCTGTGGTGTGAAAGGGT
AGGAGAGATGCTGGAGGTGGGGGTGCTGTGTTCTA

Sequence ID 1239
15 TTTCTCGGGAAGCGCGCCATTGTGTTGGTACCCGGAATTCGCGGCCGCGTCGAC
ATTTTTTTTTTTTTTTTTTTTTTTAGAAATGATTAACAATTTATTGAGTTTTATTTATC
TACAAAATATAGCAATACAGNGAACTTCACCAAATCCTAAATATTCAGTACCTGA
ACTGGCTACAACACCGNGTGCACACCCAGTTCTTGCAGAATCTCTTGCAGATATGG
GAGAGTCAGCCAGTGAAAAGATCCATTTCTTGGGAATCCTTGTCAACAAGACCAGT
20 TCAGAAATCCAGGATATATAGAAGCCTACTGTAATTTAAAAACAGTAACAAAAACC
CCAACAAAACCAAATCAACAAAGACCAAGATAAAGNGTGATAAACATTAATTGT
AATGGTTTTCTTTACATGCAATACATGCATTTTAAAATCACTAAGAAACACGAAA
TTTTGTAGAGCAAAGTTTGNGTTTTCACGTAAGTGCAAATGAATATATATTTTATTT
TTTATACTATTAAATTATATATATTTTTTTCCATACAAAAGCACACAGTGTTAATCT
25 ATAAATGACATCCAAGTGGATGATGATTGTTTTTGCATGTCCCCCTGCTTAGATT
TTTTTAAATATATAGTCAAAAATTAACATCCTTCTTTAAAATAACAGAAGGGAAA
AANGGGCAAAAAAAAAAATCTAGACTCGAGCAAGCTTATGCATGCATGCGGCCGCA
ATTCGANCTCGGNCGACTTGGCCAATTGCGCCTATAGNGAGTCGNATTACAATTCA
CTGGGCCGNCGNTTTACAACGTCGNGACTGGGAAAACCCTGGCGTTACCCNNCTNA
30 TCGNCTTGNAACAATNCCNTTTTNGCCAGNGGGG

Sequence ID 1255
TCACTTCGTATNGAANCTGTTTGGACTTGCTCAACAAGACCTTATCTTAACAAAAA
GTAACCTATAGAAAAGGGAGACATTCATTTAACTTCAAGCCCATATTATTCTTAAA
35 AGCTGACTCTTGAAATAGTATTTATTGAGTCATAGTGGAGTCATGGGACTTTTTAA
GGGCCGGAAGGGACTATTTAGATCATCCAGTCCCACCCTGTCATTTTATGGAGGAG
GAAACTGAGGCCTAGATAAGATAACCAGTTAGTGGGTCCACTGACCTTTAGGACAG

- 258 -

TAGTCTATCCGTAAGAGACAACATGGAGAAAGAAATACAACGTTTTTATAGTGAAT
TATCATCTTACAAAGAATATTCTTCCCATATCGCACTTTTAAAAAGTGGGTACCTT
AGTCAAATAGGAGAAAAAACCACCTTGAGTAGTTTCATCCTCAGGTTTTAGGTGAGG
AAACTGATACTCAGATTAAATAACTTTAAGCACACAGAGCCTGAATGATAGTCTTA
5 TTTGAGCTCATCTGTGCTTTTAATGTGTACTACGTTAGGTGTTTTCACTTGCATTT
CCTTTAGTCTTATTTGAGCTCATCTGTGCTTTTAATGTGTACTACGTTAGGTGTTTT
TCACTTGCATTTTCTTGTTTGACGTTGACAATAAATCGTGAAGCTGCCTTATCTAA
GNAGTCCTAAAGTAAATCATTGGAACACATGTANCCAGTTTGTGTTTTTATTGTC
CAGGTNTCAAATATAACTGAAAACCCATGCTAACTGACTNATTTTAAAAGNTGTNT
10 GGGGCATGAAANGATTGCTCTGCCTGGGCGGGNGGTTNANCCTGNGTCCCCCNTTT
NGGAGNCCACCCANGANGCGATATTTNAGGGNNGATTTCNAAACCCCTGGCACGNGN
NAACCCCNNTTTTAAANANAAAAANANCGGNNG

Sequence 1256

15 TTGTGTTGGTACC⁶GGGAATTCGCGGCCGCGTCGACGGAGTTTTACCTTATTACAC
TTTAATCTCTGGATTTACCCCATCTCATTTCTCTTTTAGGAAA⁶CTGTTTGTATGT
GGTGGCTTTGATGGTTCTCATGCCATCAGTTGTGTGGAAATGTATGATCCAACCTAG
AAATGAATGGAAGATGATGGGAAATATGACTTCACCAAGGAGCAATGCTGGGATTG
CAACTGTAGGGAACACCATTTATGCAGTGGGAGGATTTCGATGGCAATGAATTTCTG
20 AATACGGTGGAAGTCTATAACCTTGAGTCAAATGAATGGAGCCCCTATACAAAGAT
TTTCCAGTTTTTAACAAATTTAAGACCCTCTCAAACCTAACAGGCTTAGTGATGTAAT
TATGGTTAGCAGAGGTACACTTGTGAATAAAGAGGGTGGGTGGGTATAGATGTTGC
TAACAGCAACACAAAGCTTTTGCATATTGCATACTATTAAACATGCTGTACATACT
TTTTGGGTTTATTTGGAAAGGAATGCAAAGATGAAGGTCTGTTTTGTGTACTTTTA
25 AGACTTTGGTTATTTTACTTTTTGGAAAAGAATAAACCAAGAATTGATTGGGCACA
TCATTTCAAGAAGTCCCCTCTCCTCCACATTTGTTTTGCCAATTTGCACATTAAAT
GACTCTTCCCTCAAATGTGTACTATGGGGTAAAAGGGGTAGGGNTTAAANATGTAA
ACAGTTGGGTTTTTTAAGGGNCCTTTTTTCATAACTGGAACACTCTNTACAAGGNTN
CTTNTTAAATAAATAACTTGACTTTTTTGTNTTAAANGNANCTTCNTGCTTCCA
30 TAAAAA⁶AAATTTAANTNGNCANCTNTGCTGCTGCGNCCANTTNGCTNGNCCNT
GGCATTCCTTAGGGANGNTNAATANTGGC⁶NNNTTAACNNGGCNGNAACNNNNNCCA
NT

Sequence ID 1331

35 GGGCGATGCATGCTTTATTAAGGCTCTTGTTTCACCTGGCAGTGTACTGTATCAAC
GTATAATACAGAAAAA⁶ATCTCTTTAAGGTCCTCCTTCACAAAGACATAGAGTGA
AACTCCCTTTACATGTCAGTATTTGTTCAACACTTTAGGCAACTTGACTGTCAGTG

- 259 -

TTAAAATGGAAAACAGGAAAATGGAAAATCTGACCAATTCTGCCACCTTGAGACT
TTCATATAGACCTTGCACAACAATTGTATAGATCACACACCGGCTGTATTTAATAT
GTAACATTTTTCACACATATTAAAGATACAGAAGTATTAAAAAACCCCAATGTTAA
TGTATTTGCTTAAAAGGCACAAGTTTCACATATCTGTCTAGCTATCTGTTGGTAA
5 ACAGAAAGTATACTACTTTTTTTTAAAAAGTGGGCAGAATTCTTGTGTATGTATATT
TGTGTGTACAGTATGTGTATGTGTGTATATATATATATTATATATATAGATAATAT
ATAAATATTTTTTTTAAAGGAGAACTAGAATGTTTAGCTAGAAAATTCCACAGCCT
GTGAAGAAATATTTCAAATGGCCATAAAGGAGGTAAAAATGAAAACCATAACCTA
ACTTTTATAGAGGCTTTATCTTTAATTTAACGATGTGCGGAGGACTTTCTTGCTTG
10 AATCTGTTCCGGGCTGTCTGCTCTGTCCATCAAATGGGCAGGTCTGGGAATGAGGC
ACCTTCGGCCGTTT CAGAAGTGGCCTGAACAGAATGCTGGAACCCAGGCTGGACTCG
GAC

Sequence ID 1332

CAAACCTGCATGTTCTGCACATGTATCCAGGAACCTAAAAAAGATAGTT
15 TGTGTGTCTTAATTGAATAATAGTAGATTTATAGATTAAAGATCTATGGGTTTTTA
ATATGGATTAGAAATCTGTGGGTTTTTGATATGGATTAGAAATCTGTGGGTTTTTA
ATATGGATTGAAATCTGTGGGTTTTTAATATGGATTAAAAACATCTGTGGGTTT
TTAATATGGATTAAACATCTGTGGGTTTTTAATATGGATTAAACATCTGGGTTTTT
20 AATATGGATTAAACATCTGTGGGTTTTTAATATGGGTTAAAAATCAAAGAAAATG
AACTATTTGCTCCAGTGCAGGAAAATACAGGCAATACTGGATACAATTAGATGGTC
AGGAGCGATAACCCGGTTGCCATTGTTTGAAGAAGAGAATAAGGTGCTAGCATTC
TATCCGTAGATAATTTGACAGCTAGGAAATAGGGGGAGTCTTCTATGTAGTTAGTG
AAGGCTAAATGAACTATTATATGCAGTTATCGTAGAAGAGTACTCAAAAAATCTG
25 TAAAAATAAAGAAAGGCCGGGCGCGGTGGCTCACGCCTGTAATCCCAGCACTTTG
GGAGGCCGAGGCGGTGGATCATGAGGTCAGGAGATCGAGACCATCCTGGCTACCA
NGGTGAAACCCCGTCT

Sequence ID 1335

CAAGACTCCATCTCAAAAAAAAAAAAAATCTACAGTGCTGAGTATATAAAATTAT
30 TAACACATTT CACAACAATATGTGTTTGTGGAGTTAAATATTTTTTGTCTTTAAAA
CAGGTAATTTTAGTGCATACTTAATTTGATGATTAAATATGGTAGAATTAAGCATT
TTAAATGTTAATGTTTGTACATTGTTCAAGAAATAAGTAGAAATATATTCCTTTG
TTTTTTATTTAAATTTTTGTTCTCTGTAACTAAAAGAACACGAAGTAATTGGTC
35 ACAATTACTGGTGTTTAACTGCCAAATATGGGTAAATAAGGGAAAATTTTGTTTAA
TATTTAGTCCTTCTGAGATGGCTTGAATATTTGAATTTTGTGTACGTCTATACTG
GGTAGTCACAAGTCTTATAAACACTTTAGAGGAAAGATGGATTTTCAGTCTGTATTT

- 260 -

TTAAACATCATTTATTTTAAATCTGGTGCTGAAAAATAAGAAAAAATTAAACTGC
ATTCTGCTGTTCTTCTTTAGAAGCATTCTGCGTAAATACTGCTGTAATACTGTCA
TGCAAAGTGTATCCTTTCTTGTCTGATCCTTTTGGGGCAGTGTTTTTTTGTTTTTT
TTCCTAGAAATGTTTGTCTTCCCCACCTGTTGATCCAGGTTAAGGAATACTTTT
5 TTACACTTTATTCAA

Sequence ID 1336

CTTTTCCTCCCGCTGTCCCCACGGAGGGGACTGCTCTCCCCCGCTGCATCCTTTC
TGTGAGGTACCTTACCCACCTCAGCACCTGAGAGGGTGAAATAGAATTCTAACCTC
10 GACATTTCGGGAAGTGTTTTTGAGAAGTCTCGGTGCGTAAGGGAAGTCTTCCAAGTC
CGTGCAGCACTAACTGATTGGCACCTGCCTCCTCTTCGGCCACCCCCCAGATGAGG
CAGCTGTGACTGTGTCAAGGGAAGCCACGACTCTGACCATAGTCTTCTCTCAGCTT
CCACTGCCGTCTCCACAGGAAACCCAGAAGTTCTGTGAACAAGTCCATGCTGCCAT
CAAGGCATTTATTGCAGTGTACTATTTGCTTCCAAAGGATCAGGCCCTGAGAACAA
15 TGACCTTATTTCTACAACAGTGTCTGGGTGCGTGCCAGCAGATGCCTCAGATAC
CAAGAGATAACAAAGCTGCAGCTCTTTTGATGCTGACCAAGAATGTGGATTTTGTG
AAGGATGCACATGAAGAAATGGAGCAGGCTGTGGAAGAATGTGACCCTTACTCTGG
CCTCTTGAATGATACTGAGGAGAACAACTCTGACAACCAACAATCATGAGGATGATG
TGTTGGGGTTTCCCAGCAATCAGGACTTGTATTGGTCAGAGGACGATCAAGAGCTC
20 ATAATCCCATGCCTTGCGCTGGTGAGAGCATCCAAAGCCTGCCTGAAGAAAA

Sequence ID 1337

CAAGAACTCTGGGACATTTGCAAAGGGTATGGCATATGTGTAATGGGAATACCAGA
GGAGAGGAAAGACAGGAAGTCAAAAAAGAATTTTCCAAATTAATGATAGGTTC
25 AAACCACAGATGCAGGAAGCTTAAACACCAACAGGATAAATAAAACAAAATCTACG
CTTAAGCATATCATACTTAACCTGCAGAAAATTACAGACAAAGAAAAACACCAGA
GGGGAAGCTGGCAGAAACATACCACCTATAGCGGAAGAAGAATAAGAATTACATCA
GACTTCCCTTCAGAAATCTTGCAAACAAAAGATGTAGCACAATATTTAAAGTATT
AAAGGAGGCCGGGCCCCGGTGGCTCGGGCCTGTAATCCTAACACTTTGGGAGGCTGA
30 GGCAGGAGGACCATGAGGTGAGGATCGAGACCATCCTGGTGATGGTGATACCCC
ATCTCTACTAAAAATACAAAAATTAACCGGGCATGGTGACACGCACCTGTAATCC
CAGCTACTTGGGAGGCTGAAGCAGGAGAATCGTTTGAGCCCAGGAGGTGGAGGTTG
CAGTGAGCCGAGATCACATCACTGCACGCCTGGGCAACAGAGCGAGACTCCATCTC
AAAAAA

35

Sequence ID 1338

CGACCCGTTTTAGTCAGGATGGTCTCGATCTCCTGACCTCGTGATCCGCCTGCCTC

- 261 -

GGCCTCCCAAAGTGCTGGGATTACAGGCGTGAGCCACCGCGCCCGGCGTAAATCAG
GTTTTTTTAAATGTTTGCCAAACCTTATCACTGACTTTTATAACAAAATTATTTACT
ATAATCATTAGGGAATATTTAAGTTCTGCTAATACTTAAAATTGCAGAGTGCTAAA
ACCAGCAGTGAGTTTAGAATCAAGCTAAGCTTTATTGTTGCTACTATTTGAGGCAT
5 ATTAGTTGACTGGTGTTCATATGCAAGGCAGTCTACTGGGTGCAACAAGGGTTAGA
AGGATATTTTTTAAA^ΔAACTGACCCTATTCTCAGGATGAAAATAATACACTAGTAAT
AGTCTGCTCTGTTGGTTAACTCCTCGTAAGGAGGTACAATTAAAATGCTGTAGTGT
TGCAAGGGAAGGAGAGGAAGAATCATATTCCTTCACTAGCAGGATCAAGAAAGCTT
TTATAGAAATATACAAAATCTTCACTTCTTGAAGGATTGGTAAAATTTAATAGCCA
10 ACATTGGGCACCTTATTCATTCTCTGAGTAAATATTTATTGCATGCTTATCTTGTAT
CAAGCATTGTGATGAAAGCACAGAATGAAAGAGGAGGGAGAATGTTTAGAGAATA
AGGGCTGAAACACAGATTTTGTAGGGAGCGTAGGGGAGACTGANAAGACAGGTTCA
GGTTAGTAAGGGCGCTCATATTTTGACCCTGAATGTTA^ΔACTATGTGCACATCATGC
TAGCTATTCTAAATCAGGCATTTTCAAATGGAAGCAGGCACTGACATTTT

15

Sequence ID 1344

CGTGAAGGGTCTTTATGTATTAGTATTAGAGTGATCTTTTGATTATTTTCCTCACT
ATAAGGAAATTATTTCTCAGGATGAGCTGCCATAACATTCCACTGTCTGATGGCA
ATTTTAAAGCCTGAAATTGAAGCCCATGGCTAGGCTATGAGAACCCTAGTTCGTAT
20 AGTAAAGTTGATATCTTCTGGATGTATACTAATTTTAGGCTTTATTTTAA^ΔAACTGC
TGGA^ΔAACTGAACTTAGACAAAAGTATTTTCAGGACATCATTTACAATGTTTAGCC
CTAAAGAGTCAAGCTGTGGGATTCTGAGTCTTTCATATGTTACAGCAGAAACTTAA
AAGCAAGAGGAAATTGGCTGGGCACAGTGGCTCTGTAATCCCAGCACTTTGGGAGG
CTGAGGTGGGTGGATCATGAGGTCAAGAGATTGAGACCATCCTAGCCAACATGGTG
25 AAACCCCATCTCTACTAAAAATACAAAATTAGCTGGGCGTGGTGGCACACGCCTG
TAATCCCAGCTAGTCAGGAGGCTGAGGCAGGAGAATATCTTGA^ΔACTTGGGAGGCAG
AGGTTGCAGTGAGCCAAGATTACATCACTGCACTCCAGCCTGGTGACAGAGCGAGA
CTCCGACT

30 Sequence ID 1348

CTGAAACTGCACTGAACCCACAGGTAGGTTACATCACAGGACAGAAATCTGAGGAG
CTGGAGAAAGCAAAGAATAAAGGATGGGCTGACACCAGAAGGAATTAAAGGAATT
TTTATACTGAACTTCAATTACTTGTTCATTTGAAGTTTGT^ΔTTTTTTAATGAACGTT
TTTGCTGT^ΔTACTTAAATATAGTGT^ΔTTTGAAAGTGT^ΔTTCAAATGTATTCAAGTTGGG
35 ATTTTCCATATTTTACTACAGTTCTGTCTTAGTATGTTCAACATAAAACACTTATC
ATTAAAGCTCACAAAGTGCTTTTTTGTAAATATGAGGATAAAATGAAGCCATATAAG
AATTTTTTTATATCTGTACATTTA^ΔCCACATTTGAGCTTTAGCCAAAATATATAG

- 262 -

CTTTTTTTTTTCTGACCTGGCCAACGTATTATCCAGCAAACATCAACTGAAGCAAT
ATGGAAACACTTCCAAATGTTTGCCAATAATGCTATTAAGTGACTGATGTCAACAT
TAGTTACATGGCAAACATAAGAGGCATTATACATTTTAAAACACACTAACATATA
ACTGTAGATAATGTAAGGTTTATTTATATGCATATTCATAGTATATTTAAATGTT
5 TAAATATAAAAAAGGGTTTTTAAACACTTTTAATTTTATCTTTGATTTTTTTTAT
TGATATCTCTTTCCAGGCTACTAATAAAATTGCCAGAACTAACTATCAGGTAAAG
GTTAAGGCATCAATTGACAAGTAAGTTTTCTAATTTTCGTTTTGAATTACAATTCCA
AATGTAAGACTTTTAAAAATGAATGGCCTTTATTTTATAGAATAATTTTGACCTTT
TAAATTTACTTATCTAACATTATATAATGAATGTACTTCAAATATTTGACTTTGAA
10 GTCAACATTAACAAATTCATGGATCCTAATTAATAATTTACTATAAACTGGAATCA
TTTATTACTTCCTT'

Sequence ID 1351

TTTTTTTTTTTTTAAAGAGATGGGTCTCACTATGTTGCCCATATGTTTATGAG
15 ATTAAGTTCATCTTTTTTATCTGAGTAGTATTTTATTGTATGAATATACCACCAT
TATTTATCTGTTGGTTATTTCCAGTTTTTGGGCTATAATCCAAAATGCTTTTTTCAA
ACAATAGGCTATATATCATTAAATGTCCGTTTATCAGCAGTATAAAATATCTTACCA
TAAATATTAATAAAGAAGCATTATATATAAAATATAGATATTTCAAACCCTACA
GAGGGCCTTTTAATGATTAAATATTTTGTCTTACAAAAAGGTCCAGGTAATTACA
20 CCCATGAGGTTAACCTGCCTTAGTGCAGGACTTAAAATAAGGCTTCTCCTGCCATC
TCTCTCCATTTGTAGAATGTGAAATCTTTAAAATGCATCCTATATTAGGAATACT
ATAGCTGTGCACTGGTGTGTTCTCTCTTTAACTCGGGACCGTATATATCTGC
TCAAATTGCCCAAGTATACATATGCTGCACTCCATCAAGTGTGAGGCCACATTCTA
TCAGCACAGCGTGACTGCCTATCAGTGACAATATAAGTGAGCTCTATTTGGATCCC
25 TCTTACCCTACCTTTTATATTTATGACAGCATTATCATAAACTCCAATATTCTTC
AATAACTTACATGTTTGTGTTAGGATAAAATTATTACCCTCAATGAACTACAT

Sequence ID 1352

ACCAGCTTCTTCACAGGTTCCACGAGTCATGTCAACACAGCGTGTGCTAACACAT
30 CAACACAGACAATGGGTCCACGTCCTGCAGCTGCAGCCGCTGCAGCTACTCCTGCT
GTCCGCACCGTTCCACAGTATAAATATGCTGCAGGAGTTCGCAATCCTCAGCAACA
TCTTAATGCACAGCCACAAGTTACAATGCAACAGCCTGCTGTTTATGTACAAGGTC
AGGAACCTTTGACTGCTTCCATGTTGGCATCTGCCCCCTCCTCAAGAGCAAAAGCAA
ATGTTGGGTGAACGGCTGTTTCCTCTTATTCAAGCCATGCACCTACTCTTGCTGG
35 TAAATCACTGGCATGTTGTTGGAGATTGATAATTCAGAACTTCTTCATATGCTCG
AGTCTCCAGAGTCACTCCGTTCTAAGGTTGATGAAGCTGTAGCTGTACTACAAGCC
CACCAAGCTAAAGAGGCTGCCCAGAAAGCAGTTAACAGTGCCACCGGTGTTCCAAC

- 263 -

TGTTTAAAATTGATCAGGGACCATGAAAAGAACTTGTGCTTCACCGAAGAAAAAT
ATCTAAACATCGAAAACTTAAATATTATGAAAAAAACATTGCAAAATATAAAA
TAAATAAAAAAAGGAAAGGAAACTTTGAACCTTATGTACCGAGCAAATGCCAGGTC
TAGCAAACATAATGCTAGTCCTAGATTACTTATTGATTTAAAA

5

Sequence ID 1353

ACATTCTGAAAAGGCAAAAGGGAGGAAGAACTGATTAGTGGTTAGCCCAGGGTTA
GAGTTGGGGAGAGGATATAATGAGGGAACTTTTGTGGATTCTGTACCATGATTATG
ATTACACAAACCTATGCATACATTGAAACACATAGAACTATACGTTGAAAAAAGTG
10 AATCTGCCTGTATGTAAATTTAAAAGAAAAATTTTTTTTAAAAAACAGATGCTT
CTTAACACATTATCATCTATGTCAGTTTAAACAGTTAGTAGACTTAGGCCAGGTGTC
ATGGCTCACTCCTGTAATCCCAGTGCTTTGGGAGTCTGAGGTGGGACGATCTCTTG
AGACTAGGAGGGAGTTTGAGACAAACCTAGGCAATGTAATGAGACTCTTTCTCTAC
AAAAAATTTTAAAGTTATCTGGACATGGTGGTGCCTGCCTGTAGTCCCAGCTACTT
15 GGGAGGCTGAGGTGGGAGGATTCCTTGAGCCCAGAAGTTCAAGGCTACAGTGTGCT
ATGATAGAGCCACTGCACTCCAGCCTGGGCAACCAGGTGAGACCTTGTCTCTAAAA
TGAATAAATAAAT

Sequence ID 1355

20 TGGTCTTTTACCCAGCCAGGGAGAAGGTCTTCGCTCAGTATGAAGAAAAGCAACC
CAAACTCTCAATCTGATTTGTTTTTGTATTATGTCGATGCCCTGTAGTTTGAAAGT
GAAGTAAAGATTTAGAATTCACCTAAGTCCAAAGGAAAACACGTGGTTTTTTAAAGC
CATTAGGTAAAAAAGTTCTCAATAAAGGCATTACAATTTTTTTAGGTTTAGAAAGA
TGGACTTTTCTGATAAATCTTGGCAGACATCTAAAAAACCATATTTTTTCACA
25 AGAAAATGCAAGTTACTTTTTTTGGAAATAATACTCACTGATTATGGATAAAATGG
AATATTTTCAGATACTATATTGGCTGTTTCAAATAGTACTATTCTTTAAACTTGT
AATTTTTGCTAAGTTATTTGTCTTTGTTGTATCTATAAATATGTAAAAATATTTA
AATAGATGTACCTGTTTTGCTTTCACACTTAATAAAAAATTTTTTTTTGT

30

Sequence ID 1359

CGGGATCCCTAGTATAACACATTCAGTGTTCCCCTTTCAGTCTTACTACTTTGACC
GCGATGATGTGGCTTTGAAGAACTTTGCCAAATACTTTCTTCACCAATCTCATGAG
GAGAGGGAACATGCTGAGAACTGATGAAGCTGCAGAACCAACGAGGTGGCCGAAT
CTTCCTTCAGGATATCAAGAAACCAGACTGTGATGACTGGGAGAGCGGGCTGAATG
35 CAATGGAGTGTGCATTACATTGGAATAATGTGAATCAGTCACTACTGGAACCTGC
ACAACTGGCCACTGACAAAAATGACCCCATGTGAGTATTGGAACCCAGGAAAT
AAATGGAGGAAATCATTTGCCTTAGGGATTGGGAAAGCTGCCCACTAACTGTCTTC

- 264 -

CCCATTTGTTTTGCAGTTGTGTGACTTCATTGAGACACATTACCTGAATGAGCAGGT
GAAAGCCATCAAAGAATTGGGTGACCACGTGACCAACTTGCGCAAGATGGGAGCGC
CCGAATCTGGCTTGGCGGAATATCTCTTTGACAAGCACACCCTGGGAGACAGTGAT
AATGAAAGCTAAGCCTCGGGCTAATTTCCCCATAGCCGTGGGGTGACTTCCCTGGT
5 CACCAAGGCAGTGCATGCATGTTGGGGTTTCCTTTACCTTTTCTATAAGTTGTACC
AAACATCCACTTAAGTTCTTTGATTTGTACCATTCTTCAAATAAAGAAATTTGG
TACC

Sequence ID 1360

10

15

TGCGCAGACCAGACTTCGCTCGTACTCGTGCGCCTCGCTTCGCTTTTCCTCCGCAA
CCATGTCTGACAAACCCGATATGGCTGAGATCGAGAAATTCGATAAGTCGAAACTG
AAGAAGACAGAGACGCAAGAGAAAAATCCACTGCCTTCCAAAGAAACGATTGAACA
GGAGAAGCAAGCAGGCGAATCGTAATGAGGCGTGCGCCGCAATATGCACTGTACA
TTCCACAAGCATTGCCTTCTTATTTTACTTCTTTTAGCTGTTTAACTTTGTAAGAT
GCAAAGAGGTTGGATCAAGTTTAAATGACTGTGCTGCCCCTTTCACATCAAAGAAC
TACTGACAACGAAGGCCGCGCCTGCCTTTCCCATCTGTCTATCTATCTGGCTGGCA
GGGAAGGAAAGAACTTGCATGTTGGTGAAGGAAGAAGTGGGGTGGAAGAAGTGGGG
TGGGACGACAGTGAAAT

20

Sequence ID 1361

25

30

35

TATAAATACACTCCGGGATGATTTACCCCCGGAGGTCAGCTAGTAAAATACATGAG
TAGAATTCCTTAAAGTATGTGATAATTGCTCATCACTATCCAAGTGTGACATAAAT
CATAAAAAGAATTGACAAAATCAGGGTCGCAAAGAGAATTGAAAAAATCTGTCAC
AACCAAAATTTAAATTGACCTCTGTCCTAGAGTATGAGAGCCACACTGAACAGAAA
AACCAGATAAATCTTTTATAAAATATTCATTTGCAGCCCCATTAACGTTGCTTGTC
ACCCACCTCCCCATGTCCTTGGACAAACTGAATGTATAGTAACATCATCCCAGGC
CAGGCGCGGTGGCTCATGCCTGTAATCCCAGCACTTTGTGAGGCTAAGGCAGGCAG
ATCAGGAGGTCAGGAGTTCAGGACCAGCCTGGCCAAAAAGGTGAAACTCCGTCTCT
ACTAACAATACAAAATTAGCTGGGTGCGGTAGTAGGCGCCTGTAATCCCAGCTAC
TCGGGAGGCTGAGGCAGGAGAATTGCTCAAACCCGGAAGGTGGAGGTTGCAGTGAG
CTGAGATCGTGCCACTGCACTCCAGCCTGGGTGACAGAGCAAGACTCTGTCTCGGG
GAGGGGGGTGGCGGAGATAAAGAAATAACATCATCTTATACTGTCAAGCTCAAGGT
GTCTGCAGCCTTATCTTCAGGGGAAGTTGTGTCTTTCTCAGGGAAGATAAGATTT
CAATTTAGAGCAAGACAGAGAGAAGTTACATTCAGAGAGGAAAATGCAGTAGTCTA
ACTG

- 265 -

Sequence ID 1364

5 GCGGCCGCGCTCTTTTCAATTTTAAAAAGAAGTTTGTTTTCCATTTTCAGTAATTT
CTGCTTTGATCTTCCTTATGTCCTCCTATTGAGTTGATCAGCTTTCTTTATTCTTG
CCTTTTCTCCTCTGTGTGCCCTTTCTATTAACGTATTTACCCTTAGGCTGGGCACA
ATGGCTGATGCCTGTAATCCCTGCACTTTGGGAGGCCGAGGCAGGTGGATCACCTA
AGGTCAGGAGTTCAAGACCAGCCTGGCCAACATGGTGAAACCTGGTCTCTACTAAA
AACACAAAAATTAGCCAGGCATGGTGGTGTGCACCTGTAATCCCAGCTACTCAGGA
GGCTGAGGCAGGAGAATTGCTTGAACCTGGGAGGCGGAGATTGTGCCAAAGCACTC
CAGCCTGGGCAACAAAATGAGACTTTGTGTC

10

Sequence ID 1365

15 CACCAGGCTGTCTTCAGATACTTCATACAGAAATGAGCCTCCCTGTGGGGTCCCTCT
TCCCTCCTTCAGCCTGTCCATCAACACAGCATTGCGGGATCCTTACCATGGCATCC
AGCCCTGGAGATGCTTCAGGAAAGTTGCAGGTCCATGCTGCAGGACAGGCTCAGAT
CAGCAGAGACGCATCTCACATCGGGCTGTGAAATTCAAGTTGAGCTGCAATTGGCA
ATGAGAA

Sequence ID 1366

20 GTTATTTCACTGAGACCGTGCCCCGGTTATGAGGTTGTACCAGAAAGCAAGTATTCA
CTATGCACACTATTTACCGCTCACCTTAGCATTGAAGCCAGCCTGTAGCCTGAAAG
CCTTTGCTTTGAGGGCAGGTCTTTCCCCAAAATGCAGACACGAAGGTGCAAAGTGA
AGCTGCCAGTCTTGCAAAAGATGTAAC TTGTACGAAGGCCACGAGTGGCAGGGAG
AGCTGTCCACATTTGCGGAAGTGGCTATGTGAGGACGGGGGAGGCGGGTCCCTTA
GAGATGAGACAATCATAAGGGGAGATATCAGAGAAAATCGTAAGGGGAGCAGATGG
25 TTGTCAAGAGAATAGGCTGACCATCGAAGGACTGGCAGAAGCTTTAGAAAACCAC
TGGACGGCTGGGCACAGTGGCTTAGGCCTGTAATCCCAGCACTTTGGGAGGCTGAC
GCAGGTGAATCACTTGAGGTCAGGAGTTCCAGACCAGCCTGGCCAACATGGTGAAA
CCCCATCTCTACAGAAAATATAAAAATTAGCCAGGCGTGGTGGCACAAGCCTAGAA
TCCCAGCTACTTGGGAGGCTGAGGCAGGCGAATGGCTTGAACCCAGGAGTCAGAGG
30 CTGCAGTGAGTCGAGATTGTTCCACTGCACTCCAGCCTGGGTGACAGTGCAAGACT
CCTTCCAAAAA

Sequence ID 1367

35 TTCGTGAGTGATGGCGTCCCGGGTTGCTTGCCGGTGCTGGCCGCCGCCGGGAGAGC
CCGGGGCAGAGCAGAGGTGCTCATCAGCACTGTAGGCCCGGAAGATTGTGTGGTCC
CGTTCTTGACCCGGCCTAAGGTCCCTGTCTTGACAGCTGGATAGCGGCAACTACCTC
TTCTCCACTAGTGCAATCTGCCGATATTTTTTTTTTGTATCTGGCTGGGAGCAAGA

- 266 -

TGACCTCACTAACCAGTGGCTGGAATGGGAAGCGACAGAGCTGCAGCCAGCTTTGT
CTGCTGCCCTGTACTATTTAGTGGTCCAAGGCAAGAAGGGGGAAGATGTTCTTGGT
TCAGTGCGGAGAGCCCTGACTCACATTGACCACAGCTTGAGTCGTCAGAACTGTCC
TTTCCTGGCTGGGGAGACAGAATCTCTAGCCGACATTGTTTTGTGGGGAGCCCTAT
5 ACCCATTACTGCAAGATCCCGCCTACCTCCCTGAGGAGCTGAGTGCCCTGCACAGC
TGGTTCAGACACTGAGTACCCAGGAACCATGTCAGCGAGCTGCAGAGACTGTACT
GAAACAGCAAGGTGTCCTGGCTCTCCGGCCTTACCTCCAAAAGCAGCCCCAGCCCA
GCCCCGCTGAGGGGAAGGGCTGTACCAATGAGCCTGAGGAGGAGGAGCTGGCTACC
CTATCTGAGGAGGAGATTGCTATGGCTGTTACTGCTTGGGAGAANGGCCTAGAAAG
10 TTTTGCCCCCGCTGCGGCCCCAGCANAATCCAGTGTTGCCTGTGGCTGGAGAAAGG
AATGTGCTCATCACAGTGCCCTCCNTTACGTCAACAATGTCCCCCACCTTGGGAA
CATCATTGGTTGTGTGCTCAGTGCCCGATGTCTT

Sequence ID 1368

15 CAGTGAGCCAAGATCACACCACTGCACTCCAGCCTGGACAACAGAACGAGACTCCA
TATCAAAAAAATTAAATTAAATATAATAAATTTCTTGCCGGGCGCAGTGGCTCAC
ACCTGTAATCCCAGCACTTTGGGAGGCCGAGGTGGGCGGATCACGAAGTCAGGAGA
TTGAGACCATCCTGGCTAATACAGTGAAACCCCGTCTCTACTATAAATACAAAAA
TTAGCTGGGCATGGTGGCGGGCGTCTGTAGTCCCAGCTACTCAGGAGTCTGAGGCA
20 GGAGAATGGTGTGAACCCGGGAGGCGGAGCTTGCAGTGAGCCGAGATCGTGCCACT
GCAATCCAGCCTGGGCAGCAGAACGAGACTCCATCTCAAATAAATAAATAAATAA
ATGAATTTTCAGCTAGAAGAGCCTTATTCCATTTTCCTTTTTATTAAACATCTGGCA
TAAGTTGGTAAGTATGTGAAGTTTATCATATATTCTTATGCGAATTATTATTTTCG
CCTTTTTTTTTTATAATTCTGTCTGGGATTTGAATAGTAGAGTTTGAATTCAGGAAG
25 GACACCTGTGATAGGACAATAAAAT

Sequence ID 1369

CTGATTGCAAAACATTACAACCTCAGTACTGCGGCTTTCATTCAAATAGGTTATAT
GTATAAACTGAGGTTCAACAATATTGTATTTGAGATGGGAAAGTTAAAGAAATGCA
30 ATAATGTAAATAATACTTAAGAAAATAAGATCTCAGGAACTGTGTATACTCTGTA
CTTTTATGCAACTTTATCAGATCATTTTCAGTATATGCATCAAGGATATAGTGTATA
TGACATGAACTTTGAGTGCAAAACTGTACTATGTACCTTTTGTTTATTTTGCTGT
CAACATCTAAATAAAGGTTTTTTTG

35 Sequence ID 1370

CGAAAGGACTACAGAGCCCCGAATTAATACCAATAGAAGGGCAATGCTTTTAGATT
AAAATGAAGGTGACTTAAACAGCTTAAAGTTTAGTTTAAAGTTGTAGGTGATTAA

- 267 -

AATAATTTGAAGGCGATCTTTTAAAAAGAGATTAAACCGAAGGTGATTAAAAGACC
TTGAAATCCATGACGCAGGGAGAATTGCGTCATTTAAAGCCTAGTTAACGCATTTA
CTAAACGCAGACGAAAATGGAAAGATTAATTGGGAGTGGTAGGATGAAACAATTTG
GAGAAGATAGAAGTTTGAAGTGGAAAACCTGGAAGACAGAAGTACGGGAAGGCGAAG
5 AAAAGAATAGATAAGATAGGGAAATTAGAAGATAAAAAACATACTTTTAGAAGAAAA
AAGATAAATTTAAACCTGAAAAGTAGGAAG

Sequence ID 1371

GTCCAGNAGAAAGTTCAGTGACTTGTCAGAGCTGCAGGTCTTAAGAGGCTGAAAT
10 CTCGCCTCTGCCTCGAGGCTGCGGTTCCTGACCCATACTACTTGCCTTCAGGAA
AGAGAAATGGTGTAGGAAGGCTGTGGATGAAGACGCTTACATTCATGAAGGATTTG
GATAGGCGAACATGAGCTTTTCCACCAAATTTTCAGAATTTTAAGAAATGCCTTAAA
TTATTTCTTAAAAATCAATTTGGGGCAGACGAGAAGTTCTGATAATAGTTTTTAGG
GAACATGATAAAATTCTGACCTTAGAAGTGGTATACCAGTTTGAGAAGAAGAACAA
15 GCTATAAACGGTGTAGATAACATTCACGGCTATTTAAGAAAGAGTTACTAAGGGAA
ACCAGAATGACTTAAGAGTGTTACTCTTCTTTTCTGAGAGAACAATAGCATCATC
TCAGAAAGCCTTTTCATGCCATTAATAGGTAAGAATCTGGGCTTCTTGGACCATGGG
TTAGACTTTTCTTACAAAACCATAATATGCATTTCTTAGCAAAATTTATGCTATTAC
ATTTCTTTATCTCAACAAAGACTGGTAAATTCAGTACTTATTCCTCAATTTTCTTA
20 CCCTTAAATGGGGATATTCTGCCTCTCCAAGGAATGCTGGGAACAAGCAAGTCCT
CATGTTAGGGGTCTTTGAGTTTTTCATGGAAGTTTAGGTTATTTATATGATGACATA
GTTGTCAACTTACTTTTCAGGATGGACTTTTCTTTTGTGAGTTTGTGACCTAAATAC
AATAGTTGTTATGCATGTCCAGTTTATGGAAGTACCACTGCAATANCAG

25 Sequence ID 1372

CAGTGCAGCCAAGTATCACACCACTGCACTCCAGTCCTGGACAACAGAAACGANTA
CTCCATATCAAAAAAATTAAATTAAANGATAATAAATTTCTTGCCGGGCGCAGTGG
CTCACACCTGTAATCCCAGCACTTTGGGAGGCCGAGGTGGGCGGATCACGAAGTCA
GGAGATTGAGACCATCCTGGCTAATACAGTGAAATCCCCGTCTCTACTATAAATAC
30 AAAAAATTAGCTGGGCATGGTGGCGGGCGTCTGTAGTCCCAGCTACTCAGGAGTCT
GAGGCAGGAGAATGGTGTGAACCCGGGAGGCGGAGCTTGCAGTGAGCCGAGATCGT
GCCACTGCAATCCAGCCTGGGCAGCAGAACGAGACTCCATCTCAAATAAATAAATA
AATAAAATGAATTTTCAGCTAGAAGAGCCTTATTCCATTTTCTTTTATTAAACAT
CTGGCATAAGTTGGTAAGTATGTGAAGTTTATCATATATTCTTATGCGAATTATTA
35 TTTTCGCCTTTTTTTTTTATAATTCTGTCTGGGATTTGAATAGTAGAGTTTGAATTC
AGGAAGGACACCTGTGATAGGACAATAAAATCTA

- 268 -

Sequence ID 1374

GAAAGCACATATGATATACATGTGTGTCATATGTATTATTTTGTGGCCATCTGAG
TCTTCAAAATTTGTTACAGAATACCTGCATATTAATATTTCAAGGTATGGATTAAT

5 Sequence ID 1378

CTGAGTATTAATAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Sequence ID 1380

CCAAACCCAACTGGTCCAGTAGGATACTCACCTTACAGGGGGCGTCTCAAGAGTCT
10 CACAGTTCCCTTGGGTCTTAAGAGACTCACTGTTGGACCAGGCGTGGTGACTCACG
CCTGTAAAACCAAGCACTTTGGGAGGCCGAGGCGGGCGGATCAGTTGAGGTCAAGAG
TTCAAGACCAGCCTGACCAAGGTGCTGAAACCCCGTCTCTACTAAAAATACAAAAA
TTAGCCAGGCATGGTGGTGTGCGCCTGTAATCCAGCTACTCCAGAGGCTGAGGCA
GGAGAATCTCTTGAACCCAGGAGGTGGAGGTTGCAGTGAGTCGAGATCATGCCACT
15 GCACTCCAGCCTGGGTGACAGAGCGAGACTCCGTCTTAGAAAAAAAAAAAAAAAAAA
AAAAGAACCTCACAGTTCAGCAGGGTTCTAGCATGAGACAATGAGGACAAGGGTAG
GTGAGCAGGTGGAAAGAGTGAGAACAGGTCAATTGTGATGGAGAAAATAATAAAGA
CAGAAAAGGCAGAAGACTGCCTGGCAGAAGACCTGTCCCAGCAGATACAAAATAC
AGACAACAGGAGCCAGCATAGACCCTTGACCTGTGTAAGTCTTTCTCAGGCCTTCT
20 TTAAAGTAGAAACATGCCTTTGAAAAAAAAAGTTTTAATAAACAGGAAATCATAAAT
CCCTATTTACATAAATAATATATCCTGGTCTTATTCTTAAACCATTTGATTTTTCA
CGGCTCATTAANAAAGCTGGGCGAGGTGGCTCACGCCCCGTATCCTAGCACTTTGG
GAGGCCGAGGCGGGCANATCACAAGGTGAGGAGTTGGGAGACCAGCCTGACCAACA
CGGTGAAACCCAGTCTCTACTAAAAATACAAAATTANCTGGGGGTGGTGGTGTGT
25 GCCTGTAATCCAAGCTACTCGGGAGGCTGAGGCAGGA

Sequence ID 1382

CTTACTACCTCCAACATGAAACAAGCAGCCCCGCACTTCTCGAAGGTCTGAGTTAC
TTGGAATCGTTTTACCACATGATGGACAGAAGGAATATTTAGATATCTCTGAAAA
30 CCTCAAGCGTTACCTTCTTCAGTATTTTAAGCCAGTGATTGACAGGCAAAGCTGGA
GTGACAAGGGCTCAGTCTGGGACAGGATGCTCCGCTCGGCTCTCTTGAAGCTGGCC
TGTGACCTGAACCATGCTCCTTGCATCCAGAAAGCTGCTGAACTCTTCTCCAGTG
GATGGAATCCAGTGGAATAATTAATATACCAACAGATGTTTTAAAGATTGTGTATT
CTGTGGGTGCTCAGACAACAGCAGGATGGAATTACCTTTTAGAGCAATATGAACTG
35 TCAATGTCAAGTGCTGAACAAAACAAAATTCTGTATGCTTTGTCAACGAGCAAGCA
TCAGGAAAAGTTACTGAAGTTAATTGAACTAGGAATGGAAGGAAAGGTTATCAAGA
CACAGAACTTGGCAGCTCTCCTTCATGCGATTGCCAGACGTCCAAAGGGGCAGCAA

- 269 -

CTAGCATGGGATTTTGTAAAGAGAAAATTGGACCCATCTTCTGAAAAAATTTGACTT
GGGCTCATATGACATAAGGATGATCATCTCTGGCACAACAGCTCACTTTTCTTCCA
AGGATAAGTTGCAAGAGGTGAAACTATTTTTTTGAATCTCTTGAGGCTCAAGGATCA
CATCTGGATATTTTTCAAACGTCTGGAAACGATAACCAAAAATATAAAATGGCT
5 GGAGAAGAATCTTCCGACTCTGAGGACTTGGCTAATGGTTAATACTTAAATGGTCA
ATAGAAAAAGTAGGCTGGGCGCGGTGGCTCACGCCTGTAATCCCAGCACTTTGGGA

Sequence ID 1387

AAAATT
10 TTTCAGT
GTTAAAGTAGGTTTGTGCGACGCGGCCACGAATTTCCCGGGGACCAA

Sequence ID 1389

TTTTTTTTTTTTTTTTTGGGAGTCAGTTTCTTTTCTTTTCTTTCTTTTTTTTTTTTT
15 GNTTTTCGGAAACGGAGTCTCGCTTTCTCGCCCACTCTGGAGTGGNGCAGTGGGGN
GGTCTCAGCTCACACAGCCTCCACCTCCTGGGCCCAAGCGATCCTNTCACCTCAG
CCTCCTGCGTAGCTGGGACTACAGGCGTGCACCACCATTCCCAGGTAATTTTGTGTA
TTTTTTGTANANACAGGGTTTCACTGTTGTTGCCAGGCTGGTCTCGAACTCCTGC
TTCAGTCTGCCANAATGCTGGATTCTAGGCGTGAGCCACCGNGCCTGGCCCCAAAAG
20 TTAATTTTCTTACAGAAGCAAAGCTTTAATGCATTTTACTGAATGCTTATAGCTTT
GTAGATACTGAAAAGAGTATGAGCGTCACATACAGACACATNTAACAGCACTGCCT
CCAACCAGCCCCCTACCCACTGGTCAGGNGAGTAANAATCAAAATTTCTTTTCTGNGA
GTGGAACGGAAATTTTCTCTCTCCTCCTCAGGCAAGTAGTTAANAGGCTGGNGGGA
GTCATGGCCCCATTTTGTTCAAAATACAAGCTCCACAGGAACAAAAGGCTGAACTG
25 CTCACCTCCCAACTGATGAACCTCGTCTTTGTTCCATGTCAAAGGGGCCTTTGTGT
TACTGCAGCAGAACTCCAGCTATCAAACCATCAGGCACCAAAAGTAAAACTCCTT
TCTCTAAAAAGACCTCTCTTTACCTGAGCCTTTCAATGCATCTTTGCCCCCANATA
ATCCTGGATGAGATAATCCCCAGAGGAANACCAGCGCTTGCCTAGTGAAATTATAC
TATGAGACAAGGGTAAAAGACCTCAAANACCGGGTTGGCAGGTAAGGGAGTAGGGN
30

Sequence ID 1390

TCNGTGGCACCCGTTTCCGGCACCTTCAGACTCTGAAGAGCCACCTGCGAATCCAC
ACAGGAGAGAAACCTTACCATGTACGTAAGCCTCTTGAGGCCGCTCTCTGACCTGC
GGGGATGTGGAGGGCAGGGAAGGAGGTGGAGCGCAGGGAAGGAGGTGGAGCAGGGA
35 GGCAGTGGAACTGTTTGCTCCCATCTCAAGCACACAGTGGGGCAACCACTACGCTA
ATGGTTGGAAGACCTAGATCTGGGCCCAATGGCCAGACACCCTGCTTGACCTTGGC
CCAAGCATTAGGGGACTCATCTTTAAATGAGGGTATGGGACTAGATGATCTGGGC

- 270 -

CTTAGGAGAGGAGT

Sequence ID 1391

5 CGGCTNCTACCCTGCGGAGATCACACTGACCTGGCAGTGGGATGGGGAGGACCAAA
CTCAGGACACCGAGCTTGTGGAGACCAGGCCAGCAGGAGATGGAACCTTCCAGAAG
TGGGCAGCTGTGGTGGTGCCTTCTGGAGAAGAGCAGAGATACACGTGCCATGTTCA
GCACGAGGGGCTGCCGGAGCCCCCTACCCTGAGATGGAAGCCGTCTTCCCAGCCCA
CCATCCCCATCGTGGGCATCGTTGCTGGCCTGGCTGTCCTGGCTGTCTAGCTGTC
CTAGGAGCTATGGTGGCTGTTGTGATGTGTAGGAGGAAGAGCTCAGGTGGAAAAGG
10 AGGGAGCTGCTCTCAGGCTGCGTCCAGCAACAGTGCCCAGGGCTCTGATGAGTCTC
TCATCGCTTGTAAGCCTGAGACAGCTGCCTGTGTGGGACTGAGATGCAGGATTTTC
TTCACACCTCTCCTTTGTGACTTCAAGAGCCTCTGGCATCTCTTTCTGCAAAGGCA
TCTGAATGTGTCTGCGTTCCTGTTAGCATAATGTGAGGAGGTGGAGAGACAGCCCA
CCCCCGTGTCCACCGTGACCCCTGTCCCCACACTGACCTGTGTTCCCTCCCCGATC
15 ATCTTTCTGTTCCAGAGAAGTGGGCTGGATGTCTCCATCTCTGTCTCAACTTCAT
GGTGCCTGAGCTGCAACTTCTTACTTCCCTAATGAAGTTAAGAACCTGAATATAA
ATTTGTTTTCTCAAATATTTGCTATGAAGGGTTGATGGATTAATTAAATAAGTCAA
TTCCTGGAAGTTGAGAGAGCAAATAAAGACCTGAGAACCTTCCANAATCCG

20 Sequence ID 1392

TGAAACAAAATGAATTTNTATGGGTAAGAGAGGGTAATATTTTAGAGTTGTGTTAC
AAACTACAAATTTTTATTAAATTAATAAATCAGAATACTAAATCCATGTGTTTTT
TTCTTTCTTAAAAAATATCTTTTGGCTGGGCACGGTAGCTCATGGCTGTAATCCCA
GCACTTTGGGAGGCTGAGGTGGGTGGATCGCCTGATGTCAGGAGTTCAAGACCAGC
25 CTGGTCAACATGTTGAAACCCCATCTCTACTAAAAATATAAAAAATTAGCCGGTGTG
GTGGTGGGCGCCTGTAATCCCAGCTACTCAGGAGGCTAAGGCAGGAGAATTGCGTG
AACCAGGAGTTCAGTGATGTAGCGGGGAGCTGAGATTGTGCCACTACACTCCAGC
CTGGATGACAGAGTGAGACTCCATCTCAAAAAAAAAAAAAAAAAAAAAA

30 Sequence ID 1394

GCATAATGTGAGGAGGTGGAGAGACAGCCCACCCCGTGTCCACCGTGACCCCTGT
TCCCATGCTGACTTGTGTTTCCTCCCCAGTCATCTTTCCTGTTCCAGAGAGGTGGG
GCTGGATGTCTCCATCTCTGTCTCAACTTTATGTGCACTGAGCTGCAACTTCTTAC
TTCCCTACTGAAAATAAGAATCTGAATATAAATTTGTTTTCTCAAATATTTGCTAT
35 GAGAGGTTGATGGATTAATTAAATAAGTCAATTCCTGGAATTTGAGAGAGCAAATA
AAGACCTGAGAACCTTCCAGAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

- 271 -

Sequence ID 1395

CTTACCATGTCAGTGCACAGAAATGCTGTCTTGGGATGTAGGAAAAATAAATCCAC
AAAAGCTACCAAGTTTGAAGGGGACCATGAGTCTTCAGGCTGGAGCTTCCAAACCA
GATGAAAACCCACAATTAACCTGCAGTTTAAGATCCAGCAGCTGGCCATTTCTGG
5 ACTCAAGGTGAATCGTCTGGATATGTATGGAGAAAAGTACAAACCTTTAAGGGCA
TAAAATACATGACCAAAGCTGGGAAGTTCGAAGTTCGAACCTGAAGGGAGCATTG
CTGAGGGGAATAGTCTTGACATTTTTTTCATTTCTTACTTGTCTAAAAGTAAAAAA
AATATCAGCCTGTCTCCTAGGTCAGTCCCCTCCTGGACCCACCCGCTCCCTTTTTT
CCTTAGCCTTCAGTGCCATGGAACATAATCAAGGGAGGAAAAGGTCACCAGGGAGAA
10 CTGGACAGAACTGAAACACAGCAACACCAGTTCTCAAGGACAAGGTGTGTGATGGG
GGTAGGAAGCTTGGTGCTTATGTAACCATTTTAAACGTGGTTTCTATAGGAAAGAC
CAACATTTGTTTAGCTTGCTTGGCTTTAATTATCTAAAGCCAATGAAAGACTTCTT
TGTTGATTTTTTAAGATAGAAAGATT

15 Sequence ID 1396

CAAACACTATGTTATTTTATGAANAAGACTTGAACATCTATGGATTTTGGTATTTG
CAAGGGGTGAATGGGGTATTTGCAAGCAGTGAATGAGGAGGCCTGGAACCAATCTT
CTGCTGATATTGAGGCACAACCTGAAAAAGGTATATTACTTAAATCTCTTATTGTAT
TGTAAACTGTATAAGTAATGAAATTAAAAGGCAGAAATTGTCAGACTGAATAAAAT
20 GAAAAGACCAAACAATATGCTGCTTACAAGAAACACAATTCAAATATAAGGACACA
ATTAGTTTAAAGGAAAAGAACTGGAAAAGATATACCATGATAACACAAGTCAGAAG
AAAGCTGCTGTGGATATATTAATATGAGATGTAGATTTTCAGAGCAGTGAATATTGC
CAGGCATAAAGAAAGTTATTACATAATAATTAAGGTATCAGTTCATCAAGAAGATG
TAATAACCCTAAGTATTTATACAACATAATCAGAGCTTCAAAATACATGAAGCAA
25 AAACCAGTGGAATTGATAGGAGAAACACACAATTACACAATTATAGTCAGAATTTT
CAACATATCTTTCTCAATGGAGAAAACAACCTAGACAGGAAATCATTAAGGATATAG
ATGATTTAAATTATATGATCAACTACCTGGACGTAATTGGCATTATGGAACACTG
CACCACCAACAGCAGAGTACATATTATTTTCAAGTACACAGAAAACAGTTACCAAT
ATAGACCATTTTCTGGGTCATAAAACACATCTCAATAAATGTAAAACAATTAATGT
30 TATATAAAGTATGTGCTCTGACCNCAAAGGAATTAGAGATCAATAAAAGAACATCT
TTGAAAAATCTCACNTATTTAAAACTAATAACTCACTTCTAAATAACTCCTGTNT
CAAGAGAATNAAANGG

Sequence ID 1397

35

CCCAGCCTCACTGCGCCCCGTCAGGCCAGGCAGCTGCCCTCAGGGTCTGCCAAGGT
GGGGGTCAAGGGCCATGGGGGCAGGTAGCTCTGCCTGCAAAGCCCACAAGCATGTC

- 272 -

AGATCACCTGGGCTGCAGACAGACAAACACCTGAGCTGTTCTGAATACCTTCAGGT
TCCTGGCCTCGCTGAGCAAGTGCAGAAATTTTACCTTCAAGGATCAGGGTTTTTC
TGTTTGTTTGTTTTTTAACACACACATATGTGAACAAAGAGTATGCGTTTGTA
GCAGAAGAAGCGTCTGGTAAGACAACCAGCAAGTTAACAATGGTCACCTCCAGAAA
5 TGGGCTGGGTAAACCAAAGAATTTTTTGTGTTTTTGTGTTTTTTTGAGTCAGGGTCTA
GCTCTGTCACCCAGGCTGGAACGCACTGGTGTGATCACGGCTCACTGCAGCCTTGA
CCTCCCTGGCTCAAGCAATCCTCCCAGCTCAGCCTCCTGAGTCGTTGGGACTACAG
GCACGTGCCACCACGCCTGACACATTTTTTAAATTTTTGTAGAGACAGTGTTTCAC
CATGTTGCCCAGGCAGGTCTCAAACCTGGGCTCAAGTGGTCCTCCAGCTTCAGC
10 CTCCCAAAGTGCTAGGATTATAGGTGTGAGCCACAGTGCCAGCCCCGTAGTGGAG
AATTTCTGTTGAATGAACCAAAGCAACTGCCAACCTCTCCATGCACCATGTGTTT
CAGAGGAGAAAGCACAGTGAAGAATGCAGTGTGTTCTGAGGTCCTGTCACCCCTGA
GGCTGTGTGTGTCCTTTGCCAAATTAAAGAGTCTTACTGAATGCGGTGCATCCAGG
AGACAGGCCNAGGTTTGGACTGGTAAAAA

15

Sequence ID 1399

CAGACACCTGGNAGAACGGGAAGGAGACGCTGCAGCGCGCGGACCCCCAAAGACA
CATGTGACCCACCACCCCATCTNTGACCATGAGGCCACCCTGAGGTGCTGGGCCCT
GGGCTTCTACCCTGCGGAGATCACACTGACCTGGCAGCGGGATGGCGAGGACCAAA
20 CTCAGGACACCGAGCTTGTGGAGACCAGACCAGCAGGAGACAGAACCTTCCAGAAG
TGGGCAGCTGTGGTGGTGCCTTCTGGAGAAGAGCAGAGATACACATGCCATGTACA
GCATGAGGGGCTGCCGAAGCCCCCTCACCTGAGATGGGAGCCATCTTCCCAGTCCA
CCGTCCCCATCGTGGGCATTGTTGCTGGCCTGGCTGTCCTAGCAGTTGTGGTCATC
GGAGCTGTGGTCGCTGCTGTGATGTGTAGGAGGAAGAGTTCAGGTGGAAGGAGG
25 GAGCTACTCTCAGGCTGCGTCCAGCGACAGTGCCAGGGCTCTGATGTGTCTCTCA
CAGCTTGAAAAGCCTGAGACAGCTGTNTTGTGAGGGACTGAGATGCAGGATTTCTT
CACGCCTCCCCTTTGTGACTTCAAGAGCCTCTGGCATCTCTTTCTGCAAAGGCACC
TGAATGTGTCTGCGTCCTTGTTAGCATAATGTGAGGAGGTGGAGAGACAGCCCACC
CTTGTGTCAACTGTGACCCCTGTTCCCATGCTGACCTGTGTTTCTCCCCAGTCA
30 TCTTTTTTGTTCNCAATAGGTGGGGCCTGGATGTCTCCATCTCTGTNTCA

Sequence ID 1440

TTATAAGTACTTTTAAGGTATTTTAGTTGTCTTAGTCTATATTTCTGTACTCACC
35 TTTCTTTATCCACTCATCAGTTGATGGGCATGTAGGTTGGTTCCATATCTTTGCAA
TTCTGAATTGTGCTGTGATCAGGTGTCTTTTGTAGTATAATGATTTACTCTCCTTTG
GGTAGATACCCAGTAGTGGGATTGCTGGATCGAATGGTTTTTATAATTTCTATTT

- 273 -

TACCACAGTTTCTCTCTGCATTTTTCCTCTTTGACCACTAACCATGTGAAATTCTC
ATATTGACCTTTATAATGATCATGAACTCTTAGTATCATTGGGAAGGCCACATTTG
CCACTTATGATTGTAAACCTTATCCTCCATTTTTCCTGTTATTGTTGGTGCAAAA
GCACCTATTATACCAGGACTTTAAAAATCAGTCTGATAAGTCTTTGATAAGTCTAA
5 TAATAATAACTGATAAGTCCATTGAATTTGCTTCTGATTACTTTTTCTTTAGTAGC
TAAACATGTATGTACTCCTATGATTACAATGAACACTCCTCTCCATTTAAATTAAT
TATTTACATTGATGAAATAGCAAAATGTTAATGACTAAATACTGTCTTGGTTTTTT
CGTTCCAGGTCAGTCAATATTAACCTTCTTATAATTTTCTTTTTTTCTTT

10 Sequence ID 1447

GCAAGGACTAACCCCTATACCTTCTGCATAATGAATTAAGTAGAAATAACTTTGCA
AGGAGAGCCCAAAGCTAAGACCCCCGAAACCAGACGAGCTACCTAAGAACAGCTAAA
AGAGCACACCCGCTCTATGTAGCAAAATAGTGGGAAGATTTATAGGTAGAGGCGACA
AACCTACCGAGCCTGGTGATAGCTGGTTGTCCAAGATAGAATCTTAGTTCAACTTT
15 AAATTTGCCACAGAACCCCTCTAAATCCCCTTGTAATTTAACTGTTAGTCCAAAG
AGGAACAGCTCTTTGGACACTAGGAAAAACCTTGTAAGAGAGAGTAAAAAATTTAA
CACCCATAGTAGGCCTAAAAGCAGCCACCAATTAAGAAAGCGTTCAGCTCAACAC
CCACTACCTAAAAAATCCCAAACATATAACTGAACTCCTCACACCCAATTGGACCA
ATCTATCACCCCTATAGAAGAACTAATGTTAGTATAAGTAACATGAAAACATTCTCC
20 TCCGCATAAGCCTGCGTCAGATTAAACACTGAACTGACAATTAACAGCCCAATAT
CTACAATCAACCAACAAGTCATTATTACCCTCACTGTCAACCCAACACAGGCATGC
TCATAAGGAAAGGT

Sequence ID 1448

25 GGCCACCGGGTGCAAGGTCAGGGCTGGGGTGAGGCTGGGAAGCCCAGGGCTTGGC
CCACTGTGGCCGCCTTGTGTGGTCACTGCTTTCCTGGGCCTGCTGTGAGCTCCCTC
TAGGACCCCAGGCCTGTCTGGTGGGTCACTGTGACCACCACCTTGCACAGCACCTG
GCGCGTGGCAGGTGCTCAAACATTACTTGTTCGGAATGAACTTCATCTTGCTCTT
GGCTTTTGGACTAATGCTGTGGAACATCTGACTAATTAGTGACTCTTTGGGGCCCC
30 CAGTTTCCCAGCTATAAAGTGGTAATATTAAGATAATAATTCGGCCGGGCGCGGTG
GCTCACGCCTGTAATCCCAGCAGCACTTTGGGAGGCCGAGGTGGGCAGATCACGAG
GTCAGAAGATCGAGACCATCCTGGCTAACACGGTGAAACCCCATCTCTACTAAAAA
TACAAAAAATTANCCGGGCGTGGTGGCGGGCGCCTGTAGTCCCAGCTACTCANGAG
GCTGANGCAGGAGAATGGTGTGAACCCGGGAGGCAGAGGTTGCAGTGAACCAAGAT
35 CGNNCCACTGCACTCCAGCCTGGGCAACAGAGCGAGACTCCATCTTAAAAA

Sequence ID 1449

- 274 -

AATCAGGGCCGCAGTGTGTTCTGCGCCTGCCAGAGCTGACTCCTGATTTAACCGC
TGGCGTAACCGCGGGTTGCACGCATGCGTGCTGAAAAGCCTTTCACCCTCACGTGG
TTTCTTTTTTAACCAAGTCATCAAGCGAGGCTCGCGCGCAGGCCCCGCGTTGGAAAA
TGGCGGGGAAGCTGAAACCTCTGAATGTGGAGGCGCCAGAAGCTGCTGAGGAGGCT
5 GAAGGTAGTGAGGGCAAGTGGGCTGCACTCCTTTCTCTCCAACCAGGGCAGAAAGG
AGGGAGGATTTCGTCCCATTTACAATAATGAAATAATGATATTCTAATTTTTTTTAAAT
AAAATGTTAAGCCTTTTGTTATTGAA

Sequence ID 1450

10 GGAAANCATGAGGCTTCGGGAGCCGCTCCTGAGCGGCAGCGCCGCGATGCCAGGCG
CGTCCCTACAGCGGGCCTGCCGCTGCTCGTGGCCGTCTGCGCTCTGCACCTTGCG
GTCACCCTCGTTTACTACCTGGCTGGCCGCGACCTGAGCCGCCTGCCCCAACTGGT
CGGAGTCTCCACACCGCTGCAGGGCGGCTCGAACAGTGCCGCCGCCATCGGGCAGT
CCTCCGGGGAGCTCCGGACCGGAGGGGCCCCGGCCGCCGCTCCTNTAGGCGCCTCC
15 TCCCAGCCGCGCCCGGGTGGCGACTCCAGCCCAGTCGTGGATTCTGGCCCTGGCCC
CGCTAGCAACTTGACCTCGGTCCCAGTGCCCCACACCACCGCACTGTCGCTGCCCG
CCTGCCCTGAGGAGTCCCCGCTGCTTGTAAGGACTCGGGTCGGCGCCAGTCGGAG
GATTGGGACCCCCCGGATTTCCCCGACAGGGTCCCCCANACATTCCTCAGGCTG
GCTCTTCTACGACAGCCAGCCTCCCTCTTCTGGATCAGAGTTTTAAATCCCANACA
20 GAGGCTTGGGACTGGATGGGAGAGAAGGTTTGCGAGGTGGGTCCCTGGGGAGTCCT
GTTGGAGGCGTGGGGCCGGGACCGCACAGGGAAGTCCCGAGGCCCTCTAGCCCCA
AAACCANAGAAGGCCTTGGGAGACTTCCCTGCTGTGGCCCGAGGCTNAGGAAGTTTT
GGAGTTTTGGGTCTGCTTANGGCTTCNAGCAGCCTTGCACTGAGAACTTTGGTAGG
GACCTCGAGTAATCCACTCCNTTTTNGGGACTGACGTGAGGCTCCCGGTGGGGAAA
25 GANACTGACCTNTC

Sequence ID 1453

CCGACCTGTCTCGCTCCGTGGCCTTAGCTGTGCTCGCGCTACTCTCTCTTTCTGGC
CTGGAGGCTATCCAGCGTACTCCAAAGATTCAGGTTTACTCACGTCATCCAGCAGA
30 GAATGGAAAGTCAAATTTCTGAATTGCTATGTGTCTGGGTTTCATCCATCCGACA
TTGAAGTTGACTTACTGAAGAATGGAGAGAGAATTGAAAAAGTGGAGCATTCAGAC
TTGTCTTTTACGAAGGACTGGTCTTTCTATCTCTTGTACTACACTGAATTCACCCC
CACTGAAAAAGATGAGTATGCCTGCCGTGTGAACCATGTGACTTTGTACAGCCCA
AGATAGTTAAGTGGGATCGAGACATGTAAGCAGCATCATGGAGGTTTGAAGATGCC
35 GCATTTGGATTGGATGAATTCCAAATTCGTCTTGCTTGCTTTTTTAATATTGATATG
CTTATACACTTACACTTTATGCACAAAATGTAGGGTTATAATAATGTTAACATGGA
CATGATCTTCTTTATAATTCTACTTTGAGTGCTGTCTCCATGTTTGATGTATCTGA

- 275 -

GCAGGTTGCTCCACAGGTAGCTCTAGGAGGGCTGGCAACTTAGAGGTGGGGAGCAG
AGAATTCTCTTATCCAACATCAACATCTTGGTCAGATTTGAACTCTTCAATCTCTT
GCACTCAAAGCTTGTTAAGATAGTTAAGCGTGCATAAGTTAACTTCCAATTTACAT
ACTCTGCTTAGAATTTGGGGGAAAATTTAGAAATATAATTGACAGGATTATTGGAA
5 ATTTGTTATAATGAATGAAACATTTTGTGCATATAAGATTCATATTTACTTCTTAT
ACA

Sequence ID 1454

TAAATAGGGAATCCTTTCCCATTTGCTTGTTTTCTCAGGTTTGTCAAAGATCAGA
10 TAGTTGTAGATATGCGACGTTATTTCTGAGGGCTCTGTTCTGTTCCATTGATCTAT
ATCTCTGTCACATGCACACGTATGTTTGTGTTGGCACTATTCACAGTGGCAAAGAC
TTGGAACCAACCCAAATGTCCAACAATGATAGACCGGGTTAAGAAAATGCGGCACA
TATACACCATGGAATACTATGTAGCCATAAAAAATGATGAGTTCGTGTCCTTTGTA
GGGACATGGATGAAATTGGAAATCATCATTCTCAGTAACTATCGCAGGAACAAAA
15 AACCAAACACTGCATATTCTCACTCATAGGTGGGAATTGAACAGTGGGAACACATG
GACACAGGAAGGGGAACATCACACTCTGAGGACTGTTGTGGGGTGGGGGGAGGGAG
GAGGGATAGCATTGGGAGATATACCTAGTGCTGGATGACGAGTTAGTGGGTGCAGC
GCACCAGCATGTCACATGTATACATATGTAACCTGACATTTGTGCACATGTA
CCCTAAAACTTAAGGTAT

20

Sequence ID 1456

CCGCAACAAACACGGGAGTGCAGATATCGCTGCGATGGGCTGATTTCTTTTATTTG
GGTATATACCCAGCAGTGGGATTGCTGGATTGTATGGTAGCTCTATTAGTTTTTTG
25 AGGAACCTCCAACTGTTCTNCATAGTGGTTGTACTCATTTACATTCCCACTGTGA
ACCCTGAAAATTTGAGGCAGGTCTCAGTTAAATTAGAAAGTTGATTTTGCCAAGTT
GGGGACACGCACTCGTGACACAGCCTCAGGAGGAACTGATGACATGTGCCCAGGTG
GTCAGAGCACAGCTTGTTTTATACATTTTAGGGAAACCTGAGCCATCAATCAACA
TACGTAAAATGGGCCGGGCACAGCAGCTCAAGCTGTAATCCCAGCACTCTGGGAGG
30 CCGAGGCGGGTGGATCACTTGAGGTGAGGAGTTCGAGACCAGCCTGGCCAACATGG
TGAAACCCCGTCTCTATTAAAAATACAAAGCTTAGCTGGATGTGGTGGCGCATGCC
TGTAGTCCCAGCTGCTCTAGGAGGCTGAGGCATGAGAATTGCTTGAACCTGGGAGG
CAGAGGCTGCAGTGAGCCGAGATCGAGCCACTATACTCCAGCCTGGTCAACAGAGT
GAGACCCTGTCT

35

Sequence ID 1460

CCCAACTGTGTTCACTAGCAACCTCAAACAGACACCATGGTGCACCTGACTCCTG

- 276 -

AGGAGAAGTCTGCCGTTACTGCCCTGTGGGGCAAGGTGAACGTGGATGAAGTTGGT
GGTGAGGCCCTGGGCAGGCTGCTGGTGGTCTACCCTTGGACCCAGAGGTTCTTTGA
GTCCTTTGGGGATCTGTCCACTCCTGATGCTGTTATGGGCAACCCTAAGGTGAAGG
CTCATGGCAAGAAAGTGCTCGGTGCCTTTAGTGATGGCCTGGCTCACCTGGACAAC
5 CTCAAGGGGCACCTTTGCCCACTGAGTGAGCTGCACTGTGACAAGCTGCACGTGGA
TCCTGAGAACTTCAGGCTCCTGGGCAACGTGCTGGTCTGTGTGCTGGCCCATCACT
TTGGCAAAGAATTCACCCACCAGTGCAGGCTGCCTATCAGAAAGTGGTGGCTGGT
GTGGCTAATGCCCTGGGCCACAAGTATCACTAAGCTCGCTTTCTTGCTGTCCAATT
TCTATTAAAGGTTCTTTGTTCCCTAAGTCCAATACTAAACTGGGGGATATTATG
10 AAGGGCCTTGAGCATCTGGATTCTGCCTAATAAAAAACATTTATTTTCATTG

Sequence ID 1490

ATGGGCATCTCTCGGGACAACCTGGCACAAGCGCCGCAAAACCGGGGGCAAGAGAAA
GCCCTACCACAAGAAGCGGAAGTATGAGTTGGGGCGCCAGCTGCCAACACCAAGA
15 TTGGCCCCCGCCGCATCCACACAGTCCGTGTGCGGGGAGGTAAACAAGAAATACCGT
GCCCTGAGGTTGGACGTGGGGAATTTCTCCTGGGGCTCANAGTGTTGTACTCGTAA
AACAAGGATCATCGATGTTGTCTACAATGCATCTAATAACGAGCTGGTTCGTACCA
AGACCCTGGTGAAGAATTGCATCGTGCTCATCGACAGCACACCGTACCGACAGTGG
TACGAGTCCCACTATGCGCTGCCCCCTGGGCCGCAAGAAGGGAGCCAAGCTGACTCC
20 TGAGGAAGAAGAGATTTTAAACAAAAACGATCTAAAAAATTCAGAAGAAATATG
ATGAAAGGAAAAAGAATGCCAAAATCAGCAGTCTCCTGGAGGAGCAGTTCAGCAG
GGCAAGCTTCTTGCGTGATCGCTTCAAGGCCGGGACAGTGTGGCCGAGCAGATGG
CTATGTGCTAGAGGGCAAAGAGTTGGAGTTCTATCTTAGGAAAATCAAGGCCCGCA
AAGGCAAATAAATCCTTGTTTTGTCTTCACCCATGTAATAAAGGTGTTTATTGTTT
25 TTGTT

Sequence ID 1491

CTTNCACATACTGATTGATGTCTCATGTCTCTCTAAAATGTGTAAAACCAAGCTGT
GCCCCAACCACCTTGGGNACATGTGGNGAGGACCTCCTGAGGCTGTGTCATGGGCA
30 CACCTTAACCCTGGGAAAATAAACTTTCTAACTGACTTGAGAGCTGTCTCAGATA
TTCTGAGCTTACAGTTATTGTGAAATCATTTTAATTATAAATTAAGTGGAGATTTA
CTTAAATCATGTGTAGAAGTAGCCTGTGATATAGTCCTAGATACATACATTATCA
TCTTATGTATCTTCCCTCCCTCTTCCAGGTTCTGATAAAAACAGATGAAATCTGAA
AGACCATGACAGTAGTATTTTGAAAATGACAGTATTTGAAATTAAAAAATTGTAAA
35 AGTGTCTGTCTATCACTGCCAAAGGATAAGTTACAAATTGGTTCTTGGAACGTA
ATATGTACTATGTGCTTGCTATTTAATAATTTACCAGTCTTAGTCTTTTTTATTCA
GACTAATTTTACCTTTTTTTTAACTATGACTCTTTAGTTATAGTAGTACAAAAAAG

- 277 -

TAGTTTTAGTTATAGTTTTAGTTGTAGTACAAAAAGCATTTTCTGTAAGCTTAAT
TTCTTTCCCCTTCCCGCTTTCCAGTCAGATGACTTTAGTGATTTGGAGTTGTGTG
CTTTATAAGTGCATTCCTCAGAGGACTTAATATTACTAAGATTTTAGCAACNCTGA
AATATGTT

5

Sequence ID 1492

TGTNCCTGTAGTCCTGTGTGGGAGGATTGCCTGAGCCTAGGAGCTCAAAGTTGCAG
TGAGCCCAGATCGNGNCATTGCAGTCCAGCCTGGGTGACAGAGTGAGACCCCATGT
CAAAAAAAAAAAAAACAAAAACAGGGGCCTGCCTCANCCAGCAGGTGAGGTCTGCC
10 ACTGAGAGCACTTCTAGCAGCAGGAACAGCCTCCACCCCCACACTGCAATCAAGTT
TTTTGGGTGAGCCTTAGGAGCTAANAAAGGGCCTAGTTTGNCTAAATAGCAGGAGT
TATATCCAGGGATCTT[•]CAGGCCCAGGAATGCTAATGAGTAGGCATTCCATGGGCCC
TGGGAATGGCTTTGTGTGCCANAAATGATGGCCACAAAGGCCTTGCTGCCTTTTTTT
CAAATGGCTGCATCCAGCTGAGTGCTCTCTGCCAAAGGGGANAANAATAAGTC
15 TCCAGTGCATTTAGATTGGTCTCTCATCATCTCTCTCCTTTTTTGTTTTTATTAGTC
TCCTTAACCAA[•]AACTGCCAAGAAAGGCTTGGAAATTGAAACAAAACCTGATANAANA
GGTAAGAGGTTGTTCTTTT

Sequence ID 1493

TGTNTCAAAAAAAAAAAAAAGAACGGNAATGTACTGGAGATGTATTTGATAACCAA
20 GGNTTTAGGTAAATTTTCACCAGTATTAGTTNTATTTGCAAACCTGAAAAATGTTGT
AGGCTTAATATAAAATAACCACATTAGTGAACATTATATCTCTTAGAAGAAAGGCC
ATATTTTGCTCCTGCTTCTGTAAAAATATTATTTGTTTGAAGGGGAAATAATGGTA
GTGTGACCTTTCACTTAATTCCTACTCCCTTAATGTGAGAGAGACAAAATGAGCTG
25 AAGAAGGAA[•]AATCTGGAGTTACACTCCACAACCTTGAACATACTGACGGACATCT
CTGTTTTGACAACGATTTCTCCATGCCACCCATGCTNTAATGCCTTGTGGATCACG
GACAACCCTCTTTGCACAAGCTACAGCATCAGCGATGTTATCTTGCAGCAAAGCAC
TGCAGGATAAATGACAGGCATTAAC[•]TGCTCCTGGGGTTTTGCCATCATTACACCAG
TAGCGGCTATTGATCTGAAATATCCATAATCAGTGCTTCTGTCTCCAGCATTGTA
30 GTTTGTAGCTCGTGTGTTGTAACCACTCTCCCATTTGGCCAAACACATCCAGTTTG
CTAGGCTGATTTCCCCTGTAGCCATCCATTCCCAATCTTTTCAGAGTTCTGGCCAAC
TCACACCTTTCAAAGACCTTGCCCTGGACCGTAACAGAAAGGAGGACAAGCCCCAG
AACAAATGAGAGCCTTCATGTTGAC

35

Sequence ID 1494

TTGGTACCCGGGAAATTTCTTTGCCGCGTCGACGGCCGGTGAGGCAGATCACCTGAG
CCCAGGAGTTCAGGACCAGCCTGGGCAGCATACCGGGATTCCATCTNNACTAAAAA

- 278 -

CAGTAGGCTGGGTGTGGTGGCTCATGTCTGTAAGCTCAGGACTTTGGAAGGCCAAG
ATGGGAGGATCACTTGAGCCTGGGAGTTTGACACCAGCTTGAGCATCGTAGCCAGG
CCCTGACTCTACAAAAAAGTGAAATAATTAGCCGAGTGTGGTGGTTCACACCTGTA
ATCCCAGCTGCTCAGGAGGCTGAGGTAGGAGAATCATTTGAACCCGGGAGGTGGAG
5 GTTGCAGTTAGCCGAGATCACGCCATTGCACTCCGGCCTGGGCGATAAAGCGAGAC
TCTGTCTCAAAAAAAAAAAAAA

Sequence ID 1495

ATTCGGGCCGAGATGTCTCGCTCCGTGGCCTTAGCTGTGCTCGCGCTACTCTCTCT
10 TTCTGGCCTGGAGGCTATCCAGCGTACTCCAAAGATTGAGTTTACTCACGTCATC
CAGCAGAGAATGGAAAGTCAAATTTCTGAAATTGCTATGTGTCTGGGTTTCATCCA
TCCGACATTGAAGTTGACTTACTGAAGAATGGAGAGAGAATTGAAAAAGTGGAGCA
TTCAGACTTGTCTTTTCAAGCAAGGACTGGTCTTTCTATCTCTTGTACTACACTGAAT
TCACCCCCACTGAAAAGATGAGTATGCCTGCCGTGTGAACCATGTGACTTTGTCA
15 CAGCCCAAGATAGTTAAGTGGGATCGAGACATGTAAGCAGCATCATGGAGGTTTGA
AGATGCCGCATTTGGATTGGATGAATTCCAAATTCGTCTTGCTTTTAAATAT
TGATATGCTTATACACTTACACTTTATGCACAAAATGTAGGGTTATAATAATGTTA
ACATGGACATGATCTTCTTTATAATTCTACTTTGAGTGCTGTCTCCATGTTTGATG
TATCTGAGCAGGTTGCTCCACAGGTAGCTCTAGGAGGGCTGGCACCTTAGAGGTGG
20 GGAGCAGAGAATTCTCTTATCCAACATCAACATCTTGGTCAGATTTGAACTCTT

Sequence ID G6

GGATTTTTGGTCCGCACGCTCCTGCTCCTGACTCACCGCTGTTGCTCTCGCCGAG
GAACAAGTCGGTCAGGAAGCCCGCGCGCAACAGCCATGGCTTTTAAGGATACCGGA
25 AAAACACCCGTGGAGCCGGAGGTGGCAATTCACCGAATTCGAATCACCTAACAAG
CCGCAACGTAAAATCCTTGGAAGGTGTGTGCTGACTTGATAAGAGGCGCAAAAG
AAAAGAATCTCAAAGTGAAAGGACCAGTTCGAATGCCTACCAAGACTTTGAGAATC
ACTACAAGAAAACTCCTTGTGGTGAAGGTTCTAAGACGTGGGATCGTTTCCAGAT
GAGAATTCACAAGCGACTCATTGACTTGACACAGTCCTTCTGAGATTGTTAAGCAGA
30 TTACTTCCATCAGTATTGAGCCAGGAGTTGAGGTGGAAGTCACCATTCAGATGCT
TAAGTCAACTATTTTAATAAATTGATGACCAGTTGTTAAAA

Sequence ID - 61 nt: 362

CTTATTGAAAATTTTACTAATTTCTTACTTTTTAGGTTTTAGGAGAATACTTTTGGA
TAATTGACTAGCCTCACATTATATTGATAGAGGTTCTTGAAAACTTTAATGCCAAT
TCATGTATCTTATGACTAAAATAGATAATCCATTTAGAAATTTAAGTCATTCTTGC
GTGCTTGATATGTGTCAGCACTATCCAAGTTGCTAGGGGATAACAATGGTGAAGTG
AAAATATCAGCTAGGTGCCGGTGGCTCACACCTGTTATCCCAACAGTTTGGGAGG
CCAGGGTGGGAGGATCACTCAAGCACANGCGTTTCACACCAGCCTGGACAACAT
ACAAGACCCCATCTTTACCAAAAGTTAAG

Sequence ID - 490 nt: 382

TTTTCTTAGAACTTTATTTTTTCTGGCCAGGCGCAGTGGCTCACACCTGTAATCCC
AGCACTTTGGGAGGCCAAGGCAGGTCGATCACCTGAGGTCAGGAGCTCAAGACC
AGCCTGGCCAACATGGTGAAACCCTGTCTCTACTAAAAATACAAAAATTAGCTGG
GCGTGGTGGCGCATGCCTGTAATCCCANCTACTCAGGAGGCTGAGGCAGGAGAA
TTGTTTGAACCCGGGAGGCGGAGGTTGCANTGAGCCGAGATTGCGCCACTGCACT
CCAGCCTGGGCAACAGAGCGAACTCCATCTCAAAAAAAAAAAAAAAAAAAAAAAC
CTTTATTTTTTCTGATTTTAAAAGTAATAACTAGTTTGTAGAAACATTAAAAGT

Sequence ID - 892 nt: 559

TCTTTCGGAAGCGCGCCTTGTGTTGGTACCCGGGAATTCGCGGCCGCGTCGACGC
GGTCGTAAGGGCTGAGGATTTTTGGTCCGCACGCTCCTGCTCCTGACTACCGCT
GTTCGCTCTCGCCGAGGAACAAGTCGGTCAGGAAGCCCGCGCGCAACAGCCATG
GCTTTTAAGGATACCGGAAAAACACCCGTGGAGCCGGAGGTGGCAATTCACCGA
ATTCGAATCACCTAACAAGCCGCAACGTAAAAATCCTTGGAAGAGGTGTGTGCTG
ACTTGATAAGAGGCGCAAAAGAAAAGAATCTCAAAGTGAAAGGACCAGTTCGAA
TGCCTACCAAGACTTTGAGAATCACTACAAGAAAACTCCTTGTTGGTGAAGGTTT
TAAGACGTGGGATCGTTTCCAGATGAGAATTCACAAGCGACTCATTGACTTGCAC
AGTCCCTTCTGAGATTGTTAAGCAGATTACTTCCATCAGTATTGAGCCAGGAGTTG
AGGTGGAAGTCACCATTGCAGATGCTTAAGTCAACTATTTTAATAAATTGATGAC
CAGTTGTTT

Sequence ID - 77 nt: 464

GCGGCTGCTGTTGGTTGGGGGCCGTCCCGCTCCTAAGGCAGGAAGATGGTGGCCG
CAAAGAAGACGAAAAAGTCGCTGGAGTCGATCAACTCTAGGCTCCAACCTCGTTAT
GAAAAGTGGGAAGTACGTCTTGGGGTACAAGCAGACTCTGAAGATGATCAGACA
AGGCAAAGCGAAATTGGTCATTCTCGCTAACAACCTGCCCAGCTTTGAGGAAATCT
GAAATAGAGTACTATGCTATGTTGGCTAAAACTGGTGTCCATCACTACAGTGGCA
ATAATATTGAACTGGGCACAGCATGCGGAAAAATACTACAGAGTGTGCACACTGG
CTATCATTGATCCAGGTGACTCTGACATCATTAGAAGCATGCCAGAACAGACTGG
TGAAAAGTAAACCTTTTCACCTACAAAATTTACCTGCAAACCTTAAACCTGCAA
AATTTTCCTTTAATAAAATTTGCTTG